# AutoRep™ E Repeating Dispenser

## Operating Manual

Please read carefully before use, and follow all operating and safety instructions!

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## Contents

	Page
Safety considerations	4
Purpose and limitations of use	7
Operating exclusions	8
Operating elements	9
Initial operation	10
Insert battery pack	10
Charge battery pack	11
Replace battery pack	14
Switching on/off	15
Inserting tips	16
Removing tips	18
Functions of the AutoRep E	19
Aspirating sample	20
DISP mode	22
AUTO DISP mode	24
PIP mode	26
Adjusting STEP volume	27
Adjusting aspiration and dispensing speed	28
Changing mode	30
Use of compatible tips from other manufacturers	31
Residual stroke	33
Non-contact dispensing	33
Maintenance and cleaning	33
Default settings	33
Troubleshooting	34, 35
Technical data	42
Accuracy table	43
Ordering data	44
Return for repair	46
Disposal	47
Warranty	48

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## Safety considerations

This instrument may sometimes be used with hazardous materials, operations, and equipment. To address all of the safety problems associated with its use in such applications is beyond the scope of this manual. It is the responsibility of the user of this instrument to consult and establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

1. Every user must read and understand this operating manual prior to using the instrument and observe these instructions during use.
2. Observe general instructions for hazard prevention and safety instructions; e.g. wear protective clothing, goggles and gloves. When working with infectious or other hazardous samples, all appropriate regulations and precautions must be followed.
3. Observe all safety precautions provided by reagent manufacturers.
4. Never use the instrument in an atmosphere with danger of explosion. Highly flammable liquids must not be pipetted.
5. Use the instrument only for dispensing liquids, with strict regard to the defined Operating

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## Safety considerations

- continued -

14. To charge the NiMH battery pack use only the original AC adapter in combination with the original charging dock.
15. The AC adapter has to be protected against moisture and must be used only for this instrument.
16. Dispose of battery pack only when completely discharged, and according to applicable regulations.

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### Note:

Attention: Risk of injury!  
Do not touch the plunger seat inside the instrument!

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### Warning!

Improper use of the instrument or the batteries (short circuit, mechanical damage, overheating, incorrect AC adapter etc.) may, in extreme cases, lead to explosion of the batteries.

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## Purpose

The AutoRep E is a microprocessor-controlled battery-operated repetitive pipette. The instrument automatically recognizes new Encode™ Syringe Tips, which are encoded for tip volume. This allows fastest possible volume selection. If you want to use compatible tips from other manufacturer, you must set the volume manually.

**From now on, the term "tips" refers to Encode Syringe Tips and compatible tips.**

## Limitations of use

The instrument is designed for repetitive dispensing of samples, observing the following limitations

- +15 °C to +40 °C (59 °F to 104 °F) (instrument and reagents)
  - other temperatures on request
- vapor pressure up to 500 mbar
- viscosity:
  - 20 mPa s with 50 mL Encode tips
  - 260 mPa s with 5 mL Encode tips
  - 977 mPa s with 1.25 mL Encode tips

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## Operating exclusions

During proper operation, the dispensed sample only comes in contact with the tip and not with the AutoRep E.

The user has to ensure the compatibility of the instrument with the intended application.

Never dispense samples which attack polypropylene, polyethylene or polycarbonate (housing). Avoid aggressive fumes (danger of corrosion)!

The device may not be used for oxidizing acids, since they can attack metallic parts and the electronics.

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### Note:

The instrument and the tips are not autoclavable (121 °C / 250 °F). See sterile tips on page 45.  
Encode tips are disposables. For disinfection you may use commercially available disinfecting solutions as long as they do not attack PC/PBT or damage the display.

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## Initial operation

### Insert battery pack

1. Open the battery compartment cover and insert battery pack, ensuring that the electrical contacts on the battery pack connect with those of the AutoRep E.
2. Replace the battery compartment cover.

#### Warning!

The instrument must not be used or charged in potentially explosive environments.

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
## Charge battery pack

The battery pack can be charged either inside the AutoRep E or in the charging dock.

**Before first use charge the battery pack inside the AutoRep E in the original charging dock for at least 3.5 hours. After 2 - 3 further charging cycles inside the instrument, the full capacity is reached.**

**Charge indicator is most accurate when battery is charged 3.5 hours in the instrument.**

**Removal of battery pack resets indicator. Upon reinsertion battery-charge frame flashes to indicate that charge status is approximate**

### Note:

Stored data and settings are re-tained when the battery pack is completely discharged or during replacement. Remove the battery pack when not in regular or daily use.

- One single tone indicates that the battery pack has to be charged or replaced. The battery symbol will flash on the display.

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## Charging inside the AutoRep™ E

### Attention:

Before charging, make sure that the AC adapter is appropriate for your line voltage.

1. Plug the AC adapter connector into the side of the charging dock.
2. Place the instrument with the battery pack into the charging dock. A tone indicates proper positioning of the instrument.
3. A green light on the charging dock indicates that the instrument is charging.

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The charging time for a nearly empty battery pack mounted in the instrument is approximately 3.5 hours.

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### Attention:

Upon each battery insertion, charge-indicator frame flashes until battery is charged for at least 3.5 hours inside instrument.

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## Charging the battery pack in the charging dock

The charging dock allows you to charge a separate battery pack while the AutoRep E is in use.

1. Plug the AC adapter connector into the side of the charging dock.
2. Insert the battery pack into the central opening of the charging dock, making sure the battery contacts connect with the charging dock contacts.
3. A green light on the charging dock indicates that the battery pack is charging.

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The charging time for a nearly empty battery pack in the charging dock is approximately 2.5 hours.

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### Replace battery pack

A fully charged battery pack provides enough power for approximately 2.5 hours of continuous use with a 12.5 mL Encode Tip filled with a sample with similar viscosity and density to water.

- A tone indicates the need for recharging.
- The battery symbol will flash on the display.

When indicated, the battery should then be recharged or replaced. Before removing the battery, make sure the tip is completely empty.

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## Switching on/off

- The AutoRep E can be switched on by inserting and locking a tip, or by pressing the ENTER key.
- When the locking lever is in the open position, the display shows **CLO** (Close).
- Display: **rEF** (reference run): Press the STEP-key to start reference run. **rEF** is displayed after every replacement of the battery pack.

### - Auto Power Off


- The AutoRep E is factory-set to switch off after 3 minutes of inactivity.
- To switch off manually, press the ENTER key briefly

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
## Inserting tips

### Attention:

**Risk of injury! Do not touch the plunger seat inside the instrument!**

1. Push the locking lever to its upper stop.
2. Slide a tip laterally into the AutoRep E.
3. Push the locking lever down until you hear a distinct clicking (2 steps).
4. If using Encode tips, the tip capacity is recognized by the AutoRep E and its capacity is shown briefly on the display.
5. Then the partial volume to be dispensed in each step is displayed. The default setting is 10 % of the maximum tip capacity.

If using compatible tips (not Encode tips) you must enter the capacity of the tip yourself (see page 31).

### Note:

When using 25 mL and 50 mL Encode Tip, the supplied adapter must be used. It can be connected to the Encode tips using a twist lock and can be removed after use.

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


- If the locking lever is not fully engaged, horizontal bars are shown on the display.
- If the locking lever was fully closed but no tip was inserted, **"no TIP"** appears on the display.

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## Removing tips

1. Hold the AutoRep E over a suitable container.
2. Completely discharge the tip by continuously pressing the STEP key (releasing the residual stroke lock), or by using the "E" mode (see page 23).
3. Push the locking lever to its upper stop.
4. Remove the tip.

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### Attention:

Don't remove Encode-tip before piston of the tip is in the lowest position (tip completely discharged!).

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
### Note:

After dispensing highly viscous or crystallizing sample, the Encode tip may wear and require replacement

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## Functions of the AutoRep™ E

The instrument performs the following functions:

1. Dispensing **DISP** (see page 22): Default mode after switching on. A user-defined volume can be dispensed repeatedly, depending on the capacity of the tip and the filling volume.
2. Automatic Dispensing **AUTO DISP** (see page 24): Identical to DISP mode, but with automatic dispensing by continuously pressing the STEP key. Releasing the STEP key interrupts the process.  
In the AUTO DISP mode the instrument has a "**Learn**" function: After pipetting at least three steps, the AutoRep E calculates the average dispensing interval, and will repeatedly dispense at that interval as long as the STEP key is pressed.
3. Pipetting **PIP** (see page 26): In "PIP" mode, the AutoRep E functions as a positive displacement non-repeating pipette, particularly useful for volatile, foamy or high viscous materials. A preselected volume is aspirated and dispensed.

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### Note:

Default setting for the AutoRep E is DISP mode (see page 22).

1. The arrow on the display points upward ( $\Delta$  "Aspirate").
2. Push the STEP key once. The tip will be filled completely.

### Note:

To interrupt the filling process, press the STEP key. To continue filling, press the  $\Delta$  key first and hold down the STEP key in succession. The step count on the display will increase.

3. During filling, the number of remaining steps is continuously displayed. The arrow on the display points downward ( $\nabla$  "Dispense") when the aspirating is stopped.

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
**Attention:**

After aspirating liquid, the Auto-Rep E will automatically dispense a small amount of sample in order to compensate for gear backlash.

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## DISP mode

1. The arrow on the display points downward (▽ "Dispense").
2. Press the STEP key to dispense the selected step volume (see "Adjusting STEP volume" page 27).
3. The display shows the remaining number of steps.
4. After dispensing the last full step, a small amount of sample will be left in the tip. To discharge this volume, hold down the STEP key.
5. Remove tip (see page 18).

### Note:

The tip must be completely empty before you try to remove it.

### Attention:

Dispensing can be interrupted at any time by pressing the Δ key.

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
#### Note:

To empty the tip at any time, press the ▽ key and then hold down the STEP key. The letter E appears on the display and the arrow pointing downward flashes. This process can be interrupted at any time by pressing the Δ key.

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A partially filled tip can be completely filled by pressing and holding, in sequence, the Δ key and the STEP key.

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⇒ siehe Seite 26.

## AUTO DISP mode

1. To activate the AUTO DISP mode, briefly press the MENU key twice.

### Note:

The AUTO DISP mode must be actively selected.

2. Select the AUTO DISP mode with the  $\nabla\Delta$  keys (ON/OFF).
3. Confirm by pressing ENTER.
4. **"Learn"** feature: After pipetting at least three steps, the instrument automatically determines the average time interval between steps (see page 19).
5. Hold down the STEP key to dispense the selected volume (see page 27) at the calculated intervals.

### Note:

Intervals can range from 0.1 to 5 seconds.

6. To interrupt dispensing, release the STEP key. The calculated step interval remains in memory.
7. In AUTO DISP mode, the "Learn" process can be restarted at any time.

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
## PIP mode

(Reverse pipetting, positive displacement principle)

1. Select the PIP mode by holding down the MENU key (> 3s), and press it briefly (see page 30).

**Make sure that plunger is at lowest position!**

2. Select a volume (page 27).
3. To aspirate the selected volume, press the STEP key. A volume slightly larger than the selected volume will be drawn into the tip.
4. To dispense the volume, press the STEP key again. Some residual volume will remain in the tip.
5. To aspirate a new sample, press the STEP key again.
6. When you have completed pipetting, discharge the residual stroke by holding down the STEP key.

### Attention!

To interrupt the dispensing process at any time, press the  $\Delta$  key.

### Note:

To change from PIP mode to AUTO DISP mode, you have to change to DISP mode first.

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## Adjusting STEP volume

### Note:

After the Encode tip size has been recognized, the dispensing volume is shown on the display. The default step volume is 10 % of the nominal capacity.

1. Press the MENU key once. The volume (VOL) is shown on the display.
2. Use the  $\nabla\Delta$  keys to select the appropriate volume. Holding down the  $\nabla\Delta$  keys will accelerate the rate of change.
3. Confirm the volume selection by pressing ENTER.

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

## Adjusting aspiration and dispensing speed

1. Hold down the MENU key (> 3s) to adjust the aspiration speed. Choose between three different aspiration speeds.

### Note:

#### Lowest speed:

Display shows outline of arrow with one bar.

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
2. Use the  $\nabla\Delta$  keys to select aspiration speed.
3. Confirm aspiration speed setting by pressing ENTER.
4. The down arrow will now flash.
5. Use the  $\nabla\Delta$  keys to select dispensing speed.
6. Confirm settings by pressing ENTER. You are returned to the main menu.

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## Changement de mode

## Cambio de modo

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1. Hold down the MENU key (> 3s) and then press it again briefly. The mode can then be selected.

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- Modes available are PIP an DISP.
- DISP is the default mode.

2. Use the  $\nabla\Delta$  keys to select desired mode.
3. Confirm setting by pressing ENTER.

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


## Use of compatible tips from other manufacturers

Compatible dispenser tips that are not encoded for size (see page 32) can be used with the AutoRep E (see "Insert tips" page 16). After insertion of one of these tips, the display will automatically go to tip selection mode. By default, the 5 mL tip is displayed flashing regardless of the inserted tip. 12 possible tip sizes are available.

1. Use the  $\nabla\Delta$  keys to select appropriate tip size.
2. Confirm setting by pressing ENTER.

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– Fortsetzung –

### Attention!

For convenience, when inserting the next unencoded tip, the previously selected tip size will come up first on the menu. Confirm tip size by pressing ENTER.

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At the time of printing, the following dispenser tips were compatible with the AutoRep E:

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## Residual stroke

A lock prevents the residual stroke from being dispensed as a regular dispensed or pipetted volume. The symbol  $\Delta$  indicates this remaining volume. It can be discharged by holding the STEP key.

## Non-contact dispensing

When tips with volumes of 5 mL or larger are selected, the dispensing speed adjusts to allow accurate dispensing without contacting the walls of vessel (free-jet dispensing). This can be overridden by selecting a different speed (see p. 28).

## Maintenance and cleaning

The AutoRep E is calibrated at the factory and is maintenance-free. If required, clean the instrument with a cloth dampened with water or a mild soap solution.

**The AutoRep E must not be dismantled.**

## Default settings

To revert the AutoRep E to its factory default settings, simultaneously press the MENU and ENTER keys for > 3 seconds, with the plunger in the lowest position.

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## Troubleshooting

Problem	Error message displayed	Possible cause	Corrective action
Tip dripping	No message	Leaking tip	Replace tip
Faint display; reduced speed	Battery symbol flashing	Battery low	Recharge battery pack or replace with charged battery pack.
Instrument not responding	ERR 001	<ul style="list-style-type: none"> <li>- Locking lever not engaged during aspirating or dispensing</li> <li>- Locking lever open at lower third of the tip</li> </ul>	<p>1. Press <b>ENTER</b> to confirm error.</p> <p><b>Follow instructions to avoid instrument damage!</b></p> <p>Display: <b>CLO</b> (Close)</p> <p>2. Engage the locking lever</p> <p>Display: <b>rEF</b> (Reference run)</p> <p>3. Press the STEP key</p> <p><b>Attention: Liquid remaining in tip will be discharged during the reference run (rEF). Discharge sample into a suitable vessel</b> (see pages 4-6).</p>
Instrument not responding	ERR 002	<p>Motor blocked</p> <ul style="list-style-type: none"> <li>- Battery low</li> <li>- Viscosity is too high (see page 7)</li> <li>- Tip orifice obstructed</li> </ul>	<p>1. Press <b>ENTER</b> to confirm error.</p> <p><b>Follow instructions to avoid instrument damage!</b></p> <ul style="list-style-type: none"> <li>- Recharge or replace battery pack</li> <li>- Reduce intake/dispensing speed</li> </ul> <p>Display: <b>rEF</b> (Reference run)</p> <p>2. Press the STEP key</p> <p><b>Attention: Liquid remaining in tip will be discharged during the reference run (rEF). Discharge sample into a suitable vessel</b> (see pages 4-6).</p>
Instrument not responding	ERR TIP	Locking lever open at lower part of the tip	<p>1. Press <b>ENTER</b> to confirm error.</p> <p><b>Follow instructions to avoid instrument damage!</b></p> <p>2. Display: <b>OPE</b> (Open) Push the locking lever to its upper stop.</p>

## Troubleshooting - continued -

Problem	Error message displayed	Possible cause	Corrective action
Instrument not responding - continued -			3. Display alternating between <b>out TIP</b> and <b>CLO</b> . Remove tip and close locking lever (without tip). 4. AutoRep E carries out a reference run. 5. Display: <b>no TIP</b> Insert new tip (see page 16).
Instrument not responding	ERR 003	Electronic fault <b>or</b> tip connector not at lowest position when inserting new tip	1. Press <b>ENTER</b> to confirm error. <b>2. Follow MENU-instructions:</b> - Display rEF (reference run) - Press STEP-Key <b>Attention: Liquid remaining in tip will be discharged during reference run (rEF). Discharge sample into suitable vessel (see page 5, point 10).</b> If error message is still shown on the display remove battery pack and insert it again (reset). If this does not solve the problem send instrument back for repair.
No tone, no light when placing instrument into charging dock	ACC	- Instrument contains no battery - Battery pack inserted improperly - Battery pack faulty	- Insert battery pack. - Make sure battery contacts are aligned with contacts in the instrument. - Replace battery pack.
Charge-indicator is flashing	No message	Removal or insertion of a new battery pack	- Charge battery pack for at least 3.5 hours inside the AutoRep E.
Dispensing process interrupted. <b>Attention!</b> Volume not completely discharged	ERR VOL	Dispensing process interrupted by - Pressing the $\Delta$ key during dispensing - electronic faults	1. Press <b>ENTER</b> to confirm error. 2. Restart dispensing by pressing the <b>STEP</b> key.

### Note:

If other error messages are displayed, remove battery pack and insert it again.  
 If this does not solve the problem, send instrument in for repair.

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## Technical data

### Battery pack:

Rechargeable Nickel Metal Hydride battery pack:  
4 NiMH, AAA cells,  
4,8 V/700 mAh

### Charging times for battery:

in battery compartment of charging dock: approx. 2.5 h

inside the AutoRep E: approx. 3.5 h

When battery capacity starts to diminish, recharge battery pack in the battery compartment of the charging dock.

### Battery capacity:

A full charge gives approx. 2.5 hours of continuous work with a 12.5 mL Encode Tip.

This sign certifies that the product meets the requirements of the EEC directive and has been tested according to the specified test methods.

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### Important:

For safety reasons, only clean and decontaminated instruments will be accepted for inspection or repair.

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- Thoroughly clean and decontaminate the instrument.
- Obtain a Return Goods Authorization number by calling the appropriate Service number listed in page 53.
- Send the AutoRep E in a suitable packing container to the appropriate Service address shown in page 53, describing:
  - The nature of the problem
  - What samples have been used

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## Disposal

The adjoining symbol means that storage batteries and electronic devices must be disposed of separately from household trash (mixed municipal waste) at the end of their service life.


- According to the Directive 2002/96/EC of the European Parliament and of the Council on Waste Electrical and Electronic Equipment (WEEE) of 27 January 2003, electronic equipment requires disposal according to the relevant national disposal regulations.
- Batteries contain substances that can have harmful effects on the environment and human health. Therefore according to the Directive 2006/66/EC of the European Parliament and the Council on Waste Batteries of 6 September 2006 batteries require disposal according to the relevant national disposal regulations. Dispose of batteries only when completely discharged.

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### Warning!

Do not short-circuit the battery to discharge it!

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## Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operation or unauthorized repairs of the instrument or the consequences of normal wear and tear as well as the failure to follow the instructions of the operating manual. We are not liable for damage resulting from any actions not described in the operating manual or if non-original spare parts or components have been used.

Subject to technical modification without notice. Errors excepted.

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



**Autorized Pipette Repair and Calibration facilities:**

California: 7500 Edgewater Drive, Oakland CA 94621  
Tel. 800-662-7027, Fax 510-564-1683  
Massachusetts: Rainin Road, Woburn, MA 01801  
Tel. 800-662-7027, Fax 781-935-7631  
Japan: 4-1-11, Bunkyo-Ku, Tokyo 113-0033  
Tel. (03) 5689-8311, Fax (03) 5689-2670  
Belgium: N.V. Mettler-Toledo s.a., Leuvensesteenweg 384, B -1932 Zaventem  
Tel. (02) 334 02 11, Fax (02) 334 03 34  
Denmark: Mettler-Toledo A/S, Naverland 8, DK - 2600 Glostrup  
Tel. (43) 27 08 00, Fax (43) 27 08 28  
France: HTS, 14, allée Prométhée, ZAC Propylées III, F - 28000 Chartres  
Tel. (02) 37 88 31 00, Fax (02) 37 88 31 09  
Germany: Mettler-Toledo GmbH, Ockerweg 3, D - 35396 Giessen  
Tel. (0) 641 507 307, Fax (0) 641 507 397  
Italy: Mettler-Toledo S.p.A., I-20026 Novate Milanese  
Tel. (02) 333 321, Fax (02) 356 29 73  
Netherlands: Mettler-Toledo B.V., Postbus 6006, NL - 4000 HA Tiel  
Tel. (0) 344 63 81 38, Fax (0) 344 63 83 95  
Norway: Mettler-Toledo AS, Prof. Bikelandsvei 24c, N - 1008 Oslo  
Tlf: 22 30 44 90, Fax: 22 32 70 02  
Spain: Mettler-Toledo S.A.E., L'Hospitalet del Llobregat, E - 08038 Barcelona  
Tel. (93) 223 76 00, Fax (93) 223 67 01  
Sweden: Mettler-Toledo AB, Vikesvägen 10, Södra Hammarbyhamnen, S - 12008  
Stockholm, Tel. (08) 702 50 00, Fax (08) 642 45 62  
Switzerland: Metter-Toledo (Schweiz) AG, Im Langacher, 8606 Greifensee  
Tel. +41 (61) 716 9001, Fax +41 (61) 716 9002  
UK: Anachem Ltd. Anachem House, Charles Street, Luton, Bedfordshire LU2 0EB.  
Tel. +44 (0)1582 747500, Fax +44 (0)1582 745105

Printed in Germany