

**Operating Instructions**  
**Bedienungsanleitung**  
**Mode d'emploi**  
**Instrucciones de manejo**

**METTLER TOLEDO**  
**DP550**  
**DP660**

Operating Instructions .....	4
Bedienungsanleitung .....	8
Mode d'emploi .....	12
Instrucciones de manejo .....	16

## DP550/DP660 Phototrodes

### Order numbers of accessories

Pos. No.	Designation	Order No.
1	Adjustment knob	----
2	Removable cone	22986
3	Window	----
4	Mirror	----
5	Screw-on mirror assembly	22684
6	Signal cable with DIN plug *	----
7	Signal cable with banana plugs *	23557
8	Signal cable with USA plug *	23554
9	Power supply plug	----
10	Power unit, Europe (230 V) *	23555
11	Power unit, USA (120 V) *	23556
	Power unit, UK (240 V) *	23563
	Adapter cable (DIN-LEMO)	89600
	Steel probe special version	999186
	2 mirrors, sealings	999192

\* Choice of equipment

---

## 1. General information

The photoelectric probe “DP550 Phototrode” operates with green light (555 nm), and the “DP660 Phototrode” operates with red light (660 nm). The output signal is 1 V at 100% transmittance. The high-frequency light modulation permits significant suppression of external disturbances caused by sunlight or artificial light. Quite aside from their different designations, the only difference in their external appearance is a red (DP660) or a green (DP550) cap on the adjustment knob.

## 2. Setting up the instrument

- Plug power unit (10 or 11) into an outlet (230 V European plug; 120 V USA plug).
- Insert the small black plug (9) from the power supply cable into the socket of the power unit.
- Connect the signal cable (6, 7 or 8) to the titrator.

## 3. Adjustment

- Connect the Phototrode as described in Paragraph 2.
- Dip the probe into water or in the solvent to be used.
- Lightly depress the adjustment knob (1) and turn it until it snaps into position.
- Adjust the output potential to  $1000 \pm 5$  mV (100% transmittance) by turning the knob which has been snapped in place.
- For precise measurements, it is recommended that the warm-up time of 30 min. be observed.

## 4. Maintenance

Make absolutely sure that the mirror (4) and the window (3) are always clean. It is best to clean them with a cotton swab. The probe is resistant to most acids and solvents. Make sure that the probe is rinsed thoroughly with water after use, especially if alkaline mixtures which contain chloroform are used. If there is any dulling of the mirror, the mirror assembly (5) can be replaced easily (Order No. 22684).

---


## 5. Malfunctions and possible causes

No signal:	<ul style="list-style-type: none"><li>– Mirror (4) or window (3) is dirty.</li><li>– Unscrew mirror assembly, connect Phototrode to power supply and check whether the light is coming through the window (3). If no light is visible, send the phototrode along with the power unit to the nearest METTLER TOLEDO Service facility.</li></ul>
Signal fluctuations:	<ul style="list-style-type: none"><li>– These can be caused by air bubbles or solid particles in the sample solution. Vary the stirrer speed until the air bubbles disappear.</li></ul>

## 6. Technical specifications

Measuring principle:	Single-beam measuring process; chopper-modulated
Light path in sampler:	2 cm
Light wavelength:	DP550: $\lambda = 555$ nm; DP660: $\lambda = 660$ nm
Spectral bandwidth:	Approx. 30 nm (at 50% relative intensity)
Ambient conditions:	Temperature 15° C...35 C Maximum relative atmospheric humidity up to 80%
Maximum sample temperature:	75 °C
Measuring range:	0...100 % transmittance
Measuring accuracy:	Better than 1.5% transmittance
Stability:	After 30 min, warm-up time: $\leq 0.5\%$ transmittance/10 min, (at constant ambient temp.)
Signal output:	100 % transmittance = 1000 mV DC
Time constant:	0.7 s typically
Probe material:	Polypropylene and glass
Power supply:	29 V DC/70 mA from standard equipment power unit
Weight:	Approx. 250 g (not including power unit)

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---

To protect your METTLER TOLEDO product's future:  
METTLER TOLEDO Service assures the quality, measuring accuracy and  
preservation of value of all METTLER TOLEDO products for years to come.  
Please send for full details about our attractive terms of service.

Thank you.

\*P701158\*

Subject to technical changes and to the availability  
of the accessories supplied with the instruments.

© Mettler Toledo GmbH 1990, 1993, 1994, 1996, 1997, 1999, 2000  
ME-701158C Printed in Switzerland 0001/2.33

**Mettler-Toledo GmbH, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach,**

Tel. (01) 806 77 11

Fax (01) 806 7350

Internet: <http://www.mt.com>