


Handheld Density Meter

Densito/DensitoPro



METTLER TOLEDO

Table of Contents

1	Introduction	5
1.1	Further information.....	5
1.2	Explanation of conventions and symbols.....	5
1.3	Compliance information.....	5
2	Safety information	7
2.1	Definitions of signal words and warning symbols.....	7
2.2	Product specific safety notes.....	7
3	Design and function	9
3.1	Instrument overview.....	9
3.2	Overview of functions.....	9
3.3	User interface.....	11
3.3.1	Home screen.....	11
3.3.2	Icons on the screen.....	11
4	Installation and commissioning	12
4.1	Scope of delivery.....	12
4.2	Download the Reference Manual.....	12
4.3	Unpack the density meter.....	13
4.4	Charge the density meter.....	13
4.4.1	Charge with the AC/DC adapter.....	13
4.4.1.1	Assemble the AC/DC adapter.....	14
4.4.1.2	Connect the power supply.....	14
4.4.2	Charge with a computer.....	14
4.4.3	Charge with a charging stand.....	15
4.5	Setup for working with the sample pump.....	15
4.5.1	Activate the sample pump mode.....	15
4.5.2	Install the filling tube.....	16
4.6	Activate and deactivate the RFID reader (DensitoPro only).....	16
4.7	Activate and deactivate the barcode reader (DensitoPro only).....	16
4.8	Install accessories.....	17
4.8.1	Setup for working with a syringe.....	17
4.8.1.1	Activate the syringe mode.....	17
4.8.1.2	Install the connector for the syringe.....	18
4.8.2	Install the charging stand.....	18
4.8.2.1	Assemble the AC/DC adapter.....	19
4.8.2.2	Position the charging stand.....	19
4.8.2.3	Connect the power supply.....	20
4.8.3	Install and remove the protective cover.....	21
4.8.3.1	Install the protective cover.....	21
4.8.3.2	Remove the protective cover.....	21
4.8.4	Connect a printer.....	22
4.8.4.1	Connect a USB printer.....	22
4.8.4.2	Connect a bluetooth printer.....	23
4.8.4.3	Print a test page.....	23
4.8.5	Connect to EasyDirect Density & Refractometry.....	23
4.9	Configure the density meter.....	24
4.9.1	Change the date and time format.....	24
4.9.2	Change the temperature unit.....	24
4.9.3	Change the density unit.....	24
4.9.4	Change energy saver settings.....	25
4.9.5	Protect actions with a password.....	25
4.9.5.1	Enable the password protection.....	26
4.9.5.2	Disable the password protection.....	26

4.9.5.3	Change the password	26
4.9.6	Activate and deactivate the auto rotation of the screen	26
4.9.7	Activate and deactivate the audio signal	26
4.9.8	Change the keyboard layout	27
4.9.9	Change the language	27
4.9.10	Activate and deactivate a cleaning reminder	27
5	Configuration of density determinations	28
5.1	Create, change or delete a method	28
5.1.1	Create a method	28
5.1.2	Change a method	28
5.1.3	Delete a method	29
5.2	Change a test	29
5.3	Change an adjustment	29
5.4	Configure the user guidance	29
5.4.1	Define the level of user guidance	29
5.4.2	Define the pump speed	30
5.4.3	Define the fill cycle (only guided workflow)	30
5.4.4	Define the cleaning type (only guided workflow)	30
5.4.5	Activate and deactivate the confirmation of results	30
5.5	Configure the measurement reliability	31
5.6	Configure the result conversion	31
5.6.1	Density and specific gravity	31
5.6.2	Alcohol concentration	32
5.6.3	American Petroleum Institute (API)	33
5.6.4	Baumé scale	33
5.6.5	Sulfuric acid concentration	33
5.6.6	Sugar concentration	33
5.6.7	User defined concentration	34
5.6.7.1	Define the density or specific gravity	34
5.6.7.2	Define the function	35
5.6.7.3	Define the format of the result	35
5.6.8	Define the temperature-compensation coefficient α	35
5.6.8.1	Enter a known value for α	36
5.6.8.2	The density meter calculates a fixed value for α	36
5.6.8.3	Define the second order polynomial used to calculate α	36
5.7	Configure the automatic documentation of results	37
5.7.1	Configure the sample or standard identifier	37
5.7.2	Configure the automatic printing of results	38
5.7.3	Configure the automatic export of results	38
5.8	Configure result limits	38
5.9	Configure a standard	39
5.10	Configure the tolerance for a test	39
6	Operation	40
6.1	Start up the density meter	40
6.2	Shut down the density meter	40
6.3	Navigate and enter information	41
6.3.1	Navigate through menus and confirm settings	41
6.3.2	Enter text and numbers	41
6.4	Typical phases of density determinations	42
6.4.1	Fill the measuring cell and measure the density	42
6.4.2	Prepare the measuring cell for the next sample	42
6.5	Example: Density determination using the sample pump	42
6.5.1	Configure the method	43
6.5.2	Perform the density determination	44
6.6	Example: Density determination using a syringe	45
6.6.1	Configure the method	45

6.6.2	Perform the density determination	46
6.7	Example: Density determination of cold or hot samples	48
6.7.1	Configure the method	48
6.7.2	Perform the density determination	49
6.8	Stop a method	50
6.9	Read and write information with the RFID reader (DensitoPro only)	51
6.9.1	Write information to a Smart Tag	51
6.9.2	Read information from a Smart Tag	51
6.10	Read information with the barcode reader (DensitoPro only)	52
6.11	View, export, print and delete results	52
6.11.1	View results	52
6.11.2	Delete all results	52
6.11.3	Export and print results	53
6.11.3.1	Export results to EasyDirect Density & Refractometry	53
6.11.3.2	Export results to a USB flash drive	53
6.11.3.3	Print results	54
6.11.3.4	Configure the filter	55
7	Maintenance	56
7.1	Maintenance schedule	56
7.2	Clean the density meter	56
7.2.1	Clean the housing	56
7.2.2	Clean the measuring cell	56
7.2.2.1	Clean using the sample pump	57
7.2.2.2	Clean using a syringe	57
7.2.3	Clean the cylinder and the piston of the sample pump	58
7.2.3.1	Disassemble the cylinder and remove the piston	58
7.2.3.2	Clean the piston and the cylinder sections	59
7.2.3.3	Assemble the cylinder and install the piston	59
7.2.4	Clean the protective cover	59
7.3	Check the measurement accuracy	60
7.3.1	Example: Test with water	60
7.3.1.1	Configure the test	60
7.3.1.2	Perform the test	61
7.3.2	Measures if a test fails	61
7.4	Replace the battery	62
7.5	Replace the filling tube	63
7.6	Replace the cylinder and the piston of the sample pump	63
7.6.1	Remove the cylinder and the piston	63
7.6.2	Install the cylinder and the piston	64
7.7	Update the firmware	64
7.8	Prepare the density meter for storage	64
7.9	Ship the density meter	64
7.10	Dispose of the density meter	64
8	Troubleshooting	65
8.1	List of errors and problems	65
8.2	Drain the measuring cell	66
8.3	Restore the factory-adjustment data	66
8.4	Adjust the measurement accuracy	66
8.4.1	Example: Adjustment with water	67
8.4.1.1	Configure the adjustment	67
8.4.1.2	Perform the adjustment	67
8.5	Restore the factory settings	69
8.6	View the firmware version and other system information	69
9	Technical data	70
9.1	Density meter	70


9.2	Measurement	71
9.3	Charging stand.....	72
9.4	Protective cover	73
10	Accessories	74
10.1	Sampling.....	74
10.2	Printers.....	75
10.3	Miscellaneous accessories.....	75
11	Appendix	77
11.1	Density of pure water (0...40 °C).....	77
11.2	Temperature-compensation coefficients α	77
	Index	79

1 Introduction

Thank you for choosing a METTLER TOLEDO density meter. The handheld density meters Densito and DensitoPro are handheld measuring instruments used to determine the density of liquids.

The instructions in this document refer to the handheld density meters Densito and DensitoPro running firmware version V1.1.0 or higher.

See also

 View the firmware version and other system information ▶ Page 69

1.1 Further information

For applications notes, see the following links:

▶ www.mt.com/analytical-application-library

▶ www.mt.com/library

For third party licenses and open source attribution files, see the following link:

▶ www.mt.com/licenses

If you have any additional questions, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

See also

 Download the Reference Manual ▶ Page 12

1.2 Explanation of conventions and symbols

Note For useful information about the product.



Refers to an external document.

Elements of instructions

Instructions always contain action steps and can contain prerequisites, intermediate results and results. If an instruction contains more than one action step, the action steps are numbered.

- Prerequisites that must be fulfilled before the individual action steps can be executed.

- 1 Action step 1
 - ⇒ Intermediate result
- 2 Action step 2
 - ⇒ Result

1.3 Compliance information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ <http://www.mt.com/ComplianceSearch>

Application-relevant standards and norms are listed on the internet.

▶ www.mt.com/dere-norms

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

▶ www.mt.com/contact

European Union

The instrument complies with the directives and standards listed on the EU Declaration of Conformity. SVHC candidate substances according to Article 33 of the EU regulation no. 1907/2006 (REACH)

Material	CAS No.
1,3-Propanesultone	1120-71-4
PZT (Lead Zirconate Titanate)	12626-81-2

United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The following section only applies to DensitoPro density meters.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Canada

The following section only applies to DensitoPro density meters.

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

2 Safety information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.
- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.



User Manual and Reference Manual are available online. See [Download the Reference Manual ▶ Page 12].

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

WARNING	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
CAUTION	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
NOTICE	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols



General hazard: read the User Manual or the Reference Manual for information about the hazards and the resulting measures.

2.2 Product specific safety notes

Intended use

The density meter is designed to be used by trained staff and is intended for measuring the density of liquid samples. The samples must be compatible with the materials they come into contact with.

The density meter is designed for indoor and outdoor operation in dry conditions. The following site requirements apply:

- Ambient conditions within the limits specified in the technical data
- No corrosive gas atmosphere
- No explosive atmosphere
- No powerful electric or magnetic fields

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes



WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO AC/DC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the plugs for damage and replace damaged cables and plugs.



CAUTION

Injury or damage due to incorrect handling of the battery

The battery can explode or ignite if it is not handled correctly.

- 1 Do not store or use the battery at temperatures below -20 °C or above 60 °C.
- 2 Do not subject the battery to mechanical stresses like pressure, bending or impacts.



NOTICE

Damage to the measuring cell due to solidifying samples

Hot samples cool down during the measurement. If the samples solidify in the measuring cell, you can no longer empty the measuring cell.

- Only measure samples with a viscosity below 1000 mPa*s at ambient temperature.




NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

3 Design and function

3.1 Instrument overview



No.	Name	Function
1	Screen	To display settings and results.
2	Keys	To start up the density meter, select menu items and to navigate.
3	Navigation joystick	To navigate through the menus.
4	Measuring cell window	To view the measuring cell and to check if the measuring cell is filled or empty.
5	Filling opening with locking screw	To fill and drain the measuring cell using a syringe. The locking screw closes off the filling opening.
6	Filling tube	To fill and drain the measuring cell.
7	USB-C socket	To connect the density meter with the AC adapter or a computer.
8	RFID reader and barcode reader (DensitoPro only)	To scan a linear barcode or read and write information to Smart Tags.
9	USB-A socket	To connect a USB flash drive, a USB printer or a bluetooth dongle for a bluetooth printer.
10	Charging contacts	Contacts to charge the density meter with the optional charging stand.
11	Sample pump cover	To access the sample pump.
12	Sample pump window	To view the sample pump and to check if the sample pump is filled or empty.
13	Safety sign	Warns that incorrect handling of the battery can cause injury or damage the density meter.
14	Battery cover	To access the battery.


See also

📖 Accessories ▶ Page 74

3.2 Overview of functions

The density meter uses the oscillating body method as measurement technology. The density meter has no active temperature control.

If you measure samples with low viscosity, you can use the sample pump (1) and the filling tube (2) to fill the measuring cell (5). For samples with high viscosity you can use a syringe (3) and the connector for the syringe (4) to fill the measuring cell (5). The user can observe the filling through the window of the measuring cell.



For repeated measurements, you can define and save workflows as methods. Two workflow types are available:

- Customizable guided workflow: guides the user step by step through the density determination.
- Unguided workflow: the user decides which steps are needed.

The results are automatically converted into one of the following predefined units or a unit that the user defines.

- Density
- Specific gravity
- Alcohol concentration
- American Petroleum Institute (API density, API specific gravity, API degree)
- Baumé scale
- Sulfuric acid concentration
- Sucrose concentration

3.3 User interface

3.3.1 Home screen













No.	Name	Explanation
1	Status ribbon	Displays the current date, time, battery charge status and connection to accessories like EasyDirect Density & Refractometry.
2	Method ribbon	Displays the type of density determination and the method identifier or the selected standard. <ul style="list-style-type: none"> • MS: measurement using a method • TE: test • ADJ: adjustment
3	Sample ID ribbon	Displays the sample ID.
4	Measuring ribbon	Displays the unit, the measured value or result and the temperature of the measuring cell.
5	Sampling ribbon	Displays the current status of the sample pump or syringe.
6	Key functions ribbon	Displays the function of the keys.

3.3.2 Icons on the screen

Icon	Location	Explanation
	Status ribbon	Barcode reader is scanning.
	Status ribbon	RFID reader is reading or writing.
	Status ribbon	EasyDirect Density & Refractometry is connected.
	Status ribbon	Density meter is sending data to the printer.
	Status ribbon	Shows the charge of the battery.
	Status ribbon	Battery is charging.
	Sampling ribbon	Syringe mode activated. The sample pump is deactivated.
	Sampling ribbon	Sample pump activated and piston in the lowermost position.
	Sampling ribbon	Sample pump activated and piston in the uppermost position.

4 Installation and commissioning

4.1 Scope of delivery

Part	Order number	Densito	DensitoPro
 Handheld density meter	–	•	•
 Filling Tube 190 mm • Tube • Connector syringe/filling tube • Washer	30330847	•	•
 Battery Li-Ion 2400 mAh	30330855	•	•
 Power Supply and Worldwide Adapters	30449255	•	•
 USB-C Cable	30449253	•	•
 Density standards (3 pcs) 6 mL	–	•	•
 Tag SmartSample (10 pcs)	30449268	–	•
 User Manual	–	•	•
 Declaration of conformity	–	•	•
 Test report	–	•	•

4.2 Download the Reference Manual

- 1 Go to the website www.mt.com/library.
- 2 Select the **Technical Documentation** tab.
- 3 Enter the product type of your density meter in the search field and start the search.
- 4 Select the Reference Manual from the result list.

- 5 Select the link.
 - ⇒ The Reference Manual is either opened or downloaded depending on the browser settings.
- 6 Check which firmware version is installed on your density meter.
- 7 If the Reference Manual is not written for the installed firmware version, contact your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact

See also

📖 View the firmware version and other system information ► Page 69

4.3 Unpack the density meter

- 1 Unpack the density meter.
- 2 Store the protective packaging for later transport over long distances.
- 3 Check that you have received all parts listed in the scope of delivery.
- 4 Inspect the parts visually for flaws or damage.
- 5 If parts are missing or damaged, report it to your authorized METTLER TOLEDO service representative or dealer.


► www.mt.com/contact

4.4 Charge the density meter

You can either charge the density meter using the supplied AC/DC adapter, a computer or the optional charging stand. Recharging an empty battery takes approximately three hours.

Charging stand

The charging stand can be used on a level surface or mounted on a wall.



4.4.1 Charge with the AC/DC adapter



WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO AC/DC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the plugs for damage and replace damaged cables and plugs.



NOTICE

Damage to the AC/DC adapter due to overheating

An AC/DC adapter that does not have adequate air circulation around it, cannot cool sufficiently and overheats.


- Do not cover the AC/DC adapter.

The AC/DC adapter is suitable for all supply line voltages ranging from 100...240 V AC and 50/60 Hz.

4.4.1.1 Assemble the AC/DC adapter


Install the prongs

- 1 Slide the desired set of prongs (1) into the AC/DC adapter (2).
- 2 Insert the USB-A plug of the USB-A-C cable in the USB-A socket of the AC/DC adapter.




Change the prongs


- 1 Pull the plug of the AC/DC adapter out of the power outlet.
- 2 Push the button (2) and slide the set of prongs (1) out.
- 3 Slide the desired set of prongs into the AC/DC adapter.



4.4.1.2 Connect the power supply


- 1 Insert the USB-C plug of the USB-A-C cable into the USB-C socket on the density meter.
 - 2 Install the cables in such a way that they cannot be damaged or interfere with operation.
 - 3 Insert the plug of the AC/DC adapter into a power outlet that is easily accessible.
- ⇒ The density meter starts charging and the  icon is displayed.


4.4.2 Charge with a computer

- A running computer with USB-A sockets (USB 2.0 or higher) is available.
- 1 Insert the USB-C plug of the USB-A-C cable into the USB-C socket on the density meter.
 - 2 Install the cables in such a way that they cannot be damaged or interfere with operation.
 - 3 Insert the USB-A plug in the USB-A socket of the computer.
- ⇒ The density meter starts charging and the  icon is displayed.


4.4.3 Charge with a charging stand


Place the density meter on a level charging stand

- The charging stand is installed on a level surface.
 - 1 Align the rear of the USB-A socket (1) of the density meter with the knob (2) on the charging stand.
 - 2 Lower the density meter onto the charging stand.
 - 3 Make sure that the charging pins (4) touch the charging contacts (3).
- ⇒ The density meter starts charging and the  icon is displayed.
- ⇒ If the density meter is shutdown, it beeps when it starts charging.



Place the density meter on a wall-mounted charging stand

- The charging stand is mounted on a wall.
 - 1 Align the rear of the USB-A socket (1) of the density meter with the knob (2) on the charging stand.
 - 2 Push the density meter onto the charging stand.
 - 3 Make sure that the charging pins (4) touch the charging contacts (3).
- ⇒ The density meter starts charging and the  icon is displayed.
- ⇒ If the density meter is shutdown, it beeps when it starts charging.




See also

- 📖 Install the charging stand ▶ Page 18


4.5 Setup for working with the sample pump

4.5.1 Activate the sample pump mode


- 1 Press the **Menu** key.
 - 2 Navigate to **Settings > Measurement > Sampling**.
 - ⇒ The **Sampling** window opens.
 - 3 Navigate to **Sample pump**.
 - 4 Press the **OK** key.
 - 5 Navigate back to the home screen.
- ⇒ The icon  is displayed in the sampling ribbon.

4.5.2 Install the filling tube

- The filling tube (tube (1), connector (2) and washer (3)) is assembled.



- The locking screw (2) is installed and tight.
- Screw the connector (1) of the filling tube into the opening (3) and tighten it.



4.6 Activate and deactivate the RFID reader (DensitoPro only)

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > RFID / Barcode > Type**.
⇒ The **Type** window opens.
- 3 To activate the RFID reader, navigate to **RFID**.
- 4 To deactivate the RFID reader, navigate to **None**.
- 5 Press the **OK** key.

See also

- 📖 Read and write information with the RFID reader (DensitoPro only) ▶ Page 51

4.7 Activate and deactivate the barcode reader (DensitoPro only)

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > RFID / Barcode > Type**.
⇒ The **Type** window opens.
- 3 To activate the barcode reader, navigate to **Barcode**.
- 4 To deactivate the barcode reader, navigate to **None**.
- 5 Press the **OK** key.


See also

- 📖 Read information with the barcode reader (DensitoPro only) ▶ Page 52


4.8 Install accessories

4.8.1 Setup for working with a syringe

When you work with a syringe you can place the density meter either on the bench or hold it upright.



4.8.1.1 Activate the syringe mode


- The piston is in the lowermost position.
- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Measurement > Sampling**.
 - ⇒ The **Sampling** window opens.
- 3 Navigate to **Syringe**.
- 4 Press the **OK** key.
- 5 Navigate back to the home screen.
 - ⇒ The icon  is displayed in the sampling ribbon.
- 6 If you want place the density meter on the bench, activate the auto rotation of the screen.

See also

- 📖 Activate and deactivate the auto rotation of the screen ▶ Page 26

4.8.1.2 Install the connector for the syringe


- The connector for the syringe is available.
 - The filling tube is installed.
- 1 Unscrew the locking screw (2) counterclockwise with a flat screwdriver (1).
 - 2 Remove the locking screw.



- 3 Screw the connector (1) clockwise into the filling opening (2) and tighten it with your fingers.

See also

- 📄 Accessories ▶ Page 74
- 📄 Install the filling tube ▶ Page 16



4.8.2 Install the charging stand



⚠ WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO AC/DC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the plugs for damage and replace damaged cables and plugs.



NOTICE

Damage to the AC/DC adapter due to overheating

An AC/DC adapter that does not have adequate air circulation around it, cannot cool sufficiently and overheats.

- Do not cover the AC/DC adapter.

The AC/DC adapter is suitable for all supply line voltages ranging from 100...240 V AC and 50/60 Hz.


See also

 Charge with a charging stand ▶ Page 15

4.8.2.1 Assemble the AC/DC adapter


Install the prongs

- 1 Slide the desired set of prongs (1) into the AC/DC adapter (2).
- 2 Insert the USB-A plug of the USB-A-C cable in the USB-A socket of the AC/DC adapter.



Change the prongs

- 1 Pull the plug of the AC/DC adapter out of the power outlet.
- 2 Push the button (2) and slide the set of prongs (1) out.
- 3 Slide the desired set of prongs into the A/DC adapter.




4.8.2.2 Position the charging stand

The charging stand has been developed for indoor operation.

The following site requirements apply:


- Ambient conditions within the limits specified in the technical data
- No powerful vibrations
- No direct sunlight
- No corrosive gas atmosphere
- No explosive atmosphere
- No powerful electric or magnetic fields

The charging stand can be placed on a level surface or mounted on a wall.





Mount the charging stand on a wall



- The wall and the screws can hold the weight of the charging stand and the density meter.
- There is enough space above the charging stand to connect the USB-A-C cable.



- The screws fit into the mounting holes on the charging stand.
- The screws are mounted on the wall.



- 1 Position the wider part of the mounting holes over the screws (1) and push the charging stand toward the wall.
- 2 Slide the charging stand downwards.




4.8.2.3 Connect the power supply

- The charging stand is installed on a level surface or mounted on a wall.
 - The AC/DC adapter is assembled.
- 1 Insert the USB-C plug of the USB-A-C cable into the USB-C socket (1) on the charging stand.
 - 2 Install the cables in such a way that they cannot be damaged or interfere with operation.
 - 3 Insert the plug of the AC/DC adapter into a power outlet that is easily accessible.

See also

- 📖 Position the charging stand ▶ Page 19
- 📖 Assemble the AC/DC adapter ▶ Page 19




4.8.3 Install and remove the protective cover

The protective cover provides a better grip and protects the density meter against the following:

- Dust
- Dirt
- Moisture
- Bumps
- Scratches

4.8.3.1 Install the protective cover


- The density meter is shut down.
- 1 Remove the filling tube.
 - 2 Remove the connector for the syringe if it is installed.
 - 3 Pull the flaps in the middle apart and slide the density meter into the cover.



- 4 Insert the USB-C seal (1) into the USB-C socket.
- 5 Insert the USB-A seal (2) into the USB-A (3) socket.
- 6 Install the filling tube.
- 7 If needed, install the connector for the syringe.

See also


- 📖 Replace the filling tube ▶ Page 63
- 📖 Install the connector for the syringe ▶ Page 18



4.8.3.2 Remove the protective cover

- The density meter is shut down.
- 1 Remove the filling tube.
 - 2 Remove the connector for the syringe if it is installed.

- 3 Pull the USB-C seal (1) out of the USB-C socket.
- 4 Pull the USB-A seal (2) out of the USB-A socket (3) .
- 5 Push the top of the protective cover (4) off the density meter and slide the density meter out of the protective cover.
- 6 Install the filling tube.
- 7 If needed, install the connector for the syringe.



See also

- 📖 Replace the filling tube ▶ Page 63
- 📖 Install the connector for the syringe ▶ Page 18

4.8.4 Connect a printer

In order to print the results right after the measuring procedure, the density meter has to be connected to a USB or bluetooth printer. Alternatively the results can be printed at a later time from the **Results** menu.

The following printer types are compatible with the density meter:

Connection type	Printer
Bluetooth	Godex MX20
USB	P25

Note

The bluetooth function is only available in selected countries with the necessary approvals. For more information, contact an authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact



Read the user documentation of the USB printer for more information on how to set up and operate the USB printer.

See also

- 📖 Accessories ▶ Page 74

4.8.4.1 Connect a USB printer


- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Printer > Printer type**.
⇒ The **Printer type** window opens.
- 3 Navigate to **USB**.
- 4 Press the **OK** key.
- 5 Insert the USB-B plug of USB-A-B cable into the USB-B socket of the printer.
- 6 Insert the USB-A plug into the USB-A socket of the density meter.
- 7 Switch on the printer.
⇒ The connection is established.

See also

- 📖 Configure the automatic printing of results ▶ Page 38
- 📖 Print results ▶ Page 54

4.8.4.2 Connect a bluetooth printer

- The bluetooth printer is set up and running.
- 1 Insert the bluetooth dongle into the USB-A socket of the density meter.
- 2 Press the **Menu** key.
- 3 Navigate to **Settings > Printer > Printer type**.
 - ⇒ The **Printer type** window opens.
- 4 Navigate to **Bluetooth**.
- 5 Press the **OK** key.
- 6 Navigate to **Printer selection**.
 - ⇒ The density meter searches for bluetooth printers.
 - ⇒ A list of available bluetooth printers is displayed.
- 7 Navigate to the bluetooth printer you want to connect.
- 8 Press the **OK** key.



See also

- 📖 Configure the automatic printing of results ▶ Page 38
- 📖 Print results ▶ Page 54

4.8.4.3 Print a test page

- The printer is connected and running.
- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Printer > Test printer**.
 - ⇒ If there is no connection yet, the notification **Connecting to printer** is displayed.
 - ⇒ When the connection is established, a test page is printed.

See also

- 📖 Connect a printer ▶ Page 22


4.8.5 Connect to EasyDirect Density & Refractometry

The EasyDirect Density & Refractometry computer software is a data management software to manage a database of measurement results. The data can be directly exported to the computer software if a computer is connected, on which the EasyDirect Density & Refractometry computer software is installed.



Read the electronic help of the EasyDirect Density & Refractometry software for more information on how to operate the computer software.

Procedure

- EasyDirect Density & Refractometry is installed on the computer.
- A USB-A-C cable is available.
- 1 Insert the USB-C plug of the USB-A-C cable into the USB-C socket on the density meter.
- 2 Insert the USB-A plug into a USB-A socket on the computer.
- 3 Add the density meter to the list of instruments on EasyDirect Density & Refractometry.
 - ⇒ The icon  is displayed in the status ribbon.

See also

- 📖 Export results to EasyDirect Density & Refractometry ▶ Page 53

4.9 Configure the density meter

4.9.1 Change the date and time format

The following date formats are available:

- **dd/mm/yyyy**: day/month/year, e.g. 30/11/2020
- **mm/dd/yyyy**: month/day/year, e.g. 11/30/2020
- **yyyy/mm/dd**: year/month/day, e.g. 2020/11/30
- **dd.mm.yyyy**: day.month.year, e.g. 30.11.2020

The following two time formats are available.

- **24 h**: 24-hour clock, e.g. 6:30 and 18:30
- **am/pm**: 12-hour clock, e.g. 6:30 am and 6:30 pm

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Date / Time**.
⇒ The **Date / Time** window opens.
- 3 Navigate to **Date** or **Time**.
- 4 To change the date or time, move the navigation joystick up to raise or down to lower the number.
- 5 To switch between day, month and year or hours and minutes, move the navigation joystick to the left or right.
- 6 Press the **OK** key.
- 7 Navigate to **Date format** or **Time format**.
- 8 Navigate to the format you want to use.
- 9 Press the **OK** key.

4.9.2 Change the temperature unit

Two temperature units are available.

- **°C**
- **°F**

When you change the temperature unit, temperatures are displayed in the selected unit. The temperature unit of results that were determined before you changed the unit remain unchanged.

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Global units > Temperature unit**.
⇒ The **Temperature unit** window opens.
- 3 Navigate to the desired **Temperature unit**.
- 4 Press the **OK** key.
⇒ The temperatures are displayed in the selected unit.

4.9.3 Change the density unit

Three density units are available.

- **g/cm³**
- **kg/m³**
- **lb/gal (US)**

When you change the density unit, density values are displayed in the selected unit. There are two exceptions.

- The densities of results that were determined before you changed the unit remain unchanged.

- If you select the unit [lb/gal], the unit [g/cm³] is still used for adjustments and tests.

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Global units > Density unit**.
⇒ The **Density unit** window opens.
- 3 Navigate to the desired **Density unit**.
- 4 Press the **OK** key.
⇒ The densities are displayed in the selected unit

4.9.4 Change energy saver settings

The settings of the parameters **Brightness**, **Energy saver** and **Auto shutdown** influence the power consumption and thus the battery life.

Parameter	Description	Values
Brightness	Defines the brightness of the screen.	10...100 %
Energy saver	Defines if the screen is dimmed when no action has been performed for a given time.	Off 30 60 90 120 seconds
Auto shutdown	Defines if the density meter is shut down when no action has been performed for a given time.	Off 2 5 10 minutes

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Screen**.
⇒ The **Screen** window opens.
- 3 Navigate to the desired parameter and change the settings.

4.9.5 Protect actions with a password

With the parameter **Password** you can prevent users from performing certain actions without a password. The following table shows which actions are password-protected.

Action	Password protected
Start methods.	No
Start tests.	No
Start adjustments.	Yes
View the configuration of methods, tests and adjustments.	No
Change the configuration of methods, tests and adjustments.	Yes
Create and delete methods.	Yes
View, print and export results.	No
Delete all results.	Yes
View the settings of the density meter.	No
Change the settings of the density meter.	Yes
Update the firmware.	Yes
Restore factory adjustment.	Yes
Restore factory settings.	Yes

The password has to be entered once before a user performs a password-protected action and does not have to be reentered until the density meter is restarted.

If the user forgets the password, the universal password 8606 unlocks all actions and settings.

4.9.5.1 Enable the password protection

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Password > Password protection**.
⇒ The **Password protection** window opens.
- 3 Navigate to **Yes**.
- 4 Press the **OK** key.
⇒ The **Password** window opens.
- 5 Enter a password with the navigation joystick.
- 6 Press the **OK** key.
- 7 Reenter the password.
- 8 Press the **OK** key.
⇒ The password has to be entered before the user can perform a protected action.

4.9.5.2 Disable the password protection

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Password > Password protection**.
⇒ The **Password protection** window opens.
- 3 Navigate to **No**.
- 4 Press the **OK** key.
⇒ The **Password** window opens.
- 5 Enter the password with the navigation joystick.
- 6 Press the **OK** key.
⇒ Users can perform all actions.

4.9.5.3 Change the password

- 1 Disable the password protection.
- 2 Enable the password protection and enter the new password.

4.9.6 Activate and deactivate the auto rotation of the screen

The **Auto rotate** parameter defines whether the screen rotates by 180° when the density meter is placed on a bench.

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Screen > Auto rotate**.
⇒ The **Auto rotate** window opens.
- 3 Navigate to **Yes** to switch on or **No** to switch off the **Auto rotate** mode.
- 4 Press the **OK** key.
⇒ The screen rotates when the density meter is placed on a bench.

4.9.7 Activate and deactivate the audio signal

The parameter **Audio signal** defines whether an audio signal is played when a key or the navigation joystick is pressed and when measurements are finished.

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Audio signal**.
⇒ The **Audio signal** window opens.
- 3 Navigate to **Yes** to switch on or **No** to switch off the **Audio signal**.

- 4 Press the **OK** key.

4.9.8 Change the keyboard layout

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Language / Keyboard > Keyboard layout**.
⇒ The **Keyboard layout** window opens.
- 3 Navigate to the desired keyboard layout.
- 4 Press the **OK** key.

4.9.9 Change the language

The system language defines the language of the user interface and printouts.

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Language / Keyboard > Language**.
⇒ The **Language** window opens.
- 3 Navigate to the desired language.
- 4 Press the **OK** key.
⇒ The user interface is changed to the selected language.

4.9.10 Activate and deactivate a cleaning reminder

Because the measuring cell always contains residue from the last sample, there is a risk that the measuring cell can be damaged by the sample. To prevent such damages, METTLER TOLEDO recommends that you clean the measuring cell at the end of each sample series.

To remind users to clean the measuring cell, you can activate a reminder. This reminder is displayed before manual shutdown of the instrument, and users have to confirm or cancel the reminder before shutting down the density meter.

Action of the system	Cleaning reminder	Auto shutdown
When users shut down the density meter, a reminder is displayed.	Yes	No
If the density meter automatically shuts down, no reminder is displayed.	Yes	Yes
No reminder is displayed when the density meter shuts down.	No	Yes/No

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Settings > Measurement > Cleaning reminder**.
⇒ The **Cleaning reminder** window opens.
- 3 Navigate to **Yes** to activate or **No** to deactivate the **Cleaning reminder**.
- 4 Press the **OK** key.

See also

- 📖 Change energy saver settings ▶ Page 25

5 Configuration of density determinations

Three different types of density determinations are available.

- **Methods:** determine the density of a sample.
- **Test:** determine the density of a standard to check that the measurement cell is clean and measuring accurately.
- **Adjustment:** determine the density of a standard and adjust the density meter based on the results.

5.1 Create, change or delete a method

5.1.1 Create a method

The new method is listed in the **Methods** window with the name you enter in the parameter **Method ID**. If you do not change the parameter **Method ID**, it is set by default to (A{X}{X}), e.g. A03.

- Less than 30 methods are saved.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
⇒ The **Methods** window opens.
 - 2 Press the **New** key.
⇒ The **Method settings** window opens.
 - 3 Navigate to **Method ID** and enter a name with up to 25 characters.
 - 4 Press the **OK** key.
⇒ The method is listed in the **Methods** window.
 - 5 Change the parameters as needed.

See also

- 📖 Configure the user guidance ▶ Page 29
- 📖 Configure the measurement reliability ▶ Page 31
- 📖 Configure the result conversion ▶ Page 31
- 📖 Configure the automatic documentation of results ▶ Page 37
- 📖 Configure result limits ▶ Page 38

5.1.2 Change a method

- The home screen is open.
- 1 Move the navigation joystick to the left.
⇒ The **Methods** window opens.
 - 2 Navigate to the method you want to edit.
⇒ The **Method settings** window opens.
 - 3 Change the parameters as needed.

See also

- 📖 Configure the user guidance ▶ Page 29
- 📖 Configure the measurement reliability ▶ Page 31
- 📖 Configure the result conversion ▶ Page 31
- 📖 Configure the automatic documentation of results ▶ Page 37
- 📖 Configure result limits ▶ Page 38

5.1.3 Delete a method

- More than one method is created.
 - The method you want to delete is not active. The active method is marked by a square.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
⇒ The **Methods** window opens.
 - 2 Navigate to the method you want to delete.
⇒ The **Method settings** window opens.
 - 3 Press the **Delete** key.
⇒ The message **Do you want to delete the method?** is displayed.
 - 4 Navigate to **Yes**.
 - 5 Press the **OK** key.
⇒ The method is no longer listed in the **Methods** window.

5.2 Change a test

- 1 Press the **Menu** key.
- 2 Navigate to **Test**.
⇒ The **Test** window opens.
- 3 Change the parameters as needed.

See also

- 📖 Configure the user guidance ▶ Page 29
- 📖 Configure the automatic documentation of results ▶ Page 37
- 📖 Configure a standard ▶ Page 39
- 📖 Configure the tolerance for a test ▶ Page 39

5.3 Change an adjustment

- 1 Press the **Menu** key.
- 2 Navigate to **Adjustment**.
⇒ The **Adjustment** window opens.
- 3 Change the parameters as needed.

See also

- 📖 Configure the user guidance ▶ Page 29
- 📖 Configure the automatic documentation of results ▶ Page 37
- 📖 Configure a standard ▶ Page 39

5.4 Configure the user guidance

The parameters that are available to configure the user guidance are described in the following chapters.

5.4.1 Define the level of user guidance

The **Workflow** parameter defines the level of user guidance during the density determination. Two options are available.

- **Manual:** users decide how and when to fill the measuring cell, start the measurement and drain or clean the measuring cell.
- **Guided:** the density meter prompts the user to fill, drain or clean the measuring cell.

5.4.2 Define the pump speed

The parameter **Fill speed** defines the speed of the pump for filling the measuring cell. The sample viscosity places limits on the fill speed.

Speed level	% of maximum pump speed	Fill time for water	Viscosity
High	80 %	Approx. 4 s	Up to 40 mPa*s
Medium	40 %	Approx. 8 s	Up to 75 mPa*s
Low	10 %	Approx. 30 s	Up to 100 mPa*s

5.4.3 Define the fill cycle (only guided workflow)

The parameter **Fill cycle** defines, how many times the measuring cell is filled and emptied before the density is measured. The sample that is in the measuring cell after the last filling is used for the measurement.

Recommendation:

- Density determination with the sample pump: set **Fill cycle** to 3.
- Density determination with a syringe: set **Fill cycle** to 2.

5.4.4 Define the cleaning type (only guided workflow)

Two types of cleaning methods are available, **Drain** and **Rinse**.

Drain

The user is prompted to drain the measuring cell after the measurement.

This setting is useful if the samples that are measured dissolve in each other.

Rinse

The user is prompted to fill and drain the measuring cell with a cleaning solution after the measurement. The parameter **Rinse cycle** defines how many times the user is prompted to fill and drain the measuring cell.

This setting is only useful if the samples that are measured do not dissolve in each other.

See also

 Typical phases of density determinations ▶ Page 42

5.4.5 Activate and deactivate the confirmation of results

After the measurement reliability is achieved, the results are frozen. The parameter **Confirmation of results** defines whether users have to confirm the results before they are replaced by the current density value.

Regardless of the setting, all results are stored.

The setting is not method specific but valid for all density determinations.

Action of the system	Confirmation of results	Workflow
The screen is frozen and the sample pump is locked until users confirm the result.	Yes	Guided
The screen is frozen until users confirm the result but the sample pump is not locked.	Yes	Manual
Users are prompted to place the filling tube into the waste beaker. As soon as users press OK , the result is replaced by the current density value.	No	Guided
As soon as the user starts draining the cell, the result is replaced by the current density value.	No	Manual

Procedure

- 1 Press the **Menu** key.

- 2 Navigate to **Settings > Measurement > Confirmation of results**.
⇒ The **Confirmation of results** window opens.
- 3 Navigate to **Yes** to activate or **No** to deactivate the **Confirmation of results**.
- 4 Press the **OK** key.

See also

 Define the level of user guidance ▶ Page 29

5.5 Configure the measurement reliability

Measurement reliability is only available for methods.

The parameter **Measurement reliability** defines how the stability of a measured value is evaluated. Once a measured value is determined to be stable, it is saved as a result.

- **Medium**: The result is saved as soon as the measured temperature value stays within 0.4 °C (0.72 °F) for 10 s. Results are delivered faster, but with less accuracy compared to the **High** mode.
- **High**: The result is saved as soon as the measured temperature value stays within 0.2 °C (0.36 °F) for 10 s. Results are very accurate, but may take a longer time compared to the **Medium** mode in cases where the sample temperature differs from the ambient temperature.
- **Manual** (only for unguided density determinations): users decide when the result is stable and save it.

See also

 Define the level of user guidance ▶ Page 29

5.6 Configure the result conversion

The density meter converts the measured oscillation to the unit you define. The conversions are grouped in the following categories.

- **Density**: Density and specific gravity
- **Alcohol**: Alcohol concentration
- **API**: American Petroleum Institute (API density, API specific gravity, API degree)
- **Baumé**: Baumé scale
- **H2SO4**: Sulfuric acid concentration
- **Sugar**: Sucrose concentration
- **User-defined**: Conversion as defined by the user

Select the category and the conversion type

- 1 Press the **Menu** key.
 - 2 Navigate to **Methods > Method settings > Calculation > Category**
 - 3 Navigate to the desired category.
 - 4 Press the **OK** key.
⇒ The **Calculation** window opens.
 - 5 Navigate to the conversion type you want to use.
 - 6 Press the **OK** key.
- ⇒ The parameters that you need to configure for the selected conversion type are displayed.

You can find more information about the conversion types and the parameters in the following chapters.

5.6.1 Density and specific gravity

Available conversions:

- **d**: density at the measurement temperature
- **d comp.**: temperature compensated density

- **SG**: temperature compensated specific gravity
- **SG (t/t)**: specific gravity at the measurement temperature

Density at the measurement temperature

Unit: [g/cm³], [kg/m³] or [lb/gal (US)] as defined in **Settings > Global units**.

The calculation is based on a reference for the density of water.

► www.mt.com/dere-norms

Temperature compensated density

Unit: [g/cm³], [kg/m³] or [lb/gal (US)] as defined in **Settings > Global units**.

The measured density is converted to the value it would have at a compensation temperature. The following formula is used to compensate the temperature influence.

$$d_{\text{comp}} = d_{\text{meas}} * (1 + \alpha * (T_{\text{meas}} - T_{\text{comp}}))$$

- d_{comp} : temperature compensated density
- d_{meas} : density at the measurement temperature
- α : temperature coefficient as defined in **Alpha*1000 > Alpha entry**, see [Define the temperature-compensation coefficient α ► Page 35]
- T_{meas} : measurement temperature, temperature of the measuring cell during the measurement
- T_{comp} : compensation temperature as defined in **Compensation temp.**

Temperature compensated specific gravity

$$SG = \frac{d_{\text{comp}}(\text{Sample})}{d_{\text{comp}}(\text{H}_2\text{O})}$$

- SG: specific gravity
- $d_{\text{comp}}(\text{Sample})$: density of the sample compensated to the temperature defined in **Compensation temp. 1**. The temperature influence is compensated by the temperature coefficient α as defined in **Alpha*1000 > Alpha entry**. See [Define the temperature-compensation coefficient α ► Page 35]
- $d_{\text{comp}}(\text{H}_2\text{O})$: density of water compensated to the temperature defined in **Compensation temp. 2**

Specific gravity at the measurement temperature

$$SG(t/t) = \frac{d(\text{Sample})}{d(\text{H}_2\text{O})}$$

- SG(t/t): specific gravity
- $d(\text{Sample})$: density of the sample at the measurement temperature.
- $d(\text{H}_2\text{O})$: density of water compensated to the measurement temperature.

See also

- 📖 Define the temperature-compensation coefficient α ► Page 35
- 📖 Configure the result conversion ► Page 31
- 📖 Density of pure water (0...40 °C) ► Page 77
- 📖 Change the density unit ► Page 24

5.6.2 Alcohol concentration

Available calculations:

- **Alcohol [% w/w]**: mass percent [%w/w]
- **Alcohol 60°F [% v/v]**: volume percent [%v/v], compensated to 60 °F
- **Alcohol 20°C [% v/v]**: volume percent [%v/v], compensated to 20 °C
- **Proof US 60°F [Proof]**: Proof degree with US unit [Proof], compensated to 60 °F

- **Proof IP 60°F [Proof]**: Proof degree with IP unit [Proof], compensated to 60 °F

See also

 Configure the result conversion ▶ Page 31

5.6.3 American Petroleum Institute (API)

The conversion of the measured value to the selected unit is based on the tables published by the American Petroleum institute (API).

The following tables are available.

- **Crude oil (A)**
- **Refined products (B)**
- **Lubricants (D)**

The **Compensation temp.** parameter defines the compensation temperature.

The **Output format** defines to which of the following units the result is converted:

- **d**: API density
- **SG**: API specific gravity
- **°API**: API degree

See also

 Configure the result conversion ▶ Page 31

5.6.4 Baumé scale

Unit: [°Bé]

Available conversions:

- **H. Bé 15°C [°Bé]**: for densities > 1 g/cm³, compensated to 15 °C
- **L. Bé 15°C [°Bé]**: for densities < 1 g/cm³, compensated to 15 °C
- **H. Bé 60°F [°Bé]**: for densities > 1 g/cm³, compensated to 60 °F
- **L. Bé 60°F [°Bé]**: for densities < 1 g/cm³, compensated to 60 °F

The temperature influence is compensated by the temperature coefficient α as defined in **Alpha*1000 > Alpha entry**.

See also

 Define the temperature-compensation coefficient α ▶ Page 35

 Configure the result conversion ▶ Page 31

5.6.5 Sulfuric acid concentration

Available conversions:

- **H2SO4 [% w/w]**: mass percent [%w/w]
- **H2SO4 20°C [% v/v]**: volume percent [%v/v], compensated to 20 °C

See also

 Configure the result conversion ▶ Page 31

5.6.6 Sugar concentration

Available conversions:

- **Brix [% w/w]**: mass percent [w/w %], compensated to 20 °C
- **Plato [°P]**: degree Plato [°P], compensated to 20 °C

See also

📖 Configure the result conversion ▶ Page 31

5.6.7 User defined concentration

You can define a function to convert a density or specific gravity to a concentration. The following steps are involved:

1. Define the conversion of the measured value to a density or specific gravity.
2. Define the function to convert the density or specific gravity to a concentration.
3. Define the unit and the number of decimal places of the concentration.

See also

📖 Define the density or specific gravity ▶ Page 34

📖 Define the function ▶ Page 35

📖 Define the format of the result ▶ Page 35

📖 Configure the result conversion ▶ Page 31

5.6.7.1 Define the density or specific gravity

The conversions are either based on temperature compensated densities (**d comp.**) or specific gravities (**SG**).

Define a temperature-compensated density

The following formula is used to calculate the temperature-compensated density.

$$d_{\text{comp}} = d_{\text{meas}} * (1 + \alpha * (T_{\text{meas}} - T_{\text{comp}}))$$

- d_{comp} : temperature compensated density
- d_{meas} : density at the measurement temperature
- α : temperature coefficient as defined in **Alpha*1000 > Alpha entry**, see [Define the temperature-compensation coefficient α ▶ Page 35]
- T_{meas} : measurement temperature, temperature of the measuring cell during the measurement
- T_{comp} : compensation temperature as defined in **Compensation temp.**

Steps to define a temperature-compensated density:

- The conversion category is set to **User-defined**.

- 1 Navigate to **X** and select one of the following options: **d comp.**, **d comp. - 1** or **(1/d comp.) - 1**.
- 2 Navigate to **Compensation temp.** and define the compensation temperature.
- 3 Navigate to **Alpha*1000 > Alpha entry** and define α . See [Define the temperature-compensation coefficient α ▶ Page 35]

Define a temperature-compensated specific gravity

The following formula is used to calculate the temperature-compensated specific gravity.

$$SG = \frac{d_{\text{comp}}(\text{Sample})}{d_{\text{comp}}(\text{H}_2\text{O})}$$

- SG: specific gravity
- $d_{\text{comp}}(\text{Sample})$: density of the sample compensated to the temperature defined in **Compensation temp. 1**. The temperature influence is compensated by the temperature coefficient α as defined in **Alpha*1000 > Alpha entry**. See [Define the temperature-compensation coefficient α ▶ Page 35]
- $d_{\text{comp}}(\text{H}_2\text{O})$: density of water compensated to the temperature defined in **Compensation temp. 2**

Steps to define the temperature-compensated specific gravity:

- The conversion category is set to **User-defined**.

- 1 Navigate to **X** and select one of the following options: **SG**, **SG-1** or **(1/SG)-1**.
- 2 Navigate to **Compensation temp. 1** and define the compensation temperature for the sample.

- 3 Navigate to **Compensation temp. 2** and define the compensation temperature for water.
- 4 Navigate to **Alpha*1000 > Alpha entry** and define α . See [Define the temperature-compensation coefficient α ▶ Page 35]

See also

 Define the function ▶ Page 35

5.6.7.2 Define the function

You can either use a linear function (**Linear (Ax+B)**) or a second order polynomial function (**2nd order pol. (Ax²+Bx+C)**). The functions are based on a graph with the density or specific gravity on the x-axis and the concentration on the y-axis.

The following formulae are used.

- Linear function: $Ax + B$
- Second order polynomial function: $Ax^2 + Bx + C$


Procedure

- The coefficients A and B are known for a linear function.
 - The coefficients A, B and C are known for a second order polynomial function.
- 1 Navigate to **Fit type** and select **Linear (Ax+B)** or **2nd order pol. (Ax²+Bx+C)**.
 - 2 Press the **OK** key.
 - 3 Enter the values for the coefficients.

Calculate the coefficients A, B and C

- 1 If you want to use a second order polynomial function, prepare at least 4 solutions with different concentrations.
- 2 If you want to use a linear function, prepare at least 2 solutions with different concentrations.
- 3 Determine the density or specific gravity you want to use for each concentration.
- 4 Plot the concentration (y-axis) against the density or specific gravity (x-axis).
- 5 Fit a line or second order polynomial to the curve.

See also

 Define the density or specific gravity ▶ Page 34

5.6.7.3 Define the format of the result




Using the parameters **Unit** and **Decimal places**, you can define the unit and the number of decimal places for the display of the current value and the result.

5.6.8 Define the temperature-compensation coefficient α

The temperature-compensation coefficient α is used to convert a measured value to a reference temperature. You have three options to define α .

- **Manual:** users enter the value for α that they calculated or read from a table.
- **Calculated:** the density meter calculates α based on a linear equation and densities and temperatures that users enter.
- **2nd order pol. (ax²+bx+c):** the density meter calculates α based on a second order polynomial using coefficients that users enter.

See also

-  Density and specific gravity ▶ Page 31
-  Baumé scale ▶ Page 33
-  Define the density or specific gravity ▶ Page 34

5.6.8.1 Enter a known value for α

- The result of $\alpha \cdot 1000$ is known.
- 1 Navigate to **Alpha*1000** > **Alpha entry** and select **Manual**.
 - 2 Press the **OK** key.
 - 3 Navigate to **Alpha*1000**.
 - 4 Enter the value for $\alpha \cdot 1000$.
 - 5 Press the **OK** key.

Calculate α

- 1 Measure the density (d_1) at a temperature (T_1) above the normal measuring temperature.
- 2 Measure the density (d_2) at a temperature (T_2) below the normal measuring temperature.
- 3 Use the measured densities and the temperatures in [$^{\circ}\text{C}$] to calculate α according to the following formula.
$$\alpha = \left(\frac{d_2}{d_1} - 1 \right) / (T_1 - T_2)$$
- 4 Calculate the result of $\alpha \cdot 1000$.

Example

Density at 26 $^{\circ}\text{C}$ (T_1) = 0.7844 g/cm³

Density at 15 $^{\circ}\text{C}$ (T_2) = 0.7937 g/cm³

$$\alpha = \left(\frac{0.7937}{0.7844} - 1 \right) / (26 - 15)$$

$$\alpha = ((1.011856 - 1)/11) = 0.011856/11 = 0.001078$$

$$\alpha \cdot 1000 = 1.078$$

See also

📖 Temperature-compensation coefficients α ▶ Page 77

5.6.8.2 The density meter calculates a fixed value for α

- The density (d_1) at a temperature (T_1) above the normal measuring temperature is known.
 - The density (d_2) at a temperature (T_2) below the normal measuring temperature is known.
- 1 Navigate to **Alpha*1000** > **Alpha entry** and select **Calculated**.
 - 2 Navigate to **Density 1** and enter the value for d_1 .
 - 3 Navigate to **Temperature 1** and enter the value for T_1 .
 - 4 Navigate to **Density 2** and enter the value for d_2 .
 - 5 Navigate to **Temperature 2** and enter the value for T_2 .
- ⇒ The density meter calculates the value for **Alpha*1000** and displays it.

See also

📖 List of errors and problems ▶ Page 65

5.6.8.3 Define the second order polynomial used to calculate α

The density meter can use the following second order polynomial to calculate α .

$$\alpha = (a \cdot d_{\text{meas}}^2 + b \cdot d_{\text{meas}} + c)$$

Procedure

- The three coefficients a , b and c are known.
- 1 Navigate to **Alpha*1000** > **Alpha entry** and select **2nd order pol. ($\alpha x^2 + bx + c$)**.
 - 2 Press the **OK** key.
 - 3 Enter the values for the coefficients a , b and c .

Calculate the coefficients a, b and c

- 1 Prepare at least 4 solutions with different concentrations.
- 2 Define two temperatures, one above and one below the normal measurement temperature.
- 3 Measure the density of each solution at both temperatures.
- 4 Calculate the temperature-compensation coefficient α for each solution.
- 5 Plot the values for α (y-axis) against the densities (x-axis) measured at one of the temperatures.
- 6 Fit a second order polynomial to the curve.

See also

 List of errors and problems ▶ Page 65

5.7 Configure the automatic documentation of results



NOTICE

Loss of data

A maximum of 1100 results is saved in the density meter. When the result list contains 1100 results and a new results is saved, the oldest result is overwritten.

- Print or export the results before the maximum number of results is reached.

See also

 Change a method ▶ Page 28

 Change a test ▶ Page 29

 Change an adjustment ▶ Page 29

5.7.1 Configure the sample or standard identifier

You can define what kind of identifier is assigned to a sample or to a standard.

Automatic numbering is only available for methods.

Action of the system	Sample ID entry Standard ID entry	Sample ID Standard ID
No identifier is assigned to the sample or the standard.	None	Not available
<ul style="list-style-type: none">• The identifier you define in the Sample ID parameter is assigned to all samples measured with the same method.• The identifier you define in the Standard ID parameter is assigned to all standards measured with the same test.	Fixed	0...25 characters
Users enter the identifier before the density measurement starts.	Variable	Not available
<ul style="list-style-type: none">• Guided density determination: users are prompted to enter the identifier when the method starts. The density determination continues when the user has entered a valid identifier.• Unguided density determination: Users have to enter the identifier before the measurement starts but are not prompted to do so. If they do not enter an identifier, the identifier of the previous sample is used.		

Action of the system	Sample ID entry Standard ID entry	Sample ID Standard ID
An automatically generated identifier is assigned to the sample. The identifier automatically increases by 1, when a new measurement is started with the same method.	Automatic numbering	Not available

5.7.2 Configure the automatic printing of results

The **Printout** parameter defines whether a result is automatically printed when it is saved and which information is printed.

- **None:** No report is printed.
- **Short:** A short report is printed.
- **Long:** A long report is printed

Information	Methods	Test	Adjustment
Date / Time	Short/Long	Short/Long	Short/Long
Method ID	Short/Long	–	–
Sample ID	Short/Long	–	–
Application name with standard	–	Short/Long	Short/Long
Standard ID	–	Short/Long	Short/Long
Status	Short/Long	Short/Long	–
Calculation	Short/Long	–	–
Density	Short/Long	Short/Long	Short/Long
Dev. d	–	Short/Long	Short/Long
Temp.	Short/Long	Short/Long	Short/Long
d nominal	–	Short/Long	Short/Long
SNR	Long	Long	Long
Additional settings like the compensation temperature and the temperature-compensation coefficient α	Long	Long	Long

See also

 Connect a printer ▶ Page 22

5.7.3 Configure the automatic export of results

Results are automatically exported to EasyDirect Density & Refractometry if the density meter is connected to a computer with the software installed.

See also

 Connect to EasyDirect Density & Refractometry ▶ Page 23

5.8 Configure result limits

Result limits are only available for methods.

The **Result limits** parameter defines whether the density meter evaluates if a result lies within a defined range. The parameters **Lower limit** and **Upper limit** define the range.

- Result within the limits: the background of the measuring ribbon turns green.

- Result outside the limits: the background of the measuring ribbon turns red and the text **Out of limits** is displayed in the sampling ribbon

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Methods > Method settings > Calculation > Result limits**.
- 3 To deactivate **Result limits**, navigate to **No**.
- 4 To activate **Result limits**, navigate to **Yes**.
- 5 Press the **OK** key.
 - ⇒ **Lower limit** and **Upper limit** are displayed.
- 6 Set the **Lower limit** and **Upper limit**.

5.9 Configure a standard

Use water as standard

- The **Test** window or the **Adjustment** window are open.
- 1 Navigate to **Standard type** and select **Water**.
 - 2 Press the **OK** key.
- ⇒ The density meter uses predefined values for the **Test** or the **Adjustment**.

Use another substance as standard

- The **Test** window or the **Adjustment** window are open.
- 1 Navigate to **Standard type** and select **Other standard**.
 - 2 Press the **OK** key.
 - 3 Navigate to **Compensation temp.** and enter the compensation temperature.
 - 4 Navigate to **d nominal** and enter the nominal density of the standard at the temperature you entered as compensation temperature.
 - 5 Navigate to **Alpha*1000** and define the temperature-compensation coefficient α .

See also

- 📖 Define the temperature-compensation coefficient α ▶ Page 35

5.10 Configure the tolerance for a test

The **Tolerance d** parameter is used to evaluate the measured density and decide if the test passed or failed. Recommendations for defining the value for **Tolerance d**:


- Deionized water: $2 * (\text{instrument resolution} + \text{operator repeatability})$
- Certified organic standards with a high temperature-compensation coefficient: $\text{uncertainty of the standard} + \text{limit of error of instrument} + \text{temperature error} + 2 * \text{operator repeatability}$

The operator repeatability is the standard deviation of ten density determinations performed by the operator.

6 Operation

6.1 Start up the density meter

- Press one of the keys (1) or the navigation joystick (2).
 - ⇒ A short beep is played.
 - ⇒ A welcome message is displayed.
- ⇒ The home screen opens. The density meter is ready to use.



6.2 Shut down the density meter

- Press and hold any of the keys or the navigation joystick for more than 3 s.
 - ⇒ A beep is played and the screen goes dark.
- ⇒ The control circuits for the keys and the navigation joystick are energized. The rest of the density meter is no longer energized.

Shut down using the Menu

- The home screen is open.
 - 1 Press the **Menu** key.
 - 2 Navigate to **Shut down**.
 - ⇒ The dialog **Do you want to shut down the instrument?** opens.
 - 3 Press the **Yes** key.
 - ⇒ A beep is played and the screen goes dark.
- ⇒ The control circuits for the keys and the navigation joystick are energized. The rest of the density meter is no longer energized.

See also

- 📖 Activate and deactivate a cleaning reminder ▶ Page 27

6.3 Navigate and enter information

6.3.1 Navigate through menus and confirm settings

- 1 Press the **Menu** key.
- 2 To navigate to a desired menu, move the navigation joystick up or down.
- 3 To navigate to a submenu of the selected menu, move the navigation joystick to the right.
- 4 To go back, move the navigation joystick to the left.



Confirm settings

In every menu, where the function of key 3 is marked with a green bar, you can press the navigation joystick to confirm the selection.


Shortcut to the method list






- The home screen is open.
- Move the navigation joystick to the left.
 - ⇒ The **Methods** window opens.

6.3.2 Enter text and numbers

- The keyboard window is open.
- 1 Navigate to the character, number, symbol or icon.
 - 2 Press the navigation joystick to confirm the selection of the character, number, symbol or icon.
 - 3 To confirm the entered information and leave the keyboard window, press the **OK** key.

Keyboard overview



No.	Icon	Function
1	  	Cycles between the screens with small, capital or special characters.
2		Deletes all entered characters.
3		Deletes the character to the left of the cursor.

No.	Icon	Function
4	◀ ▶	Moves the cursor.
5	–	Displays available characters, numbers and symbols.
6	–	Maximum number of characters you can enter.
7	–	Minimum number of characters you need to enter.
8	–	Displays the entered characters.

See also

 Change the keyboard layout ▶ Page 27

6.4 Typical phases of density determinations

Density determinations include the following two phases.

- Fill the measuring cell with sample and perform the measurement.
- Prepare the measuring cell for the next sample.

6.4.1 Fill the measuring cell and measure the density

Because the measuring cell always contains residue from the previous sample or the cleaning solution, it is important that you remove the residue using the new sample. Make sure that the residue is soluble in the new sample. To remove any residue from the previous sample, fill and empty the measuring cell three times with the new sample.

To achieve good results, make sure that the measuring cell contains only the sample you want to measure and that there are no bubbles in the measuring cell.

You can use the sample pump or syringes to fill, drain and rinse the measuring cell. The sample pump is used for samples with a viscosity up to 100 mPa*s. Syringes are typically used if the viscosity of the samples is above 100 mPa*s.

See also

 Example: Density determination using the sample pump ▶ Page 42

 Example: Density determination using a syringe ▶ Page 45

 Example: Density determination of cold or hot samples ▶ Page 48

6.4.2 Prepare the measuring cell for the next sample

At the end of this phase, the residue in the measuring cell must be soluble in the next sample. The solubility of the sample defines if you can just drain the measuring cell or if you need to rinse it as well.

- The next sample dissolves the sample you just measured: you can just drain the measuring cell.
- The next sample does not dissolve the sample you just measured: you have to drain the measuring cell and rinse it with a cleaning solution. The cleaning solution must dissolve the sample that you just measured and the next sample. When you rinse the measuring cell, the cleaning solution removes the residue from the original sample.

See also

 Clean the measuring cell ▶ Page 56

6.5 Example: Density determination using the sample pump

The following chapters show you how to configure and perform a density determination of deionized water at room temperature (23 °C or 73.4 °F).

Summary of the configuration

- Guided density determination
- Filling and draining with the sample pump

- Evaluation if result lies within specified limits

For additional information about the configuration of a method, read the chapter [Configuration of density determinations ▶ Page 28].

6.5.1 Configure the method

- The density unit is set to [g/cm³].
 - **Confirmation of results** in **Settings** is activated.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
 - ⇒ The **Methods** window opens.
 - 2 Press the **New** key.
 - ⇒ The **Method settings** window opens.
 - 3 Navigate to **Method ID** and enter a name with up to 25 characters.
 - 4 Press the **OK** key.
 - ⇒ The method is listed in the **Methods** window.
 - 5 Change the parameter settings to the values shown in the following table.

Parameter	Setting	Explanation
Workflow	Guided	The density meter prompts users to fill and drain the measuring cell.
Sample ID entry	Variable	The density meter prompts users to enter the identifier before the density is measured.
Fill cycle	3	The density meter prompts users to fill and drain the measuring cell twice and then fill the measuring cell for the measurement.
Fill speed	High	The pump fills the measuring cell with 80 % of the maximum pump speed.
Measurement reliability	Medium	The result is saved as soon as the measured temperature value stays within 0.4 °C (0.72 °F) for 10 s.
Calculation > Category	Density	A density is calculated.
Calculation > Calculation d		The density is not compensated to a certain temperature.
Calculation > Result limits	Yes	The density meter evaluates if the result is within a defined range.
Lower limit	0.9972 g/cm ³	Defines the lower limit of the range.
Upper limit	0.9977 g/cm ³	Defines the upper limit of the range.
Clean	Drain	Users are prompted to drain the measuring cell.
Printout	None	The result is not printed.

See also

- 📖 Change the density unit ▶ Page 24
- 📖 Activate and deactivate the confirmation of results ▶ Page 30
- 📖 Configuration of density determinations ▶ Page 28


6.5.2 Perform the density determination

Start the method

- The density meter is set up to work with the sample pump.
 - The piston is in the lowermost position.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
 - 2 Navigate to the method you want to use and press the **Start** key.
 - ⇒ The home screen opens.
 - 3 Press the **Start** key.
 - ⇒ The **Sample ID** window opens.
 - 4 Enter the identifier.
 - 5 Press the **OK** key.


Rinse the measuring cell

- 1 Immerse the filling tube into the sample (1).
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Hold the filling tube over the waste beaker (2).
- 4 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 5 Clean the end of the filling tube with a clean tissue.
- 6 Repeat the steps once.




Fill the measuring cell and measure the density

- 1 Immerse the filling tube into the sample (1).
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Make sure there are no bubbles in the measuring cell.
 - ⇒ The measurement starts.
 - ⇒ When the background of the measuring ribbon changes color, the measurement is finished.
 - ⇒ The measured value is saved in the results.
- 4 Press the **OK** key.



Drain the measuring cell

- 1 Hold the filling tube over the waste beaker (1).
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 3 Clean the end of the filling tube with a clean tissue.



6.6 Example: Density determination using a syringe

The following chapters show you how to configure and perform a density determination of a viscous sample such as 40 % w/w sugar solution at room temperature (23 °C or 73.4 °F).

Summary of the configuration

- Guided density determination
- Filling and draining with a syringe

For additional information about the configuration of a method, read the chapter [Configuration of density determinations ▶ Page 28].

6.6.1 Configure the method

- **Confirmation of results** in **Settings** is activated.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
 - ⇒ The **Methods** window opens.
 - 2 Press the **New** key.
 - ⇒ The **Method settings** window opens.
 - 3 Navigate to **Method ID** and enter a name with up to 25 characters.
 - 4 Press the **OK** key.
 - ⇒ The method is listed in the **Methods** window.
 - 5 Change the parameter settings to the values shown in the following table.

Parameter	Setting	Explanation
Workflow	Guided	The density meter prompts users to fill and drain the measuring cell.
Sample ID entry	Variable	The density meter prompts users to enter the identifier before the density is measured.
Fill cycle	2	The density meter prompts users to rinse the measuring cell and then fill the measuring cell for the measurement.
Fill speed	High	Has no effect because the syringe mode is activated.
Measurement reliability	Medium	The result is saved as soon as the measured temperature value stays within 0.4 °C (0.72 °F) for 10 s.

Parameter	Setting	Explanation
Calculation > Category	Sugar	A sugar concentration is calculated.
Calculation > Calculation	Brix	A Brix value is calculated.
Calculation > Result limits	None	The density meter does not evaluate if the result is within a defined range.
Clean	Drain	Users are prompted to drain the measuring cell.
Printout	None	The result is not printed.

See also

📖 Activate and deactivate the confirmation of results ▶ Page 30


6.6.2 Perform the density determination

Start the method


- The density meter is set up to work with a syringe.
 - The measuring cell is drained.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
 - 2 Navigate to the method you want to use and press the **Start** key.
 - ⇒ The home screen opens.
 - 3 Press the **Start** key.
 - ⇒ The **Sample ID** window opens.
 - 4 Enter the identifier.
 - 5 Press the **OK** key.

Rinse the measuring cell

- 1 Fill a new syringe with sample.
- 2 Insert the syringe into the connector for the syringe.




- 3 Hold the filling tube over the waste beaker (1).
- 4 Press the syringe plunger continuously into the syringe.
 - ⇒ The sample flows through the measuring cell and removes the residue from the previous sample.
- 5 Pull the plunger continuously out of the syringe.
 - ⇒ The sample is pulled back into the syringe and the measuring cell fills with air.
- 6 Remove the syringe and empty it into a suitable waste beaker.




Fill the measuring cell and measure the density

- 1 Fill a new syringe with sample.
- 2 Insert the syringe into the connector for the syringe.




- 3 Hold the filling tube over the waste beaker (1).
- 4 Press the syringe plunger slowly (5...10 cm/s) and continuously into the syringe until only a small amount of sample is left in the syringe.
 - ⇒ The sample flows with a laminar flow into the measuring cell.
- 5 Leave the syringe in the connector.
- 6 Make sure there are no bubbles in the measuring cell.
- 7 Press the **OK** key.
 - ⇒ The measurement starts.
 - ⇒ When the background of the measuring ribbon changes color, the measurement is finished.
 - ⇒ The measured value is saved in the results.
- 8 Press the **OK** key.



Drain the measuring cell

- 1 Pull the plunger continuously out of the syringe.
⇒ The sample is pulled back into the syringe and the measuring cell fills with air.
- 2 Remove the syringe and empty it into a suitable waste beaker.



6.7 Example: Density determination of cold or hot samples

To achieve good results, the measuring cell and the sample must have the same temperature. If you fill the measuring cell with a sample that is more than 5 °C (9 °F) above or below the temperature displayed on the home screen, this condition is not met. Because the density meter has no active temperature control, you must wait until the temperature of the measuring cell and the temperature of the sample are the same, before you make the measurement.

The following chapters show you how to configure and perform a density determination of deionized water that is more than 5 °C (9 °F) above or below the temperature displayed on the home screen.

Summary of the configuration

- Unguided density determination
- Filling and draining with the sample pump

For additional information about the configuration of a method, read the chapter [Configuration of density determinations ▶ Page 28].

6.7.1 Configure the method

- The density unit is set to [g/cm³].
 - **Confirmation of results** in **Settings** is activated.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
⇒ The **Methods** window opens.
 - 2 Press the **New** key.
⇒ The **Method settings** window opens.
 - 3 Navigate to **Method ID** and enter a name with up to 25 characters.
 - 4 Press the **OK** key.
⇒ The method is listed in the **Methods** window.
 - 5 Change the parameter settings to the values shown in the following table.

Parameter	Setting	Explanation
Workflow	Manual	Users define their individual measuring workflows.
Sample ID entry	Variable	Users need to enter the identifier before they start the measurement.
Fill speed	High	The pump fills the measuring cell with 80 % of the maximum pump speed.
Measurement reliability	Medium	The result is saved as soon as the measured temperature value stays within 0.4 °C (0.72 °F) for 10 s.

Parameter	Setting	Explanation
Calculation > Category	Density	A density is calculated.
Calculation > Calculation d		The density is not compensated to a certain temperature.
Calculation > Result limits	None	The density meter does not evaluate if the result is within a defined range.
Clean	Drain	Users are prompted to drain the measuring cell.
Printout	None	The result is not printed.

See also

📖 Activate and deactivate the confirmation of results ▶ Page 30


6.7.2 Perform the density determination

Start the method

- The density meter is set up to work with the sample pump.
 - The piston is in the lowermost position.
 - The home screen is open.
- 1 Move the navigation joystick to the left.
 - 2 Navigate to the method you want to use and press the **Start** key.
 - ⇒ The home screen opens.
 - 3 To enter a **Sample ID**, move the navigation joystick to the right.
 - ⇒ The **Sample ID** window opens.
 - 4 Enter the identifier.
 - 5 Press the **OK** key.


Rinse the measuring cell

- 1 Immerse the filling tube into the sample (1).
- 2 Move the navigation joystick upwards and let go.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Hold the filling tube over the waste beaker (2).
- 4 Move the navigation joystick downwards and let go.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 5 Clean the end of the filling tube with a clean tissue.
- 6 Repeat the steps once.




Fill the measuring cell and measure the density

- 1 Immerse the filling tube into the sample (1).
- 2 Move the navigation joystick upwards and let go.
⇒ The sample pump moves the piston to its uppermost position.
- 3 Make sure there are no bubbles in the measuring cell.
- 4 Press the **Start** key.
⇒ The measurement starts.
⇒ When the background of the measuring ribbon changes color, the measurement is finished.
⇒ The measured value is saved in the results.
- 5 Press the **OK** key.
- 6 Wait for about 20 seconds and press the **Start** key again.
⇒ The sample that is still in the measuring cell is measured again.
- 7 Tap **Results** and compare the last two results.
- 8 Tap **Back**.
- 9 If the results are not the same, continue to measure the same sample until the temperatures and the results of two measurements in a row are the same.



Drain the measuring cell

- 1 Hold the filling tube over the waste beaker (1).
- 2 Move the navigation joystick downwards and let go.
⇒ The sample pump moves the piston to its lowermost position.
- 3 Clean the end of the filling tube with a clean tissue.



6.8 Stop a method

Stop a method when working with the sample pump

- 1 Press the **Stop** key.
- 2 If you stopped the method before you started to fill the measuring cell, you do not need to do anything else.
- 3 If you have started to fill the measuring cell before you stopped the method, drain the measuring cell. See [Drain the measuring cell ▶ Page 66]
- 4 Clean the measuring cell. See [Clean using the sample pump ▶ Page 57]

Stop a method when working with a syringe

- 1 Press the **Stop** key.
- 2 If you stopped the method before you started to fill the measuring cell, you do not need to do anything else.
- 3 If you have started to fill the measuring cell before you stopped the method, drain the measuring cell. See [Drain the measuring cell ▶ Page 66]
- 4 Clean the measuring cell. See [Clean using a syringe ▶ Page 57]


6.9 Read and write information with the RFID reader (DensitoPro only)

With the RFID reader, you can read or write identifiers for samples or standards to Smart Tags.

The RFID reader only recognizes the following characters.

- A-Z
- a-z
- 0-9


Metal containers interfere with reading and writing to Smart Tags. For this reason you need to place the Smart Tag (1) on a SmartSample film (2).



See also


📖 Miscellaneous accessories ▶ Page 75

6.9.1 Write information to a Smart Tag

- The RFID reader is activated.
 - METTLER TOLEDO Smart Tags are available.
- 1 Press the **Menu** key.
 - 2 Navigate to **Settings > RFID / Barcode > Write ID to RFID tag**.
 - ⇒ The icon  is displayed in the status ribbon.
 - 3 Enter the identifier.
 - 4 Hold the RFID reader over a Smart Tag and press the **Write** key.
 - ⇒ An audio signal is played and the notification **RFID tag is written** is displayed.
- ⇒ The identifier is saved on the Smart Tag.

See also


📖 Activate and deactivate the RFID reader (DensitoPro only) ▶ Page 16



6.9.2 Read information from a Smart Tag

Enter the identifier of a sample or a standard

- The RFID reader is activated.
 - The parameter **Sample ID entry** or **Standard ID entry** is set to **Variable**.
- When you are prompted to enter an identifier, press the **RFID** key and hold the RFID reader over the Smart Tag.
- ⇒ An audio signal is played and the identifier is saved.




Read information from a Smart Tag

- The RFID reader is activated.
- 1 Press the **Menu** key.

- 2 Navigate to **Settings > RFID / Barcode > Scan RFID tag**.
 - ⇒ The notification **Scan RFID tag** is displayed.
- 3 Hold the RFID reader over a Smart Tag and wait until an audio signal is played.
 - ⇒ The read information is displayed.

See also

- 📖 Activate and deactivate the RFID reader (DensitoPro only) ▶ Page 16
- 📖 Configuration of density determinations ▶ Page 28



6.10 Read information with the barcode reader (DensitoPro only)

With the barcode reader, you can read identifiers for samples and standards from linear barcodes. The barcode reader only recognizes the following characters.

- A-Z
- a-z
- 0-9

Procedure

- The barcode reader is activated.
- The parameter **Sample ID entry** or **Standard ID entry** is set to **Variable**.
- When you are prompted to enter an identifier, press the **Barcode** key and hold the barcode reader over the barcode.
 - ⇒ An audio signal is played and the identifier is saved.

See also

- 📖 Activate and deactivate the barcode reader (DensitoPro only) ▶ Page 16
- 📖 Configuration of density determinations ▶ Page 28

6.11 View, export, print and delete results

6.11.1 View results

- The home screen is open.
 - 1 Press the **Results** key.
 - ⇒ A list with all results is displayed.
 - 2 To view an individual result, navigate to the result and move the navigation joystick to the right.

6.11.2 Delete all results

You can only delete all results but not individual results or a subgroup of results.

- The home screen is open.
 - 1 Press the **Results** key.
 - 2 Press the **Delete** key.
 - ⇒ The message **Do you want to delete all results?** is displayed.
 - 3 Press the **Yes** key.
 - ⇒ An empty list is displayed.

6.11.3 Export and print results



NOTICE

Loss of data

A maximum of 1 100 results is saved in the density meter. When the result list contains 1 100 results and a new results is saved, the oldest result is overwritten.

- Print or export the results before the maximum number of results is reached.

You can export or print results if you have the required accessories.

- To export results you need the EasyDirect Density & Refractometry software or a USB flash drive.
- To print results, you need a compatible USB printer or bluetooth printer.

See also

- 📖 Accessories ▶ Page 74
- 📖 Connect to EasyDirect Density & Refractometry ▶ Page 23
- 📖 Connect a printer ▶ Page 22

6.11.3.1 Export results to EasyDirect Density & Refractometry

Export an individual result

- EasyDirect Density & Refractometry is connected.
 - An individual result is displayed.
- 1 Press the **Report** key.
 - 2 Navigate to **Option** and select **EasyDirect**.
 - 3 Press the **OK** key.
 - 4 Press the **Start** key.
- ⇒ The export is completed and the result is displayed.

Export all results or a subset of results

- EasyDirect Density & Refractometry is connected.
 - The result list is displayed.
- 1 Press the **Report** key.
 - 2 Navigate to **Option** and select **EasyDirect**.
 - 3 Press the **OK** key.
 - 4 To export all results, deactivate the filter.
 - 5 To export a subset of results, configure the filter.
 - 6 Press the **Start** key.
- ⇒ The export is completed and the result list is displayed.

See also

- 📖 Configure the filter ▶ Page 55
- 📖 Connect to EasyDirect Density & Refractometry ▶ Page 23

6.11.3.2 Export results to a USB flash drive

Individual results, subsets of results or all results can be exported as CSV files to a USB flash drive. Each time an export is performed, the data is stored in a new CSV file.

Characteristics of the CSV files

- A CSV file only contains data of one method type. Example: If results from test methods and measurement methods are exported, they are saved to two separate CSV files. One file contains the results from the test methods. The other file contains the results from the measurement methods.
- File names are constructed according to the following pattern: Method type_product type_serial number_date_time.csv. Example: MeasureExport_DensitoPro_A001711238_20200519_073724.csv
- Data delimiter: tab

Export an individual result

- A USB flash drive is connected.
 - An individual result is displayed.
- 1 Press the **Report** key.
 - 2 Navigate to **Option** and select **USB stick**.
 - 3 Press the **OK** key.
 - 4 Press the **Start** key.
- ⇒ The result is exported.

Export all results or a subset of results

- A USB flash drive is connected.
 - The result list is displayed.
- 1 Press the **Report** key.
 - 2 Navigate to **Option** and select **USB stick**.
 - 3 Press the **OK** key.
 - 4 To export all results, deactivate the filter.
 - 5 To export a subset of results, configure the filter.
 - 6 Press the **Start** key.
- ⇒ The results are exported.

See also

📖 Configure the filter ▶ Page 55

6.11.3.3 Print results

The **Printout** parameter defines whether a result is automatically printed when it is saved and which information is printed.

Information	Methods	Test	Adjustment
Date / Time	Short/Long	Short/Long	Short/Long
Method ID	Short/Long	–	–
Sample ID	Short/Long	–	–
Application name with standard	–	Short/Long	Short/Long
Standard ID	–	Short/Long	Short/Long
Status	Short/Long	Short/Long	–
Calculation	Short/Long	–	–
Density	Short/Long	Short/Long	Short/Long
Dev. d	–	Short/Long	Short/Long
Temp.	Short/Long	Short/Long	Short/Long
d nominal	–	Short/Long	Short/Long
SNR	Long	Long	Long

Information	Methods	Test	Adjustment
Additional settings like the compensation temperature and the temperature-compensation coefficient α	Long	Long	Long

Print an individual result

- A printer is connected.
 - An individual result is displayed.
- 1 Press the **Report** key.
 - 2 Navigate to **Option** and select **Print**.
 - 3 Press the **OK** key.
 - 4 Navigate to **Printout** and set it to **Short** or **Long** as needed.
 - 5 Press the **OK** key.
- ⇒ The result is printed.

Print all results or a subset of results

- A printer is connected.
 - The result list is displayed.
- 1 Press the **Report** key.
 - 2 Navigate to **Option** and select **Print**.
 - 3 Press the **OK** key.
 - 4 Navigate to **Printout** and set it to **Short** or **Long** as needed.
 - 5 Press the **OK** key.
 - 6 To print all results, deactivate the filter.
 - 7 To print a subset of results, configure the filter.
 - 8 Press the **Start** key.
- ⇒ The results are printed.

See also

- 📖 Configure the filter ▶ Page 55
- 📖 Connect a printer ▶ Page 22

6.11.3.4 Configure the filter

If you activate the filter, you can print or export a subset of the results. The following filters are available.

- **Off**: Filter is deactivated and all results are printed or exported.
- **Today**: Only the results from current day are printed or exported.
- **Range**: Only the results within the set date range are printed or exported.

Procedure

- 1 Press the **Menu** key.
- 2 Navigate to **Results > Report > Filter**.
- 3 Select the desired **Filter** type.
- 4 Press the **OK** key.
- 5 If you selected **Range**, define the values for **Start date** and **End date**.

See also

- 📖 Export results to a USB flash drive ▶ Page 53
- 📖 Export results to EasyDirect Density & Refractometry ▶ Page 53
- 📖 Print results ▶ Page 54

7 Maintenance

In this chapter you find descriptions of the maintenance tasks you should perform on your density meter. Any other maintenance tasks need to be performed by a service technician that has been qualified by METTLER TOLEDO.

If you experience problems with your density meter, contact your authorized METTLER TOLEDO service representative or dealer.

METTLER TOLEDO recommends that a preventive maintenance and calibration certification is done at least once a year through your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact

7.1 Maintenance schedule

If the standard operating procedures of your company require other maintenance intervals, use the intervals listed in the standard operating procedures.

Frequency	Task	See
Daily	Clean the measuring cell at the end of the work day.	[Clean the measuring cell ► Page 56]
	Perform a test with deionized water.	[Check the measurement accuracy ► Page 60]

See also

📖 Activate and deactivate a cleaning reminder ► Page 27

7.2 Clean the density meter



NOTICE

Damage to the density meter due to inappropriate cleaning methods

Inappropriate cleaning agents can damage the housing or other parts of the density meter. If liquids enter the housing, they can damage the density meter.

- 1 Make sure the cleaning agent is compatible with the material of the part you want to clean.
- 2 Make sure that no liquid enters the interior of the density meter.

If you have questions about the compatibility of cleaning agents, contact your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact

7.2.1 Clean the housing

METTLER TOLEDO recommends the following cleaning agents:

- Water
- Water with a mild detergent

Procedure

- The density meter is shut down.
- Wipe the housing with a cloth moistened with the cleaning agent.

7.2.2 Clean the measuring cell

Because the measuring cell always contains residue from the last sample, there is a risk that the measuring cell can be damaged by the sample. To prevent such damages, METTLER TOLEDO recommends that you clean the measuring cell at the end of each sample series.

To clean the measuring cell, you need to rinse the measuring cell with a cleaning solution. The cleaning solution must have the following properties.

- Does not chemically interact with the material of the measuring cell.
- Dissolves the sample that was just measured.
- Evaporates without leaving incrustations.

METTLER TOLEDO recommends the following cleaning solutions.

Sample	Cleaning solution
Water Water based	Deionized water
Acids, concentrated	Water (flush the measuring cell with plenty of water to remove the heat from the reaction of water and acid)
Alkaline, concentrated	0.3...0.5 % deconex solution
Samples with fats or oily components	0.3...0.5 % deconex solution
Petrochemical samples Edible oils and fats	Toluene, xylene or petrol ether mixtures


See also

📖 Activate and deactivate a cleaning reminder ▶ Page 27

7.2.2.1 Clean using the sample pump

- The density meter is set up to work with the sample pump.
- The cleaning solution is defined and available.
- The piston is in the lowermost position.


- 1 Immerse the filling tube in the cleaning solution (1).
- 2 Move the navigation joystick upwards and let go.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Lift the filling tube out of the cleaning solution and hold it over the waste beaker (2).
- 4 Move the navigation joystick downwards and let go.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 5 Clean the end of the filling tube with a clean tissue.
- 6 Repeat the steps until all residue from the sample is removed.
- 7 Move the navigation joystick upwards and let go.
 - ⇒ The sample pump fills the measuring cell with air.
- 8 Move the navigation joystick downwards and let go.
 - ⇒ Most of the cleaning solution is removed from the measuring cell.
- 9 Clean the end of the filling tube with a clean tissue.




7.2.2.2 Clean using a syringe

- The density meter is set up to work with a syringe.
 - The cleaning solution is defined and available.
 - The measuring cell is drained.
- 1 Fill a syringe with the cleaning solution.


- 2 Insert the syringe into the connector for the syringe.



- 3 Hold the filling tube over the waste beaker (1).
- 4 Press the plunger continuously into the syringe until the syringe is empty.



- 5 Pull the plunger continuously out of the syringe.
⇒ The cleaning solution is pulled back into the syringe and the measuring cell fills with air.
- 6 Remove the syringe and empty it into a suitable waste beaker.
- 7 Repeat the steps until all residue from the sample is removed.




7.2.3 Clean the cylinder and the piston of the sample pump

If rinsing is not enough to clean the cylinder and the piston, you can remove both for more thorough cleaning.


7.2.3.1 Disassemble the cylinder and remove the piston

- 1 Remove the cylinder and the piston from the density meter. See [Remove the cylinder and the piston ► Page 63].

- 2 Push the piston (1) into the cylinder until it touches the bottom section (2) of the cylinder.
- 3 Carefully apply pressure to the piston (1) and push the bottom section (2) partway out of the glass section of the cylinder (3).
- 4 Pull the bottom section (2) out of the glass section (3).



- 5 Pull the piston (1) out of the glass section (2).



7.2.3.2 Clean the piston and the cylinder sections

METTLER TOLEDO recommends the following cleaning agents:

- Water
- Water with a mild detergent

Procedure


- 1 Wipe the piston and the cylinder sections with a cloth moistened with the cleaning agent.
- 2 If needed, soak the piston and the cylinder sections in the cleaning agent.
- 3 If needed, scrub the piston and the cylinder sections with a non-abrasive brush.
- 4 Rinse the piston and the cylinder sections with deionized water.
- 5 Dry the piston and the cylinder sections.

See also

📖 Technical data ▶ Page 70

7.2.3.3 Assemble the cylinder and install the piston

- 1 Insert the bottom section (3) into the glass section (1).
- 2 Push the bottom section (3) into the glass section (1) until the ledge (2) touches the glass.
- 3 Insert the piston (4) into the glass section (1).
- 4 Install the cylinder and the piston into the density meter. See [Install the cylinder and the piston ▶ Page 64].



7.2.4 Clean the protective cover

The protective cover is made of EPDM (ethylene propylene diene monomer (M-class) rubber).

METTLER TOLEDO recommends the following cleaning agents:

- Water
- Water with a mild detergent

Procedure

- The density meter is shut down.
- 1 Remove the protective cover.
 - 2 Wipe the outside and the inside of the protective cover with a cloth moistened with the cleaning agent.

See also

- 📖 Install and remove the protective cover ▶ Page 21
- 📖 Clean the housing ▶ Page 56

7.3 Check the measurement accuracy

To check the measurement accuracy of the density meter, you need to configure and perform a test.

For additional information about the configuration of a test, read the chapter [Configuration of density determinations ▶ Page 28].

7.3.1 Example: Test with water

The following chapters show you how to configure and perform a guided test with deionized water at room temperature (23 °C or 73.4 °F).

7.3.1.1 Configure the test

- The density unit is set to [g/cm³].
 - **Confirmation of results** in **Settings** is activated.
- 1 Press the **Menu** key.
 - 2 Navigate to **Test**.
 - ⇒ The **Test** window opens.
 - 3 Change the parameter settings to the values shown in the following table.

Parameter	Setting	Explanation
Workflow	Guided	The density meter prompts users to fill and drain the measuring cell.
Standard type	Water	The density meter uses predefined values for the water measurement.
Standard ID entry	Variable	The density meter prompts users to enter the identifier before the density is measured.
Fill cycle	3	The density meter prompts users to fill and drain the measuring cell twice and then fill the measuring cell for the measurement.
Fill speed	High	The pump fills the measuring cell with 80 % of the maximum pump speed.
Tolerance d	0.0002 g/cm ³	The result is marked as passed if it falls within the following range. d nominal ± Tolerance d
Clean	Drain	Users are prompted to drain the measuring cell.
Printout	None	The result is not printed.

See also

- 📖 Activate and deactivate the confirmation of results ▶ Page 30
- 📖 Change the density unit ▶ Page 24

7.3.1.2 Perform the test

Start the test

- The density meter is set up to work with the sample pump.
 - The measuring cell has been cleaned.
 - The piston is in the lowermost position.
- 1 Press the **Menu** key.
 - 2 Navigate to **Test**.
 - ⇒ The **Test** window opens.
 - 3 Press the **Start** key.
 - ⇒ The home screen opens.
 - 4 Press the **Start** key.
 - ⇒ The **Standard ID** window opens.
 - 5 Enter the identifier.
 - 6 Press the **OK** key.

Rinse the measuring cell

- 1 Immerse the filling tube in the water.
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Hold the filling tube over the waste beaker.
- 4 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 5 Clean the end of the filling tube with a clean tissue.
- 6 Repeat the steps once.

Fill the measuring cell and measure the density

- 1 Immerse the filling tube in the water.
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Make sure there are no bubbles in the measuring cell.
 - ⇒ The measurement starts.
 - ⇒ When the background of the measuring ribbon changes color, the measurement is finished.
 - ⇒ The measured value is saved in the results.
- 4 Press the **OK** key.

Drain the measuring cell

- 1 Hold the filling tube over the waste beaker.
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 3 Clean the end of the filling tube with a clean tissue.

7.3.2 Measures if a test fails

- 1 Check if the correct standard has been used, and if needed repeat the test with the correct standard.

- 2 If the test continues to fail, clean the measuring cell with a cleaning solution that dissolves your samples and dissolves in the standard.
- 3 Repeat the test.
- 4 If the test continues to fail, repeat the test twice and compare all three results.
- 5 If the results are not the same, clean the measuring cell and perform tests until the test passes or the results of three consecutive tests are the same.
- 6 If the test continues to fail and the results of three consecutive tests are the same, restore the factory adjustment data and repeat the test.
- 7 If the test continues to fail, perform an adjustment.

See also

- 📖 Clean the measuring cell ▶ Page 56
- 📖 Restore the factory-adjustment data ▶ Page 66
- 📖 Adjust the measurement accuracy ▶ Page 66

7.4 Replace the battery



⚠ CAUTION



Injury or damage due to incorrect handling of the battery

The battery can explode or ignite if it is not handled correctly.

- 1 Do not store or use the battery at temperatures below $-20\text{ }^{\circ}\text{C}$ or above $60\text{ }^{\circ}\text{C}$.
- 2 Do not subject the battery to mechanical stresses like pressure, bending or impacts.

The density meter works with a rechargeable 2400 mAh Li-Ion battery that you can order from METTLER TOLEDO.


- 1 Press on the lower part (1) of the battery cover and push the battery cover up.
- 2 Pull the upper part of the battery out of the battery compartment and lift the battery out.



- 3 Insert the lower part of the battery at an angle into the battery compartment.
- 4 Push the upper part of the battery into the battery compartment.
- 5 Slide the battery cover from above onto the density meter.

See also

- 📖 Accessories ▶ Page 74
- 📖 Technical data ▶ Page 70




7.5 Replace the filling tube

- 1 Unscrew the filling tube by hand and remove it.
- 2 Install the new filling tube.

See also


- 📖 Accessories ▶ Page 74
- 📖 Install the filling tube ▶ Page 16



7.6 Replace the cylinder and the piston of the sample pump

7.6.1 Remove the cylinder and the piston

- The measuring cell is drained.
- 1 Move the navigation joystick upwards and let go.
⇒ The sample pump moves the piston to its uppermost position.
 - 2 Shut down the density meter.
 - 3 Press on the grip recesses (4) of the cover and lift the lower part of the cover out of the housing.
 - 4 Pull the top of the cover out of the housing.
 - 5 Lift the piston (1) out of the guide system.
 - 6 Pull the cylinder (2) up and out of its guide system (3).




7.6.2 Install the cylinder and the piston

- 1 Insert the cylinder (3) of the sample pump into its guide system (4).
- 2 Insert the piston (1) into its guide system (2).
- 3 Reinstall the cover of the sample pump.
- 4 Start up the density meter.
- 5 Move the navigation joystick downwards and let go.
⇒ The sample pump moves the piston to its lowermost position.

See also

📖 Accessories ▶ Page 74



7.7 Update the firmware

The firmware of your density meter can be updated. If you want to update the firmware of your density meter, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

7.8 Prepare the density meter for storage

The density meter can be stored at temperatures of -20...+50 °C (-4...+122 °F). But the lifespan of the battery is reduced when it is stored at high temperatures.

- 1 Clean the measuring cell.
- 2 Shut down the density meter.
- 3 Clean the housing.
- 4 Store the density meter in a dry and clean place.

See also

📖 Clean the density meter ▶ Page 56

📖 Technical data ▶ Page 70

7.9 Ship the density meter

If you have questions about transporting your density meter, contact your authorized METTLER TOLEDO dealer or service representative.

▶ www.mt.com/contact

- 1 Shut down the density meter.
- 2 Clean the density meter.
- 3 Ship the density meter using the original packaging.

7.10 Dispose of the density meter

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.



8 Troubleshooting

8.1 List of errors and problems

Problem	Possible cause	Measure
The sample pump does not aspirate any liquid.	The locking screw is not installed.	– Install the locking screw.
	The locking screw is not tight.	– Tighten the locking screw.
	The O-rings of the cylinder are damaged.	– Replace the cylinder and the piston. See [Replace the cylinder and the piston of the sample pump ▶ Page 63].
Test of measurement accuracy fails.	The measuring cell is not clean.	See [Measures if a test fails ▶ Page 61].
	Mismatch between measured and configured standard	See [Measures if a test fails ▶ Page 61].
	The tolerance range is too small.	See [Configure the tolerance for a test ▶ Page 39].
	The density meter is not adjusted correctly.	See [Measures if a test fails ▶ Page 61].
The message Drain cell over waste beaker. is displayed when you try to start a method.	The piston is not in the lowermost position.	See [Drain the measuring cell ▶ Page 66]
There are bubbles in the measuring cell.	Fill speed is too high.	1 If you use the sample pump, lower the setting for fill speed. See [Define the pump speed ▶ Page 30].
		2 If you use a syringe, press the syringe plunger more slowly into the syringe.
	The liquid contains dissolved gases.	– Degas the sample before you fill it into the measuring cell.
	Connector for syringe or the filling tube is not tight.	1 Tighten the connector of the filling tube.
		2 Tighten the locking screw or the connector for the syringe.
The message Pump position unknown. is displayed.	The saved pump position was lost because the battery was removed too long.	The density meter is set up to work with a syringe. – You can ignore the message and continue to work with the density meter. The density meter is set up to work with the sample pump. See [Drain the measuring cell ▶ Page 66]
"--" is displayed, instead of a value for Alpha*1000.	The calculated value is out of range.	– Check and correct the density values and the values for the compensation temperatures you entered.

Problem	Possible cause	Measure
"--" is displayed, instead of a measured value or a result.	The measured value or the result is out of range.	<ol style="list-style-type: none"> 1 Check that the unit defined in the method matches the sample you want to measure. 2 If you use a user-defined conversion, check that all the settings are correct.


8.2 Drain the measuring cell

Drain using the sample pump

- 1 Hold the filling tube over the waste beaker.
- 2 Move the navigation joystick downwards and let go.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 3 Clean the end of the filling tube with a clean tissue.

Drain using a syringe

- 1 Pull the plunger continuously out of the syringe.
 - ⇒ The sample is pulled back into the syringe and the measuring cell fills with air.
- 2 Remove the syringe and empty it into a suitable waste beaker.



8.3 Restore the factory-adjustment data

You can delete the current adjustment data and restore the factory-adjustment data.

- 1 Press the **Menu** key.
- 2 Navigate to **Maintenance & Service > Reset to factory adjust..**
 - ⇒ The message **Do you want to reset to factory adjustment?** opens.
- 3 To overwrite the current adjustment data with the factory adjustment data, press **Yes**.
 - ⇒ The factory adjustment data are restored.

8.4 Adjust the measurement accuracy



NOTICE

Inaccurate results due to incorrect adjustment

Adjustments that are not performed properly can lead to inaccurate results.

- Only perform an adjustment if you are qualified to do it.

The following chapters show you, how to configure and perform an exemplary adjustment with water.

For additional information about the configuration of an adjustment, read the chapter [Configuration of density determinations ▶ Page 28].

See also

- 📖 Measures if a test fails ▶ Page 61

8.4.1 Example: Adjustment with water

The following chapters show you how to configure and perform a guided adjustment with deionized water at room temperature (23 °C or 73.4 °F).

8.4.1.1 Configure the adjustment

- **Confirmation of results** in **Settings** is activated.

- 1 Press the **Menu** key.
- 2 Navigate to **Adjustment**.
 - ⇒ The **Adjustment** window opens.
- 3 Change the parameter settings to the values shown in the following table.

Parameter	Setting	Explanation
Workflow	Guided	The density meter prompts users to fill and drain the measuring cell.
Standard type	Water	The density meter uses predefined values for the water measurement.
Standard ID entry	Variable	The density meter prompts users to enter the identifier before the density is measured.
Fill cycle	3	The density meter prompts users to fill and drain the measuring cell twice and then fill the measuring cell for the measurement.
Fill speed	High	The pump fills the measuring cell with 80 % of the maximum pump speed.
Clean	Drain	Users are prompted to drain the measuring cell.
Printout	None	The result is not printed.


8.4.1.2 Perform the adjustment

Start the adjustment

- The density meter is set up to work with the sample pump.
 - The measuring cell has been cleaned.
 - The piston is in the lowermost position.
- 1 Press the **Menu** key.
 - 2 Navigate to **Adjustment**.
 - ⇒ The **Adjustment** window opens.
 - 3 Press the **Start** key.
 - ⇒ The home screen opens.
 - 4 Press the **Start** key.
 - ⇒ The **Standard ID** window opens.
 - 5 Enter the identifier.
 - 6 Press the **OK** key.


Rinse the measuring cell

- 1 Immerse the filling tube in the water (1).
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Hold the filling tube over the waste beaker (2).
- 4 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 5 Clean the end of the filling tube with a clean tissue.
- 6 Repeat the steps once.



Fill the measuring cell and measure the density

- 1 Immerse the filling tube in the water (1).
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its uppermost position.
- 3 Make sure there are no bubbles in the measuring cell.
 - ⇒ The measurement starts.
 - ⇒ When the background of the measuring ribbon changes color, the measurement is finished.
 - ⇒ The measured value is saved in the results.
- 4 Accept or reject the adjustment.




Drain the measuring cell

- 1 Hold the filling tube over the waste beaker (1).
- 2 Press the **OK** key.
 - ⇒ The sample pump moves the piston to its lowermost position.
- 3 Clean the end of the filling tube with a clean tissue.

See also

- 📄 Activate and deactivate the confirmation of results ▶ Page 30



8.5 Restore the factory settings



NOTICE

Loss of data

With a factory reset all user-specific settings are set to the default values and all saved data (e.g. sample IDs, results, Method IDs and adjustment data) are deleted.

- Document the data you need before you restore the factory settings.

- 1 Press the **Menu** key.
- 2 Navigate to **Maintenance & Service > Reset to factory settings**.
 - ⇒ The dialog **Do you want to reset all values to factory settings?** opens.
- 3 To reset the density meter to factory settings, press the **Yes** key.
 - ⇒ The density meter restarts.

8.6 View the firmware version and other system information

- The home screen is open.
- 1 Press the **Menu** key.
 - 2 Navigate to **Information**.
 - ⇒ The firmware version and other system information are displayed.

Parameter	Description
Instrument type	Shows if the density meter is a Densito or DensitoPro.
SNR instrument	Shows the serial number of the instrument.
Firmware version	Shows the firmware version of the density meter.
Licenses	Opens a window with information on the End User License Agreement, third party licenses and open source attribution files.

9 Technical data

9.1 Density meter

Characteristic		Value
Power rating instrument	Input values	5 V DC, 1.0 A
	Connector type	USB C
Power rating AC/DC adapter	Input values	100...240 V AC, 0.3 A
	Input frequency	50 - 60 Hz
	Output values	5 V DC, 1.0 A
Battery	Type	Rechargeable lithium-ion battery
	Capacity	2400 mAh
	Temperature range for charging	0...40 °C
Dimensions	Width	56 mm
	Length	288 mm
	Depth	80 mm
	Weight	355 g
Screen	Technology	Color screen
	Size	2.4"
	Resolution	240 x 320 pixels
USB type A	Data rate	USB 2.0, host fullspeed
	Output voltage	5 V DC \pm 5 %
	Max. output current	250 mA
USB type C	Data rate	USB 2.0, device fullspeed
	Input voltage	5 V DC \pm 10 %
	Max. input current	400 mA, max. 1 A with AC adapter
RFID reader (DensitoPro only)	Compatibility	METTLER TOLEDO Smart Tags
Barcode reader (DensitoPro only)	Format	Linear barcodes
Materials	Housing	PBT (polybutylene terephthalate)
	Measuring cell	Borosilicate glass
	Piston	PTFE (polytetrafluorethylene) Stainless steel (X8CrNiS18-9 (1.4305), X10CrNi18-8 (1.4310) and X5CrNi18-10 (1.4301))
	Cylinder	PVDF (polyvinylidene difluoride) FKM (fluoroelastomer) FFKM (perfluoro-elastomers) Borosilicate glass
	Filling tube	PTFE (polytetrafluorethylene)
	Connector for filling tube or syringe	PP (polypropylene)
	Locking screw	PVDF (polyvinylidene difluoride)
Ambient conditions		
Characteristic		Value
Ambient temperature (Densito)		-10...+50 °C

Characteristic	Value
Ambient temperature (DensitoPro)	0...+50 °C
Relative humidity	Max. 80 % (non-condensing) at 31 °C, linear fall to 50 % at 40 °C
Altitude	Up to 2000 m above sea level
Pressure	Atmospheric pressure
Use	In interior spaces
Overvoltage category	II
Pollution degree	2
IP code	IP5X

Storage temperatures

Characteristic	Value	
Instrument	-20 °C...+50 °C	
Battery	Storage for up to 1 month	-20 °C...+60 °C
	Storage for up to 3 months	-20 °C...+45 °C
	Storage for up to 1 year	-20 °C...+30 °C

9.2 Measurement

Characteristic	Value	
Density	Measuring range	0.000...3.000 g/cm ³
	Accuracy ¹⁾	± 0.001 g/cm ³
	Repeatability	± 0.0005 g/cm ³
	Resolution	0.0001 g/cm ³
Temperature	Temperature range for sample	0...50 °C
	Resolution	0.1 °C
	Accuracy	±0.2 °C
Viscosity	Using sample pump	0...100 mPa*s
	Using syringe	0...1000 mPa*s
Volumes	Sample pump	5 mL
	Min. sample volume	2 mL
Temperature compensation	Automatic or by using a user-defined temperature correction coefficient	
Application	Max. number of methods	30
	Stored results	1100

¹⁾ For viscosity range of 0...200 mPa*s

Measurement scales

Characteristic	Value	
Ethanol (d) % w/w	Range	0.0...100.0 % w/w
	Accuracy	±0.5 % w/w
Ethanol (d) % v/v	Range	0.0...100.0 % v/v
	Accuracy	±0.5 % v/v
Proof (IP)	Range	0.0...175.0 Proof
	Accuracy	±1.75 Proof

Characteristic		Value
Proof (US)	Range	0.0...200.0 Proof
	Accuracy	±2.0 Proof
Light Baumé	Range	10...100 °Bé
	Accuracy	±0.4 °Bé
Heavy Baumé	Range	0...72 °Bé
	Accuracy	±0.1 °Bé
H₂SO₄ % w/w	Temperature range	10...50 °C
	Concentration range	0.0...90.0 % w/w
	Accuracy	±1.0 % w/w
H₂SO₄ % v/v	Temperature range	10...50 °C
	Concentration range	0.0...90.0 % v/v
	Accuracy	±1.0 % v/v
Brix (d) % w/w	Range	0.0...85.0 % w/w
	Accuracy	±0.3 % w/w
Plato	Range	0.0...85.0 °P
	Accuracy	±0.3 °P

9.3 Charging stand

Characteristic		Value
Power rating charging stand	Input values	5 V DC, 1.0 A
	Connector type	USB C
Power rating AC/DC adapter	Input values	100...240 V AC, 0.3 A
	Input frequency	50 - 60 Hz
	Output values	5 V DC, 1.0 A
Dimensions	Width	71 mm
	Length	250 mm
	Depth	117 mm
	Weight	390 g
Materials	Housing	ABS (acrylonitrile butadiene styrene) POM-C (polyacetal copolymer) Stainless steel (X5CrNi18-10 (1.4301))

Ambient and storage conditions

Ambient temperature	-10...+50 °C
Storage temperature	-20...+50 °C
Relative humidity	Max. 80 % (non-condensing) at 31 °C, linear fall to 50 % at 40 °C
Altitude	Up to 2000 m above sea level
Pressure	Atmospheric pressure
Use	In interior spaces
Overvoltage category	II
Pollution degree	2

9.4 Protective cover

Characteristic		Value
Dimensions	Width	63 mm
	Length	277 mm
	Depth	122 mm
	Weight	75 g
Materials		EPDM (ethylene propylene diene monomer (M-class) rubber)

10 Accessories


The following accessories and parts are available for the density meters Densito and DensitoPro.

All parts and accessories are specified with their order number, and quantity where more than one part is included.

If you have any questions, contact your authorized METTLER TOLEDO service representative or dealer.

► www.mt.com/contact


10.1 Sampling



Filling Tube
190 mm

- Tube
- Connector syringe/filling tube
- Washer


30330847



Filling tube
600 mm


- Tube
- Connector syringe/filling tube
- Washer

30330848




Locking Screw, Densito

30330852




Connector syringe/filling tube

30330853



Disposable syringe (100 pcs)
10 mL


51338100



Sample pump

- Piston
- Glass cylinder


30330849



Sample pump cover


30330851

10.2 Printers




Bluetooth Printer Godex MX20 with Dongle

30330864




Bluetooth Printer Paper (5 rolls)

30330865



Bluetooth Dongle Godex MX20

30330863



USB-P25 printer


11124301
US/CA 11124321

10.3 Miscellaneous accessories




EasyDirect Density & Refractometry computer software

30451628




Density standards (10 pcs)
6 mL

51325005




Protective Cover

30330860




Charging stand

30330842




Carrying Case

30330861




Battery Li-Ion 2400mAh

30330855




Battery cover

30330854



USB-C Cable

30449253



Power Supply and Worldwide Adapters


30449255

Accessories for DensitoPro



Tag SmartSample (10 pcs)

30449268



Tag SmartSample Film for Tanks (10 pcs)

30449269

11 Appendix

11.1 Density of pure water (0...40 °C)

Temperature [°C]	Density [g/cm ³]	Temperature [°C]	Density [g/cm ³]
0	0.99984		
1	0.99990	21	0.99799
2	0.99994	22	0.99777
3	0.99996	23	0.99754
4	0.99997	24	0.99730
5	0.99996	25	0.99705
6	0.99994	26	0.99679
7	0.99990	27	0.99652
8	0.99985	28	0.99624
9	0.99978	29	0.99595
10	0.99970	30	0.99565
11	0.99961	31	0.99534
12	0.99950	32	0.99503
13	0.99938	33	0.99471
14	0.99925	34	0.99438
15	0.99910	35	0.99404
16	0.99894	36	0.99369
17	0.99878	37	0.99333
18	0.99860	38	0.99297
19	0.99841	39	0.99260
20	0.99821	40	0.99222

[Chemical Handbook Fundamental Version, Rev. 3, Table 5.2 (1984)]

11.2 Temperature-compensation coefficients α

Substance	Temperature range [°C]	$\alpha \cdot 10^3 / ^\circ\text{C}$
Water	15...30	0.23
Ethanol	0...30	1.09
m-xylene	0...30	0.99
p-xylene	15...30	1.02
Glycerine	15...30	0.49
Chloroform	0...30	1.26
Carbon tetrachloride	0...30	1.22
Toluene	0...30	1.07
Benzene	6...30	1.21
Methanol	6...30	1.18
Acetone	0...30	1.42
Bromobenzene	0...30	0.91
Cyclohexane	0...30	1.20
Isopropanol	0...30	1.06
n-nonane	0...30	1.08

Index

A

AC/DC adapter	
assemble	14, 19
change prongs	14, 19
connect charging stand	20
connect density meter	14
install prongs	14, 19
power rating	70, 72
accessories	74
activate	
audio signal	26
auto rotation	26
barcode reader	16
cleaning reminder	27
confirmation of results	30
password protection	26
RFID reader	16
sample pump mode	15
syringe mode	17
adjustment	
change	29
configure	29, 35, 37, 38, 39
configure example	67
density unit	24
perform	66
perform example	67
alcohol concentration	31
mass percent	32
proof (IP)	32
proof (US)	32
volume percent	32
alpha	
configure	35
table	77
ambient conditions	
charging stand	72
density meter	70
American Petroleum Institute	31
API	31
API degree	33
density	33
specific gravity	33
assemble	
AC/DC adapter	14, 19

filling tube	16
audio signal	
activate	26
deactivate	26
auto rotation	
activate	26
deactivate	26

B

barcode reader	
activate	16
deactivate	16
read	52
battery	
install	62
remove	62
replace	62
storage temperature	71
technical data	70
Baumé	31, 33
beep	
activate	26
deactivate	26
bluetooth printer	
connect	22, 23
brightness	
change	25
dimming	25
Brix	33

C

celsius	24
change	
adjustment	29
date	24
date format	24
density unit	24
keyboard layout	27
language	27
method	28
password	26
prongs AC/DC adapter	14, 19
screen brightness	25
temperature unit	24
test	29
time	24

time format	24	confirmation of results	
charge		activate	30
density meter	13	deactivate	30
charging stand		connect	
connect AC/DC adapter	20	AC/DC adapter	14, 20
dimensions	72	bluetooth printer	22, 23
function	13	EasyDirect Density & Refractometry	23
material	72	Godex MX20	22, 23
mount on wall	20	power supply	14, 20
power rating	72	USB printer	22
site requirements	19, 72	connector for syringe	
storage temperature	72	install	18
weight	72	material	70
clean		consumables	74
housing	56	contact	5
measuring cell	56	create	
protective cover	59	method	28
cleaning reminder		D	
activate	27	date	
deactivate	27	change	24
concentration		date format	
alcohol	31	change	24
sugar	31	deactivate	
sulfuric acid	31	audio signal	26
user-defined	31, 34	auto rotation	26
configure		barcode reader	16
adjustment	29, 35, 37, 38, 39	cleaning reminder	27
alpha	35	confirmation of results	30
drain	30	password protection	26
measurement reliability	31	RFID reader	16
method	28, 29, 31, 35, 37, 38	sample pump	17
result limits	38	sample pump mode	17
rinse	30	syringe mode	15
standard	39	delete	
test	29, 35, 37, 38, 39	method	29
tolerance	39	results	52
user-guidance	29	density	31, 32
configure example		temperature compensated	32
adjustment	67	density determination	
guided	42, 45	typical phases	42
sample pump	42, 48	density meter	
syringe	45	charge	13
test	60	functions	9
unguided	48	overview	9
confirm settings	41	power rating	70
		ship	64

shut down	40	fahrenheit	24
site requirements	7, 70	fill	
start up	40	background	42
store	64	guided	44, 47
unpack	13	sample pump	50
density of pure water		syringe	47
table	77	unguided	50
density unit		fill speed	30
adjustment	24	filling tube	
change	24	assemble	16
test	24	install	16
dimensions		material	70
charging stand	72	remove	63
protective cover	73	replace	63
dimming		firmware	
automatic	25	open source attribution files	5
Disposal	64	third party licenses	5
download		update	64
Reference Manual	12	version	5, 69
drain		G	
background	42	Godex MX20	
configure	30	connect	22, 23
guided	45, 48	H	
sample pump	50	housing	
syringe	48	clean	56
unguided	50	material	70
E		I	
EasyDirect Density & Refractometry		icons	11
connect	23	identifier	
export	38	sample	37
energy saver	25	standard	37
enter		install	
number	41	battery	62
special character	41	charging stand	20
symbol	41	connector for syringe	18
text	41	filling tube	16
errors		protective cover	21
list	65	sample pump	64
export		syringe adapter	18
results	38, 53	intended use	7
F		K	
factory adjustment		key functions ribbon	11
restore	66	keyboard	
factory settings		change layout	27
restore	69		

overview	41
L	
language	
change	27
locking screw	
material	70
M	
maintenance	
schedule	56
material	
charging stand	72
connector for filling tube or syringe	70
filling tube	70
housing	70
locking screw	70
measuring cell	70
protective cover	59, 73
sample pump	70
measurement accuracy	
adjust	66
check	60
measurement reliability	
configure	31
measuring cell	10
clean	56
material	70
measuring ribbon	11
method	
change	28
configure	29, 31, 35, 37, 38
create	28
delete	29
stop	50
method ribbon	11
N	
navigate	41
number	
enter	41
O	
open source attribution files	5
overview	
functions	9
instrument	9
keyboard	41

P	
parts	74
password	
change	26
password protection	25
activate	26
deactivate	26
perform	
adjustment	66
test	60
perform example	
adjustment	67
method	44, 46, 49
test	61
Plato	33
power rating	
AC/DC adapter	70, 72
charging stand	72
density meter	70
power supply	
connect charging stand	20
connect density meter	14
print	
results	38, 54
printer	
connect	22
test	23
printout	38
protective cover	
clean	59
dimensions	73
install	21
material	59, 73
remove	21
weight	73
R	
read	
barcode reader	52
RFID reader	51
Reference Manual	
download	12
remove	
battery	62
filling tube	63
protective cover	21

sample pump	63	rinse	49
replace		set up	15
battery	62	sample pump mode	
filling tube	63	activate	15
sample pump	63	deactivate	17
reset		sampling ribbon	11
factory adjustment	66	scope of delivery	12
factory settings	69	screen	
restore		auto rotation	26
factory settings	69	brightness	25
factory adjustment	66	dimming	25
result limits		serial number	69
configure	38	set up	
results		sample pump	15
delete	52	syringe	17
export	38, 53	ship	
print	38, 54	density meter	64
view	52	shortcut method list	41
RFID reader		shut down	
activate	16	automatic	25
deactivate	16	density meter	40
metal container	51	Signal words	7
read	51	site requirements	
write	51	charging stand	19, 72
rinse		density meter	7, 70
background	42	Smart Tag	
configure	30	metal container	51
guided	44, 46	read	51
sample pump	49	write	51
syringe	46	spare parts	74
unguided	49	special character	
S		enter	41
Safety information		specific gravity	31, 32
Signal words	7	temperature compensated	32
Warning symbols	7	standard	
sample		configure	39
identifier	37	identifier	37
sample ID ribbon	11	standard delivery	12
sample pump		start adjustment	
deactivate	17	guided	67
drain	50	start method	
fill	50	guided	44, 46
install	64	sample pump	44
material	70	syringe	46
remove	63	unguided	49
replace	63		

start test		configure	29, 35, 37, 38, 39
guided	61	configure example	60
start up		density unit	24
density meter	40	failed	61
status ribbon	11	perform	60
stop		perform example	61
method	50	printer	23
store		start guided	61
battery	71	text	
charging stand	72	enter	41
density meter	64	third party licenses	5
sucrose concentration	31	time	
sugar concentration	31	change	24
Brix	33	time format	
Plato	33	change	24
sulfuric acid concentration	31	tolerance	
mass percent	33	configure	39
volume percent	33	troubleshooting	
switch off		list of errors	65
density meter	40	turn off	
switch on		density meter	40
density meter	40	turn on	
symbol	11	density meter	40
enter	41	U	
syringe		unit	
drain	48	density	24
fill	47	temperature	24
rinse	46	unpack	13
set up	17	USB printer	
syringe adapter		connect	22
install	18	user guidance	
syringe mode		configure	29
activate	17	V	
deactivate	15	view	
T		firmware version	69
technical data	70	results	52
temperature compensated		W	
density	32	Warning symbols	7
specific gravity	32	weight	
temperature unit		charging stand	72
change	24	protective cover	73
temperature-compensation coefficient	35	write	
configure	35	RFID reader	51
values	77		
test			
change	29		

To protect your product's future:
METTLER TOLEDO Service assures
the quality, measuring accuracy and
preservation of value of this product
for years to come.

Please request full details about our
attractive terms of service.

www.mt.com/density2go

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contact

Subject to technical changes.
© Mettler-Toledo GmbH 08/2020
30451668B



30451668