


SevenExcellence™

pH/ORP/Ion/Conductivity/DO



METTLER TOLEDO

Table of Contents

1	Introduction	5
2	Safety information	6
2.1	Definitions of signal words and warning symbols	6
2.2	Product specific safety notes	6
3	Design and Function	8
3.1	Description of the instrument	8
3.2	Overview	9
3.3	Module connections	10
3.4	User interface	13
3.4.1	Homescreen	13
4	Putting into Operation	16
4.1	Scope of delivery	16
4.2	Connecting and removing modules	16
4.3	Mounting uPlace™ electrode arm	18
4.4	Installing power supply	19
4.5	Switching the instrument on and off	20
4.6	Initial startup wizard	20
4.7	Connecting sensors	20
4.8	Connecting autosampler	20
4.9	Connecting stirrers	21
4.10	Connecting PnP devices	21
4.11	Connecting computers	21
5	Configuring the instrument	22
5.1	Buffers & standards	22
5.1.1	Predefined buffers and standards	22
5.1.2	User-Defined buffers and standards	23
5.2	Hardware	26
5.2.1	Sensors	27
5.2.2	Automation	30
5.2.2.1	Rondolino	30
5.2.2.2	InMotion	30
5.2.3	Peripherals	30
5.2.3.1	USB stick	31
5.2.3.2	Fingerprint reader	31
5.2.3.3	Printer	31
5.2.3.4	Barcode reader	32
5.2.3.5	Stirrer	32
5.2.4	Communication settings	33
5.2.5	Network settings	33
5.2.6	Auxiliary instruments	33
5.3	User settings	34
5.3.1	Language	34
5.3.2	Screen	35
5.3.3	Beep	35
5.3.4	Shortcuts	35
5.3.5	Keyboards	36
5.4	Global settings	36
5.4.1	System	36
5.4.1.1	Identification	37
5.4.1.2	Date/Time	37
5.4.1.3	Header and footer	38

5.4.2	User management.....	39
5.4.2.1	User	39
5.4.2.2	Account policies.....	40
5.4.2.3	User management settings and login screen.....	41
5.4.2.4	User groups and user rights.....	42
5.4.3	Analysis and resources behavior.....	42
5.4.3.1	Analysis sequence settings.....	43
5.4.3.2	Actions when sensors expire.....	43
5.4.4	Physical properties.....	43
5.4.5	Operation mode.....	44
5.5	Tables.....	44
5.5.1	METTLER TOLEDO tables.....	44
5.5.2	User-defined tables.....	45
5.6	Maintenance & service.....	47
5.6.1	MT service.....	47
5.6.2	Import/export.....	47
5.6.3	Reset to factory settings.....	48
5.6.4	Firmware.....	48
5.6.5	Update.....	48
6	Sensors	49
6.1	Available sensors.....	49
6.2	Sensor parameters.....	49
6.3	Calibration history.....	49
7	Starting Analysis	50
7.1	Starting direct calibration.....	50
7.2	Starting direct measurement.....	53
7.2.1	Starting from module settings.....	53
7.2.2	Starting from home screen.....	56
7.3	Starting methods/series.....	56
7.3.1	Start directly after creating methods/series.....	56
7.3.2	Starting from methods/series list.....	56
7.3.3	Starting from shortcut/direct shortcut.....	57
7.4	Interrupting analysis.....	57
7.5	Errors in the analysis sequences.....	57
7.5.1	Malfunction Types: Error.....	57
7.5.2	Malfunction types: Terminate error.....	57
7.5.3	Malfunction types: Critical error.....	57
7.5.4	Malfunction types: Failed.....	58
8	Methods	59
8.1	Method types.....	59
8.2	METTLER TOLEDO methods.....	60
8.3	Creating methods.....	61
8.3.1	Copying an existing method.....	61
8.3.2	Creating methods based on templates.....	61
8.4	Creating shortcuts for methods.....	62
8.5	Editing methods.....	62
8.6	Deleting methods.....	69
8.7	Method functions.....	71
8.7.1	Overview.....	71
8.7.2	Title.....	72
8.7.3	Configuration.....	72
8.7.4	Sample.....	76
8.7.5	Sample (Calibration).....	76
8.7.6	Sample (Incremental).....	77
8.7.7	Sample (Sensor test).....	78

8.7.8	Blank (BOD)	79
8.7.9	Seeded blank (BOD)	79
8.7.10	Standard (BOD).....	80
8.7.11	Sample (BOD).....	81
8.7.12	Sample verification.....	82
8.7.13	Sensor check.....	82
8.7.14	Measure.....	83
8.7.15	Measure (Calibration).....	86
8.7.16	Measure (Temperature)	87
8.7.17	Measure (Interval).....	88
8.7.18	Measure (Incremental)	91
8.7.19	Measure (Sensor test)	92
8.7.20	Measure (Blank).....	93
8.7.21	Measure (Seeded blank).....	94
8.7.22	Measure (Standard)	94
8.7.23	Measure (BOD)	95
8.7.24	Calibration analysis.....	96
8.7.25	Sensor evaluation	97
8.7.26	Analysis (Blank)	97
8.7.27	Analysis (Seeded blank)	99
8.7.28	Analysis (Standard)	101
8.7.29	Analysis (BOD).....	103
8.7.30	Instruction	106
8.7.31	Wait/Stir.....	106
8.7.32	Calculation.....	107
8.7.33	Report.....	107
8.7.34	Auxiliary instrument.....	108
8.8	Endpoint criteria.....	110
9	Series	111
9.1	Creating series.....	111
9.2	Creating shortcuts for series	111
9.3	Modifying series	112
9.3.1	Change single sample ID	112
9.3.2	Insert samples	112
9.3.3	Delete samples	113
9.4	Delete series.....	113
10	Results	114
10.1	Statuses of measurements	114
10.2	Results of analysis.....	115
10.3	Statistics	115
10.4	Deleting all analysis.....	115
10.5	Deleting single analysis.....	115
10.6	Viewing single analysis data	116
10.7	Printing analysis data.....	116
11	Formula Syntax	117
11.1	Formula symbols.....	117
11.2	Creating formulas	119
11.3	Examples.....	119
11.3.1	Formula in method function calculation	119
11.3.2	Formula in conditions.....	119
11.4	Formula symbols in text	120
12	News, Tasks and Online Display	121
12.1	News	121
12.2	Tasks	121

12.3	Online display	122
13	Maintenance and Care	123
13.1	Cleaning the Instrument	123
13.2	Maintenance of electrodes.....	123
13.3	Transporting the instrument.....	124
13.4	Disposal	124
14	Accessories	125
15	Technical Data	128
15.1	SevenExcellence™	128
15.2	pH/mV module	129
15.3	Conductivity module.....	130
15.4	pH/Ion module	131
15.5	DO/BOD module.....	132
16	Appendix	133
16.1	Predefined buffers & standards	133
16.2	METTLER TOLEDO tables	139
	Index	143

1 Introduction

The METTLER TOLEDO SevenExcellence™ is a modern, professional instrument suitable for use in a wide variety of application areas and segments. It can be used, for example, in quality control as well as in research and development and meets the most demanding requirements.

This flexible meter perfectly combines easy-to-understand operation with a high level of measurement accuracy and outstanding reliability. Thanks to the plug & play capability, the detection of external devices and ISM® sensors occurs automatically, in this way supporting the comfortable handling of SevenExcellence™.

Straightforward user guidance on the large color touch screen makes the operation particularly intuitive. User-definable shortcuts allow methods to be activated directly from the home screen with a single tap.

About this document

The instructions in this document refer to running firmware version **4.3.0** or higher.

The firmware license is subject to the End User License Agreement EULA version 3.0. See the following link for the license text:

▶ www.mt.com/EULA

For third party licenses and open source attribution files, see the following link:

▶ www.mt.com/licenses

If you have any additional questions, contact your authorized METTLER TOLEDO service representative or dealer.

▶ www.mt.com/contact

Conventions and symbols

Note For useful information about the product.

Elements of instructions

Instructions always contain action steps and can contain prerequisites, intermediate results and results. If an instruction contains more than one action step, the action steps are numbered.

- Prerequisites that must be fulfilled before the individual action steps can be executed.

1 Action step 1

➔ Intermediate result

2 Action step 2

➔ Result

2 Safety information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.
- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.



User Manual and Reference Manual are available online.

► www.mt.com/library

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

WARNING A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.

NOTICE A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols



Electrical shock

2.2 Product specific safety notes

Intended use

This instrument is designed to be used by trained staff. The SevenExcellence™ is intended for measuring pH, mV, ion, conductivity, DO/BOD.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes



WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO AC/DC adapter designed for your instrument.
- 2 Keep all electrical cables and connections away from liquids and moisture.
- 3 Check the cables and the plugs for damage and replace damaged cables and plugs.



NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

FCC Rules

This device complies with Part 15 of the FCC Rules and Radio Interference Requirements of the Canadian Department of Communications. Operation is subject to the following conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

3 Design and Function

3.1 Description of the instrument

The SevenExcellence™ bench meters are based on precise electrochemical measuring technique. Modular expansion units, automatically recognized ISM® sensors and the plug & play function for external devices make SevenExcellence a comfortable instrument for reliable measurement.

For the optimal use of your SevenExcellence™ bench meters, the following devices can be connected:

- **Modules**


- METTLER TOLEDO pH/mV module
- METTLER TOLEDO pH/Ion module
- METTLER TOLEDO conductivity module
- METTLER TOLEDO DO/BOD module

- **Sensors**

- Select from the wide range of METTLER TOLEDO electrodes, the electrode value box supports you to find the right tool. www.electrodes.net.
- ISM® sensors are detected automatically.
- Temperature probes can be connected separately.
- Rondolino sample changer allows fully automatic analysis of up to 9 samples. The optional Power-Shower™ perfectly rinses contaminated electrodes.
- InMotion sample changer allows fully automatic analysis of up to 303 samples in different beaker sizes (only if operated by **LabX**).
- Computer for operation with **LabX** or **EasyDirect pH** software.
- uMix™ magnetic stirrer to optimize workflows and to enhance the reproducibility of results.
- Barcode reader for scanning sample data.
- Logstraight™ fingerprint reader for user identification.
- Compact printer (USB-P25), to print out results.
- USB stick to store and transfer data.
- LAN for printing with network printer.
- Additional compact printers, including RS-P25, RS-P26, RS-P28.


3.2 Overview

Front view



1	On/Off button	2	USB A interface
3	USB A interface	4	Touch screen
5	Info Opens a screen with information about the instrument and the attached modules.	6	Home Brings you back to the homescreen.


Rear view



1	LAN	Network connection
2	USB	USB B interface for PC connection (LabX or EasyDirect pH software)
3	Data	RS232 interface
4	Stir	Socket for METTLER TOLEDO stirrer (Mini-DIN)
5	DC	Power supply socket

Pin Assignments


Below the PIN assignments for the RS-232 interface are shown. METTLER TOLEDO printers such as RS-P25 can be connected to this interface.



Pin 1	NC	Pin 6	NC
Pin 2	TxD (out)	Pin 7	NC
Pin 3	RxD (in)	Pin 8	NC
Pin 4	NC	Pin 9	NC
Pin 5	RSGND		

3.3 Module connections


Instrument connection



In the lower area you find the socket for connecting the module to the instrument. A guide rail on the backside of the housing supports attaching the module.

Connections on the pH/mV module


The type of module is identified by a colored tag. Yellow for pH/mV.



Digi	Mini-LTW socket for digital sensors
ISFET	Mini-DIN socket for ISFET sensors
ATC int	RCA (Cinch) socket for internal temperature input (NTC30k)
pH	BNC socket for mV/pH sensors

Connections on the conductivity module


The type of module is identified by a colored tag. Purple for conductivity.



- Digi** Mini-LTW socket for digital sensors
- ATC ext** RCA (Cinch) socket for external temperature input (NTC30k or PT1000)
- Cond** Mini-DIN socket for conductivity sensors

Connections on the pH/Ion module


The type of module is identified by a colored tag. Orange for pH/Ion.



- Digi** Mini-LTW socket for digital sensors
- ATC ext** RCA (Cinch) socket for external temperature input (NTC30k or PT1000)
- ref** 2 mm banana socket for reference electrode
- ATC int** RCA (Cinch) socket for internal temperature input (NTC30k)
- pH** BNC socket for mV/pH sensors

Connections on the DO/BOD module

The type of module is identified by a colored tag. Blue for DO/BOD.



Digi	Mini-LTW socket for digital sensors
ATC ext	RCA (Cinch) socket for external temperature input (NTC22k)
ATC int	RCA (Cinch) socket for internal temperature input (NTC22k)
DO	BNC socket for the DO sensorss

3.4 User interface




NOTICE

Danger of damaging the touch screen with pointed or sharp objects!

Pressing on the touch screen with pointed or sharp objects may be damaged it.

- Operate the touch screen by applying gentle pressure with the pad of your finger.

3.4.1 Homescreen



Name	Explanation
1 Shortcuts	User-specific shortcuts for frequently used methods. Shortcuts are saved in the user profile and can be defined, changed and deleted by the user.
2 Modules	Attached modules and connected electrodes are shown.
3 Status bar	The status bar contains the current menu item, user name as well as date and time.
4 News	Indicates important changes of the system like expiration of resources etc. or disconnected PnP devices.
5 Tasks	Shows the current working status of the instrument. The list of tasks can be accessed.
6 Read	Initiate a direct measurement from the Home Screen, without the need to first set up a measurement method.
7 User data	Access an overview of the currently logged in user, the user group, and the time since the current user has been logged in.
8 Log out	Allows you to log out of the instrument, change user levels, or lock the instrument with a personalized password.
9 Menus	<p>Methods Create and handle methods for every measurement type.</p> <p>Series Create and manage series of individual samples</p> <p>Results Display all measurement results, print out or export them. Find detail information about every single result.</p> <p>Setup Define all system settings in this menu, e.g. hardware settings, user management or user preferences. These settings are usually made during installation of the instrument.</p> <p>Sensors Modify and manage sensors.</p>

Footer buttons

Depending on the selected submenu there are specific buttons in the footer.

AddToHome	Create a shortcut of a quick analysis.
Back	Go back one step in the menu structure.
Cancel	Cancel the current input without saving.
Delete method	Delete the selected method. Shortcuts and series which reference to this method will also be deleted.
Expired resources	Shows the expired resources connected to this instrument.
Insert	Insert a method function in an existing method.
Log out	Log out the current user.
New	Creates a new method, sample series, sensor etc.
OK	Confirm the entered settings.
Preview	When entering a text with an alphanumeric keypad, you can preview your input.
Print	Print out the currently shown screen in tabular form. A printer must be connected to the instrument.
Samples	Open the list with your defined samples.
Read	Start a quick analysis directly from the homescreen.
Start	Start a selected method.
User data	Access an overview of the currently logged in user

Input field types

In the single menus you have different kinds of parameter-fields to enter information, values or names or to choose an option in a drop-down list. Depending on the input field (value or name) you will get a numeric or alphabetic keypad for entering your values.



Text input field

Any text comprised of letters (up to 30 characters), numbers and symbols can be entered into these fields.



Text input field (extended)

Any text comprised of letters (up to 500 characters), numbers and symbols can be entered into these fields.



Number input field

Numbers and formulas can be entered into these fields.



Drop-down list

A drop-down list opens from which you can select an entry.



Shortlist

A shortlist opens from which you can select an entry.



Formula field

A formula must be entered in these fields.




Info field

The displayed content is only for information (read only).

Keypads

Alphabetic keypad



- Tap (1) to preview the entry.
- Tap (2) to enter capital letters.
- Tap (3) to enter lowercase letters.
- Tap (4) to switch to a numeric keypad and (2) to turn back to alphanumeric.
- Tap (5) to delete all entered letters or numbers.
- Tap (6) to delete the last entered letter or number.

Numeric keypad



- Tap (1) to delete all entered numbers.
- Tap (2) to delete the last entered number.

If the user language is set to Chinese, an additional button **CN** appears in the alphabetic keypad at the end of the input field. Tap it to enter Chinese character in Pinyin and select the character from the list below the input field. Be aware that Chinese characters cannot be printed on compact printers.

4 Putting into Operation

4.1 Scope of delivery

Unpack the instrument and check the scope of delivery. Keep the calibration certificate in a safe place.

SevenExcellence™ is delivered with:

- uPlace™ electrode arm
- Measurement modules and blank modules (according to your order)
- Sensors (kit version or according to your order)
- Universal AC adapter
- User Manual
- Declaration of conformity
- Calibration certificate
- EasyDirect pH PC software

4.2 Connecting and removing modules

Note

- Before connecting or removing modules switch off the instrument.


METTLER TOLEDO offers different types of modules. The type of module is identified by a colored tag.

The instrument can be equipped with up to three modules. The modules can be combined in any composition and order according to your requirements.


With less than three modules installed, protect the open sockets against moisture or contamination by installing a blank module.

Connecting modules

- 1 Insert the expansion unit into the slot. The guide rails simplify the connecting.




- 2 Attach the cover.




Removing modules

- 1 To remove the cover, press the clips on both sides and lift it up.




- 2 To remove a module, pull it out of the slot.




4.3 Mounting uPlace™ electrode arm

The electrode arm can be used as stand alone or it can be attached to the instrument on the left or right side, according to your preferences. The height of the electrode arm can be varied by using the extension shaft part. Use the wrench to attach the extension part .





Assembly of the electrode arm


- 1 Use the wrench to attach the base to the electrode arm by tightening the screws. The electrode arm can now be used in the stand alone mode.



- Then insert the foot of the meter to the arm base and shift the meter in the direction of the arrow to make the foot fit.



- Use the lock screw to attach the meter to the base of the arm.



4.4 Installing power supply



WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- Only use the METTLER TOLEDO AC/DC adapter designed for your instrument.
- Keep all electrical cables and connections away from liquids and moisture.
- Check the cables and the plugs for damage and replace damaged cables and plugs.



NOTICE


Danger of damage to the AC adapter due to overheating!

If the AC adapter is covered or in a container, it is not sufficiently cooled and overheats.

- Do not cover the AC adapter.
- Do not put the AC adapter in a container.

The instrument is operated using an AC adapter. The AC adapter is suitable for all supply line voltages ranging from 100...240 V AC \pm 10% and 50-60 Hz.

- Insert the correct connector plug into the AC adapter until it is completely inserted.
- Connect the cable of the AC adapter with the DC socket of the instrument.
- Install the cables in such a way that they cannot be damaged or interfere with operation.



- 4 Insert the plug of the AC/DC adapter into a power outlet that is easily accessible.
 - ➔ To remove the connector plug, push the release button and withdraw the connector plug.

4.5 Switching the instrument on and off

Switching the instrument on

- 1 Press the **On/Off** button.
 - ➔ The power button lights up green.
- 2 The instrument is initialized and a system test is performed. This process will take some time.
 - ➔ If only one user is defined, the instrument starts directly with the home screen without showing the login screen.

Switching the instrument off

- 1 Tap **Home > Log out > Shut down**.
 - or -
 - Press the **On/Off** button.
- 2 The instrument stops running tasks and shuts down. This process will take some time.
 - ➔ When the light goes out, the instrument is in the standby mode.

Note

- In the standby mode, the control circuit for the **On/Off** switch is energized. The rest of the instrument is no longer energized.

4.6 Initial startup wizard

At first startup of the instrument or after a factory reset an intuitive wizard guides you through the most important settings of the instrument. In three steps you can select a language, date, time and region as well as if you would like to have predefined shortcuts on the Home screen. These shortcuts refer to predefined METTLER TOLEDO methods and enable an immediate start with one click. If you enable predefined shortcuts the Home screen shows a calibration and a measurement method shortcut per measurement parameter. Depending on the selected region, the default pH buffer group for calibration is defined.

4.7 Connecting sensors

Connect the sensor and make sure that the plugs are properly inserted. For sensors with a built-in temperature probe, connect the second cable to the appropriate ATC socket (**ATC int**). When using separate temperature sensors, connect the cable to the appropriate ATC socket (**ATC ext**). Twist the RCA (Cinch) plug to ease the attachment of the sensor.

ISM[®] sensor

ISM[®] sensors are automatically recognized upon connection.

When connecting an ISM[®] sensor to the instrument, verify that no analysis is running. In case an analysis is running, the sensor identification will be displayed after the endpoint of the analysis is reached.

4.8 Connecting autosampler

InMotion autosampler

With the InMotion autosampler the automatic measurement of up to 303 samples is possible. InMotion is only supported when the instrument is connected to the **LabX** software. The InMotion autosampler is connected to the instrument with a USB cable and can only be used in **LabX** connected mode. Several different configurations are available for the InMotion autosampler depending on your requirements.

Rondolino sample changer

With the Rondolino sample changer the automatic measurement of up to 9 samples is possible. To connect the Rondolino sample changer to the instrument, a "Rondolino USB-TTL box" is required. Use the USB cable to connect the USB A interface of the instrument to the USB-TTL box. Use the 9-pin cable to connect the the Rondolino sample changer with the TTL box. Connect the compact stirrer to the stirrer adapter cable and then to the instrument. The compact stirrer is an overhead stirrer.

See also

 Rondolino ▶ Page 30

 InMotion ▶ Page 30

4.9 Connecting stirrers

Connect uMix™, the METTLER TOLEDO external magnetic stirrer, or the METTLER TOLEDO Compact Stirrer (overhead stirrer, adapter cable needed) to the instrument. The stirrers are powered by the instrument and will be automatically switched on/off during an analysis according to the settings. Before connecting a stirrer, switch off the instrument. Use the stirrer socket (**Stir**) on the rear side, note the designation. Note the corresponding chapter of the operating instructions concerning the settings of peripheral devices.

4.10 Connecting PnP devices

USB sticks, barcode reader, US-P25 printer and fingerprint reader can be connected to the instrument via the USB connections on the left side of the instrument for immediate use. The instrument has a plug & play detection to support the comfortable handling of peripheral devices.

4.11 Connecting computers

The instrument can be connected to a PC with **LabX** or **EasyDirect pH** software.

Connected to the PC software **LabX** via USB or Ethernet, SevenExcellence profits from additional functionalities such as data and resource management, report creation and ERP/LIMS integration. The instrument can be controlled via touchscreen and PC software. In addition, **LabX** offers regulation support (e.g. 21 CFR part 11) and allows combined analysis with other METTLER TOLEDO laboratory instruments.

Find out more and visit us on:

▶ www.mt.com/LabX

The PC software **EasyDirect pH** allows for easy and efficient transfer of results to a PC. SevenExcellence is connected via USB, one instrument per PC.

For simplified data transmission to a PC or network via RS232, the RS232 Compact Printer can be selected in the printer settings.

5 Configuring the instrument

Navigation: Home > Setup

This chapter provides you information on setting up the instrument according to your requirements.

The following summary shows the buttons available for the various setting options.

Buffers & Standards	Buffers & Standards
Hardware	Sensors
	Automation
	Peripherals
	Communication settings
	Network settings
	Auxiliary instruments
User settings	Language
	Screen
	Beep
	Shortcuts
	Keyboards
Global settings	System
	User management
	Analysis and resources behavior
	Physical properties
	Operation mode
Tables	Tables
Maintenance & Service	MT-Service
	Import / Export
	Reset to factory settings
	Firmware
	Update

5.1 Buffers & standards

This section describes how to manage buffers and standards. You can use predefined buffers and standards or you can create, edit or delete user-defined buffers and standards.

A maximum of 20 user-defined buffer sets and standard groups can be created with up to 20 temperature-dependent values per buffer or standard. A total of 13 buffer sets and standard groups for pH, ion and conductivity calibration are predefined.

The predefined pH buffer sets offer the option to create individual buffer sets based on a mixture of different commercial buffers.

Note

- Before you perform a calibration, buffers or standards must be defined.

5.1.1 Predefined buffers and standards

The button [**Buffers & Standards**] opens the list of predefined buffer sets and standards. The predefined buffer sets and standards cannot be modified or deleted. The parameters of the following buffers and standards are stored, see in the Appendix [Predefined buffers & standards ▶ Page 133].

Type	Name
pH	METTLER TOLEDO USA (Ref. 25 °C)
	METTLER TOLEDO Europe (Ref. 25 °C)
	METTLER TOLEDO Verification Buffers (Ref. 25°C)
	MERCK (Ref. 20 °C)
	DIN(19266)/NIST (Ref. 25 °C)
	DIN(19267) (Ref. 25 °C)
	JJG119 (Chinese) (Ref. 25 °C)
	Technical (Ref. 25 °C)
Conductivity	International (Ref. 25 °C)
	Saturated NaCl (Ref. 25 °C)
	Chinese (Ref. 25 °C)
	Japanese (Ref. 20 °C)
Ion	METTLER TOLEDO ION (Ref. 25 °C)

- 1 Tap a buffer set or standard group.
 - ➔ **Name of the buffer set or standard** is displayed with more detailed information about the stored values.
- 2 Tap a value.
 - ➔ **Temp. related value** is displayed.
The parameters **Temperature** and **Value** are shown.

Note

This list has information character concerning the predefined buffers and standards. When user-defined buffers and standards are defined, they are shown in this list too.

5.1.2 User-Defined buffers and standards

The button [**Buffers & Standards**] leads to the list of the stored buffer sets and standard groups. When creating a user-defined buffer set or standard group you can make use of the already existing buffer sets and standard groups. In this way you can create new individual sets and groups with buffers and standards from different groups. User-defined buffers and standards will be added to the list of predefined buffers and standards.

Creating user-defined buffer and standard group

Navigation: Home > Setup > Buffers & Standards > New

Parameter	Description	Values
Type	Defines the type of the buffer or standard for the calibration.	pH Ion Conductivity
Name	Defines the name of the buffer or standard group. The name must be unique and must not contain any special characters or spaces.	-
Unit	For preparing a pH buffer group, pH as standard value is shown as information. For preparing an ion buffer group, select a unit.	pH mmol/L mol/L mg/L ppm % pX
Reference temperature	Defines the temperature of the nominal buffer or standard value (usually 20 °C or 25 °C).	15.0...35.0 °C
Use predefined buffers	Active: This group is used to mix buffers from existing groups. Inactive: User-defined pH buffers with temperature dependent values are entered.	Active Inactive

- 1 In **Type** select the type of the buffer or standard.
- 2 Give the buffer or standard group a meaningful name.

➔ A name **Calibration Standard** with consecutive number is entered automatically, it can be overwritten.

3 In **Reference temperature** enter a temperature.

4 If type **Ion** is selected, select a unit.

- or -

If type **pH** is selected and you want to assemble a buffer set from stored buffer sets, activate **Use predefined buffers**.

5 Tap [**Save**].

➔ Now you have created a user-defined buffer or standard group. The name of the user-defined buffer or standard group is displayed.

In the next step values need to be determined for:

- User-defined pH buffer groups based on predefined buffer sets
- User-defined pH buffer groups
- User-defined ion standards
- User-defined conductivity standards

Parameter	Description	Values
Reference temperature	Defines the temperature of the nominal buffer or standard value (usually 20 °C or 25 °C).	15.0...35.0 °C

Add values for user-defined pH buffer groups based on predefined buffers

Navigation: Home > Setup > Buffers & Standards > My group > New

Parameter	Description	Values
Buffer group	Opens the list of predefined buffers and standards. Displayed if Use predefined buffers is activated.	-
pH value	Opens the list of predefined values. Displayed if Use predefined buffers is selected.	-
Reference temperature	Information on the defined reference temperature.	-

1 Tap the list field **Buffer group**.

➔ **Select buffer / standard group** is displayed.

2 Select a buffer group.

3 Tap the list field **pH value**.

➔ **Select buffer / standard** is displayed.

4 Select a pH value.

5 Tap [**OK**]

➔ The pH value at reference temperature is displayed in the list.

6 To enter additional values, tap [**New**] and repeat the previous steps.

- or -

To finish and save the entries tap [**Save**].

Add values for user-defined pH buffer groups

Navigation: Home > Setup > Buffers & Standards > My group > New

Parameter	Description	Values
Reference temperature	Information on the defined reference temperature.	-
pH value	Defines the pH value of the calibration standard to be used.	-

1 Enter a pH value and confirm with [**OK**].

- ➔ The nominal value of the user-defined buffer or standard is displayed. The value is set.
- 2 To enter temperature related values for this standard, select it.
- 3 Tap **[New]** and enter the values in **Reference temperature** and **pH value**, then tap **[OK]**. Repeat this step for every temperature, then tap **[OK]**.
- 4 To enter additional values, tap **[New]** and repeat the previous steps.
 - or -
 - To finish and save the entries tap **[Save]**.

Add values for user-defined ion standards

Navigation: Home > Setup > Buffers & Standards > My group > New

Parameter	Description	Values
Reference temperature	Information on the defined reference temperature.	-
Standard value	Defines the value of the calibration standard to be used.	-

- 1 Enter a standard value and confirm with **[OK]**.
 - ➔ The nominal value of the standard is displayed. The value is set.
- 2 To enter temperature related values for this standard, select it.
- 3 Tap **[New]** and enter the values in **Reference temperature** and **Standard value**, then tap **[OK]**. Repeat this step for every temperature, then tap **[OK]**.
- 4 To enter additional values, tap **[New]** and repeat the previous steps.
 - or -
 - To finish and save the entries tap **[Save]**.

Add values for user-defined conductivity standard groups

Navigation: Home > Setup > Buffers & Standards > My group > New

Parameter	Description	Values
Unit	Defines the measurement unit for conductivity.	μS/cm mS/cm S/m μS/m mS/m
Reference temperature	Information on the defined reference temperature.	-
Conductivity	Defines the value for the calibration standard to be used in the unit selected in Unit .	-

- 1 In **Unit** select a type to which your calibration refers to.
- 2 Enter a conductivity value and confirm with **[OK]**.
 - ➔ The nominal value of the standard is displayed. The value is set.
- 3 To enter temperature related values for this standard, select it.
- 4 Tap **[New]** and enter the values in **Reference temperature** and **Conductivity**, then tap **[OK]**. Repeat this step for every temperature, then tap **[OK]**.
- 5 To enter additional values, tap **[New]** and repeat the previous steps.
 - or -
 - To finish and save the entries tap **[Save]**.

Modify user-defined buffers and standards

Navigation: Home > Setup > Buffers & Standards

- 1 Select the appropriate buffer set or standard group from the **Buffer / Standard group** list.
 - ➔ Name of the user-defined buffer or standard group is displayed.
- 2 Tap the standard to be modified.

- ➔ **Temp. related value** is displayed.
- 3 Tap the temperature related value to be modified.
 - ➔ The value, depending on the selected buffer or standard is displayed.
- 4 Modify the value and confirm with [**OK**].
- 5 Confirm with [**OK**].
 - ➔ **Temp. related value** is displayed.
- 6 To modify additional values repeat previous steps.
- 7 To finish, tap [**Back**].
 - ➔ Name of the user-defined buffer or standard group is displayed.
- 8 To store the entries tap [**Save**].

Note

Values of user-defined buffer groups created from predefined buffer groups can not be modified.

Delete user-defined buffers and standards

Navigation: Home > Setup > Buffers & Standards

- 1 Select the appropriate buffer set or standard group.
 - ➔ Name of the user-defined buffer or standard group is displayed.
- 2 Tap the value to be deleted.
 - ➔ **Temp. related value** is displayed.
- 3 Tap [**Delete**].
- 4 To finish and store entries tap [**Save**].

Note

When a value of a buffer set or standard group is deleted, methods that refer to the deleted buffer set or standard group are no longer executable.

Delete user-defined buffer and standard groups

Navigation: Home > Setup > Buffers & Standards

- 1 Select the appropriate buffer set or standard group.
 - ➔ Name of the user-defined buffer or standard group is displayed.
- 2 Tap [**Parameters**] to open the parameters.
- 3 Tap [**Delete**].

Note

When buffer sets or standard groups are deleted, methods that refer to the deleted buffer sets or standard groups are no longer executable.

Modify the name of a user-defined buffer or standard group

Navigation: Home > Setup > Buffers & Standards

- 1 Select the appropriate buffer set or standard group.
 - ➔ Name of the user-defined buffer or standard group is displayed.
- 2 Tap [**Parameters**] to open the parameters.
- 3 Overwrite the name for the user-defined buffer set or standard group.
- 4 To finish and store entries, tap [**Save**].

5.2 Hardware

This section describes how to configure the hardware components connected to the instrument, such as:

- **Sensors**

- **Automation**
- **Peripherals**
- **Communication settings**
- **Network settings**

5.2.1 Sensors

Sensors for pH, redox, ISFET, conductivity, dissolved oxygen, ion and temperature measurement can be connected to the instrument. One default sensor per sensor type is predefined.

Sensors can be added, modified and deleted. The usable life and life span of a sensor can be set and monitored (for redox and temperature sensors only the life span can be monitored). You can display or print out the calibration history.

Note

Digital ISM[®] sensors are detected automatically. The factory calibration will be transferred when ISM[®] sensors are used.

Navigation: Home > Setup > Hardware > Sensors

Parameter	Description	Values
Type	Information on the type of measurement.	-
Name	Information on the name of the sensor.	-
Serial number	Information on the serial number of the sensor.	-
Module	Module to which the sensor is attached. ISM [®] sensors are detected automatically. For ISM [®] sensors A-ISM , B-ISM or C-ISM is displayed. For non ISM sensors A , B or C is displayed. Displayed if a module is connected to the instrument and a sensor is connected to a module.	-

Adding sensors

Navigation: Home > Setup > Hardware > Sensors > New

Parameter	Description	Values
Sensor type	Select a sensor type according to the measurement type.	pH Redox ISFET Conductivity DO Ion
Temperature signal	Select a type according to the specifications of the temperature sensor. Displayed if Sensor type = Temperature .	NTC30kOhm Pt1000
Sensor name	Define a name for the sensor.	-
Serial number	Each sensor has a serial number. The serial number of ISM [®] sensors are detected automatically.	-
Ion	The ion charge and molar mass values of predefined ions are stored in the instrument. Displayed if Sensor type = Ion .	Ag ⁺ Ca ²⁺ Cl ⁻ CN ⁻ Cu ²⁺ F ⁻ I ⁻ K ⁺ Na ⁺ NO ₃ ⁻ Pb ²⁺ S ²⁻ Other
Molar mass	Molar mass of an ion is calculated by adding the atomic weight of the individual atoms the ion is composed of. This parameter is needed for calculations of the units mol/L and mmol/L. Displayed if Ion = Other .	-
Ion charge	Ion charge of the ion to be detected. Displayed if Ion = Other .	-3 -2 -1 +1 +2 +3

Detection	Information about the detection mode of a DO sensor. Manually added sensors always have Polarographic detection. Displayed if Sensor type = DO .	-
Cell constant	Conductivity measurement depends on the cell constant of the probe. The precise cell constant is noted in the certificate of the probe. Displayed if Sensor type = Conductivity .	-
Cal. mode	Information on calibration mode. Displayed if Sensor type = pH or ISFET or Ion .	-
Offset	Information on offset. Displayed if Sensor type = pH or ISFET or Ion .	-
Slope	Information on the slope. Displayed if Sensor type = pH or ISFET or Ion or DO .	-
Calibration procedure	Information on calibration procedure. Displayed if Sensor type = pH or ISFET or Conductivity or Ion or DO or Temperature .	-
Date/Time of calibration	Information on date and time of calibration. Displayed if Sensor type = pH or ISFET or Conductivity or Ion or DO or Temperature .	-
Calibrated by	Name of the user performing the calibration. Displayed if Sensor type = pH or ISFET or Conductivity or Ion or DO or Temperature .	-
Monitoring usable life	Specifies whether the usable life is to be monitored. Displayed if Sensor type = pH or ISFET or Conductivity or Ion or DO .	Active Inactive
Time period	Defines the unit for the period of time. Displayed if Monitoring usable life is activated.	Days Hours
Usable life	Defines the period in days or hours up to next calibration. Displayed if Monitoring usable life is activated.	-
Expiry date	Information on the sensors lifetime, calculated from the predefined parameters. Displayed if Monitoring usable life is activated.	-
Monitoring life span	Based on your requirements a period of time for usable life can be defined. Not displayed for optical DO sensor because its sensor cap anyway has a limited life span of one year.	-
Initial operation	Here you can enter the date of initial commissioning of the sensor. Displayed if Monitoring life span is activated.	-
Life span	Describes the period of time after which a resource is consumed and should be replaced. Displayed if Monitoring life span is activated.	-
Expiry date	Information on the end of the life span, calculated on base of the predefined parameters. Displayed if Monitoring life span is activated.	-

- 1 Select **Sensor type**. Different parameters and information are displayed depending on the sensor type.
 - ➔ If **Temperature** is defined, select a unit in **Temperature signal**.
 - ➔ If **Ion** is defined, select a predefined parameter or **Other** in **Ion type**.
 - ➔ If **Other** is defined, enter a value in **Molar mass** and select a value in **Ion charge** of the sample to be analyzed.
 - ➔ If **Conductivity** is selected, enter the value in **Cell constant** of your conductivity probe (found in the certificate of the probe). The unit is predefined.

- 2 Enter a sensor name. The name **Sensor** with consecutive number is entered automatically. The name in **Sensor name** is unique. If an already existing name is used, a warning will be displayed. Tap **[OK]** and change the name.
- 3 Enter the serial number.
- 4 To monitor the usable life of the sensor activate **Monitoring usable life**. Enter days or hours in **Time period** and **Usable life**. **Expiry date** is shown as information.
- 5 To monitor the life span of the sensor, activate **Monitoring life span**. Enter the date and time in **Initial operation**. In **Life span**, enter the number of months until the sensor must be replaced. **Expiry date** is shown as information.
- 6 To finish and store the entries tap **[Save]**.
 - ➔ The sensor was added to the list.

Modifying sensors

Navigation: Home > Setup > Hardware > Sensors

The parameters that can be modified are:

- **Sensor name**
- **Monitoring usable life**
- **Monitoring life span**
- **Cell constant** (conductivity only)

- 1 Tap the sensor to be modified.
 - ➔ **Parameters** is displayed.
- 2 Modify the required parameters.
- 3 To finish and store the entries tap **[Save]**.

Note

- The type of the sensor cannot be modified
- A sensor that is currently in use cannot be modified.
- Changing the sensor name does not create a new sensor in the list.

Deleting sensors

Navigation: Home > Setup > Hardware > Sensors

- 1 Tap the sensor to be deleted.
 - ➔ **Parameters** is displayed.
- 2 Tap **[Delete]**.
 - ➔ The sensor has been deleted without warning.

Note

- The last sensor of a type cannot be deleted.
- A sensor that is currently in use cannot be deleted.
- Predefined sensors cannot be deleted.

Show and print calibration history

Navigation: Home > Setup > Hardware > Sensors

- 1 Tap the required sensor.
 - ➔ **Parameters** is displayed.
- 2 Tap **[Calibration history]** to display the history.
- 3 To print the history, tap **[Print]**.
 - ➔ The history is printed.

5.2.2 Automation

The following automation units can be installed.

5.2.2.1 Rondolino

The instrument has PnP (Plug & Play) identification. As soon as Rondolino is connected, the appropriate information is displayed.

Navigation: Home > Setup > Hardware > Automation > Rondolino

Parameter	Description	Values
Type	Shows the instrument type.	-
Status	Information if Rondolino is installed or not.	Installed Not installed

5.2.2.2 InMotion

The instrument has PnP (Plug & Play) identification. As soon as InMotion autosampler is connected, the appropriate information is displayed.

The InMotion Autosamplers can hold multiple samples depending on the base unit (**Flex**, **Pro** or **Max**) and the accompanying sample rack. The InMotion Autosampler connects to a USB interface on the instrument. InMotion autosamplers and attached towers are PnP devices that are automatically recognized and installed upon connection to the instrument.

InMotion autosamplers can only be used if the instrument is connected to the **LabX** PC software. Only one tower can be used.

Navigation: Home > Setup > Hardware > Automation > InMotion

Parameter	Description	Values
Type	Shows the instrument type.	-
Base type	Indicates the sample changer type.	Flex Pro Max
Status	Information if an InMotion autosampler is installed.	Installed Not installed Installed, not supported
Stirrer output	Defines the stirrer output on the InMotion where the stirrer is connected.	InMotion/Stirrer1 ... InMotion/Stirrer6
PowerShower output	Specifies the pump connection on the InMotion where the PowerShower pump is connected.	InMotion/Pump1 ... InMotion/Pump6
Rate	Specifies the pump rate for the PowerShower pump [mL/min].	0.1 ... 1000
Drain pump output	Specifies the pump connection on the InMotion where the drain pump is connected.	InMotion/Pump1 ... InMotion/Pump6
Rate	Specifies the pump rate for the drain pump [mL/min].	0.1 ... 1000
CoverUp	Information if lid handling is installed.	Installed Not installed
Rack	Indicates the type of the installed rack. Standard: Standard rack. Water bath: Rack including a wather bath.	Standard Water bath
Rack size	Indicates the size of the installed rack.	Number of positions on the rack
Beaker height	Defines the beaker height [mm].	65...215
Barcode reader	Information if a barcode reader is installed.	Installed Not installed

5.2.3 Peripherals

Navigation: Home > Setup > Hardware > Peripherals

In the dialog **Peripherals**, the following devices and settings can be configured:

- **USB stick**

- **Fingerprint reader**
- **Printer**
- **Barcode reader**
- **Stirrer**

The following chapter describes the management of the various peripherals which can be connected to the instrument. This includes the instruction how to add peripherals and how to change the parameters.

5.2.3.1 USB stick

Navigation: Home > Setup > Hardware > Peripherals > USB stick

Commercially available USB sticks of USB Version 1.1 are supported with FAT12, FAT16, or FAT32, or exFAT (but not NTFS) file system. Two USB sockets are located on the right side of the instrument. As soon as the USB stick is connected, the instrument recognizes the peripheral device. The status field informs if the external memory is installed or not.

5.2.3.2 Fingerprint reader

Connect the fingerprint reader to the USB output of the instrument. Two USB sockets are located on the right side of the instrument. Before the fingerprint reader can be used, the corresponding functions need to be configured, see [Account policies ▶ Page 40].

Navigation: Home > Setup > Hardware > Peripherals > Fingerprint reader

Parameter	Description	Values
Activate fingerprint reader	Activate the fingerprint reader.	Active Inactive
Status	Information that the peripheral device is installed. Displayed if check box Activate fingerprint reader is selected.	Installed Not installed

– To activate the fingerprint reader, select **Activate fingerprint reader**.

➔ A status field is shown and informs that the peripheral device is connected.

5.2.3.3 Printer

The instrument supports different printers and a PDF writer to store the data on a USB stick. The selected type of printer is used for all kind of print-outs, e.g. during execution of a direct measurement or a method with method function Report. The language for print-outs is defined in the user setting as report language.

Navigation: Home > Setup > Hardware > Peripherals > Printer

Parameter	Description	Values
Printer type	Several types of printers and pdf writers are offered to store or print data.	PDF Writer USB Compact Printer RS232 Compact Printer Network printer
Auto baudrate	If activated, the settings for the RS232 connection will be determined automatically.	Active Inactive
Baud rate	Defines the baud rate for data transmission via the RS232 interface.	1200 2400 4800 9600 19200
Number of bits	Defines the number of data bits for RS232 connections.	7 8
Stop bits	Defines the stop bit used for RS232 connections. Read-only if Auto baudrate is enabled.	1 1.5 2
Parity	Defines the parity method for RS232 connections.	Even Odd None
Handshake	Defines the handshake used for RS232 connections.	None Xon/Xoff

Type	Offers different types of command language for the connected network printer.	HP PLC Epson ESC/P2
IP address	Enter the IP address of the network printer.	-
Port number	Enter the port number of the network printer.	-
Paper size	Defines the paper size suitable for the selected network printer.	A4 Letter
PDF storage location	Information on the data storage. Displayed if PDF Writer is selected.	USB stick

- 1 Select a type in **Printer type**
- 2 If **RS232 Compact Printer** is selected and **Auto baudrate** not enabled, enter the corresponding parameters.
- or -
If **Network printer** is selected, enter the corresponding parameters.
- 3 To print out a test page, tap [**Test page**].
- 4 To finish and save the entries tap [**Save**].

See also

- 📖 Language ▶ Page 34
- 📖 Calibration history ▶ Page 49
- 📖 Starting direct calibration ▶ Page 50
- 📖 Starting from module settings ▶ Page 53
- 📖 Measure (Interval) ▶ Page 88
- 📖 Printing analysis data ▶ Page 116

5.2.3.4 Barcode reader

Navigation: Home > Setup > Hardware > Peripherals > Barcode reader

Connect the barcode reader to a USB socket. Two USB sockets are located on the right side of the instrument. The PnP function of the instrument recognizes if a barcode reader is connected. The status field informs that a barcode reader is connected.

5.2.3.5 Stirrer

Connect the stirrer to the corresponding socket at the backside of the housing. The PnP function of the instrument recognizes if a stirrer is connected.

Navigation: Home > Setup > Hardware > Peripherals > [Stirrer]

Parameter	Description	Values
Status	Information if a stirrer is connected.	Installed Not installed
Stirrer type	Select a stirrer from different standard types or define individual parameters for a stirrer.	Magnetic stirrer Overhead stirrer User-defined stirrer
10%	Defines the minimum voltage for the defined stirrer.	-
100%	Defines the maximum voltage for the defined stirrer.	-

- 1 In **Stirrer type** select the type of stirrer to be used.
- 2 Enter a value for the minimum voltage for the current stirrer.
- 3 Enter a value for the maximum voltage for the current stirrer.
- 4 To check the settings of the stirring speed, tap [**Test**].
- 5 To finish and store the entries tap [**Save**].

Note

The stirrer type can only be selected, if a stirrer is connected to the instrument.

5.2.4 Communication settings

SevenExcellence can be connected to a PC for communication with the laboratory software **LabX** or **EasyDirect pH** for remote control.

Navigation: Home > Setup > Hardware > Communication settings

Parameter	Description	Values
Communication	Defines the type of communication. Transfer results to EasyDirect: To transfer results to EasyDirect pH . Remote control at start-up: Activates remote interface via Ethernet and USB. Connect to LabX at start-up: A connection to LabX will be established on startup.	Off Transfer results to EasyDirect Remote control at start-up Connect to LabX at start-up
Connection type	Defines how the instrument is connected to the PC, either via the network connection or via the USB connection. Displayed if Communication = Remote control at start-up or Connect to LabX at start-up is selected.	Ethernet USB
Port number	Defines the port for a network connection of the instrument. Displayed if Connection type = Ethernet is selected.	1024...65535
Host status	Information if a host is connected or not. Displayed if Communication = Transfer results to EasyDirect is selected.	Installed Not installed

5.2.5 Network settings

Configure these settings if you have your instrument connected to a network.

Navigation: Home > Setup > Hardware > Network settings

Parameter	Description	Values
Type	Information on the type of network connection.	Ethernet
Obtain IP address automatically	If activated, the device automatically obtains an IP address.	Active Inactive
IP address	Defines the IP address of the instrument. Only if Obtain IP address automatically is not activated.	000.000.000.000 ... 255.255.255.255
Subnet mask	Defines the subnet mask to link the subnet's IP address, to run the instrument on a local subnetwork. Only if Obtain IP address automatically is not activated.	000.000.000.000 ... 255.255.255.255
Standard gateway	Defines the address of the standard gateway for communication between the various networks. Only if Obtain IP address automatically is not activated.	000.000.000.000 ... 255.255.255.255

5.2.6 Auxiliary instruments

During method execution, SevenExcellence can send data to and get data from an auxiliary instrument which is connected via RS232 cable using a specific USB-RS232 adapter. The communication settings can be defined per auxiliary instrument.

Navigation: Home > Setup > Hardware > Auxiliary instruments

Select an existing auxiliary instrument from the list or tab **New** to create a new one. To delete an auxiliary instrument from the list, select it and tab **Delete**. One item has to remain in the list and cannot be deleted.

Parameter	Description	Values
Control type	Defines the way the auxiliary instrument is connected to the instrument.	USB-RS232
Name	Specify a descriptive name of your choice.	Arbitrary
Baud rate	Defines the baud rate for data transmission via the RS232 interface.	1200 2400 4800 9600 19200
Number of bits	Defines the number of data bits for RS232 connections.	7 8
Stop bits	Defines the stop bit used for RS232 connections.	1
Parity	Defines the parity protocol.	Even Odd None
Handshake	Defines the handshake used for RS232 connections.	None Xon/Xoff

See also

 Accessories ▶ Page 125

5.3 User settings

Navigation: Home > Setup > User settings

The dialog **User settings** contains the settings that can be made specifically for each currently logged in user.

User settings include the following settings:

- **Language**
- **Screen**
- **Beep** for signal tones
- **Shortcuts**
- **Keyboards** (alphanumeric and numeric)

5.3.1 Language

The language can be defined separately for the operation of the user interface as well as for the reports that are to be printed out.

Navigation: Home > Setup > User settings > Language

Parameter	Description	Values
Screen	Defines the language setting for user interface.	English French German Italian Spanish Portuguese Russian Chinese Japanese Korean
Report	Defines the language setting for printout.	English French German Italian Spanish Portuguese Russian

- 1 Select a language in **Screen**.
- 2 Select a language in **Report**.
- 3 To finish and store the entries tap [**Save**].

5.3.2 Screen

The screen can be adapted to individual needs. Activate the screen saver if it is desired and set the time period up to displaying.

Navigation: Home > Setup > User settings > Screen

Parameter	Description	Values
Primary color	Information on the color scheme of the user interface.	Blue
Brightness	Defines the brightness of the display.	50...100%
Screen saver	Activates the screen saver.	Active Inactive
Wait time	Defines how long in [min] the system should wait after the user's last action on the terminal before activating the screen saver.	1...120

- 1 Select a percentage value for **Brightness**.
- 2 To use the screen saver, activate **Screen saver** and set a time for **Wait time**.
- 3 To finish and store the entries tap [**Save**].

5.3.3 Beep

A signal tone can be activated. With the push of a button, in case of error messages, stability of measurement results, additional news or if user interaction is required you hear a beep.

Navigation: Home > Setup > User settings > Beep

Parameter	Description	Values
Push of a button	Enables the audio signal when tapping on the touch screen.	Active Inactive
Error message	Activates a beep in case of error messages.	Active Inactive
Stability signal	Enables the audio signal when a measurement result is stable.	Active Inactive
News	Enables the audio signal when events occur that appear under News .	Active Inactive
User interaction required	Activates a beep when user interaction is required.	Active Inactive

- 1 To hear a signal in certain cases, activate the check boxes according to your requirements.
- 2 To finish and store the entries tap [**Save**].

5.3.4 Shortcuts

Navigation: Home > Setup > User settings > Shortcuts

Shortcuts can be created after the parameters for methods and series have been set, see [Creating shortcuts for methods ▶ Page 62], [Creating shortcuts for series ▶ Page 111]. The administration of shortcuts takes place in the menu **Setup**. This chapter describes how to modify or delete shortcuts.

Modify shortcuts

Shortcuts allow you to start methods, series, and manual operations directly from the home screen. You can place shortcuts on the home screen by pressing the **AddToHome** button. **AddToHome** is located in the start dialog of each method, series and manual operation.

Navigation: Home > Setup > User settings > Shortcuts

- Select a method or series from the list.
 - ➔ **Shortcut parameters** is displayed.

Parameter	Description	Values
Type	If the shortcut refers to a series, Series as type is shown, if it refers to a method, Method as type is displayed as information.	-

Description	A meaningful name for the series or method can be entered and will be displayed on the home screen.	-
Series ID or Method ID	If the shortcut refers to a series the Series ID is show, if it refers to a method, the Method ID is shown.	-
Immediate start	If activated, the method or series starts can be started immediately.	Active Inactive
Homescreen position	Shows the position on the home screen. Nine positions are offered, they are numbered 1-9 from left to right in three rows.	-
Created by	Information about the administrator who created the shortcut.	-

- 1 To change the entry in **Description**, overwrite it.
- 2 To change the options for **Immediate start**, activate or deactivate the check box.
- 3 To change the position, tap the list field in **Homescreen position**.
Touching a free area in **Select location for shortcut** will assign the new home screen position.
- 4 To finish and store the entries tap [**Save**].

Delete shortcuts

Navigation: Home > Setup > User settings > Shortcuts

- 1 Select the relevant shortcut from the list.
➔ **Shortcut parameters** is displayed.
- 2 Tap [**Delete**].
➔ The shortcut has been deleted without warning.

5.3.5 Keyboards

The layout for the alphanumeric and the numeric input fields can be set here.

Navigation: Home > Setup > User settings > Keyboards

Parameter	Description	Values
ABC keyboard	Defines the layout of the alphanumeric input field.	English French German
123 keyboard	Defines the layout of the numeric input field.	Calculator Phone

- 1 Select a layout for **ABC keyboard**.
- 2 Select a layout for **123 keyboard**.
- 3 To finish and store the entries tap [**Save**].

5.4 Global settings

Navigation: Home > Setup > Global settings

In **Global settings**, you can make general instrument settings that apply to all users. The settings in this dialog can only be changed by users with the appropriate authorizations.

Global settings include settings for:

- **System**
- **User management**: Create user accounts and set assigning rights.
- **Analysis and resources behavior**
- **Physical properties** include **Temperature unit** and **Barometric pressure unit**.
- **Operation mode**

5.4.1 System

Instrument identification, date and time format and current date and time can be set in this dialog.

Navigation: Home > Setup > Global settings > System

Parameter	Description	Values
Identification	Information on identification and software versions of the instrument.	-
Date/Time	Settings for date and time.	-
Header and footer	Settings for header, footer, and signature lines of print-outs.	-

5.4.1.1 Identification

You can enter your own code to identify the instrument.

Navigation: Home > Setup > Global settings > System > Identification

Parameter	Description	Values
Instrument ID	Define the instrument identification.	-
Firmware version	Information on the firmware version of the instrument.	-
Instrument SN	Information on the serial number of the instrument.	-
Module A	Information on the type of module currently plugged in.	pH/mV pH/Ion Conductivity DO/BOD
Module ID	Define the module identification of the module currently plugged in.	-
Module serial number	Information on the serial number of the module currently plugged in.	-
Module firmware version	Information on the firmware of the module currently plugged in.	-
Module B	Information on the type of module currently plugged in.	pH/mV pH/Ion Conductivity DO/BOD
Module ID	Define the module identification of the module currently plugged in.	-
Module serial number	Information on the serial number of the module currently plugged in.	-
Module firmware version	Information on the firmware of the module currently plugged in.	-
Module C	Information on the type of module currently plugged in.	pH/mV pH/Ion Conductivity DO/BOD
Module ID	Define the module identification of the module currently plugged in.	-
Module serial number	Information on the serial number of the module currently plugged in.	-
Module firmware version	Information on the firmware of the module currently plugged in.	-
Radio clock firmware version	Information on the firmware of the radio controlled clock.	-

- 1 Enter the ID in **Instrument ID**.
- 2 To enter an ID for **Module ID**, connect a module to the instrument.
- 3 To finish and store the entries tap [**Save**].

5.4.1.2 Date/Time

You can enter a display format for date and time and set the current date and time.

Navigation: Home > Setup > Global settings > System > Date/Time

Parameter	Description	Values
Date format	Defines the format for displaying the date.	dd/mm/yyyy mm/dd/yyyy yy/dd/mm yyyy/mm/dd dd-mmm-yy mmm-dd-yy yy-dd-mmm yy-mmm-dd mmm-dd-yyyy dd-mmm-yy yyyy-dd-mmm yyyy-mmm-dd
Time format	Defines the format for displaying the time.	24h a.m./p.m.
Enable radio clock	Activate the radio clock for time synchronization. Only available on instrument hardware with radio clock.	Active Inactive
Sender	Define the transmitter for the reception. Only if Enable radio clock is activated.	Auto List of transmitter
Date/Time	Enter the current date and time. Only if Enable radio clock is not activated.	-
Time zone	Define the time zone of your location.	UTC+ UTC-
Daylight saving time	Activates to switch automatically to daylight saving time.	Active Inactive

- 1 Select a format in **Date format**.
- 2 Select a format in **Time format**.
- 3 To use the radio clock, activate **Enable radio clock**.
- or -
Tap list field [**Date/Time**], enter the current day and time and confirm with [**OK**].
- 4 Select your location in **Time zone**.
- 5 To switch automatically to daylight saving time, activate **Daylight saving time**.
- 6 To finish and store the entries tap [**Save**].

Note

- Some instrument versions are equipped with a radio clock. The radio clock will attempt to synchronize every 7 hours. To avoid interruption of measurements during synchronization of the radio clock, the synchronization will only occur when the instrument is switched off. For reduction of interference during synchronization of the radio clock it is recommended to disconnect the power supply overnight or during long periods of non-use.
- Automated Daylight saving time and the radio clock should not be used at the same time.

5.4.1.3 Header and footer

You can define the header, footer, and signature lines of print-outs.

Navigation: Home > Setup > Global settings > System > Header and footer

Parameter	Description	Values
Header text	Defines a text which is shown in the header lines, e.g. company name or instrument location.	-
End of report	If activated, signature lines are added at the end of the report.	Active Inactive
Created by	If activated, a signature line Created by will be printed if End of report is activated.	Active Inactive
Modified by	If activated, a signature line Modified by will be printed if End of report is activated.	Active Inactive

Checked by	If activated, a signature line Checked by will be printed if End of report is activated	Active Inactive
Approved by	If activated, a signature line Approved by will be printed if End of report is activated.	Active Inactive
Header/footer on compact print-out	If activated, header and footer lines are printed on: - RS232 Compact Printer - USB Compact Printer . Header and footer lines are always printed on: - PDF Writer and - Network printer .	Active Inactive

- 1 Enter a text in **Header text**.
- 2 To add signature lines at the end of a report, activate **End of report** and activate which signature lines to print
- 3 To have the header and footer printed on a compact printer, activate **Header/footer on compact print-out**.
- 4 To finish and store the entries tap [**Save**].

5.4.2 User management

Navigation: Home > Setup > Global settings > User management

In the dialog **User management** you can administer the settings for **User** and **Account policies**. A maximum of 30 different users can be defined, whereby only one can be logged in (1 user operation). User accounts can be deleted and edited. There is a default user with a predefined user name (user name: **Administrator**, user group: Administrator) that cannot be deleted.

5.4.2.1 User

In this dialog you have access to the list of users. You can enter user names and assign the users to a user group. The user groups have different user rights, see [User groups and user rights ► Page 42]. User accounts can be locked, edited or deleted.

User list

All registered users are listed in the user list. The user list gives a quick overview on the users and their assigned user groups.

Navigation: Home > Setup > Global settings > User management > User

Enter user and assign user rights

New user can be added at any time. A maximum of 30 user can be added to the list. If the maximum is reached, delete a user to create a new entry.

The number of parameters that are displayed depend on the settings in **Account policies**.

Navigation: Home > Setup > Global settings > User management > User > New

Parameter	Description	Values
User name	Defines the user name which uniquely identifies him to the system.	-
Full name	Defines the full name of the user.	-
User group	Assignment of the user to a user group. Depending on the user group, the user has various rights.	Administrator Expert Technician Operator
Reset password	If activated, the password for the user will be reset to 123456 . Displayed if Enforce password/fingerprint is activated.	Active Inactive

Lock user	If activated, the user will be locked and cannot access the instrument. Displayed if Enforce password/fingerprint is activated.	Active Inactive
Enforce password change	If activated, the entry of a new password will be enforced with the next login. Displayed if Enforce password/fingerprint is activated.	Active Inactive
Created by	Information on the administrator logged in at time of account creation.	-
Created on	Information on date and time of account creation.	-
Modified by	Information on the administrator logged in at time of account modification.	-
Modified on	Information on date and time of account modification.	-

- 1 Enter a user name in **User name**.
- 2 Enter the full name of the user in **Full name**.
- 3 Assign a group in **User group**.
- 4 To lock a user, activate **Lock user**.
- 5 To finish and store entries tap [**Save**].

Note

Reset password and **Enforce password change** are displayed when the option **Enforce password/fingerprint** is activated. No changes are possible in this submenu. To edit the settings, see [Account policies ▶ Page 40].

Edit user

Navigation: Home > Setup > Global settings > User management > User

The settings for access permissions and user data can be changed.

- 1 Select a user from the list **User**.
➔ **User parameters** is displayed.
- 2 Edit the entries.
- 3 To finish and store entries tap [**Save**].

Delete user

Navigation: Home > Setup > Global settings > User management > User

User entries can be deleted.

- 1 Select a user from the list **User**.
➔ **User parameters** is displayed.
- 2 Tap [**Delete**].
➔ The user has been deleted without warning.

5.4.2.2 Account policies

In the dialog **Account policies** you can set the options for password or fingerprint reader, also see [Peripherals ▶ Page 30].

Navigation: Home > Setup > Global settings > User management > Account policies

Parameter	Description	Values
Enforce password/fingerprint	Activates to log in by entering a password or via the fingerprint reader.	Active Inactive
Min. password length	Defines the minimum length of the user passwords Displayed if Enforce password/fingerprint is activated.	0...20

Password complexity required	Enforces the user to define a password that must contain number, Active Inactive capital letters and small letters. Displayed if Enforce password/fingerprint is activated.
------------------------------	---

- 1 To activate the login screen when turning on the instrument, activate **Enforce password/fingerprint**. In **Min. password length**, enter the number of digits required.
- 2 To finish and store the entries tap [**Save**].

5.4.2.3 User management settings and login screen

Depending on the settings in **User** and **Account policies**, the following login options are available:

Enforce password/fingerprint is deactivated

- If only one user is defined, the instrument starts directly with the home screen without showing the login screen.
- If several users are defined, the instrument starts with the login screen, the user can select an entry from the list. The input of a password is not required.

Enforce password/fingerprint is deactivated, fingerprint reader is connected and activated

- If only one user is defined, the instrument starts directly with the home screen without showing the fingerprint logon screen.
- If several users are defined, the instrument starts with the fingerprint logon screen. The identification by fingerprint leads to the home screen. In case of problems with the fingerprint identification the user can change to the manual login screen and selects an entry from the list. The input of a password is not required.

Enforce password/fingerprint is activated, fingerprint reader is not activated

- If the user logs in for the first time or **Enforce password change** is deactivated, the instrument starts with the login screen, enter username and password.
- If the user logs in for the first time or **Enforce password change** is active, the instrument starts with the login screen, enter username and password. Change password screen is displayed, enter a new password, repeat the password and confirm with [**OK**].

Note

- Login resets the parameter in **Enforce password change**.

Enforce password/fingerprint is activated, fingerprint reader is activated

- If **Enforce password change** is deactivated, the instrument starts with the fingerprint logon screen. The identification by fingerprint leads to the home screen. In case of problems with the fingerprint logon the user can change to the manual login screen to enter user name and password.
- If user logs in for the first time or **Enforce password change** is activated, the instrument starts with the fingerprint logon screen. The identification by fingerprint leads to the Change password screen, enter new password, repeat the password and confirm with [**OK**].
- If user logs in for the first time or **Enforce password change** is activated, the instrument starts with the fingerprint logon screen. The user changes to the password login and enters user name and password. Change password screen is displayed, enter a new password, repeat the password and confirm with [**OK**].

5.4.2.4 User groups and user rights

Every user is assigned to a user group.

The following table shows the user rights that are assigned to the corresponding user group:

User right	Operator	Technician	Expert	Administrator
Performance of analysis				
Start method or series via shortcut	•	•	•	•
Start method from method list		•	•	•
Start series from series list		•	•	•
Start direct measurement from home screen	•	•	•	•
Start direct measurement from module screen	•	•	•	•
Start direct calibration from module screen	•	•	•	•
Review results	•	•	•	•
Delete results			•	•
Modification of analysis parameters				
Change sample ID in start analysis screen or module screen	•	•	•	•
Create, edit and delete methods			•	•
Create, edit and delete series			•	•
Import / export of methods and series			•	•
Change module settings		•	•	•
User interface settings				
Modify user-specific screen language	•	•	•	•
Modify user-specific notification signal (beep)	•	•	•	•
Modify user-specific keyboard layout	•	•	•	•
Modify user-specific screen brightness	•	•	•	•
Create, edit and delete shortcuts		•	•	•
Instrument settings				
Review sensor data and calibration history (last 5 calibrations)	•	•	•	•
Create, edit and delete sensors			•	•
Manually modify cell constant of conductivity sensors			•	•
Create, edit and delete buffer and standard solution lists			•	•
Import / export of sensors, buffers, tables			•	•
Modify peripherals settings			•	•
Modify analysis sequence settings			•	•
Modify actions when sensor expire			•	•
Modify physical properties (e.g. temperature unit)				•
Modify operation mode			•	•
Modify report language			•	•
Edit system settings (instrument ID, date/time, header and footer)				•
Create, edit and delete users				•
Reset to factory settings				•
Update instrument and module firmware				•

5.4.3 Analysis and resources behavior

Navigation: Home > Setup > Global settings > Analysis and resources behavior

The following settings can be defined:

- **Analysis sequence settings**
- **Actions when sensors expire**

5.4.3.1 Analysis sequence settings

With these settings the instrument shows all required resources at the start of a method or series. Additionally you can specify whether the analysis ends automatically or manually.

Navigation: Home > Setup > Global settings > Analysis and resources behavior > Analysis sequence settings

Parameter	Description	Values
Show required resources at start	All resources that are required to perform the analysis are shown.	Active Inactive
Allow taking of manual endpoint	A button Take manual endpoint is shown on the online display to finish the measurement manually. If enabled this button is shown independent of the defined endpoint type. If disabled this button is only shown if Endpoint type = Manual .	Active Inactive
Confirm end of analysis	Instrument waits for confirmation after every measurement.	Active Inactive
Show calculated results after analysis	Before the end of a method, all calculated results are shown and have to be confirmed.	Active Inactive
Suppress limits-/threshold popupsConfirm end of analysis	Hinders the instrument to show messages if limits or thresholds are exceeded.	Active Inactive
Save results on instrument	Results of up to 250 analysis are stored in the instrument. Please note that the method step Follow for analysis type BOD and BOD cannot be started if activated.	Active Inactive

5.4.3.2 Actions when sensors expire

With these settings you specify whether the instrument warns, warns and blocks or does not react, when usable life and life span expired.

Navigation: Home > Setup > Global settings > Analysis and resources behavior > Actions when sensors expire

Parameter	Description	Values
Usable life	Select the action when exceeding usable life.	Warn Warn and block None
Life span	Select the action when exceeding life span.	Warn Warn and block None

- 1 To be warned and/or blocked, when **Usable life** expires, activate **Usable life**.
- 2 To be warned and/or blocked, when **Life span** expires, activate **Life span**.
- 3 To finish and store the entries tap [**Save**].

5.4.4 Physical properties

In **Physical properties** dialog you can define the parameters for **Temperature unit**. The setting will be applied after restarting the instrument.

Navigation: Home > Setup > Global settings > Physical properties

Parameter	Description	Values
Temperature unit	Defines the temperature unit applicable for all measurements. All entries and representations will be provided in the selected unit.	°C °F

Barometric pressure unit	Defines the barometric pressure unit applicable for dissolved oxygen measurements. All entries and representations will be provided in the selected unit.	mbar hPa mmHg atm
--------------------------	--	-------------------------

5.4.5 Operation mode

You can define if direct measurements and direct calibrations can be started beside methods and series. The settings defined for the module are used. If direct measurements are not allowed, the **Read** button in the home screen and the **Read** and **Calibrate** button in the module settings are hidden.

Navigation: Home > Setup > Global settings > Operation mode

Parameter	Description	Values
Allow direct measurements	If enabled, starting direct measurements and calibrations is possible.	Active Inactive

See also

📖 Starting direct measurement ▶ Page 53

5.5 Tables

Navigation: Home > Setup > [Tables]

This instrument offers the comfort to enter data and calculate results, based on these data, see [Formula Syntax ▶ Page 117] and [Results ▶ Page 114]. Tables can be used in calculations.

Tables are always assigned to an application. The list can be sorted according to the application.

There are two types of tables:

- **METTLER TOLEDO tables:**
These tables are included in the factory settings; they are listed and neither can be modified nor deleted.
- **User-defined tables:**
The layout of the table is a value table (x-y). It can be created, edited and deleted. User-defined tables are added to the list of METTLER TOLEDO tables.

5.5.1 METTLER TOLEDO tables

Navigation: Home > Setup > Tables

In **Tables** you find a list of predefined buffer and standard tables. User-defined buffer and standard tables can be added.

The following buffer and standard tables with their temperature related values are available:

- METTLER TOLEDO Europe pH 2.00
- METTLER TOLEDO Europe pH 4.01
- METTLER TOLEDO Europe pH 7.00
- METTLER TOLEDO Europe pH 9.21
- METTLER TOLEDO Europe pH 11.00
- **Ultra pure water** (Ultra-pure water (USP/EP))
- **Purified water** (Purified water (EP))
- Conductivity Standard 10 µS/cm
- Conductivity Standard 84 µS/cm
- Conductivity Standard 1413 µS/cm
- Conductivity Standard 12.88 mS/cm
- Oxygen solubility in water
- USP645 Stage 3 pH and conductivity requirements

For detailed information concerning the temperature related values:

- Tap a table from the list.
- ➔ **Tables** is displayed and shows the selected pH value related to the temperature, starting at 5 °C to 50 °C in steps of 5 °C.
The conductivity values for ultra-pure water and purified water are displayed in $\mu\text{S}/\text{cm}$.
The table with the temperature related values of ultra-pure water contains all data between 0 °C and 100 °C in steps of 5 °C.
The table with the temperature related values of purified water contains all data between 0 °C and 100 °C in steps of 10 °C.

For more information about the conditions for the calculation of the values:

- Tap **Parameter**.
 - ➔ **Parameters** is displayed and provides information on:
 - **Name**
 - **Input value**
 - **Output value**
 - **Fit type**
 - **Comment**

See also

 METTLER TOLEDO tables ▶ Page 139

5.5.2 User-defined tables

Navigation: Home > Setup > Tables

The following section describes how to create and administer a user-defined table, how to enter, modify and delete values and how to delete tables.

Note

- A total of 10 user-defined tables can be created. If the limit is exceeded delete a table before you create a new one.

Creating tables

Navigation: Home > Setup > Tables > New

Parameter	Description	Values
Name	Enter a name that uniquely identifies the table in the system. The name of the table can be entered in the method function Calculation > Formula .	-
Input value	Defines the table heading for the input value.	-
Output value	Defines the table heading for the output value.	-
Fit type	Information on the definition of the curve type for the calculation. Linear interpolation corresponds to a segmented curve.	-
Comment	Defines a short comment that will be displayed.	-

- 1 Enter a name in **Name**. The new table name must be unique.
 - 2 Enter a table heading value for **Input value**.
 - 3 Enter a table heading value for **Output value**.
 - 4 Enter a comment in **Comment**.
 - 5 Tap [**Save**].
- ➔ **Table** is displayed.

Entering values in tables

- Tap [**New**].

➔ **Values** is displayed.

Parameter	Description	Values
x	Defines the input value of the value pair.	-
y	Defines the output value of the value pair.	-

- 1 Enter value for **x**.
- 2 Enter value for **y**.
- 3 Confirm with [**OK**].
- 4 To enter additional values, tap [**New**] and repeat previous steps.
- 5 To finish and store the entries tap [**Save**].

Note

The maximum number of value entries per table is 25. If you enter by mistake identical values for **Input value** and **Output value**, a prompt informs that his value already exists.

Modifying table parameters

Navigation: Home > Setup > Tables

- 1 Select a table from **Tables**.
➔ **Table** is displayed.
- 2 Tap [**Parameter**].
➔ **Parameters** is displayed.
- 3 To Following parameters can be modified:
 - **Name**
 - **Input value**
 - **Output value**
 - **Comment**
- 4 To finish and store the entries tap [**Save**].

Modifying values

Navigation: Home > Setup > Tables

- 1 Select a table from **Tables**.
➔ **Table** is displayed.
- 2 Tap a value from the list.
➔ **Values** is displayed.
- 3 To modify the value, overwrite it and confirm with [**OK**].
- 4 To modify additional values, repeat previous steps.
- 5 To finish and store the entries tap [**Save**].

Deleting values

Navigation: Home > Setup > Tables

- 1 Select a table from **Tables**.
➔ **Table** is displayed.
- 2 Tap the value on the list.
➔ **Values** is displayed.
- 3 Tap [**Delete**].
- 4 To delete additional values, repeat previous steps.
- 5 To finish and store the entries tap [**Save**].
➔ The values have been deleted.

Deleting tables

Navigation: Home > Setup > Tables

- 1 Select a table from **Tables**.
 - ➔ **Table** is displayed.
- 2 Tap [**Parameter**].
 - ➔ **Parameters** is displayed.
- 3 Tap [**Delete**].
 - ➔ The table has been deleted.

Note

Calculations which refer to deleted tables can no longer be carried out.

5.6 Maintenance & service

Navigation: Home > Setup > Mainten. & Service

The following functions are available:

- **MT-Service**
- **Import / Export**
- **Reset to factory settings**
- **Firmware**
- **Update**

5.6.1 MT service

Navigation: Home > Setup > Mainten. & Service > MT-Service

In this dialog a list is shown with (max. 10 entries) METTLER TOLEDO maintenance information. Each entry contains information about the date of last maintenance and the executive service technician. The most recently performed maintenance appears at the bottom of the list.

Set service life

Navigation: Home > Setup > Mainten. & Service > MT-Service > Settings

Define the parameter for the service life to determine the date of the next maintenance interval. You can set a reminder to be informed on the upcoming maintenance. The following parameters can be set:

- **Service life:** Set the maintenance intervals (in days).
- **Reminder:** Activates whether a reminder informs about the days before the next maintenance takes place.
- **Days before expiration:** Define the number of days before the reminder informs about expiration of the maintenance interval.
The entered number of days must be smaller than the number of days you set for the service life (displayed if **Reminder** is activated).

5.6.2 Import/export

Navigation: Home > Setup > Mainten. & Service > Import / Export

With the aid of this function you can either transfer data to or from USB sticks. Using a **Memory copy** is an easy way to transfer the settings from one instrument to another. Quickly you have installed two instruments with identical functions.

The following data can be imported or exported:

- **Methods**
- **Series**
- **Buffers & Standards**
- **Sensors**

- **Memory copy**
- **User management**
- **Log file**

Note

- The function **Import / Export** is only possible when a USB stick is connected and when no task is running.
 - All the saved results, analysis data, shortcuts and METTLER TOLEDO tables, METTLER TOLEDO methods as well as predefined buffer sets and standard groups are **not** contained in a memory copy.
 - When you import or export user management settings, the entire user management settings, with all users and their properties, will be imported or exported.
 - You need to have administrator rights to create and re-import a backup copy.
- 1 In **Action**, select **Import** or **Export**.
 - 2 Tap [**Data**] and select an item from the list.
 - 3 To transfer data tap [**Start**].

5.6.3 Reset to factory settings

Navigation: Home > Setup > Mainten. & Service > Reset to factory settings

You can reset all instrument settings to factory settings. A pop-up informs you that all existing data and settings will be deleted.

Note

If reset is performed all created data, amendments, settings, setup entries and results will be lost. Please create a **Memory copy** before you reset to factory settings.

- 1 Tap [**Reset to factory settings**].
 - ➔ An information is displayed.
- 2 Tap [**Start**] to proceed.
 - ➔ A second information is displayed.
- 3 Tap [**Continue**] to start the procedure.
 - ➔ The instrument is shutting down and all the data will be deleted.

5.6.4 Firmware

Navigation: Home > Setup > Mainten. & Service > Firmware

The list shows all the firmware updates. The first entry on the list corresponds to initial operation.

All list entries show the following information:

- **Date:** The date of the installation.
- **FW Version:** The software version installed.
- **User name:** Name of technician.

5.6.5 Update

Navigation: Home > Setup > Mainten. & Service > Update

Firmware updates of the instrument and modules can be managed from this screen. Data can be transferred from USB stick.

Note

- You need administrator rights to implement a firmware update.
- METTLER TOLEDO provides additional information for updating instruments. These instructions will be part of the firmware update.

6 Sensors

Navigation: Home > Sensors

Additional to the sensor management in the setup menu, SevenExcellence™ allows managing the sensors and sensor settings directly from the home screen.

6.1 Available sensors

Navigation: Home > Sensors

Sensors leads you to the list of available sensors according to the sensor settings, see [Sensors ▶ Page 27]. You have a quick view over the defined sensors and the modules to which they are connected. Furthermore you have access to the calibration history.

6.2 Sensor parameters

Navigation: Home > Sensors

To get detailed information concerning the sensor parameters, to change the sensor name and the conditions for usable life and life span, proceed as follows:

- 1 Select a sensor.
 - ➔ **Parameters** is displayed with all information concerning sensor settings.
- 2 Define the parameters according to your requirements.
- 3 To finish and store the entries, tap [**Save**].

6.3 Calibration history

To get reliable and exact measurement values, check the calibration history before starting a measurement. Calibrate the sensor or replace it if necessary. You can print the calibration history for your documentation.

- 1 Tap [**Sensors**].
 - ➔ **Sensors** is displayed.
- 2 Select a sensor.
 - ➔ **Parameters** is displayed.
- 3 Tap **Calibration history**.
 - ➔ **Calibration history** is displayed.
- 4 If desired, tap [**Print**].
 - ➔ The calibration history is being printed.

7 Starting Analysis

An analysis can be started in various ways:

- Start by tapping [**Read**] on the **Home screen**
- Start by tapping [**Read**] from **Module settings**
- Start by tapping [**Calibrate**] from **Module settings**
- Start by tapping [**Shortcuts**] on the **Home screen**
- Start from the Methods/Series editor

The different options that are offered depend on the settings you have done before and on the kind of analysis. The following chapters show how to start:

- **Direct calibration**
- **Direct measurement**
- **Methods**
- **Series**

Note

- Before you start a calibration or measurement using a METTLER TOLEDO method, connect a module to the instrument and connect a sensor to the module suitable to the method.
- Add the sensor to the **Sensors** list, see [Sensors ▶ Page 27].
ISM® sensors are recognized automatically and enable to start the method immediately.

7.1 Starting direct calibration

You can start a direct calibration from the module displayed on the home screen. Before you start a direct calibration you have to set the following parameters:

- **Sensor**
- **Temperature sensor** (optional)
- **General settings**
- **Calibration settings**

The following section guides you through this process.

- Tap the colored tag of the relevant module displayed on the top of the **Home screen**.
➔ **Module settings** is displayed.

Parameter	Description	Values
Sensor name	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. A selected sensor is referenced by its name.	List of defined sensors
Temperature sensor	Opens the list of temperature sensors, according to the sensor settings. Displayed for pH/Ion , Conductivity , and DO/BOD module.	-
General settings	All parameters concerning endpoint and temperature for direct calibration and direct measurement.	-
Measurement settings	All parameters concerning the measurement and the result (only for direct measurement).	-
Calibration settings	All parameters concerning buffers or standards (only for direct calibration).	-

- 1 Tap list field **Sensor name**.
➔ **Select sensor** is displayed.
- 2 Select a type.
- 3 To use an external temperature sensor, tap list field **Temperature sensor**.

➔ **Select sensor** is displayed.

4 Select a type.

General settings

In **General settings** you set the conditions for **Endpoint type**, **Endpoint criteria** and **Temperature capture**. Additionally you can set whether the results of the calibration should be printed or not.

– Tap [**General settings**].

➔ **General settings** is displayed.

Parameter	Description	Values
Measurement type	Information on the measurement type.	-
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
Temperature capture	Select the temperature capture mode. Internal : The used temperature sensor is integrated in the measurement sensor. External : A separate temperature sensor is used. Manual : The temperature is entered manually.	Internal External Manual
Temperature	Defines the temperature for the measurement. Displayed if Temperature capture = Manual is selected.	-30...130 °C / -22...266 °F
Barometric pressure capture	Select the Barometric pressure capture mode. Displayed in case of Measurement type = Dissolved Oxygen	Automatic Manual
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed in case of Measurement type = Dissolved Oxygen and Barometric pressure capture = Manual .	500.0...1100.0 mbar
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
Print	Defines if data will be printed to the connected printer.	Active Inactive
Print format	Summary : Covers most important data concerning date, time, user and parameters according to the settings of the measurement type.	Summary

1 Select **Endpoint type**.

2 Select **Endpoint criteria**.

3 Select **Temperature capture**.

If **Manual** is selected, enter a value for **Temperature**.

4 If desired, activate check box **Print**.

5 Confirm with [**OK**].

Calibration settings

– Tap [**Calibration settings**].

➔ **Calibration settings** is displayed.

Parameter	Description	Values
Measurement type	Information on the measurement type.	-
Buffer group	Select a buffer group from the list of predefined and user-defined buffer groups. Displayed if Measurement type = pH .	-
Calibration standard group	Opens the list of predefined and user-defined standard groups. Displayed if Measurement type = Ion or Conductivity .	-
Unit	Information on the preselected unit. Displayed if Measurement type = Ion	-
Cal. mode	Specifies whether calibration data will be approximated with one straight line segment or left as a series of segments. Displayed if Measurement type = pH or Ion .	Linear Segmented
Number of buffers	Defines the number of buffers for the calibration. Displayed if Measurement type = pH .	1 2 3 4 5
Number of standards	Defines the number of calibration standards to be used. For Measurement type = Ion up to 5, for Conductivity up to 2.	1 2 3 4 5
Standard 1 - Standard 5	Number of fields depends on the number of defined standards. Up to 5 standards for ion sensors and up to 2 standards for conductivity sensors are displayed with consecutive numbers. Each field opens the list of preselected standard groups. Displayed if Measurement type = Ion or Conductivity .	-
Automatic buffer recognition	Activates the automatic recognition of buffers. Only available for predefined pH buffer groups. Displayed if Measurement type = pH .	Active Inactive
Buffer 1 - Buffer 5	The number of fields depends on the number of defined buffers. Up to 5 buffers are displayed with consecutive numbers. Select a buffer for each calibration point from the list. Displayed if Measurement type = pH and Automatic buffer recognition is not activated.	-
Calibration points	Option to choose between different number of calibration points. Displayed if Measurement type = Dissolved Oxygen .	1 2
Calibration standard 1	Select the standard for the first calibration point. Displayed if Measurement type = Dissolved Oxygen .	Saturated air
Calibration standard 2	Defines the standard for the second calibration point. Cannot be edited. Displayed if Measurement type = Dissolved Oxygen and Calibration points = 2 .	Zero point standard

- 1 Tap the list field **Buffer group** or **Calibration standard group**.
➔ **Buffer / Standard group** is displayed.
- 2 Select a buffer or standard from the list.
- 3 If **Measurement type** pH or ion is selected, select a **Cal. mode**.
- 4 Select the number of buffers you want to use from the list **Number of buffers**.
- or -
Select **Automatic buffer recognition**.
- 5 If **Number of buffers** is selected, select specific values for **Buffer**.
- or -
Select **Standard**.
- 6 Confirm with [**OK**].

- 7 Tap [**Save**].
- 8 Immerse the sensor into the first buffer.
- 9 Tap [**Calibrate**].

Note

- The calibration settings are saved even if you log out.

See also

 Endpoint criteria ▶ Page 110

7.2 Starting direct measurement

A direct measurement can be started in two ways. You can start the analysis via the **Module settings** or directly from the **Home screen** via the button [**Read**].

7.2.1 Starting from module settings

You can start a direct measurement from the module displayed on the **Home screen**. Before you start a direct measurement you have to set the following parameters:

- Select a sensor
- Select a temperature sensor (optional)
- Adjust the **General settings**
- Adjust the **Measurement settings**

The following section guides you through this process.

- Tap the colored tag of the relevant module displayed on the top of the **Home screen**.
 - ➔ **Module settings** is displayed.

Parameter	Description	Values
Sensor name	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. A selected sensor is referenced by its name.	List of defined sensors
Temperature sensor	Opens the list of temperature sensors, according to the sensor settings. Displayed for pH/Ion , Conductivity , and DO/BOD module.	-
General settings	All parameters concerning endpoint and temperature for direct calibration and direct measurement.	-
Measurement settings	All parameters concerning the measurement and the result (only for direct measurement).	-
Calibration settings	All parameters concerning buffers or standards (only for direct calibration).	-

- 1 Tap list field **Sensor name**.
 - ➔ **Select sensor** is displayed.
- 2 Select a type.
- 3 To use an external temperature sensor, tap list field **Temperature sensor**.
 - ➔ **Select sensor** is displayed.
- 4 Select a type.

General settings

- Tap **General settings**.
 - ➔ **General settings** is displayed.

Parameter	Description	Values
Measurement type	Information on the measurement type.	-
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
Temperature capture	Select the temperature capture mode. Internal: The used temperature sensor is integrated in the measurement sensor. External: A separate temperature sensor is used. Manual: The temperature is entered manually.	Internal External Manual
Temperature	Defines the temperature for the measurement. Displayed if Temperature capture = Manual is selected.	-30...130 °C / -22...266 °F
Barometric pressure capture	Select the Barometric pressure capture mode. Displayed in case of Measurement type = Dissolved Oxygen	Automatic Manual
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed in case of Measurement type = Dissolved Oxygen and Barometric pressure capture = Manual .	500.0...1100.0 mbar
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
Print	Defines if data will be printed to the connected printer.	Active Inactive
Print format	Summary: Covers most important data concerning date, time, user and parameters according to the settings of the measurement type.	Summary

- 1 Select **Endpoint type**.
- 2 Select **Endpoint criteria**.
- 3 Select **Temperature capture**.
If **Manual** is selected, enter a value for **Temperature**.
- 4 If desired, activate check box **Print**.
- 5 Confirm with [**OK**].

Measurement settings

- Tap **Measurement settings**.
➔ **Measurement settings** is displayed.

Parameter	Description	Values
Measurement type	Information on the measurement type.	-
Sample ID	Defines the sample ID.	-
Auto sequential sample ID	Activates auto sequential sample ID.	Active Inactive
Conductivity mode	Defines the mode of the method. Displayed if Measurement type = Conductivity is selected.	Conductivity TDS Salinity Resistivity

Unit	The unit to be used for the measurement. The displayed units depend on the selected measurement types.	pH mV μS/cm mS/cm S/m μS/m mS/m ppm ppt (‰) mg/L g/L psu ppt (‰) MΩ.cm Ω.cm mV Rel.mV mmol/L mol/L ppm % pX mV
Offset	Defines the offset [mV]. Displayed if Measurement type = Redox and Unit = Rel.mV is selected.	-2000...2000 mV
Temperature correction	Defines the relationship between conductivity, temperature and ion concentration. Linear: Use for the temperature correction of medium and highly conductive solutions. Non-linear: Use for natural water (only for temperature between 0...36 °C). The measured conductivity at the sample temperature is corrected to the defined reference temperature (20 °C or 25 °C). Off: The conductivity value at the current temperature is displayed. Pure water: An optimized type of temperature algorithm is used for Measurement type = Conductivity or Resistivity . Displayed if Measurement type = Conductivity, TDS or Resistivity is selected.	Linear Non-linear Off Pure water
α-coefficient	Defines the factor for the linear dependency. Displayed if Cond.mode = Conductivity, TDS or Resistivity in combination with Temperature correction = Linear is selected.	0.00...10.00
Reference temperature	The conductivity reading will be directly corrected to the set reference temperature. If Measurement type = Conductivity or Resistivity in combination with Temperature correction = Pure water is selected a reference temperature of 25 °C is automatically set. Displayed if Measurement type = Conductivity, TDS or Resistivity in combination with Temperature correction = Linear is selected.	20°C 25°C
TDS factor	To calculate the TDS value, the conductivity will be multiplied with this factor. Displayed if Conductivity = TDS is selected.	0.00...10.00
Decimal places	Defines the number of digits for the displayed measurement result. The displayed decimal places depend on the selected unit.	1 2 3 4
Interval	Activates the interval measurement. The measurement data will be stored after each period specified in Interval time .	Active Inactive
Interval time	Defines the time period between successive readings [s].	1...10 ⁶
Print after each interval	Activates printing out the result after each interval. Only if method type = Interval .	Active Inactive
Limits	Activates to define the limits for the measurement.	Active Inactive
Lower limit	Defines a value for the deviation of the measured value downward.	-
Upper limit	Defines a value for the deviation of the measured value upward.	-

- 1 Enter a sample ID.
- 2 If desired, activate **Auto sequential sample ID**.
- 3 Select a unit.
- 4 Select the resolution decimals to be displayed.
- 5 Enter measurement parameters according to your measurement type and requirements.
- 6 If desired:
 - Activate **Interval** and enter a value for **Interval time**.
 - Activate **Limits** and enter values for **Lower limit** and **Upper limit**.
- 7 Confirm with [**OK**].
- 8 To finish and store entries, tap [**Save**].
- 9 Immerse the sensor into the sample.
- 10 Tap [**Read**].

See also

 Endpoint criteria ▶ Page 110

7.2.2 Starting from home screen

A measurement can also be started directly from the home screen. In this case the last saved settings of the module are used for the measurement.

Note

- Verify that the sensor according to the module settings is connected.

- 1 Activate the check box of the module.
- 2 Tap the colored tag of the module.
- 3 Immerse the sensor into the sample.
- 4 Tap [**Read**].

Note

- If stirring is enabled in the settings of one or several selected modules, the stirrer is activated and uses the lowest of all stirring speeds.

7.3 Starting methods/series

Methods and **Series** can be started in different ways.

- Start directly after creating a method or series.
- Start from the method an series list.
- Start from the **Home screen** via shortcut or direct shortcut.

7.3.1 Start directly after creating methods/meries

When you have created a method or series and saved all settings you can immediately start the analysis.

- 1 Tap [**Start**].
 - ➔ The start analysis dialog is displayed.
- 2 Tap [**Start**].

7.3.2 Starting from methods/series list

You can start predefined and user-defined methods and series directly from the method or series lists.

- 1 Tap [**Methods**] or [**Series**].
 - ➔ **Methods** or **Series list** is shown.
- 2 Select the method or series you want to run.
- 3 Tap [**Start**].

- **Start analysis** dialog is displayed.
- 4 Tap [**Start**].

7.3.3 Starting from shortcut/direct shortcut

You can start methods and series via shortcuts and direct shortcuts displayed on the home screen.

- 1 Tap the shortcut on the home screen.
 - **Start analysis** is displayed.
- 2 Tap [**Start**].

Note

- Tapping a direct shortcut triggers an immediate start of an analysis. The on-line display is shown forthwith.

7.4 Interrupting analysis

Ongoing analyses can be interrupted in the on-line display via [**Terminate**]. To continue the task tap [**Cancel**].

7.5 Errors in the analysis sequences

Analyses may show errors. Four types of errors can be identified during an analysis sequence:

- **Error**
- **Terminate Error**
- **Critical error**
- **Failed**

7.5.1 Malfunction Types: Error

Malfunctions of the **Error** type are:

- Interrupting the analysis.

The malfunction **Error** triggers the following behavior:

- A message appears with information and hints on the error
- All remaining method functions will not be executed
- The processing of further tasks is interrupted
- The analysis with the status Error is listed in the analysis list.

See also

 Results ▶ Page 114

7.5.2 Malfunction types: Terminate error

Malfunctions of the type **Terminate Error** are:

- Disconnected ISM[®] sensor
- Connecting ISM[®] sensor during measurement

The malfunction **Terminate Error** triggers the following response:

- An information message for the relevant error appears
- The task is terminated immediately
- The processing of further tasks is interrupted
- No results are generated

7.5.3 Malfunction types: Critical error

Malfunctions of the type **Critical error** are:

- Disconnected module

- Connecting the module during measurement

The malfunction **Critical error** triggers the following response:

- A message appears with information on the error
- The task is terminated immediately
- The processing of further tasks is interrupted
- The instrument is shut down automatically to prevent further damage
- No results are generated
- After the restart of the instrument a popup is shown, indicating which critical error has occurred

7.5.4 Malfunction types: Failed

Malfunctions of the type **Failed** are:

- Calibration analysis failed
- Method function **Sensor check** failed
- Method function **Calculation** failed

The malfunction **Failed** triggers the following response:

- A message appears with information on the error
- The task is terminated immediately
- The processing of further tasks is interrupted
- All remaining method functions will not be executed, except **Report**
- An entry is generated in **Results** with information according to the failure

8 Methods

Navigation: Home > Methods

To perform a calibration, sensor test or measurement you can create methods. Methods can also be created for interval measurements or incremental ion measurements. A method is an analysis program and consists of a sequence of method functions. Method functions encompass the type of method, its configuration, the selection of expansion units, the definition of measurement criteria desired and ultimately the creation of a report for the measurements recorded. Parameters of method functions can be defined depending on the requirements of the method.

8.1 Method types

SevenExcellence offers the following seven method types.

Method type	Description
Measurement	This is the method type for a normal measurement. Up to three (parallel) measurement types can be defined in one method. Multiple measurements and several calculations can be defined.
Calibration	In this method type, the number of buffers or standards to perform a calibration is defined. The method function Measure (Calibration) repeats for the defined number of buffers/standards. The calibration data of the sensor is calculated automatically and will be transferred to the sensor setup.
Interval	This method type is very similar to method type Measurement . The difference is that data points during the measurement are recorded and stored in a value table or printed out during execution.
Incremental	Incremental methods are exclusively for the measurement type Ion. Only one module can be used. The method function Measure (Incremental) repeats for the defined number of samples/standard additions all other method functions are not repeated. The concentration of the sample is calculated automatically.
Sensor Test	This method type is related to method type Calibration and gives information about the condition of the sensor. The sensor calibration data is not transferred to the sensor setup, it serves only informative purposes.
BOD	Method type BOD (Biological Oxygen Demand) is only for the measurement type Dissolved Oxygen . Only one module can be used. The method function Measure (BOD) repeats for the defined number of bottles; all other method functions are not repeated. The BOD value of the sample is calculated automatically. The same method is used for the DO measurement before and after the incubation time. The measurement at start is called Base , the measurement afterwards (e.g. 5 days for BOD5) is called Follow . With one Base measurement several measurements of type Follow can be done (e.g. after 5 days and after 10 days).
BCV	Method type BCV (BOD Check Values) is only for the measurement type Dissolved Oxygen. Only one module can be used. The method functions Measure (Blank) , Measure (Seeded blank) and Measure (Standard) are repeated for the defined number of bottles in method types Blank (BOD) , Seeded blank (BOD) and Standard (BOD) ; all other method functions are not repeated. A BCV method determines the BOD of different solutions needed for measurement correction or validation of the application. Every BCV method supports a BOD method; without relation to a BOD method it is useless. At BOD method start the relation between BCV and BOD method is built. The same BCV method is used for the DO measurement before and after the incubation time of the check solutions. The measurement at start is called Base , the measurement afterwards (e.g. 5 days for BOD5) is called Follow . With one Base measurement several measurements of type Follow can be done (e.g. after 5 days and after 10 days). It is needed to run the BCV method always before running the related BOD method.

8.2 METTLER TOLEDO methods

Navigation: Home > Methods

A wide range of standard calibration and measurement methods, used in daily laboratory practice, are covered by predefined METTLER TOLEDO methods. This way the instrument allows the immediate use of methods to carry out calibrations and measurements or to create a series in combination with a METTLER TOLEDO method.

Predefined METTLER TOLEDO methods

The following METTLER TOLEDO methods for different applications are available:

Method	Used for
M001	pH calibration with METTLER TOLEDO Europe buffer group
M001-US	pH calibration with METTLER TOLEDO USA buffer group
M002	conductivity calibration
M003	ion calibration
M004	pH measurement
M005	pH measurement of pure water according to USP/EP/Ch.P.
M006	conductivity measurement
M007	conductivity measurement with α -coefficient determination
M008	USP645 Bulk Water Stage 1
M011	conductivity of bioethanol
M012	Rel.mV measurement with offset determination
M013	incremental measurement with single standard addition
M014	incremental measurement with multiple standard additions
M015	incremental measurement with single sample addition
M016	interval pH measurement
M017	sensor test of pH sensor
M018	DO calibration
M019	DO measurement
M020	biochemical oxygen demand (BOD) analysis
M021	BOD check values (BCV) analysis
M027	pH Calibration in pH 4.01 and 7.00
M028	USP645 sterile water analysis for <10 mL volume
M029	USP645 sterile water analysis for >10 mL volume
M030	USP645 bulk water stage 2
M031	USP645 bulk water stage 3
M032	advanced sensor test for pH sensor
M033	pH calibration verification in pH 4.01, 7.00, and 9.21
M034	conductivity calibration verification in 84 μ S/cm
M035	Purified Water EP/ChP
M036	Water Injection EP/ChP Stage 1
M037	Water Injection EP/ChP Stage 2
M038	Water Injection EP/ChP Stage 3

8.3 Creating methods

Various ways of generating a method are available to you in the instrument. Either you create a new method by changing the parameters of a delivered METTLER TOLEDO method and saving it under a new method ID or you select an appropriate method template from the list of proposals, modify it and save it under a new method ID. Lastly, you can also create a method by copying and saving an existing method.

A maximum number of 50 user-defined methods can be stored.

8.3.1 Copying an existing method

Navigation: Home > Methods

The easiest way to create a new method, is copying an existing method by changing the method ID.

- 1 Select a method with the parameters that are nearest to the calibration or measurement you wish to carry out.
 - ➔ The method ID of the method is displayed.
- 2 Tap **[Title]**.
- 3 Enter a new ID in **Method ID** and confirm with **[OK]**.
- 4 Tap **[OK]** to accept the entries.
 - ➔ A message box appears with the information that the new method has been created as copy of the preselected method.
- 5 Confirm the message box with **[OK]**.
- 6 Adapt the method according to your requirements, e.g. insert additional method functions and adjust method parameters.
- 7 Tap **[Save]** to save your entries.
 - ➔ The new method has been saved.

8.3.2 Creating methods based on templates

Using method templates is another way to create a method. The parameters can be modified according to your requirements.

Navigation: Home > Methods

Type	ID	Title	Description
MS	T0001	Measure	Method type for normal measurements.
CAL	T0002	Calibration	Method type for calibrating sensors.
INT	T0003	Interval	Method type that record data points during the measurements.
INC	T0004	Incremental	Incremental methods for the measurement type Ion.
ST	T0005	Sensor Test	Method type that provide information about the condition of sensors.
BOD	T0007	BOD	Method type for the measurement type Dissolved Oxygen.
BCV	T0006	BOD Check Values	BCV methods determine the BOD of different solutions needed for measurement correction or validation of applications.

- 1 In **Methods**, tap **[New]** and select a template from the list.
 - ➔ The method function **Configuration** of the new method is displayed.
- 2 Define the measurement type(s) and tap **[OK]** to proceed.
Note: Once confirmed, the measurement type(s) cannot be modified anymore.
- 3 Adapt the method according to your needs, e.g. insert additional method functions and adjust method parameters.
- 4 To finish and store the method tap **[Save]**.
 - ➔ The new method has been saved.

See also

- 📖 Editing methods ▶ Page 62
- 📖 Method functions ▶ Page 71

8.4 Creating shortcuts for methods

Navigation: Home > Methods

Shortcuts can be placed on the **Home screen** and constitute links to executable methods. Shortcuts can be created from the **Start analysis** screen using the **AddToHome** button. The administration of shortcuts takes place in the setup of the instrument. This chapter describes how you can create shortcuts to start a method from the **Home screen**.

- 1 Select a method from the list.
 - ➔ The name of the method is displayed with all parameters.
- 2 Tap [**Start**].
 - ➔ **Start analysis** is displayed.
- 3 Tap [**AddToHome**].
 - ➔ **Shortcut parameters** is displayed.
- 4 Enter a meaningful text in **Description** that will be displayed in the shortcut.
- 5 Select **Immediate start**, to create a direct shortcut.
- 6 Tap the list field **Homescreen position**.
 - ➔ **Select location for shortcut** is displayed.
- 7 Select a free position.
- 8 Tap [**Save**] to save your entries.
 - ➔ The shortcut is displayed on the selected position of the home screen.

Note

- With direct shortcuts a series can be started immediately. Shortcuts lead to the **Start analysis** screen from where you can start the serial analysis.

See also

- 📖 Shortcuts ▶ Page 35

8.5 Editing methods

Navigation: Home > Methods

- 1 Select the method to be modified.
- 2 Select and edit the parameters to be modified.
- 3 Adapt the method according to your requirements, e.g. insert additional method functions and adjust method parameters.
- 4 Tap [**Save**] to save your entries.
 - ➔ The modified method has been saved.

Note

- Predefined METTLER TOLEDO methods cannot be modified. Changing the method ID will create a method copy, which can be modified according to your requirements.

The following tables show the predefined method functions and the method functions that can be added per method type.

Measure

Predefined Method Functions	Additional Method Functions
Title	< Instruction
Configuration	< Instruction
Sample	< Instruction < Wait/Stir < Sample verification < Sensor check < Measure (Temperature) < Auxiliary instrument
Measure	< Instruction < Measure < Calculation < Wait/Stir < Measure (Temperature) < Auxiliary instrument
Report	< Instruction < Measure < Calculation < Wait/Stir < Measure (Temperature) < Auxiliary instrument

Calibration

Predefined Method Functions	Additional Method Functions
Title	< Instruction
Configuration	< Instruction
Sample (Calibration)	< Instruction < Wait/Stir < Auxiliary instrument
Measure (Calibration)	< Instruction < Wait/Stir < Calibration analysis < Auxiliary instrument
Report	< Instruction < Wait/Stir < Calibration analysis < Auxiliary instrument

Interval

Predefined Method Functions	Additional Method Functions
Title	< Instruction
Configuration	< Instruction
Sample	< Instruction < Wait/Stir < Sample verification < Sensor check < Measure (Temperature) < Auxiliary instrument
Measure (Interval)	< Instruction < Measure (Interval) < Calculation < Wait/Stir < Measure (Temperature) < Auxiliary instrument
Report	< Instruction < Measure (Interval) < Calculation < Wait/Stir < Measure (Temperature) < Auxiliary instrument

Incremental

Predefined Method Functions	Additional Method Functions
Title	< Instruction
Configuration	< Instruction
Sample (Incremental)	< Instruction < Wait/Stir < Sample verification < Sensor check < Measure (Temperature) < Auxiliary instrument
Measure (Incremental)	< Instruction < Calculation < Wait/Stir < Measure (Temperature) < Auxiliary instrument
Report	< Instruction < Calculation < Wait/Stir < Measure (Temperature) < Auxiliary instrument

Sensor Test

Predefined Method Functions	Additional Method Functions
Title	< Instruction
Configuration	< Instruction
Sample (Sensor test)	< Instruction < Wait/Stir < Auxiliary instrument
Measure (Sensor test)	< Instruction < Wait/Stir < Auxiliary instrument
Sensor evaluation	< Instruction < Wait/Stir < Auxiliary instrument
Report	< Instruction < Wait/Stir < Auxiliary instrument

BOD

Predefined Method Functions	Additional Method Functions
Title	< Instruction
Configuration	< Instruction
Sample (BOD)	< Instruction < Wait/Stir < Sensor check < Auxiliary instrument
Measure (BOD)	< Instruction < Analysis (BOD) < Calculation < Wait/Stir < Auxiliary instrument
Report	< Instruction < Analysis (BOD) < Calculation < Wait/Stir < Auxiliary instrument

BCV

Predefined Method Functions	Additional Method Functions
Title	Instruction
Configuration	Instruction
Blank (BOD)	< Instruction < Wait/Stir < Sensor check < Auxiliary instrument
Measure (Blank)	< Instruction < Analysis (Blank) < Calculation < Wait/Stir < Auxiliary instrument
Seeded blank (BOD)	< Instruction < Wait/Stir < Sensor check < Auxiliary instrument
Measure (Seeded blank)	< Instruction < Analysis (Seeded blank) < Calculation < Wait/Stir < Auxiliary instrument
Standard (BOD)	< Instruction < Wait/Stir < Sensor check < Auxiliary instrument
Measure (Standard)	< Instruction < Analysis (BOD) < Calculation < Wait/Stir < Auxiliary instrument
Report	

8.6 Deleting methods

Navigation: Home > Methods

- 1 Select the method you want to delete.
- 2 Tap [**Delete method**].

➔ A message box appears with the information that references to the method will also be deleted.

3 Tap [**Delete**].

➔ The method has been deleted.

Note

- METTLER TOLEDO methods cannot be deleted.
- Referred shortcuts and sample series will also be deleted.

8.7 Method functions

8.7.1 Overview

Up to three measurement types can be defined in one method of the type measurement. The method functions will run sequentially from top to bottom. You can select the following functions to carry out the method. The total number of method functions is limited to 100.

Method function	Details	Used in method types	Maximum occurrences per method
Title	Method name, ID, author and SOP text	All	1
Configuration	Measurement types, sensors, temperature capture, barometric pressure capture, basic settings for BOD	All	1
Sample	Sample ID and other sample behaviors	MS, INT	1
Sample (Calibration)	Calibration standards, calibration mode	CAL	1
Sample (Incremental)	Type of incremental analysis, volume and concentration of added solutions	INC	1
Sample (Sensor test)	Calibration standards, drift test	ST	1
Blank (BOD)	BOD bottles used for blank value determination	BCV	1
Seeded blank (BOD)	BOD bottles used for seeded blank value determination	BCV	1
Standard (BOD)	BOD bottles used for standard value determination	BCV	1
Sample (BOD)	BOD bottles used for BOD measurement	BOD	1
Sample verification	Request to confirm of sample ID	MS, INT, INC	1
Sensor check	Check of calibration values, calibration date and calibration range of used sensor	MS, INT, INC, BCV, BOD	1
Measure	Measurement parameters like unit, resolution, endpoint type, stirring	MS	6
Measure (Calibration)	Measurement parameters like endpoint type and stirring	CAL	1
Measure (Interval)	Measurement parameters like unit, resolution, endpoint type, stirring	INT	6
Measure (Temperature)	Measurement parameters like temperature source, unit, resolution, endpoint type, stirring	MS, INT, INC	6
Measure (Incremental)	Measurement parameters like unit, resolution, endpoint type, stirring	INC	1
Measure (Sensor test)	Measurement parameters like endpoint type and stirring	ST	1
Measure (Blank)	Measurement parameters like unit, resolution, endpoint type, stirring for blank value	BCV	1
Measure (Seeded blank)	Measurement parameters like unit, resolution, endpoint type, stirring for seeded blank value	BCV	1
Measure (Standard)	Measurement parameters like unit, resolution, endpoint type, stirring for standard value	BCV	1
Measure (BOD)	Measurement parameters like unit, resolution, endpoint type, stirring	BOD	1
Calibration analysis	Check of calibration results like slope, offset or cell constant	CAL	1

Method function	Details	Used in method types	Maximum occurrences per method
Sensor evaluation	Check of sensor test results like slope, offset or drift	ST	1
Analysis (Blank)	Check for blank value results	BCV	1
Analysis (Seeded blank)	Check for seeded blank value results	BCV	1
Analysis (Standard)	Check for standard value results	BCV	1
Analysis (BOD)	Check for BOD results	BOD	1
Instruction	Message in a pop-up window	All	no limit
Wait/Stir	Waiting time or stirring duration	All	no limit
Calculation	Name, unit, formula und resolution to calculate a result	MS, INT, INC, BCV, BOD	20
Report	Content of the report	All	1
Auxiliary instrument	Sends data to or gets data from an auxiliary instrument	All	20

8.7.2 Title

This method function describes the contents and context of the method.

Parameter	Description	Values
Method type	Information on type of measurement.	Measurement Calibration Interval Incremental Sensor Test BCV BOD
Method ID	The method is uniquely identified in the system via the method ID. - Once the method has been saved, the method ID cannot be changed. A change of the method ID results in the creation of a copy of the method. METTLER TOLEDO methods begin with "M" followed by numbers (M is reserved for METTLER TOLEDO methods).	-
Title	Defines the name of the method.	-
Author	Information on the author (logged in user) of this method.	-
Created on	Information on the date of creation.	-
Modified on	Information on the date of modification.	-
Modified by	Information on the author (logged in user) of modification.	-
Protect	Protects the method against deletion or modification by other users than the author (logged in user) or administrator.	Active Inactive
SOP	Activates a text to be displayed before the start of the analysis.	Active Inactive
SOP text	Defines the SOP text. Displayed if check box SOP is activated.	-
Method type	Information on type of measurement.	Measurement Calibration Interval Incremental Sensor Test BCV BOD

8.7.3 Configuration

This method function defines the sensors to be used, how to capture temperature and barometric pressure, and basic parameters for BOD analysis.

Method types Measurement, Calibration, Interval, Incremental, Sensor Test

Parameter	Description	Values
Measurement type 1	Information on the measurement type 1.	pH Ion Conductivity Dissolved Oxygen Redox
Sensor name	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. A selected sensor is referenced by its name.	List of defined sensors
Check value ID	Defines the ID for all check values.	-
Temperature capture	Select the temperature capture mode. Internal: The used temperature sensor is integrated in the measurement sensor. External: A separate temperature sensor is used. Manual: The temperature is entered manually.	Internal External Manual
Barometric pressure capture	Select the barometric pressure capture mode. Displayed in case of Automatic Manual Measurement type 1 = Dissolved Oxygen	
Temperature sensor	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. Displayed if Temperature capture = External is selected.	List of defined sensors
Measurement type 2	Information on measurement type 2. Displayed if Measurement type 2 is defined.	pH Ion Conductivity Dissolved Oxygen Redox
Sensor name	Opens the sensor list, according to the sensor settings. Displayed if Measurement type 2 is selected.	-
Temperature capture	Different kinds of temperature captures can be selected. Displayed if Measurement type 2 is selected.	Internal External Manual
Barometric pressure capture	Select the barometric pressure capture mode. Displayed if Measurement type 2 = Dissolved Oxygen	Automatic Manual
Temperature sensor	Opens the list of temperature sensors, according to the sensor settings. Displayed if Measurement type 2 is selected in combination with Temperature capture = External .	List of defined sensors
Measurement type 3	Information on measurement type 3. Displayed if Measurement type 3 is defined.	pH Ion Conductivity Dissolved Oxygen Redox
Sensor name	Opens the sensor list, according to the sensor settings. Displayed if Measurement type 3 is selected.	List of defined sensors
Temperature capture	Different kinds of temperature captures can be selected. Displayed if Measurement type 3 is selected.	Internal External Manual
Barometric pressure capture	Select the barometric pressure capture mode. Displayed if Measurement type 3 = Dissolved Oxygen	Automatic Manual
Temperature sensor	Opens the list of temperature sensors, according to the sensor settings. Displayed if Measurement type 3 is selected in combination with Temperature capture = External .	List of defined sensors
Shared temperature value	If parallel measurements have been defined, one temperature sensor can be used for all measurements.	Active Inactive

Used temperature	Select the channel to be used for temperature capture. Displayed if Shared temperature value is activated.	Measurement type 1 Measurement type 2 Measurement type 3
------------------	--	--

Method type BCV

Parameter	Description	Values
Measurement type	Define the measurement type.	Dissolved Oxygen
Sensor name	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. A selected sensor is referenced by its name.	List of defined sensors
Check value ID	Defines the ID for all check values.	-
Blank	Activates the measurement of a blank value, which is the BOD in the dilution water including inorganic nutrients.	Active Inactive
Seeded blank	Activates the measurement of a seeded blank value, which is the BOD of a solution including all ingredients except the sample.	Active Inactive
Standard	Activates the measurement of a standard value, which is the BOD of a solution with a well-defined BOD value (typically 2% glucose glutamine acid).	Active Inactive
Salinity correction	Activates salinity corrections of several solutions. The reading of dissolved oxygen in the samples is corrected in regards of the entered Salinity of dilution water , Enter salinity of seed solution and Enter salinity of standard solution .	Active Inactive
Salinity of dilution water	Defines the salinity of the water used to dilute the samples. Displayed if check box Salinity correction is activated and any of the check boxes Blank , Seeded blank or Standard is activated.	0.0...70.0 ppt
Enter salinity of seed solution	Defines the salinity of the seed solution which is added to either the dilution water or the BOD bottle. Displayed if check box Salinity correction is activated and check box Seeded blank or Standard is activated.	0.0...70.0 ppt
Enter salinity of standard solution	Defines the salinity of the standard stock solution used to prepare the standards. Displayed if check boxes Salinity correction and Standard are activated.	0.0...70.0 ppt
Seed added	Defines to whether the seed is added to the bottle directly or to the dilution water first. If this selection is changed, all set volumes for sample, seed, and dilution water are reset to default values.	To bottle To dilution water
Seed dilution factor	Ratio of seed solution to dilution water. Displayed if Seed added = To dilution water is selected.	1.0...999.9
Blank correction	Activates the blank correction for BOD measurements. When activated, the BOD of dilution water including the inorganic nutrients used to dilute the samples is determined to check the presence of oxidizable compounds in the dilution water. The result will be subtracted from the BOD of the sample.	Active Inactive
Bottle volume	Defines the volume of used BOD bottle. If this value is changed, all set volumes for sample, seed, and dilution water are adapted proportionally.	10.0 - 1000.0 mL

Temperature capture	Select the temperature capture mode. Internal: The used temperature sensor is integrated in the measurement sensor. External: A separate temperature sensor is used. Manual: The temperature is entered manually.	Internal External Manual
Temperature sensor	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. Displayed if Temperature capture = External is selected.	List of defined sensors
Barometric pressure capture	Select the barometric pressure capture mode.	Automatic Manual

Method type BOD

Parameter	Description	Values
Measurement type	Define the measurement type.	Dissolved Oxygen
Sensor name	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. A selected sensor is referenced by its name.	List of defined sensors
Salinity correction	Activates salinity corrections of several solutions. The reading of dissolved oxygen in the samples is corrected in regards of the entered Salinity of dilution water, Enter salinity of seed solution and Enter salinity of standard solution .	Active Inactive
Enter salinity of seed solution	Defines the salinity of the seed solution which is added to either the dilution water or the BOD bottle. Displayed if check box Salinity correction is activated and check box Seeded blank or Standard is activated.	0.0...70.0 ppt
Salinity of dilution water	Defines the salinity of the water used to dilute the samples. Displayed if check box Salinity correction is activated and any of the check boxes Blank, Seeded blank or Standard is activated.	0.0...70.0 ppt
Salinity of undiluted sample	Defines the salinity of the sample before addition of seed or dilution water. Displayed if check box Salinity correction is activated.	0.0...70.0 ppt
Seed added	Defines to whether the seed is added to the bottle directly or to the dilution water first. If this selection is changed, all set volumes for sample, seed, and dilution water are reset to default values.	To bottle To dilution water
Seed dilution factor	Ratio of seed solution to dilution water. Displayed if Seed added = To dilution water is selected.	1.0...999.9
Blank correction	Activates the blank correction for BOD measurements. When activated, the BOD of dilution water including the inorganic nutrients used to dilute the samples is determined to check the presence of oxidizable compounds in the dilution water. The result will be subtracted from the BOD of the sample.	Active Inactive
Bottle volume	Defines the volume of used BOD bottle. If this value is changed, all set volumes for sample, seed, and dilution water are adapted proportionally.	10.0 - 1000.0 mL

Temperature capture	Select the temperature capture mode. Internal: The used temperature sensor is integrated in the measurement sensor. External: A separate temperature sensor is used. Manual: The temperature is entered manually.	Internal External Manual
Temperature sensor	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. Displayed if Temperature capture = External is selected.	List of defined sensors
Barometric pressure capture	Select the barometric pressure capture mode.	Automatic Manual

8.7.4 Sample

In this method function you can enter the sample ID and an optional comment that will be displayed in the Start analysis screen.

Parameter	Description	Values
Sample ID	Defines the sample ID.	-
Temperature	Defines the temperature for the measurement. Displayed if Temperature capture = Manual is selected.	-30...130 °C / -22...266 °F
Salinity correction factor	Defines the sample's salinity. The oxygen concentration depends on temperature and salinity of the sample. Displayed if Measurement type = Dissolved Oxygen .	-
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed if Measurement type = Dissolved Oxygen and Barometric pressure capture = Manual	500.0...1100.0 mbar
Comment	Defines a short comment that will be displayed.	-

8.7.5 Sample (Calibration)

In this method function buffer sets and standard groups and calibration modes can be defined.

Parameter	Description	Values
Measurement type	Information on the measurement type.	-
Sensor name	Information on the sensor name selected for the method.	-
Buffer group	Select a buffer group from the list of predefined and user-defined buffer groups. Displayed if Measurement type = pH .	-
Calibration standard group	Opens the list of predefined and user-defined standard groups. Displayed if Measurement type = Ion or Conductivity .	-
Unit	Information on the preselected unit. Displayed if Measurement type = Ion	-
Cal. mode	Specifies whether calibration data will be approximated with one straight line segment or left as a series of segments. Displayed if Measurement type = pH or Ion .	Linear Segmented
Number of buffers	Defines the number of buffers for the calibration. Displayed if Measurement type = pH .	1 2 3 4 5
Automatic buffer recognition	Activates the automatic recognition of buffers. Only available for predefined pH buffer groups. Displayed if Measurement type = pH .	Active Inactive

Buffer 1 - Buffer 5	The number of fields depends on the number of defined buffers. Up to 5 buffers are displayed with consecutive numbers. Select a buffer for each calibration point from the list. Displayed if Measurement type = pH and Automatic buffer recognition is not activated.	-
Number of standards	Defines the number of standards for the calibration (up to 5 for ion sensors, up to 2 for conductivity sensors). Displayed if Measurement type = Ion or Conductivity .	1 2 3 4 5
Standard 1 - Standard 5	Number of fields depends on the number of defined standards. Up to 5 standards for ion sensors and up to 2 standards for conductivity sensors are displayed with consecutive numbers. Each field opens the list of preselected standard groups. Displayed if Measurement type = Ion or Conductivity .	-
Calibration points	Option to choose between different number of calibration points. Displayed if Measurement type = Dissolved Oxygen .	1 2
Calibration standard 1	Select the standard for the first calibration point. Displayed if Measurement type = Dissolved Oxygen .	Saturated air
Calibration standard 2	Defines the standard for the second calibration point. Cannot be edited. Displayed if Measurement type = Dissolved Oxygen and Calibration points = 2 .	Zero point standard
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed if Barometric pressure capture = Manual is selected.	500.0...1100.0 mbar 500...1100 hPa 375...825 mmHg 0.493...1.086 atm
Temperature	Defines the temperature for the measurement. Displayed if Temperature capture = Manual is selected.	-30...130 °C / -22...266 °F
Comment	Defines a short comment that will be displayed.	-

8.7.6 Sample (Incremental)

In this method function you can define the basic parameters to perform a concentration determination. You can select the type of incremental measurement, enter the sample parameters and the concentration of the standard, the number of standard additions and the added volume.

Parameter	Description	Values
Sample ID	Defines the sample ID.	-
Measurement type	Information on the measurement type.	Ion
Incremental type	Select the incremental type.	Standard addition Standard subtraction Sample addition Sample subtraction
Sample volume	Defines the volume of the sample (mL is set). Displayed if Incremental type = Standard addition or Standard subtraction is selected.	-
Standard volume	Defines the volume of the standard (unit mL is set). Displayed if Incremental type = Sample addition or Sample subtraction is selected.	-
Standard unit	Defines the unit of the ion standard to be used.	mmol/L mol/L mg/L ppm % pX
Standard concentration	Defines the concentration of the added sample.	-

Ion ratio	Defines the ion ratio of the precipitated salts. Displayed if Incremental type = Standard subtraction or Sample subtraction is selected.	-
Number of standard additions	Defines the number of additions. Displayed if Incremental type = Standard addition is selected.	1 2 3 4 5
Added volume per addition	Defines the value for the volume to be added (unit mL is set).	0.010...1000 mL
Used slope	Defines the kind of slope to be detected.	From calibration Theoretical Enter known slope
Enter known slope	Defines the value of the known slope to be used (mV/pX is set). Displayed if Used slope Enter known slope is selected.	-
Comment	Defines a short comment that will be displayed.	-
Temperature	Defines the temperature for the measurement. Displayed if Temperature capture = Manual is selected.	-30...130 °C / -22...266 °F

8.7.7 Sample (Sensor test)

This method function gives an overview of the conditions used for the sensor test such as buffer sets and standard groups and drift.

Parameter	Description	Values
Measurement type	Information on the measurement type.	pH
Sensor name	Select a sensor from the list of defined sensors, if the method always runs with the same sensor. If not selected, you can choose any suitable sensor when you start the method. A selected sensor is referenced by its name.	List of defined sensors
Buffer group	Select a buffer group from the list of predefined and user-defined buffer groups.	-
Calibration	Activates using sensor test in context with calibration.	Active Inactive
Cal. mode	Specifies whether calibration data will be approximated with one straight line segment or left as a series of segments.	Linear Segmented
Number of buffers	Defines the number of buffers for the calibration.	1 2 3 4 5
Automatic buffer recognition	Activates the automatic recognition of buffers. Only available for predefined pH buffer groups.	Active Inactive
Buffer 1 - Buffer 5	The number of fields depends on the number of defined buffers. Up to 5 buffers are displayed with consecutive numbers. Select a buffer for each calibration point from the list of predefined and user defined buffer groups. Displayed if Automatic buffer recognition is not activated.	-
Temperature	Defines the temperature for the measurement. Displayed if Temperature capture = Manual is selected.	-30...130 °C / -22...266 °F
Drift test	Activate to start a drift test.	Active Inactive
Test measurements	Activates test measurements.	Active Inactive
No. of test measurements	Defines the number of measurements. Displayed if Test measurements is activated.	1 2 3 4 5

Test buffer 1 - Test buffer 5	Number of fields depends on the number of defined buffers. Up to 5 test buffers are displayed with consecutive numbers. Each field opens the list of predefined and user-defined buffer groups. Displayed if Test measurements is activated.	-
Comment	Defines a short comment that will be displayed.	-

8.7.8 Blank (BOD)

In this method function you can define the basic parameters to perform a BOD determination of a blank solution (dilution water including inorganic nutrients). You can define a comment and the number of bottles. In the bottle list you can enter the bottle ID and a comment.

Parameter	Description	Values
Blank ID	Information of the Blank ID entered in method function Configuration.	-
Comment	Defines a short comment that will be displayed.	-
Sample type	Information on the sample type.	Blank
Same bottles used for base and follow	Activates that DO is measured before (Base) and after (Follow) incubation time in the same bottles.	Active Inactive
Number of bottles (base)	Defines the number of bottles of this sample ID for base measurement (before incubation).	1...10
Number of bottles (follow)	Defines the number of bottles of this sample ID for follow measurement (after incubation). Displayed if check box Same bottles used for base and follow is de-activated.	1...10
Temperature	Defines the sample temperature. Displayed if Temperature capture = Manual is selected in method function Configuration .	0...60 °C / 32...140 °F
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed if Barometric pressure capture = Manual is selected.	500.0...1100.0 mbar 500...1100 hPa 375...825 mmHg 0.493...1.086 atm

Bottle settings

Parameter	Description	Values
Blank ID	Information of the Blank ID entered in method function Configuration.	-
Bottle ID	Defines the bottle ID.	-
Comment	Defines a short comment that will be displayed.	-
Seed volume	Defines the seed volume in the bottle, which equals always 0 mL.	-
Dilution volume	Information on the dilution volume in the bottle. Value is calculated from other volumes above.	-

8.7.9 Seeded blank (BOD)

In this method function you can define the basic parameters to perform a BOD determination of a seeded blank solution (dilution water including all ingredients except the sample). You can define a comment and the number of bottles. In the bottle list you can enter the bottle ID and a comment.

Parameter	Description	Values
Blank ID	Information of the Blank ID entered in method function Configuration.	-

Comment	Defines a short comment that will be displayed.	-
Sample type	Information on the sample type.	Blank
Same bottles used for base and follow	Activates that DO is measured before (Base) and after (Follow) incubation time in the same bottles.	Active Inactive
Number of bottles (base)	Defines the number of bottles of this sample ID for base measurement (before incubation).	1...10
Number of bottles (follow)	Defines the number of bottles of this sample ID for follow measurement (after incubation). Displayed if check box Same bottles used for base and follow is de-activated.	1...10
Temperature	Defines the sample temperature. Displayed if Temperature capture = Manual is selected in method function Configuration .	0...60 °C / 32...140 °F
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed if Barometric pressure capture = Manual is selected.	500.0...1100.0 mbar 500...1100 hPa 375...825 mmHg 0.493...1.086 atm

Bottle settings

Parameter	Description	Values
Seeded blank ID	Information on given seeded blank ID.	-
Bottle ID	Defines the bottle ID.	-
Comment	Defines a short comment that will be displayed.	-
Seed volume	Defines the seed volume in the bottle. The seed volume equals the bottle volume if Seed added = To dilution water in Configuration method function.	-
Dilution volume	Information on the dilution volume in the bottle. Value is calculated from other volumes above.	-

8.7.10 Standard (BOD)

In this method function you can define the basic parameters to perform a BOD determination of a seeded blank solution (dilution water including all ingredients except the sample). You can define a comment and the number of bottles. In the bottle list you can enter the bottle ID and a comment.

Parameter	Description	Values
Blank ID	Information of the Blank ID entered in method function Configuration.	-
Comment	Defines a short comment that will be displayed.	-
Sample type	Information on the sample type.	Blank
Same bottles used for base and follow	Activates that DO is measured before (Base) and after (Follow) incubation time in the same bottles.	Active Inactive
Number of bottles (base)	Defines the number of bottles of this sample ID for base measurement (before incubation).	1...10
Number of bottles (follow)	Defines the number of bottles of this sample ID for follow measurement (after incubation). Displayed if check box Same bottles used for base and follow is de-activated.	1...10

Temperature	Defines the sample temperature. Displayed if Temperature capture = Manual is selected in method function Configuration .	0...60 °C / 32...140 °F
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed if Barometric pressure capture = Manual is selected.	500.0...1100.0 mbar 500...1100 hPa 375...825 mmHg 0.493...1.086 atm

Bottle settings

Parameter	Description	Values
Standard ID	Information on given standard ID.	-
Bottle ID	Defines the bottle ID.	-
Comment	Defines a short comment that will be displayed.	-
Standard volume	Defines the standard volume in the bottle.	-
Seed volume	Defines the seed volume in the bottle. The seed volume equals the residual bottle volume if Seed added = To dilution water in Configuration method function.	-
Dilution volume	Information on the dilution volume in the bottle. Value is calculated from other volumes above.	-

8.7.11 Sample (BOD)

In this method function you can define the basic parameters to perform a BOD determination. You can the sample ID, the number of bottles, a corresponding comment and the seed addition unit. In the bottle list the you can enter the bottle ID, the sample volume, the volume of seed addition and the salinity correction factor.

Parameter	Description	Values
Sample ID	Defines the sample ID.	-
Comment	Defines a short comment that will be displayed.	-
Sample type	Information on the sample type.	Sample
Number of bottles (base)	Defines the number of bottles of this sample ID for base measurement (before incubation).	1...10
Number of bottles (follow)	Defines the number of bottles of this sample ID for follow measurement (after incubation). Displayed if check box Same bottles used for base and follow is de-activated.	1...10
Temperature	Defines the sample temperature. Displayed if Temperature capture = Manual is selected in method function Configuration .	0...60 °C / 32...140 °F
Barometric pressure	Defines the barometric (atmospheric) pressure. Displayed if Barometric pressure capture = Manual is selected.	500.0...1100.0 mbar 500...1100 hPa 375...825 mmHg 0.493...1.086 atm

Bottle settings

Parameter	Description	Values
Sample ID	Defines the sample ID.	-
Bottle ID	Defines the bottle ID.	-
Comment	Defines a short comment that will be displayed.	-

Sample volume	Defines the volume of sample in the bottle. Editable for base bottle. Editable for follow bottle only if checkbox Same bottles used for base and follow is de-activated.	-
Seed volume	Defines the seed volume in the bottle. Editable if checkbox Same bottles used for base and follow is de-activated and Seed added = To bottle is selected in method function Configuration.	-
Dilution volume	Information on the dilution volume in the bottle. Value is calculated from other volumes above.	-

8.7.12 Sample verification

This method function is mainly used in remote control mode. It requests the operator to verify if the correct sample is taken for the subsequent measurement. The sample ID defined in the method or at method start must be retyped or scanned with a barcode reader.

Parameter	Description	Values
Instruction	Text that will be displayed on the screen. Formula symbols can be used.	-
Interrupt after time span	The sample ID request will be interrupted when exceeding the time limit, if activated.	Active Inactive
Time	Defines the period of time for the method to be interrupted. Displayed if Interrupt after time span is activated.	-

See also

 Formula Syntax ► Page 117

8.7.13 Sensor check

This method function allows you to define special criteria for a sensor so that only sensors with a satisfying last calibration will be used. Calibration data can be checked and optionally the measurement can be interrupted if the sensor criteria are outside limits.

Parameter	Description	Values
Check offset and slope	Activates checking the offset and slope value of the last calibration. Displayed if Measurement type = pH or Ion .	Active Inactive
Min. slope	Defines the minimum slope value, in [%]. Displayed if Check offset and slope is enabled.	10...200
Max. slope	Defines the maximum slope value, in [%]. Displayed if Check offset and slope is enabled.	10...200
Min. offset	Defines the minimum offset value, in [mV]. Displayed if Check offset and slope is enabled.	-2000...2000
Max. offset	Defines the maximum offset value, in [mV]. Displayed if Check offset and slope is enabled.	-2000...2000
Check slope	Activates checking the offset and slope value of the last calibration. Displayed if Measurement type is Dissolved Oxygen .	Active Inactive
Min. slope	Defines the minimum slope value, in [%]. Displayed if Check slope is enabled.	10...200
Max. slope	Defines the maximum slope value, in [%]. Displayed if Check slope is enabled.	10...200

Check cell constant	Activates checking the offset and slope value of the last calibration. Displayed if Measurement type is Conductivity .	Active Inactive
Min. cell constant	Defines the minimum cell constant, in [°cm]. Displayed if Check cell constant is enabled.	0...100
Max. cell constant	Defines the maximum cell constant, in [°cm]. Displayed if Check cell constant is enabled.	0...100
Check calibration date	Activates checking the calibration date of the sensor.	Active Inactive
Monitoring period	Defines the unit for the monitoring period. Displayed if Check calibration date is enabled.	Days Hours
Max. elapsed period	Defines the maximum time period between calibration date and method execution, in days or hours.	1...100
Check calibration range	Activates checking the calibration ranges of the sensor.	Active Inactive
Min. calibration points	Defines the minimum number of calibration points within the below defined calibration range. Displayed if Check calibration range is enabled.	-
Unit for calibration range	Defines the unit for the below values of the calibration range. Displayed if Check calibration range is enabled.	-
Lower limit for calibration range	Defines the lower limit for the calibration range in the unit defined before. Displayed if Check calibration range is enabled.	-
Upper limit for calibration range	Defines the upper limit for the calibration range in the unit defined before. Displayed if Check calibration range is enabled.	-
Accept calibration points outside range	Defines if further calibration points outside the above defines calibration range can be accepted. Displayed if Check calibration range is enabled.	Active Inactive
Interrupt outside limits	Activates to interrupt the measurement when exceeding the limits.	Active Inactive

8.7.14 Measure

This method function defines all measurement parameters like unit, decimal places, endpoint type and temperature correction.

Measurement types pH, Redox, Ion, Dissolved oxygen

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Unit	The unit to be used for the measurement. The available units depend on the selected measurement type.	pH mV mg/L ppm mol/L mmol/L % pX Rel.mV
Offset	Defines the offset [mV]. Displayed if Measurement type = Redox and Unit = Rel.mV is selected.	-2000...2000 mV
Decimal places	Defines the number of digits for the displayed measurement result. The displayed decimal places depend on the selected unit.	1 2 3 4
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed

Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
dE	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	-
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

 Endpoint criteria ▶ Page 110

 Formula Syntax ▶ Page 117

Measurement type Conductivity

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Unit	Defines the measurement unit for conductivity.	µS/cm mS/cm S/m µS/m mS/m
Decimal places	Defines the number of digits for the displayed measurement result. The displayed decimal places depend on the selected unit.	1 2 3 4

Temperature correction	<p>Defines the relationship between conductivity, temperature and ion concentration.</p> <p>Linear: Use for the temperature correction of medium and highly conductive solutions.</p> <p>Non-linear: Use for natural water (only for temperature between 0...36 °C). The measured conductivity at the sample temperature is corrected to the defined reference temperature (20 °C or 25 °C).</p> <p>Off: The conductivity value at the current temperature is displayed.</p> <p>Pure water: An optimized type of temperature algorithm is used for Measurement type = Conductivity or Resistivity.</p> <p>Displayed if Measurement type = Conductivity, TDS or Resistivity is selected.</p>	Linear Non-linear Off Pure water
α-coefficient	<p>Defines the factor for the linear dependency.</p> <p>Displayed if Cond.mode = Conductivity, TDS or Resistivity in combination with Temperature correction = Linear is selected.</p>	0.00...10.00
Reference temperature	<p>The conductivity reading will be directly corrected to the set reference temperature.</p> <p>If Measurement type = Conductivity or Resistivity in combination with Temperature correction = Pure water is selected a reference temperature of 25 °C is automatically set.</p> <p>Displayed if Measurement type = Conductivity, TDS or Resistivity in combination with Temperature correction = Linear is selected.</p>	20°C 25°C
TDS factor	<p>To calculate the TDS value, the conductivity will be multiplied with this factor.</p> <p>Displayed if Conductivity = TDS is selected.</p>	0.00...10.00
Endpoint type	<p>Defines how the endpoint of the measurement is to be terminated.</p>	Automatic Manual Timed
Endpoint criteria	<p>Defines the parameters for the endpoint criteria.</p> <p>Displayed if Endpoint type = Automatic is selected.</p>	Strict Standard Fast User-defined
dE	<p>Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	
dt	<p>Defines the time component for dE. $dt > t_{min}$ and $t_{max} > dt$.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	1...600
tmin	<p>Earliest possible time for the measured value acquisition.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	1...100000
tmax	<p>Latest possible time for the measured value acquisition.</p> <p>Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	1...100000
Endpoint time	<p>Period of time [s] until the endpoint of the measurement is reached.</p> <p>Displayed if Endpoint type = Timed.</p>	5...1000000
Stir	<p>Activates the stirrer.</p>	Active Inactive

Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

 Endpoint criteria ▶ Page 110

 Formula Syntax ▶ Page 117

8.7.15 Measure (Calibration)

In this method function endpoint type and endpoint criteria are selected. You can also define whether stirring during calibration takes place or not.

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
dE	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	-
dt	Defines the time component for dE. $dt > t_{min}$ and $t_{max} > dt$. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%

See also

 Endpoint criteria ▶ Page 110

8.7.16 Measure (Temperature)

This method function defines all measurement parameters like temperature source, unit, resolution, endpoint type and stirring.

Parameter	Description	Values
Temperature signal	Select the channel to be used for the temperature capture.	Measurement type 1 Measurement type 2 Measurement type 3
Sensor name	Information on the sensor name selected for the method.	-
Unit	Information on the temperature unit.	°C °F
Decimal places	Defines the number of digits for the displayed measurement result.	0 1
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed Set value
Endpoint criteria	Defines the parameters for the endpoint criteria. Standard: Value varies less than 0.1 °C for the last 6 seconds. User-defined: The relevant settings are displayed. Displayed if Endpoint type = Automatic is selected.	Standard User-defined
dT	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dT, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	0.1 ... 1.0
dt	Defines the time component for dE. dt > tmin and tmax > dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1 ... 600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1 ... 100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1 ... 100000
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5 ... 1000000
Mode	Defines how the endpoint is reached when a set value is defined. Displayed if Endpoint type = Set value .	T > set value T < set value T within range
Set value	Defines the temperature for the set value. Displayed if Endpoint type = Set value and Mode = T > set value or T < set value .	pH/Redox Ion Conductivity: -30... 130 °C / -22... 266 °F DO: 0... 60 °C / 32... 140 °F
Lower limit	Defines the temperature for the lower limit. Displayed if Endpoint type = Set value and Mode = T within range .	See Set value

Upper limit	Defines the temperature for the upper limit. Displayed if Endpoint type = Set value and Mode = T within range .	See Set value
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the criteria for the set value of the endpoint type are not fulfilled. Displayed if Endpoint type = Set value .	1...100000
Mean value	The mean value of all measured values over a defined time span must exceed the set value before the measurement stops. Displayed if Endpoint type = Set value .	Active Inactive
Time span	Defines the time span to calculate the mean value.	1...60
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

 InMotion ▶ Page 30

 Formula Syntax ▶ Page 117

8.7.17 Measure (Interval)

This method function defines all measurement parameters like unit, decimal places, endpoint type and temperature correction.

Measurement types pH, Redox, Ion, Dissolved oxygen

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Unit	The unit to be used for the measurement. The available units depend on the selected measurement type.	pH mV mg/L ppm mol/L mmol/L % pX Rel.mV
Offset	Defines the offset [mV]. Displayed if Measurement type = Redox and Unit = Rel.mV is selected.	-2000...2000 mV
Decimal places	Defines the number of digits for the displayed measurement result. The displayed decimal places depend on the selected unit.	1 2 3 4
Interval time	Time period between starting and storing measurement data [s]. Only if method type = Interval .	-
Print after each interval	Activates printing out the result after each interval. Only if method type = Interval .	Active Inactive
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed

Endpoint criteria	<p>Defines the parameters for the endpoint criteria.</p> <p>Strict: Value varies less than 0.03 mV for the last 8 seconds or 0.1 mV for the last 20 seconds.</p> <p>Standard: Value varies less than 0.1 mV for the last 6 seconds.</p> <p>Fast: Value varies less than 0.6 mV for the last 4 seconds.</p> <p>User-defined: The relevant settings are displayed.</p> <p>Displayed if Endpoint type = Automatic is selected.</p>	Strict Standard Fast User-defined
dE	<p>Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	-
dt	<p>Defines the time component for dE. $dt > t_{min}$ and $t_{max} > dt$.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	1...600
tmin	<p>Earliest possible time for the measured value acquisition.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	1...100000
tmax	<p>Latest possible time for the measured value acquisition.</p> <p>Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled.</p> <p>Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined.</p>	1...100000
Endpoint time	<p>Period of time [s] until the endpoint of the measurement is reached.</p> <p>Displayed if Endpoint type = Timed.</p>	5...1000000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	<p>Defines the stirring speed.</p> <p>Displayed if Stir is activated.</p>	10%...100%
Threshold	Activates thresholds to be monitored, with optionally finding an endpoint of the measurement upon exceeding a threshold.	Active Inactive
Lower threshold	<p>Defines the value for the threshold.</p> <p>Displayed if Threshold is activated.</p>	-
Upper threshold	<p>Defines the value for the threshold.</p> <p>Displayed if Threshold is activated.</p>	-
Endpoint when threshold exceeded	<p>Defines that the endpoint is reached when the threshold value is exceeded.</p> <p>Displayed if Threshold is activated.</p>	Active Inactive
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	<p>Here you can enter a formula whose result (true or false) will determine the execution of the method function.</p> <p>Displayed if Condition is activated.</p>	-

See also

 [Formula Syntax](#) ▶ Page 117

Measurement type Conductivity

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Unit	Defines the measurement unit for conductivity.	μS/cm mS/cm S/m μS/m mS/m
Decimal places	Defines the number of digits for the displayed measurement result. The displayed decimal places depend on the selected unit.	1 2 3 4
Interval time	Time period between starting and storing measurement data [s]. Only if method type = Interval .	-
Print after each interval	Activates printing out the result after each interval. Only if method type = Interval .	Active Inactive
Temperature correction	Defines the relationship between conductivity, temperature and ion concentration. Linear : Use for the temperature correction of medium and highly conductive solutions. Non-linear : Use for natural water (only for temperature between 0...36 °C). The measured conductivity at the sample temperature is corrected to the defined reference temperature (20 °C or 25 °C). Off : The conductivity value at the current temperature is displayed. Pure water : An optimized type of temperature algorithm is used for Measurement type = Conductivity or Resistivity . Displayed if Measurement type = Conductivity, TDS or Resistivity is selected.	Linear Non-linear Off Pure water
α-coefficient	Defines the factor for the linear dependency. Displayed if Cond.mode = Conductivity, TDS or Resistivity in combination with Temperature correction = Linear is selected.	0.00...10.00
Reference temperature	The conductivity reading will be directly corrected to the set reference temperature. If Measurement type = Conductivity or Resistivity in combination with Temperature correction = Pure water is selected a reference temperature of 25 °C is automatically set. Displayed if Measurement type = Conductivity, TDS or Resistivity in combination with Temperature correction = Linear is selected.	20°C 25°C
TDS factor	To calculate the TDS value, the conductivity will be multiplied with this factor. Displayed if Conductivity = TDS is selected.	0.00...10.00
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
dE	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600

tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
Threshold	Activates thresholds to be monitored, with optionally finding an endpoint of the measurement upon exceeding a threshold.	Active Inactive
Lower threshold	Defines the value for the threshold. Displayed if Threshold is activated.	-
Upper threshold	Defines the value for the threshold. Displayed if Threshold is activated.	-
Endpoint when threshold exceeded	Defines that the endpoint is reached when the threshold value is exceeded. Displayed if Threshold is activated.	Active Inactive
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

 Endpoint criteria ▶ Page 110

 Formula Syntax ▶ Page 117

8.7.18 Measure (Incremental)

In this method function the endpoint type, endpoint criteria and other important parameters are determined. You can also define whether stirring takes place during the measurement or not.

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Unit	Defines the unit for ion.	mmol/L mg/L ppm % pX
Decimal places	Defines the number of digits for the displayed measurement result. The displayed decimal places depend on the selected unit.	1 2 3 4
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000

dE	Defines the measured value interval. As soon as the change in the - measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%

See also

 Endpoint criteria ▶ Page 110

8.7.19 Measure (Sensor test)

In this method function endpoint type and endpoint criteria are determined. You can also define whether stirring takes place during the measurement or not.

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
dE	Defines the measured value interval. As soon as the change in the - measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000

tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%

See also

 Endpoint criteria ▶ Page 110

8.7.20 Measure (Blank)

In this method function you can define the measuring units and resolution as well as endpoint type and criteria for BOD determination in blank. You can also define whether stirring takes place during the measurement or not.

Parameter	Description	Values
Sensor name	Information on the sensor selected for this method.	-
DO unit	Defines the unit for single DO measurement.	mg/L ppm
BOD unit	Defines the unit for BOD result.	mg/L
DO resolution	Defines the number of digits for the displayed DO result.	1 2 3
BOD resolution	Defines the number of digits for the calculated BOD value.	1 2 3
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
dE	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	-
dt	Defines the time component for dE. dt > tmin and tmax > dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive

Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%
----------------	---	------------

See also

 Endpoint criteria ▶ Page 110

8.7.21 Measure (Seeded blank)

In this method function you can define the measuring units and resolution as well as endpoint type and criteria for BOD determination in seeded blank. You can also define whether stirring takes place during the measurement or not.

Parameter	Description	Values
DO unit	Defines the unit for single DO measurement.	mg/L ppm
BOD unit	Defines the unit for BOD result.	mg/L
DO resolution	Defines the number of digits for the displayed DO result.	1 2 3
BOD resolution	Defines the number of digits for the calculated BOD value.	1 2 3
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
dE	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%

See also

 Endpoint criteria ▶ Page 110

8.7.22 Measure (Standard)

In this method function you can define the measuring units and resolution as well as endpoint type and criteria for BOD determination in the standard solution (solution with well-defined BOD value). You can also define whether stirring takes place during the measurement or not.

Parameter	Description	Values
DO unit	Defines the unit for single DO measurement.	mg/L ppm
BOD unit	Defines the unit for BOD result.	mg/L
DO resolution	Defines the number of digits for the displayed DO result.	1 2 3
BOD resolution	Defines the number of digits for the calculated BOD value.	1 2
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed
Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
dE	Defines the measured value interval. As soon as the change in the - measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%

See also

 Endpoint criteria ▶ Page 110

8.7.23 Measure (BOD)

In this method function you can define the measuring units and resolution as well as endpoint type and criteria. You can also define whether stirring takes place during the measurement or not.

Parameter	Description	Values
Sensor name	Information on the sensor selected for this method.	-
DO unit	Defines the unit for single DO measurement.	mg/L ppm
BOD unit	Defines the unit for BOD result.	mg/L
DO resolution	Defines the number of digits for the displayed DO result.	1 2 3
BOD resolution	Defines the number of digits for the calculated BOD value.	1 2 3
Endpoint type	Defines how the endpoint of the measurement is to be terminated.	Automatic Manual Timed

Endpoint criteria	Defines the parameters for the endpoint criteria. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Endpoint time	Period of time [s] until the endpoint of the measurement is reached. Displayed if Endpoint type = Timed .	5...1000000
dE	Defines the measured value interval. As soon as the change in the measured value over the time period dt is less than dE, the measured value will be acquired. This occurs within the defined time interval. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	
dt	Defines the time component for dE. dt>tmin and tmax>dt. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...600
tmin	Earliest possible time for the measured value acquisition. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
tmax	Latest possible time for the measured value acquisition. Note The measurement ends after the time defined, even if the stability criteria of dE and dt are not fulfilled. Displayed if Endpoint type = Automatic and Endpoint criteria = User-defined .	1...100000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10%...100%

See also

 Endpoint criteria ▶ Page 110

8.7.24 Calibration analysis

In this method function the limits for accepting a calibration can be entered.

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Min. slope	Defines the lower limit for the slope, in [%]. Displayed if Measurement type = pH, Ion or Dissolved Oxygen .	10...200
Max. slope	Defines the upper limit for the slope, in [%]. Displayed if Measurement type = pH, Ion or Dissolved Oxygen .	10...200
Min. offset	Defines the lower limit for the offset, in [mV]. Displayed if Measurement type = pH or Ion .	-2000...2000
Max. offset	Defines the upper limit for the offset, in [mV]. Displayed if Measurement type = pH or Ion .	-2000...2000
Min. cell constant	Defines the lower limit for the cell constant [cm ⁻¹]. Displayed if Measurement type = Conductivity .	0...100
Max. cell constant	Defines the upper limit for the cell constant [cm ⁻¹]. Displayed if Measurement type = Conductivity .	0...100
Interrupt outside limits	Activates to interrupt the measurement when exceeding the limits.	Active Inactive
Display actual pH and temp. values	Activates the option to display additional values in the results. Displayed if Measurement type = pH .	Active Inactive

8.7.25 Sensor evaluation

In this method function you can set calibration limits and test measurement limits for sensors in the method type sensor test.

Parameter	Description	Values
Sensor name	Information on the sensor name selected for the method.	-
Calibration limits	Activates the parameters to set the limits.	Active Inactive
Min. slope	Defines the lower limit for the slope, in [%]. Displayed if Calibration limits is activated.	10...200
Max. slope	Defines the upper limit for the slope, in [%]. Displayed if Calibration limits is activated.	10...200
Min. offset	Defines the lower limit for the offset, in [mV]. Displayed if Calibration limits is activated.	-2000...2000
Max. offset	Defines the upper limit for the offset, in [mV]. Displayed if Calibration limits is activated.	-2000...2000
Max. drift	Define values for maximum drift during the 5 minutes drift test, in [mV]. Displayed if Calibration limits is activated.	0...2000
Test meas. limits	Activate to set a tolerance value and to determine interruption of the method if the sensor is outside limits.	Active Inactive
Tolerance	Define values for maximum difference between measured and theoretical value, in [pH]. Displayed if Test meas. limits is activated.	0.01...1.00

8.7.26 Analysis (Blank)

In this method function different limiting parameters for the BOD blank determination can be selected and modified in order to create related warnings, entries in the report, suspension of measurement or even interruption of measurement. The method function has three parts: **Analysis (base)**, **Analysis (follow)**, and **Analysis results**.

Analysis (base)

Parameter	Description	Values
Temperature limits	Defines if temperature limits are applied.	Active Inactive
Max. temperature	Defines the upper temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Min. temperature	Defines the lower temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Action when outside limits	Defines the behavior in case the temperature is outside the limits. Displayed if Temperature limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the temperature is outside the limits. Displayed if Temperature limits is activated.	Active Inactive
Max. DO limit	Defines if an upper limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Max. DO	Defines the upper oxygen content limit, in [%]. Displayed if Max. DO limit is activated.	90...200

Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Active Inactive
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1 ... 7.0
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive

Analysis (follow)

Parameter	Description	Values
Time tolerance limit	Defines if a time tolerance between base and follow measurement is applied	Active Inactive
Time tolerance	Defines the time tolerance to full days between measurement Base and Follow , in [h]. Example: a time tolerance of 3 hours is applied. The measurement Base is done at 10:00 a.m. The measurement Follow must be started at any later day between 7:00 a.m. and 1:00 p.m. Displayed if Time tolerance limit is activated.	0.1 ... 12.0
Action when outside limits	Defines the behavior in case the time tolerance exceeds the limit. Displayed if Time tolerance limit is activated.	Disable start Save and warn

Analysis results

Parameter	Description	Values
Max. BOD limit of bottle	Defines if an upper limit for the BOD value per bottle is applied.	Active Inactive
Max. BOD	Defines the upper oxygen content limit, in [mg/L]. Displayed if Max. BOD limit of bottle is activated.	0.1 ... 15.0
Action when outside limits	Defines the behavior in case the BOD value is below the limit. Displayed if Max. BOD limit of bottle is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the BOD value drops below the limit. Displayed if Max. BOD limit of bottle is activated.	Active Inactive

Action when outside limits

The following table explains the different actions when limits are exceeded.

Save and report	The measured data is saved and marked as outside limit; the method proceeds with the next bottle.
Repeat	The measured data is discarded and the last measurement must be repeated. This happens for infinite times if the limits are not met, except the task is interrupted.
Skip bottle	The measured data is discarded; the method proceeds with the next bottle.
Interrupt	The running task is terminated.

8.7.27 Analysis (Seeded blank)

In this method function different limiting parameters for the BOD seeded blank determination can be selected and modified in order to create related warnings, entries in the report, suspension of measurement or even interruption of measurement. The method function has three parts: **Analysis (base)**, **Analysis (follow)**, and **Analysis results**.

Analysis (base)

Parameter	Description	Values
Temperature limits	Defines if temperature limits are applied.	Active Inactive
Max. temperature	Defines the upper temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Min. temperature	Defines the lower temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Action when outside limits	Defines the behavior in case the temperature is outside the limits. Displayed if Temperature limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the temperature is outside the limits. Displayed if Temperature limits is activated.	Active Inactive
Max. DO limit	Defines if an upper limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Max. DO	Defines the upper oxygen content limit, in [%]. Displayed if Max. DO limit is activated.	90...200
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Active Inactive
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1...7.0
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive

Analysis (follow)

Parameter	Description	Values
Time tolerance limit	Defines if a time tolerance between base and follow measurement is applied	Active Inactive
Time tolerance	Defines the time tolerance to full days between measurement Base and Follow , in [h]. Example: a time tolerance of 3 hours is applied. The measurement Base is done at 10:00 a.m. The measurement Follow must be started at any later day between 7:00 a.m. and 1:00 p.m. Displayed if Time tolerance limit is activated.	0.1 ... 12.0
Action when outside limits	Defines the behavior in case the time tolerance exceeds the limit. Displayed if Time tolerance limit is activated.	Disable start Save and warn
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Follow) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1 ... 15.0
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive

Analysis results

Parameter	Description	Values
Min. BOD limit of bottle	Defines if a lower limit for the BOD value per bottle is applied	Active Inactive
Min. BOD	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. BOD limit of bottle is activated.	0.3 ... 15.0
Limit applied to	Defines if the min. BOD limit is applied to the BOD calculated from measurement (Base and Follow) or correction coming from a BCV method are included. Displayed if Min. BOD limit of bottle is activated.	Corrected BOD value Uncorrected BOD value
Action when outside limits	Defines the behavior in case the BOD value drops below the limit. Displayed if Min. BOD limit of bottle is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the BOD value drops below the limit. Displayed if Min. BOD limit of bottle is activated.	Active Inactive
O ₂ depletion limits	Defines if limits oxygen depletion are applied. The oxygen depletion is the ratio of oxygen that is gone between measurement Base and Follow .	Active Inactive
Min. O ₂ depletion	Defines the lower oxygen depletion limit, in [%] Displayed if O₂ depletion limits is activated.	0 ... 100
Max. O ₂ depletion	Defines the upper oxygen depletion limit, in [%]. Displayed if O₂ depletion limits is activated.	0 ... 100

Action when outside limits	Defines the behavior in case the correction factor exceeds the limits. Displayed if O₂ depletion limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the correction factor exceeds the limits. Displayed if O₂ depletion limits is activated.	Active Inactive

Action when outside limits

The following table explains the different actions when limits are exceeded.

Save and report	The measured data is saved and marked as outside limit; the method proceeds with the next bottle.
Repeat	The measured data is discarded and the last measurement must be repeated. This happens for infinite times if the limits are not met, except the task is interrupted.
Skip bottle	The measured data is discarded; the method proceeds with the next bottle.
Interrupt	The running task is terminated.

8.7.28 Analysis (Standard)

In this method function different limiting parameters for the BOD standard determination can be selected and modified in order to create related warnings, entries in the report, suspension of measurement or even interruption of measurement. The method function has three parts: **Analysis (base)**, **Analysis (follow)**, and **Analysis results**.

Parameter	Description	Values
Temperature limits	Defines if temperature limits are applied.	Active Inactive
Max. temperature	Defines the upper temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Min. temperature	Defines the lower temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Action when outside limits	Defines the behavior in case the temperature is outside limits. Displayed if Temperature limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the temperature is outside the limits. Displayed if Temperature limits is activated.	Active Inactive
Max. DO limit	Defines if an upper limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Max. DO	Defines the upper oxygen content limit, in [%]. Displayed if Max. DO limit is activated.	90...200
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Active Inactive
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1...7.0

Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive

Analysis (follow)

Parameter	Description	Values
Time tolerance limit	Defines if a time tolerance between base and follow measurement is applied	Active Inactive
Time tolerance	Defines the time tolerance to full days between measurement Base and Follow , in [h]. Example: a time tolerance of 3 hours is applied. The measurement Base is done at 10:00 a.m. The measurement Follow must be started at any later day between 7:00 a.m. and 1:00 p.m. Displayed if Time tolerance limit is activated.	0.1 ... 12.0
Action when outside limits	Defines the behavior in case the time tolerance exceeds the limit. Displayed if Time tolerance limit is activated.	Disable start Save and warn
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Follow) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1 ... 15.0
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive

Analysis results

Parameter	Description	Values
Min. BOD limit of bottle	Defines if a lower limit for the BOD value per bottle is applied	Active Inactive
Min. BOD	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. BOD limit of bottle is activated.	0.3 ... 15.0
Limit applied to	Defines if the min. BOD limit is applied to the BOD calculated from measurement (Base and Follow) or correction coming from a BCV method are included. Displayed if Min. BOD limit of bottle is activated.	Corrected BOD value Uncorrected BOD value
Action when outside limits	Defines the behavior in case the BOD value drops below the limit. Displayed if Min. BOD limit of bottle is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the BOD value drops below the limit. Displayed if Min. BOD limit of bottle is activated.	Active Inactive

Seed correction limits	Defines if limits for the seed correction factor are applied. The seed correction factor is that ratio of the calculated BOD that originates from the added seed and not from the standard. This option is only useful if values from a seeded blank determination are available.	Active Inactive
Min. correction factor	Defines the lower correction factor limit. Displayed if Seed correction limits is activated.	0.1...9.9
Max. correction factor	Defines the upper correction factor limit. Displayed if Seed correction limits is activated.	0.1...9.9
Action when outside limits	Defines the behavior in case the correction factor exceeds the limits. Displayed if Seed correction limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the correction factor exceeds the limits. Displayed if Seed correction limits is activated.	Active Inactive
BOD limits of standard	Defines if limits for the BOD of the sample (calculated over all bottles) are applied.	Active Inactive
Max. BOD	Defines the upper BOD limit, in [mg/L]. Displayed if BOD limits of standard is activated.	0.1...1000
Min. BOD	Defines the lower BOD limit, in [mg/L]. Displayed if BOD limits of standard is activated.	0.1...1000
Action when outside limits	Defines the behavior in case the sample BOD exceeds the limits. Displayed if BOD limits of standard is activated.	Save and report Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the sample BOD exceeds the limits. Displayed if BOD limits of standard is activated.	Active Inactive

Action when outside limits

The following table explains the different actions when limits are exceeded.

Save and report	The measured data is saved and marked as outside limit; the method proceeds with the next bottle.
Repeat	The measured data is discarded and the last measurement must be repeated. This happens for infinite times if the limits are not met, except the task is interrupted.
Skip bottle	The measured data is discarded; the method proceeds with the next bottle.
Interrupt	The running task is terminated.

8.7.29 Analysis (BOD)

In this method function different limiting parameters can be selected and modified in order to create related warnings, entries in the report, suspension of measurement or even interruption of measurement. The method function has three parts: **Analysis (base)**, **Analysis (follow)**, and **Analysis results**.

Analysis (base)

Parameter	Description	Values
Temperature limits	Defines if temperature limits are applied.	Active Inactive
Max. temperature	Defines the upper temperature limit. Displayed if Temperature limits is activated.	5...40 °C
Min. temperature	Defines the lower temperature limit. Displayed if Temperature limits is activated.	5...40 °C

Action when outside limits	Defines the behavior in case the temperature is outside limits. Displayed if Temperature limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the temperature is outside the limits. Displayed if Temperature limits is activated.	Active Inactive
Max. DO limit	Defines if an upper limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Max. DO	Defines the upper oxygen content limit, in [%]. Displayed if Max. DO limit is activated.	90...200
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Max. DO limit is activated.	Active Inactive
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Base) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1...7.0
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive

Analysis (follow)

Parameter	Description	Values
Time tolerance limit	Defines if a time tolerance between base and follow measurement is applied	Active Inactive
Time tolerance	Defines the time tolerance to full days between measurement Base and Follow , in [h]. Example: a time tolerance of 3 hours is applied. The measurement Base is done at 10:00 a.m. The measurement Follow must be started at any later day between 7:00 a.m. and 1:00 p.m. Displayed if Time tolerance limit is activated.	0.1...12.0
Action when outside limits	Defines the behavior in case the time tolerance exceeds the limit. Displayed if Time tolerance limit is activated.	Disable start Save and warn
Min. DO limit	Defines if a lower limit for oxygen content in the measurement (Follow) is applied.	Active Inactive
Min. DO	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. DO limit is activated.	0.1...15.0
Action when outside limits	Defines the behavior in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Save and report Repeat Skip bottle Interrupt

Show instruction	Defines if an automatically generated instruction will be displayed in case the oxygen content exceeds the limit. Displayed if Min. DO limit is activated.	Active Inactive
------------------	--	-------------------

Analysis results

Parameter	Description	Values
Min. BOD limit of bottle	Defines if a lower limit for the BOD value per bottle is applied	Active Inactive
Min. BOD	Defines the lower oxygen content limit, in [mg/L]. Displayed if Min. BOD limit of bottle is activated.	0.1 ... 15.0
Limit applied to	Defines if the min. BOD limit is applied to the BOD calculated from measurement (Base and Follow) or correction coming from a BCV method are included. Displayed if Min. BOD limit of bottle is activated.	Corrected BOD value Uncorrected BOD value
Action when outside limits	Defines the behavior in case the BOD value drops below the limit. Displayed if Min. BOD limit of bottle is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the BOD value drops below the limit. Displayed if Min. BOD limit of bottle is activated.	Active Inactive
Seed correction limits	Defines if limits for the seed correction factor are applied. The seed correction factor is that ratio of the calculated BOD that originates from the added seed and not from the sample. This option is only useful if check values from a BCV method are available.	Active Inactive
Min. correction factor	Defines the lower correction factor limit. Displayed if Seed correction limits is activated.	0.1 ... 9.9
Max. correction factor	Defines the upper correction factor limit. Displayed if Seed correction limits is activated.	0.1 ... 9.9
Action when outside limits	Defines the behavior in case the correction factor exceeds the limits. Find more details further below. Displayed if Seed correction limits is activated.	Save and report Repeat Skip bottle Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the correction factor exceeds the limits. Displayed if Seed correction limits is activated.	Active Inactive
BOD limits of sample	Defines if limits for the BOD of the sample (calculated over all bottles) are applied.	Active Inactive
Max. BOD	Defines the lower BOD limit, in [mg/L]. Displayed if BOD limits of sample is activated.	0.01 ... 1000000
Min. BOD	Defines the upper BOD limit, in [mg/L]. Displayed if BOD limits of sample is activated.	0.01 ... 1000000
Action when outside limits	Defines the behavior in case the sample BOD exceeds the limits. Displayed if BOD limits of sample is activated.	Save and report Interrupt
Show instruction	Defines if an automatically generated instruction will be displayed in case the sample BOD exceeds the limits. Displayed if BOD limits of sample is activated.	Active Inactive

Action when outside limits

The following table explains the different actions when limits are exceeded.

Save and report	The measured data is saved and marked as outside limit; the method proceeds with the next bottle.
Repeat	The measured data is discarded and the last measurement must be repeated. This happens for infinite times if the limits are not met, except the task is interrupted.
Skip bottle	The measured data is discarded; the method proceeds with the next bottle.
Interrupt	The running task is terminated.

8.7.30 Instruction

In this method function you can enter a text which will be displayed on the screen and you can set the conditions when the text will disappear. There are two possibilities either the text disappears after a predefined period of time or after confirmation.

Parameter	Description	Values
Instruction	Text that will be displayed on the screen. Formula symbols can be used.	-
Continue after	There are two possibilities either the text disappears after a predefined period of time or after confirmation.	Confirmation Time span
Time	Defines the period of time for the text to disappear. Displayed if Continue after = Time span .	-
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

 Formula Syntax ▶ Page 117

8.7.31 Wait/Stir

In this method function you can set a period of time for a pause before the next method function starts. You can define to stir during the waiting period.

Parameter	Description	Values
Wait time	Time [s] the method waits or the stirrer is activated.	1 ... 1000000
Stir	Activates the stirrer.	Active Inactive
Stirring speed	Defines the stirring speed. Displayed if Stir is activated.	10% ... 100%
Instruction	Activates the option to display a text on the screen, after the waiting / stirring time has elapsed.	Active Inactive
Text	Enter a text which will be displayed on the screen. Formula symbols can be used. Displayed if Instruction is activated.	-
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

 Formula Syntax ▶ Page 117

8.7.32 Calculation

In this method function you can enter a calculation based on the results of your measurement. You can also set result limits and determine to interrupt the measurement when exceeding the limits.

Parameter	Description	Values
Name	Defines the name of the calculation.	-
Unit	Enter the unit that will be displayed for the calculation.	-
Formula	Enter a calculation based on the results of your measurement.	-
Decimal places	Defines the number of digits for the displayed measurement result.	1 ... 6
Result limits	Defines whether limits should be observed for the result. If this function is activated, there will be a message in the record if the result falls outside the defined limits.	Active Inactive
Lower limit	Defines the lower result limit. Appears only if Result limits is activated.	-10^8 ... 10^8
Upper limit	Defines the upper result limit. Appears only if Result limits is activated.	-10^8 ... 10^8
Interrupt outside limits	Activates to interrupt the measurement when exceeding the limits.	Active Inactive

See also

 Formula Syntax ▶ Page 117

8.7.33 Report

The details for creating a report, printing or exporting data can be defined here.




Parameter	Description	Values
Print	Defines if data will be printed to the connected printer.	Active Inactive
Print format	Summary: Covers most important data related to date, time, user name, sample ID, sensor name, value, temperature, results, endpoint type and most important parameters according to the settings of the measurement type. User-defined: Allows to define which information to include.	Summary User-defined

Selectable parameters if User-defined is activated

Parameter	Description	Values
Measured values	Interval and endpoint values or calibration results can be exported or printed.	Active Inactive
Raw values	Raw values of the measurements can be exported or printed.	Active Inactive
Calculated results	Results of method function Calculation can be exported or printed. For method type Sensor Test the findings of method function Sensor evaluation is exported or printed; for method type Incremental it is the final sample concentration.	Active Inactive
Status	The overall status of the analysis can be exported or printed.	Active Inactive
Date/Time	The date and time of analysis execution can be exported or printed.	Active Inactive
User name	The name of the user who performed the analysis can be exported or printed. Avoid names longer than 10 characters in case of printing with a compact printer.	Active Inactive

Sensor name	The name of the used sensor can be exported or printed. Avoid names longer than 10 characters in case of printing with a compact printer.	Active Inactive
Sensor details	The sensor serial number, last calibration date, temperature sensor name, and further sensor details can be exported or printed.	Active Inactive
Sample ID	The sample ID can be exported or printed. For BOD and BCV method type the bottle ID is used. Avoid IDs longer than 10 characters in case of printing with a compact printer.	Active Inactive
Sample details	Further sample details can be exported or printed. For Calibration and Sensor Test the buffers or standards are used; for method type Incremental the ion ratio, standard additions and standard volume. For method types BOD and BCV the sample volume, seed volume, salinity, and further values are used.	Active Inactive
Method data	The method ID and measurement types can be exported or printed.	Active Inactive
Measurement details	The measurement parameters like temperature capture, temperature correction, interval time, stirrer speed, calibration mode, number of BOD bottles, and so on can be exported or printed.	Active Inactive
Endpoint	The measurement endpoint settings can be exported or printed.	Active Inactive
Endpoint criteria	Defines the parameters for the endpoint criteria. Strict: Value varies less than 0.03 mg/L during the last 20 seconds. Standard: Value varies less than 0.08 mg/L during the last 20 seconds. Fast: Value varies less than 0.08 mg/L during the last 10 seconds. Displayed if Endpoint type = Automatic is selected.	Strict Standard Fast User-defined
Instrument data	The instrument ID and serial number, module type and serial number, and last radio clock synchronization time can be exported or printed.	Active Inactive

See also

-  Printer ▶ Page 31
-  Header and footer ▶ Page 38
-  Printing analysis data ▶ Page 116




8.7.34 Auxiliary instrument

In this method function you can define data which is sent to an auxiliary instrument or to wait for incoming data from an auxiliary instrument and how to treat that data.

Parameter	Description	Values
Control type	Type of connection to the auxiliary instrument.	USB-RS232
Name	Defines which of the auxiliary instruments defined in the setup is used.	List of auxiliary instruments
Send out sequence	Defines if the method sends data to the auxiliary instrument.	Active Inactive

Output sequence	Defines the string which is sent out. Besides text also raw values - embracketed with % signs can be used. To send a specific ASCII character, the format \xxx must be used with x = digits, e.g. \010 for a line break. Displayed if Send out sequence is activated.	-
Wait for input sequence	Defines if the method waits for incoming data from an auxiliary instrument.	Active Inactive
Maximum time	Defines the duration which the methods waits at maximum to get an input sequence. Displayed if Wait for input sequence is activated.	0...1000000 s infinite
Input sequence	Defines the exact input sequence for which the method is waiting. - Displayed if Wait for input sequence is activated and Input sequence with results is deactivated.	-
Input sequence with results	Defines if the input sequences contains information that is to be stored in raw value AuxInstr. If activated, for example measurement results of the auxiliary instrument can be used later in the method. Displayed if Wait for input sequence is activated.	Active Inactive
Start sequence	Defines the beginning of the incoming data and hence allows to start the separation into results at a specific position. Displayed if Wait for input sequence and Input sequence with results are activated.	-
Total length	Defines the total length of the incoming data. The separation into results does not start before this number of characters are received. Additional characters are cut off. The total length must be at least the sum of the lengths of each result. Displayed if Wait for input sequence and Input sequence with results are activated.	1...1000
Number of results	Defines how many results are extracted from the incoming data string. Displayed if Wait for input sequences and Input sequence with results are activated. For each result, a "Start position" and "Max. length" has to be defined. The instrument tries to identify a number in that section and stores it in the respective AuxInstr value; other characters are ignored.	1...10
Condition	A logical condition can be defined. The method function is executed or not on the result (true or false) of the calculation.	Active Inactive
Formula	Here you can enter a formula whose result (true or false) will determine the execution of the method function. Displayed if Condition is activated.	-

See also

-  Configuration ▶ Page 72
-  Accessories ▶ Page 125
-  Formula Syntax ▶ Page 117

8.8 Endpoint criteria

Measurement type	Endpoint criteria		
	Strict	Standard	Fast
pH or Redox	Value varies less than 0.03 mV during the last 8 seconds or 0.1 mV for the last 20 seconds.	Value varies less than 0.1 mV during the last 6 seconds.	Value varies less than 0.6 mV during the last 4 seconds.
Ion	Value varies less than 0.03 mV for the last 8 seconds or 0.08 mV for the last 20 seconds.	Value varies less than 0.08 mV for the last 8 seconds.	Value varies less than 0.3 mV for the last 4 seconds.
Conductivity	Value varies less than 0.4% during the last 8 seconds.	Value varies less than 0.6% during the last 6 seconds.	Value varies less than 0.8% during the last 4 seconds.
Dissolved Oxygen	Value varies less than 0.03 mg/L during the last 20 seconds.	Value varies less than 0.08 mg/L during the last 20 seconds.	Value varies less than 0.08 mg/L during the last 10 seconds.

9 Series

Navigation: Home > Series

Defining and using series support you in performing identical analysis sequences for several samples. Serial measurement can be carried out with METTLER TOLEDO methods or with user-defined methods. Verify that a user-defined method is composed before you set the parameters for a series. Series can be defined in combination with the following method types.

- **Measurement**
- **Interval**

A maximum number of 9 samples per series can be defined, Rondolino can be used to run the series. A maximum of 60 series can be stored in the instrument. You can create shortcuts for series. Series can be created, modified and deleted.

9.1 Creating series

Note

Verify that a method is defined suitable to the series you want to create.

Navigation: Home > Series > New

Parameter	Description	Values
Series ID	According to the headline of the screen, the Series ID is entered automatically starting with S and consecutive number.	-
Method ID	Opens the list of METTLER TOLEDO Methods and user-defined methods.	-
Method type	Information on type of measurement.	Measurement Interval
Number of samples	Defines the number of samples for the series.	-
Default sample ID	Defines an ID for the default sample.	-

- 1 Enter a series ID and confirm with **[OK]**.

Note

If you enter the **Series ID** of an existing series, a pop-up window opens with the information that the series ID already exists!

Tap **[Overwrite]** to use the identical series ID for a new series.

- or -

Tap **[Cancel]** and change the series ID.

- 2 Select a method ID.
- 3 In **Number of samples**, select the number you want to use and confirm with **[OK]**.
- 4 If desired, enter the ID in **Default sample ID**.
- 5 Tap **[Save]**.

➔ Now you have created a series. A dialog with the series ID as title of the newly created series is displayed.

Note

- A maximum of 60 series can be stored at the instrument. If the maximum number of series is reached the **[New]** button is disabled. You have to delete minimum one series before you can create new series

9.2 Creating shortcuts for series

Navigation: Home > Series

Shortcuts can be placed on the **Home screen** and constitute links to executable series. Shortcuts can only be created from the **Start analysis** screen using the **[AddToHome]** button. The administration of shortcuts takes place in the **Setup** menu, see Shortcuts. This chapter describes how a user can create shortcuts to start a series from the home screen.

- 1 Select a series from the list.
 - ➔ **Series ID** is displayed.
- 2 Tap [**Start**].
 - ➔ **Start analysis** is displayed.
- 3 Tap [**AddToHome**].
 - ➔ **Shortcut parameters** is displayed.
- 4 If desired, enter a meaningful **Description**, that will be displayed in the shortcut.
- 5 Select **Immediate start**, to create a direct shortcut.
- 6 Tap the list field **Homescreen position**.
 - ➔ **Select location for shortcut** is displayed.
- 7 Select a free position.
- 8 Tap [**Save**]
 - ➔ The shortcut is displayed on the selected position of the home screen.

Note

With direct shortcuts a series can be started immediately. Shortcuts lead to the **Start analysis** screen from where you can start the serial analysis.

9.3 Modifying series

You can modify series by changing the sdample ID. You can insert or delete samples.

9.3.1 Change single sample ID

Navigation: Home > Series

- 1 Select a series.
 - ➔ **Series ID** is displayed.
- 2 Tap the sample to be modified.
 - ➔ **Series item** is displayed.
- 3 Change the sample ID and confirm with [**OK**].
- 4 Confirm with [**OK**].
- 5 To modify additional sample ID's, repeat previous steps.
- 6 To finish and store entries, tap [**Save**].

9.3.2 Insert samples

Navigation: Home > Series

- 1 Select a series.
 - ➔ **Series ID** is displayed.
- 2 Tap [**Insert**].
 - ➔ Arrow-shaped **Insert** buttons are displayed.
- 3 Tap [**Insert**] at the position where you want to insert one or more samples.
 - ➔ **Series item** is displayed.
- 4 Enter a sample ID and confirm with [**OK**].
- 5 Enter the number to be inserted in **Number of samples** and confirm with [**OK**].
- 6 Confirm with [**OK**].
- 7 To finish and store entries, tap [**Save**].

Note

- A maximum of 9 samples per series can be stored.

9.3.3 Delete samples

Navigation: Home > Series

- 1 Select a series.
 - ➔ **Series ID** is displayed.
- 2 Tap the sample to be deleted.
 - ➔ **Series item** is displayed.
- 3 Tap [**Delete**].
 - ➔ The sample has been deleted without confirmation or warning.
- 4 To delete additional samples, repeat previous steps.

9.4 Delete series

Navigation: Home > Series

- 1 Tap the series to be deleted.
 - ➔ **Series ID** is displayed.
- 2 Tap [**Delete**].
 - ➔ An Information screen opens with the warning that shortcuts which refer to the series will also be deleted.
- 3 Tap [**Delete**].
 - ➔ The series has been deleted.
- 4 To delete additional series, repeat previous steps.

10 Results

Navigation: Home > Results

Results on the **Home screen** opens the analysis list. The results of the most recent 250 analyses consisting of analysis status, analysis data, user-defined calculations, information on the analysis and statistical data of series are found. The performed analyses are listed in chronological order, with the newest analysis at the top of the list. If the maximum number of possible analysis entries is reached, the oldest entry will be deleted.

The analysis list as a whole can be deleted, single analysis can be deleted as well. You can print or transfer data of single entries if the corresponding printer settings are done. Statistics are available only for series.

Note

- Discontinued analyses are not included in the list.

Opening the analysis list leads to the status view automatically. You can switch between the analysis list including data concerning the status of the measurements and the analysis list concerning the results of the measurements. Tapping the [**Results**] button switches immediately to the result view. Return to the status view with **Status**.

See also

- 📖 Peripherals ▶ Page 30
- 📖 Errors in the analysis sequences ▶ Page 57

10.1 Statuses of measurements

Navigation: Home > Results

The status view of the analysis list is displayed and contains the following information:

- **Date**
- **Type**
- **Method/Series ID**
- **Status**

The following designations characterize the type of measurement.

- DM: **Direct measurement**
- DC: **Direct calibration**
- MS: method type **Measurement**
- CAL: method type **Calibrate**
- INC: method type **Incremental**
- INT: method type **Interval**
- S: **Series**
- ST: method type **Sensor Test**
- BCV: Method type BCV
- BOD: Method type BOD

Note

- For **Direct measurement** and **Direct calibration** no method ID is displayed.
- For sample series, a single entry is shown in the analysis list.

Executed analyses can have the following status.

- **OK**
 - Measurement proceeded flawlessly.

- **OK***

Corresponds to the status **OK** with one of the following restrictions:

- Sensor has expired.

- Limits were set and exceeded but **Interrupt outside limits** was disabled.
The task was not interrupted.
- **Error**
 - The user tapped [**Terminate**] and finished the analysis before the end of the method was reached.
The task was interrupted.
- **Failed**
 - Limits were set and exceeded and **Interrupt outside limits** was enabled.
The task was interrupted.

10.2 Results of analysis

Navigation: Home > Results > Results

The result view of the analysis list is displayed and contains the following information:

- **Date**
- **Type**
- **Sample**
- **Result**

Note

- Listed series do not show **Result** or **Sample** in this view.
- For additional information on the results of series, see [Viewing single analysis data ▶ Page 116]

10.3 Statistics

[**Results**] includes statistical calculations based on the measurement results of series.

Navigation: Home > Results > Series > Statistics

Statistics are displayed with the following parameters:

- **Results**, measurement results
- **Mean**, average of all measurement results
- **SD**, standard deviation of all measurement results
- **Min.** , **Max.**, minimum and maximum measurement result

10.4 Deleting all analysis

Navigation: Home > Results

You can delete the complete content of the analysis list.

- Tap [**Delete all**].
- ➔ The analysis list is emptied.

10.5 Deleting single analysis

Navigation: Home > Results

- 1 Tap the relevant analysis on the analysis list.
 - 2 Tap [**Delete**].
 - 3 Tap again [**Delete**] to confirm the deletion.
- ➔ The analysis has been deleted, the entry is deleted from the analysis list.

10.6 Viewing single analysis data

Navigation: Home > Results

Single analysis data or series data can be reviewed.

- 1 Tap the relevant analysis on the analysis list.
- 2 Tap **Measurement data**, **Settings**, **Sample** or **Resources**.
 - ➔ **Measurement data** includes endpoint values of all measurements as well as calculated values including limits. For calibrations it includes the calibration results. Click on any measurement result to see the status, endpoint information, and raw value.
 - ➔ Click on **Data (Interval)** to see the timed interval readings for interval measurement analysis.
 - ➔ **Settings** includes all important method parameters like method ID, temperature capture, calibration mode and so on.
 - ➔ **Sample** includes the sample IDs and comments. For calibrations it includes the buffer/standard group.
 - ➔ **Resources** includes date and time, user name, sensor name, module type and so on.

See also

📖 Peripherals ▶ Page 30

10.7 Printing analysis data

Navigation: Home > Results

Single analysis data or series data can be printed or transferred to an USB stick. The printer used is defined in **Setup > Hardware > Peripherals > Printer**.

- 1 Tap the relevant analysis on the analysis list.
- 2 Tap [**Print**].
- 3 Select the **Print format**.
- 4 If **User-defined** is selected in **Print format**, activate the types of data to be printed.
- 5 Tab [**OK**] to print.

Print format

- **Summary** prints only the most important data.
- **User-defined** provides the selection what data to print.
- **According to method** will re-print exactly the data that was printed during method execution.

Note

- For additional information on the content of print-outs, see method function **Report**.

See also

📖 Report ▶ Page 107

11 Formula Syntax

The formula syntax enables you to formulate criteria and conditions to evaluate the results of your analysis. Formulas can be entered under the method function Calculation or whenever the tickbox Condition is selected within the appropriate method function:

- **Measure**
- **Measure (Incremental)**
- **Measure (Interval)**
- **Measure (Temperature)**
- **Calculation**
- **Wait/Stir**

11.1 Formula symbols

The following symbols are available under the measure function

Symbol	Description
U	May be used for method function Measure , Measure (Incremental) , Calculation and Measure (Interval) . Analysis value at the time when the endpoint is taken, all value corrections are included. The unit depends on the settings in method function Measure , Measure (Interval) , Measure (Temperature) or Measure (Incremental) .
E	May be used for method function Measure , Measure (Interval) and Measure (Temperature) . Analysis value at the time when the endpoint is taken without value corrections. Units are mV for pH, ion and Redox and Ω for conductivity.
UST	May be used for method function Measure , Measure (Interval) and Measure (Temperature) . Analysis values at the start of an analysis, all value corrections are included. The unit depends on the settings in method function Measure or Measure (Interval) .
T	May be used for method function Measure , Measure (Incremental) , Measure (Interval) and Measure (Temperature) . Temperature at the time when the endpoint is taken.
t	May be used for method function Measure . Time period from the start of an analysis until the endpoint is taken. Unit is "s".
P	May be used for method function Measure and Measure (Interval) . Barometric pressure at the time when the endpoint is taken. Only available for the Dissolved Oxygen parameter

The following symbols for raw results are generally available

Symbol	Description
E1-Ex	The number following the result refers to a specific measurement type in method function Configuration .
R1...Rx	Fix marker for all results.
AuxInstr	Results created by method function Auxiliary instrument .
TAB[Tablename()]	Fix marker to calculate tables. The value in round brackets is the input value 'x' of the table, and the completed and calculated fix marker is the output value 'y' of the table.

Square brackets listed in the above tables refer to table names and for indexing the method function **Measure** with a method (e.g. E1[1], E1[2]...).

The following symbols for raw results are available for BOD calculations

Symbol	Description
U	Average BOD value generated after the follow measurement of the last bottle of a sample or check value is taken, including also seed or blank correction. The unit is always mg/L
EB	BOD value of a bottle generated after each follow measurement by subtracting base value and seeded blank or blank. The unit is always mg/L
E	DO value generated after a measurement has been end-pointed. Temperature, salinity and barometric pressure corrections included. The unit is always mg/L.
UST	DO values at the start of an analysis, temperature, salinity and barometric pressure corrections included. The unit is mg/L.
T	Temperature at the time when the endpoint is taken during the DO measurement.
t	Time period from the start of a DO measurement until the endpoint is taken. Unit is seconds.
P	Barometric pressure at the time when the endpoint is taken during the DO measurement.
e	Elapsed days. Time interval between endpoint of base and follow measurement The unit is days. First bottle of base measurement to first bottle follow measurement.

The number following a one of the above listed raw results (U1 = Blank, U2 = seeded Blank, U3 = Standard) refers to a specific check value in MF_Configuration. The number in brackets following the result EB[1], EB[2] refers to the bottle numbers. The index 1 (E[1]1, E[2]1 ...) after the bracket indicates a base measurement, whereas the index 2 (E[1]2, E[2]2 ...) stands for a follow measurement.

The following mathematical operators are available

Symbol	Description
+	Addition
-	Subtraction
*	Multiplication
/	Division

The following comparison operators are available

Symbol	Description
>	Greater than
>=	Greater than or equal to
=	Numerical equality
<=	Less than or equal to
<	Less than
..<..<..	In the range of
<>	Not equal to

The following logical operators are available

Symbol	Description
and	Conjunction
or	Disjunction
true/false	Condition

The following mathematical formulas are available

Symbol	Description
lg()	Logarithm base 10
ln()	Logarithm base e

Symbol	Description
pw()	Exponentials base 10
ex()	Exponential base e
sq()	Square
sr()	Square root

11.2 Creating formulas

Raw results, symbols and mathematical formulas can either be entered directly or taken from the proposal list. If a value is selected from the proposal list, it will be entered at the cursor position.

11.3 Examples

11.3.1 Formula in method function calculation

Formula	Description
U1[2], also possible: R1 = U[2]	Second result of a single channel measurement
E3[1], also possible: R1 = E3	First non-corrected result of third channel
EB1[2]	BOD value in second blank bottle (BCV method type)
P2[3]2	Barometric pressure of third bottle of seeded blank in follow step (BCV method type)
UST1[1]-U1[1] - also possible - UST-U	Difference between initial and final reading of single channel measurement
(T1+T2)/2	Average temperature of channel 1 and 2 in first measurement
First calculation: t[1]+t[2] Second calculation: R1+t[3]	Total measuring time of two and three measurements of single channel measurement
U-EB[5]	Difference between average BOD value and BOD value in fifth sample bottle
U-TAB[Buffer pH 7.00(T)]	Difference between measured and theoretical pH value at measured temperature in pH buffer 7.0
lg(U)	Logarithm (base 10) of result
AuxInstr2	Second result of third method function Auxiliary instrument in the method.

11.3.2 Formula in conditions

The corresponding method function is only executed if formula is true.

Examples for formula in Conditions

Formula	Description
U[3]>=100	Third result is greater or equal to 100
1.0<R1<1.2	Result R1 is between 1.0 and 1.2
T[1]<>T[2], also possible: T-T[2]<>0	First and second temperature are different
e<7	Less than 7 days between base and follow measurement of BOD analysis

11.4 Formula symbols in text

Formula symbols in text In the instruction texts of method functions **Sample verification**, **Instruction** and **Wait/Stir**, formula symbols can be used as placeholders for values. Each symbol has to be bracketed with the % signs.

Example: %U% is used to insert the results of the first measurement in the method. %R1% will show the result of the first calculation. %t[2]% will show the duration of the second measurement.

Beside the formula symbols above, the following placeholders can be used in method functions **Sample verification** and **Instruction**.

Examples for formula in Conditions

Symbol	Description
MethodID	Method ID of the running method.
SampleID	Sample ID from either method function Sample, the start analysis screen or the series parameters.
Comment	Sample comments from either method function Sample or the start analysis screen.

12 News, Tasks and Online Display

If an analysis is started, the online display is shown. If a task is already running and a new analysis is started, the task list with the waiting tasks can be displayed. Tapping the **[News]** button provides information concerning the peripheral devices connected to the instrument.

12.1 News

In case of disconnected PnP resources, exceeded calibration limits etc. the **[News]** button opens a screen with additional information about these issues.

The **[News]** button is located in the top left-hand corner of the **Home screen**. The color of the symbol indicates if the list is empty (grey) or not (white).

Tap the **[News]** button to see the content of news. The list contains entries about connected and disconnected external devices, sensors with expired life span and if the radio clock synchronization process is finished. All information is listed in chronological order.

When a new device is added, the life span of a sensor exceeds or the radio clock finished synchronization, the **[News]** button flashes.

Note

- The maximum number of news entries in the list is 30. If the capacity of the news list is fully exploited, the oldest message will be overwritten. If you switch the instrument off and on again, all old news will be cleared.
- The news list can be cleared manually with **[Remove all]**.

12.2 Tasks

Each execution of a calibration, measurement, method or series is a task. Tasks are always processed in sequence. Each task is listed in the task list and receives a number, based on the chronological order of its start.

For tasks, the following rules must be observed:

- It is allowed to start several identical measurements.
- Methods can also be started when a calibration or direct measurement is running. They will be queued.
- Starting a direct measurement in which more than one module is involved will be listed as one task in the task list.
- It is not possible to start several direct calibrations via **Calibrate**.
- It is not possible to start several direct measurements via **Read**.

Tasks is located in the top right-hand corner of the **Home screen**. It is activated, as soon as at least one task is lined up.

By tapping **Tasks** you switch to the **Online screen** or, if several tasks were started, the task list.

Status display of the **Tasks** button:

Blue	No task is lined up.
Yellow	A task is running right now.
Yellow/blue blinking	A task is waiting for user interaction.
Orange	The task list is interrupted and no task is running.

Note

- Direct calibration via **[Calibrate]** and direct measurement via **[Read]** can only be started when no other task is running. In case a task is running, **[Read]** and **[Calibrate]** are deactivated.

The task list offers the following options:

Show **Online display** By tapping the list entry of the ongoing task, the online display is shown.

[Interrupt]	By tapping the button Interrupt , the processing of the tasks is interrupted. The ongoing task is carried out to completion.
[Resume]	By tapping Resume you can continue the processing of the task.
[Remove all]	This button is only visible, when the task list is interrupted and no tasks are running. By tapping this button all lined up tasks will be removed.
Move tasks	Tap a task. By changing the number, you can move the task.
Delete individual tasks	Tap a task. By tapping Remove the task is deleted from the list.

12.3 Online display

The online display shows the status and readings of the currently running analysis. The data is updated every 0.5 seconds, for conductivity and dissolved oxygen every second.

The **Title bar** at the top of the online screen shows:

- **Method ID** and name
 - If you started a direct measurement: **DM: Direct Measure**
 - If you started a direct calibration: **DC: Direct Calibration**
- **News** for opening a screen including information about connected and disconnected devices and exceeded life span.
- **Tasks** for showing an ongoing task and to open a screen including the task list.

The **Status bar**, located below the **Title bar**, shows the following information:

- Three squares indicating the used modules.
- Name of the method function that is running.
- A timer to indicate how long the method function is running.

A maximum of three **Data fields** can be shown. The number depends on the number of used modules. If two or three Data fields are shown, **Switch view** allows to toggle between full information and large digits.

Note

- Depending on the analysis different parameters are shown in the **Data field**, the commands on the button below the data field vary.

13 Maintenance and Care

Do not open the housing of the instrument; it does not contain any parts that can be maintained, repaired or replaced by the user. If you experience problems with your instrument, contact your authorized METTLER TOLEDO dealer or service representative.

► www.mt.com/contact

13.1 Cleaning the Instrument



NOTICE

Danger of damage to the instrument due to inappropriate cleaning agents!

The housing is made of acrylonitrile butadiene styrene/polycarbonate (ABS/PC). This material is sensitive to some organic solvents, such as toluene, xylene and methyl ethyl ketone (MEK). If liquids enter the housing they can damage the instrument.

- 1 Use only water and a mild detergent to clean the housing.
- 2 Wipe off any spills immediately.
- 3 The instrument is IP54 splash water proof: Do not immerse the instrument in liquid.

- The instrument is turned off and disconnected from the electrical outlet.
- Clean the housing of the instrument using a cloth dampened with water and a mild detergent.

13.2 Maintenance of electrodes

The instrument monitors the condition of the attached pH electrodes.



Slope: 95-105%
and offset: \pm (0-20) mV
Electrode is in good condition



Slope: 90-94%
or offset: \pm (20-35) mV
Electrode needs cleaning



Slope: 85-89%
or offset: \pm (>35) mV
Electrode is defective or too old

When cleaning, always follow the instructions in the manual of the electrodes used. Make sure the pH electrode is always kept filled with the appropriate filling solution. For maximum accuracy, any filling solution that may have "crept" and encrusted the outside of the electrode should be removed with deionized water. Always store the electrode according to the manufacturer's instructions and do not allow it to dry out.

If the electrode slope falls rapidly, or if the response becomes sluggish, the following procedures may help. Try one of the following, depending on your sample.

Problem	Action
Fat or oil build-up	Either rinse the membrane with soap solution or acetone/ethanol or shortly soak the tip of the electrode in hot water. When rinsed with organic solvent, place the membrane overnight in 0.1 mol/L HCl.
pH electrode membrane has dried out	Soak the tip of the electrode overnight in 0.1 mol/L HCl. If this procedure has no effect, soak the tip of the electrode for a few minutes in reactivation solution for pH electrodes.
Protein build-up in the diaphragm of a pH electrode	Remove deposits by soaking the electrode in an HCl/pepsin solution.
Silver sulfide contamination of pH electrode	Remove deposits by soaking the electrode in a thiourea solution.

Run a new calibration after treatment.

Note

- Cleaning and filling solutions should be handled with the same care as that given to toxic or corrosive substances.

- The condition of the pH electrode can also be checked, using the provided METTLER TOLEDO method Sensor Test.

13.3 Transporting the instrument

Note the following instructions when transporting the instrument to a new location:

- Transport the instrument with care to avoid damage! The instrument may be damaged if not transported correctly.
- Unplug the instrument and remove all connected cables.
- Remove the electrode arm.
- To avoid damage to the instrument when transporting it over long distances, please use the original packaging.
- If the original packaging is no longer available, choose packaging that will ensure safe handling.

13.4 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.



14 Accessories

Accessories	Order No.
SevenExcellence™ pH/mV module	30034472
SevenExcellence™ conductivity module	30034473
SevenExcellence™ pH/Ion module	30034471
SevenExcellence™ DO/BOD module	30034474
SevenExcellence™ blank module	30034475
SevenExcellence™ protective film (2 pcs.)	30041155
uPlace™ electrode arm, complete with base and extension	30019823
SevenExcellence™ semi-transparent cover	30041154
InMotion sample changer bundle: Flex base and 100 mL rack (18 samples)	30094120
Rondolino sample changer, complete for SevenExcellence™	51302888
Rondolino USB-TTL box	30046261
Magnetic stirrer uMix™	30040002
Compact Stirrer Kit (stirrer incl. 2 propeller stirrer and adapter cable)	30115728
Compact Stirrer adapter cable	30098212
Disposable beakers for InMotion and Rondolino (PP, 100 mL, 1400 pcs.)	00101974
PowerShower™	51108219
LogStraight™ fingerprint reader	51192107
Barcode reader	21901297
USB cable for barcode reader	21901309
Printer USB-P25	11124301
EasyDirect pH PC software	-
USB A-B cable 1.8 m for EasyDirect pH software (not included with EasyDirect pH software)	51191926
Adapter for auxiliary instruments (USB-RS232 adapter)	51105856

Buffers and standards	Order No.
pH 4.01 buffer sachets, 30 x 20mL	51302069
pH 4.01 buffer solution, 6 x 250mL	51350018
pH 7.00 buffer sachets, 30 x 20mL	51302047
pH 7.00 buffer solution, 6 x 250mL	51350020
pH 9.21 buffer sachets, 30 x 20mL	51302070
pH 9.21 buffer solution, 6 x 250mL	51350022
pH 10.01 buffer sachets, 30 x 20mL	51302079
pH 10.00 buffer solution, 6 x 250mL	51350024
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/9.21)	51302068
Rainbow (box of each, 10 x 20 mL, 4.01/7.00/10.01)	51302080
10 µS/cm conductivity standard solution, 250 mL	51300169
84 µS/cm conductivity standard solution, 250 mL	51302153
500 µS/cm conductivity standard solution, 250 mL	51300170
1413 µS/cm conductivity standard solution, 30 x 20 mL	51302049
1413 µS/cm conductivity standard solution, 6 x 250 mL	51350096
12.88 mS/cm conductivity standard solution, 30 x 20 mL	51302050
12.88 mS/cm conductivity standard solution, 6 x 250 mL	51350098
Zero oxygen tablets (24 pcs.)	51300140
pH sensors with fixed cable	Order No.
InLab® Expert Pro-ISM, rugged 3-in-1 pH sensor, PEEK shaft, ATC	30014096
ISM® sensors with MultiPin™ head	Order No.
InLab® Routine Pro-ISM, 3-in-1 pH sensor, glass shaft, ATC, refillable	51344055
InLab® Micro Pro-ISM, 3-in-1 pH sensor, glass shaft, 5 mm shaft diameter, ATC, refillable	51344163
InLab® Power Pro-ISM, 3-in-1 pH sensor, glass shaft, ATC, pressurized SteadyForce™ reference system	51344211
InLab® Pure Pro-ISM, 3-in-1 pH sensor, glass shaft, immovable glass sleeve, ATC, refillable	51344172
InLab® Science Pro-ISM, 3-in-1 pH sensor, glass shaft, movable glass sleeve, ATC, refillable	51344072
InLab® Solids Pro-ISM, 3-in-1 pH sensor, glass shaft, open junction, sharp membrane, ATC	51344155
Electrode cable 1.2 m, BNC/RCA-MultiPin™	30281896
Electrode cable 3 m, BNC/RCA-MultiPin™	30281897
ATC probe, temperature sensor	12997876
Digital ISM® pH sensors with diagnostic functions	Order No.
InLab® Smart Pro-ISM, 3-in-1 pH sensor, glass shaft, ATC, pressurized SteadyForce™ reference system	30027775
Cable InLab® Digital, 1.2 m	30027776

Conductivity sensors with fixed cable	Order No.
InLab® 731-ISM, 4-electrode graphite conductivity sensor, ATC	30014092
InLab® 741-ISM, 2-electrode steel conductivity sensor, ATC	30014094
InLab® Trace, high-accuracy sensor for low conductivities, ATC	30014097
Flow-Cell for InLab® Trace	30014098
InLab® Trace Kit (Sensor and Flow-Cell)	30014099
Polarographic Dissolved Oxygen sensors with fixed cable	Order No.
InLab® 605-ISM-2 m	51344611
InLab® 605-ISM-5 m	51344612
InLab® 605-ISM-10 m	51344613
Optical Dissolved Oxygen sensors with fixed cable	Order No.
InLab® OptiOx, 1.8 m	51344621
InLab® OptiOx, 5 m	51344622
InLab® OptiOx, 10 m	51344623
OptiOx parts	Order No.
OptiOx replacement cap	51344630
OptiOx calibration tube	51344631
OptiOx protective guard	51344632
OptiOx BOD adapter	51344633
Solutions	Order No.
HCl/Pepsin solution (removes protein contamination), 1 x 250 mL	51350100
Thiourea solution (remove silver sulfide contamination), 1 x 250 mL	51350102
Reactivation solution for pH electrodes, 1 x 25 mL	51350104
InLab Storage solution, 1 x 250 mL	30111142

15 Technical Data

15.1 SevenExcellence™

Screen	Color TFT	
Interfaces	RS232, USB A, USB B, Ethernet	
Ambient conditions	Ambient temperature	5 °C to 40 °C
	Relative humidity	5% to 80% (non-condensing)
	Overtoltage category	Class II
	Pollution degree	2
	Range of application	For indoor use only
	Maximum operating altitude	Up to 2000 m
Dimensions	Width	235 mm
	Depth	188 mm
	Height	75 mm
Weight	Basic device	1120 g
	1 module	111 - 130 g
Power rating instrument	Input voltage	12 V $\overline{=}$
	Power consumption	10 W
Power rating AC adapter	Line voltage	100 - 240 V $\sim \pm 10\%$
	Input frequency	50/60 Hz
	Input current	0.3 A
	Output voltage	12 V $\overline{=}$
	Output current	0.84 A
Materials	Housing	ABS/PC
	Expansion units	ABS/PC
	Electrode arm	ABS/PC
	Touch panel	Tempered glass

15.2 pH/mV module

Sensor inputs	Digi	Mini-LTW digital sensors	
	ISFET	Mini-DIN ISFET sensor	
	ATC int	RCA (Cinch) NTC30k	
	pH	BNC mV/pH sensors, Impedance $>3 \cdot 10^{12} \Omega$	
pH mode			
	pH	mV	Temperature °C
Measurement range	-2.000...20.000	± 2000.0	-30.0...130.0
Resolution	0.001/0.01/0.1	0.1	0.1
Limits of error	± 0.002	± 0.1 mV (-1000...+1000 mV) ± 0.2 mV ($> \pm 1000$ mV)	0.0...100.0 °C: ± 0.1 -30.0...0.0 °C: ± 0.3 100.0...130.0 °C: ± 0.3
Relative mV	-	Yes	-
Temperature compensation	Automatic	-30.0 °C...130.0 °C	
	Manual	-30.0 °C...130.0 °C	
ISFET mode			
	pH	Temperature °C	
Measurement range	0.000...14.000	-30.0...130.0	
Resolution	0.001/0.01/0.1 pH	0.1	
Limits of errors	± 0.05 pH	0.0...100.0 °C: ± 0.1 -30.0...0.0 °C: ± 0.3 100.0...130.0 °C: ± 0.3	

15.3 Conductivity module

Sensor inputs Digi	Mini-LTW digital sensors	
ATC ext	RCA (Cinch) NTC30k or PT1000	
Cond	Mini-DIN conductivity sensors	
Conductivity mode	Measurement range	
	0.001...999999 $\mu\text{S}/\text{cm}$	
	0.001...2000 mS/cm	
	0.001...200 S/m	
	0.001...200000 mS/m	
	0.1...999999 $\mu\text{S}/\text{m}$	
Resolution	0.001...1 $\mu\text{S}/\text{cm}$, mS/cm	
	0.0001...0.1 S/m , $\mu\text{S}/\text{m}$, mS/m	
Limits of errors	$\pm 0.5\%$ of measured value	
TDS mode	Measurement range	
	0.001...1000 ppt, g/L	
	0.001...999999 mg/L, ppm	
	TDS Factor	0.00...10.00
	Resolution	0.0001...1 ppt, g/L
	0.001...1 mg/L, ppm	
Limits of errors	$\pm 0.5\%$ of measured value	
Salinity mode	Measurement range	0.0...80 psu, ppt
	Resolution	0.01...0.1 psu, ppt
	Limits of errors	$\pm 0.5\%$ of measured value
Resistivity mode	Measurement range	
	0.01...100 $\text{M}\Omega\cdot\text{cm}$	
	0.01...999999 $\Omega\cdot\text{cm}$	
	Resolution	0.0001...1 $\text{M}\Omega\cdot\text{cm}$
		0.01...1 $\Omega\cdot\text{cm}$
Limits of errors	$\pm 0.5\%$ of measured value	
Temperature capture	Temperature compensation	Automatic: -30 °C...130 °C
		Manual: -30 °C...130 °C
	Limits of errors	0.0...100.0 °C: ± 0.1 -30.0...0.0 °C: ± 0.3 100.0...130.0 °C: ± 0.3

15.4 pH/Ion module

Senor inputs	Digi	Mini-LTW digital sensors	
	ATC ext	RCA (Cinch) NTC30k or PT1000	
	ref	Reference electrode	
	ATC int	RCA (Cinch) NTC30k	
	pH	BNC mV/pH sensors, Impedance $>3 \cdot 10^{12} \Omega$	
Ion mode			
	Ion	Temperature °C	
Measurement range	0...999 999 mg/L, ppm	-30.0...130.0	
	0...100 mol/L, %		
	0...100000 mmol/L		
	-2.000...20.000 pX		
Resolution	0.001...1 mg/L, ppm, mmol/L, pX	0.1	
	0.0001...100 mol/L, %		
Limits of error	± 0.5 % of measured value	0.0...100.0 °C: ±0.1	
		-30.0...0.0 °C: ±0.3	
		100.0...130.0 °C: ±0.3	
pH mode			
	pH	mV	Temperature °C
Measurement range	-2.000...20.000	± 2000.0	-30.0...130.0
Resolution	0.001/0.01/0.1	0.1	0.1
Limits of error	± 0.002	± 0.1	0.0...100.0 °C: ±0.1
			-30.0...0.0 °C: ±0.3
			100.0...130.0 °C: ±0.3
Relative mV	-	yes	-
Temperature compensation	Automatic	-30.0 °C...130.0 °C	
	Manual	-30.0 °C...130.0 °C	

15.5 DO/BOD module

Sensor inputs	Digi	Mini-LTW digital sensors
	ATC ext	RCA (Cinch) NTC30k or PT1000
	ATC int	RCA (Cinch) NTC22k
	DO	BNC DO sensors
Dissolved oxygen		
	Optical sensor (digital)	Polarographic sensor (analog)
Measuring range	0.000...50 mg/L (ppm)	0.000...99 mg/L (ppm)
Resolution	0.001/0.01/0.1	0.001/0.01/0.1
Limits of error	± 0.1 mg/L from 0...8 ± 0.2 mg/L from 8...20 ± 10% from 20...50	± 0.5%
Units	mg/L, ppm	mg/L, ppm
DO Saturation		
	Optical sensor (digital)	Polarographic sensor (analog)
Measuring range	0.0...500%	0.0...600%
Resolution	0.0001/0.001/0.01/0.1	0.0001/0.001/0.01/0.1
Temperature		
	Optical sensor (digital)	Polarographic sensor (analog)
Measuring range	0.0 °C...50.0 °C	0.0 °C...60.0 °C
Resolution	0.1 °C	0.1 °C
Limits of error	± 0.1 °C	± 0.1 °C
Pressure		
	Optical sensor (digital)	Polarographic sensor (analog)
Measuring range	500 to 1100 mbar	500 to 1100 mbar
Resolution	1	1
Limits of error	± 1 mbar	± 1 mbar
General		
	Optical sensor (digital)	Polarographic sensor (analog)
Calibration points	2	2
Pressure compensation	automatic / manual	automatic / manual
Pressure units	mbar, hPa, mmHg, atm	mbar, hPa, mmHg, atm

16 Appendix

16.1 Predefined buffers & standards

Buffer Sets

METTLER TOLEDO Europe (Ref. 25°C)

T [°C]	2.00	4.01	7.00	9.21	11.00
0.0	2.03	4.01	7.12	9.52	11.90
5.0	2.02	4.01	7.09	9.45	11.72
10.0	2.01	4.00	7.06	9.38	11.54
15.0	2.00	4.00	7.04	9.32	11.36
20.0	2.00	4.00	7.02	9.26	11.18
25.0	2.00	4.01	7.00	9.21	11.00
30.0	1.99	4.01	6.99	9.16	10.82
35.0	1.99	4.02	6.98	9.11	10.64
40.0	1.98	4.03	6.97	9.06	10.46
45.0	1.98	4.04	6.97	9.03	10.28
50.0	1.98	4.06	6.97	8.99	10.10
55.0	1.98	4.08	6.98	8.96	-
60.0	1.98	4.10	6.98	8.93	-
65.0	1.98	4.13	6.99	-	-
70.0	1.99	4.16	7.00	8.88	-
75.0	1.99	4.19	7.02	-	-
80.0	2.00	4.22	7.04	8.83	-
85.0	2.00	4.26	7.06	-	-
90.0	2.00	4.30	7.09	8.79	-
95.0	2.00	4.35	7.12	8.77	-

METTLER TOLEDO USA (Ref. 25°C)

T [°C]	1.68	4.01	7.00	10.01
0.0	1.67	4.01	7.12	10.32
5.0	1.67	4.01	7.09	10.25
10.0	1.67	4.00	7.06	10.18
15.0	1.67	4.00	7.04	10.12
20.0	1.68	4.00	7.02	10.06
25.0	1.68	4.01	7.00	10.01
30.0	1.68	4.01	6.99	9.97
35.0	1.69	4.02	6.98	9.93
40.0	1.69	4.03	6.97	9.89
45.0	1.70	4.04	6.97	9.86
50.0	1.71	4.06	6.97	9.83
55.0	1.72	4.08	6.98	-
60.0	1.72	4.10	6.98	-
65.0	-	4.13	6.99	-
70.0	1.74	4.16	7.00	-
75.0	-	4.19	7.02	-
80.0	1.77	4.22	7.04	-
85.0	-	4.26	7.06	-
90.0	1.79	4.30	7.09	-
95.0	1.81	4.35	7.12	-

METTLER TOLEDO Verification buffers (Ref. 25°C)

T [°C]	5.00	8.00
0.0	5.04	8.07
5.0	5.03	8.06
10.0	5.02	8.07
15.0	5.01	8.04
20.0	5.00	8.02
25.0	5.00	8.00
30.0	5.01	7.98
35.0	5.01	7.95
40.0	5.03	7.94
45.0	5.05	7.91
50.0	5.06	7.90
55.0	5.08	7.89
60.0	5.11	7.86
65.0	5.14	7.88
70.0	5.17	7.87
75.0	5.20	7.86
80.0	5.23	7.85
85.0	5.26	7.86
90.0	5.29	7.87
95.0	5.32	7.87

MERCK (Ref. 20°C)

T [°C]	2.00	4.00	7.00	9.00	12.00
5.0	2.01	4.05	7.07	9.16	12.41
10.0	2.01	4.03	7.05	9.11	12.26
15.0	2.00	4.02	7.02	9.05	12.10
20.0	2.00	4.00	7.00	9.00	12.00
25.0	2.00	3.99	6.98	8.95	11.88
30.0	2.00	3.98	6.98	8.91	11.72
35.0	2.00	3.98	6.96	8.88	11.67
40.0	2.00	3.98	6.95	8.85	11.54
45.0	2.00	3.98	6.95	8.82	11.44
50.0	2.00	3.98	6.95	8.79	11.33

DIN(19266)/NIST (Ref. 25°C)

T [°C]	1.679	4.005	6.865	9.180	12.454
0.0	1.666	4.000	6.984	9.464	-
5.0	1.668	3.998	6.951	9.395	13.207
10.0	1.670	3.997	6.923	9.332	13.003
15.0	1.672	3.998	6.900	9.276	12.810
20.0	1.675	4.000	6.881	9.225	12.627
25.0	1.679	4.005	6.865	9.180	12.454
30.0	1.683	4.011	6.853	9.139	12.289
35.0	1.688	4.018	6.844	9.102	12.133
37.0	-	4.022	6.841	9.088	-
38.0	1.691	-	-	-	12.043
40.0	1.694	4.027	6.838	9.068	11.984
45.0	-	-	-	-	11.841
50.0	1.707	4.050	6.833	9.011	11.705

DIN(19267) (Ref. 25°C)

T [°C]	1.09	4.65	6.79	9.23	12.75
0.0	1.08	4.67	6.86	9.48	-
10.0	1.09	4.66	6.84	9.37	13.37
20.0	1.09	4.65	6.80	9.27	12.96
25.0	1.09	4.65	6.79	9.23	12.75
30.0	1.10	4.65	6.78	9.18	12.61
40.0	1.10	4.66	6.76	9.09	12.29
50.0	1.11	4.68	6.76	9.00	11.98
60.0	1.11	4.70	6.76	8.92	11.69
70.0	1.11	4.72	6.76	8.88	11.43
80.0	1.12	4.75	6.78	8.85	11.19
90.0	1.13	4.79	6.80	8.82	10.99

JJG119 (Chinese) (Ref. 25°C)

T [°C]	1.680	4.003	6.864	9.182	12.460
0.0	1.668	4.006	6.981	-	13.416
5.0	1.669	3.999	6.949	9.391	13.210
10.0	1.671	3.996	6.921	9.330	13.011
15.0	1.673	3.996	6.898	9.276	12.820
20.0	1.676	3.998	6.879	9.226	12.637
25.0	1.680	4.003	6.864	9.182	12.460
30.0	1.684	4.010	6.852	9.142	12.292
35.0	1.688	4.019	6.844	9.105	12.130
37.0	1.694	4.022	6.839	-	12.069
40.0	1.694	4.029	6.838	9.072	11.975
45.0	1.700	4.042	6.834	9.042	11.828
50.0	1.706	4.055	6.833	9.015	11.697
55.0	1.713	4.070	6.834	8.990	11.553
60.0	1.721	4.087	6.837	8.968	11.426
70.0	1.739	4.122	6.847	8.926	-
80.0	1.759	4.161	6.862	8.890	-
90.0	1.782	4.203	6.881	8.856	-
95.0	1.795	4.224	6.891	8.839	-

Technical (Ref. 25°C)

T [°C]	2.00	4.01	7.00	10.00
0.0	2.03	4.01	7.12	-
5.0	2.02	4.01	7.09	10.65
10.0	2.01	4.00	7.06	10.39
15.0	2.00	4.00	7.04	10.26
20.0	2.00	4.00	7.02	10.13
25.0	2.00	4.01	7.00	10.00
30.0	1.99	4.01	6.99	9.87
35.0	1.99	4.02	6.98	9.74
40.0	1.98	4.03	6.97	9.61
45.0	1.98	4.04	6.97	9.48
50.0	1.98	4.06	6.97	9.35
55.0	1.98	4.08	6.98	-
60.0	1.98	4.10	6.98	-
65.0	1.98	4.13	6.99	-
70.0	1.99	4.16	7.00	-
75.0	1.99	4.19	7.02	-
80.0	2.00	4.22	7.04	-
85.0	2.00	4.26	7.06	-
90.0	2.00	4.30	7.09	-
95.0	2.00	4.35	7.12	-

JIS Z 8802 (Japanese) (Ref. 25°C)

T [°C]	1.679	4.008	6.865	9.180
5.0	1.668	3.999	6.951	9.395
10.0	1.670	3.998	6.923	9.332
15.0	1.672	3.999	6.900	9.276
20.0	1.675	4.002	6.881	9.225
25.0	1.679	4.008	6.865	9.180
30.0	1.683	4.015	6.853	9.139
35.0	1.688	4.024	6.844	9.102
40.0	1.694	4.035	6.838	9.068
45.0	1.700	4.047	6.834	9.038
50.0	1.707	4.060	6.833	9.011

Standard Groups**International (Ref. 25°C)**

T [°C]	10 µS/cm	84 µS/cm	500 µS/cm	1413 µS/cm	12.88 mS/cm
5.0	6.13	53.02	315.3	896	8.22
10.0	7.10	60.34	359.6	1020	9.33
15.0	7.95	67.61	402.9	1147	10.48
16.0	8.15	69.25	-	-	-
17.0	8.36	70.89	-	-	-
18.0	8.56	72.52	-	-	-
19.0	8.77	74.16	-	-	-
20.0	8.97	75.80	451.5	1278	11.67
21.0	9.18	77.44	-	-	-
22.0	9.38	79.08	-	-	-
23.0	9.59	80.72	-	-	-
24.0	9.79	82.36	-	-	-
25.0	10.00	84.00	500.0	1413	12.88
30.0	11.03	92.19	548.5	1552	14.12
35.0	12.14	100.92	602.5	1696	15.39
40.0	13.29	109.21	-	-	-
45.0	14.44	118.05	-	-	-
50.0	15.55	126.80	-	-	-

Chinese (Ref. 25°C)

T [°C]	146.5 µS/cm	1408 µS/cm	12.85 mS/cm	111.3 mS/cm
15.0	118.5	1141.4	10.455	92.12
18.0	126.7	1220.0	11.163	97.80
20.0	132.2	1273.7	11.644	101.70
25.0	146.5	1408.3	12.852	111.31
35.0	176.5	1687.6	15.353	131.10

Japanese Standards (Ref. 20°C)

T [°C]	1330 µS/cm	133.0 µS/cm	26.6 µS/cm
0.0	771.40	77.14	15.428
5.0	911.05	91.11	18.221
10.0	1050.70	105.07	21.014
15.0	1190.35	119.04	23.807
20.0	1330.00	133.00	26.600
25.0	1469.65	146.97	29.393
30.0	1609.30	160.93	32.186
35.0	1748.95	174.90	34.979

Saturated NaCl (Ref. 25°C)

T [°C]	251.3 mS/cm
5.0	155.5
10.0	177.9
15.0	201.5
20.0	226.0
25.0	251.3
30.0	277.4
35.0	304.1

METTLER TOLEDO ION (Ref. 25°C)

T [°C]	0.1 mg/L	1 mg/L	10 mg/L	100 mg/L	1000 mg/L
20.0	0.1	1	10	100	1000
25.0	0.1	1	10	100	1000
30.0	0.1	1	10	100	1000

16.2 METTLER TOLEDO tables

METTLER TOLEDO Europe (Ref. 25°C)

T [°C]	2.00	4.01	7.00	9.21	11.00
0.0	2.03	4.01	7.12	9.52	11.90
5.0	2.02	4.01	7.09	9.45	11.72
10.0	2.01	4.00	7.06	9.38	11.54
15.0	2.00	4.00	7.04	9.32	11.36
20.0	2.00	4.00	7.02	9.26	11.18
25.0	2.00	4.01	7.00	9.21	11.00
30.0	1.99	4.01	6.99	9.16	10.82
35.0	1.99	4.02	6.98	9.11	10.64
40.0	1.98	4.03	6.97	9.06	10.46
45.0	1.98	4.04	6.97	9.03	10.28
50.0	1.98	4.06	6.97	8.99	10.10
55.0	1.98	4.08	6.98	8.96	-
60.0	1.98	4.10	6.98	8.93	-
65.0	1.98	4.13	6.99	-	-
70.0	1.99	4.16	7.00	8.88	-
75.0	1.99	4.19	7.02	-	-
80.0	2.00	4.22	7.04	8.83	-
85.0	2.00	4.26	7.06	-	-
90.0	2.00	4.30	7.09	8.79	-
95.0	2.00	4.35	7.12	8.77	-

Purified Water (EP/Ch.P.)

T [°C]	µS/cm	T [°C]	µS/cm
0	2.4	60	8.1
10	3.6	70	9.1
20	4.3	75	9.7
25	5.1	80	9.7
30	5.4	90	9.7
40	6.5	100	10.2
50	7.1		

USP645 Stage 1

T [°C]	µS/cm	T [°C]	µS/cm
0	0.6	54.99	1.9
4.999	0.6	55	2.1
5	0.8	59.99	2.1
9.999	0.8	60	2.2
10	0.9	64.99	2.2
14.99	0.9	64	2.4
15	1.0	69.99	2.4
19.99	1.0	70	2.5
20	1.1	74.99	2.5
24.99	1.1	75	2.7
25	1.3	79.99	2.7
29.99	1.3	80	2.7
30	1.4	84.99	2.7
34.99	1.4	85	2.7
35	1.5	89.99	2.7
39.99	1.5	90	2.7
40	1.7	94.99	2.7
44.99	1.7	95	2.9
45	1.8	99.99	2.9
49.99	1.8	100	3.1
50	1.9		

Water for Injection Stage 1

T [°C]	µS/cm	T [°C]	µS/cm
0	0.6	54.99	1.9
4.999	0.6	55	2.1
5	0.8	59.99	2.1
9.999	0.8	60	2.2
10	0.9	64.99	2.2
14.99	0.9	64	2.4
15	1.0	69.99	2.4
19.99	1.0	70	2.5
20	1.1	74.99	2.5
24.99	1.1	75	2.7
25	1.3	79.99	2.7
29.99	1.3	80	2.7
30	1.4	84.99	2.7
34.99	1.4	85	2.7
35	1.5	89.99	2.7
39.99	1.5	90	2.7
40	1.7	94.99	2.7
44.99	1.7	95	2.9
45	1.8	99.99	2.9
49.99	1.8	100	3.1
50	1.9		

Water for Injection Stage 3

T [°C]	µS/cm	T [°C]	µS/cm
5.0	4.7	6.1	2.4
5.1	4.1	6.2	2.5
5.2	3.6	6.3	2.4
5.3	3.3	6.4	2.3
5.4	3.0	6.5	2.2
5.5	2.8	6.6	2.1
5.6	2.6	6.7	2.6
5.7	2.5	6.8	3.1
5.8	2.4	6.9	3.8
5.9	2.4	7.0	4.6
6.0	2.4		

Index

A

Account policies	40
Actions when sensors expire	43
AddToHome	62, 111
Analysis	
Start	50
Terminate	57
Analysis and resources behaviour	
Actions when sensors expire	43
Analysis sequence settings	43
Analysis data	
Delete	115
Print	116
View	116
Analysis list	114
Delete	115
Analysis sequence settings	43
Analysis sequences	
Errors	57
Automation	20, 30

B

Barcode reader	32
Barometric pressure	50, 71
Beep	35
Buffers and standards	
Create sets and groups	23
Delete sets and groups	26
Delete values	26
Modify name	26
Modify sets and groups	25
Predefined	22
Select value, conductivity	25
Select value, ion	25
Select value, pH	24
Select values predefined buffer sets	24
Setup	22
User-defined	23
Button	
AddToHome	62, 111
Calibrate	50
News	121
Read	50
Tasks	121

C

Calibration	
Template calibration	86
Calibration history	49
Print	29
Show	29
Calibration settings	51
Condition	
Formula syntax	117
Method functions	117
Critical error	57

D

Data	
Analysis	115
Import/Export Buffers and standards	47
Import/Export Log file	47
Import/Export Memory copy	47
Import/Export Methods	47
Import/Export Series	47
Import/Export User management	47
Data fields	122
Date/Time	38
Direct calibration	50
Calibration settings	51
General settings	51
Module settings	50
Direct measurement	53
General settings	53
Measurement settings	54
Module settings	53
Direct shortcuts	35
Create, methods	62
Create, series	111
Disposal	124
Dissolved oxygen	50, 71

E

end user license agreement	5
Error	57
Errors	57

F

Failed	58
Fingerprint reader	31

Firmware	48	Malfunction type	
end user license agreement	5	Critical error	57
open source attribution files	5	Error	57
third party licenses	5	Failed	58
Formula syntax	117	Terminate error	57
Available symbols	117	Measurement	
BOD calculations	118	Results	115
Comparison operators	118	Status	114
General	117	Measurement settings	54
Logical operators	118	Memory copy	47
Mathematical formulas	118	Method types	59
Mathematical operators	118	Methods	59
G		Create	61
General settings	51, 53	Create, with method templates	61
Global settings		Create, with METTLER TOLEDO methods	61
Analysis and resources behaviour	42	Delete	69
Physical properties	43	METTLER TOLEDO methods	60
Setup	36	METTLER TOLEDO methods, use	60
System	36	Modify	62
H		Setup	59
Hardware		Template Measure	71
Setup	26	METTLER TOLEDO tables	44
Header and footer	38	Module settings	50, 53
I		MT service	47
Identification	37	Service life, set up	47
Incremental		N	
Template incremental	91	Navigation bar	122
InMotion	30	Network settings	33
InMotion autosampler	20	News	121
K		News button	121
Keyboards	36	O	
L		Online display	122
LabX	21, 33	Data fields	122
Language	34	Navigation bar	122
M		Title bar	122
Maintenance & service settings	22	open source attribution files	5
Maintenance and service	47	Operation mode	44
Export, Data	47	P	
Firmware	48	Peripherals	
Import, Data	47	Barcode reader	32
MT Service	47	Fingerprint reader	31
Reset to factory settings	48	Printer	31
Update	48	Setup	30
		Stirrer	32
		USB stick	31

Physical properties	43	Hardware	22
Plug and Play		Maintenance and service	22, 47
Rondolino	30	Methods	59
Stirrer	32	Network settings	33
USB stick	31	Series	111
Pressure	50, 71	Tables	22, 44
Printer	31	User settings	22
R		Shortcuts	35
Registration screen	41	Create, methods	62
Reset to factory settings	48	Create, series	111
Results	114	Delete	36
Measurement results	115	Modify	35
Measurement, status	114	Start analysis	
Rondolino	30	Calibrate	50
Rondolino sample changer	21	Direct calibration	50
S		Direct measurement	53
safety information	6	From home screen	56
Sample		From module, calibration	50
Template calibration	76	From module, measurement	53
Template incremental	77	Method/Series	56
Template sensor test	78	Method/series editor	50
Screen	35	Read	50
Sensor test		Shortcut	50
Template sensor test	92	Start Method/Series	
Sensors	49	After editing	56
Actions when sensors expire, set up	43	From method/series list	56
Add	27	Via shortcut/direct shortcut	57
Calibration history	49	Statistics	115
Delete	29	Stirrer	32
Modify	29	Symbol	5
Sensor list	49	warning	6
Sensor parameters	49	System	36
Setup	27	Date/Time	38
Series	111	Identification	37
Change sample ID	112	User management	39
Delete sample	113	T	
Delete series	113	Tables	44
Insert sample	112	Predefined	44
Modify	112	User-defined	45
Setup	111	User-defined, create	45
Statistics	115	User-defined, delete	47
Setup		User-defined, delete values	46
Buffers and standards	22	User-defined, enter values	45
Communication settings	33	User-defined, modify	46
Global settings	22	User-defined, modify values	46

Task	
Display method function	121
Interrupt	122
Resume	122
Task button	121
Status	121
Task list	
Delete task	122
Move task	122
Remove tasks	122
Tasks	121
Template calibration	
Calibration	86
Sample	76
Template incremental	
Incremental	91
Sample	77
Template sensor test	
Sample	78
Sensor test	92
Terminate error	57
third party licenses	5
Title bar	122

U

Update	48
USB stick	31
User	39
Assign rights	39
Delete	40
Edit	40
Enter	39
User groups	42
User list	39
User rights	42
User management	39
Account policies	40
User	39
User management settings	41
User setting	
Keyboards	36
User settings	
Beep	35
Language	34
Screen	35
Setup	34

User settings	
Shortcuts	35

W

warning symbol	6
----------------	---

To protect your product's future:
METTLER TOLEDO Service assures
the quality, measuring accuracy and
preservation of value of this product
for years to come.

Please request full details about our
attractive terms of service.

www.mt.com/phlab

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
Tel. +41 22 567 53 22
Fax +41 22 567 53 23
www.mt.com/contact

Subject to technical changes.
© Mettler-Toledo GmbH 03/2021
30034459H en



30034459