
METTLER TOLEDO
InLab® 605 O₂-Sensor
InLab® 605-ISM-Sensor

General instructions
Betriebsanleitung
Instructions générales
Istruzioni d'uso generali
Instrucciones generales
一般的な手順
般指示
Руководство

METTLER TOLEDO

The logo graphic consists of a diamond shape filled with parallel diagonal lines. A stylized 'M' is formed by the intersection of these lines, with the right side of the 'M' being more densely shaded.

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Intelligent Sensor Management (ISM®)

Intelligent Sensor Management (ISM®)

The METTLER TOLEDO dual channel portable meters SevenGo Duo (pro) SG23, 68 and 78 have the ability to automatically recognize ISM sensors upon connection, transfer calibration data to and from the sensor etc. For more information regarding ISM® please read the respective chapter of the operating instructions supplied with the meter.

Intelligent Sensor Management (ISM®)

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ISM® (Intelligent Sensor Management)

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Intelligent Sensor Management (ISM™)

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Instruction manual InLab® 605 O₂-Sensor

1 Product description

1.1 Utilization, conformity and identification

The InLab 605 is a polarographic oxygen sensor (Clark principle) designed for the simultaneous measurement of dissolved oxygen and temperature in water applications. The robust IP67 connectors provide waterproof performance in outdoor or industrial applications. The InLab 605 represents a state-of-the-art product and is built under the strict quality guidelines for ISO 9001. The serial number is engraved on the sensor head for traceability.

1.2 Key data and characteristics

Length	120 mm
Diameter	12 mm
Connectors	IP67 BNC/Cinch
Temperature sensor	NTC
Shaft material	PPS, glass fiber reinforced
Membrane material	Silicone/Teflon (reinforced with steel mesh)
O-Ring material	Viton, Silicone
Temperature range	0...60 °C
Humidity	0...100% rel.
Media	The sensor is designed for use in aqueous samples.

1.3 Equipment features and scope of delivery

The oxygen sensor is shipped with a silicone membrane body in place and filled with electrolyte, thus it is ready for use. The InLab 605 features a built-in temperature device (NTC) that allows for automatic compensation of temperature changes.

2 Safety

Please read this instruction manual entirely before using the sensor.

3 Measuring dissolved oxygen

3.1 Initial start-up

For instructions concerning the connection of the IP67 cable to the instrument, please refer to the instrument instruction manual.

Notice: The protective cap at the tip of the sensor should be removed before putting the sensor into operation. When the sensor is installed for the first time or has been disconnected from the instrument for longer than 5 minutes, the sensor must be polarized before calibration or the first measurements. This can be achieved by connecting the sensor to the instrument for at least six hours. If the sensor has been disconnected for less than 5 minutes a polarization time of 10 minutes suffices.

3.2 General hints

There should be enough fluid circulation at the head of the sensor. At least 1 liter per hour should pass by the membrane. Either the sample needs to be stirred or the sensor should be used in a place with natural flow (e.g. river). Placing the sensor in a dead spot would lead to false results (too low). It is preferable to hold it at an angle against the flow to achieve the best possible measurement. An angle of 15 to 75° against the flow prevents the accumulation of air bubbles on the sensor head.

4 Storage and maintenance

The sensor and the membrane must be cleaned before storing the sensor. The sensor can be stored for several months, provided it is filled with O₂ electrolyte and the protective cap is placed over the membrane. To avoid the 6 hour polarization requirement after storage, the sensor can be stored connected to the instrument. After storage > 3 months the electrolyte should be replaced. If storage > 6 months is intended the sensor should be stored dry, without electrolyte.

When used in water applications, the sensor is designed to require minimal service. However in some cases it can be necessary to recalibrate the sensor or to change the electrolyte or the membrane body. The calibration interval depends on the required accuracy and individual use of the sensor. Therefore, a general recommendation cannot be made.

4.1 Safety precautions

Danger: the electrolyte has an alkaline pH value of 13. Contact of the electrolyte with the skin, especially mucous membrane or eyes, should be avoided. If such contact occurs, the affected area should be well rinsed with running water. Get medical attention if adverse

signs appear. As contact with the electrolyte is very likely during the exchange of electrolyte or membrane body, the use of protective gloves is recommended.

4.2 Cleaning and care

The outside of the sensor can be gently cleaned with soapy water (or with a mild bleach) and rinsed thoroughly. Caution: Do not clean the interior body with cleaning agents or ethanol. This might damage the sensor and lead to erroneous results. Inspect the membrane. If it has cracks or if the sensor shows longer response time or if the signal drifts significantly, then the membrane should be replaced as follows:

4.3 Replacing the membrane and the electrolyte

Please follow the following instructions to change the membrane body and the electrolyte.

- Unscrew the cap sleeve from the shaft and carefully pull it off the sensor.
- Pull off the membrane body from the interior body. If it remains inside the cap sleeve, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve.
- Clean the interior body using de-ionized water and a soft tissue. Make sure it is dry before proceeding.
- Check the O-rings for mechanical defects and replace if necessary.
- Half fill the membrane body with O₂ electrolyte and make sure it is bubble-free. Air bubbles can be removed by carefully tipping on the side of the membrane body.
- Slowly slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and should be wiped off with a paper tissue.
- Carefully slip the cap sleeve over the fitted membrane body and screw it down. The cap sleeve must be clean and dry on the inside.
- After each replacement of the electrolyte or of the membrane, the sensor must be polarized over a period of 6 hours.
- Recalibrate before measuring.

4.4 Troubleshooting and rectification work

On a regular basis, you can take the sensor out of the water, clean and dry it. If it gives a reading of 100% saturation in air, the sensor does not need to be recalibrated. If the sensor gives too high or too low values, it should be recalibrated. After a long storage time (more than 3 months) or after a certain period of time in operation (typically one year for water applications), the electrolyte should be replaced. Replace damaged membranes as necessary. If after a calibration the sensor still gives values that are too high or too low or unstable, you should change the electrolyte and the membrane.

5 Waste disposal

 In accordance with the requirements of European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), this device must not be disposed of with household waste. This also applies to countries outside the EU, per their specific requirements. Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related. Thank you for your contribution to environmental protection.

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Subject to technical changes. 02/09 © Mettler-Toledo AG. Printed in Switzerland.
Mettler-Toledo AG