METTLER TOLEDO

**pH Sensor
OptiMax SG/325**

Instruction Manual

de	Betriebsanleitung	2
en	Instruction Manual	4
es	Manual de instrucciones	6
fr	Instructions d'utilisation	8
it	Istruzioni per l'uso	10
ja	取扱説明書	12
nl	Gebruikershandleiding	14
zh	说明书	16

OptiMax SG/325

pH-Elektroden

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1. Einleitung

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Instruction Manual

1 Introduction

The METTLER TOLEDO OptiMax SG/325 pH single rod measuring cell is a low maintenance, autoclavable and pressure resistant pH electrode with a pre-pressurised liquid electrolyte and a reference electrolyte together with a silver ions barrier. Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use. We recommend that you only operate the electrode in conjunction with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

2 Safety Instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis both protective spectacles and protective gloves should be worn. Do not touch the electrode at the plug-in contacts, since electrostatic discharges can damage the electronics.

3 Product Description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
pH Sensor OptiMax SG/325	Type Designation: SG = auxillary electrode/ 325 = rod length in mm
pH 0 ... 14	pH-measurement range
0 ... 140 °C	Temperature range for operation (for sterilisation)
Order No. 30934431	Order number
IS/I, II, III/1/ABCDEFG/T6	FM tested
See instruction manual	Note

In addition each electrode is provided with a serial number on the black part of the plug-in head in order to enable identification.

4 Installation and Commissioning

1. While unpacking check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the protective cap and rinse the electrode for a short time with deionised water. Carefully remove the silicon protection from the diaphragm with the aid of the knife supplied (cf. green instruction images on the packaging). If the electrode has already been used the sensitive part can be cleaned if necessary with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionised water. After rinsing the electrode should only be dabbed. Any rubbing of the pH-sensitive glass can increase the response time as a result of electrostatic charge accumulation.
3. Check the space behind the pH sensitive glass for air bubbles and remove any air bubbles that may be present by shaking the electrode lightly in a vertical direction.
4. Connect the electrode and the OptiMax with the corresponding cable (VP6, VP8, AK9). Observe the connections scheme that accompanies the cable or the transmitter.

5 Operation:

Calibration of the electrode and the OptiMax

Before a calibration immerse the electrode for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connecting it to the OptiMax. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. For further details please consult the operating instructions for the OptiMax.

6 Maintenance

1. After each working cycle the electrode tip and diaphragm(s) are to be carefully rinsed with deionised water. Any drying out of residues from the solution being measured on these parts must be avoided at all costs!
2. If the electrode is not in use it is stored with the electrode tip and the diaphragm immersed in a tank of electrolyte (9823/3M KCl or 9848 FriscoIyt™).
3. If the electrode is stored in its armature the rules described under Item 2 apply, although the electrolyte has to be slightly modified for storage purposes (2 parts of buffer solution pH 9.2 to 10 parts of electrolyte) in order to avoid any corrosion of the armature parts.
4. If an electrode is inadvertently stored in a dry condition for several days it must be rinsed for several hours in a standard storage electrolyte.
5. The plug-in head should occasionally be checked for possible traces of moisture. If necessary clean it thoroughly with deionised water or alcohol and then dry it off carefully.

7 Removal of Sources of Malfunction

With a soft toothbrush and water any residues can easily be removed from the electrode. Mild washing agents can also be used. Protein contaminations on the diaphragm can be most likely removed with the cleaning solution (order no. 51340070).

⚠ Warning: Solutions containing hydrofluoric acid damage the pH sensitive part of the electrode.

8 Environmental Protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9 Warranty

On manufacturing defects 12 months after delivery.

10 Storage Conditions

We recommend that the electrodes are stored at room temperature.

Friscolyt is a trademark of the METTLER TOLEDO Group.

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All rights reserved