English

Reference Manual **Pipette Calibration Module** MCP-R

Deutsch

Referenzhandbuch **Pipettenkalibriermodul** MCP-R

Español

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Table of Contents

1	Introduction	3
1.1	Document purpose	3
1.2	Further documents and information	3
1.3	Explanation of conventions and symbols used	3
1.4	Acronyms and abbreviations	3
1.5	Product compatibility	4
2	Safety Information	4
2.1	Definitions of signal words and warning symbols	4
2.2	Product-specific safety notes	4
3	Function	5
3.1	Function description	5
4	Installation	5
4.1	Selecting the location	6
4.2	Scope of delivery	7
4.2.1	MCP-R pipette calibration module	7
4.2.2	Documentation	7
4.3	Installing the MCP-R pipette calibration module on the balance	8
5	Operation	8
5.1	Filling the evaporation trap	8
5.2	Performing a pipette calibration	8
5.3	Performing a test	9
5.3.1	Performance test or sensitivity test	9
5.3.2	Positioning test	10
5.3.3	MCP repeatability test	11
6	Maintenance	11
6.1	Maintenance tasks	12
6.2	Emptying the containers	12
6.2.1	Emptying the containers manually	12
6.2.1.1	Emptying the pipetting tanks	12
6.2.1.2	Emptying the evaporation trap	12
6.2.2	Emptying the containers with the pump	13
6.2.2.1	Assembling the pump	13
6.2.2.2	Emptying the pipetting tanks	13
6.2.2.3	Emptying the evaporation trap	13
6.3	Disassembling for cleaning	14
6.4	Cleaning	15
6.4.1	Cleaning the instrument	16
6.4.2	Cleaning after an overflow	16
6.5	Putting into operation after cleaning	16
6.6	Adjusting the holder for pipetting tanks	17
6.6.1	Horizontal adjustment	17
6.6.2	Vertical adjustment	17
7	Troubleshooting	17
7.1	Error symptoms	18
8	Technical Data	18
8.1	General data	18
8.2	Dimensions	19

9	Accessories and Spare Parts	20
9.1	Accessories.....	20
9.2	Spare parts.....	21
10	Disposal	21
11	Compliance Information	22

1 Introduction

Thank you for choosing a METTLER TOLEDO instrument. The instrument combines high performance with ease of use.

1.1 Document purpose

This manual is intended to support trained staff when using the product described herein.

Restriction of content

This manual does not describe the following items:

- The balance connected to the product described herein.
- Accessories used with the product described herein. Exceptions might apply.

1.2 Further documents and information

This document is available in other languages online.

▶ www.mt.com/MCP-R-RM

Instructions for cleaning a balance, "8 Steps to a Clean Balance":

▶ www.mt.com/lab-cleaning-guide

Search for software:

▶ www.mt.com/labweighing-software-download

Search for documents:

▶ www.mt.com/library

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

▶ www.mt.com/contact

1.3 Explanation of conventions and symbols used

Conventions and symbols

Key and/or button designations and display texts are shown in graphic or bold text, e.g.,  **Edit**.

 **Note** For useful information about the product.

Refers to an external document.

Elements of instructions

In this manual, step-by-step instructions are presented as follows. The action steps are numbered and can contain prerequisites, intermediate results and results, as shown in the example. Sequences with less than two steps are not numbered.

- Prerequisites that must be fulfilled before the individual steps can be executed.
- 1 Step 1
 - ⇒ Intermediate result
- 2 Step 2
 - ⇒ Result

1.4 Acronyms and abbreviations

Original term	Explanation
AC	Alternating Current
ASTM	American Society for Testing and Materials
DC	Direct Current
EMC	Electromagnetic Compatibility
OIML	Organisation Internationale de Métrologie Légale

1.5 Product compatibility

The MCP-R pipette calibration module is only compatible with the listed balance models.

Balance	Models designation
	<ul style="list-style-type: none"> • XPR105DUHR • XPR106DUHR

2 Safety Information

A document named "Reference Manual" is available for this instrument.

- This Reference Manual contains a full description of the instrument and its use.
- Keep the Reference Manual for future reference.
- Include the Reference Manual if you transfer the instrument to other parties.

Only use the instrument according to the Reference Manual. If you do not use the instrument according to the Reference Manual or if it is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

DANGER	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
WARNING	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
CAUTION	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
NOTICE	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Notice

2.2 Product-specific safety notes

Intended use

This instrument is designed to be used by trained staff. The instrument is intended for volumetric calibration of pipettes only.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your accessory.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

3 Function

3.1 Function description

The MCP-R pipette calibration module can be added to compatible balances to perform volumetric calibration of single and multichannel pipettes. It can contain 12 pipetting tanks that are placed inside the module. The MCP-R pipette calibration module has a removable evaporation trap door on top that keeps liquids from evaporating during measurements and reduces influences from the outside.

The balance with the installed MCP-R pipette calibration module is connected to a computer with the specific Calibry software installed that initiates the calibration process and collects all the data. When a calibration is initiated, the user opens the evaporation trap door and pipettes the required volume into each tank. The tanks are weighed successively and the data is collected for the evaluation of the pipettes.

Note

The Calibry software is not part of the delivery and can be purchased optionally, see [Accessories ▶ Page 20].

See also

 Product compatibility ▶ Page 4

4 Installation

This instrument must be installed by a METTLER TOLEDO service technician.

4.1 Selecting the location

The MCP-R pipette calibration module is a sensitive precision instrument. The location of use has a profound effect on the accuracy of the weighing results.


Requirements of the location

Place indoors on stable

Ensure sufficient spacing

Level the instrument

Provide adequate lighting


Avoid direct sunlight

Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations

Take into account the environmental conditions. See "Technical Data".


Sufficient spacing for balances: > 15 cm all around the instrument

See also

[Technical Data](#) ▶ Page 18

4.2 Scope of delivery

4.2.1 MCP-R pipette calibration module

1	Evaporation trap door	7	Position pin, 2 pcs
2	Cover	8	Base plate
3	Pipetting tanks, 12 pcs	9	Centering gauge
4	Holder for pipetting tanks	10	Height checking tool
5	MCP-R casing	11	Tweezers
6	Holder support	12	Test weighing pan

4.2.2 Documentation

- Reference Manual
- Production certificate
- Declaration of Conformity

4.3 Installing the MCP-R pipette calibration module on the balance


This instrument must be installed by a METTLER TOLEDO service technician.

5 Operation

5.1 Filling the evaporation trap

The evaporation trap includes two elongated recesses that should be filled with water. When the evaporation trap is closed, the humidity inside increases until the air is saturated. This prevents evaporation of the liquid in the pipetting tanks and ensures accurate weighing results.

- The evaporation trap door is closed.
 - 1 Press the button to remove the cover.
 - 2 Fill both recesses, each with a maximum of 8 mL of distilled water.
 - **Note**
In the event of an overflow, the water will drain outside the instrument onto the underlying surface.
 - 3 Place the cover and slightly press it down until it snaps in.
 - 4 Wait at least **120 minutes** to allow the air in the evaporation trap to reach saturation.

5.2 Performing a pipette calibration

For further information, consult the Reference Manual of the connected balance and the Operating Instructions of the Calibry software.

► www.mt.com/MCP-R-RM

According to ISO 8655-6:2022, the readability of the balance should correspond to the nominal pipette volume to be tested. The readability can be set in the tolerance profile of the method.

For more details about specific balance settings for pipette calibration, refer to the Reference Manual of the balance.

The 12 pipetting tanks have a capacity of 17 mL each. They must be emptied regularly to avoid an overflow, especially when calibrating pipettes with a larger volume (> 200 μ L).

For performing a pipette calibration, a computer with the Calibry software installed is needed. The Calibry software leads the user through the pipette calibration. The following steps explain the handling of the MCP-R pipette calibration module while performing a pipette calibration.

- The evaporation trap is filled.
- The evaporation trap door has been closed for at least 120 minutes, allowing the air in the evaporation trap to reach saturation.
- The balance is connected to a computer with the Calibry software installed.
- An appropriate tolerance profile is selected.

- 1 Open the evaporation trap door (1).
- 2 Carefully pipette the required volume into the pipetting tanks (2).


Note

In the event of an overflow, remove the excess liquid immediately, see [Cleaning after an overflow ▶ Page 16].

- 3 Close the evaporation trap door (3).
- 4 Follow the instructions in the Calibry software.
- 5 Wait until all the pipetting tanks are weighed.
 - ➔ The weighing results are shown on the terminal display.
- 6 Repeat the procedure until the Calibry software indicates that the calibration process is finished.

Note

The results can be evaluated in the Calibry software.

5.3 Performing a test

For further information, consult the Reference Manual of the connected balance and the Operating Instructions of the Calibry software.

5.3.1 Performance test or sensitivity test


The balance software leads the user through the test. The following steps explain the preparation of the MCP-R pipette calibration module for the test.

Note

Performance test: METTLER TOLEDO recommends using OIML/ASTM test weights between 1 g and 20 g.

Sensitivity test: METTLER TOLEDO recommends using an OIML/ASTM test weight of 100 g.

- An appropriate tolerance profile is selected.
 - 1 Press the button to remove the cover.
 - 2 Remove the pipetting tanks one by one with tweezers.
 - 3 Vertically lift out the holder for pipetting tanks (1).
 - 4 Place the test weighing pan (2) on the holder support (3).
- ➔ The instrument is ready for testing.
- 5 While performing the test, close the cover to avoid draft (4).


5.3.2 Positioning test

For performing the positioning test, a computer with the Calibry software installed is needed. The Calibry software leads the user through the test. The following steps show how to use the height checking tool.


Testing the pipetting position

- 1 Place the height checking tool (1) with the pipetting position side onto the border of the inner housing (2).
 - ➔ The upper end of the spring blades must be between the two tolerance lines (3). This applies to both sides.
- 2 Check all 12 positions.

Testing the weighing position

- 1 Place the height checker tool (1) with the weighing position side onto the border of the inner housing (2).
 - ➔ The upper end of the spring blades must be between the two tolerance lines (3). This applies to both sides.
- 2 Check all 12 positions.

Failed test

If the test fails, see "Troubleshooting".

See also

🔗 [Troubleshooting](#) ▶ Page 17

5.3.3 MCP repeatability test

For performing the test, a computer with the Calibry software installed is needed.

- The pipetting tanks are empty.
- All 12 pipetting tanks are correctly placed in the holder for pipetting tanks.
- The balance is connected to a computer with the Calibry software installed.
 - Follow the instructions in the Calibry software.

Failed test

If the test fails, see "Troubleshooting".

See also

🔗 [Troubleshooting](#) ▶ Page 17

6 Maintenance

To guarantee the functionality of the instrument and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

The appropriate maintenance interval depends on your standard operating procedure (SOP).

6.1 Maintenance tasks

Maintenance action	Recommended interval	Remarks
Cleaning	<ul style="list-style-type: none">• Depending on the degree of pollution• In the event of an overflow• Depending on your internal regulations (SOP)	see "Cleaning"
Emptying the containers	<ul style="list-style-type: none">• Depending on the frequency of usage of the instrument• Depending on your internal regulations (SOP)	see "Emptying the containers"
Performing routine tests (performance test, positioning test, MCP repeatability test, sensitivity test).	<ul style="list-style-type: none">• After cleaning• After assembling the instrument• After a software update• Depending on your internal regulations (SOP)	see "Performing a test" see the Reference Manual of the balance

See also

- 🔗 [Cleaning](#) ▶ Page 15
- 🔗 [Emptying the containers](#) ▶ Page 12
- 🔗 [Performing a test](#) ▶ Page 9

6.2 Emptying the containers

The pipetting tanks and the evaporation trap can be emptied manually or with the help of a pump.

Note

The pump is not part of the delivery and can be purchased optionally, see [Accessories ▶ Page 20].

6.2.1 Emptying the containers manually

6.2.1.1 Emptying the pipetting tanks

- 1 Remove the evaporation trap door.
 - 2 Press the button to remove the cover.
 - 3 Remove the pipetting tanks one by one with tweezers.
 - 4 With tweezers, carefully pull out the glass tubes from the pipetting tanks.
 - 5 Over a container suitable for liquid waste, turn the pipetting tanks upside down and drain out the water.
 - 6 Dry the glass tubes and the pipetting tanks.
 - 7 Insert the glass tubes into the pipetting tanks.
 - 8 Insert the pipetting tanks into the holder for pipetting tanks.
- ➔ The pipetting tanks are empty.

6.2.1.2 Emptying the evaporation trap

An absorbing sponge is usually sufficient to empty the evaporation trap.

- 1 Remove the evaporation trap door.
 - 2 Press the button to remove the cover.
 - 3 With a sponge, remove the liquid in the evaporation trap recesses.
- ➔ The evaporation trap is empty.

6.2.2 Emptying the containers with the pump

6.2.2.1 Assembling the pump

⚠ WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your accessory.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

NOTICE


Damage to the pump due to inappropriate handling

If the pump is switched on when it is not in use, it can get damaged.

- Switch off the pump when it is not in use.

During normal usage, the pipetting tanks and the evaporation trap recesses contain water. The pump serves to remove the water in a safe and efficient manner, without having to disassemble the instrument. The inlet tubing has a metallic tube attached for easier handling. The outlet tubing serves to dispose of the aspirated liquid.

- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
- 2 Insert the plug of the AC/DC adapter (1) into the power socket.
- 3 Secure the plug by firmly tightening the knurled nut.
- 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
- 5 Connect the tubing with the metallic tube (2) to the nozzle IN (3).
- 6 Connect the other tubing to the nozzle OUT (4).
- 7 Place the free end of the tubing into a container suitable for liquid waste (5).

6.2.2.2 Emptying the pipetting tanks

- The pump is assembled.
 - The outlet tubing is placed in a container suitable for liquid waste.
- 1 Open the evaporation trap door.
 - 2 Place the metallic tube inside one pipetting tank.
 - 3 Switch on the pump.
 - ➔ The liquid is pumped out and collected in the container for liquid waste.
 - 4 Switch off the pump.
 - 5 Repeat this procedure for the other pipetting tanks.
 - ➔ The pipetting tanks are empty.

6.2.2.3 Emptying the evaporation trap


- The pump is assembled.
 - The outlet tubing is placed in a container suitable for liquid waste.
- 1 Remove the evaporation trap door.
 - 2 Press the button to remove the cover.

- 3 Place the metallic tube into one elongated recess of the evaporation trap.
- 4 Switch on the pump.
 - ➔ The liquid is pumped out and collected in the container for liquid waste.
- 5 Switch off the pump.
- 6 Repeat this procedure for the other elongated recess of the evaporation trap.
 - ➔ The evaporation trap is empty.


6.3 Disassembling for cleaning

- The containers are empty, see [Emptying the containers ▶ Page 12].
- The connected balance is disconnected from the power supply.


- 1 Remove the evaporation trap door (1).
- 2 Press the button (2) and remove the cover (3).
- 3 Press the button to remove the cover.


- 4 Remove the pipetting tanks one by one with tweezers (4).


- 5 With tweezers, carefully pull out the glass tubes from the pipetting tanks (5).


6 Vertically lift out the holder for pipetting tanks (6).

7 Lift off the MCP-R casing from the balance (7).

8 Remove the holder support (8).

See also

 Emptying the containers ▶ Page 12

6.4 Cleaning

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

▶ www.mt.com/lab-cleaning-guide

NOTICE

Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

6.4.1 Cleaning the instrument

Cleaning around the instrument

- Remove any dirt or dust around the instrument and avoid further contaminations.

Cleaning removable parts

- Clean the removed part with a damp cloth or a tissue and a mild cleaning agent.

Cleaning not removable parts

- 1 Disconnect the connected balance from the AC/DC adapter.
- 2 Use a lint-free cloth moistened with a mild cleaning agent to clean the surface of the instrument.
- 3 Remove powder or dust with a disposable tissue first.
- 4 Remove sticky substances with a damp lint-free cloth and a mild solvent, e.g., isopropanol or ethanol 70%.

6.4.2 Cleaning after an overflow

In the event of an overflow, for example when overfilling the pipetting tanks, excess liquid must be removed immediately.

- 1 If needed, empty the pipetting tanks, see Emptying the pipetting tanks.
- 2 Disassemble the remaining items.
- 3 Clean the items with a lint-free cloth or tissue, see [Disassembling for cleaning ▶ Page 14].
- 4 Reassemble the instrument.

6.5 Putting into operation after cleaning

For further information, consult the Reference Manual of the connected balance.


▶ www.mt.com/XPR-analytical-RM

- 1 Reassemble the instrument in reverse order.
 - 2 Respect the warm-up time of the connected balance.
 - 3 If needed, adjust the holder for pipetting tanks, see [Adjusting the holder for pipetting tanks ▶ Page 17].
 - 4 Perform a positioning test, see [Positioning test ▶ Page 10].
 - 5 Perform an MCP repeatability test, see [MCP repeatability test ▶ Page 11].
- ⇒ The instrument is ready to be used.

6.6 Adjusting the holder for pipetting tanks


6.6.1 Horizontal adjustment

- A screwdriver of type T8 is available.
- 1 Place the centering gauge (1) from underneath (2) into the holder for pipetting tanks.
 - 2 Make sure the pin holes (3) are orientated towards the balance and slowly place the holder for pipetting tanks onto the holder support (4).
 - 3 Enter the screwdriver through the centering gauge to the screw (5).
 - 4 Turn the screw 180° in counterclockwise direction to release the tension (6).
 - 5 Turn the screw approximately 180° in clockwise direction to tighten it again (7).
 - ➔ The holder for pipetting tanks is now correctly positioned horizontally.
 - 6 Vertically lift out the holder for pipetting tanks.
 - 7 Remove the centering gauge from the holder for pipetting tanks.
 - 8 Place the holder for pipetting tanks.

6.6.2 Vertical adjustment

- A screwdriver of type Hex 1.27 is available.
- 1 With tweezers, place the two outer pipetting tanks (1) onto the holder for pipetting tanks.
 - 2 Adjust the central screw of the holder for pipetting tanks (2) until the holder for pipetting tanks lies at all four edges approximately 1 to 2 mm below the pipetting tank (3).

7 Troubleshooting

For further information, consult the Reference Manual of the connected balance.

► www.mt.com/XPR-analytical-RM

7.1 Error symptoms

Error symptom	Possible cause	Diagnostic	Remedy
The time to reach weight value stability is longer than usual.	The environment is unstable.	-	Make sure the environment is stable, see [Selecting the location ▶ Page 6].
	The pipetting tanks are touching the holder for pipetting tanks when the instrument is not in use.	Check for touching points.	Adjust the position of the holder for pipetting tanks, see [Adjusting the holder for pipetting tanks ▶ Page 17].
		Check if one or more spring blades at the support for pipetting tanks are bent, see [Positioning test ▶ Page 10].	Contact a METTLER TOLEDO service technician.
The weight value is unstable at certain tank positions.	The outer surface of one or more pipetting tanks is wet.	Check for wet surfaces.	Dry the outer surface of the pipetting tanks.
	The pipetting tanks are touching the holder for pipetting tanks when the instrument is not in use.	Check if the holder for pipetting tanks is adjusted.	Adjust the position of the holder for pipetting tanks, see [Adjusting the holder for pipetting tanks ▶ Page 17].
		Check if one or more spring blades at the support for pipetting tanks are bent, see [Positioning test ▶ Page 10].	Contact a METTLER TOLEDO service technician.
The positioning test failed.	The support for pipetting tanks is bent.	Check if one or more spring blades at the support for pipetting tanks are bent, see [Positioning test ▶ Page 10].	Contact a METTLER TOLEDO service technician.
The MCP repeatability test failed.	The pipetting tanks are touching each other or are not placed properly.	Check visually.	Perform a positioning test.
	The balance is defective.	-	Contact a METTLER TOLEDO service technician.

8 Technical Data

8.1 General data

Weight:	2.5 kg
Dimensions (W × D × H):	145 × 160 × 116 mm
Power consumption:	12 V DC ± 10%, 1 A

Protection and standards

Degree of pollution:	2
Overvoltage category:	II
Range of application:	Use only indoors in dry locations

Operating conditions

Temperature:	+15 – +30 °C
Humidity:	40 – 70 %, non-condensing

The instrument can be used under the following operating conditions. However, the performance of the instrument may be outside the limit values:

Temperature:	+5 – +40 °C
Humidity:	80% for temperatures up to 31 °C, decreasing linearly to temperature 5–40 °C non-condensing


Environmental conditions

Height above mean sea level:	Up to 5000 m
Ambient temperature:	+5 – +40 °C
Relative air humidity:	20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

Storage conditions (in packaging)

Ambient temperature:	-25 – +70 °C
Relative air humidity:	10 – 90%, non-condensing

8.2 Dimensions

↔	Clear dimensions [mm]
↔	Outer dimensions [mm]


9 Accessories and Spare Parts

9.1 Accessories

Accessories are additional components that could help you in your workflow.

	Description	Order no.
Software		
	Calibry Single workstation; for calibration with one system MPC	11138419
	Calibry Network; Installation on several PC of the network accessing the same database	11138420
Various		
	Reagent reservoirs, 5 pcs	17007886
	Suction pump	30676828

9.2 Spare parts

	Order no.	Designation	Remarks
1	30480598	Evaporation trap door	–
2	30480599	Cover	–
3	30480596	Pipetting tank	Set of 2 pcs
4	30480602	MCP-R casing	Excluding: type label, cover, evaporation trap door, pipetting tanks, holder pipetting tanks, weighing pan, base plate
5	30480594	Support, holder	–
6	30480595	Holder, pipetting tanks	–

10 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

11 Compliance Information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

► www.mt.com/ComplianceSearch

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

► www.mt.com/contact

United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

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© 10/2022 METTLER TOLEDO. All rights reserved.