
English	User Manual Microbalances XPR Essential
Deutsch	Benutzerhandbuch Mikrowaagen XPR Essential
Español	Manual de usuario Microbalanzas XPR Essential
Français	Guide de l'utilisateur Microbalances XPR Essential
Italiano	Manuale per l'utente Microbilance XPR Essential
简体中文	用户手册 微量天平 XPR Essential
日本語	ユーザマニュアル マイクロ天びん XPR Essential
Čeština	Návod k použití Mikrováhy XPR Essential
Dansk	Brugervejledning Mikrovægte XPR Essential
Hrvatski	Korisnički priručnik Mikro vage XPR Essential

METTLER TOLEDO

► www.mt.com/XPR-Essential-micro-RM

en For full information, always refer to the Reference Manual (RM).

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Table of Contents

1	Introduction	3
1.1	Document purpose	3
1.2	Further documents and information	3
1.3	Acronyms and abbreviations	3
2	Safety Information	4
2.1	Definitions of signal words and warning symbols	4
2.2	Product-specific safety information	4
3	Design and Function	5
3.1	Overview balance	6
3.2	Overview interface connections	6
3.3	Components description	7
3.3.1	Draft shield	7
3.3.2	Door-release bar	7
3.3.3	Handles	7
3.3.4	Weighing pan	8
3.3.5	Drip tray	8
3.3.6	Leveling feet	8
3.3.7	Terminal	8
3.4	Overview terminal	9
3.5	User interface	10
3.5.1	Main sections at a glance	10
3.5.2	Main weighing screen	10
4	Installation and Putting into Operation	11
4.1	Selecting the location	11
4.2	Unpacking the balance	12
4.3	Installation	13
4.3.1	Attaching the terminal	13
4.3.2	Assembling the balance	15
4.4	Putting into operation	16
4.4.1	Connecting the balance	16
4.4.2	Switching on the balance	16
4.4.3	Leveling the balance	17
4.4.4	Performing an internal adjustment	17
4.4.5	Entering / Exiting standby mode	17
4.4.6	Switching off the balance	17
4.5	Performing a simple weighing	18
4.5.1	Opening and closing the draft shield doors	18
4.5.2	Zeroing the balance	18
4.5.3	Taring the balance	18
4.5.4	Performing a weighing	18
4.5.5	Completing the weighing	18
4.6	Transporting, packing, and storing	19
4.6.1	Transporting the balance over short distances	19
4.6.2	Transporting the balance over long distances	19
4.6.3	Packing and storing	19
5	Maintenance	20
5.1	Maintenance tasks	20
5.2	Cleaning	20
5.2.1	Disassembling for cleaning	20
5.2.2	Cleaning the balance	22

5.2.3	Putting into operation after cleaning	22
5.3	Service.....	23
6	Technical Data	23
6.1	General data.....	23
7	Disposal	24
8	Compliance Information	24

1 Introduction

Thank you for choosing a METTLER TOLEDO balance. The balance combines high performance with ease of use.

EULA

The software in this product is licensed under the METTLER TOLEDO End User License Agreement (EULA) for Software.

When using this product you agree to the terms of the EULA.

▶ www.mt.com/EULA

1.1 Document purpose

This User Manual provides brief instructions about the first steps to take with the instrument. This ensures safe and efficient handling. Personnel must have carefully read and understood this manual before performing any task.

1.2 Further documents and information

This document is available in other languages online.

▶ www.mt.com/XPR-Essential-micro-UM

Product page:

▶ www.mt.com/XPR-Essential-micro

Instructions for cleaning a balance, "8 Steps to a Clean Balance":

▶ www.mt.com/lab-cleaning-guide

Search for software:

▶ www.mt.com/labweighing-software-download

Search for documents:

▶ www.mt.com/library

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

▶ www.mt.com/contact

1.3 Acronyms and abbreviations

Original term	Explanation
AC	Alternating Current
ASTM	American Society for Testing and Materials
DC	Direct Current
EMC	Electromagnetic Compatibility
FCC	Federal Communications Commission
GWP	Good Weighing Practice
HID	Human Interaction Device
ID	Identification
LED	Light-Emitting Diode
LPS	Limited Power Source
MAC	Media Access Control
MT-SICS	METTLER TOLEDO Standard Interface Command Set
NA	Not Applicable

OIML	Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology)
RAM	Random Access Memory
RFID	Radio-frequency identification
RM	Reference Manual
SELV	Safety Extra Low Voltage
SOP	Standard Operating Procedure
SQC	Statistical Quality Control
UM	User Manual
USB	Universal Serial Bus
USP	United States Pharmacopeia

2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is available online in various languages.
- A printed version of the User Manual is delivered with the instrument.
- The Reference Manual is available online. This manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

DANGER	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
WARNING	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
CAUTION	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
NOTICE	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Notice

2.2 Product-specific safety information

Intended use

This instrument is designed to be used by trained staff. The instrument is intended for weighing purposes.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

A list of spare parts and accessories can be found in the Reference Manual.


3 Design and Function

For further information, consult the Reference Manual (RM).


► www.mt.com/XPR-Essential-micro-RM

3.1 Overview balance

1	StatusLight	7	Model label
2	Terminal	8	Side door, draft shield (right/left)
3	Door-release bar	9	Handle, side door
4	Weighing pan	10	Drip tray
5	Front panel, draft shield	11	Leveling feet
6	Lid, draft shield		

3.2 Overview interface connections

1	Ethernet port	4	Socket for cable to terminal
2	USB-A ports (to device)	5	USB-B port (to host)
3	Service seal	6	Socket for AC/DC adapter

NOTICE

Possible electromagnetic interference with other devices


If the Ethernet cable is longer than 30 meters, electromagnetic interference with other devices may be possible.

- Use an Ethernet cable that is shorter than 30 meters.

3.3 Components description


3.3.1 Draft shield

The draft shield is a housing device that protects the weighing area against environmental impacts like drafts or moisture. The side doors are wing doors. They rotate open when the user presses the door-release bar. The lid and the front panel are to be removed manually.


3.3.2 Door-release bar

The door-release bar is located in front of the draft shield. The bar ends are used to open the corresponding side door.


3.3.3 Handles

The handles on the side doors are used to close the side doors. The handle on the lid is used to remove or insert the lid.


3.3.4 Weighing pan

The SmartPan weighing pan is the load receptor that serves directly to accommodate the weighing item.


3.3.5 Drip tray

The drip tray is positioned below the weighing pan on the weighing chamber base plate. The primary purpose of the drip tray is to ensure quick cleaning of the balance.


3.3.6 Leveling feet

The balance stands on two height-adjustable feet. These feet are used to level the balance.


3.3.7 Terminal

The 7-inch balance terminal has a touch-sensitive display. Further, on the front side of the terminal is a StatusLight LED strip that indicates the current status of the balance.

3.4 Overview terminal


	Name	Description
	Standby	By tapping  , the balance is not completely switched off but goes into standby mode. To switch the balance completely off, it must be unplugged from the power supply. Note Do not disconnect the balance from the power supply unless the balance is not used for an extended period of time. After switching on the instrument, it must warm up before giving accurate results.
	Tare	Tares the balance. This function is used when the weighing process involves containers. After taring the balance, the screen shows Net which indicates that all displayed weights are net.
	Zero	Zeroes the balance. The balance must always be zeroed before starting the weighing process. After zeroing, the balance sets a new zero point.
	Home	To return from any menu level to the main weighing screen.

3.5 User interface

3.5.1 Main sections at a glance

The main weighing screen (1) is the central navigation point where all the menus and settings can be found. The **Balance menu** (2), **Methods** (3) and **Results** (4) open when tapping the tabs along the sides of the main weighing screen.

See also

[Main weighing screen](#) ▶ Page 10

3.5.2 Main weighing screen

Name	Description
1 User name	Shows the name of the current user.
2 Weighing value field	Shows the current weighing value.

	Name	Description
3	Level indicator	Indicates if the balance is leveled (green) or not (red).
4	Methods menu	Accesses the user-defined list of methods, tests, and alignments.
5	Info weight	Shows the current weighing value in another unit.
6	Warning and error message area	Shows current warning and/or error messages.
7	Results list	Shows the weighing results saved for this task.
8	Sample status OK	Result status indicator green: indicates that the result fulfills a set of criteria. For example: <ul style="list-style-type: none"> • The balance is in level. • The internal adjustment was performed and ok. • The weighing result is within the defined tolerance range (only if tolerance is defined).
9	Sample status Excluded	Result status indicator black: indicates that the result was excluded from the Results list .
10	Sample status Not OK	Result status indicator red: indicates that the result criteria are not fulfilled, e.g., "The weighing result was out of the defined tolerances".
11	Button Add result	Adds the result to the Results list . Depending on the selected method, the button can have different functions.
12	Action bar	Contains actions referring to the current task.
13	Balance menu	Accesses the balance properties.
14	Method information area	Contains information about the sample, method or task IDs.
15	SmartTrac	Used as a weighing aid to define a target weight with upper and lower tolerances.
16	Weighing value area	Shows the results of the current weighing process.
17	Method name	Shows the name of the current method.

4 Installation and Putting into Operation

4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.


Requirements of the location

Place indoors on stable table

Ensure sufficient spacing

Level the instrument

Provide adequate lighting



Avoid direct sunlight

Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations


Sufficient spacing for balances: > 15 cm all around the instrument

Take into account the environmental conditions. See "Technical Data".


4.2 Unpacking the balance

Check the package, the packaging elements and the delivered components for damages. If any components are damaged, please contact your METTLER TOLEDO service representative.


- 1 Open the box and lift the package out using the lifting strap (1).


- 2 Open the lifting strap (1) and remove the User Manual (2).


- 3 Remove the upper part of the package and remove the set with the AC adapter and power cable (3), the box containing several accessories (4), and the weighing pan (5).

4 Carefully remove the terminal (6).

5 Carefully remove the package set with the draft shield doors and the display holder (7).


6 Carefully remove the weighing unit (8) from the bottom packaging.

7 Remove the protective bag.

8 Store all parts of the packaging in a safe place for future use.

➔ The weighing unit is ready for assembling.

4.3 Installation

4.3.1 Attaching the terminal

NOTICE


Damage to the cables due to careless handling

- Do not kink or twist the cables.


Note

Attaching the terminal to the balance is optional. The terminal can also be placed next to the balance.


1 Attach the front cover (1) to the front of the weighing unit (2).


- 2 Insert the slides of the display holder (3) into the front of the weighing unit (2).


- 3 Connect the terminal cable (4) with the terminal (5).
NOTICE: The mark on the cable plug must point upwards in order to insert the plug in the correct orientation.


- 4 Place the terminal (5) onto the terminal holder (6).

- 5 Carefully tilt the balance to its side.
- 6 Lead the cable (7) through the cable channel (8).
- 7 Carefully put the balance back on its feet.

- 8 Insert the plug of the terminal cable (9) into the socket of the balance (10).
➔ The terminal is ready.

4.3.2 Assembling the balance


CAUTION

Injury due to sharp objects or broken glass


Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.

- 1 Insert the cover (1) and push it down on the weighing unit (2).

- 2 Insert the drip tray (3).
- 3 Tilt the weighing pan (4) and carefully mount it on the weighing hooks (5).


- 4 Slide the side door (6) onto the pivot (7).


Note

The screw hole (8) must face towards the front.


- 5 Tighten the screw (8) with the supplied 1.5 mm hex key (9).
- 6 Repeat the procedure to attach the other side door.

- 7 Push the front panel (10) firmly into the groove (11).

- 8 Push down the handle (12) to close the side door (13, right/left).
 - 9 Slide in the lid (14) until it engages.
NOTICE: Be careful not to drop the lid into the weighing area.
- ➔ The balance is assembled and ready to be put into operation.

4.4 Putting into operation

4.4.1 Connecting the balance


WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
 - 2 Insert the plug of the AC/DC adapter (1) into the power socket of the instrument.
 - 3 Secure the plug by firmly tightening the knurled nut.
 - 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
- ➔ The balance automatically switches on.

Note

Do not connect the instrument to a power outlet controlled by a switch. After switching on the instrument, it must warm up before giving accurate results.

See also

🔗 General data ▶ Page 23

4.4.2 Switching on the balance

When connected to the power supply, the balance automatically switches on.

EULA (End User License Agreement)

When the balance is switched on the first time, the EULA (End User License Agreement) appears on the screen.

- 1 Read the conditions.
- 2 Tap **I accept the terms in the license agreement.** and confirm with **✓ OK.**

- ➔ The main weighing screen appears.

Acclimatization and warm up

Before the balance gives reliable results, it must:

- acclimatize to the room temperature
- warm up by being connected to the power supply

The acclimatization time and warm-up time for balances are available in "General data".

Note

When the balance is exiting standby, it is ready immediately.

See also

-  General data ▶ Page 23
-  Switching off the balance ▶ Page 17
-  Entering / Exiting standby mode ▶ Page 17

4.4.3 Leveling the balance


Exact horizontal and stable positioning are essential for repeatable and accurate weighing results.

If the message **Balance is out of level** appears:

- 1 Tap ▶ **Level the balance**.
 - ➔ The **Leveling aid** opens.
- 2 Turn both leveling feet (1) as instructed on the display until the dot is in the center of the level indicator.

The leveling aid can also be accessed through the **Balance menu**:

☰ **Navigation:** ▶ **Balance menu** > ⦿ **Leveling aid**

4.4.4 Performing an internal adjustment

☰ **Navigation:** ▼ **Methods** > ⚙ **Adjustments**

- The adjustment **Strategy** is set to **Internal adjustment**.
- 1 Open the **Methods** section, tap ⚙ **Adjustments**, select the adjustment, and tap ▶ **Start** - or - from the main weighing screen, tap ... **More** and tap **Start adjustment**.
 - ➔ **Internal adjustment** is being executed.
 - ➔ When the adjustment has been completed, an overview of the adjustment results appears.
 - 2 Tap 🖨 **Print** if you want to print the results.
 - 3 Tap ✓ **Finish adjustment**.
 - ➔ The balance is ready.

4.4.5 Entering / Exiting standby mode

- 1 To enter standby mode, hold ⏻.
 - ➔ The display is dark. The balance is still switched on.
- 2 To exit standby mode, press ⏻.
 - ➔ The display is turned on.

4.4.6 Switching off the balance

To completely switch off the balance, it must be disconnected from the power supply. By holding ⏻, the balance goes only into standby mode.

Note

When the balance was completely switched off for some time, it must warm up before it can be used.


See also

 Switching on the balance ▶ Page 16

4.5 Performing a simple weighing

4.5.1 Opening and closing the draft shield doors

- 1 Press the end of the door-release bar (1) to open the side door (2).
- 2 Push down the handle (3) to close the side door (2).


4.5.2 Zeroing the balance

- 1 Open the draft shield.
- 2 Clear the weighing pan.
- 3 Close the draft shield.
- 4 Press $\rightarrow 0 \leftarrow$ to zero the balance.
➔ The balance is zeroed.

4.5.3 Taring the balance

If a sample vessel is used, the balance must be tared.

- 1 Open the draft shield.
- 2 Clear the weighing pan.
- 3 Close the draft shield.
- 4 Press $\rightarrow 0 \leftarrow$ to zero the balance.
- 5 Open the draft shield.
- 6 Place the sample vessel (1) on the weighing pan (2).
- 7 Close the draft shield.
- 8 Press $\rightarrow T \leftarrow$ to tare the balance.
➔ The balance is tared. The icon  Net appears.

4.5.4 Performing a weighing

- 1 Open the draft shield.
- 2 Place the weighing object into the sample vessel.
- 3 Close the draft shield.
- 4 Tap $\rightarrow + \text{Add result} \leftarrow$ if you want to report the weighing result.
➔ The result is added to the **Results list**.

4.5.5 Completing the weighing

- 1 To save the **Results list**, tap  **Complete**.
➔ The window **Complete task** opens.
- 2 Select an option to save or print the **Results list**.
➔ The respective dialog opens.

- 3 Follow the instructions from the wizard.
- 4 Tap **✓ Complete**.
 - ➔ The **Results list** is saved/printed and then cleared.


4.6 Transporting, packing, and storing

4.6.1 Transporting the balance over short distances

- 1 Disconnect the AC/DC adapter and unplug all interface cables.
- 2 Hold the weighing platform with both hands and carry the balance in horizontal position to the target location. Consider the requirements of the location.

If you want put the balance into operation, proceed as follows:

- 1 Connect in reverse order.
- 2 Level the balance.
- 3 Perform an internal adjustment.

See also

- 🔗 [Selecting the location](#) ▶ Page 11
- 🔗 [Switching on the balance](#) ▶ Page 16
- 🔗 [Leveling the balance](#) ▶ Page 17
- 🔗 [Performing an internal adjustment](#) ▶ Page 17
- 🔗 [Switching on the balance](#) ▶ Page 16

4.6.2 Transporting the balance over long distances

METTLER TOLEDO recommends using the original packaging for transportation or shipment of the balance or balance components over long distances. The elements of the original packaging are developed specifically for the balance and its components and ensure maximum protection during transportation.

See also

- 🔗 [Unpacking the balance](#) ▶ Page 12

4.6.3 Packing and storing

Packing the balance

Store all parts of packaging in a safe place. The elements of the original packaging are developed specifically for the balance and its components, and ensures maximum protection during transportation and storage.

Storing the balance

Only store the balance under the following conditions:

- Indoor and in the original packaging
- According to the environmental conditions, see "Technical Data"

i Note

When storing for longer than 6 months, the rechargeable battery may become empty (only date and time get lost).

See also

- 🔗 [Technical Data](#) ▶ Page 23

5 Maintenance

To guarantee the functionality of the balance and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

For further information, consult the Reference Manual (RM).

► www.mt.com/XPR-Essential-micro-RM

5.1 Maintenance tasks

Maintenance action	Recommended interval	Remarks
Performing an internal adjustment	<ul style="list-style-type: none">• Daily• After cleaning• After leveling• After changing the location	see "Performing an internal adjustment"
Performing routine tests (eccentricity test, repeatability test, sensitivity test). METTLER TOLEDO recommends to at least perform a sensitivity test.	<ul style="list-style-type: none">• After cleaning• After assembling the balance• After a software update• Depending on your internal regulations (SOP)	see "Tests" in the Reference Manual
Cleaning	<ul style="list-style-type: none">• After every use• After changing the substance• Depending on the degree of pollution• Depending on your internal regulations (SOP)	see "Cleaning"
Updating the software	<ul style="list-style-type: none">• Depending on your internal regulations (SOP).• After a new software release.	see "Software update" in the Reference Manual

See also

🔗 Performing an internal adjustment ► Page 17

🔗 Cleaning ► Page 20

5.2 Cleaning

5.2.1 Disassembling for cleaning


⚠ CAUTION

Injury due to sharp objects or broken glass


Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.


- 1 Pull the lid (1) towards the front and lift it off.
NOTICE: Be careful not to drop the lid into the weighing area.


- 2 Press the ends of the door-release bar (2) to open the side doors (3).
- 3 Tilt the front panel (4) sideways and carefully pull it out of the groove.


- 4 Optional: Hold the side door (3) and loosen the screw (5) with the supplied 1.5 mm hex key (6).
- 5 Pull the side door (3) towards the side to remove it.
- 6 Repeat the two previous steps to remove the other side door.

- 7 Carefully tilt the weighing pan (7) and lift it off the weighing hooks (8).
- 8 Remove the drip tray (9).

- 9 Slightly lift the cover (10) and pull it towards the front to remove it.
 - 10 Store all removed components in a safe place.
- ➔ The balance is ready for cleaning.

5.2.2 Cleaning the balance

NOTICE

Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► www.mt.com/lab-cleaning-guide

Cleaning around the balance

- Remove any dirt or dust around the balance and avoid further contaminations.

Cleaning the terminal

- Clean the terminal with a damp cloth or a tissue and a mild cleaning agent.

Cleaning the removable parts

- Clean the removed part with a damp cloth or a tissue and a mild cleaning agent or clean in a dishwasher up to 80 °C.

Cleaning the weighing unit

- 1 Disconnect the balance from the AC/DC adapter.
- 2 Use a lint-free cloth moistened with a mild cleaning agent to clean the surface of the balance.
- 3 Remove powder or dust with a disposable tissue first.
- 4 Remove sticky substances with a damp lint-free cloth and a mild solvent, e.g., isopropanol or ethanol 70%.

5.2.3 Putting into operation after cleaning

- 1 Reassemble the balance.
- 2 Check that the draft shield doors and the lid open and close normally.
- 3 Check if the terminal is connected to the balance.
- 4 Reconnect the balance to the AC/DC adapter.
- 5 Check the level status, level the balance if necessary.
- 6 Respect the warm-up time specified in the "Technical Data".

- 7 Perform an internal adjustment.
- 8 Perform a routine test according to the internal regulations of your company. METTLER TOLEDO recommends performing a sensitivity test after cleaning the balance.
- 9 Press **→0←** to zero the balance.
 - ➔ The balance is ready to be used.

See also

- 🔗 Leveling the balance ▶ Page 17
- 🔗 Technical Data ▶ Page 23
- 🔗 Performing an internal adjustment ▶ Page 17

5.3 Service

Regular servicing by an authorized service technician ensures reliability for years to come. Contact your METTLER TOLEDO representative for details about the available service options.

6 Technical Data

6.1 General data

Power supply

AC/DC adapter (model no. FSP060-DHAN3):	Input: 100 – 240 V AC ± 10%, 50 – 60 Hz, 1.8 A Output: 12 V DC, 5 A, LPS, SELV
AC/DC adapter (model no. FSP060-DIBAN2):	Input: 100 – 240 V AC ± 10%, 50 – 60 Hz, 1.5 A Output: 12 V DC, 5 A, LPS, SELV
Cable for AC/DC adapter:	3-core, with country-specific plug
Balance power consumption:	12 V DC ± 10%, 2.25 A
Polarity:	⊖—●—⊕

Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Standards for safety and EMC:	See Declaration of Conformity
Range of application:	Use only indoors in dry locations

Environmental conditions

The limit values apply when the balance is used under the following environmental conditions:

Height above mean sea level:	Up to 5000 m
Ambient temperature:	+10 – +30 °C
Temperature change, max.:	5 °C/h
Relative humidity:	30 – 70%, non-condensing
Acclimatization time:	At least 12 hours after placing the instrument in the same location where it will be put into operation.
Warm-up time:	At least 120 minutes after connecting the balance to the power supply. When switched on from standby, the instrument is ready for operation immediately.

The balance can be used under the following environmental conditions. However, the weighing performances of the balance may be outside the limit values:

Ambient temperature:	+5 °C – +40 °C
Relative humidity:	20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

The balance can be disconnected and stored in its packaging under the following conditions:

Ambient temperature: -25 – +70 °C

Relative humidity: 10 – 90%, non-condensing

7 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

8 Compliance Information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ www.mt.com/ComplianceSearch

For further information, consult the Reference Manual (RM).

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► www.mt.com/contact

1.3 Akronymmer og forkortelser

Originalt udtryk	Oversat udtryk	Forklaring
AC		Alternating Current
ASTM		American Society for Testing and Materials
DC		Direct Current
EMC		Electromagnetic Compatibility
FCC		Federal Communications Commission
GWP		Good Weighing Practice
HID		Human Interaction Device
ID		Identification
LED		Light-Emitting Diode
LPS		Limited Power Source
MAC		Media Access Control
MT-SICS		METTLER TOLEDO Standard Interface Command Set
NA		Not Applicable

OIML	Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology)
RAM	Random Access Memory
RFID	Radio-frequency identification
RM	Reference Manual
SELV	Safety Extra Low Voltage
SOP	Standard Operating Procedure
SQC	Statistical Quality Control
UM	User Manual (Brugervejledning)
USB	Universal Serial Bus
USP	United States Pharmacopeia

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© 04/2024 METTLER TOLEDO. All rights reserved.