English

Reference Manual **Q3 Dosing Module** for XPR Analytical Balances

Deutsch

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---

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# Table of Contents

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Further documents and information .....	3
1.2	Explanation of conventions and symbols used .....	3
1.3	Acronyms and Abbreviations .....	4
1.4	Compliance information .....	4
<b>2</b>	<b>Safety Information</b>	<b>6</b>
2.1	Definitions of signal words and warning symbols .....	6
2.2	Product-specific safety notes .....	6
<b>3</b>	<b>Design and Function</b>	<b>10</b>
3.1	Function description .....	10
3.2	Overview .....	10
3.3	Components description .....	11
3.4	Dosing heads and equipment .....	13
3.5	Dosing head RFID tag .....	14
<b>4</b>	<b>Installation and Putting into Operation</b>	<b>15</b>
4.1	Selecting the location .....	15
4.2	Scope of delivery .....	15
4.3	Setting up the powder dosing configuration .....	16
4.4	Setting up the liquid dosing configuration .....	16
4.4.1	Overview .....	16
4.4.2	Interface of the pump .....	16
4.4.3	Wiring the pump .....	17
4.4.4	Connecting the tubes .....	17
<b>5</b>	<b>Operation</b>	<b>22</b>
5.1	Installing and removing the dosing head .....	22
5.2	Installing the sealing insert .....	23
5.3	Selecting and using an ErgoClip .....	23
5.4	Powder dosing .....	23
5.4.1	Using dosing heads without a vial .....	23
5.4.2	Filling the vial of a powder dosing head .....	24
5.4.3	Storing powder dosing heads .....	24
5.4.4	Avoiding electrostatic charges .....	24
5.5	Liquid dosing .....	25
5.5.1	Attaching the dosing head to the bottle cap .....	25
5.5.2	Handling the bottle .....	26
5.5.3	Using the QL3 pump .....	27
5.6	Using other weighing methods .....	27
<b>6</b>	<b>Maintenance</b>	<b>29</b>
6.1	Cleaning .....	29
6.1.1	Cleaning the housing .....	29
6.1.2	Cleaning the sealing insert .....	29
6.1.3	Cleaning the weighing chamber .....	30
6.1.4	Cleaning powder dosing heads .....	30
6.1.5	Cleaning the liquid dosing elements .....	30
6.2	Replacing the sealing ring and fastening nut on the liquid tube .....	31
<b>7</b>	<b>Troubleshooting</b>	<b>32</b>
7.1	Error symptoms .....	32
<b>8</b>	<b>Technical Data</b>	<b>34</b>
8.1	General data .....	34

8.2	Model-specific data.....	35
8.3	Material specifications.....	36
8.3.1	Glossary of material acronyms.....	36
8.3.2	Powder dosing head QH008.....	36
8.3.3	Powder dosing head QH012.....	37
8.3.4	Powder dosing head QH002-CNMW.....	38
8.3.5	Powder dosing head QH012-LNJW.....	39
8.3.6	Liquid dosing head QL001.....	40
8.3.7	Liquid dosing head QL003.....	41
8.3.8	QLL standard/advanced kit for bottles.....	41
8.4	Dimensions.....	42
<b>9</b>	<b>Disposal</b>	<b>44</b>
<b>10</b>	<b>Accessories and Spare Parts</b>	<b>45</b>
10.1	Accessories.....	45
10.2	Spare parts.....	52
10.2.1	Q3 dosing module.....	52
10.2.2	QL3 pump.....	52
10.2.3	QLL spare parts kit for bottles.....	53
	<b>Index</b>	<b>55</b>

# 1 Introduction

Thank you for choosing a METTLER TOLEDO instrument. The instrument combines high performance with ease of use.

## 1.1 Further documents and information

► [www.mt.com/XPR-automatic](http://www.mt.com/XPR-automatic)

This document is available in other languages online.

► [www.mt.com/Q3-RM](http://www.mt.com/Q3-RM)

Instructions for cleaning a balance: "8 Steps to a Clean Balance"

► [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

Search for software downloads

► [www.mt.com/labweighing-software-download](http://www.mt.com/labweighing-software-download)

Search for documents

► [www.mt.com/library](http://www.mt.com/library)

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

► [www.mt.com/contact](http://www.mt.com/contact)

## 1.2 Explanation of conventions and symbols used

### Conventions and symbols

Key and/or button designations and display texts are shown in graphic or bold text, e.g., , **Edit**.

#### Note

For useful information about the product.

Refers to an external document.

### Elements of instructions

In this manual, step-by-step instructions are presented as follows. The action steps are numbered and can contain prerequisites, intermediate results and results, as shown in the example. Sequences with less than two steps are not numbered.

- Prerequisites that must be fulfilled before the individual steps can be executed.
- 1 Step 1
  - ➔ Intermediate result
- 2 Step 2
  - ➔ Result

## 1.3 Acronyms and Abbreviations

Original term	Explanation
EMC	Electromagnetic Compatibility
FCC	Federal Communications Commission
LPS	Limited Power Source
POM	Polyoxymethylene
RFID	Radio-frequency identification
RM	Reference Manual
sd	Standard deviation
SELV	Safety Extra Low Voltage
SOP	Standard Operating Procedure
UM	User Manual
USB	Universal Serial Bus

## 1.4 Compliance information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

► <http://www.mt.com/ComplianceSearch>

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

► [www.mt.com/contact](http://www.mt.com/contact)

### United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

### Canada

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

### Singapore

Complies with  
IMDA Standards  
DA108390

### South Africa

### Thailand

 The image contains the nanb. logo at the top center, which features a red bird-like icon above the text 'nanb.'. Below the logo is a block of Thai text:
 

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## 2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.
- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

### 2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

#### Signal words

<b>DANGER</b>	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
<b>WARNING</b>	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
<b>CAUTION</b>	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
<b>NOTICE</b>	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

#### Warning symbols

General hazard

Notice

### 2.2 Product-specific safety notes

#### Intended use

This dosing system is designed to be used in analytical laboratories by trained staff. The dosing system is intended for weighing and dosing powder or liquid samples.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

#### Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

#### Protective equipment

Chemical-resistant gloves

Goggles

Lab coat

## Safety notes

### **WARNING**

#### **Death or serious injury due to electric shock**

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

### **WARNING**

#### **Injury and/or damage due to hazardous substances**

Chemical, biological or radioactive hazards can be associated with the substances processed by the instrument. During dosing procedures, small amounts of the dosed substance may become airborne and penetrate the instrument or contaminate its surroundings.

The substance characteristics and related hazards is the full responsibility of the instrument owner.

- 1 Be aware of possible hazards associated with the substance and take adequate safety measures, e.g., those stated on the safety data sheet provided by the manufacturer.
- 2 Make sure that every instrument part in contact with the substance will not get altered or damaged by the substance.

### **WARNING**

#### **Injury or damage due to powder handling**

Powders may get compacted in the dosing head and block it. Too much force applied by the dosing head mechanism may lead to a break and potentially hazardous substances may get airborne.

- 1 Handle the dosing heads with care.
- 2 If the dosing head seems blocked, avoid further dosing. Remove the head from the instrument and turn it upside down to loosen the powder.
- 3 Cease work immediately in the event of leak or break.

### **WARNING**

#### **Injury and/or damage due to reacting, flammable, or explosive substances**

During the dosing procedure, substances could be combined and cause an exothermic reaction or explosion. This includes powders, liquids, and gases.

The sample characteristics and related hazards is the full responsibility of the instrument owner.

- 1 Be aware of possible hazards associated with reacting, flammable, or explosive substances.
- 2 Ensure a working temperature low enough to prevent the formation of flames or an explosion.

### **⚠ WARNING**

#### **Injury or death due to toxic, explosive, or flammable substances**

If you use toxic, explosive, or flammable liquids with the pump, the exhaust air will be contaminated.

- Connect a tube to the exhaust air outlet to collect the contaminated air.

### **⚠ WARNING**

#### **Injury and/or damage due to reacting substances**

When pressure is released from the bottle, the air/gas in the bottle moves back towards the pump. The air/gas coming from the coupled outlets mixes in the pump. Molecules of the substances in the various bottles can get in contact through this contaminated air/gas.

- 1 Do not connect bottles with incompatible liquids to the same pump simultaneously.
- 2 Before connecting a second, incompatible liquid to the pump, disconnect the first bottle and purge the pump cavity with clean air/gas.

### **⚠ WARNING**

#### **Injury and damage to the pump or the bottle due to high pressure**

High pressure from external gas can damage the pump or the bottle.

- 1 Use a regulator on the external gas line.
- 2 Ensure that the pressure of the external gas does not exceed 0.5 bar (7.2 psi).

### **⚠ CAUTION**

#### **Injury due to splashing liquids**

If the pressure in the bottle is not released, liquid might splash when removing the micro dosing valve, opening the bottle, or removing the liquid tube.

- Always release pressure before removing the micro dosing valve, opening the bottle, or removing the liquid tube.

### **⚠ CAUTION**

#### **Injury due to leaking liquids**

Wrongly cut tubing can result in leaking connections.

- Cut the tubes with a tube cutter or a sharp knife.

### **⚠ CAUTION**

#### **Injury due to moving parts**

- Do not reach into the working area while parts of the instrument are moving.

### **⚠ CAUTION**

#### **Injury due to sharp objects or broken glass**

Instrument components, e.g., glass, can break and lead to injuries.

- Always proceed with focus and care.

**NOTICE****Damage to the instrument or malfunction due to the use of unsuitable parts**

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

**NOTICE****Damage to the instrument**

The instrument contains no user-serviceable parts.

- 1 Do not open the instrument.
- 2 In the event of problems, please contact a METTLER TOLEDO representative.

**NOTICE****Damage to the instrument due to inappropriate cleaning methods**

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

## 3 Design and Function


### 3.1 Function description

The **Q3 dosing module** can be added to any **XPR analytical balance** to perform automatic dosing of powders and/or liquids. The powder or liquid dosing head is attached to the dosing module and stays outside of the weighing chamber during the whole dosing procedure. The top door of the balance is replaced by a top panel with an opening through which the powder or liquid is dosed. The top panel moves down with the dosing module and the dosing head until the distance to the sample vessel is appropriate for dosing. The weighing chamber is therefore closed during the whole dosing procedure, providing the best possible dosing performances.

The **Q3 dosing module** is delivered with three ErgoClips: the **ErgoClip adapter holder**, the **ErgoClip adapter stabilizer**, and the **ErgoClip height extension**. The different designs can accommodate vessels of different sizes and shapes. All three ErgoClips are compatible with XPR analytical balances with tall draft shield.

When liquid is dosed, the **QL3 pump** and **QLL kit for bottle** are used together with the **Q3 dosing module**. The **QL3 pump** is used to build pressure in the bottle. As soon as the pressure is sufficiently high, the micro dispensing valve in the liquid dosing head opens and liquid ascends the liquid tube. Several pumps can be connected to the balance simultaneously and each pump can accommodate up to three bottles.

### 3.2 Overview


1	Dosing lift	6	Powder dosing head
2	Dosing module	7	Sealing insert
3	Top panel	8	Dosing head holder
4	Vial adapter	9	Dosing head release button
5	ErgoClip adapter holder		

### 3.3 Components description


#### Dosing module

The dosing module is designed to hold the dosing head and can automatically move up and down during the dosing procedure.

#### Dosing lift


The dosing lift is fixed to the back of the weighing chamber. It contains the mechanism that moves the dosing module up and down during the dosing procedure and is controlled through the terminal.

#### Top panel


The Q3 dosing module includes a top panel especially designed to minimize the influence of the environment on the dosing process. An opening in the top panel allows the tip of the dosing head to penetrate the weighing chamber while dosing. The top panel is attached to the bottom part of the dosing module, moving up and down together with it.

While the dosing module is at its uppermost position, other weighing methods and accessories can be used normally.


#### Dosing head holder

The dosing head holder consists of four pins that hold the dosing head in place. Once installed, the dosing head is secured and can only be released by pressing one of the release buttons on either side of the holder, while the dosing module is at its uppermost position. The holder includes an RFID reader to identify the attached dosing head.

### Optical sensor for HeightDetect


The dosing module is equipped with an optical sensor to detect the height of the sample vessel opening (Height-Detect). This allows automatic vertical positioning of the dosing module for fast and safe dosing.

### ErgoClip adapter holder and adapters


The Q3 dosing module includes the ErgoClip adapter holder, which can be installed on the SmartGrid weighing pan. It can hold adapters of various sizes allowing perfect positioning of different sample vessels. Additional adapters are available as accessories. See [Accessories ▶ Page 45].

The Q3 dosing module is delivered with a variety of ErgoClips for various usage. See [Selecting and using an ErgoClip ▶ Page 23]


### ErgoClip adapter stabilizer

The Q3 dosing module includes the ErgoClip adapter stabilizer, which can be installed on the SmartGrid weighing pan. The ErgoClip adapter stabilizer is especially designed to provide more stability to the adapter. It can therefore be used for tall and narrow vessels.

### ErgoClip height extension

The Q3 dosing module includes the ErgoClip height extension, which can be installed on the SmartGrid weighing pan. This ErgoClip has a flat surface and can therefore accommodate any sample vessel. It is especially useful for short sample vessels, such as weighing boats.

### 3.4 Dosing heads and equipment

#### Powder dosing head

Powder dosing heads are used for automated powder dosing. If a small quantity of powder is required, the dosing head can be used without a vial attached to it. In this case, the opening can then be closed using the cap included in the delivery. The dosing heads can be stored in appropriate storage containers. Vials and storage containers need to be purchased separately.

If a printer is connected, you can use a pre-defined template to print a label with the dosing-head data from the RFID tag and affix this label to the dosing head.

#### Powder test head




The powder test head is used to test a powder dosing system. The dosing head is filled with calcium carbonate ( $\text{CaCO}_3$ ) and uses a predefined test program to check the dosing functionality of the system. The powder test head is programmed to perform 15 tests. Each test weighs 10 doses of a pre-defined amount of powder. If your instrument fails the tests, inform a METTLER TOLEDO service technician.

#### Weight test head

The weight test head performs an automated repeatability test by placing a small test weight on given tare weights. The repeatability of the system is determined based on 10 consecutive measurements of the test weight. If your instrument fails the tests, inform a METTLER TOLEDO service technician.

#### Liquid dosing head

Liquid dosing heads are used for automated liquid dosing. They are used together with a pump and bottle. The liquid dosing head contains an RFID tag which stores information about the dosing head and the used substance.

### 3.5 Dosing head RFID tag

Each powder and liquid dosing head is equipped with an integrated RFID tag (1) that stores and exchanges data with the instrument.


Various data is stored in the RFID tag of the dosing heads, such as the substance name, lot ID, filling date, expiry date, etc. It also includes customized data fields.

This data can be edited on the terminal and should be set before using a new dosing head such that the data is available for reports and labels.

For powder dosing heads, the counter for the remaining powder dosing cycles is based on the **Dose limit** of the dosing head, which can be defined by the customer. With every dosing cycle started, the counter is decreased by 1.

If the counter drops to zero, consider replacing the dosing head. If the vial of the old dosing head still contains a considerable amount of powder, you may remove the vial from the old dosing head and screw it onto the new dosing head. Copy the user data and the powder content value from the old dosing head to the new dosing head.

The quantity of powder added to the dosing head can be stored in the RFID tag during the filling procedure. The RFID information is updated after each dosing operation such that the powder quantity in the dosing head is up to date. If the remaining quantity of powder is insufficient for the next dosing cycle, a warning message will appear on the terminal.

For further information, consult the Reference Manual (RM) of your XPR balance.

► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

## 4 Installation and Putting into Operation

This instrument must be installed by a METTLER TOLEDO service technician.  
The Q3 dosing module is compatible with all XPR analytical balances.

### 4.1 Selecting the location


A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.

#### Requirements of the location


Place indoors on stable table


Avoid direct sunlight

Ensure sufficient spacing

Avoid vibrations


Level the instrument

Avoid strong drafts

Provide adequate lighting

Avoid temperature fluctuations

Sufficient spacing for balances: > 15 cm all around the instrument

Take into account the environmental conditions. See "Technical Data".

#### Note

The powder delivery mechanism of the Q3 dosing module might cause the instrument to vibrate. Do not place on the same working surface as other instruments that are sensitive to vibrations.

### 4.2 Scope of delivery

#### Q3 dosing module

- Dosing module
- Dosing lift
- Dosing lift back cover with screws
- Top panel dosing module
- Sealing inserts, 5 pcs
- Powder dosing head
- ErgoClip adapter holder
- ErgoClip adapter stabilizer
- ErgoClip height extension
- Various vial adapters, 4 pcs
- Declaration of Conformity
- User Manual

#### Recommended options

- EasyHub USB
- Integrable antistatic kit


The scopes of delivery for the **QL3 pump** and the **QLL kits for bottles** are in the documentation delivered with these products.

### 4.3 Setting up the powder dosing configuration

After the Q3 dosing module has been installed on your XPR analytical balance by a METTLER TOLEDO service technician, simply install a powder dosing head and start a **Automated dosing** or **Automated solution prep.** method.


### 4.4 Setting up the liquid dosing configuration

#### 4.4.1 Overview

1	Dosing lift	7	ErgoClip adapter holder
2	Dosing module	8	QL3 pump (available separately)
3	Liquid dosing head	9	QLL kit for bottles (available separately)
4	Sealing insert	10	Air tube (included in QLL kit for bottles)
5	Top panel	11	Liquid tube (included in QLL kit for bottles)
6	Vial adapter		

#### 4.4.2 Interface of the pump

1	Air outlets (to bottle)	4	Socket for AC/DC adapter
2	Pressure release button and indicator light	5	USB-B port (to host)
3	Air exhaust outlet	6	Air inlet

### 4.4.3 Wiring the pump


- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
- 2 Insert the plug of the AC/DC adapter in the power inlet of the pump.
- 3 Secure the plug by firmly tightening the knurled nut.
- 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
- 5 Using the USB cable, connect the USB-B port of the pump to one of the USB-A ports of the balance.

### 4.4.4 Connecting the tubes

#### Tubes definition

The liquid tube is the thinner tube (1) used for transporting liquid from the bottle to the liquid dosing head. The air tube is the slightly bigger tube (2) used for pumping air into the bottle. By adding air through the air tube, pressure rises in the bottle. When the pressure reaches the target pressure, between 0.3 and 0.5 bar (4.4 to 7.2 psi), the micro dispensing valve in the dosing head opens and liquid can ascend the liquid tube.

- 1 Liquid tube
- 2 Air tube

#### Preparing the liquid tube

#### CAUTION

##### Injury due to leaking liquids

Wrongly cut tubing can result in leaking connections.

- Cut the tubes with a tube cutter or a sharp knife.

- The liquid dosing head is inserted in the liquid dosing head support of the bottle.

- 1 Using a tube cutter or a sharp knife, cut a sufficient amount of tubing. The appropriate length mainly depends on the distance between the balance and the bottle while dosing.


Recommended length: about 0.9 m

- 2 **NOTICE: Leaking liquid due to wrong assembly. Pay attention to the orientation of the sealing ring when threading it on the tube.**

Place the sealing ring (1) on a flat stable surface, e.g., a table or a workbench, with the wider end downwards.

- 3 Take the end of the liquid tube and press it into the sealing ring.

- ➔ This is the dosing-head end of the tube. The opposite end is the bottle end.


- 4 Thread the two fastening nuts (2), paying attention to the orientation.
- 5 **NOTICE: Leaking liquid due to wrong assembly. Pay attention to the orientation of the sealing ring when threading it on the tube.**  
Thread the sealing ring (3) from the bottle end of the tube.
- 6 Slide the sealing ring until the distance (D) is enough for the tube to reach the bottom of the bottle. Recommended distances for typical bottle volumes are listed next.

**Typical distances (D) between the sealing ring and the bottle-end of the tube**

Thread	Bottle volume	Distance (D)
GL45	1000 ml	220 mm
	500 ml	170 mm
	250 ml	135 mm
GL25	25 ml	90 mm

**Connecting the liquid tube**

- 1 Insert the dosing-head end of the tube (1) in the dosing head (2).
- 2 Tightly fasten the fastening nut to the dosing head.
- 3 Insert the bottle end of the tube through the corresponding hole in the bottle cap (3). The tube should reach the bottom of the bottle.
- 4 If needed, attach the suction filter (4) to the bottle end of the tube.
- 5 Tightly fasten the fastening nut to the bottle cap.
- 6 Screw the cap to the bottle (5).

The suction filter is used to ensure that no particles or impurities are carried through the liquid dosing head. Using the suction filter will extend the lifetime of the dosing head. However, when dosing solutions, molecules of one of the substances might be absorbed by the suction filter, altering the concentration of the solution. The suction filter should only be used when dosing pure solvents.

**Connecting the air tube**


**WARNING**

**Injury and/or damage due to reacting substances**

When pressure is released from the bottle, the air/gas in the bottle moves back towards the pump. The air/gas coming from the coupled outlets mixes in the pump. Molecules of the substances in the various bottles can get in contact through this contaminated air/gas.


- 1 Do not connect bottles with incompatible liquids to the same pump simultaneously.
- 2 Before connecting a second, incompatible liquid to the pump, disconnect the first bottle and purge the pump cavity with clean air/gas.

- 1 Using a tube cutter or a sharp knife, cut a sufficient amount of tubing. The appropriate length mainly depends on the distance between the bottle and the pump while dosing.  
Recommended length: about 0.7 m
- 2 Connect the air tube (1) to the air inlet of the bottle (2)


- 3 Insert the other end of the air tube (1) in a tube adapter (3) and press firmly.
- 4 Connect the tube adapter (3) to one of the air outlets (4) of the pump. Press until you hear a click.

When a tube is connected to the air outlet of the pump, the valve of the air outlet opens. Never leave a tube that is connected to the air outlet unconnected at the other end because pressure cannot build up.  
You can connect up to three bottles to each pump.

### Connecting the muffler

- Insert the muffler (1) into the air inlet (2) to absorb the noise.

### Removing the air tube


#### NOTICE

##### Damage to tube connectors due to mishandling


If the tubes are not removed correctly, the connectors of the pump and the bottle cap can be damaged.

- To remove the tubes, press down the ring on the connector and pull out the tube carefully.

- 1 Remove the tube adapter (1) from the pump by pushing the connector (2) of the air outlet towards the pump.  
➔ The adapter is released and can be removed.

- 2 Remove the air tube from the bottle by firmly pressing down the ring (3) and pulling the tube (4) simultaneously.
- 3 If necessary, e.g., for maintenance purposes, use the same procedure to remove the air tube from the tube adapter: firmly press the ring and pull the tube simultaneously.

### Using the pump with external gas

The liquid can be protected by feeding an external gas, e.g., nitrogen, to the pump. Make sure that the pressure of the external gas does not exceed 0.5 bar (7.2 psi).

### **WARNING**

#### **Injury and damage to the pump or the bottle due to high pressure**

High pressure from external gas can damage the pump or the bottle.


- 1 Use a regulator on the external gas line.
- 2 Ensure that the pressure of the external gas does not exceed 0.5 bar (7.2 psi).

- A regulator is connected to the external gas line.
- 1 Remove the muffler from the air inlet (1).
  - 2 Connect the tube of external gas (2) to the air inlet (1).

#### **Note**

Outer tube diameter: 6 mm

Pressure in external gas line: The pressure must be at least 0.1 bar (1.5 psi). The pressure must not exceed the configured dosing pressure set on the terminal. The external gas tube is not provided by METTLER TOLEDO.

## Collecting contaminated air

### **WARNING**

#### **Injury or death due to toxic, explosive, or flammable substances**

If you use toxic, explosive, or flammable liquids with the pump, the exhaust air will be contaminated.


- Connect a tube to the exhaust air outlet to collect the contaminated air.

- Connect a tube (1) to the exhaust air outlet (2) to collect the contaminated air into a safe container (3).

#### **Note**

Outer tube diameter: 6 mm

The exhaust air tube and the container are not provided by METTLER TOLEDO.

## 5 Operation

For further information, consult the Reference Manual (RM) of your XPR balance.


► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

### 5.1 Installing and removing the dosing head

When the dosing module is in its uppermost position (**Home position**), the dosing head is unlocked. That means that it can be released by pressing the release button of the holder. During dosing, the dosing head is locked and cannot be removed until it returns to its uppermost position.


#### Installing the dosing head

- The dosing module is in its uppermost position.
- 1 Slide the dosing head (1) onto the dosing head holder (2) until it comes to a stop.
- 2 Slightly press down the dosing head until it properly sits in the dosing head holder (2).
  - ➔ The dosing head is locked in the holder.
- ➔ The dosing head is ready to dose.


#### Removing the dosing head

- The dosing module is in its uppermost position.
- 1 Press the release button (1) on either side of the dosing head holder to release the dosing head (2).
- 2 Simultaneously, pull the dosing head (2) up and outwards.

## 5.2 Installing the sealing insert


Always install a sealing insert (1) to avoid air drafts in the weighing chamber.

## 5.3 Selecting and using an ErgoClip

While dosing in vials, capsules, or tubes, specific adapters can be used. To ensure stability and ease the alignment of the sample vessel with the tip of the dosing head, three ErgoClips are delivered with the Q3 dosing module.

- When using sample vessels taller than 75 mm, place directly on the weighing pan or use the ErgoClip adapter holder (1) with an adapter. The nose on the ErgoClip adapter holder is designed to make sure that the adapter (and therefore the vessel opening) is centered on the weighing pan.
- When using sample vessels taller than 75 mm that are very narrow and need more stability, use the ErgoClip adapter stabilizer (2) with an adapter.
- When using sample vessels shorter than 75 mm, use the taller ErgoClip height extension (3). The surface of this ErgoClip is completely flat, allowing to use vessels without adapters, for example, weighing boats.


The ErgoClips delivered with the Q3 dosing module are compatible with all XPR analytical balances with tall draft shield. If using an XPR analytical balance with short draft shield, place your sample vessel and/or adapter directly on the weighing pan. For more stability, use the ErgoClip vial, available separately. See [Accessories ▶ Page 45].

## 5.4 Powder dosing

### 5.4.1 Using dosing heads without a vial

Most powder dosing heads can be used with or without a vial, since the neck of the dosing head can accommodate a certain amount of powder. Using dosing heads without vials is particularly relevant when dosing very small quantities. In this case, the dosing head needs to be closed using the provided spare cap (1). The neck (2) of the QH002, QH008, QH010, and QH012 dosing heads have a capacity of 2 ml, 8 ml, 10 ml, and 12 ml, respectively.

## 5.4.2 Filling the vial of a powder dosing head

Software functions to easily fill or refill the dosing head are available through the terminal of the balance. Consult the Reference Manual of the balance for more information.

- A method **General weighing** is running on the terminal.
- The dosing head is removed.
  - 1 Turn the dosing head upside down.
  - 2 If the attached vial already contained powder, gently tap on it to let the remaining powder flow back into the vial.
  - 3 **⚠ WARNING: Injury or contamination from hazardous substances. Take the appropriate measures to avoid injuries and contamination.**  
Unscrew the vial.
  - 4 Place the vial on the weighing pan.
  - 5 Press **→T←** to tare the balance.
  - 6 Fill the vial with the desired amount of the appropriate substance.
  - 7 Note down the weight of added substance.
  - 8 Remove the vial from the weighing pan.
  - 9 Holding the vial upright and the dosing head upside down, screw the dosing head onto the vial.
  - 10 Turn the dosing head upright and gently tap on it to let the powder flow towards the dosing head.
  - 11 Install the dosing head on the dosing module.
  - 12 On the terminal, edit the settings of the dosing head to include the added amount of powder.

## 5.4.3 Storing powder dosing heads


For short-term and long-term storage of the powder dosing heads, METTLER TOLEDO recommends using storage containers to:

- reduce the risk of contamination
- reduce the quantity of humidity absorbed by the powder

You can store the dosing heads right-side up (1) or upside-down (2).

### See also

📖 Accessories ▶ Page 45

## 5.4.4 Avoiding electrostatic charges

When preparing sample vessels, electrostatic charge may build up. Electrostatic charges may make correct dosing impossible or adversely affect the dosing result.

Factors increasing the amount of electric charges:

- using sample vessels made of plastic
- wearing latex gloves

The integrable antistatic kit is used to remove electrostatic charges from sample vessels. This antistatic kit is especially designed for dosing purposes and consists of two ionizing electrodes, one placed on each side of the weighing chamber. The electrodes can be activated automatically when dosing or manually. Using those two electrodes helps removing electrostatic charges from sample vessels.


Make sure to install a dosing head before placing the sample vessel on the weighing pan. This way, ionizing is active while placing the sample vessel on the weighing pan and electric charges are neutralized. In addition, try not to touch the upper edge of the sample vessel, e.g., close to a vial opening, when handling it.


## 5.5 Liquid dosing

### 5.5.1 Attaching the dosing head to the bottle cap

- 1 Insert the liquid dosing head in the liquid dosing head support.


- 2 To remove the liquid dosing head from the liquid dosing head support, pull the lever (1) towards the dosing head and remove the liquid dosing head (2).


## 5.5.2 Handling the bottle

### Filling the bottle

The suction filter (1) has to be covered with liquid at all times. Before the suction filter gets dry, refill the bottle.

- The pressure is released.
- 1 **⚠ CAUTION: Injury due to splashing liquids. Make sure that the pressure in the bottle is released.**  
Unscrew the cap.
- 2 Fill in the liquid. Do not exceed the maximum (max. is shown on bottle, e.g., 1000 ml). The air above the liquid is necessary for dosing.
- 3 Screw the cap on tightly.

### Changing the bottle content


This procedure is used if you have a single cap with dosing head and want to change the liquid to be dosed.

- The pressure is released.
- 1 **⚠ CAUTION: Injury due to splashing liquids. Make sure that the pressure in the bottle is released.**  
Unscrew the cap.
- 2 Remove the suction filter, if applicable.
- 3 If the cap needs to be cleaned, remove the air tube.  
Unscrew the fastening nut from the cap.  
Rinse the cap with the appropriate solvent or liquid.  
Insert the liquid tube into the cap.
- 4 If the liquid tube needs to be cleaned with a solvent, fill the bottle with the appropriate solvent.  
Screw the cap on the bottle.  
Insert the air tube on the cap.  
Purge using the **Purge** function on the terminal.  
Unscrew the cap.  
Dispose of the remainder of solvent.
- 5 Attach a new suction filter, if applicable.
- 6 Screw the cap on the bottle containing the new liquid.
- 7 Check that the cap is tight.
- 8 Connect the air tube to the new bottle.
- 9 Purge using the **Purge** function.

### Connecting the air tube to a different bottle

If you have more than one bottle equipped with cap and dosing head, and you want to use the same pump outlet and air tube to dose from another bottle:

- The pressure is released.
- 1 Install the dosing head on the dosing head support of the bottle.
- 2 Unplug the air tube by pressing down the ring (1) and simultaneously pulling out the tube (2).
- 3 To seal the bottle, insert the pin delivered with the QLL kit (3) into the air tube connector.
- 4 Take the new bottle.
- 5 Connect the air tube to the new bottle.
- 6 To continue dosing with the new bottle, install the dosing head.

### 5.5.3 Using the QL3 pump


#### Pressure indicator light

The pressure indicator light displays the status of the pump:

- Light on: pressure is building up or pressure is established
- Light off: there is no pressure and no pressure is building up
- Light blinking: pump error and/or warning

#### Releasing the pressure

- The pressure indicator light (1) is on.
- Press the pressure release button (2) to release the pressure.
- ➔ The status light (1) turns off when pressure is released.

#### Purging the pump cavity

If bottles containing incompatible liquids (and for which the fumes should not mix) are connected to the pump one after the other, it is recommended to purge the pump cavity before connecting the second bottle to the pump.

- A dosing head is attached to the dosing or liquid module. The bottle of this dosing head is not connected to any pump.
- A method **Automated dosing** or **Automated solution prep.** is running on the terminal.
- 1 Disconnect all tube adapters from the pump.
- 2 Connect an empty tube adapter to the right-most air outlet at the front of the pump.
  - ➔ The pump is trying to build pressure and air flows through the cavity of the pump, purging it.
  - ➔ The pump cavity has been purged and bottles can safely be re-connected to the air outlets.

### 5.6 Using other weighing methods

To use other methods with the balance, move the dosing module to its uppermost position.

For further information, consult the Reference Manual (RM) of your XPR balance.

▶ [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

## 6 Maintenance

The appropriate maintenance interval depends on your standard operating procedure (SOP).

Please contact your METTLER TOLEDO representative for details about the available service options. Regular servicing by an authorized service technician ensures constant accuracy for years to come and prolongs the service life of your instrument.

### 6.1 Cleaning

#### NOTICE

##### Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

Detailed information about the compatibility of cleaning agents can be found in the Reference Manual (RM) of your XPR balance.

► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)


#### 6.1.1 Cleaning the housing

The housing material of the dosing module and the dosing lift is the same as the one of the balance. All surfaces can therefore be cleaned with a commercially available, mild cleaning agent.

#### 6.1.2 Cleaning the sealing insert


The sealing insert is located right under the tip of the dosing head. It is therefore the element that is the most likely to be contaminated by the substances dosed with the dosing module. Clean or replace the sealing insert regularly, depending on the substances used and your SOP (Standard Operating Procedure).

- The dosing module is in its uppermost position.
  - 1 Remove the sealing insert (1).
  - 2 Clean the sealing insert with a wet cloth and a mild cleaning agent or dispose of it. Do not wash the sealing insert in the dishwasher.
  - 3 Install the cleaned sealing insert or a new one.

### 6.1.3 Cleaning the weighing chamber

- The dosing module is in its uppermost position.
  - The dosing head is removed.
  - The sealing insert is removed.
  - The draft shield is removed.
- 1 Tilt the top panel (1) upwards and remove.
  - 2 Remove any adapter holder or ErgoClip, the weighing pan, and the drip tray.
  - 3 Clean the weighing chamber. Wipe the covering strip (2) gently using vertical movements.
    - ⚠ **WARNING: Damage to the covering strip. The covering strip is delicate and must be handled with care.**
  - 4 Reassemble all parts in reverse order.
  - 5 Put the balance into operation.

For more details about cleaning and putting into operation after cleaning, see "Maintenance" in the Reference Manual of the balance.

### 6.1.4 Cleaning powder dosing heads

Powder dosing heads are designed to be used with a unique substance. METTLER TOLEDO does not recommend cleaning them. When dosing a different powder, use a new dosing head to avoid contamination.

### 6.1.5 Cleaning the liquid dosing elements

#### Cleaning the suction filter

- 1 **Once a week**, check visually that the suction filter is clean. If necessary, rinse the suction filter with solvent using the **Purge** function on the terminal or replace the suction filter.
- 2 At least **once a year**, change suction filter. Maintenance interval depends on the liquid used.

#### Purging the liquid dosing head

The liquid dosing head can be rinsed/purged by letting a large quantity of solvent (or another liquid) run through it. Use the **Purge** function. See the Reference Manual of the XPR balance.

- The bottle is filled with enough solvent to purge the dosing head.
- A big sample vessel is on the weighing pan to collect the solvent used to purge the dosing head.
- The liquid dosing head is installed on the dosing or liquid module.
- Use the **Purge** function on the terminal to rinse the dosing head.

#### Cleaning the micro dosing valve


The liquid dosing heads QL001 have a micro dosing valve that can be disassembled and cleaned, for example, in an ultrasonic bath.

#### Note

The QL003 dosing head does not contain a micro dosing valve. Cleaning this dosing head using the **Purge** function is sufficient (see above).

- The liquid dosing head is installed on the dosing or liquid module.
  - A sample vessel is on the weighing pan, big enough for the quantity of liquid in the liquid tube.
  - The pressure is released.
- 1 ⚠ **CAUTION: Injury due to splashing liquids. Make sure that the pressure in the bottle is released.** Empty the liquid tube by replacing the bottle by an empty bottle and using the **Purge** function.
    - ➔ The liquid tube is empty.


- 2 Release pressure.
- 3 **⚠ CAUTION: Injury due to splashing liquids. Make sure that the pressure in the bottle is released.**  
Remove the dosing head from the dosing or liquid module.
- 4 **⚠ WARNING: Injury or contamination from hazardous substances. Be aware of any liquid draining from the tubing and the valve.**  
Open the dosing head with the micro dosing valve tool (1), if applicable.
- 5 Take out the micro dosing valve and clean it, e.g., in an ultrasonic bath.
- 6 After the cleaning, reinstall the micro dosing valve and refill/replace the bottle.


## 6.2 Replacing the sealing ring and fastening nut on the liquid tube

### Removing the sealing ring and the fastening nut

If liquid is leaking from the dosing head nut, replace the sealing ring and fastening nut at the dosing head (1). If pressure cannot be built in the bottle, replace the sealing ring and fastening nut at the bottle cap (2).

- The pressure is released.
- 1 **⚠ CAUTION: Injury due to splashing liquids. Make sure that the pressure in the bottle is released.**  
Unscrew the fastening nut at the dosing head or the bottle cap.
  - 2 To access the sealing ring (1), slide back the fastening nut (2).
  - 3 If changing the nut and ring at the bottle cap, remove the suction filter and slide the sealing ring all the way to the bottle-end of the tube.
  - 4 Using a tube cutter or a sharp knife, cut off the tube above the sealing ring (1).
  - 5 Remove the fastening nut.

### Reconnecting the liquid tube

Insert the new fastening nut sealing ring, and reconnect the liquid tube.

#### See also

Connecting the tubes ▶ Page 17

## 7 Troubleshooting

For further information, consult the Reference Manual (RM) of your XPR balance.

► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

### 7.1 Error symptoms

Error symptom	Possible cause	Diagnostic	Remedy
The powder dosing head does not dispense any powder.	The powder does not reach the dispensing mechanism of the dosing head.	If the dosing head is transparent, visually inspect through the housing if the powder is clogged up.	Remove the dosing head and shake it to loosen the powder.
The liquid dosing head does not dispense any liquid.	The suction filter is clogged.	Check if there is liquid in the liquid tube.	Clean or replace the suction filter.
	The dosing head is not installed properly.	–	Remove and reinstall the dosing head in the holder. Make sure to press until you hear a click.
Liquid leaks from the bottle cap and/or the dosing head.	The tube is not properly connected.	–	Tighten the fastening nut at the bottle cap and/or the dosing head. Make sure that the end of the tube is cut straight.
	The sealing ring and/or fastening nut are damaged.	–	Change the sealing ring and the fastening nut at the bottle cap and/or the dosing head. See "Maintenance". Make sure that the end of the tube is cut straight.
Liquid drips from the QL003 liquid dosing head.	Impurities have entered the QL003 liquid dosing head.	–	Use the <b>Purge</b> function for at least 10 seconds to rinse the dosing head. Check if the dosing head is still dripping. Repeat if necessary. If the dosing head cannot be cleaned properly, replace it.
The status light of the pump does not blink when the device is connected to power.	The pump is disconnected from power.	Disconnect and reconnect the AC/DC adapter from the pump. The indicator light should blink once when the device is connected to power. Check that the AC/DC adapter and the power cable are not damaged.	Replace the AC/DC adapter and power cable.
	The pump is damaged.	Check with another pump, if available.	Replace the pump.

Error symptom	Possible cause	Diagnostic	Remedy
			Contact your METTLER TOLEDO service representative.
On the display, the QL3 pump does not appear in the list of devices connected to the balance.	The pump is disconnected from power.	Disconnect and reconnect the AC/DC adapter from the pump. The indicator light should blink once when the device is connected to power. Check that the AC/DC adapter and the power cable are not damaged.	Replace the AC/DC adapter and power cable.
	The USB cable is not properly connected.	Check that the USB cable is connected properly.	Connect the USB cable properly.
	The USB cable is damaged.	Check that the USB cable is not damaged.	Replace the USB cable.
	The USB-A port of the balance is damaged.	Disconnect the pump from the USB-A port of the balance. Connect a USB mouse to the same USB-A port. Verify that a pointer (arrow) appears on the terminal and can be moved by moving the mouse.	If the mouse pointer does not appear, contact your METTLER TOLEDO service representative.
	The pump is damaged.	Check with another pump, if available.	Replace the pump. Contact your METTLER TOLEDO service representative.
The value on the display is unstable when dosing powder.	The weighing sample is electrostatically charged.  <b>Note</b> Additional possible causes can be found in the Reference Manual of your XPR balance.	Check if the weighing result is stable when using a test weight.	Increase the air humidity in the weighing chamber. Use an antistatic kit. See "Accessories" in the Reference Manual.

## 8 Technical Data

### 8.1 General data

Dosing lift, weight: 1250 g

Dosing module, weight: 450 g

#### Power consumption

Dosing lift: 12 V DC  $\pm$  6%, 1 A

Dosing module: 12 V DC  $\pm$  6%, 1 A

#### Protection and standards

Overvoltage category: II

Degree of pollution: 2

Range of application: Use only indoors in dry locations

#### Environmental conditions

Height above mean sea level: Up to 5000 m

Ambient temperature: +5 – +40 °C

Relative air humidity: 20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

#### Storage conditions (in packaging)

Ambient temperature: -25 – +70 °C

Relative air humidity: 10 – 90%, non-condensing

## 8.2 Model-specific data

### Liquid dosing heads

	QL001	QL003
<b>Dosing properties</b>		
Suitable liquids: max. viscosity <sup>1)</sup>	20 mPa·s	20 mPa·s
<b>Typical values<sup>2)</sup></b>		
Repeatability (sd)	1 mg	2 mg
Dosing time <sup>3)</sup>	40 s	25 s

<sup>1)</sup> 1 mPa·s = 1 cP = 1 cSt = 1 mm<sup>2</sup>/s

<sup>2)</sup> for H<sub>2</sub>O

<sup>3)</sup> for 10 g

### Tubing

	Outer diameter	Inner diameter	Recommended length
Liquid tube, for GL45 bottles (dosing head to bottle)	3.2 mm	1.6 mm	0.9 m
Liquid tube, for GL25 bottles (dosing head to bottle)	1.6 mm	0.8 mm	0.9 m
Air tube (pump to bottle)	4.0 mm	2.4 mm	0.7 m
Tube for exhaust air	6 mm	–	–
Tube for external gas	6 mm	–	–

### Suction filter

Tubing outer diameter	3.2 mm
Filter pore size	10 µm

### Bottle

Pressure resistance, min. <sup>1)</sup>	1.5 bar
Burst pressure, min.	3 bar
Volume, max.	2 l

<sup>1)</sup> According to DIN EN 1595: Pressure Equipment made from Borosilicate Glass 3.3 – General Rules for Design, Manufacture and Testing

### QL3 pump

Maximum pressure <sup>1)</sup>	1.5 bar
--------------------------------	---------

<sup>1)</sup> In the unlikely event of a malfunction, the overpressure valve opens if the maximum pressure is reached, to preserve the integrity of the system.

### 8.3 Material specifications

Materials in contact with the substance to be dosed. Please note that traces of all contact material could be transferred to the sample.


#### 8.3.1 Glossary of material acronyms

Acronyms defined in ISO 1043: Plastics – Symbols and abbreviated terms.

<b>ETFE</b>	=	Ethylene tetrafluoroethylene
<b>FEP</b>	=	Fluorinated ethylene propylene
<b>FFKM</b>	=	Perfluoroelastomer
<b>PE</b>	=	Polyethylene
<b>PEEK</b>	=	Polyetheretherketone
<b>PMP</b>	=	Polymethylpentene
<b>POM</b>	=	Polyoxymethylene
<b>PP</b>	=	Polypropylene
<b>PP (ESD)</b>	=	Polypropylene (electrostatic dissipative)
<b>PTFE</b>	=	Polytetrafluoroethylene


#### 8.3.2 Powder dosing head QH008

	<b>Designation</b>	<b>Material</b>
<b>1</b>	Dosing head body	PMP
<b>2</b>	Bushing	POM
<b>3</b>	Dosing pin	Stainless steel 1.4404
<b>4</b>	Scraper	Stainless steel 1.4310
<b>5</b>	Cone	Stainless steel 1.4404
<b>6</b>	Insert	POM
<b>7</b>	Adapter cap	POM


### 8.3.3 Powder dosing head QH012

	Designation	Material
1	Dosing head body	PMP
2	Bushing	POM
3	Dosing pin	Stainless steel 1.4404
4	Scraper	Stainless steel 1.4310
5	Cone	Stainless steel 1.4404
6	Insert	POM
7	Adapter cap	POM


### 8.3.4 Powder dosing head QH002-CNMW

	Designation	Material
1	Body	PP
2	Closing cap	PP
3	Dosing pin	PP (ESD)
4	Cone body	PP (ESD)


### 8.3.5 Powder dosing head QH012-LNJW

	Designation	Material
1	Adapter cap	POM
2	Dosing head body	PMP
3	Dosing pin	PP (ESD)
4	Cone body	PP (ESD)
5	Insert	POM

### 8.3.6 Liquid dosing head QL001


	Designation	Material
1	Fastening nut	PEEK
2	Sealing ring	ETFE
3	Tube connector	Stainless steel 1.4404
4	O-Ring	FFKM
5	Solenoid valve	PEEK Sapphire-Ruby Stainless steel 1.4105IL Stainless steel 1.14301 / 1.4306 Stainless steel 1.4305

### 8.3.7 Liquid dosing head QL003

	Designation	Material
1	Fastening nut	PEEK
2	Sealing ring	ETFE
3	Nozzle	Ceramic Al <sub>2</sub> O <sub>3</sub> > 99.7%
4	Valve block	PTFE
5	Gasket	EPDM
6	Valve	PEEK


### 8.3.8 QLL standard/advanced kit for bottles

**Note**

The QLL kits for bottle include a liquid dosing head, see [Liquid dosing head QL001 ▶ Page 40], [Liquid dosing head QL003 ▶ Page 41].


	Designation	Material
1	Cap insert	PE
2	Bottle	Borosilicate glass
3	Tubing	FEP
4	Suction filter	PP

## 8.4 Dimensions

### Q3 dosing module and XPR balance

Dimensions in mm.


Dimensions	Q3 on XPR with tall draft shield	Q3 on XPR with half-height draft shield
Maximum height of sample vessel, without adapter <sup>1)</sup>	212 mm	135 mm
Minimum opening of sample vessel (diameter)	6 mm	6 mm
Vertical range of movement of dosing module	150 mm	73 mm
Maximum height of Q3 on XPR balance with powder dosing head and 125-ml vial	453 mm	426 mm

<sup>1)</sup> If the vial opening is higher than the bottom of the optical sensor, the HeightDetect feature will be unavailable.


### Dosing module

Dimensions in mm.

### Dosing lift

Dimensions in mm.

## 9 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

## 10 Accessories and Spare Parts

### 10.1 Accessories

Accessories are additional components that could help you in your workflow.

Please contact your METTLER TOLEDO sales representative for more details about using accessories with your product.

	Description	Order no.
<b>Dosing heads</b>		
	QH008-BNMW dosing head 10 pcs.	11141533
	QH012-LNMW dosing head 10 pcs.	11141532
	QH012-LNLW dosing head 10 pcs.	11150145
	QH012-LNJW dosing head 10 pcs.	30366317
	QH002-CNMW dosing head 30 pcs. 90 pcs. 500 pcs.	30083440 30098264 30244518
	QH010-CNMW dosing head 10 pcs. 40 pcs. 150 pcs.	30132790 30132791 30046355
	QH008-BNMP dosing head 10 pcs.	11150120

	QH012-LNMP dosing head 10 pcs.	11150118
	QH012-LNCT dosing head 10 pcs.	11150171
	QH012-LNLT dosing head 10 pcs.	11150170
	QH012-LNMX dosing head 10 pcs.	30112276
	QH012-LNLX dosing head 10 pcs.	11150155

Additional dosing heads for specific applications, for example, sterile dosing heads, are available online or through your METTLER TOLEDO sales representative.

	Storage containers with vials 16 ml (for QH008), 10 pcs.	30139824
	Storage containers with vials 125 ml (for QH012), 10 pcs.	30036965
	QA075-P powder test head Filled with $\text{CaCO}_3$ for 15 tests of 10 doses each.	11141506

	QA000-W weight test head	11141507
---	--------------------------	----------

	Dosing head starterkit Set of various dosing heads	30132792
---	---	----------

## Printers

	CLS-631 label printer (RS232C/USB-A) Label and ink ribbon kit The NetCom Kit might be needed for this printer. Contact your local METTLER TOLEDO service representative.	11141820 30004309
---	--	----------------------

	P-52RUE dot matrix printer RS232C, USB and Ethernet connections, simple print-outs Paper roll (length: 20 m), set of 5 pcs Paper roll (length: 13 m), self-adhesive, set of 3 pcs Ribbon cartridge, black, set of 2 pcs	30237290 00072456 11600388 00065975
---	--	--

	P-56RUE thermal printer with RS232C, USB and Ethernet connections, simple print-outs, date and time Paper roll, white (length: 27 m), set of 10 pcs Paper roll, white, self-adhesive (length: 13 m), set of 10 pcs	30094673 30094723 30094724
---	--	----------------------------------

	P-58RUE thermal printer with RS232C, USB and Ethernet connections, simple print-outs, date and time, label printing, balance applications, e.g., statistics, formulation, totaling Paper roll, white (length: 27 m), set of 10 pcs Paper roll, white, self-adhesive (length: 13 m), set of 10 pcs Paper roll, white, self-adhesive labels (550 labels), set of 6 pcs Dimension of the label 56×18 mm	30094674 30094723 30094724 30094725
---	--	--

## ErgoClips

	ErgoClip adapter holder	30521809
---	-------------------------	----------

ErgoClip height extension

30542824

ErgoClip adapter stabilizer

30542831

ErgoClip vial

30521808

### ErgoDiscs

ErgoDisc support plate

30604881

ErgoDisc plate (12 pos)

- capsule size 00
- capsule size 0
- capsule size 1
- capsule size 2
- capsule size 3
- capsule size 4

11141754  
11141751  
11141755  
11141758  
11141752  
30336822

### Antistatic kits

Integrable antistatic kit standard

Including one pair of multiple-point electrodes and power supply.

Must be installed by a METTLER TOLEDO service technician.

30521821

Integrable antistatic kit small

Including one pair of multiple-point electrodes and power supply.

Must be installed by a METTLER TOLEDO service technician.

30521822

**RFID readers / writers / cards**

EasyScan USB  
Reads and writes RFID tags.

30416173

Smart Tag  
Set of 50 pieces  
Set of 200 pieces

30101517

30101518

**Barcode readers**

Corded USB barcode reader

30417466

**Adapters**

**Vial adapters (POM), diameter < 15 cm**

- 8.5 mm × 15 mm (5 pcs.) 30428901
- 9 mm × 25 mm (5 pcs.) 30428902
- 9.5 mm × 25 mm (5 pcs.) 30428903
- 10 mm × 25 mm (5 pcs.) 30428904
- 10.5 mm × 25 mm (5 pcs.) 30428905
- 11 mm × 30 mm (5 pcs.) 30428906
- 11.5 mm × 30 mm (5 pcs.) 30428907
- 12 mm × 20 mm (5 pcs.) 30428908
- 12.5 mm × 30 mm (5 pcs.) 30428909
- 13.5 mm × 20 mm (5 pcs.) 30428910
- 14.5 mm × 25 mm (5 pcs.) 30428911

**Vial adapters (POM), diameter > 15 cm**

- 15.5 mm × 25 mm (5 pcs.) 30428912
- 16.5 mm × 25 mm (5 pcs.) 30428913
- 17.5 mm × 25 mm (5 pcs.) 30428914
- 18.5 mm × 25 mm (5 pcs.) 30428915
- 19.5 mm × 25 mm (5 pcs.) 30428916
- 20.5 mm × 25 mm (5 pcs.) 30428917
- 21.5 mm × 25 mm (5 pcs.) 30428918
- 22.5 mm × 25 mm (5 pcs.) 30428919
- 23.5 mm × 25 mm (5 pcs.) 30428920
- 26 mm × 25 mm (5 pcs.) 30428926
- 27.5 mm × 25 mm (5 pcs.) 30428921
- 28.5 mm × 25 mm (5 pcs.) 30428922
- 29.5 mm × 35 mm (5 pcs.) 30428923
- 31.5 mm × 30 mm (5 pcs.) 30428924

**Bottle adapters (POM)**

33.5 mm × 15 mm (5 pcs.)	30459921
35.5 mm × 15 mm (5 pcs.)	30459922
37.5 mm × 15 mm (5 pcs.)	30459923
40.5 mm × 15 mm (5 pcs.)	30459924
44.5 mm × 15 mm (5 pcs.)	30459925
48 mm × 15 mm (5 pcs.)	30459926
52 mm × 15 mm (5 pcs.)	30459927

**Tube adapters (POM)**

Eppendorf 1.5 ml (5 pcs.)	30306209
Falcon 16.4 mm × 60 mm (5 pcs.)	30459992
Falcon 29.3 mm × 60 mm (5 pcs.)	30459991

**Capsule adapters (stainless steel)**

size 000 (5 pcs.)	30006416
size 00 (5 pcs.)	30006417
size 0 (5 pcs.)	30006418
size 1 (5 pcs.)	30006419
size 2 (5 pcs.)	30006430
size 3 (5 pcs.)	30006431
size 4 (5 pcs.)	30006432

**Dosing components**

QL3 pump	30418660
----------	----------



QS3 sample changer	30418662
Must be installed by a METTLER TOLEDO service technician.	

QLL standard kit for bottles	30008318
<ul style="list-style-type: none"><li>• QL001 liquid dosing head</li><li>• 250 ml bottle, pressure resistant</li><li>• Bottle cap (GL45) with tube connector and dosing head support</li><li>• Spare parts for QLL standard kit</li></ul>	


QLL standard kit for small bottles	30237340
<ul style="list-style-type: none"><li>• QL001 liquid dosing head</li><li>• 25 ml bottle, pressure resistant</li><li>• Stand for bottle</li><li>• Bottle cap (GL25) with tube connector</li><li>• Spare parts for QLL standard kit small</li></ul>	

	<p>QLL advanced kit for bottles</p> <ul style="list-style-type: none"> <li>• QL003 liquid dosing head</li> <li>• 250 ml bottle, pressure resistant</li> <li>• Bottle cap (GL45) with tube connector and dosing head support</li> <li>• Spare parts for QLL standard kit</li> </ul>	30521817
	Holder for 250 ml bottles	30542827
<b>Software</b>		
	<p>LabX Balance Express</p> <p>Stand-alone system, includes one balance license.</p>	11153120
	<p>LabX Balance Server</p> <p>Client server system, includes one balance license.</p>	11153121
	<p>LabX 1 Balance instrument license</p> <p>Single additional instrument license for Express or Server edition.</p>	11153220
<b>Various</b>		
	Cable box	11141845
	EasyHub USB	30468768
	Drip pan, gray	30460856


## 10.2 Spare parts

### 10.2.1 Q3 dosing module


	Order no.	Designation	Remarks
1	30542819	Sealing insert top panel	Including: 5 sealing inserts
2	30525850	Panel top draft shield dosing module	Material: glass
3	30521809	ErgoClip adapter holder	—
4	30542824	ErgoClip height extension	—
5	30542831	ErgoClip adapter stabilizer	—

### 10.2.2 QL3 pump

	Order no.	Designation	Remarks
1	30378287	Muffler for pump	Including: muffler, adapter

### 10.2.3 QLL spare parts kit for bottles

	Order no.	Designation	Remarks
1	30080473	Dosing head liquid QL001	–
2	30348196	Micro dosing valve	Including: micro dosing valve tool
3	30542814	Dosing head liquid QL003	–
4	30080708	Bottle cap with tube connector and support	Thread: GL45
5	30025649	Spare parts for QLL kit	Including: tubing OD 4 mm, ID 2.4 mm (2 m), tubing OD 3.2 mm, ID 1.6 mm (2.5 m), 10 sealing rings for tubing OD 3.2 mm, 5 fastening nuts for tubing OD 3.2 mm, 2 closure pins ISO 2338 4H8x20 A1, 2 suction filters for tubing OD 3.2 mm, tube adapter
6	30460218	Spare parts for QLL kit small	Including: tubing OD 4 mm, ID 2.4 mm (2 m), tubing OD 1.6 mm, ID 0.8 mm (2.5 m), 10 sealing rings for tubing OD 1.6 mm, 5 fastening nuts for tubing OD 1.6 mm, 2 closure pins ISO 2338 4H8x20 A1, tube adapter
7	30542832	Suction filter for tubing OD 3.2 mm	Including: 2 suction filters

## Index

### A

altitude	34
antistatic kit	24

### B

bottle	16
change	26
change content	26
dosing head	25
fill	26

### C

cap, dosing head	13
change	
bottle	26
liquid	26
clean	
housing	29
liquid dosing head	30
sealing insert	29
suction filter	30
weighing chamber	30
compliance information	4
connect	
pump	16
tube	17
contaminated air	21
convention	3
counter dosing cycles	14

### D

dimension	42, 43
disposal	44
dosing cycle	14
dosing head	22, 23
cap	13
clean	30
counter	14
data	14
fill	24
holder	11
liquid	13, 25, 30, 35
material	36
powder	13
powder test head	13
valve	30

dosing lift	11
dosing module	11
draft shield	11

### E

electrostatic charge	24
environmental condition	15, 34
ErgoClip	23
exhaust air	21
external gas	20

### F

fastening nut	31
fill	
bottle	26
dosing head	24

### H

HeightDetect	12
humidity	34

### I

install	
dosing head	22
fastening nut	31
muffler	19
pump	16
sealing insert	23
sealing ring	31
site	15
tube	17
ionizer	24

### L

light	
pressure indicator	27
liquid dosing	16, 25, 27
location	15

### M

manual weighing	27
material	36
micro dosing valve	30
muffler	19

### O

optical sensor	12
options	15

**P**

---

powder dosing	10, 23, 24
power consumption	34
pressure	
indicator light	27
security valve	35
pump	16, 27
contaminated air	21
external gas	20
muffler	19

**R**

---

recommended options	15
remove	
dosing head	22
replace	
fastening nut	31
sealing ring	31
suction filter	30
RFID	13, 14

**S**

---

safety information	6
scope of delivery	15
sealing insert	23, 29
sealing ring	31
spare cap	23
storage	
condition	34
storage container	13, 24
suction filter	30
symbol	3
warning	6

**T**

---

technical data	34
temperature	34
test	
powder test head	13
weight test head	13
top panel	11
toxic substance	21
tube	16, 17, 31

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© Mettler-Toledo GmbH 01/2022