English

Reference Manual **QL3 Pump** for XPR Analytical Balances

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
---

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## Table of Contents

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Further documents and information .....	3
1.2	Explanation of conventions and symbols used .....	3
1.3	Acronyms and Abbreviations .....	3
1.4	Compliance information .....	4
<b>2</b>	<b>Safety Information</b>	<b>4</b>
2.1	Definitions of signal words and warning symbols .....	4
2.2	Product-specific safety notes .....	5
<b>3</b>	<b>Design and Function</b>	<b>7</b>
3.1	Function description .....	7
3.2	Interface of the pump .....	7
<b>4</b>	<b>Installation and Putting into Operation</b>	<b>8</b>
4.1	Scope of delivery .....	8
4.2	Wiring the pump .....	8
4.3	Connecting the tubes .....	8
<b>5</b>	<b>Operation</b>	<b>12</b>
5.1	Operating elements .....	12
5.2	Purging the pump cavity .....	12
<b>6</b>	<b>Maintenance</b>	<b>13</b>
6.1	Cleaning .....	13
6.2	Service .....	13
<b>7</b>	<b>Troubleshooting</b>	<b>13</b>
7.1	Error symptoms .....	13
<b>8</b>	<b>Technical Data</b>	<b>14</b>
8.1	General data .....	14
8.2	Explanatory notes for the METTLER TOLEDO AC/DC adapter .....	15
8.3	Model-specific data .....	15
8.4	Dimensions .....	16
<b>9</b>	<b>Disposal</b>	<b>16</b>
<b>10</b>	<b>Spare Parts</b>	<b>17</b>

## 1 Introduction

Thank you for choosing a METTLER TOLEDO instrument. The instrument combines high performance with ease of use.

### 1.1 Further documents and information

► [www.mt.com/XPR-automatic](http://www.mt.com/XPR-automatic)

This document is available in other languages online.

► [www.mt.com/QL3-RM](http://www.mt.com/QL3-RM)

Search for software downloads

► [www.mt.com/labweighing-software-download](http://www.mt.com/labweighing-software-download)

Search for documents

► [www.mt.com/library](http://www.mt.com/library)

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

► [www.mt.com/contact](http://www.mt.com/contact)

### 1.2 Explanation of conventions and symbols used

#### Conventions and symbols

Key and/or button designations and display texts are shown in graphic or bold text, e.g., , **Edit**.

For useful information about the product.

Refers to an external document.

#### Elements of instructions

In this manual, step-by-step instructions are presented as follows. The action steps are numbered and can contain prerequisites, intermediate results and results, as shown in the example. Sequences with less than two steps are not numbered.

- Prerequisites that must be fulfilled before the individual steps can be executed.

1 Step 1

⇒ Intermediate result

2 Step 2

⇒ Result

### 1.3 Acronyms and Abbreviations

Original term	Explanation
EMC	Electromagnetic Compatibility
FCC	Federal Communications Commission
LPS	Limited Power Source
POM	Polyoxymethylene
RFID	Radio-frequency identification
RM	Reference Manual
sd	Standard deviation
SELV	Safety Extra Low Voltage
SOP	Standard Operating Procedure
UM	User Manual
USB	Universal Serial Bus

## 1.4 Compliance information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

► <http://www.mt.com/ComplianceSearch>

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

► [www.mt.com/contact](http://www.mt.com/contact)

### United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## 2 Safety Information

- Read and understand the instructions in this manual before using the device.
- Keep this manual for future reference.
- Include this manual if you pass on the device to other parties.

If the device is not used according to the instructions in this manual or if it is modified, the safety of the device may be impaired and Mettler-Toledo GmbH assumes no liability.

### 2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

#### Signal words

<b>DANGER</b>	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
<b>WARNING</b>	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
<b>CAUTION</b>	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
<b>NOTICE</b>	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

#### Warning symbols

General hazard

Notice

## 2.2 Product-specific safety notes

### Intended use

This dosing system is designed to be used in analytical laboratories by trained staff. The dosing system is intended for weighing and dosing powder or liquid samples.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

### Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

### Protective equipment

Chemical-resistant gloves

Goggles

Lab coat

### Safety notes

#### **⚠ WARNING**

##### **Death or serious injury due to electric shock**

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

#### **⚠ WARNING**

##### **Injury and/or damage due to hazardous substances**

Chemical, biological or radioactive hazards can be associated with the substances processed by the instrument. During dosing procedures, small amounts of the dosed substance may become airborne and penetrate the instrument or contaminate its surroundings.

The substance characteristics and related hazards is the full responsibility of the instrument owner.

- 1 Be aware of possible hazards associated with the substance and take adequate safety measures, e.g., those stated on the safety data sheet provided by the manufacturer.
- 2 Make sure that every instrument part in contact with the substance will not get altered or damaged by the substance.

### **⚠ WARNING**

#### **Injury and/or damage due to reacting, flammable, or explosive substances**

During the dosing procedure, substances could be combined and cause an exothermic reaction or explosion. This includes powders, liquids, and gases.

The sample characteristics and related hazards is the full responsibility of the instrument owner.

- 1 Be aware of possible hazards associated with reacting, flammable, or explosive substances.
- 2 Ensure a working temperature low enough to prevent the formation of flames or an explosion.

### **⚠ WARNING**

#### **Injury or death due to toxic, explosive, or flammable substances**

If you use toxic, explosive, or flammable liquids with the pump, the exhaust air will be contaminated.

- Connect a tube to the exhaust air outlet to collect the contaminated air.

### **⚠ WARNING**

#### **Injury and/or damage due to reacting substances**

When pressure is released from the bottle, the air/gas in the bottle moves back towards the pump. The air/gas coming from the coupled outlets mixes in the pump. Molecules of the substances in the various bottles can get in contact through this contaminated air/gas.

- 1 Do not connect bottles with incompatible liquids to the same pump simultaneously.
- 2 Before connecting a second, incompatible liquid to the pump, disconnect the first bottle and purge the pump cavity with clean air/gas.

### **⚠ WARNING**

#### **Injury and damage to the pump or the bottle due to high pressure**

High pressure from external gas can damage the pump or the bottle.

- 1 Use a regulator on the external gas line.
- 2 Ensure that the pressure of the external gas does not exceed 0.5 bar (7.2 psi).

### **⚠ CAUTION**

#### **Injury due to splashing liquids**

If the pressure in the bottle is not released, liquid might splash when removing the micro dosing valve, opening the bottle, or removing the liquid tube.

- Always release pressure before removing the micro dosing valve, opening the bottle, or removing the liquid tube.

### **⚠ CAUTION**

#### **Injury due to leaking liquids**

Wrongly cut tubing can result in leaking connections.

- Cut the tubes with a tube cutter or a sharp knife.

### **NOTICE**

#### **Damage to the instrument or malfunction due to the use of unsuitable parts**

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

## NOTICE

### Damage to the instrument

The instrument contains no user-serviceable parts.

- 1 Do not open the instrument.
- 2 In the event of problems, please contact a METTLER TOLEDO representative.

## NOTICE

### Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.


## 3 Design and Function

### 3.1 Function description

The **QL3 pump** is used together with the **Q3 dosing module** or the **QLX3 liquid module** to dose liquids. The setup also involves the **QLL kit for bottle**, which includes the liquid dosing head and all required tubing. The pump builds pressure in the bottle through the air tube. When the pressure is high enough, the micro dispensing valve in the liquid dosing head opens and liquid can ascend the liquid tube. Up to 3 bottles can be connected to the pump simultaneously, and several pumps can be connected to the balance.

The contaminated exhaust air can be collected from an air outlet. Pressure can also be build up by connecting a gas container to the air inlet. Otherwise, the air inlet is connected to a muffler, which is used to absorb the noise.

### 3.2 Interface of the pump

1	Air outlets (to bottle)	4	Socket for AC/DC adapter
2	Pressure release button and indicator light	5	USB-B port (to host)
3	Air exhaust outlet	6	Air inlet

## 4 Installation and Putting into Operation

### 4.1 Scope of delivery

#### QL3 pump

- Pump
- Muffler
- AC/DC adapter with country-specific cable
- USB cable A – B
- Declaration of Conformity
- Reference Manual

#### Recommended options

- EasyHub USB

### 4.2 Wiring the pump


- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
- 2 Insert the plug of the AC/DC adapter in the power inlet of the pump.
- 3 Secure the plug by firmly tightening the knurled nut.
- 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
- 5 Using the USB cable, connect the USB-B port of the pump to one of the USB-A ports of the balance.

### 4.3 Connecting the tubes

#### Tubes definition

The liquid tube is the thinner tube (1) used for transporting liquid from the bottle to the liquid dosing head. The air tube is the slightly bigger tube (2) used for pumping air into the bottle. By adding air through the air tube, pressure rises in the bottle. When the pressure reaches the target pressure, between 0.3 and 0.5 bar (4.4 to 7.2 psi), the micro dispensing valve in the dosing head opens and liquid can ascend the liquid tube.

- 1 Liquid tube
- 2 Air tube

#### Preparing the liquid tube


#### CAUTION

##### Injury due to leaking liquids

- Wrongly cut tubing can result in leaking connections.
- Cut the tubes with a tube cutter or a sharp knife.

- The liquid dosing head is inserted in the liquid dosing head support of the bottle.
- 1 Using a tube cutter or a sharp knife, cut a sufficient amount of tubing. The appropriate length mainly depends on the distance between the balance and the bottle while dosing.  
Recommended length: about 0.9 m

- 2 **NOTICE: Leaking liquid due to wrong assembly. Pay attention to the orientation of the sealing ring when threading it on the tube.**

Place the sealing ring (1) on a flat stable surface, e.g., a table or a workbench, with the wider end downwards.


- 3 Take the end of the liquid tube and press it into the sealing ring.  
 ⇒ This is the dosing-head end of the tube. The opposite end is the bottle end.

- 4 Thread the two fastening nuts (2), paying attention to the orientation.

- 5 **NOTICE: Leaking liquid due to wrong assembly. Pay attention to the orientation of the sealing ring when threading it on the tube.**

Thread the sealing ring (3) from the bottle end of the tube.

- 6 Slide the sealing ring until the distance (D) is enough for the tube to reach the bottom of the bottle. Recommended distances for typical bottle volumes are listed next.


#### Typical distances (D) between the sealing ring and the bottle-end of the tube

Thread	Bottle volume	Distance (D)
GL45	1000 ml	220 mm
	500 ml	170 mm
	250 ml	135 mm
GL25	25 ml	90 mm

#### Connecting the liquid tube

- 1 Insert the dosing-head end of the tube (1) in the dosing head (2).
- 2 Tightly fasten the fastening nut to the dosing head.
- 3 Insert the bottle end of the tube through the corresponding hole in the bottle cap (3). The tube should reach the bottom of the bottle.
- 4 If needed, attach the suction filter (4) to the bottle end of the tube.
- 5 Tightly fasten the fastening nut to the bottle cap.
- 6 Screw the cap to the bottle (5).

The suction filter is used to ensure that no particles or impurities are carried through the liquid dosing head. Using the suction filter will extend the lifetime of the dosing head. However, when dosing solutions, molecules of one of the substances might be absorbed by the suction filter, altering the concentration of the solution. The suction filter should only be used when dosing pure solvents.

#### Connecting the air tube


#### ⚠ WARNING

##### Injury and/or damage due to reacting substances

When pressure is released from the bottle, the air/gas in the bottle moves back towards the pump. The air/gas coming from the coupled outlets mixes in the pump. Molecules of the substances in the various bottles can get in contact through this contaminated air/gas.

- 1 Do not connect bottles with incompatible liquids to the same pump simultaneously.
- 2 Before connecting a second, incompatible liquid to the pump, disconnect the first bottle and purge the pump cavity with clean air/gas.


- 1 Using a tube cutter or a sharp knife, cut a sufficient amount of tubing. The appropriate length mainly depends on the distance between the bottle and the pump while dosing.  
Recommended length: about 0.7 m
- 2 Connect the air tube (1) to the air inlet of the bottle (2)

- 3 Insert the other end of the air tube (1) in a tube adapter (3) and press firmly.
- 4 Connect the tube adapter (3) to one of the air outlets (4) of the pump. Press until you hear a click.


When a tube is connected to the air outlet of the pump, the valve of the air outlet opens. Never leave a tube that is connected to the air outlet unconnected at the other end because pressure cannot build up.

You can connect up to three bottles to each pump.

### Connecting the muffler

- Insert the muffler (1) into the air inlet (2) to absorb the noise.

### Removing the air tube


#### NOTICE

##### Damage to tube connectors due to mishandling


If the tubes are not removed correctly, the connectors of the pump and the bottle cap can be damaged.

- To remove the tubes, press down the ring on the connector and pull out the tube carefully.

- 1 Remove the tube adapter (1) from the pump by pushing the connector (2) of the air outlet towards the pump.  
➔ The adapter is released and can be removed.

- 2 Remove the air tube from the bottle by firmly pressing down the ring (3) and pulling the tube (4) simultaneously.
- 3 If necessary, e.g., for maintenance purposes, use the same procedure to remove the air tube from the tube adapter: firmly press the ring and pull the tube simultaneously.

### Using the pump with external gas

The liquid can be protected by feeding an external gas, e.g., nitrogen, to the pump. Make sure that the pressure of the external gas does not exceed 0.5 bar (7.2 psi).

#### **WARNING**

##### **Injury and damage to the pump or the bottle due to high pressure**

High pressure from external gas can damage the pump or the bottle.

- 1 Use a regulator on the external gas line.
- 2 Ensure that the pressure of the external gas does not exceed 0.5 bar (7.2 psi).


- A regulator is connected to the external gas line.
- 1 Remove the muffler from the air inlet (1).
  - 2 Connect the tube of external gas (2) to the air inlet (1).

#### **Note**

Outer tube diameter: 6 mm

Pressure in external gas line: The pressure must be at least 0.1 bar (1.5 psi). The pressure must not exceed the configured dosing pressure set on the terminal.

The external gas tube is not provided by METTLER TOLEDO.

### Collecting contaminated air

#### **WARNING**

##### **Injury or death due to toxic, explosive, or flammable substances**


If you use toxic, explosive, or flammable liquids with the pump, the exhaust air will be contaminated.

- Connect a tube to the exhaust air outlet to collect the contaminated air.

- Connect a tube (1) to the exhaust air outlet (2) to collect the contaminated air into a safe container (3).

**Note**

Outer tube diameter: 6 mm  
The exhaust air tube and the container are not provided by METTLER TOLEDO.

## 5 Operation


For further information, consult the Reference Manual (RM) of your XPR balance.

► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

### 5.1 Operating elements

#### Releasing the pressure

- The pressure indicator light (1) is on.
- Press the pressure release button (2) to release the pressure.
- ➔ The status light (1) turns off when pressure is released.

#### Pressure indicator light

The pressure indicator light displays the status of the pump:

- Light on: pressure is building up or pressure is established
- Light off: there is no pressure and no pressure is building up
- Light blinking: pump error and/or warning

### 5.2 Purging the pump cavity

If bottles containing incompatible liquids (and for which the fumes should not mix) are connected to the pump one after the other, it is recommended to purge the pump cavity before connecting the second bottle to the pump.

- A dosing head is attached to the dosing or liquid module. The bottle of this dosing head is not connected to any pump.
- A method **Automated dosing** or **Automated solution prep.** is running on the terminal.
  - 1 Disconnect all tube adapters from the pump.
  - 2 Connect an empty tube adapter to the right-most air outlet at the front of the pump.

- ➔ The pump is trying to build pressure and air flows through the cavity of the pump, purging it.
- ➔ The pump cavity has been purged and bottles can safely be re-connected to the air outlets.

## 6 Maintenance

### 6.1 Cleaning

#### NOTICE

##### Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► [www.mt.com/lab-cleaning-guide](http://www.mt.com/lab-cleaning-guide)

The housing material of this accessory is the same as the one of the balance. All surfaces can therefore be cleaned with a commercially available, mild cleaning agent.

Detailed information about the compatibility of cleaning agents can be found in the Reference Manual (RM) of your XPR balance.

► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

### 6.2 Service

There are no service options available for this instrument.

## 7 Troubleshooting

### 7.1 Error symptoms

Error symptom	Possible cause	Diagnostic	Remedy
The status light of the pump does not blink when the device is connected to power.	The pump is disconnected from power.	Disconnect and reconnect the AC/DC adapter from the pump. The indicator light should blink once when the device is connected to power.  Check that the AC/DC adapter and the power cable are not damaged.	Replace the AC/DC adapter and power cable.
	The pump is damaged.	Check with another pump, if available.	Replace the pump.  Contact your METTLER TOLEDO service representative.

<b>Error symptom</b>	<b>Possible cause</b>	<b>Diagnostic</b>	<b>Remedy</b>
On the display, the QL3 pump does not appear in the list of devices connected to the balance.	The pump is disconnected from power.	Disconnect and reconnect the AC/DC adapter from the pump. The indicator light should blink once when the device is connected to power.  Check that the AC/DC adapter and the power cable are not damaged.	Replace the AC/DC adapter and power cable.
	The USB cable is not properly connected.	Check that the USB cable is connected properly.	Connect the USB cable properly.
	The USB cable is damaged.	Check that the USB cable is not damaged.	Replace the USB cable.
	The USB-A port of the balance is damaged.	Disconnect the pump from the USB-A port of the balance. Connect a USB mouse to the same USB-A port. Verify that a pointer (arrow) appears on the terminal and can be moved by moving the mouse.	If the mouse pointer does not appear, contact your METTLER TOLEDO service representative.
	The pump is damaged.	Check with another pump, if available.	Replace the pump.  Contact your METTLER TOLEDO service representative.

## 8 Technical Data

### 8.1 General data

Weight (without packaging): 1600 g

#### Power supply

AC/DC adapter (model no. FSPO60-DHAN3): Input: 100 – 240 V AC  $\pm$  10%, 50 – 60 Hz, 1.8 A  
Output: 12 V DC, 5 A, LPS, SELV

AC/DC adapter (model no. FSPO60-DIBAN2): Input: 100 – 240 V AC  $\pm$  10%, 50 – 60 Hz, 1.5 A  
Output: 12 V DC, 5 A, LPS, SELV

Cable for AC/DC adapter: 3-core, with country-specific plug

Power consumption: 12 V DC  $\pm$  6%, 1 A

Polarity: 

#### Protection and standards

Overvoltage category: II

Degree of pollution: 2

Range of application: Use only indoors in dry locations

#### Environmental conditions

Height above mean sea level: Up to 5000 m

Ambient temperature: +5 – +40 °C

Relative air humidity: 20% to max. 80% at 31 °C, decreasing linearly to 50% at 40 °C, non-condensing

### Storage conditions (in packaging)

Ambient temperature:	-25 – +70 °C
Relative air humidity:	10 – 90%, non-condensing


## 8.2 Explanatory notes for the METTLER TOLEDO AC/DC adapter

The certified external AC/DC adapter complies to the requirements for Class II double insulated equipment. It is not provided with a protective earth connection but with a functional earth connection for EMC purposes. This earth connection **is not** a safety feature. Further information about the compliance of our products can be found in the "Declaration of Conformity" delivered with every product.

In case of testing with regard to the European Directive 2001/95/EC, the AC/DC adapter and the instrument have to be handled as Class II double insulated equipment.

Consequently, a grounding test is not required. It is not necessary to carry out a grounding test between the earth connector of the power plug and any exposed part of the metallic housing of the instrument.

Because the instrument is sensitive to static charges, a leakage resistor of 10 kΩ is connected between the earth connector (1) and the negative pole (2) of the AC/DC adapter. The arrangement is shown in the equivalent circuit diagram. This resistor is not part of the electrical safety arrangement and does not require testing at regular intervals.

## 8.3 Model-specific data

### Liquid dosing heads

For further information about the performance of the liquid dosing heads (QL001, QL003) on your **Q3 dosing module** or **QLX3 liquid module**, consult the corresponding Reference Manual (RM).

► [www.mt.com/Q3-RM](http://www.mt.com/Q3-RM)

► [www.mt.com/QLX3-RM](http://www.mt.com/QLX3-RM)

### Tubing

	Outer diameter	Inner diameter	Recommended length
Liquid tube, for GL45 bottles (dosing head to bottle)	3.2 mm	1.6 mm	0.9 m
Liquid tube, for GL25 bottles (dosing head to bottle)	1.6 mm	0.8 mm	0.9 m
Air tube (pump to bottle)	4.0 mm	2.4 mm	0.7 m
Tube for exhaust air	6 mm	–	–
Tube for external gas	6 mm	–	–

### Suction filter

Tubing outer diameter	3.2 mm
Filter pore size	10 µm

## Bottle

Pressure resistance, min. <sup>1)</sup>	1.5 bar
Burst pressure, min.	3 bar
Volume, max.	2 l

<sup>1)</sup> According to DIN EN 1595: Pressure Equipment made from Borosilicate Glass 3.3 – General Rules for Design, Manufacture and Testing


## QL3 pump

Maximum pressure <sup>1)</sup>	1.5 bar
--------------------------------	---------

<sup>1)</sup> In the unlikely event of a malfunction, the overpressure valve opens if the maximum pressure is reached, to preserve the integrity of the system.

## 8.4 Dimensions

Dimensions in mm.

## 9 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

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▶ [www.mt.com/Q3-RM](http://www.mt.com/Q3-RM)

▶ [www.mt.com/QLX3-RM](http://www.mt.com/QLX3-RM)

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<sup>1)</sup> 準拠：DIN EN 1595: Pressure Equipment made from Borosilicate Glass 3.3 – General Rules for Design, Manufacture and Testing

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

Good Weighing Practice™

---

GWP® is the global weighing standard, ensuring consistent accuracy of weighing processes, applicable to all equipment from any manufacturer. It helps to:

- Choose the appropriate balance or scale
- Calibrate and operate your weighing equipment with security
- Comply with quality and compliance standards in laboratory and manufacturing

 [www.mt.com/GWP](http://www.mt.com/GWP)

[www.mt.com/XPR-automatic](http://www.mt.com/XPR-automatic)

For more information

**Mettler-Toledo GmbH**

Im Langacher 44  
8606 Greifensee, Switzerland  
[www.mt.com/contact](http://www.mt.com/contact)

Subject to technical changes.  
© Mettler-Toledo GmbH 02/2021  
30491781B en, de, es, fr, zh, ja

30491781