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
METTLER TOLEDO

en

This User Manual provides brief instructions about the first steps to take with the instrument. This ensures safe and efficient handling. Personnel must have carefully read and understood this manual before performing any task.

For full information, always refer to the Reference Manual (RM).

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▶ www.mt.com/e10-RM

User Manual **Robotic Mass Comparator**

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Table of Contents

1	Introduction	3
1.1	Further documents and information	3
1.2	Acronyms and abbreviations	3
2	Safety Information	3
2.1	Definition of signal words and warning symbols	4
2.2	Product-specific safety notes	4
2.3	Warning symbols on the robot arm	5
2.4	Switching off in case of an emergency	6
3	Design and Function	6
3.1	Overview instrument	7
3.2	Components description	8
3.2.1	Overview XPR10U balance	8
3.2.2	Overview terminal	9
3.2.3	Overview robot system	10
3.2.4	Overview control unit	10
3.2.5	Control software	10
3.3	Sensors and status indicator	11
4	Installation and Putting into Operation	12
4.1	Selecting the location	12
4.2	Scope of delivery	12
4.3	Installation	13
4.4	Putting into operation	13
4.4.1	Connecting the instrument	13
4.4.2	Switching on the instrument	14
4.4.3	Starting e10control	14
4.4.3.1	Main menus at a glance	15
4.4.3.2	File menu	15
4.4.4	Preparing the weight magazine	15
4.4.4.1	Available weight carriers	16
4.4.4.2	Selecting a suitable weight carrier	17
4.4.4.3	Loading the weight magazine	19
4.4.5	Switching off the instrument	20
4.5	Performing a simple weighing	20
5	Maintenance	21
5.1	Maintenance of the balance	21
5.1.1	Maintenance tasks	21
5.2	Maintenance of the robot system	22
5.3	Cleaning	22
5.3.1	Cleaning the components	22
5.4	Service	23
6	Technical Data	23
6.1	General data	23
7	Disposal	24
8	Compliance Information	24

1 Introduction

Thank you for choosing a METTLER TOLEDO instrument. The instrument combines high performance with ease of use.

Disclaimer for comparators

In this document, the term "balance" is used to describe comparators.

Comparators are characterized by their higher resolution compared to balances. They are mainly used for differential weighing applications, such as the calibration of standard weights. Beside standard balance tests, comparators have also been tested with differential repeatability (ABA repeatability) during production.

EULA

The software in this product is licensed under the METTLER TOLEDO End User License Agreement (EULA) for Software.

When using this product you agree to the terms of the EULA.

▶ www.mt.com/EULA

1.1 Further documents and information

This document is available in other languages online.

Product page:

▶ <http://www.mt.com/lab-robotic-MC>

Instructions for cleaning a balance, "8 Steps to a Clean Balance":

▶ www.mt.com/lab-cleaning-guide

Search for documents:

▶ www.mt.com/library

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

▶ www.mt.com/contact

1.2 Acronyms and abbreviations

Original term	Explanation
AC	Alternating Current
ASTM	American Society for Testing and Materials
DC	Direct Current
EMC	Electromagnetic Compatibility
FACT	Fully automatic time- and temperature-controlled internal adjustment
FCC	Federal Communications Commission
GWP	Good Weighing Practice
ID	Identification
NA	Not Applicable
OIML	Organisation Internationale de Métrologie Légale (International Organization of Legal Metrology)
RM	Reference Manual
SELV	Safety Extra Low Voltage
SOP	Standard Operating Procedure
UM	User Manual
USB	Universal Serial Bus

2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.

- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definition of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

DANGER	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
WARNING	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
CAUTION	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
NOTICE	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Heavy object

Bruising

Notice

2.2 Product-specific safety notes

Intended use

This instrument is designed to be used by trained staff. The automated mass comparator is intended for measuring calibration weights using direct comparison or down-/upward calibration.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

⚠ WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

⚠ WARNING

Injury from lifting heavy objects

The instrument weighs more than what should be lifted by a single person.

- Do not move or lift this equipment without assistance.

NOTICE

Damage to the instrument due to incorrect installation or incorrect repair

- 1 Installation and repairs must be carried out by specially trained METTLER TOLEDO personnel.
- 2 Do not open the balance, the control unit, or the robot system.

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

2.3 Warning symbols on the robot arm


⚠ CAUTION

Injury due to moving parts

The robot arm may move unexpectedly.

- Do not reach into the working area while parts of the instrument are moving.

2.4 Switching off in case of an emergency

- 1 Unplug the power cable to switch off the instrument.
- 2 Contact your METTLER TOLEDO representative.


3 Design and Function

For further information, consult the Reference Manual (RM).

► www.mt.com/e10-RM


3.1 Overview instrument

1	Balance (XPR10U mass comparator)	5	Weight magazine, with weight carriers (60 pcs)
2	Robot hand, with light barrier	6	Terminal for balance
3	Weighing chamber, with sliding door	7	Control unit for robot system
4	Robot system (3 axes)		

3.2 Components description

3.2.1 Overview XPR10U balance


1	Cooling element	8	Door handle
2	USB-B port (to host)	9	Weighing pan
3	USB-A ports (to device)	10	Weighing chamber
4	Weighing unit	11	Draft shield
5	Leveling foot	12	Draft shield cover
6	Weighing chamber plate	13	Model plate
7	Drip tray	14	Weighing display (SmartView)

For further information, consult the Reference Manual (RM).

► www.mt.com/XPR-micro-RM

3.2.2 Overview terminal


	Standby		Tare
	Home screen		Zero
	Open/close door		

 **Note**


No changes must be made on the terminal of the balance. During weighing, the terminal of the balance is disabled.

3.2.3 Overview robot system

1	Robot arm	3	Robot hand
2	Status indicator	4	Light barrier

3.2.4 Overview control unit

1	RS232C serial port (to balance)	4	Socket for interface cable to robot system
2	RS232C serial port (to computer)	5	Socket for AC/DC adapter
3	Ethernet port		

3.2.5 Control software


e10control is a software that serves to operate the instrument. It displays measuring results and other data, and it serves to manage the instrument settings.

The software is installed on a computer provided by METTLER TOLEDO.

3.3 Sensors and status indicator


Light barrier

The light barrier (1) checks the weight carrier before each measurement. It is located on the robot hand. If an unsuitable weight carrier is detected, the measuring process stops.

Balance position sensor

The balance position sensor (1) monitors the position of the balance. The exact position is defined during the installation of the instrument. Do not change the position of the balance, for example, during cleaning. In the event of a position error, contact a METTLER TOLEDO representative.

Status indicator

The robot system is equipped with a status indicator. It indicates whether the robot system is switched off or switched on.

Status indicator	Description
	<p>Status indicator is off</p> <ul style="list-style-type: none"> • The robot system is switched off. • The robot arm does not move automatically.
	<p>Status indicator is on</p> <ul style="list-style-type: none"> • The robot system is switched on. • The robot arm can move automatically. <p>– ⚠ CAUTION: When the status indicator is on, the robot arm can move unexpectedly. Do not reach into the working area while parts of the instrument are moving.</p>

4 Installation and Putting into Operation

4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.

Requirements of the location

Place indoors on stable table


Ensure sufficient spacing

Level the instrument

Provide adequate lighting table


Avoid direct sunlight



Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations

Sufficient spacing: > 30 cm all around the instrument

Note

Place the computer on a separate table to avoid interferences due to vibrations.

Take into account the environmental conditions. See "Technical Data".

4.2 Scope of delivery

Instrument and accessories

- Robot system
- Balance (XPR10U micro comparator)
- Draft shield
- Weight magazine (60 positions)
- Weight carrier, design 1, 30 pcs
- Weight carrier, design 2, 12 pcs
- Weight carrier, design 3, 12 pcs
- Weight carrier, design 4, 6 pcs
- Tweezers, for weights of 1 mg to 50 g
- Tweezers, for weights of 1 g to 200 g
- Air bellow

Instrument control

- Computer
- **e10control** software

Documentation

- User Manual
- Declaration of Conformity
- Production certificate

4.3 Installation

NOTICE

Damage to the instrument and property due to incorrect installation

Incorrect installation and commissioning can lead to damage of the instrument and property.

- Installation and commissioning must be carried out by METTLER TOLEDO specialists or authorized personnel.

See also

 Product-specific safety notes ▶ Page 4

4.4 Putting into operation

After switching on the instrument, it must warm up before giving accurate results.

See also

 General data ▶ Page 23

4.4.1 Connecting the instrument

⚠ WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable and AC/DC adapter designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

NOTICE

Damage to the AC/DC adapter due to overheating

If the AC/DC adapter is covered or in a container, it is not sufficiently cooled and will overheat.

- 1 Do not cover the AC/DC adapter.
- 2 Do not put the AC/DC adapter in a container.

NOTICE


Damage to the balance during startup due to weight carriers or weights

If the robot hand or the balance are equipped with a weight carrier or a weight during the startup procedure, the balance can get damaged.


- Remove weight carriers and weights from the robot hand and the balance before switching on the instrument.

- The instrument is installed by the manufacturer.
 - The instrument is not yet connected to the power supply.
 - The robot hand and the weighing pan of the balance are free of weight carriers or weights.
- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.

- 2 Insert the plug of the AC/DC adapter in the power inlet of the balance.
- 3 Secure the plug by firmly tightening the knurled nut.
- 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
 - ➔ The balance is switched on.

- 5 Check that the robot system and the control unit are connected (1).
- 6 Insert the plug of the AC/DC adapter in the power inlet of the control unit (2).
- 7 Insert the plug of the power cable into a grounded power outlet that is easily accessible.
 - ➔ The robot system is switched on.

Note

Do not connect the instrument to a power outlet controlled by a switch. After switching on the instrument, it must warm up before giving accurate results.

To ensure optimum weighing conditions, leave the robot system on all the time.

Note

If any of the following situations is detected, each robot axis moves to its home position:

- balance is turned off
- interface connection is inactive
- balance was moved from its original position

4.4.2 Switching on the instrument

When the instrument is connected to the power supply, it automatically switches on.

Acclimatization and warm up

Before the balance gives reliable results, it must:

- acclimatize to the room temperature
- warm up by being connected to the power supply

The acclimatization time and warm-up time for the balance are available in "General data".

4.4.3 Starting e10control

- Double-click the **e10control** icon (**e10control.exe**).
 - ➔ The software opens.
 - ➔ A blank settings file named **Untitled.e10** is displayed.

4.4.3.1 Main menus at a glance

Menu	Description
File	Contains file-related topics, for example creating a new file.
Edit	Serves to edit a file.
View	Serves to adapt the software view.
Weights	Gives access to the weights database, including all relevant data on your standards and test weights.
Magazine	Serves to identify and register the weights placed on the weight magazine.
Process	Serves to set up the weighing process.
Report	Serves to define the content of the report file.
Adjustment	Serves to start the adjustment procedure, using the internal adjustment weights of the balance.
System	Serves to adjust system settings.
Start	Serves to start a weighing process.
Help	Contains the help file and further information about the software.

4.4.3.2 File menu

Command	Description
New	Closes the current settings file and opens a new, blank file named Untitled.e10 .
Open...	Serves to open an existing file. Once a particular file is selected, the currently open settings file closes.
Import...	Serves to import an existing text file (extension .imp or .txt). During the import, the text file is converted into a new standard settings file (extension: .e10). The new settings file is immediately loaded.
Save	Saves the changes made to the current settings file under the current file name (extension: .e10).
Save As...	Serves to save a settings file under a new name (extension: .e10).
Save as Text...	Serves to save the settings as a text file (current name with extension .txt). The following options are available: <ul style="list-style-type: none"> • Standards data for selected sets • Test weights data • Magazine places allocation • Weighing process settings • Series scheme • Report heading
Exit	Saves any changes and quits e10control .

For further information, consult the Reference Manual (RM).

► www.mt.com/e10-RM

4.4.4 Preparing the weight magazine

Each test weight or standard used during the weighing process needs to be placed onto one weight carrier. The selection of an adequate weight carrier is determined by the weight geometry.


In order to ensure a trouble-free operation of the comparator and to minimize corner load errors, strict rules must be followed when choosing the carrier type.

Note

Do not touch the weight carriers or the weights with bare hands. Use the provided tweezers or powder-free gloves.


4.4.4.1 Available weight carriers

Design 1


1	Hook for wire weights	3	Handle (grey) for design 1
2	Indent for cylindrical and sheet weights		

Design 2


1	Indent for cylindrical and sheet weights	2	Handle (grey) for design 2
----------	--	----------	----------------------------

Design 3

1	Indent for cylindrical and sheet weights	2	Handle (red) for design 3 and design 4
----------	--	----------	--

Design 4

1	Indent for cylindrical and sheet weights	2	Handle (red) for design 3 and design 4
----------	--	----------	--

4.4.4.2 Selecting a suitable weight carrier

Note

Only use weights that fit into one of the listed categories.


Note

Make sure to place the weights correctly on the corresponding carrier. Weights must never be placed over the edge of the indent.

In the mode **down-/upward calibration**, combinations of up to three weights can be weighed:


- Design 1 or 2: up to three weights, each placed on its own carrier
- Design 1 or 2, in combination with design 3 or 4: up to two weights, each placed on its own carrier

Cylindrical weight with knob


Weight shape	Weight dimension	Suitable weight carrier
	diameter: $4 \text{ mm} \leq d_c \leq 8.2 \text{ mm}$ height: $h_c \leq 16 \text{ mm}$	Design 1 Design 2
	diameter: $8.2 \text{ mm} \leq d_c \leq 14 \text{ mm}$ height: $h_c \leq 19 \text{ mm}$	Design 3
	diameter: $8.2 \text{ mm} \leq d_c \leq 10 \text{ mm}$ height: $h_c \leq 19 \text{ mm}$	Design 4


Wire weight




Weight shape	Weight dimension	Suitable weight carrier
triangular 	triangle inner side: $5.5 \text{ mm} \leq s_w \leq 18 \text{ mm}$ wire diameter: $0.05 \text{ mm} \leq d_w \leq 1.5 \text{ mm}$ height: $h_w \leq 6 \text{ mm}$	Design 1
square 	square inner side: $5.5 \text{ mm} \leq s_w \leq 12 \text{ mm}$ wire diameter: $0.05 \text{ mm} \leq d_w \leq 1.5 \text{ mm}$ height: $h_w \leq 6 \text{ mm}$	Design 1

Weight shape	Weight dimension	Suitable weight carrier
pentagonal 	pentagon inner side: $5.5 \text{ mm} \leq s_w \leq 12 \text{ mm}$ wire diameter: $0.05 \text{ mm} \leq d_w \leq 1.5 \text{ mm}$ height: $h_w \leq 6 \text{ mm}$	Design 1

Sheet weight, polygonal

Weight shape	Weight dimension	Suitable weight carrier
	diameter circumscribed circle: $d_s \leq 4 \text{ mm}$ distance l_s: $l_s \geq 3 \text{ mm}$ angle sheet-handle: $\alpha \geq 90^\circ$	Design 1 Design 2
	diameter circumscribed circle: $4 \text{ mm} \leq d_s \leq 8.2 \text{ mm}$	Design 1 Design 2
	diameter circumscribed circle: $8.2 \text{ mm} \leq d_s \leq 14 \text{ mm}$	Design 3

See also


 Available weight carriers ▶ Page 16

4.4.4.3 Loading the weight magazine


To place the weights onto their respective weight carrier in the magazine, proceed as follows:

- The instrument is switched on.
- 1 Open the doors of the weighing chamber.

- 2 Insert the required weight carriers into the weight magazine.
Note
 The carriers must be placed with their handle on the side opposite to the position numbers.
Note
 User tweezers or wear gloves to handle the weights and the weight carriers. Do not touch the weights or the weight carriers with bare hands.
- 3 Center the weight carriers in the provided space on the magazine.

- 4 Correctly place the weights on the appropriate weight carriers.
- 5 Once the magazine is loaded, close the weighing chamber doors.

See also

Switching on the instrument ▶ Page 14

4.4.5 Switching off the instrument

To completely switch off the instrument, it must be disconnected from the power supply.

4.5 Performing a simple weighing

Connecting software with system

- 1 Select **System > Comparator serial port... > Serial port** to connect the software with the system.
- 2 Enter the port to which the interface cable is connected.
- 3 Fill in the required information for **Robot system controller type** and for **Balance type**.

Starting the weighing process

- 1 Select **Start > Start measurement** to start the weighing process.
 - ➔ A window to enter a name for the report file opens.
- 2 Enter the name and the target location of the report file.
- 3 Click **Save**.
 - ➔ The window **Weighing process settings** opens.
- 4 Click **OK**.
 - ➔ The weighing process starts.
 - ➔ The weighing process monitor opens.

Adjust settings

- 1 Start **e10control**.
 - ➔ A new blank setting file opens.
- 2 Select **Weights** to enter and edit standards data and test weights data. The default password to change the currently stored data is "" (empty character string = no character).

- 3 Select **Magazine > Places allocation... > Allocation of weight magazine places** to identify and register the position of the test weights on the weight magazine.
- 4 Select **Process > Settings... > Weighing process settings** to set the comparisons for this weighing process, the precise timing, and the sequence.
- 5 Set the following parameters:
 - **Weighing mode**
 - **Pre-run requested**
 - **History-specific pause requested**
 - **Start delay**
 - **No. of nonreported preweighings per group (0-5)**
 - **No. of reported comparisons per group (1-20)**
 - **No. of series (1-20)**
 - **Stabilisation time (10-60 s)**
 - **Integration time (0-60 s)**
 - **Comparison scheme**
 - **Sensitivity check**
 - **Sensitivity check standard**
- 6 Select **Process > Settings... > Weighing process settings > Series scheme** to determine the series scheme.
- 7 Select **Report > Contents...** to define the contents of the report file.

For further information, consult the Reference Manual (RM).

► www.mt.com/e10-RM

5 Maintenance

To guarantee the functionality of the balance and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

For further information, consult the Reference Manual (RM).

► www.mt.com/e10-RM

5.1 Maintenance of the balance

5.1.1 Maintenance tasks

Maintenance action	Recommended interval	Remarks
Performing an internal adjustment	<ul style="list-style-type: none"> • Daily • After cleaning • After leveling • After changing the location 	see "Performing an internal adjustment" in the Reference Manual for the balance
Performing routine tests (eccentricity test, repeatability test, sensitivity test). METTLER TOLEDO recommends to at least perform a sensitivity test.	<ul style="list-style-type: none"> • After cleaning • After assembling the balance • After a software update • Depending on your internal regulations (SOP) 	see "Tests" in the Reference Manual for the balance
Cleaning	<ul style="list-style-type: none"> • After every use • Depending on the degree of pollution • Depending on your internal regulations (SOP) 	see "Cleaning"

For further information, consult the Reference Manual (RM).

► www.mt.com/XPR-micro-RM

5.2 Maintenance of the robot system

Apart from cleaning no regular maintenance is required by the instrument owner.

5.3 Cleaning

5.3.1 Cleaning the components

Overview

Periodically, clean the following parts of the instrument:

System part	Task	Tool	Notes
Weighing chamber	Remove dust.	Air bellows	
Balance: <ul style="list-style-type: none"> • Weighing pan • Drip tray • Housing • Terminal 	Remove dust.	Air bellows	The position of the balance must remain absolutely unchanged.
Weight magazine: <ul style="list-style-type: none"> • Weights • Weight carriers 	Remove dust.	Air bellows	The centering holes underneath the carriers should also be clean. Store weight carriers that are not in use in a dust-free environment. Do not leave these weight carriers in the weight magazine.
Robot system: <ul style="list-style-type: none"> • Robot hand • 3 carrier-centering cones 	Remove dust.	Air bellows	Do not use compressed air or petroleum-based solvents.
<ul style="list-style-type: none"> • Light barrier 	If a problem with the light barrier on the robot hand occurs, clean the light barrier.	Air bellows	Do not use any solvent or ethanol!

NOTICE

Damage to the instrument due to inappropriate cleaning methods

If liquid enters the housing, it can damage the instrument. The surface of the instrument can be damaged by certain cleaning agents, solvents, or abrasives.

- 1 Do not spray or pour liquid on the instrument.
- 2 Only use the cleaning agents specified in the Reference Manual (RM) of the instrument or the guide "8 Steps to a Clean Balance".
- 3 Only use a lightly moistened, lint-free cloth or a tissue to clean the instrument.
- 4 Wipe off any spills immediately.

For further information on cleaning a balance, consult "8 Steps to a Clean Balance".

► www.mt.com/lab-cleaning-guide

Cleaning around the balance

- Remove any dirt or dust around the balance and avoid further contaminations.

Cleaning the terminal

- Clean the terminal with a damp cloth or a tissue and a mild cleaning agent.

Cleaning the removable parts

- Clean the removed part with a damp cloth or a tissue and a mild cleaning agent or clean in a dishwasher up to 80 °C.

Cleaning the weighing unit

- 1 Disconnect the balance from the AC/DC adapter.
- 2 Use a lint-free cloth moistened with a mild cleaning agent to clean the surface of the balance.
- 3 Remove powder or dust with a disposable tissue first.
- 4 Remove sticky substances with a damp lint-free cloth and a mild solvent, e.g., isopropanol or ethanol 70%.

5.4 Service

Regular servicing by an authorized service technician ensures reliability for years to come. Contact your METTLER TOLEDO representative for details about the available service options.

6 Technical Data

6.1 General data

Automated weight handler

Weight handler

For automatic determination of test weights. This can be done by comparison of one test weight against one standard weight. Comparisons of weight combinations, consisting of up to three weights, are also possible.

Measuring time (typical)

15 min for a series of five comparisons of one weight against another weight (**One-vs.-one comparisons**). 30 min for a series of five comparisons of one weight against a weight combination. Typically, a 5 × **A-B-A** comparative weighing applies.

Test weights / standards

Knob weights, wire weights, and sheet weights in various shapes and sizes. Nominal value: 1 mg up to 5 g

Weight magazine

60 places

Control software

e10control

Data interface

RS232C to controller

Balance - XPR10U Comparator

Readability:

0.1 µg

Maximum capacity:

10.1 g

Electrical weighing range:

10.1 g

Repeatability:

Determined as standard deviation of 5 × **A-B-A** comparative weighing:

- 0-1 g: $s \leq 0.15 \mu\text{g}$
- 1-2 g: $s \leq 0.25 \mu\text{g}$
- 2-6 g: $s \leq 0.40 \mu\text{g}$
- 6-10 g: $s \leq 0.60 \mu\text{g}$

Linearity:

± 2 µg

Setting time (typical):

20 s

Adjustment:

Automated internal adjustment using built-in weights of the balance, or external adjustment using external weights.

Power supply

Balance AC/DC adapter:	Primary: 100 – 240 V AC, -15%/+10%, 50/60 Hz Secondary: 12 V DC $\pm 3\%$, 2.5 A (with electronic overload protection)
Balance power consumption:	12 V DC $\pm 3\%$, 2.25 A, maximum ripple: 80 mVpp
Polarity:	⊖ ⊕ with a current limited SELV output
Robot system AC/DC adapter:	Primary: 100 – 240 V AC, $\pm 10\%$, 50/60 Hz Secondary: 24 V DC, $\pm 5\%$, 2.1 A (with electronic overload protection)
Robot power consumption:	24 V DC $\pm 5\%$, 1.5 A
Cable for AC/DC adapter:	3-core, with country-specific plug
Power consumption robot system:	36 VA max. (24 V DC 1.5 A)

Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Protection:	Protected against dust and water
Standards for safety and EMC:	See Declaration of Conformity
Range of application:	For use in closed interior rooms only

Environmental conditions

Height above mean sea level:	Up to 2000 m
Ambient temperature:	17 – 27 °C (± 0.5 °C / 12 hour)
Relative air humidity:	45 – 60 %, non-condensing
Vibrations:	Set up in a room free of vibrations
Acclimatization time:	At least 24 hours after placing the instrument in the same location where it will be put into operation.
Warm-up time:	At least 4 hours after connecting the balance to the power supply. When switched on from standby, the instrument is ready for operation immediately.

Instrument (robot system and balance)

Weight:	50 kg
---------	-------

See also

🔗 Selecting a suitable weight carrier ▶ Page 17

7 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

8 Compliance Information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

▶ www.mt.com/ComplianceSearch

For further information, consult the Reference Manual (RM).

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▶ www.mt.com/contact

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原文	译文	说明
AC		Alternating Current (交流电)
ASTM		American Society for Testing and Materials (美国试验与材料协会)
DC		Direct Current (直流电)
EMC		Electromagnetic Compatibility (电磁兼容)
FACT		Fully automatic time- and temperature-controlled internal adjustment (全自动时间和温度控制的内部校正)
FCC		Federal Communications Commission (美国联邦通讯委员会)
GWP		Good Weighing Practice
ID		Identification (标识)
NA		Not Applicable

	(不适用)
OIML	Organisation Internationale de Métrologie Légale (国际法制计量组织)
RM	Reference Manual (参考手册)
SELV	Safety Extra Low Voltage (额定安全低电压)
SOP	Standard Operating Procedure (标准操作程序)
UM	User Manual (简明用户手册)
USB	Universal Serial Bus

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GWP®

Good Weighing Practice™

GWP® is the global weighing standard, ensuring consistent accuracy of weighing processes, applicable to all equipment from any manufacturer. It helps to:

- Choose the appropriate balance or scale
- Calibrate and operate your weighing equipment with security
- Comply with quality and compliance standards in laboratory and manufacturing

 www.mt.com/GWP

www.mt.com/comparator

For more information

Mettler-Toledo GmbH

Im Langacher 44
8606 Greifensee, Switzerland
www.mt.com/contact

Subject to technical changes.

© 08/2022 METTLER TOLEDO. All rights reserved.
30572159A en, zh

30572159