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METTLER TOLEDO

This User Manual provides brief information about XK-KC comparator models. Personnel must have carefully read and understood this manual before performing any tasks.

For full information, refer to the manuals provided with the terminal and the weighing platform.

► www.mt.com/IND900

► www.mt.com/PBK9

► www.mt.com/PFK9

All XK-KC comparator models are delivered with an IND930 terminal. Refer to the table below to see which weighing platform is included (PBK or PFK).

XK-KC comparator model	Weighing platform
XK155KSD5C	PBK987-CC150
XK604KMC	PFK988-C600
XK1003KMD2C	PFK988-C600
XK1004KMD5C	PFK988-C600
XK2003KLC	PFK988-E1500
XK6002KLC	PFK988-E1500

The availability of the models is country-specific. For information about the availability of the models in your country, please contact your authorized METTLER TOLEDO dealer or service representative.

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
Overview terminal

en

1	Menu and messages	4	Selected weighing platform
2	Weighing information	5	Keys
3	Access level, region, time, date		

Key	Name	Description
	Clear	Clears the tare value.
	Tare	Tares the balance. This function is used when the weighing process involves containers.
	Zero	Zeroes the balance. The balance must always be zeroed before starting the weighing process.
	Print / Transfer	Sends or transfers data, e.g., to a connected printer.
	Preset Tare	Opens a dialog to manually set a tare value.
	Switch balance	Switches between connected weighing platforms (if available).
	Switch units	Switches the weighing unit.
	Tare table	Lists tare values of defined containers.

This is only a selection of all available keys. For further information, consult the manual provided with the terminal.

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Table of Contents

1	Introduction	3
1.1	Further documents and information	3
1.2	Acronyms and abbreviations	3
1.3	Compliance information	3
2	Safety Information	4
2.1	Definitions of signal words and warning symbols	4
2.2	Product-specific safety information	4
3	Design and Function	6
3.1	Overview	6
3.2	User interface	7
4	Installation and Putting into Operation	7
4.1	Selecting the location	7
4.2	Unpacking	8
4.3	Scope of delivery	8
4.4	Installation	8
4.4.1	Installing the terminal	8
4.5	Putting into operation	9
4.5.1	Connecting the balance	9
4.5.2	Switching on the balance	9
4.5.3	Leveling the balance	9
4.5.4	Adjusting the balance	10
4.5.5	Switching off the balance	10
4.6	Performing a simple weighing	10
4.6.1	Zeroing the balance	10
4.6.2	Taring the balance	11
4.6.3	Performing a weighing	11
4.7	Transporting, packing, and storing	11
4.7.1	Transporting over short distances	11
4.7.2	Transporting over long distances	11
4.7.3	Packing and storing	12
5	Maintenance	12
5.1	Maintenance tasks	12
5.2	Cleaning	12
6	Technical Data	13
6.1	General data	13
6.2	Model-specific data	14
6.2.1	XK-KC weighing platforms	14
6.2.2	Terminal	16
7	Accessories	17
8	Disposal	18

1 Introduction

Thank you for choosing a METTLER TOLEDO balance. The balance combines high performance with ease of use.

Disclaimer for comparators

In this document, the term "balance" is used to describe comparators.

Comparators are characterized by their higher resolution compared to balances and are mainly used for differential weighing application, such as the calibration of standard weights. Beside standard balance tests, comparators have also been tested with differential repeatability (ABA repeatability) during production.

1.1 Further documents and information

This document is available online.

► www.mt.com/XK-KC-UM

► www.mt.com/comparators

Search for documents

► www.mt.com/IND900

► www.mt.com/PBK9

► www.mt.com/PFK9

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

► www.mt.com/contact

1.2 Acronyms and abbreviations

Original term	Explanation
EMC	Electromagnetic Compatibility
FACT	Fully automatic time- and temperature-controlled internal adjustment
FCC	Federal Communications Commission
GWP	Good Weighing Practice
UM	User Manual
USB	Universal Serial Bus
USP	United States Pharmacopeia

1.3 Compliance information

National approval documents, e.g., the FCC Supplier Declaration of Conformity, are available online and/or included in the packaging.

► <http://www.mt.com/ComplianceSearch>

Contact METTLER TOLEDO for questions about the country-specific compliance of your instrument.

► www.mt.com/contact

United States of America

This equipment has been tested and found to comply with the limits for a **Class A** digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.

- This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

2 Safety Information

The following documents are available for this device.

- The User Manual at hand is printed and delivered with XK-KC comparator models.
- The other manuals delivered with the device contain a full description of the terminal (IND930) and the weighing platform (PBK or PFK models), see table at the very beginning of this manual.
- Keep all documents for future reference.
- Include all documents if you transfer the device to other parties.

Only use the device according to the provided manuals. If you do not use the device according to these documents or if the device is modified, the safety of the device may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

DANGER A hazardous situation with high risk, resulting in death or severe injury if not avoided.

WARNING A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.

CAUTION A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.

NOTICE A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

General hazard

Heavy object

Inflammable or explosive substance

Bruising

Notice

2.2 Product-specific safety information

Intended use

This device is designed to be used by trained staff. The XK-KC comparators are intended for weighing.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

Safety notes

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

WARNING

Injury due to lifting heavy objects

Attempts to lift the instrument by a single person could cause injuries.

- Do not move or lift this equipment without assistance.

WARNING

Death or serious injury due to toxic substances and vapors

Handling hazardous substances, and possible resulting vapors, can cause injuries if they come in contact with the skin or the eyes or if they are inhaled.

- 1 When using chemicals and solvents, comply with the instructions of the manufacturer and general safety rules.
- 2 Set up the instrument in a well-ventilated location.

WARNING

Death or serious injury due to flammable substances

Flammable substances can ignite and lead to fire and explosions.

- 1 Keep flammable substances away from naked flames.
- 2 When using chemicals and solvents, comply with the instructions of the manufacturer and general safety rules.

CAUTION

Injury due to moving parts

Adjusting the position of the terminal screen on its stand can bruise your hands.

- Make sure not to squeeze your hands between the screen and its stand.

NOTICE

IP23 becomes void if weighing pan not installed

According to the ingress protection grade IP23, the balance is protected, to a certain extent, from water sprays. The protection IP23 becomes void if the weighing pan is not installed on the weighing platform.

- Do not spray liquid on the balance if the weighing pan is not installed on the weighing platform.

NOTICE

Damage to the touch screen due to pointed or sharp objects

- Operate the touch screen with your fingers.

NOTICE

Damage to the balance

The balance contains no user-serviceable parts.

- 1 Do not open the balance.
- 2 In the event of problems, please contact a METTLER TOLEDO representative.

NOTICE

Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

Your balance represents state-of-the-art technology and complies with all recognized safety rules, however, certain hazards may arise in extraneous circumstances. Do not open the housing of the balance; it does not contain any parts that can be maintained, repaired or replaced by the user. If you experience problems with your balance, contact your authorized METTLER TOLEDO dealer or service representative.

The balance has been tested for the experiments and intended purposes documented in the appropriate manual. However, this does not absolve you from the responsibility of performing your own tests of the products supplied by us regarding their suitability for the methods and purposes you intend to use them for.

3 Design and Function

For further information, consult the manuals provided with the terminal and the weighing platform.

▶ www.mt.com/IND900

▶ www.mt.com/PBK9

▶ www.mt.com/PFK9

3.1 Overview

The XK-KC comparators are modular weighing systems consisting of a terminal and, by default, one weighing platform. Alternatively, up to four weighing platforms can be connected to one terminal. This requires additional accessories. For further information, please contact your authorized METTLER TOLEDO dealer or sales representative.


See the sections "Overview" (graphics and legend) at the very beginning of this manual.

3.2 User interface

The terminal serves to operate the weighing platform.

Main weighing screen

The main weighing screen is the central navigation point where all the menus and settings can be found. The functions open when tapping the corresponding icon on the touch screen.

For further information, consult the manual provided with the terminal.

► www.mt.com/IND900

4 Installation and Putting into Operation

4.1 Selecting the location

A balance is a sensitive precision instrument. The location where it is placed will have a profound effect on the accuracy of the weighing results.



Requirements of the location

Place indoors on stable floor


Ensure sufficient spacing

Level the instrument

Provide adequate lighting

Avoid direct sunlight


Avoid vibrations

Avoid strong drafts

Avoid temperature fluctuations

Take into account the environmental conditions. See "Technical Data".

4.2 Unpacking

Check the packaging and the delivered components for transportation damage or missing parts. Please inform a METTLER TOLEDO service representative in the event of missing or defective parts.

4.3 Scope of delivery

Components	XK-KC comparators
Weighing platform	✓
Terminal	✓
Terminal stand	✓
Connection cable	✓
Power cable (country-specific)	✓
MC Link software	✓
User Manual	✓
Production certificate	✓
Declaration of Conformity	✓

4.4 Installation

The XK-KC comparator must be installed by a METTLER TOLEDO service technician.

4.4.1 Installing the terminal

NOTICE

Danger of damage to the device

The connection cable must be mechanically secured against unplugging at the terminal.

The terminal is delivered with a stand. Alternatively, the terminal can be attached to a wall using a wall fixture. For information regarding the wall fixture, please contact your authorized METTLER TOLEDO dealer or service representative.

4.5 Putting into operation

4.5.1 Connecting the balance

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the METTLER TOLEDO power cable designed for your instrument.
- 2 Connect the power cable to a grounded power outlet.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and the power plug for damage and replace them if damaged.

- 1 Install the cables in such a way that they cannot be damaged or interfere with operation.
- 2 Plug the connection cable into the weighing platform and the terminal (interface SICS422).
- 3 Secure both plugs by firmly tightening the knurled nut.
- 4 Insert the plug of the power cable into a grounded power outlet that is easily accessible.

Note

Do not connect the instrument to a power outlet controlled by a switch. After switching on the instrument, it must warm up before giving accurate results.

4.5.2 Switching on the balance

- The weighing platform is connected to the terminal.
- Insert the plug of the power cable into a power outlet.
- ➔ The terminal starts up automatically.
- ➔ The balance is switched on.


4.5.3 Leveling the balance

Precise leveling is essential for repeatable and accurate weighing results.

The balance must be leveled and adjusted each time it is moved to a new location.


XK155KSD5C, XK604KMC, XK1003KMD2C, XK1004KMD5C

- 1 Level the weighing platform by turning the 4 leveling screws (1).
 - 2 Check if all leveling screws sit flat on the surface.
 - 3 Check the stability of the weighing platform by pressing down on each corner.
- ➔ If leveling is successful, the air bubble is in the center of the level indicator (2).

XK2003KLC, XK6002KLC

- 1 Level the weighing platform by turning the 4 leveling screws (1).
 - 2 The leveling screws can be adjusted with a size 30 engineer's wrench.
 - 3 Ensure that the leveling screws all sit flat on the floor evenly.
 - 4 Test the stability of the weighing platform by pressing down on or rocking the corners.
- ➔ The level is correct when the air bubble lies in the center of the level indicator (2).

4.5.4 Adjusting the balance

- The balance is switched on.
- 1 Tap to open the system menu.
 - 2 Tap Terminal, then Metrology.
 - 3 Tap to start the adjustment (FACT).
 - ➔ The adjustment screen opens.
 - 4 Tap to start the adjustment.
 - ➔ The result of the adjustment for all connected weighing platforms is displayed.

4.5.5 Switching off the balance

- The balance is switched on.
- Unplug the power cable from the power outlet.
 - ➔ The balance is switched off.

4.6 Performing a simple weighing

For further information, consult the manuals provided with the terminal and the weighing platform.

► www.mt.com/IND900

► www.mt.com/PBK9

► www.mt.com/PFK9

4.6.1 Zeroing the balance

- The balance is switched on.
- 1 Clear the weighing pan.
 - 2 Press to zero the balance.
- ➔ The balance is zeroed.

4.6.2 Taring the balance

If a weighing container is used, the balance must be tared.

- 1 Clear the weighing pan.
 - 2 Press [**→0←**] to zero the balance.
 - 3 Place the container on the weighing pan.
 - 4 Press [**→T←**] to tare the balance.
- ⇒ The balance is tared. The icon **Net** appears.

4.6.3 Performing a weighing

- The weighing pan is clear and the display shows 0.
- 1 Load a weight onto the weighing pan.
 - ⇒ The display shows the corresponding weight.
 - 2 Remove the weight from the weighing pan.
 - ⇒ The display shows 0.

4.7 Transporting, packing, and storing

⚠ WARNING

Injury due to lifting heavy objects

Attempts to lift the instrument by a single person could cause injuries.

- Do not move or lift this equipment without assistance.

Preparing the balance for transport

- The balance is switched on.
- 1 Unplug the power cable from the power outlet.
 - ⇒ The balance is switched off.
 - 2 Disconnect the connection cable from the terminal and the weighing platform.
 - ⇒ The balance is ready for transport.

4.7.1 Transporting over short distances

To move the balance over a short distance to a new location, follow the instructions below.

- The balance is prepared for transport.
- 1 Remove the weighing pan to transport it separately to the new location.
 - 2 Carefully lift the weighing platform and carry it in horizontal position to the new location.
 - 3 Carefully lift the terminal and carry it to the new location.

To put the balance into operation, proceed as follows:

- 1 Connect the cables in reverse order.
- 2 Level the balance.
- 3 Install the weighing pan.
- 4 Perform an adjustment.

See also

- 🔗 Leveling the balance ▶ Page 9
- 🔗 Adjusting the balance ▶ Page 10

4.7.2 Transporting over long distances

To transport the balance over long distances, always use the original packaging.

4.7.3 Packing and storing

Packing

Store all parts of the packaging in a safe place. The elements of the original packaging are developed specifically for the balance and its components to ensure maximum protection during transportation or storing.

Storing

Store the balance under following conditions:

- Indoor and in the original packaging.
- According to the environmental condition, see "Technical data".

5 Maintenance

To guarantee the functionality of the balance and the accuracy of the weighing results, a number of maintenance actions must be performed by the user.

Contact a METTLER TOLEDO representative to find about the service options available – regular maintenance by an authorized service engineer will ensure consistent weighing accuracy over the long term and extend the service life of the balance.

5.1 Maintenance tasks

The terminal offers various maintenance options.

For further information, consult the manual provided with the terminal.

► www.mt.com/IND900

Maintenance of the weighing platform is limited to regular cleaning.

5.2 Cleaning

NOTICE

Damage to balance due to inappropriate cleaning methods

The balance is made from high quality, resistant materials and can be damaged by certain cleaning agents, solvents or abrasives. Any liquid that enters the housing may damage the balance.

- 1 Use water and a mild detergent to clean the balance or terminal.
- 2 Wipe off any spills immediately.
- 3 Prevent liquid from entering the interior of the balance.

Periodically clean the weighing platform and the terminal using a damp cloth. Use household cleaning agents. If required, the weighing pan can be removed to clean the interior of the weighing platform.

For further information, consult the manual provided with the weighing platform.

► www.mt.com/PBK9

► www.mt.com/PFK9

6 Technical Data

6.1 General data

Power supply

Power supply:	100 – 240 V AC +10% / -15%, 50 – 60 Hz
Power cable:	3-core, with country-specific plug
Terminal current consumption:	650 mA - 275 mA

Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Protection (EN40050):	Weighing platform: IP23D Terminal: IP69k
Standards for safety and EMC:	See Declaration of Conformity
Range of application:	Use only indoors in dry locations

Environmental conditions

Height above mean sea level:	Up to 2000 m
Ambient temperature:	+10 °C – +30 °C ± 1 °C / 12 h
Storage condition:	-20 °C – +60 °C
Relative air humidity:	40 – 70 % ± 10 % / 4 h, non-condensing

Materials

Weighing platform:	XK155KSD5C, XK604KMC, XK1003KMD2C, XK1004KMD5C: powder coated XK2003KLC, XK6002KLC: hot galvanized S235JR
Terminal:	stainless steel housing AISI 304 / 1.4301, DIN X5 CrNi 1810

6.2 Model-specific data

6.2.1 XK-KC weighing platforms

	XK155KSD5C	XK604KMC	XK1003KMD2C
Limit values			
Capacity	150 kg	600 kg	1 100 kg
Readability	0.05 g	0.1 g	2 g
Repeatability (at nominal load)	0.12 g (100 kg)	0.3 g (500 kg)	6 g (1000 kg)
Repeatability (at 5% load)	0.1 g (5 kg)	0.2 g (20 kg)	4 g (50 kg)
Repeatability ABA (5 cycles at nominal load)	0.1 g (100 kg)	0.23 g (500 kg)	5 g (1000 kg)
Repeatability ABA (5 cycles at 5% load)	0.08 g (5 kg)	0.18 g (20 kg)	3 g (50 kg)
Linearity deviation	2 g	10 g	20 g
Eccentricity deviation (at test load)	3 g (50 kg)	30 g (200 kg)	40 g (500 kg)
Sensitivity offset (at nominal load) ¹⁾	1 g (100 kg)	7 g (500 kg)	10 g (1000 kg)
Sensitivity temperature drift	0.0003%/°C	0.0003%/°C	0.0003%/°C
Typical values			
Repeatability	0.06 g	0.2 g	4 g
Repeatability ABA (5 cycles at nominal load)	0.05 g	0.18 g	3.5 g
Repeatability ABA (5 cycles at 5% load)	0.03 g	0.09 g	1.2 g
Linearity deviation	1.2 g	6 g	12 g
Eccentricity deviation (at test load)	3 g (50 kg)	20 g (200 kg)	20 g (500 kg)
Sensitivity offset (at nominal load) ¹⁾	0.2 g (100 kg)	1 g (500 kg)	5 g (1000 kg)
Minimum weight (USP, tolerance = 0.10%) ²⁾	80 g	200 g	6000 g
Minimum weight (tolerance = 1%) ²⁾	8 g	20 g	600 g
Settling time	8...12 s	8...12 s	8...12 s
Dimensions & other specifications			
Platform model name	PBK987-CC150	PFK988-C600	PFK988-C600
Platform dimensions (W × D × H)	800 × 600 × 130 - 155 mm	1000 × 800 × 125 - 145 mm	1000 × 800 × 125 - 145 mm
Platform weight	52 kg	91 kg	91 kg

¹⁾ after adjustment with internal weight

²⁾ determined at 5% load, k = 2

	XK1004KMD5C	XK2003KLC	XK6002KLC
Limit values			
Capacity	1 100 kg	2500 kg	5400 kg
Readability	0.5 g	1 g	10 g
Repeatability (at nominal load)	1.7 g (1000 kg)	9 g (2000 kg)	90 g (5000 kg)
Repeatability (at 5% load)	1.5 g (50 kg)	5 g (100 kg)	10 g (200 kg)
Repeatability ABA (5 cycles at nominal load)	1.4 g (1000 kg)	7 g (2000 kg)	70 g (5000 kg)
Repeatability ABA (5 cycles at 5% load)	1 g (50 kg)	3 g (100 kg)	6 g (200 kg)
Linearity deviation	20 g	100 g	300 g
Eccentricity deviation (at test load)	40 g (500 kg)	120 g (1000 kg)	200 g (2000 kg)
Sensitivity offset (at nominal load) ¹⁾	10 g (1000 kg)	60 g (2000 kg)	50 g (5000 kg)
Sensitivity temperature drift	0.0003%/°C	0.0003%/°C	0.0003%/°C
Typical values			
Repeatability	0.9 g	7 g	40 g
Repeatability ABA (5 cycles at nominal load)	0.8 g	5 g	30 g
Repeatability ABA (5 cycles at 5% load)	0.2 g	1 g	2 g
Linearity deviation	12 g	60 g	220 g
Eccentricity deviation (at test load)	20 g (500 kg)	30 g (1000 kg)	80 g (2000 kg)
Sensitivity offset (at nominal load) ¹⁾	2 g (1000 kg)	6 g (2000 kg)	15 g (5000 kg)
Minimum weight (USP, tolerance = 0.10%) ²⁾	600 g	4000 g	6000 g
Minimum weight (tolerance = 1%) ²⁾	60 g	400 g	600 g
Settling time	8...2 s	8...12 s	8...12 s
Dimensions & other specifications			
Platform model name	PFK988-C600	PFK988-E1500	PFK988-E1500
Platform dimensions (W × D × H)	1000 × 800 × 125 - 145 mm	1500 × 1250 × 185 - 205 mm	1500 × 1250 × 185 - 205 mm
Platform weight	91 kg	345 kg	345 kg

¹⁾ after adjustment with internal weight

²⁾ determined at 5% load, k = 2

6.2.2 Terminal

IND930, desk mount	
Display	
Type	Active TFT color LED with backlight
Resolution	1280 × 800 pixels at a size of 10.1" (25.7 cm)
Dimensions and weight	
Terminal dimensions (W × D × H)	259 × 320 × 241 mm
Terminal weight	5.5 kg
Other specifications	
Maximum number of balances to be operated simultaneously	up to 4 (additional interfaces required) ¹⁾
Interfaces	Ethernet, USB ¹⁾

¹⁾ For additional interfaces, please contact your authorized METTLER TOLEDO dealer or service representative.

7 Accessories

Accessories are additional components that could help you in your workflow.

	XK155KSD5C	XK604KMC	XK1003KMD2C XK1004KMD5C	XK2003KLC XK6002KLC
Draft shield XP-WKS	✓	–	–	–
Draft shield XP-WKM	–	✓	✓	–
Draft shield XP-WKL	–	–	–	✓
LevelMatic 1000	✓	✓	–	–
LevelMatic 5000	–	–	✓	✓

Description

Order no.

Draft shields

Draft shield XP-WKS for KS models

11116556

Draft shield XP-WKM for KM models

11116557

Draft shield XP-WKL for KL models

11116558

LevelMatic

LevelMatic 1000

22001940

LevelMatic 5000

11116554

8 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

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	XK155KSD5C	XK604KMC	XK1003KMD2C XK1004KMD5C	XK2003KLC XK6002KLC
Pare-brise XP-WKS	✓	–	–	–
Pare-brise XP-WKM	–	✓	✓	–
Pare-brise XP-WKL	–	–	–	✓
LevelMatic 1000	✓	✓	–	–
LevelMatic 5000	–	–	✓	✓

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▶ www.mt.com/contact

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原文	译文	说明
EMC		Electromagnetic Compatibility (电磁兼容)
FACT		Fully automatic time- and temperature-controlled internal adjustment (全自动时间和温度控制的内部校正)
FCC		Federal Communications Commission (美国联邦通讯委员会)
GWP		Good Weighing Practice
UM		User Manual (简明用户手册)
USB		Universal Serial Bus
USP		United States Pharmacopeia (美国药典)

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© 04/2022 METTLER TOLEDO. All rights reserved.