English Reference Manual **Printer** P25, P26, P28
Deutsch Referenzhandbuch **Drucker** P25, P26, P28
Español Manual de referencia **Impresora** P25, P26 y P28
Français Manuel de référence **Imprimante** P25, P26, P28
Italiano Manuale di riferimento **Stampante** P25, P26, P28

METTLER TOLEDO

Reference Manual **Printer**

English

Referenzhandbuch **Drucker**

Deutsch

Manual de referencia **Impresora**

Español

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Français

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Italiano

1 Introduction

This printer is a simple to operate dot matrix printer for METTLER TOLEDO instruments. The printout on normal paper meets the requirements of modern quality assurance systems (GLP, GMP, ISO 9001, etc.). These high-value printers guarantee long-term traceability.

Feature	RS-P25	USB-P25	RS-P26	RS-P28
Auto Bauderate detection ¹⁾	✓	*	✓	✓
Balance Feature on/off ²⁾	✓	*	✓	✓
USB interface	–	✓	–	–
RS232 interface	✓	–	✓	✓
Zero / Tare button	–	–	✓	✓
Date / Time	–	–	✓	✓
Statistics / Totalling	–	–	–	✓

* USB has standard settings.

¹⁾ Auto Bauderate detection: the printer detects the settings of the connected device every time it is switches on/off by the key. Switch this feature off if the connected device has this feature not implemented, see chapter [Troubleshooting ▶ Page 16].

²⁾ Balance Feature on/off: this feature disables the keys on the printer, see chapter [Troubleshooting ▶ Page 16].

This document is based on the software version V 3.01.

Finding more information

▶ www.mt.com/lab-printers

Search for documents

▶ www.mt.com/library

For further questions, please contact your authorized METTLER TOLEDO dealer or service representative.

▶ www.mt.com/contact

Conventions and symbols

Key and/or button designations and display texts are shown in graphic or in bold text (e.g. ↑).

Note For useful information about the product.

This symbol indicates press key briefly (less than 1.5 s).

This symbol indicates press and hold key down (longer than 1.5 s).

Elements of instructions

- Prerequisites
- 1 Steps
- 2 ...
 - ⇒ Intermediate results
 - ⇒ Results

2 Safety Information

Two documents named "User Manual" and "Reference Manual" are available for this instrument.

- The User Manual is printed and delivered with the instrument.
- The electronic Reference Manual contains a full description of the instrument and its use.
- Keep both documents for future reference.
- Include both documents if you transfer the instrument to other parties.

Only use the instrument according to the User Manual and the Reference Manual. If you do not use the instrument according to these documents or if the instrument is modified, the safety of the instrument may be impaired and Mettler-Toledo GmbH assumes no liability.

2.1 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

Signal words

WARNING A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.

CAUTION A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.

NOTICE A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

Warning symbols

Electrical shock

General hazard: read the Operating Instructions for information about the hazards and the resulting measures.

Notice

2.2 Product specific safety notes

Intended use

This instrument is designed to be used in laboratories by trained staff. The instrument is intended for print out protocol purposes.

Any other type of use and operation beyond the limits of technical specifications without written consent from Mettler-Toledo GmbH is considered as not intended.

Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

METTLER TOLEDO assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. METTLER TOLEDO assumes that the instrument owner provides the necessary protective gear.

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the approved METTLER TOLEDO power supply cable and AC/DC adapter with a current-limited SELV output.
- 2 Connect the power cable to a grounded power outlet, ensure correct polarity.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and power plug for damage and replace damaged cables and power plugs.

NOTICE


Damage to the instrument due to the use of unsuitable parts

Using unsuitable parts with the instrument can damage the instrument or cause it to malfunction.

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

3 Design and Function

3.1 Overview printer

1	On/Off Switch	2	Pilot lamp – status indicator
3	Control panel	4	Interface connector
5	Power connection		

3.2 Overview functions

Key	P25	P26	P28
	<ul style="list-style-type: none"> Paper feed Menu option (+ On switch) 	<ul style="list-style-type: none"> Paper feed Scroll up menu option 	<ul style="list-style-type: none"> Paper feed Scroll up menu option
	–	<ul style="list-style-type: none"> Zeroing / taring the balance Abort – exit 	<ul style="list-style-type: none"> Zeroing / taring the balance Abort – exit
	–	<ul style="list-style-type: none"> Printout of date and time Scroll down menu option 	<ul style="list-style-type: none"> Printout of date and time Scroll down menu option
	–	<ul style="list-style-type: none"> Printing of stable weight values Opening the menu – select next menu option Save settings 	<ul style="list-style-type: none"> Printing of stable weight values Opening the menu – select next menu option Save settings
	–	–	<ul style="list-style-type: none"> Totaling application
	–	–	<ul style="list-style-type: none"> Statistics application

4 Installation and Putting into Operation

4.1 Scope of delivery

- AC/DC adapter
- Power cable (country specific)
- User Manual
- Declaration of conformity

4.2 Connecting the printer to the power supply

WARNING

Death or serious injury due to electric shock

Contact with parts that carry a live current can lead to death or injury.

- 1 Only use the approved METTLER TOLEDO power supply cable and AC/DC adapter with a current-limited SELV output.
- 2 Connect the power cable to a grounded power outlet, ensure correct polarity.
- 3 Keep all electrical cables and connections away from liquids and moisture.
- 4 Check the cables and power plug for damage and replace damaged cables and power plugs.

NOTICE

Damage to the AC/DC adapter due to overheating


If the AC/DC adapter is covered or in a container, it is not sufficiently cooled and will overheat.

- 1 Do not cover the AC/DC adapter.
- 2 Do not put the AC/DC adapter in a container.

The printer is supplied with a universal AC/DC adapter or an AC/DC adapter with a country-specific power cable.

- Install the cables so that they cannot be damaged or interfere with operation.
- Insert the power cable in a grounded power outlet that is easily accessible.

- 1 Connect the AC/DC adapter to the connection socket on the back of the printer.
- 2 Lock the connector with the screw nut.
- 3 Connect the 3-pin grounded power cable to the socket for the AC/DC adapter.
- 4 Connect the printer and the balance with the interface cable.
- 5 Switch the balance on.
- 6 Connect the printer to the power supply.
- 7 Switch on the printer with the **I/O** switch.


 **Note**

RS232 interface only: the pilot lamp will flash until an instrument is connected. If the lamp is still flashing check the chapter [Troubleshooting ▶ Page 16].

4.3 Inserting paper

- 1 Remove paper cover (pull on the back) (1).
- 2 Insert the paper rotary axis through the core of the paper roll (2).
- 3 Feed paper through the slot in the printing unit in the direction of the arrow (3).
- 4 Press and hold until enough paper has been fed (4).
- 5 Lead paper through the paper cover slot.
- 6 Place the paper cover (1) on.

If, after the first printouts the connection settings (Bauderate, Bit, Handshake) of the balance change, the printer must be restarted. Switch off / on to start the auto bauderate defection again. Auto bauderate can be switched off.

5 Keystroke Functions

Printer P25

Operation level	 Press briefly	 Press and hold down
	one line feed	continuous feed (until release)
  	Printer test with setting printout	open menu
Menu level	 Press briefly	 Press and hold down
	next menu item	store

Printer P26

Operation level	 Press briefly	 Press and hold down		
	zero / tare the balance	—		
	print Date and Time	—		
	one line feed	continuous feed (until release)		
	print	open menu		
Menu level	 Press briefly	 Press and hold down		
 exit	 down	 up	 next	 store

Printer P28

Operation level	 Press briefly	 Press and hold down
	zero / tare the balance	–
	print Date and Time	–
	one line feed	continuous feed (until release)
	print	open menu
Menu level	 Press briefly	 Press and hold down
 exit  down  up  next	 store	
Function level	 Press briefly	 Press and hold down
	Start totaling application	Conclude totaling application
	Start statistics application	Conclude statistics application


6 Menu

6.1 Printer P25

 Press briefly	 Press and hold down	Printout
<div style="border: 1px solid black; padding: 2px; text-align: center;">Settings</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Reset</div>	Reset the printer	Settings Reset ----- Reset done ----- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS
<div style="border: 1px solid black; padding: 2px; text-align: center;"> RS232 interface only Auto Baudrate On </div>		Store
<div style="border: 1px solid black; padding: 2px; text-align: center;">Auto Baudrate Off</div>	Store	Auto Baudrate: Off -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Char Set IBM / DOS</div>	Store	IBM/DOS -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS Cyrillic</div>	Store	IBM/DOS Cyrillic -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS Japanese</div>	Store	IBM/DOS Japanese -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Latin 1</div>	Store	ANSI/WIN Latin 1 -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Cyrillic</div>	Store	ANSI/WIN Cyrillic -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Japanese</div>	Store	ANSI/WIN Japanese -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Enable</div>	Store	Balance Feature: Enable -----STORED----- ..."current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Disable</div>	Store	Balance Feature: Disable -----STORED----- ..."current settings"

Menu Option
Factory Setting

6.2 Printer P26 / P28

Menu Option
Factory Setting

- * factory settings
- ** instruments with built-in clock: not applicable, built-in clock has priority
- *** instruments with built-in clock: format of printer has priority

Examples if balance display is used:

Settings	
Reset	rESEt
List	LISt
Header	
*No Header	no Hdr
Date/Time	d t Hdr
Date/Time/Balance	dtb Hdr
Date Day**	
18	dRy 18
Date Month**	
07	Month 07
Date Year**	
2008	YEAr2008
Time Minute**	
28	Min. 28
Time Hour**	
14	Hour 14
Date Format***	
DD.MMM.YYYY	ddMmY
MMM.DD.YYYY	Mm dd Y
DD.MM.YYYY	ddMmY
MM/DD/YYYY	Mm dd Y
YYYY-MM-DD	Y-Mm-dd
Time Format***	
24:MM:SS	24:MmSs
12:MM:SS	12:MmSs
24.MM.SS	24MmSs
12.MM.SS	12MmSs
Baud rate	
*9600	bd9600
19200	bd.19200
38400	bd.38400
600	bd600
1200	bd.1200
2400	bd.2400
4800	bd.4800
Bit/Parity	
*8 Bit No	8b-no
7 Bit Odd	7b-odd
7 Bit Even	7b-E
7 Bit No	7b-no
Handshake	
*Xon/Xoff	HS oFF
HS Soft	HS SoFt
HS Hard	HS HARd
Char Set	
*IBM/DOS	dOS
IBM/DOS Cyrillic	dOS Cyr
IBM/DOS Japanese	dOS JRP
Ansi/Win Latin 1	11 in LRt
Ansi/Win Cyrillic	11 in CRt
Ansi/Win Japanese	11 in JRP
Tare/Zero	
*Zero	ZEro
Tare	TARe

Example for setting the day (P26 / P28 only)

Open the menu

- Hold down the  key until **Reset** appears.
- Select the menu option **Day** with the  key.
- Select the day, e.g. day 21.
 Scroll up with the \uparrow key.
 Scroll down with the  key.

Save the settings

- Hold down the  key until **Stored** appears.

Statistics application (P28 only)

Sample weight must be in the range of 70% - 130% of the current average value.
 Minimum dynamic weight changed must be 100 d.

- 1 Press  to start the Statistics application (LED on).
- 2 Add a new sample and press **M+**.
- 3 Press  to print subtotal, if needed.
- 4 Press and hold  to print and end the Statistics application (LED off).


Note

- Unit change during the application is not permitted.
- Maximal 999 sample.

Totalling application (P28 only)

Minimum weights + 1 d.

- 1 Press  to start the Totalling application (LED on).
- 2 Add a new sample and press **M+**.
- 3 Press  to print subtotal, if needed.
- 4 Press and hold  to print and end the Totalling application (LED off).

Note

- Unit change during the application is not permitted.
- Maximal 999 sample.

6.3 Printer Test

The printer is equipped with a self-test that automatically prints out the character set of the printer (factory settings IBM/DOS).

Starting the test

- Press  and power on the printer.
 ⇒ The printer prints out the current settings.

Terminating the test

- Power off the printer.

7 Maintenance

WARNING

Death or serious injury due to electric shock

Contact with parts carrying a live current can lead to injury and death.

- 1 Disconnect the printer from the power supply prior to cleaning and maintenance.
- 2 Prevent liquid from entering the printer or AC/DC adapter.

7.1 Replacing the paper roll

The insertion of the new paper roll is described in chapter [Inserting paper ▶ Page 8]. The order number for the paper roll can be found in chapter [Consumables ▶ Page 21].

Note

The paper rotary axis must always be transferred from the old paper roll to the new one.

7.2 Replacing the ribbon

The order number for the ribbon cartridge can be found in chapter [Consumables ▶ Page 21].

- 1 Remove the paper cover (1).
- 2 Pull paper out of printing unit (2).

- 3 Remove ribbon cartridge, by pressing in arrow direction (3).
- 4 Insert the new ribbon cartridge. If necessary, take up slack with tensioning wheel (4).

Insertion of a new paper roll is described in chapter [Inserting paper ▶ Page 8].

7.3 Cleaning the housing

Since the printer housing is made of top grade, resistant materials, all commercially available mild cleaning agents may be used.

NOTICE

Damage to printer due to inappropriate cleaning methods

The printer housing is made from high quality, resistant materials and can be damaged by certain cleaning agents, solvents or abrasives. Any liquid that enters the housing may damage the printer.

- 1 Use water and a mild detergent to clean the printer.
- 2 Wipe off any spills immediately.
- 3 Prevent liquid from entering the interior of the printer.

Contact a METTLER TOLEDO representative to find about the service options available – regular maintenance by an authorized service engineer will ensure consistent weighing accuracy over the long term and extend the service life of the instrument.

8 Troubleshooting

Problem / Message	Possible cause	Remedy
Printout unreadable	<ul style="list-style-type: none"> Ribbon worn or entangled. Service life of printing unit at an end. 	<ul style="list-style-type: none"> Change ribbon cartridge or tense the ribbon properly. Please contact METTLER TOLEDO service.
Pilot lamp does not light	<ul style="list-style-type: none"> No power. Power supply fault. 	<ul style="list-style-type: none"> Switch printer on. Please contact METTLER TOLEDO service.
Pilot lamp flashes ---Bauderate not found--- *Please connect device and switch the printer off/on. *Alternatively you can change the settings of the printer or the connected device. Printer settings: Baudrate 9600 Bit/Parity 8 Bit None Handshake Off	<ul style="list-style-type: none"> No connection printer device. Printer did not find the settings of the connected instrument. 	<ul style="list-style-type: none"> Make sure that the printer and the instrument is connected. Switch the printer off/on. Alternatively change the settings of the printer or the connected instrument, see chapter [Recommended Settings and Printer Models ▶ Page 22].
-----Timeout -----	<ul style="list-style-type: none"> Data transfer has not completed in the last 45 s. 	<ul style="list-style-type: none"> Move balance to more stable location. Check the balance.
----- Over Load -----	<ul style="list-style-type: none"> Balance in overload range. Zeroing out of the zero range. 	<ul style="list-style-type: none"> Check the weight of sample. Change the zero range of the balance. Change the setting of the printer from zero to tare.
----- Under Load -----	<ul style="list-style-type: none"> Balance in underload range. Weighing pan is not in position. 	<ul style="list-style-type: none"> Check the balance.
-- Transmission Error --	<ul style="list-style-type: none"> Instrument does not understand printer instruction. 	<ul style="list-style-type: none"> Check interface parameter settings on printer and instrument.
---- Not Executable ----	<ul style="list-style-type: none"> Instrument can not execute printer instruction. 	<ul style="list-style-type: none"> Check the balance.
--- Date/Time Not Set -- ----- Memory Error -----	<ul style="list-style-type: none"> Printer error messages. 	<ul style="list-style-type: none"> Check the instrument. Reset the printer. Set date/time. Back-up battery depleted. Please contact METTLER TOLEDO service.
----- Memory full -----	<ul style="list-style-type: none"> Statistics or Totaling memory is full. 	<ul style="list-style-type: none"> Start a new Statistics or Totaling application.
-- Load/Unload Weight --	<ul style="list-style-type: none"> No sample or the same sample is on the weighing pan. 	<ul style="list-style-type: none"> Put a new sample or remove it and put it again on the pan.
Unit change not allowed !	<ul style="list-style-type: none"> Changed unit during Statistics or Totaling application 	<ul style="list-style-type: none"> Reset the balance. Change the unit back to the first one.

Problem / Message	Possible cause	Remedy
-- Out of range --	<ul style="list-style-type: none">• Sample weight is out of range.	<ul style="list-style-type: none">• Sample weight must be in the range of 70% – 130% of the current average value.

9 Technical Data

9.1 General technical data

Power supply

AC/DC adapter:	Primary: 100 – 240 V, 50/60 Hz Secondary: 12 V DC, 2.5 A (with electronic overload protection)
Cable for AC/DC adapter:	3-core, with country-specific plug
Polarity:	 with a current limited SELV (Safety Extra Low Voltage) output
Power line voltage:	12 V DC
Power consumption:	Max. 2.5 A

Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Standards for safety and EMC:	See Declaration of Conformity
Range of application:	For use in closed interior rooms only

Environmental conditions

Height above mean sea level:	Up to 4000 m
Ambient temperature:	0 to 40 °C
Storage condition:	-25 to 60 °C
Relative air humidity:	Max. 80% up to 31 °C, linearly decreasing to 50% at 40 °C, noncondensing

Printer features

Print technology:	Dot matrix printer 5 × 7, 24-character line length
Print speed:	2.3 lines per second
Ribbon cartridge:	Exchangeable, black
Paper roll:	Standard paper 57.5 mm × ø 50 to 60 mm, integrated in housing, commercial size
Print quality:	Light resistant and thermally stable printing (GLP, GMP, ISO 9001)
Interfaces:	RS-P25, RS-P26, RS-P28 USB-P25
Compatibility limits:	RS232C USB <ul style="list-style-type: none">• Function Date / Time: not applicable on balances with built-in clock (built-in clock has priority)• Function Zeroing / Taring: not applicable on all models of balances


Dimension / weight

Printer dimensions (W × D × H):	120 × 203 × 73 mm
Packaging dimensions (W × D × H):	255 × 205 × 150 mm
Net weight:	740 g (paper roll incl.)

9.2 RS232C Interface






RS-P25, RS-P26 and RS-P28 printers are equipped with an RS232C interface to connect METTLER TOLEDO instruments.

- 9-pin male connector.
- Matching to other device (transmission parameters), see chapter [Menu ▶ Page 11].




10 Accessories and Spare Parts


10.1 Accessories

	Description	Part No.
Cables for RS232 interface		
	RS9 – RS9 (m/f): connection cable for PC, length = 1 m	11101051
	RS9 – RS9 (f/f): connection cable, length = 1 m	51190362
	USB (A – B) connection cable for PC, length = 1 m	30241476
Various		
	AC/DC adapter (without power cable) 100–240 V AC, 0.8 A, 50/60 Hz, 12 V DC 2.5 A	11107909
	Country-specific 3-Pin power cable with grounding conductor.	
	Power cable AU	00088751
	Power cable BR	30015268
	Power cable CH	00087920
	Power cable CN	30047293
	Power cable DK	00087452
	Power cable EU	00087925
	Power cable GB	00089405
	Power cable IL	00225297
	Power cable IN	11600569
	Power cable IT	00087457
	Power cable JP	11107881
	Power cable TH, PE	11107880
	Power cable US	00088668
	Power cable ZA	00089728

10.2 Consumables

	Description	Part No.
	Paper roll (length: 20 m), set of 5 pcs	00072456
	Paper roll (length: 13 m), self-adhesive, set of 3 pcs	11600388
	Ribbon cartridge, black, set of 2 pcs	00065975

10.3 Spare parts

	Order no.	Designation	Remarks
1	12120734	Paper cover	–
2	12120735	Paper rotary axis	–
3			

11 Recommended Settings and Printer Models


Connected instrument		Instrument settings	Printer settings	Compatible printer models
MS, ML, JP, JS, PHS, PHL, AB-S/FACT, AB-S, AB-L, PB-S/FACT, PB-S, PB-L, AL, PL, PL-S, PL-L, JB-G/FACT, JB-C/FACT, JB-L-G, JB-L-C, JL-G, JL-C, JL-G/L, EL, AB-S/PH, PB-S/PH		Printer mode	Factory settings	RS-P25, RS-P26, RS-P28
PG-S, PG-S/PH		Printer mode	Factory settings	RS-P25
SB		Printer mode	Factory settings	RS-P26
XP, XS, XA		Printer mode	Factory settings	RS-P25
	Katakana	Printer mode and language ANSI/WIN Katakana	Factory settings	
	Russian	Printer mode and language ANSI/WIN Russian	CharSet IBM/DOS Cyrilic	
	Chinese, Japanese	Printer mode	Factory settings	
MJ33, HB43-S		Printer mode	Factory settings	RS-P25
	Nihongo	Printer Mode and language Nihongo	Factory settings	
	Russian	Printer Mode and language Russian	CharSet IBM/DOS Cyrilic	
HG63, HR83		Select on LocalCAN Switchbox: 1 Printer, Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff	*Reset, Auto Bauderate off, Balance Feature off	RS-P25
	Russian	Select on LocalCAN Switchbox: 1 Printer, Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff and language Russian	*Reset, Auto Bauderate off, Balance Feature off, CharSet IBM/DOS Cyrilic	
QB5		Factory settings	Factory settings	RS-P25
SQC XP		Printer mode	Factory settings	RS-P25
SQC16		Printer mode	*Reset, Auto Bauderate off	RS-P25
T50, T70, T90		USB compact printer	Factory settings	USB-P25
V20, V30, C20, C30		USB compact printer	Factory settings	USB-P25
DL15, DL22, DL28		RS-P42	*Reset, Auto Bauderate off, Balance Feature off	RS-P25 (connect with RS9 F/F cable order number: 51190362)
DL31, DL38, DL32, DL39		GA42	Auto Bauderate off, Balance Feature off, Baudrate: 1200; Bit/Parity: 8 Bit No, Handshake: None	RS-P26 (connect with RS9 F/F cable order number: 51190362)

Connected instrument	Instrument settings	Printer settings	Compatible printer models
DE40, DE45, DE51, RE40D, RE50, DR40, DR45	GA42	Auto Bauderate off, Balance Feature off, Baudrate: 1200; Bit/ Parity: 8 Bit No, Handshake: Hardware	RS-P26 (connect with RS9 F/ F cable order number: 51190362)
Refracto 30PX/GS, Densito 30PX	Printer mode and Baudrate 9600	Auto Bauderate off, Balance Feature off, Baudrate: 1200; Bit/ Parity: 8 Bit No, Handshake: Hardware	RS-P26
SevenEasy S20, S30, V 1.30 or higher	Factory settings	*Reset, Auto Bauderate off, Balance Feature off	RS-P26, RS-P25
SevenMulti S40, S47, S50, S70, S80, V 1.1 or higher	Baudrate: 9600; Bit/ Parity: 8 Bit No, Handshake: Xon/Xoff	*Reset, Auto Bauderate off, Balance Feature off	RS-P26, RS-P25
SevenGo pro SG6, SG7, SG8	Factory settings of the RS232 Infrared device	Auto Bauderate off, Baudrate: 1200, Bit/ Parity: 8 Bit No, Handshake: off	RS-P26 data transfer of SevenGo pro to the printer occurs via the Infrared-RS232 adapter (order number: 51302333)
DM40, DM45, DM50, RM40, RM50	Printer mode	Factory settings	USB-P25
MP50, MP70, MP80	Printer mode	Factory settings	USB-P25

* The printer detects automatically the settings of the connected device. Therefore disconnect the printer of the connected device and reset the printer to set the factory settings again: Baudrate: 9600; Bit/Parity: 8 Bit No, Handshake: Xon/Xoff

12 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions, please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties, the content of this regulation must also be related.

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6 Menu

6.1 Imprimante P25

Appuyez brièvement	Appuyez et maintenez la touche enfoncée.	Ticket d'impression
<div style="border: 1px solid black; padding: 2px; text-align: center;">Settings</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Reset</div>	Réinitialiser l'imprimante	Settings Reset ----- Reset done ----- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS
Interface RS232 uniquement		
<div style="border: 1px solid black; padding: 2px; text-align: center;">Auto Bauderate</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">On</div>	Enregistrer	Auto Bauderate: On -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Auto Bauderate</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Off</div>	Enregistrer	Auto Bauderate: Off -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Char Set</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS</div>	Enregistrer	IBM/DOS -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Cyrillic</div>	Enregistrer	IBM/DOS Cyrillic -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Japanese</div>	Enregistrer	IBM/DOS Japanese -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Latin 1</div>	Enregistrer	ANSI/WIN Latin 1 -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Cyrillic</div>	Enregistrer	ANSI/WIN Cyrillic -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Japanese</div>	Enregistrer	ANSI/WIN Japanese -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Enable</div>	Enregistrer	Balance Feature: Enable -----STORED----- ... "current settings"
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Disable</div>	Enregistrer	Balance Feature: Disable -----STORED----- ... "current settings"

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6 Menu

6.1 Stampante P25

 Pressione breve	 Premere e tenere premuto	Stampa
<div style="border: 1px solid black; padding: 2px; text-align: center;">Settings Reset</div>	Ripristino della stampante	Settings Reset
<div style="border: 1px solid black; padding: 2px; text-align: center;">Solo per interfaccia RS232</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Auto Bauderate On</div> <div style="border: 1px solid black; padding: 2px; text-align: center;">Auto Bauderate Off</div>		<pre>----- Reset done ----- Current settings: Baudrate: 9600 Bit/Parity: 8 Bit Non Handshake: Off Char Set: IBM/DOS</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Char Set IBM / DOS</div>	Memorizzare	<pre>Auto Bauderate: On -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS Cyrillic</div>	Memorizzare	<pre>Auto Bauderate: Off -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">IBM / DOS Japanese</div>	Memorizzare	<pre>IBM/DOS -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Latin 1</div>	Memorizzare	<pre>IBM/DOS Cyrillic -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Cyrillic</div>	Memorizzare	<pre>IBM/DOS Japanese -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">ANSI / WIN Japanese</div>	Memorizzare	<pre>ANSI/WIN Latin 1 -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Enable</div>	Memorizzare	<pre>ANSI/WIN Cyrillic -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Disable</div>	Memorizzare	<pre>ANSI/WIN Japanese -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Disable</div>	Memorizzare	<pre>Balance Feature: Enable -----STORED----- ... "current settings"</pre>
<div style="border: 1px solid black; padding: 2px; text-align: center;">Balance Feature Disable</div>	Memorizzare	<pre>Balance Feature: Disable -----STORED----- ... "current settings"</pre>

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© Mettler-Toledo GmbH 12/2018