# **InPro 7100 (i) Conductivity Sensors**

## **Instruction Manual**

cs	Návod k obsluze	2
da	Brugsanvisning	5
de	Betriebsanleitung	8
en	Instruction Manual	12
es	Manual de instrucciones	15
fi	Käyttöohje	18
fr	Instructions d'utilisation	21
hu	Használati utasítás	25
it	Istruzioni per l'uso	29
ja	取扱説明書	32
ko	사용 매뉴얼	35
nl	Gebruikershandleiding	38
pl	Instrukcja obsługi	42
pt	Manual de instruções	45
ru	Инструкция по эксплуатации	48
sv	Bruksanvisning	52
th	คู่มือคำแนะนำ	55
zh	说明书	58

# InPro 7100 (i)

## Senzory vodivosti

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#### 1 Úvod

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## Instruction Manual

### 1 Introduction

The METTLER TOLEDO InPro™ 7100 (i) sensors are low-maintenance, CIP and SIP resistant four-electrode conductivity sensors for conductivity measurements in aqueous solutions. The sensors have a Pg 13.5 thread and a diameter of 12 mm. Please read these operating instructions fully and carefully before using the sensor in order to ensure proper, trouble-free use. We recommend using the conductivity sensor only in combination with original METTLER TOLEDO components. It should only be used and serviced by trained personnel and employees who have read and understood the operating instructions.

### 2 Safety instructions

You should wear safety glasses or goggles and protective gloves when cleaning or calibrating the conductivity sensor with an acid or alkali solution.

### 3 Product description

The printed markings on each InPro 7100 (i) conductivity sensor contain the following information:

<b>METTLER TOLEDO</b>	Manufacturer of the electrode
<b>InPro 7100 (i)/12/XXX/XXXX</b>	Type designation/diameter/length/pin material
Cond. 0.02...500 mS/cm	Sensor conductivity measuring range
Order No. 52 00X XXX	Order number

Each conductivity sensor is also marked with a serial number located on the Pg 13.5 threaded fitting. The operating temperature range of the sensors in the InPro 7100 (i) series is –20 to 150 °C.

### 4 Installation and initial use

1. When you unpack the conductivity sensor, check it for possible damage. If you find any damage, report it immediately to the carrier or shipping agent and your supplier.

2. Each conductivity sensor in the InPro 7100i series has an integrated NTC temperature sensor. InPro 7100 series has an integrated Pt1000 temperature sensor (IEC 751 Class A).
3. The nominal cell constant is  $0.31 \text{ cm}^{-1}$ . The exact value is printed on the certificate.
4. An addition factor is also noted on the certificate. If a M300 or M400 transmitter is used:
  - InPro 7100i: this factor is automatically transferred to the transmitter.
  - InPro 7100: this factor must be entered under "Cell A" in the Quick Setup Menu.
5. The cell constant of each sensor is measured individually and printed on the accompanying quality certificate. All calibrations must be traceable to NIST and/or ASTM standards.
6. Install the conductivity sensor such that the medium is in direct contact with the electrodes. Any other installation is a risk of air inclusions or the formation of hard contaminant deposits. The conductivity cell should ideally be installed with a distance of at least 30 mm between the pipe and the end of the conductivity cell.
7. Vertical installation is only possible if the pipe is filled and there is no possibility of bubble formation. In case of installation at the side, the direction of flow should be upward.
8. **Note** that the measuring range depend strongly on the transmitter that is used.

## 5 Operation

For each sensor, the cell constant and temperature constant have been measured/established individually and the values documented in the accompanying Quality Certificate. In general, the cell constants and temperature constant will no change. It is not necessary to perform sensor cell constant correction or calibration with standard solutions. However, depends on plant operating procedure, cell constant can be verified/corrected after the period of initial calibration has expired by using conductivity standard solutions.

## 6 Maintenance

Contaminated conductivity sensor can give false indications. If contamination is suspected, dismount the conductivity sensor and clean the electrodes and the insulation between the electrodes with a soft cloth. A mild detergent or a highly diluted solution (less than 0.5 % by weight) of an acid, such as nitric acid, is suitable to use as a cleaning solution (always handle acids with caution!).

Only use cleaning agents that are compatible with the materials to be cleaned. In addition, do not use coarse materials such as steel wool for cleaning, as they will damage the surface.

Rinse the conductivity cell thoroughly with distilled water before installation.

## 7 Troubleshooting

1. If no signal appears at the display unit, check the electrical connections. Another possible cause of this problem is that the conductivity cell is not in contact with the medium.
2. If no temperature signal is present, the temperature sensor may not be connected properly or the device may not support NTC temperature sensor (for InPro 7100 i) or Pt1000 temperature sensor (for InPro 7100).
3. If the measuring signal is incorrect or instable, the cause may be that the conductivity cell is contaminated or fitted too close to the pipe (the distance between the pipe and the end of the conductivity cell should be 30 mm).

## 8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

## 9 Warranty

Warranted free of manufacturing defects for a term of 12 months after delivery.

## 10 Storage conditions

If the conductivity cell is not in use, it can be stored dry at room temperature.

## 11 EC declaration of conformity

The complete declaration is available in the declaration of conformity certificate.

InPro and ISM are trademarks of the METTLER TOLEDO Group.

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All rights reserved