METTLER TOLEDO

InPro 4550 VP pH Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Instruktionsvejledning	6
de	Betriebsanleitung	10
en	Instruction Manual	14
es	Manual de instrucciones	18
fr	Instructions d'utilisation	22
fi	Käyttöohje	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

InPro 4550 VP pH elektroda

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Obsah

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Instruction Manual

Contents

1	Introduction	15
2	Safety instructions	15
3	Product description.....	15
4	Installation and commissioning	16
5	Cabling.....	16
6	Calibration and measurement.....	16
7	Maintenance	16
8	Diagnostics.....	17
9	Trouble-shooting (sluggish response time, insufficient sensitivity)	17
10	Disposal	17
11	Guarantee.....	17

InPro and Xerolyt are registered trademarks
of the METTLER TOLEDO Group in Switzerland, the USA,
the European Union and a further five countries.

1 Introduction

The METTLER TOLEDO electrode versions InPro™ 4550VP are low-maintenance, temperature-compensated pH combination electrodes containing polymer electrolyte (Xerolyt® Plus). An integrated temperature sensor (RTD) provides automatic compensation of the electrode slope to ensure accurate readings.

The InPro 4550VP electrodes are suitable for pH measurements in:

- emulsions and suspensions
- sulfide-bearing media
- solutions containing heavy metals (lead, mercury, silver, etc.)
- processes with rapid pressure fluctuations (within the given range)
- heavily contaminated media.

The InPro 4550VP electrodes are unsuitable for sterilization/autoclaving procedures, for processes with rapid temperature fluctuations or food and beverage applications.

Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use. These general instructions must be readily available to operators of the electrode at all times.

2 Safety instructions

The InPro 4550VP electrodes are intended solely for the purpose described in §1 above. The electrodes consist partly of glass and should be handled carefully. If the electrode is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn. Local safety regulations are applicable to the handling of the electrode.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

3 Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
InPro™ 4550 VP/xxx/SG	type designation/xxx = type of temperature sensor (Pt 100/ Pt 1000)/SG = ancillary electrode
combination pH	Type of electrode
pH 0 ... 14	pH measurement range
0 ... 130 °C	Temperature range for operation
Order No. 52 002 40X	Order number
 SEV14 ATEX 0168 X	Identification and number of certificate
 IECEx SEV 14.0025X	Identification and number of certificate
 CML 22 UKEX 2108X	Identification and number of certificate
IS/I, II, III/1 / ABCDEFG/T6	FM approved
See instruction manual!	Instruction (Observe the operating instructions)

In addition each electrode is allocated a serial number on the plug-in head to enable identification.

4 Installation and commissioning

1. On unpacking, and again prior to installation, check the electrode visually for signs of mechanical damage to the glass membrane, the shaft and the connector. Report any faults or defects immediately to your authorized METTLER TOLEDO dealer.
2. Remove the watering cap and rinse the electrode with distilled or deionized water. Use only soft tissue to dab the electrode dry. Do not rub the glass membrane since this can lead to sluggish response behavior of the electrode.
3. Electrodes with a dry glass membrane must be conditioned for several hours in a 3 mol KCl solution before being used.
4. Connect always the solution ground (blue cable) to the transmitter.
5. Use PTFE tape to seal the NPT thread.
6. If possible use the protection sleeve (order no. 52 401 808).

5 Cabling

The InPro 4550VP is equipped with a VarioPin connector. To connect the cable to the electrode, attach and gently rotate until the connector is engaged (groove). Screw the cable cap nut onto the electrode until the connection is secure. The plug assembly should be fitted hand-tight only; no tools should be used. For details of the exact assignment of the individual cable leads to the VP pins, please refer to the instruction leaflet "Preparation of the cable for connection to the transmitter" supplied together with the cable.

6 Calibration and measurement

InPro 4550VP electrodes require two-point calibration. Before calibration, the watering cap must be removed from the pH electrode. The electrode is then dipped successively into two different buffer solutions, each of known pH value, and the transmitter calibrated to these buffer values (for example, buffers pH 4.01 and pH 7.00). Please note the instructions given in the operating manual of the transmitter. Calibration should be repeated as and when needed.

7 Maintenance

- Electrodes with Xerolyt Plus polymer electrolyte cannot be refilled.
- Electrodes should not be stored dry for any period of time.
- During storage, the watering cap should be in place, filled with Friscolyt™ or 3 mol KCl solution.
- Never store an electrode in distilled or deionized water, as this shortens the lifetime of the electrode.
- Should the electrode dry out, it must be conditioned for 24 hours in Friscolyt (ordering no. 51 340 053 for 250 ml, or 51 340 054 for 6 3 250 ml) before reuse. After conditioning, recalibration is required.
- After each workcycle, the tip of the electrode should be cleaned using deionized water.

Attention! Do **not** under any circumstance use METTLER TOLEDO reactivating solution (ordering no. 51 319 053/in USA 43 989 5125) as this contains hydrofluoric acid (HF) which would attack the Xerolyt reference system.

8 Diagnostics

The electrodes are equipped with a grounding ring which allows a diagnostic of the diaphragm condition. If this diaphragm (open junction) should become clogged or if the reference resistance exceeds a certain value (see transmitter operation manual), it is then necessary to replace the electrode. InPro 4550VP electrodes should as a general rule not be cleaned with acidic cleaning solutions, e.g. diaphragm cleaning solution, reactivation solution, etc. However, by extremely careful application, it is possible to reactivate the glass membrane, see point 9 below.

9 Trouble-shooting (sluggish response time, insufficient sensitivity)

Attention! The glass membrane must be very cautiously dabbed with HF solution using a fine brush and taking care that solution does not come into contact with the titanium pin (metallic solution ground). Also check the electrode connector (VP) for possible traces of moisture. If necessary, clean well using deionized water or alcohol, and afterwards dry carefully.

10 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

11 Guarantee

On manufacturing defects, 12 months after delivery.

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 本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○ : 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○ : Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x : 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x : Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 08/2022 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 52 002 400 E

Management System
certified according to
ISO 9001 / ISO 14001

**UK
CA**₂₅₀₃

CE₁₂₅₈

EAC

www.mt.com/pro