Type 405-DPAS-SC-VP pre-pressurized gel-reference pH Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	4
de	Betriebsanleitung	6
en	Instruction manual	8
es	Manual de instrucciones	10
fi	Käyttöohje	12
fr	Instructions d'utilisation	14
hu	Használati utasítás	16
it	Istruzioni per l'uso	18
ja	取扱説明書	20
ko	사용 매뉴얼	22
nl	Gebruikershandleiding	24
pl	Instrukcja obsługi	26
pt	Manual de instruções	28
ru	Инструкция по эксплуатации	30
sv	Bruksanvisning	32
th	คู่มือคำแนะนำ	34
zh	说明书	36

405-DPAS-SC-VP pH Electrodes
53 500 197 C

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Instructions for Type 405-DPAS-SC-VP pre-pressurized gel-reference pH Electrodes

1. Introduction

The METTLER TOLEDO Type 405-DPAS-SC-VP combination pH electrode is a low-maintenance, pre-pressurized, autoclavable electrode with a gel-type reference and a VarioPin (VP) connector. The reference element incorporates a silver-ion trap to prevent silver sulfide contamination of the reference junction.

2. Preparation

After unpacking, examine the electrode for any mechanical damage.

METTLER TOLEDO electrodes are factory tested and are shipped ready-to-use. (To maintain factory pressurization, the Type 405-DPAS-SC-VP electrodes have a silicone seal covering the ceramic diaphragm; this must be removed with the enclosed blade prior to use.) Remove the hydration cap and remove any electrolyte crystals by rinsing the electrode with distilled water. Pat dry; do not rub as this may build up a static charge and lengthen the response time.

Eliminate any air bubbles around the internal element in the pH-sensitive membrane.

3. Operation

METTLER TOLEDO pH electrodes can be used with most pH/mV meters. Follow manufacturer's instructions for use.

A two-point calibration is recommended, using fresh buffers of known pH values which bracket the sample pH; one buffer should be close to pH 7 to set the zero point. For fast routine analysis, a standardization using one buffer at a value near the sample pH can be done.

4. Sterilization

Type 405-DPAS-SC-VP electrodes are steam-sterilizable in-situ and auto-clavable. Follow operating instructions for the autoclave. It is not necessary to use a protective cap on the connector when autoclaving.

Repeated sterilization may cause the electrolyte to discolor slightly due to the elevated temperatures. This does not affect the electrode performance.

5. Maintenance

When properly cared for, METTLER TOLEDO electrodes should provide accurate results over thousands of uses.

Storage: Electrodes should be rinsed thoroughly after use and stored in 3M KCl or pH 4 buffer solution with added KCl. DO NOT store the electrode dry or in distilled or tap water. If the electrode has been stored dry or in water, it must be reconditioned by soaking in 3M KCl for several hours.

Refilling: Type 405-DPAS-SC-VP gel-reference electrodes are not refillable.

Protein Contamination: Clean the electrode with pH Electrode Cleaner for Proteins following instructions on the bottle.

Reactivating pH membrane: Use Reactivating Solution following directions CAREFULLY.

Servicing a clogged or dry reference junction: Immerse the electrode in a beaker of 3M KCl. Heat the solution to approximately 70 °C (160 °F) and allow electrode to soak in this solution as it cools to room temperature.

6. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

7. Troubleshooting

No response (symptom)

- Poor meter connection
Check connectors at the meter and at the electrode
- Cracked glass
Replace electrode

Slow response (symptom)

- Dehydrated pH membrane
Rehydrate (see Maintenance)
- Clogged reference junction
Soak in warm KCl (see Maintenance)
- Protein contamination
Use pH Electrode Cleaner (see Maintenance)
- Depleted gel layer
Use Reactivating Solution (see Maintenance)

Noise (symptom)

- Poor meter connection
Check connector
- Clogged reference junction
Soak in warm KCl (see Maintenance)

Low slope (symptom)

- Dehydrated pH membrane
Rehydrate (see Maintenance)
- Contaminated membrane
Use pH electrode cleaner (see Maintenance)

For further assistance, call our Customer Service Center at process.hotline@mt.com.

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Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○ : 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○ : Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x : 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x : Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 04/2023 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 53 500 197 C

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro