Type 465 combination pH electrodes

General instructions

cs	Návod k obsluze	2
da	Brugsanvisning	4
de	Betriebsanleitung	6
en	Instruction manual	8
es	Manual de instrucciones	10
fi	Käyttöohje	12
fr	Instructions d'utilisation	14
hu	Használati utasítás	16
it	Istruzioni per l'uso	18
ja	取扱説明書	20
ko	사용 매뉴얼	22
nl	Gebruikershandleiding	24
pl	Instrukcja obsługi	26
pt	Manual de instruções	28
ru	Инструкция по эксплуатации	30
sv	Bruksanvisning	32
th	คู่มือคำแนะนำ	34
zh	说明书	36

Type 465 pH Electrodes
59 905 937 D

Všeobecné pokyny

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General instructions

For all type 465 combination pH electrodes

For all type 465 combination pH electrodes to be mounted into the following housings: types InFit™ 763, InFit 764, and InTrac™ 776.

We recommend that the sensor may only be in operation with the original accessories from METTLER TOLEDO. The sensor may be operated and serviced only by personnel familiar with the unit, and who have read and understood this instruction manual.

1. Installation

- a When unpacking, check electrode for possible mechanical damage to the pH sensitive glass, the diaphragm or the body. Electrodes subject to complaints must be returned to us or to our agent in the original box indicating the invoice number.
- b Take off watering cap and wash the electrode tip with distilled water. After rinsing with distilled water, dab the electrode with a tissue. Do not rub as an increase of the response time might result.
- c **Important:** Reference electrolyte filler tube: The T-shaped rubber stopper **as well as the elastic securing band have to be removed completely from the electrode**, before inserting the electrode into the housing. Remove any KCl solution residing in the filling port. Re-assemble the housing as described in the instruction for the housing.
- d Eliminate any air bubbles trapped behind the pH sensitive glass by smoothly swinging the unit in the vertical plane. The lead-off element of the glass electrode must dip into the inner buffer solution without air bubbles.
- e Combination pH electrodes stored dry for several days have to be soaked in the appropriate storage solution for several hours before use. It is recommended to store combination and reference electrodes in the reference electrolyte solution.
- f After a long and dry storage period, the pH sensitive glass can be reactivated with hydrofluoric acid (order no. 51 319 053): The electrode membrane is dipped for 1 min. in this solution followed by rinsing with distilled water and soaking for 12 hours in the storage solution.
- g Checking electrolyte level: Combination electrodes of the type 465 contain various electrolytes. For refilling always use the type of electrolyte indicated on the body of the electrode. The imprint «refill 9816» means to refill with Viscolyt (order no. 51 340 235); «refill 9848» means to refill with Friscolyt-B (order no. 51 340 053).

The liquid level of a properly filled electrode is about 1 – 2 cm below the filling port.

2. Calibrating of electrode and pH transmitter

The transmitter is calibrated using the electrode either ready for mounting or already mounted into the housing. Calibration is made in two buffer solutions of a known pH value by dipping the electrode into the buffer and adjusting the pH transmitter. For further studies in pH measurement the following booklet is available free of charge: «pH Theory Guide» (order no. 30 078 149). For other language versions please contact your METTLER TOLEDO representative.

3. Maintenance

- a After each measuring cycle the electrode tip and the liquid junction must be thoroughly rinsed with distilled water. Under no circumstance must the electrode be allowed to dry with measuring solution coating it. **This is very important.**
- b On storage of an assembled probe, tip and liquid junction must be left fully immersed in an alkaline KCl solution (10 parts of 3M KCl or Viscolyt plus 2 parts of buffer solution pH 9.2). Without the addition of buffer solution pH 9.2 the housing may corrode.
- c The liquid level of the reference electrolyte must be checked regularly. If it has sunk below the bulbous shaft section of the electrode, the reference electrolyte must be replenished. **Attention is drawn to the instructions given in 1.g. above.**
- d **Rapid rise in resistance of reference electrode:** A rapid rise in resistance of the reference electrode (up to 1 M Ω) after 2 to 3 measuring cycles may occasionally be noticed in media which are highly enriched with proteins. Most transmitters are very sensitive to this rise in resistance, which manifests itself through a drift. If this happens, the electrode should be treated with our electrode cleaner (order no. 51 340 068).

A preventive action is to use combination electrodes with Friscolyt-B reference electrolyte (suffix/9848 in product designation).

- e **Black diaphragm** (liquid junction): If the solution to be measured contains sulfide compounds, the liquid junction may become black. Simultaneously, its resistance will rise significantly. Such an electrode does not work properly any longer and treating with our diaphragm cleaner (order no. 51 340 070) becomes necessary.

A preventive measure is the use of combination electrodes fitted with a silver-ion trap (code –SC– in product designation).

4. Cables

Make certain that the shield of the coaxial cable carrying **the reference electrode signal of the combined electrode is never connected to ground.**

Combination electrode with screw cap S7 suitable to cable socket AS9.

Combination electrode with screw cap K9 suitable to cable socket AK9.

5. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

InFit and InTrac are trademarks of the METTLER TOLEDO Group.

All other trademarks are the property of their respective holders.

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Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○ : 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○ : Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x : 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x : Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 04/2023 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 59 905 937 D

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro

* 5 9 9 0 5 9 3 7 0 *