METTLER TOLEDO

**Turbidity Sensors
Trübungssensoren
Sondes de turbidité**

InPro 86X0ie Series

English **page** **3**

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Turbidity Sensors

InPro 86X0ie

Instruction manual

Trademark notice:

InPro is a registered trademark of the METTLER TOLEDO Group in Switzerland, USA and a further eleven countries. All other trademarks are the property of their respective owners.

© It is forbidden to reprint this instruction manual in whole or part.

No part of this manual may be reproduced in any form, or modified, copied or distributed using electronic systems, in particular in the form of photocopies, photographs, magnetic or other recordings, without written consent of Mettler-Toledo GmbH, Process Analytics, CH - 8902 Urdorf.

All rights reserved, in particular reproduction, translation and patenting/registration.

Contents

1	Introduction	6
2	Safety	7
2.1	Introduction	7
2.2	Intended use	8
2.3	Inappropriate use	8
2.4	Basic principles	8
2.5	Responsibilities, organizational measures	9
2.5.1	Responsibilities of operator	9
2.5.2	Responsibilities of personnel	9
2.6	Product-specific hazards	9
2.6.1	Manipulation of and maintenance work on the sensor	9
2.6.2	Installation in pressurized systems	10
2.7	Residual hazards	10
2.7.1	Rupture of connections	10
2.7.2	Heat protection	10
2.7.3	External influences	10
2.8	Emergency measures	10
2.9	Safety measures	11
2.10	Modifications	11
3	Product description	12
3.1	Scope of delivery	12
3.2	Packing	12
3.3	Sensor designation	12
3.4	Product overview	13
3.5	Description of function	14
3.5.1	Measuring principle	14
3.5.2	Construction	14
4	Installation, operation and maintenance	15
4.1	Sensor installation	15
4.2	Electrical connections	17
4.3	Operation	17
4.3.1	Process conditions	17
4.4	Maintenance	18
4.4.1	Cleaning the sensor head	18
4.4.2	Calibration	18
4.4.3	Checking with an inspection unit	18
4.4.4	Tuchenhagen-VARIVENT® process adapter check and change of O-rings	18
4.4.5	Measurement window check	18
4.4.6	Humidity check	19
4.5	Storage	19
4.6	Disposal	19
5	Product specifications	20
6	Conditions of warranty	21
7	Ordering information	22

How to use this instruction manual

This instruction manual is an integral part of the METTLER TOLEDO InPro® 86X0ie sensors and contains notes and instructions that are important for safety and operation.

All persons working on or with the InPro 86X0ie sensor must have first read and understood the sections appropriate to the work in hand.

Please read this instruction manual carefully before using the InPro 86X0ie sensor. Keep this document close to the unit, so that operating personnel may easily be able to refer to it at any time.

Please first read section 1 “Introduction” and section 2 “Safety”.

Explanation of sensor designations

The general designation InPro 86X0ie used in this instruction manual refers to the following turbidity sensor models:

- InPro 8610ie: Single angle (25° scattered light)
- InPro 8630ie: Dual angle (25° and 90° scattered light) plus color

1 Introduction

Thank you for buying **the InPro 86X0ie Series sensor from METTLER TOLEDO**. The construction of the InPro 86X0ie sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and/or adverse effects on the plant or other equipment. **Therefore, the instruction manual must be read and understood by the persons involved before work is started with the sensor.**

In addition to this instruction manual please also note the following:

- all local safety regulations
- all instructions and warning remarks in the publications of the products that are used in conjunction with the sensor
- all safety precautions for the plant into which the sensor will be installed.

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

2 Safety

2.1 Introduction

This instruction manual contains all the information needed for safe and proper use of the InPro 86X0ie Series sensor. The instruction manual is intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

The instruction manual must be stored where it is constantly accessible and available to any person working with the InPro 86X0ie.

Warning notices and symbols

Warning notices and informations are marked with the following symbols:

Danger! Warning of a danger that can lead to extensive material damage, to death or grave bodily injury.

Caution! Warning of a possible dangerous situation that can lead to light bodily harm and/or material damage.

Notice: Information referring to technical requirements. Non-adherence can lead to malfunctions, uneconomic working and possibly also to loss of productivity. The manufacturer is not liable for any damages resulting from improper or inappropriate use.

2.2 Intended use

METTLER TOLEDO InPro 86X0ie Series sensors are intended solely for the measurement in liquids in industrial applications.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed
- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions
- observance of local legislation
- consultation with METTLER TOLEDO Process Analytics in the event of any uncertainties.

2.3 Inappropriate use

Any utilization other than the above mentioned, as well as any utilization with the technical data, is taken as being not in conformance with the regulations. The operator bears the sole risk for any damage caused by such utilization.

2.4 Basic principles

- The sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorized modifications. The risk is borne entirely by the user.

Notice: Manipulation of the sensors may only take place after it has been ensured that no process medium can escape through the sensor in the event of incorrect manipulation. For this reason, the complete system must be emptied and vented in advance (**safe condition**).

2.5 Responsibilities, organizational measures

2.5.1 Responsibilities of operator

- The operator is under obligation only to permit persons to work with InPro 86X0ie sensors, who are familiar with the basic requirements of work safety and accident prevention, and who have been instructed in the handling of the sensor. This instruction manual serves as the basic document.
- In addition to the instruction manual there are also generally valid legal and other binding regulations for work safety and accident prevention as well as for environmental protection, and these must be provided by the operator and instructed to personnel using the sensors.
- The operator/user must be fully aware of safety and potential danger factors during work with the sensor and this awareness must be checked by the operator at regular intervals.
- Measures must be taken to ensure that the sensor is only operated in a safe and fully functional condition.

Caution! Before the sensor is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

2.5.2 Responsibilities of personnel

- All persons whose duty it is to operate the sensor are under obligation to read chapter 1 "Introduction" and chapter 2 "Safety" as well as the warning notices in this instruction manual.
- In addition to the instruction manual, generally valid legal and other binding regulations for work safety and accident prevention must be adhered to.
- Any method of working which is doubtful from a safety perspective and which exceeds the operation according to regulations must be omitted.

Notice: Before every start-up, the sensor must be checked for:

- damage to the connections, fastenings, etc.
- leakage
- perfect functioning
- authorization for use in conjunction with other plant equipment and resources.

Caution! A defective sensor may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

2.6 Product-specific hazards

2.6.1 Manipulation of and maintenance work on the sensor

Notice: Before dismantling the sensor or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (depressurize, no explosion risk, empty, rinse, vent, etc.).

It is principally necessary to wear personal protective outfit such as protective goggles and clothing.

Only such maintenance work as is specifically described in this instruction manual may be performed on the sensor.

Only original spare parts from METTLER TOLEDO may be used for replacing defective components (see chapter 7).

Danger! Non-compliance with the prescribed maintenance instructions can endanger personnel and the environment.

2.6.2 Installation in pressurized systems

Notice: The maximum permissible temperature and pressure specifications must not be exceeded. The relative specifications are given on the type plate of the sensor.

2.7 Residual hazards

2.7.1 Rupture of connections

Sensor connections can become loosened through the effects of vibration.

Notice: All sensor interfaces including process adapter, and the interfaces between optical window and sensor head, need to be checked periodically by the customer and kept functional.

Danger! Leaky connections can allow process medium, cleaning solution or control (compressed) air to escape to the environment, presenting a hazard for persons and the environment.

2.7.2 Heat protection

Caution! The sensor is not equipped with heat protection. During steam-sterilization procedures in pipes, the surface of the sensor can reach high temperatures and cause burns.

2.7.3 External influences

Notice: Objects falling on the sensor can damage or destroy the unit, or cause leaks etc. Lateral forces may damage or destroy the unit.

2.8 Emergency measures

Notice: Always observe and comply with local regulations!

2.9 Safety measures

Notice: Always observe and comply with local laws and regulations! These are not an integral part of this instruction manual.

The operator is responsible for the instruction of personnel. Additional copies of this instruction manual can be ordered from the equipment supplier. As an integral part of the sensor, this instruction manual must at all times be readily accessible to users at the point of operation of the sensor.

The operator is obliged to inform the supplier / manufacturer of the sensor immediately about any safety-relevant incidents, or observations made, during use of the sensor.

Danger! Incorrect manipulation and/or instruction errors can lead to potential hazards for persons and for the environment.

Before every start-up the sensor must be checked for:

- damage to the connections, fastenings, etc.
- leakage
- defective cables and lines etc.
- authorization for use of the sensor in conjunction with the associated plant resources.

Caution! A defective sensor must never be installed or put into operation. Poor containment, leaky connections etc. or non-compliant installation of the sensor can lead to escape of process medium and thereby to a potential threat to life.

2.10 Modifications

Notice: No attachments or modifications to the sensor housings are allowed.

Danger! The manufacturer / supplier accepts no responsibility for any damage caused by unauthorized attachments and alterations or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

3 Product description

3.1 Scope of delivery

Standard supply of the sensor is made up of the following:

- instruction manual
- certificates in accordance with the specifications
- blanking plate

On receipt of the shipment, check immediately:

- The sensor and accessories for any sign of transport damage. Report any damage immediately to the carrier and to your supplier.
- The type designation on the sensor body.
- For completeness of the supply. Please notify your supplier immediately if the shipment is incomplete or in any way incorrect.

3.2 Packing

The packing consists of cardboard and foam plastic. Keep the packing for possible later use, such as for storage or transport of the sensor. However, if you wish to dispose of the packing, please observe your local regulations in this respect.

3.3 Sensor designation

Type designation


Sensor designation, product part number, serial number, material specifications as well as maximum permissible pressure and temperature can be noted from the product type plate.

Always include the product part number and serial number as minimum for clear identification when communicating with the manufacturer.

InPro 86X0ie

1 = single-angle version (25°)


3 = dual-angle plus color version (25° and 90° plus color)


Example sensor label (UKCA will be available from the beginning of 2022).

3.4 Product overview

The InPro 86X0ie series sensors can only be used in combination with the METTLER TOLEDO transmitter M800 1-Channel Process.

Dimension drawing InPro 86X0ie series

3.5 Description of function

3.5.1 Measuring principle

The InPro 86X0ie series of sensors are used for the determination of undissolved particles or substances in solutions. They are designed for measuring low to medium levels of turbidity and coloration.


The InPro 86X0ie sensors operate on the scattered-light principle in order to determine turbidity. Undissolved particles or substances located in the beam path of the sensor scatter the impinging light. Detectors arranged at an angle of 25° (single-angle version) or 25° and 90° (dual-angle version) measure the scattered light. The ratio of scattered light to transmitted light is output as a turbidity value. The turbidity value calculated in this manner is proportional to the concentration of undissolved particles or substances.

Additionally, the InPro 8630ie has a color measurement function (with or without turbidity compensation). Absorption of blue light is measured for the color measurement. The graduation of yellow is determined according to the Mitteleuropäische Brautechnische Analysenkommission (MEBAK) directives. This allows to measure turbidity and graduation of yellow with one instrument.

3.5.2 Construction

The sensor consists of a sensor head and a sensor body. The sensor head, which is in contact with the medium, is attached to the pipe using a Tuchenhausen-VARIVENT® process adapter.

The instrument electronics are contained in the sensor body.

Assembly of the InPro 86X0ie sensor.

- 1 Sensor head
- 2 Sensor body
- 3 5-pin connector
- 4 Cover

4 Installation, operation and maintenance

Caution! In all installation work described below, ensure that the equipment in which the sensor is to be installed is in a non-hazardous condition (depressurized, empty, rinsed, vented, etc.).


Requirements to be compliant with EHEDG: Install the sensor according to the requirements given in EHEDG Guidelines 8, 10 and 37. That is, to mount the sensor in a self-draining orientation, which is only possible in vertical pipe installations. The sensor should be installed flush to the process area. If mounting in a T-piece, the ratio between the depth of the upstand (L) and the diameter (D) of the upstand shall be $L/D < 1$. The device has been developed for Cleaning in Place (CIP) applications and must not be dismantled for cleaning.


4.1 Sensor installation

The sensor is equipped with a Tuchenhagen-VARIVENT® process adapter, which is used to attach the sensor to the process pipe in dimensions from DN40 to DN 150. Prior to the installation, the sealing surfaces (sensor and pipe adapter) need to be cleaned and checked for damages. Further, it needs to be verified that no obstacles are present within the installation cavity of the sensor.

The installation of the sensor should be in a **vertical standpipe**. If the sensor will be installed in a horizontal pipe, it must be in a horizontal position and the process pipe has to be filled completely with liquid during the measurement.

Installation on a vertical line (preferred)


Installation on a horizontal line

There is a mark on the sensor head which shows the **flow direction of the media**. The sensor needs to be installed in this direction.

Mark for the installation into a pipe

The sensor must be fitted to the pipe **at least 2 meters (6.6 ft) away from stray light sources** (such as sight glasses).

- 1 Cover
- 2 O-ring 82.00 x 3.00 (Nitril)
- 3 Desiccant bag compartment
- 4 Sensor body
- 5 O-ring 60.00 x 3.00 (EPDM)
- 6 Process Line

- 7 O-ring 60.00 x 3.00 (EPDM)
- 8 Sight-glass
- 9 Stainless steel blanking plate holder
- 10 Blanking plate
- 11 O-ring 34.60 x 2.62 (FKM)
- 12 Socket (for fixing the blanking plate to the stainless steel blanking plate holder).

Notice: All of the illustrated items must be fitted when the sensor is installed to ensure proper operation and splash protection (IP66).

4.2 Electrical connections

Caution! Some components inside the sensor are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero-potential before any operation on the wiring terminal. Be careful during preparing of operation, handling and operation of the sensor. This work is to be accomplished in principle only by qualified technical personnel..

Cable connection to M800 Process 1-Channel transmitter

The RS 485 data cable (5-pin) is available from METTLER TOLEDO in various length. Connect the data cable to the M800 Process 1-Channel transmitter as described below.

M800 Process 1-Channel	TB3	
Terminal	Function	5-pin cable wire color
1	1-Wire	–
2	GND 5V	Yellow
3	RS485 B	Blue
4	RS485 A	White
5	GND 5V	Grey
6	5V	–
7	24 V	Brown
8	GND 24 V	Black

In order to ensure optimum measurement results the following points must be considered:

- Run the connection cable through the cable gland.
- Connect the wires to the sensor.
- Hand-tighten the cable gland.

4.3 Operation

In case of any failure or malfunction, the equipment in which the sensor is installed, needs to be in a safe condition (depressurized, empty, rinse, vent, etc.).

4.3.1 Process conditions

- The plant pressure should never exceed the specification of the delivered sensor.
- The process temperature should never exceed the specification of the delivered sensor.

Caution! Exceeding the specified maximum pressure and/or the specified maximum temperature represent a very high safety risk. Please read the additional safety instructions before installation.

- Avoid air and gas bubbles inside the sensor, they cause disturbances. Noise and drift of the measurement would be the result (air bubbles are not expected at a pressures higher than 2 bar in aqueous solutions).

4.4 Maintenance

4.4.1 Cleaning the sensor head

The ratio measurement principle provides a large degree of compensation when the sensor head becomes contaminated with dirt/particles.

Depending on the operating conditions and media, the contamination can become so severe that the compensation provided by the ratio measurement principle is no longer adequate and the sensor head must be cleaned.

Window Fouling is a measurement of soiling of the measurement windows. This value is meaningful in clean water-related applications. A Window Fouling surveillance feature is implemented in the sensor to provide real-time information on sensor window contamination. A fouling alarm can be set in the range of 0...100% in the M800 Process 1-Channel transmitter relay, and to provide advance information for maintenance planning.

1. Place the system or equipment in a safe state (depressurized, empty, rinsed, vented or ventilated, etc.)
2. Disconnect electrical power
3. Dismount the sensor from the process
4. Clean the sensor head using a soft, lint-free cloth and a mild cleaning agent
5. Install the sensor

4.4.2 Calibration

The METTLER TOLEDO InPro 86X0ie turbidity sensor is factory calibrated prior to shipment. The factory calibration is done with 12-points formazin calibration over the entire measuring range. The factory calibration data is stored in the sensor and does not have to be entered manually. The LED light source in the InPro 86X0ie is very stable and there is no change in intensity over time. The ratio measuring principle compensates for sample color change, or window fouling. Hence, it is not necessary to re-calibrate the InPro 86X0ie turbidity sensor. When re-calibration with formazin is required, the InPro 86X0ie turbidity sensor should be returned to the factory. Please consult your local METTLER TOLEDO supplier for assistance.

In-situ one point process calibration (slope and offset) are possible for the InPro 86X0ie turbidity sensor. Please refer to the METTLER TOLEDO M800 1-channel transmitter instruction manual, section "Turbidity Calibration" for details.

4.4.3 Checking with an inspection unit

The sensor can be checked using an inspection unit with a built-in solid reference (turbi glass) every month or when necessary. The procedure is described in the corresponding instruction manual.

4.4.4 Tuchenhausen-VARIVENT® process adapter check and change of O-rings

The O-ring of the process adapter needs to be checked monthly and during every service call.

METTLER TOLEDO recommends that medium-wetted seals should under all circumstances be replaced every 6 months.

4.4.5 Measurement window check

The sapphire windows need to be checked monthly and during every service call. Take special attention to damage, contamination and sealing.

4.4.6 Humidity check

The sensor actively monitors the humidity level in its electronics compartment. The humidity level can be read out at the M800 transmitter. This humidity value should always be below 25% as otherwise there is a risk of condensation forming.

<Humidity in sensor body> is displayed in M800 transmitter in menu

<ISM --> ISM Diagnostics --> Sensor Monitor --> Humidity in sensor body.>

The moisture warning alarm is active when relative humidity in the sensor reaches 40%. A general alarm is active when relative humidity in the sensor reaches 50%. In such situations the desiccant must be exchanged immediately. Please contact your local METTLER TOLEDO representative for assistance.

METTLER TOLEDO recommends that the desiccant should, under all circumstances, be replaced every year.

Notice: It may take from several hours to a day for the humidity level to return to an optimal level. Removing moisture from the sensor's electronics compartment necessitates servicing. Please contact METTLER TOLEDO for assistance.

Desiccant will absorb moisture from atmosphere when it left open. Always store desiccant in a seal packaging/ container. Discard any unused desiccants when they are older than six months.

Do not change the desiccant when the pipeline is cold. There is a risk of condensate forming once the sensor cover has been unscrewed. The desiccant should be changed when the pipeline is warmer than the ambient air or when the temperature in the pipeline and the ambient temperature are approximately the same, e.g. during production downtimes.

4.5 Storage

Please inspect the instrument immediately after receiving for eventual transport damages.

The sensor should be stored in a dry, clean and protected area until time of installation.

Storage for a reinstallation

If the sensor has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

Shipment of the sensors

Please clean the instrument carefully before shipment (e.g. for revision/repair). Please use a fixed packaging to protect the instrument against transport damages. In ideal case the original packaging should be used.

Transportation and storage temperature: -10 ... +70 °C

Relative humidity: 0 ... 90% rel. (non-condensing)

4.6 Disposal

The disposal is to the responsibility of the user. Sensors contain electronic components that require special disposal without any health hazard for humans, and with no risk of harm for the environment.

Notice: Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

5 Product specifications

Specifications

Measuring principle	Turbidity: Scattered light (25°/90°**) Color**: Absorption measurement
Light source	Turbidity: 650 nm, LED (red light) Color: 430 nm, LED (blue light)
Measuring range	InPro 8610ie 0 to 250 EBC 0 to 1000 FTU InPro 8630ie 0 to 1000 EBC 0 to 4000 FTU 0 to 50 EBC Color
Units	FTU, NTU, EBC, EBCC, g/l, mg/l, ppm, ASBC, %T, ICUMSA, %Foul
Resolution	0.001 EBC
Repeatability	Range: 0 to 1 FTU (0.01 FTU) Range: 1 to 4000FTU (± 1%) Color: ± 0.2EBC color
Respond time (T90)	< 2 sec.
Process connection	Tuchenhagen-VARIVENT® Type N50/40
Pipe diameter (process)	DN40 to DN150
Surface finish (steel wetted)	N6/R _a 32 (R _a ≤ 0.8 µm/32 µin)
Sensor head (wetted)	Sapphire, Hastelloy (2.4602/C22)
Sensor housing	stainless steel (1.4404/316L)
Wetted measurement window	sapphire
O-ring (flange wetted)	EPDM
Protection degree	IP66
Admissible process temperature	-10 °C to +120 °C (14 to 248 °F) max. +120 °C/248 °F up to 2 hours max. +150 °C/302 °F short time, peaks for 15 mins
Admissible process pressure	up to 16 bar (232 psi)
Admissible ambient temperature	0 to 60 °C (32 to 140 °F)
Relative humidity	0 to 90% rel. (non-condensing)
Transportation and storage temperature	-10 to +70 °C (14 to 158 °F)
Communication	Digital (RS485)
Power supply	24 VDC (± 15%), 1.5 W, supply by transmitter
Sensor power consumption	approx. 65 mA (max.)
Certificate and Approvals	METTLER TOLEDO Quality Certificate, CE, PED, EHEDG
** InPro 8630ie only	

6 Conditions of warranty

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for the period of one year from date of delivery ex works. If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, carriage paid, to your appropriate METTLER TOLEDO agency. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the option of the Customer Service department of METTLER TOLEDO. After expiry of the period of warranty, faulty sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

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All rights reserved. 30 421 372 D