Verification Kit

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Verification Kit

for InPro 86X0i e Series Turbidity Sensors

InPro is a registered trademark of the METTLER TOLEDO Group in Switzerland, the USA and a further eleven countries.

Inspection of the turbidity sensors InPro 86X0i e


Thank you for buying the verification kit for the InPro[®] 86X0i e series turbidity sensors from METTLER TOLEDO. The kit is for periodic verification of the sensor, and calibration if required. The verification kit includes an integrated turbidglass (solid reference). The turbidglass value is determined by the user for each turbidity InPro 86X0i e series sensor the kit is used with (the verification kit is supplied without a sensor). It is highly recommended to use the verification kit to determine the turbidglass value, prior to installing a new sensor.

Procedure

Notice: The verification kit has to be matched with each InPro 86X0i e series sensor before it can be used for verification. It is necessary to record the matching value(s) for future reference. Each subsequent verification of a sensor is based on these recorded values.

It is possible that matching values may differ slightly between two new sensors. This does not mean a sensor or the verification kit is at fault. Please consult your local METTLER TOLEDO representative for further information.

Turbidity verification with turbidglass

1. Remove the sensor from the process (observe the instructions in the installation, operation and maintenance chapters in the sensor instruction manual).
2. Clean the sensor head with a soft cloth, where appropriate with a mild detergent, and finally with ultra-pure water or distilled water.
3. Place the sensor on a firm table with the sensor head upwards, and connect the sensor to the transmitter and wait for 5 minutes (in order for the electronics to warm up).
4. Clean the inside of the verification kit with a soft cloth, where appropriate with a mild detergent and finally with ultra-pure or distilled water. Clean the inside of the filler hose and discharge tube of the verification kit with a mild detergent and finally with ultra-pure or distilled water.
5. Carefully attach the verification kit to the sensor head using the process O-ring and Varivent™ clamp. Align the verification kit such that the hollow on the sensor matches with the nipple on the verification kit.
6. Fill the filler hose with ultra-pure water (with zero-FTU turbidity, bubble free). Slowly and carefully, fill it until the discharge tube is approximately 3/5 filled.
7. Expel bubbles from the unit by tilting the sensor forward and backward about 45°. Slowly tilt the sensor from side to side about 10 times.

Expel bubbles from the unit by tilting the sensor forward and backward about 45°

8. Place the sensor vertically on the table and check the measured value on the transmitter. This value should be recorded on the attached sheet.


Important: This matched measured value enables you to verify the sensor at any time in the future. Frequency of sensor verification is dependent on your plant standard operating procedures (SOPs). However, sensor verification is recommended every half-year to a year, or whenever necessary.

9. Drain the water by tilting the sensor and verification kit and empty the water through the discharge tube.
10. Remove the Varivent clamp and O-ring from the process connection.
11. Remove the verification kit carefully without damaging the turbiglass.
12. Carefully clean and dry the verification kit with a soft cloth. Keep it in its case to protect it from dirt and moisture. Store the case at room temperature.


Zero-point verification with ultra-pure water

Besides the turbidity verification with turbiglass, this verification kit can also verify the zero-point using ultra-pure water.


1. Place the sensor on a firm table with the sensor head upwards, and connect the sensor to the transmitter and wait for 5 minutes (in order for the electronics to warm up).
2. Clean the groove and wetted parts of the sensor head with ultra-pure water, rinse 3 times.
3. Completely fill the groove of the sensor head with ultra-pure water. Ensure there are no bubbles in the groove.

- Carefully attach the adapter ring to the sensor head. Align the sensor head such that the hollow on the sensor matches with the nipple on the adapter ring.

- Put the verification kit (with turbiglass) over the adapter ring. Align them such that the hollow on the adapter ring matches with the nipple on the verification kit. The purpose of the adapter ring and verification kit is to prevent light interfering with the verification.

6. Check the measured value on the transmitter. This value (zero-point) should be recorded on the attached sheet.
7. Carefully remove the verification kit and the adapter ring. Clean and dry the adapter ring. Keep the adapter ring and verification kit in their case to protect them from dirt and moisture. Store the case at room temperature.

Notice: Verification/Calibration should be done according to the appropriate instruction manuals.

For further information, please contact your local METTLER TOLEDO representative.

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© 10/2021 METTLER TOLEDO. All rights reserved.
Mettler-Toledo GmbH, Process Analytics