InPro 8000 Series Turbidity Sensors

Quick Setup Guide

Contents

cs	Průvodce rychlým nastavením	2
da	Lynvejledning	5
de	Kurzanleitung	8
en	Quick Setup Guide	11
es	Guía de configuración rápida	14
fi	Pika-asetusopas	17
fr	Guide de paramétrage rapide	20
hu	Gyorsbeállítási útmutató	23
it	Guida alla configurazione rapida	26
ja	クイック セットアップ ガイド	29
ko	사용 매뉴얼	32
nl	Beknopte handleiding	35
pl	Podręcznik szybkiej konfiguracji	38
pt	Guia de Configuração Rápida	41
ru	Руководство по быстрой настройке	44
sv	Snabbguide	47
th	คู่มือการตั้งค่าอย่างรวดเร็ว	50
zh	快速设置指南	53

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InPro 8000 Series

1 Introduction

Thank you for buying this InPro® 8000 Series turbidity sensor from METTLER TOLEDO. InPro 8000 Series are optical sensors which measure the concentration of suspended particles in a liquid medium on the basis of the principle of backscattered light. The sensor is operated with the M800 transmitter from METTLER TOLEDO.

2 Safety instructions

 **Notice:** Before every start-up, the sensor must be checked for:

- Damage to the connections, fastenings, etc.
- Leakage
- Perfect functioning
- Authorization for use in conjunction with other plant equipment and resources.

 **Notice:** The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

 **Caution:** Before the sensor is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

 **Caution:** A defective sensor may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

 **Caution:** Before installation in Ex-zones, please read and observe the guidelines (from instruction manuals) and adhere to side requirements.

 **Caution:** The fiber optic cables have a core of hard clad silica (HCS) fibers. When laying fiber optic cables, the bend radius must always exceed at least 10 cm (4"). Fiber optic cables may never be shortened independently. Fitting of the SMA plug requires special tooling. The transparent plastic cap on the sensor is to be removed before use. Any formation of deposits of solid particles on the surface is to be avoided.

 **Warning:** Before dismantling the sensor or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (depressurize, no explosion risk, empty, rinse, vent, etc.). The metal body of the turbidity sensors as well as the safety weld-in sockets METTLER TOLEDO housing InFit® type 76*-*** and InTrac® 7**-*** are, if necessary, to be included into the recurring pressure test of the unit.

3 Product description

InPro 8100/8200 series are ATEX certified to be installed in Ex-classified area. Please observe the precautions and special conditions for safe use before installing the sensor in an Ex-classified area.

Ex-classification:

⊕ II 1/2G Ex op is/op pr IIA 95 °C Ga/Gb

⊕ II 1/2G Ex op is/op pr IIB 38 °C Ga/Gb

Examination certificate number: SEV 12 ATEX 0169 X

Notified body: Electrosuisse SEV, (No.: 1258)

Technical Data

Technical data	InPro 8200/S(H)	InPro 8100/S	InPro 8050
Meas. principle	dual fiber	single fiber	single fiber
Shaft material	DIN 1.4435 or 1.4404 (316L) or Hastelloy C276 (available only with Epoxy sealed window)	1.4435 or 1.4404 (316L)	PSU (Polysulfone)
Shaft diameter	12 mm	12 mm	12 mm
Shaft length	120, 205, 297 or 407 mm	120, 205, 297 or 407 mm	120 mm
Window material	Sapphire	without window	Sapphire spigot
O-ring sealing	Viton®-FDA	Viton-FDA	Viton-FDA
Process adaption	Pg 13.5	Pg 13.5	Pg 13.5
Plug	2 × SMA	2 × SMA	2 × SMA
Meas. range	5 FTU...4000 FTU 0...30 g/l	10 FTU...4000 FTU 0...250 g/l	10 FTU...4000 FTU 0...250 g/l
Pressure range	0...6 bar (0...87 psi) with Kalrez®-FDA sealed window 0...16 bar (0...232 psi) with Epoxy bonded window	0...6 bar (0...87 psi)	0...2 bar (0...29 psi)
Temperature range	-10...130 °C (14...266 °F) with Kalrez-FDA sealed window -30...130 °C (-22...266 °F) with Epoxy bonded window	-30...130 °C (-22...266 °F)	0...60 °C (32...140 °F)
Steam-sterilizable	130 °C (only with Kalrez-FDA sealed window)	130 °C	non-sterilizable

Subject to technical changes.

4 Sensor installation

The turbidity sensors are to be installed at a distance of at least 10 cm (4") from the walls of the vessel, agitator or from the pipe walls in order to prevent disturbances through reflections. If it is not feasible in pipes of small diameter then the sensor can be installed sloped, at an angle of 15...45° or, alternatively, in a knee bend. The pipe must always be completely full of liquid.

5 Operation

If the turbidity measuring system is calibrated off-line, it is recommended to use the METTLER TOLEDO calibration tool, CaliCap®, (ordering no. 52 800 210) for fitting on the sensor shaft. If a CaliCap is not used, a minimum distance of 10 cm (4") to the vessel walls has to be observed as well. During calibration, the sample should be stirred in order to prevent the sedimentation of solid particles. In case of any failure or

malfunction, the equipment in which the sensor is installed, needs to be in a safe condition (depressurized, empty, rinse, vent, etc.).

6 Maintenance: Cleaning the sensor

Dirty or contaminated sensors can deliver incorrect measurement values. If fouling is presumed, the sensor is to be removed and cleaned. Suitable cleaning solutions are mild detergents or strongly diluted acids (< 0.5% by weight) such as hydrochloric acid.

7 Disposal

Disposal is to be carried out by the user in line with valid local regulations. The user must either deliver the sensor to an authorized private or public collection undertaking, or dispose of it himself in accordance with the regulations for such items/materials.

8 Warranty

On manufacturing defects, 12 months after delivery.

9 Storage condition

Inspect the instrument immediately after receiving for any transport damage. The sensor should be stored in a dry, clean and protected area until the time of installation. If the sensor has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

10 EC Declaration of Conformity

The complete declaration is available in the declaration of conformity certificate.

InPro, InFit, InTrac and CaliCap are registered trademarks of Mettler-Toledo GmbH in Switzerland, the USA, the European Union and a further five countries.

Viton and Kalrez are registered trademarks of DuPont Performance Elastomers LLC.

All other trademarks are the property of their respective holders.

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© 07/2019 METTLER TOLEDO. All rights reserved.
Mettler-Toledo GmbH, Process Analytics