InPro 86X0ie Turbidity Sensor

Quick Setup Guide

cs	Návod k obsluze	2
da	Instruktionsvejledning	7
de	Bedienungsanleitung	12
en	Instruction manual	17
es	Manual de instrucciones	22
fi	Käyttöohje	27
fr	Manuel d'instruction	32
hu	Használati utasítás	37
it	Manuale d'istruzioni	42
ja	取扱説明書	47
ko	사용 매뉴얼	52
nl	Gebruikershandleiding	57
pl	Instrukcja obsługi	62
pt	Manual de instruções	67
ru	Инструкция по эксплуатации	72
sv	Bruksanvisning	77
th	คู่มือคำแนะนำ	82
zh	说明书	87

Senzory InPro 86X0ie

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Obsah

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Instruction Manual

Contents

1	Introduction	18
2	Safety instructions	18
3	Product description	18
4	Sensor Installation.....	19
5	Operation	20
6	Maintenance: Cleaning the sensor head.....	20
7	Disposal	20
8	Warranty	21
9	Storage conditions	21
10	EC Declaration of Conformity	21

Appendix

Technical Drawings	92
--------------------------	----

InPro is a registered trademark of the METTLER TOLEDO Group in Switzerland, USA and a further eleven countries.

All other trademarks are the property of their respective owners.

1 Introduction

Thank you for buying this InPro® 86X0ie sensor from METTLER TOLEDO. The InPro 86X0ie sensor is used for the determination of undissolved particles or substances in solutions. They are designed for measuring low to medium levels of turbidity and coloration.

2 Safety instructions

 **Notice:** Before every start-up, the sensor must be checked for:

- Damage to the connections, fastenings, etc.
- Leakage
- Perfect functioning
- Authorization for use in conjunction with other plant equipment and resources.

 **Notice:** The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

Caution: Before the sensor is put into operation, the operator must have already clarified that use of the sensor in conjunction with the other associated equipment and resources is fully authorized.

Caution: A defective sensor may neither be installed nor put into operation. Faulty containment or installation out of conformance with regulations and instructions can lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

Caution: The sensor is not equipped with heat protection. During steam-sterilization procedures in pipes, the surface of the sensor can reach high temperatures and cause burns.

Caution: Some components inside the sensor are energized with voltages which can cause lethal shocks in case of contact. The instrument needs to be switched to zero potential before any operation on the wiring terminal.

Warning: Before dismantling the sensor or commencing any maintenance work on it, ensure that the equipment in which the sensor is installed is in a safe condition (pressurize, no explosion risk, empty, rinse, vent, etc.).

3 Product description

The label on the housing contains the following information:

InPro 86X0ie

- x: 1 = single-angle version (25°)
3 = dual-angle plus colour version (25° and 90° plus color)


Serial number

The serial number can be found on the label on the sensor.

Sensor certifications

The label located on the sensor provides the information about the sensor certifications.

Example label

(UKCA will be available from the beginning of 2022)

Cable connection to M800 Process 1-Channel transmitter

The RS-485 data cable (5-pin) is available from METTLER TOLEDO in various lengths. Connect the data cable to the M800 Process 1-Channel transmitter as described below.

M800 Process 1-Channel	TB3	
Terminal	Function	5-pin cable wire color
1	1-Wire	–
2	GND5V	Yellow
3	RS485B	Blue
4	RS485A	White
5	GND5V	Grey
6	5V	–
7	24V	Brown
8	GND24V	Black
9–16	Not used	–


In order to ensure optimum measurement results the following points must be considered:

- Run the connection cable through the cable gland.
- Connect the wires to the sensor.
- Hand-tighten the cable gland.


4 Sensor Installation

The sensor is equipped with a Tuchenhausen-VARIVENT® process adapter, which is used to attach the sensor to the process pipe in dimensions from DN 40 to DN 150. Prior to the installation, the sealing surfaces (sensor and pipe adapter) need to be cleaned and checked for damages. Further, it needs to be verified that no obstacles are present within the installation cavity of the sensor.

The installation of the sensor should be in a **vertical stand-pipe**. If the sensor will be installed in a horizontal pipe, it must be in a horizontal position and the process pipe has to be filled completely with liquid during the measurement.

Installation on a vertical line (preferred)

Installation on a horizontal line

There is a mark on the sensor head which shows the flow direction of the media. The sensor needs to be installed in this direction.

Mark for the installation into a pipe

The sensor must be fitted to the pipe at least 2 meters (6.6 ft) away from stray light sources (such as sight glasses).

5 Operation

In case of any failure or malfunction, the equipment in which the sensor is installed, needs to be in a safe condition (depressurized, empty, rinse, vent, etc.).

6 Maintenance: Cleaning the sensor head

The measurement principle provides a large degree of compensation for contamination of the sensor head. However, the contamination can become so severe, depending on the operating conditions and media, that the compensation is no longer adequate and the sensor head must be cleaned.

1. Place the system or equipment in a safe state (depressurized, empty, rinsed, vented or ventilated, etc.)
2. Disconnect electrical power
3. Dismount the sensor from process
4. Clean the sensor head using a soft, lint-free cloth and a mild cleaning agent
5. Install the sensor

7 Disposal

It is the user's responsibility to dispose of the sensor in a professional manner. This sensor contains electronic components that require special disposal without any hazard for humans, and with no risk of harm to the environment. Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

8 Warranty

On manufacturing defects, 12 months after delivery.

9 Storage conditions

Inspect the instrument immediately after receiving for any transport damage. The sensor should be stored in a dry, clean and protected area until the time of installation. If the sensor has been removed from the process, it should be thoroughly cleaned and dried. It should be stored in a dry, clean and protected area until time of installation.

10 EC Declaration of Conformity

The complete declaration is available in the declaration of conformity certificate.

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All rights reserved. 30 421 373 B
Mettler-Toledo GmbH, Process Analytics