CO₂ ISM Qualification Kit

Instruction manual
Bedienungsanleitung
Instructions d'utilisation

CO₂ ISM Qualification Kit
30 019 554 C

METTLER TOLEDO

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CO₂ ISM Qualification Kit

Instruction manual

ISM is a trademark of the METTLER TOLEDO Group.

Contents

	Page
1 Safety	5
1.1 Intended use.....	6
1.2 Environmental protection.....	6
2 Product Description	7
3 Operation	8
3.1 Connecting the transmitter to the CO ₂ ISM Qualification Kit	8
4 Product Specification	10
4.1 Simulator, Pre-registered Values	10
4.2 Sensor Wear Data	10
4.3 Warnings and Alarms.....	11
5 Installation	12
5.1 Dimensions	12
6 Ordering Information	13
7 Protocol	14

1 Safety

This instruction manual lists the main points of how to use the **CO₂ ISM™ Qualification Kit** from METTLER TOLEDO for best results and in the intended manner.

This **instruction manual**, and the **safety requirements** in particular, **must be observed** by the **operator and all other persons** working with this equipment.

The instruction manual must always be stored close at hand, in a place accessible to all people working with an ISM CO₂ sensor and the CO₂ ISM Qualification Kit.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of product.

1.1 Intended use

The CO₂ ISM Qualification Kit is intended **solely for validating and testing of the ISM transmitter family**, respectively, as described in this instruction manual.

The CO₂ ISM Qualification Kit is intended **for simulating the functions of METTLER TOLEDO ISM sensors**, as described in this instruction manual.

Any other use, or use not mentioned here, that is incompatible with the technical specifications is deemed inappropriate. The operator is solely responsible for any damage arising from such use.

Other prerequisites for appropriate use include:

- observing the instructions, notes and requirements set out in this instruction manual
- correct care and maintenance of the unit, according to the instruction manual
- observing the prescribed environmental and operational conditions
- compliance with local statutes

1.2 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

2 Product Description

The CO₂ ISM Qualification Kit is a set of 5 different sensor plug-in heads and allows the simulation of reading values of CO₂ ISM sensors with predefined values and errors (not changeable by the user) to assure safe use as qualification tool.

Each sensor plug-in head corresponds to one of METTLER TOLEDO's CO₂ ISM sensors and delivers the complete data of an ISM sensor.

The CO₂ ISM Qualification Kit has basically 2 functions:

- Tool for the qualification of measuring systems.
- Service tool for quick checks and loop trouble-shooting.

The data indicated (primary / secondary) and transmitted are as follow:

- Main measured value
- Temperature
- Error message according to indication.

All standard process and wear data of the sensor are also available in each plug-in head and are also sent to the transmitter or the PLC as well as to the iSense software.

3 Operation

3.1 Connecting the transmitter to the CO₂ ISM Qualification Kit

To qualify your measuring systems proceed as follows:

Please refer to the transmitter instruction manual for the exact description of the transmitter operation.

Step 1: Disconnect the ISM sensor from the measuring point.

Step 2: Connect one of the CO₂ ISM Qualification Kit plug-in heads to the transmitter via the standard connecting cable.

Step 3: The preregistered values then are displayed on the transmitter.

Step 4: Check if the displayed values by the transmitter or the PLC correspond to the predefined fixed values of the CO₂ ISM Qualification Kit plug-in head.

See certification for exact values and tolerances.

Step 5: In case you are using the Error simulators, check that the appropriate errors are correctly displayed by the transmitter or the PLC.

Step 6: Report results into the attached spreadsheet "Protocol for CO₂ ISM Qualification Kit".

Problem solving: If any of the tests described above failed, you should change the connecting cable and repeat the test. If the test still shows the same result, the transmitter or the cable is defective. If you are using a transmitter with current outputs, so ensure the outputs are correctly wired.

You must **never do a calibration** (the simulators are read-only).

You must never increment autoclave cycles (always set "NO" to cycle counter), this is only needed if autoclave has been set to a value higher than 0.

Different transmitters may display different values and errors.


Out of range behavior is also different depending on the transmitter in use thus **different values/readings may occur.**

Set the transmitter input filter to off (only for SIM CO₂ TGL).

Do not use the CO₂ ISM plug-in heads after the expiry date shown on the label of each plug-in head. Please send them back to METTLER TOLEDO for recertification.

Measuring point with
METTLER TOLEDO electrode

Connecting CO₂ ISM
Qualifying plug-in head

Environment operating and storing conditions:

Permissible ambient temperature (operation):

– 25 °C to 60 °C (– 13 °F to 140 °F)

Permissible ambient temperature (storage/transport):

– 25 °C to 60 °C (– 13 °F to 140 °F)

4 Product Specification

4.1 Simulator, Pre-registered Values

Order No.	Product Designation	CO ₂ (mbar)	Temp.
SIM CO₂ Kit		Order No. 30 031 035	
30 034 395	SIM CO ₂ 15 mbar	15.0	25 °C
30 038 175	SIM CO ₂ 950 mbar	950.0	40 °C
30 038 973	SIM CO ₂ TGL	15.0 & 950.0	25 °C & 40 °C
30 038 974	SIM CO ₂ ERR1	204.3	37 °C
30 038 975	SIM CO ₂ ERR2	10.1	37 °C

4.2 Sensor Wear Data

Sensor Wear Monitor CO ₂	Transmitters M400 / M800
Time to maintenance	120 d
Life time indicator	360 d
Adaptive Cal Timer	7 d
Sensor operating time	1 d
Autoclaving cycles	1
CIP cycles	3
SIP cycles	2

4.3 Warnings and Alarms

30 038 974 (SIM CO₂ ERR1) **Incident:** Slope reduced

Transmitter	Warnings / Alarms
M400	Error pH slope < 80%
M800	Error pH slope < 80%

30 038 975 (SIM CO₂ ERR2) **Incident:** Slope reduced


Transmitter	Warnings / Alarms
M400	Error pH ZeroPt > 8.0 pH
M800	Error pH Offset > 8.00 pH

Please refer to your transmitter instruction manual for specific instructions.

5 Installation

5.1 Dimensions

Markings

The markings labeled on each CO₂ ISM plug-in head state the following information:

6 Ordering Information

Order No.	Designation	Order No.	Designation
52 300 410	SIM pH Kit	52 300 428	SIM O₂ Trace Kit
52 300 411	SIM pH4	52 300 429	SIM O ₂ Trace ZERO
52 300 412	SIM pH7	52 300 430	SIM O ₂ Trace AIR
52 300 413	SIM pH TGL	52 300 431	SIM O ₂ Trace TGL
52 300 414	SIM pH ERR1	52 300 432	SIM O ₂ Trace ERR1
52 300 415	SIM pH ERR2	52 300 433	SIM O ₂ Trace ERR2
52 300 416	SIM O₂ Kit	52 300 422	SIM O₂ ppb Kit
52 300 417	SIM O ₂ ZERO	52 300 423	SIM O ₂ ppb ZERO
52 300 418	SIM O ₂ AIR	52 300 424	SIM O ₂ ppb AIR
52 300 419	SIM O ₂ TGL	52 300 425	SIM O ₂ ppb TGL
52 300 420	SIM O ₂ ERR1	52 300 426	SIM O ₂ ppb ERR1
52 300 421	SIM O ₂ ERR2	52 300 427	SIM O ₂ ppb ERR2
30 031 035	SIM CO₂ Kit		
30 034 395	SIM CO ₂ 15 mbar		
30 038 175	SIM CO ₂ 950 mbar		
30 038 973	SIM CO ₂ TGL		
30 038 974	SIM CO ₂ ERR1		
30 038 975	SIM CO ₂ ERR2		

The single simulators are not available for sale.

7 Protocol

Protocol for CO₂ ISM Qualification Kit

Facility Name _____

Instrument Manufacturer _____

Date of Purchase _____

Transmitter Model _____

Transmitter Serial Nr _____

Simulator P/N	Simulator S/N	Designation	Valid until	Parameter/s (check all that apply for this instrument)

Date / Time	Acceptable? ⇔ Yes	Acceptable? ⇔ No	Name / Signature

Protocol for CO₂ ISM Qualification Kit © Mettler-Toledo GmbH, MaCom, Im Hackacker 15, CH - 8902 Urdorf. 08/2021

A copy of this form in A4 format is delivered with the CO₂ ISM Qualification Kit.

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All rights reserved
Protocole pour set de qualification ISM CO₂ © Mettler-Toledo GmbH, MarCom, Im Hackacker 15, CH - 8902 Urdorf. 08/2021