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## Operating instructions for pH electrodes InPro® 325X, InPro® 325X (ISM), InPro® 325X i

### 1. Introduction

The METTLER TOLEDO InPro® 325X pH single rod measuring cell is a low maintenance, autoclavable and pressure resistant pH electrode with a pressurised liquid electrolyte and a reference electrolyte together with a silver ions barrier. Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use. We recommend that you only operate the electrode in conjunction with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

### 2. Safety instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis both protective spectacles and protective gloves should be worn. Do not touch the electrode at the plug-in contacts, since electrostatic discharges can damage the electronics.

### 3. Product description

The inscription on each electrode contains the following information:

<b>METTLER TOLEDO InPro® 325X i /SG /yyy /zzz /ISM</b>	Manufacturer of the electrode
	Type designation X = membrane glass, i = ISM Digital /SG = auxillary electrode / yyy = rod length in mm /zzz = type of temperature sensor / ISM (Analog)
combination <b>pH</b>	Type of electrode
pH 0 ... 14 (InPro 3250), pH 0 ... 12 (InPro 3253)	pH-measurement range
pH 1 ... 11 (InPro 3251, InPro 3252)	
0 ... 100 °C (up to 140 °C InPro 3250, InPro 3253)	Temperature range for operation (for sterilisation)
–25 ... 80 °C (InPro 3251), 0 ... 80 °C (InPro 3252)	
Order No. 30 04X XXX	Order number
Observe the operating instructions	Instruction (see: operating instructions)

In addition each electrode is provided with a serial number on the red part of the plug-in head in order to enable identification.

This product complies with the ATEX guideline 94/9/EC. The complete declaration can be downloaded under www.mt.com/pro/.

#### 4. Installation and comissioning

- While unpacking check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
- Remove the protective cap and rinse the electrode for a short time with deionised water. Carefully remove the silicon protection from the diaphragm with the aid of the knife supplied (cf. green instruction images on the packaging). If the electrode has already been used the sensitive part can be cleaned if necessary with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionised water. After rinsing the electrode should only be dabbed. Any rubbing of the pH-sensitive glass can increase the response time as a result of electrostatic charge accumulation.
- Check the space behind the pH sensitive glass for air bubbles and remove any air bubbles that may be present by shaking the electrode lightly in a vertical direction.
- Place the electrode into the armature as described in the instructions for the armature.
- Connect the electrode and the pH transmitter with the corresponding cable (VP6, VP8, AK9). Observe the connections scheme that accompanies the cable or the transmitter.
- Intelligent ISM sensors such as InPro 325X (ISM) and InPro 325X i enable «plug and measure» and extended diagnostics. Operation of the electrodes in the ISM mode requires the M700 transmitter as well as an ISM pH-module (2700 i X). For installation, commissioning and operation of the ISM system please observe the instructions for the transmitter, module and cable.

#### 5. Operation: Calibration of the electrode and pH transmitter

Before a calibration immerse the electrode for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connecting it to the transmitter. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. For further details please consult the operating instructions for the pH transmitter.

#### 6. Maintenance

- After each working cycle the electrode tip and diaphragm(s) are to be carefully rinsed with deionised water. Any drying out of residues from the solution being measured on these parts must be avoided at all costs!
- If the electrode is not in use it is stored with the electrode tip and the diaphragm immersed in a tank of electrolyte (9823/3M KCl or 9848 Friscolyt).
- If the electrode is stored in its armature the rules described under Item 2 apply, although the electrolyte has to be slightly modified for storage purposes (2 parts of buffer solution pH 9.2 to 10 parts of electrolyte) in order to avoid any corrosion of the armature parts.
- If an electrode is inadvertently stored in a dry condition for several days it must be rinsed for several hours in a standard storage electrolyte.
- The plug-in head should occasionally be checked for possible traces of moisture. If necessary clean it thoroughly with deionised water or alcohol and then dry it off carefully.

#### 7. Removal of sources of malfunction

With a soft toothbrush and water any residues can easily be removed from the electrode. Mild washing agents can also be used. Protein contaminations on the diaphragm can be most likely removed with the cleaning solution (order no. 51 340 070). Warning: Solutions containing hydrofluoric acid damage the pH sensitive part of the electrode.

#### 8. Disposal

It is the user’s responsibility to dispose of the sensor in a professional manner. The sensor may contain electronic components that require correct disposal such that no persons or the environment are endangered.

#### 9. Warranty

On manufacturing defects 12 months after delivery.

#### 10. Storage conditions

We recommend that the electrodes are stored at room temperature.

## Instructions générales de pH-électrode InPro® 325X, InPro® 325X (ISM), InPro® 325X i

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Subject to technical changes. 05./2012. © Mettler-Toledo AG. Printed in Switzerland.