InPro 4501 pH Combination Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	5
de	Betriebsanleitung	8
en	Instruction Manual	11
es	Manual de instrucciones	14
fr	Instructions d'utilisation	17
fi	Käyttöohje	20
hu	Használati utasítás	23
it	Istruzioni per l'uso	26
ja	取扱説明書	29
ko	사용 매뉴얼	32
nl	Gebruikershandleiding	35
pl	Instrukcja obsługi	38
pt	Manual de instruções	41
ru	Инструкция по эксплуатации	44
sv	Bruksanvisning	47
th	คู่มือคำแนะนำ	50
zh	说明书	53

InPro 4501 pH Combination Electrode
52 002 390 C

InPro 4501

Kombinované pH elektrody

Všeobecné pokyny

1. Úvod

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General instructions

1. Introduction

These general instructions must be readily available to operators of the electrode at all times. Before using the electrode, these instructions must be read and understood by all operating personnel.

The METTLER TOLEDO electrode versions InPro™ 4501 VP and InPro 4501 fixed cable are low-maintenance, temperature-compensated pH combination electrodes containing polymer electrolyte (Xerolyt™ Extra), and are hereafter referred to as the InPro 4501 family. An integrated temperature sensor (RTD) provides automatic compensation of the electrode slope to ensure accurate readings.

The InPro 4501 family is suitable for pH measurements in:

- emulsions and suspensions
- sulfide-bearing media
- solutions containing heavy metals (lead, mercury, silver, etc.)
- processes with rapid pressure fluctuations (within the given range)
- measurements in wastewater and in decontamination processes for industrial effluents
- heavily contaminated media

ATTENTION: The InPro 4501 family electrodes are unsuitable for sterilization/autoclaving procedures, for processes with rapid temperature fluctuations or food and beverage applications.

2. Safety remarks

The InPro 4501 family of electrodes is intended solely for the purpose described in §1 above. The electrodes consist partly of glass. If they are to be cleaned or calibrated using acid or alkaline solutions, protective clothing (goggles and gloves) must be worn. In addition, local and in-house-safety regulations must be observed.

We recommend that the electrodes only be put into operation in conjunction with original products/accessories from METTLER TOLEDO. The electrodes may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions.

3. Product description

The markings on each electrode provide the following information:

InPro 4501 VP/xxx/yyy	Designation/temperature sensor/solution ground (SG)
InPro 4501/xxx/yyy/zzz/www	Fixed cable designation/temperature sensor/solution ground (SG)/cable length/cable end
Combination electrode	type of sensor
Pt100/Pt1000	Type of built-in temperature sensor
pH 1...14/0...100 °C	pH measuring range/temperature range
Order no. 41 453 310x	Ordering number VP
Order no. 10 453 30xx	Ordering number fixed cable

4. Installation and preparation for use

1. On unpacking, and again prior to installation, check the electrode visually for signs of mechanical damage to the glass membrane, the shaft and the connector. Report any faults or defects immediately to your authorized METTLER TOLEDO dealer.
2. Remove the watering cap and rinse the electrode with distilled or de-ionized water. Use only soft tissue to dab the electrode dry. Do not rub the glass membrane since this can lead to sluggish response behavior of the electrode.
3. Electrodes with a dry glass membrane must be conditioned for several hours in a 3 mol KCl solution before being used.
4. Connect always the solution ground (blue cable) to the transmitter.
5. Use PTFE tape to seal the NPT thread.
6. If possible use the protection sleeve (order no. 52 401 808).

5. Cabling

The InPro 4501 VP is equipped with a VarioPin connector. To connect the cable to the electrode, attach and gently rotate until the connector is engaged (groove). Screw the cable cap nut onto the electrode until the connection is secure. The plug assembly should be fitted hand-tight only; no tools should be used. For details of the exact assignment of the individual cable leads to the VP pins, please refer to the instruction leaflet "Preparation of the cable for connection to the transmitter" supplied together with the cable.

InPro 4501 fixed cable. There are two options for termination for the InPro 4501 with fixed cable:

- BNC connector for pH combination electrode with stripped and tinned leads for RTD and solution ground.
- Stripped and tinned leads for pH combination electrode, for RTD and solution ground.

Refer to your transmitter instruction manual for correct wiring.

ATTENTION: The solution ground lead is designed to go to a galvanically isolated terminal on the transmitter and should not be used as an outer ground shield. If the transmitter does not have sensor diagnostics, this lead should be left unterminated.

Color	Function
black/transparent	Glass
red	Coax shield/reference
grey	RTD sensing (temperature sensor, 3-wire connection)
blue	Solution ground
white	RTD return (temperature sensor, 2-wire connection)
green	RTD in (temperature sensor, 2-wire connection)
green/yellow	Cable shield

6. Calibration and measurement

InPro 4501 VP electrodes require two-point calibration. Before calibration, the watering cap must be removed from the pH electrode. The electrode is then dipped successively into two different buffer solutions, each of known pH value, and the transmitter calibrated to these buffer values (for example, buffers pH 4.01 and pH 7.00). Please note the instructions given in the operating manual of the transmitter. Calibration should be repeated as and when needed.

7. Maintenance

- Electrodes with Xerolyt Extra polymer electrolyte cannot be refilled.
- Electrodes should not be stored dry for any period of time
- During storage, the watering cap should be in place, filled with Friscolyt™ or 3 mol KCl solution.
- Never store an electrode in distilled or deionized water, as this shortens the lifetime of the electrode.
- Should the electrode dry out, it must be conditioned for 24 hours in Friscolyt (ordering no. 51 2340 053 for 250 mL, or 51 340 054 for 6 x 250 mL) before reuse. After conditioning, recalibration is required.
- After each workcycle, the tip of the electrode should be cleaned using de-ionized water.

ATTENTION: Do not under any circumstance use METTLER TOLEDO reactivating solution (ordering no. 51 319 053/in USA 43 989 5125) as this contains hydrofluoric acid (HF) which would attack the Xerolyt reference system.

8. Diagnostics

The electrodes are equipped with a grounding ring which allows a diagnostic of the diaphragm condition. If this diaphragm (open junction) should become clogged or if the reference resistance exceeds a certain value (see transmitter operation manual), it is then necessary to replace the electrode. Electrodes of the InPro 4501 family should as a general rule not be cleaned with acidic cleaning solutions, e.g. diaphragm cleaning solution, reactivation solution, etc. However, by extremely careful application, it is possible to reactivate the glass membrane, see «point 9» below.

9. Troubleshooting (sluggish response time, insufficient sensitivity, etc.)

ATTENTION: The glass membrane must be very cautiously dabbed with HF solution using a fine brush and taking care that solution does not come into contact with the titanium pin (metallic solution ground). Also check the electrode connector (only VP version) for possible traces of moisture. If necessary, clean well using deionized water or alcohol, and afterwards dry carefully.

10. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

11. Guarantee

METTLER TOLEDO grants 12 months warranty covering manufacturing defects, valid from the date of production of the electrode.

InPro, Xerolyt and Friscolyt are trademarks of the METTLER TOLEDO Group.

All other trademarks are the property of their respective holders.

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本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 04/2023 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 52 002 390 C

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro