
InPro 6000 G Series Gas Phase Oxygen Sensors

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© It is forbidden to reprint this Instruction Manual in whole or part. No part of this manual may be reproduced in any form, or modified, copied or distributed using electronic systems, in particular in the form of photocopies, photographs, magnetic or other recordings, without written consent of Mettler-Toledo GmbH, Process Analytics, CH-8902 Urdorf, Switzerland.

All rights reserved, in particular reproduction, translation and patenting/registration.

InDip, InFit, InPro, ISM, and InTrac are trademarks of the METTLER TOLEDO Group. Kalrez is a registered trademarks of DuPont Performance Elastomers LLC.”

Contents

1	Introduction	4
2	Important notes	5
2.1	Notes on operating instructions	5
2.2	Intended use	6
2.3	Safety instructions	6
2.4	Environmental protection	9
2.5	Examples of some typical applications	9
2.6	Use in Ex-zones	9
2.7	Ex-classification ATEX	10
2.7.1	Introduction	10
2.7.2	Rated data	10
2.7.3	Special conditions	11
2.8	Ex-classification – FM approved	12
2.8.1	Introduction	12
2.8.2	Rated data	13
2.8.3	Special conditions	13
2.8.4	Applied standards:	14
3	Product description	15
3.1	General information	15
3.2	Principle	15
3.3	Enhanced diagnostics (only for ISM types available)	17
3.4	Scope of delivery	18
3.5	Equipment features	19
4	Installation	20
4.1	Mounting the sensor	20
4.2	Connection	20
4.2.1	Connecting the sensor to a AK9 cable	20
4.2.2	Connecting the cable to the transmitter	21
5	Operation	22
5.1	Start-up and polarizing	22
5.2	Calibration	23
5.2.1	Purpose of calibration	23
5.2.2	What you have to know for calibration	23
5.2.3	Single point calibration	24
5.2.4	Dual point calibration	24
6	Maintenance	25
6.1	Inspection of the sensor	25
6.1.1	Visual inspection	25
6.1.2	Using the METTLER TOLEDO O ₂ sensor master	26
6.1.3	Testing the sensor via a transmitter	27
6.1.4	ISM design	28
6.2	Changing the electrolyte, the membrane body or the interior body	29
7	Storage	32
8	Product specification	32
8.1	Certificates	32
8.2	Specifications	33
9	Ordering information	35
9.1	Sensors	35
9.2	Transmitter	35
9.3	Accessoires	35
9.4	Spare parts	36
9.5	Cables	36
9.6	Recommended housings	36
10	Theory of the polarographic sensor	37
10.1	Introduction	37
10.2	Principle of the design of an oxygen electrode	38
10.3	Parameters determining current	42
10.4	Polarization voltage	42
10.5	Temperature	43
10.6	Dependence on flow	43
10.7	Oxygen partial pressure – oxygen concentration	43

1 Introduction

Thank you for buying the **InPro™ 6800 G / 6850 i G / 6900 i G / 6950 i G sensor from METTLER TOLEDO.**

The construction of the InPro series employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and / or adverse effects on the plant or other equipment.

Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.

The instruction manual must always be stored close at hand, in a place accessible to all people working with the InPro sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. He will be glad to assist you.

2 Important notes

2.1 Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the InPro 6800 G / 6850 i G / 6900 i G / 6950 i G sensor.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and/or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

2.2 Intended use

METTLER TOLEDO InPro 6800 G / 6850 i G / 6900 i G / 6950 i G sensors are intended solely for in-line measurement of the oxygen partial pressure in liquids and gases, as described in this instruction manual.

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose.

The manufacturer/supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user/operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (housings, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.
- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

2.3 Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. housings, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant

operator. This applies in particular in the case of plants operating in hazardous zones.

- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, “9.3 Spare parts”).
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.
- When using cables for sensor connection not supplied and recommended by the manufacturer the following standards need to be considered:
 - US National Electrical Code® (ANSI/NFPA 70 [NEC®]), where applicable.
 - Canadian Electrical (CE) Code® (CEC Part 1, CAN/CSA-C22.1), where applicable.
- For guidance on US installations, see ANSI/ISA-RP12.06.01, Installation of Intrinsically Safe Systems for Hazardous (Classified) Locations.
- Control room equipment connected to intrinsically safe associated apparatus in the US should not use or generate more than the specified Um of the associated apparatus.
- Care must be taken during installation to avoid impacts or friction that could create an ignition source.
- Tampering and replacement with non-factory components may adversely affect the safe use of the system.

- Insertion or withdrawal of removable electrical connectors or modules is to be accomplished only when the area is known to be free of flammable vapors.
- **WARNING –**
INTRINSICALLY SAFE APPARATUS CAN BE A SOURCE OF IGNITION IF INTERNAL SPACINGS ARE SHORTED OR CONNECTIONS OPENED.
- **WARNING –**
DO NOT OPEN WHEN AN EXPLOSIVE ATMOSPHERE IS PRESENT.
- **WARNING –**
SUBSTITUTION OF COMPONENTS MAY IMPAIR INTRINSIC SAFETY.
- **WARNING –**
SUBSTITUTION OF COMPONENTS MAY IMPAIR SUITABILITY OF THE EQUIPMENT.
- **WARNING –**
FOR CONNECTION ONLY TO NON-FLAMMABLE PROCESSES.
- **WARNING –**
TO MAINTAIN THE ENCLOSURE IP66 RATING, THE CONNECTOR MUST BE FULLY ENGAGED.

2.4 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

2.5 Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

Measurement in gases:

- CO₂ recovery
- CO₂ purity
- Product storage
- Inert production
- Offgas monitoring

2.6 Use in Ex-zones

Attention!

For an installation in Ex-zones please read the guidelines following hereafter:

Ex-classification ATEX:

 II 1/2G Ex ia IIC T6...T3 Ga/Gb

 II 1/2D Ex ia IIIC T69° C...T161° C Da/Db

Number of the test certificate:

SEV 14 ATEX 0169X

IECEX SEV 14.0026X

CML 22 UKEX 2212X

Ex-classification FM approved:

IS / I, II, III / 1 / ABCDEFG / T6 Ta = 60° C
- 53 800 002; Entity

Project ID: 3021227

FM certificate number: FM16US0034X

FM18CA0021X

2.7 Ex-classification ATEX

2.7.1 Introduction

According to Directive 2014/34/EU (ATEX 114)¹⁾ Appendix I, InPro 6XXX*/**/* oxygen sensors are devices group II, category 1/2G and according to RL 99/92/EG²⁾ (ATEX 137) may be used in zones 0/1 or 1/2 or 1 or 2 and gas groups IIA, IIB and IIC that are potentially explosive due to combustible substances in above stated temperature classes.

For use/installation, the requirements of EN 60079-14 must be observed.

According to Directive 2014/34/EU (ATEX 114)¹⁾ Appendix I, InPro6XXX*/**/* oxygen sensors are devices group II, category 1/2D and according to RL 99/92/EG²⁾ (ATEX 137) may also be used in zones 20/21 or 21 or 22 that contain combustible dusts.

The digital sensor circuit is part of a common intrinsically safe system and is for operation connected to a separately certified transmitter.

The digital sensor circuit as part of an intrinsically safe system is isolated from the not-intrinsically safe electric circuits up to a maximum rated voltage of 375 V and from grounded parts up to a maximum rated voltage of 30 V.

2.7.2 Rated data

Measuring circuit

Intrinsically safe ignition protection class Ex ia IIC only for connection to a certified intrinsically safe circuit.

Analog O₂ oxygen sensor maximum values

$$U_i \leq 16 \text{ V}$$

$$I_i \leq 190 \text{ mA}$$

$$P_i \leq 200 \text{ mW}$$

$$L_i = 0 \text{ (effective internal inductance)}$$

$$C_i = 900 \text{ pF (effective internal capacitance)}$$

Digital O₂ oxygen sensor maximum values

Two-wire current circuit:

$$U_i \leq 16 \text{ V}$$

$$I_i \leq 30 \text{ mA}$$

$$P_i \leq 50 \text{ mW}$$

$$L_i = \text{negligible}$$

$$C_i = \text{negligible}$$

Note: The above maximum values are each the total of all individual circuits of the associated intrinsically safe power supply and transmitter.

¹⁾ For UK statutory Requirements SI 2016 No. 1107

²⁾ For UK regulations 7 and 11 of the Dangerous Substances and Explosive Atmospheres Regulations 2002 (DSEAR).

2.7.3 Special conditions

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1G applications, zone 0, is shown in the following table:

Temperature class	Max. ambient or media temperature
T 6	68 °C
T 5	80 °C
T 4	108 °C
T 3	160 °C

- The relationship between the maximum permissible ambient or media temperature and temperature class, for category 1D applications, zone 20, is shown in the following table:

Temperature class	Max. ambient or media temperature
T 69 °C	68 °C
T 81 °C	80 °C
T 109 °C	108 °C
T 161 °C	160 °C


- The capacitance and inductance of the connecting cable has to be considered.
- The O₂ Oxygen sensor type InPro 6XXX can be used in / with the fittings InFit™ 76*-*** or InTrac™ 7**-***, or in / with other suitable fittings in potentially explosive areas.
- The metal body of the O₂ Oxygen sensors, or the fittings InFit 76*-*** or InTrac 7**-***, or other appropriate fitting is optionally included in the routine pressure test of the system.
- The independent fitting used for installation of O₂ Oxygen sensor must be conductively connected to the equipotential bonding system.

2.8 Ex-classification – FM approved

2.8.1 Introduction

The following FM control drawing and the standards listed in section 2.8.4 must be observed, where applicable:

2.8.2 Rated data

For rated data please observe section 2.7.2

2.8.3 Special conditions

In type of protection intrinsically safe apparatus, the analytical probes (InPro 6000 Series with InFit 76X/Y and InTrac 7XX/YY Series housings) equipment is designated with the following specific conditions of use.

1. The oxygen sensor shall be installed in compliance with the enclosure, mounting, spacing and segregation requirements of the ultimate application, including a tool removable cover, and is suitable for use with FM Approved InFit 76X/Y and InTrac 7XX/YY Series Housings.
2. Process temperature no greater than +130 °C.
3. Maximum permissible working pressure is 12 barg (174 psig).
4. Potential Electrostatic Charging Hazard – To prevent the risk of electrostatic sparking, the non-metallic surface should only be cleaned with a damp cloth.
5. Enclosures containing titanium constitute a potential risk of ignition by impact or friction. Care must be taken into account during installation and use to prevent impact or friction.

2.8.4 Applied standards:

United States Standards

Title	Number	Issue Date
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	ANSI/UL 61010-1:2018	2018

Canadian Standards

Title	Number	Issue Date
Explosive Atmospheres – Part 0: Equipment – General Requirements.	CAN/CSA-C22.2 No. 60079-0	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”.	CAN/CSA-C22.2 No. 60079-11	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	CAN/CSA-C22.2 No. 61010-1	R2017

3 Product description

3.1 General information

The **oxygen sensor series InPro 6800 G / 6850i G / 6900i G / 6950i G** with integrated temperature probe are used for measurement of oxygen at low and medium concentrations.

InPro 6xxxi G sensors with ISM™ functionality offer Plug and Measure as well as enhanced diagnostics features.

3.2 Principle

Amperometric oxygen sensors:

The amperometric oxygen sensors of METTLER TOLEDO base all on the same measuring principle according to Clark. However different series are offered, which clearly differ in the number and arrangement of their electrodes and thus in their specifications:

- The measurement system of the InPro sensors consists of a working electrode (cathode), a counter electrode (anode), a reference electrode and a guard ring electrode. The measurement system is separated from the process medium by an oxygen permeable membrane.
- The measuring technique in principle is alike with all sensors.
- The transmitter supplies a constant voltage between cathode and anode.
- The electrolyte creates a conductive connection between the electrodes.
- The oxygen molecules migrate from the measurement medium through the membrane to the cathode to which the voltage is supplied and are reduced. At the same time, oxidation takes place at the anode.
- This causes a current to flow between the anode and cathode which is directly proportional to the partial pressure of oxygen (pO_2) in the process medium. With digital ISM sensors, the sensor itself converts the current into oxygen concentration and communicates the value to the transmitter.
- The guard ring electrode of the InPro 6900i G and 6950i G built into the sensor reduces the oxygen that migrates to the side of the cathode and can affect the measurement. The guard ring electrode therefore enables traces of oxygen to be precisely determined at even the lowest concentrations.

Note: Please refer to "10 Theory of the polarographic sensor" for further information.

ISM Sensors:

All oxygen sensors with the index "i" (6850i G, 6900i G, 6950i G) are equipped with ISM, Intelligent Sensor Management.

Principle: In the sensor head a chip is integrated, which takes over the entire monitoring and control of the sensor and, beyond that, stores all sensor data. This chip is responded via the transmitter.

The following data are available permanently in the sensor:

- type of sensor
- serial no.
- software version
- hardware version
- order no.
- operating time
- calibration time and calibration date
- calibration table

To check the system, the following indicators are monitored:

- temperature
- slope
- zero current
- air current
- polarization voltage
- cathode current and voltage

Based on these informations, the transmitter calculates a wear monitor and displays it depending upon transmitter type differently. (see the respective instruction manuals)

ISM enables the connection of the digital sensor to the iSense™ Asset Suite. The user-friendly software allows to manage all sensor information and store it in a data base. Furthermore, the sensor can be calibrated via the iSense software.

3.3 Enhanced diagnostics (only for ISM types available)

DLI Dynamic Lifetime Indicator The dynamic lifetime indication allows an estimation, when the sensor inner body is at the end of his lifetime, based on the actual stress it is exposed to. The sensor permanently takes the averaged stress of the past days into consideration and is able to increase/decrease the lifetime accordingly.

The following parameters affect the lifetime indicator:

Dynamic parameters:

- Temperature
- Oxygen value

Static parameters:

- Calibration history
- Zero and Slope

The sensor keeps the information stored in the built-in electronics and can be retrieved via a transmitter or the iSense asset management suite.

The alarm will be reset if the Lifetime Indicator is not 0 days anymore (e.g. after connecting a new sensor or changes on the measurement conditions)

For amperometric oxygen sensors, the lifetime indicator is related to the inner-body of the sensor. After exchanging the inner-body, reset the lifetime indicator in the ISM setup menu (Menu Configure).

TTM Time to Maintenance:

This timer defines when the next electrolyte or membrane replacement should be done to maintain the best possible measurement performance. The timer depends on the DLI parameters.

The sensors InPro 6850i G / 6950i G are equipped additionally with an electrolyte level monitor, ELM.

The electrolyte level monitor, ELM, automatically sets the TTM to zero as soon as a critical electrolyte level in the sensor tip is reached.

The alarm needs to be reset in the menu "ISM set-up". For oxygen sensors, the time to maintenance indicates a maintenance cycle for the membrane and electrolyte.

Note: The recommended initial value for the maintenance interval will be uploaded from the sensor to the transmitter and can be adapted according to the application experience of the user.

ACT Adaptive Cal Timer:

This timer estimates when the next calibration should be performed to keep the best possible measurement performance. The timer is influenced by significant changes on the DLI parameters.

Note: The recommended initial value for the calibration interval will be uploaded from the sensor to the transmitter and can be adapted according to the application experience of the user.

3.4 Scope of delivery

Each sensor is supplied fully assembled and factory-tested for correct function together with:


- a quality control certificate
- inspection certificates 3.1 (complying with EN 10204)

Digital sensors must be filled with electrolyte before use.

See electrolyte order details at section “9.3 Spare parts”

3.5 Equipment features

12 mm sensor

METTLER TOLEDO ISM O₂ sensors are supplied with fitted membrane body but without electrolyte and covered with the protection cap and have been checked for proper function.

4 Installation

4.1 Mounting the sensor


Important! Remove the protection cap before mounting the sensor.

Mounting the sensor in a housing

Please refer to the instruction manual of your housing explaining on how to mount the sensor in place.

Mounting the sensor directly on a pipe or a vessel

The **12 mm sensors** can be mounted directly through a socket with inside thread Pg 13.5 and securely tightened via the Pg 13.5 threaded sleeve.

4.2 Connection

4.2.1 Connecting the sensor to a AK9 cable



Note: The sensor is connected to the transmitter via a AK9 cable. The AK9 cable ensures a secure connection between the transmitter and the sensor under harsh industrial conditions. The robust watertight IP 68 connector housing guarantees maximum process safety.


Do not touch the sensor at the AK9 connector plug!

Tightly screw the plug to fasten the two parts.

4.2.2 Connecting the cable to the transmitter

Note: For connecting the AK9 or VP cable to the terminals of the transmitter, please refer to the instructions given in the METTLER TOLEDO transmitter manual.

5 Operation

Important! Before using the sensors for the first time, the electrolyte must be filled in (see Chapter "6.2").

5.1 Start-up and polarizing

Important! The protection cap must be removed before mounting the sensor in the process.

When the system is operated for the first time or if the sensor has been disconnected from the voltage source (transmitter or O₂ Sensor Master) for longer than 5 minutes, the sensor has to be polarized prior to calibration by connecting it to the operating O₂ transmitter or to a sensor master. After 6 hours, the sensor is fully polarized and ready for operation.

Attention: The InPro 6950i G should never be polarized at air!

During Polarization, we suggest to use the protection cap filled with cleaning and conditioning solution always (Chapter "7 Storage").

A shorter polarization period is sufficient if the sensor has been disconnected for only a few minutes. The following table serves to establish the correct polarization time in relation to the depolarization time.

Depolarization time ¹ t _{depol} [Min.]	Minimum required polarization time ² [Min.]
t _{depol} > 30	360
30 > t _{depol} > 15	6 × t _{depol}
15 > t _{depol} > 5	4 × t _{depol}
t _{depol} < 5	2 × t _{depol}

1 Depolarization time: Time span in which the polarization voltage is cut off from the sensor. This is the case:

- during the time the cable is disconnected or no transmitter or sensor master is connected to the cable, or the transmitter has been disconnected from the current supply.
- after changing the electrolyte or membrane body. In this case at least 6 hours of polarization must follow.

2 Polarization time: Time span during which the sensor is under a polarization voltage.

Important! Setting of the polarization voltage on the transmitter for correct measurements:

- Standard applications InPro 6800 G / 6850i G:
– **675 mV**
- Measurement of permanently low oxygen concentrations (< 500 ppb in liquids or < 10,000 ppm [vol.] in gases) in the presence of volatile acidic components (e.g. carbon dioxide during measurements in breweries) e.g. InPro 6900i G / 6950 i G:
– **500 mV**

Note: To ensure the supply of the correct polarization voltage the transmitter must be set according to the transmitter manual.

5.2 Calibration

5.2.1 Purpose of calibration

Each oxygen sensor has its own individual slope and own individual zero point. Both values are subject to change, for example, through electrolyte consumption or after exchange of electrolyte or membrane body. To ensure high measurement accuracy of the sensor, a calibration must be carried out regularly, but at least after each change of electrolyte or membrane. Prior to calibration, the sensor has to be polarized for at least 6 hours.

Important: Please remove the protection cap from the sensor, rinse the sensor with water and dry it wait at least 10 minutes before starting calibration. Store sensor upright with membrane pointing down.

Note: To check if your sensor needs a recalibration, you may dry it and take it in the air to check that the reading is close to 100 %. If not, then the sensor needs a new calibration.

With ISM sensors all calibration data are stored in the sensor.

5.2.2 What you have to know for calibration

General remarks:

- **For calibration in gas (e.g. air), the sensor membrane must be dry outside**, since adhering water drops can falsify the measured oxygen value.
- Make sure that the **oxygen saturation index** of the calibration medium is **correct and remains constant** during calibration.
- For correct calibration, a minimum flow rate of the calibration medium is necessary.
- Make sure that all other parameters, such as temperature and pressure, are constant.

For continuous applications, we recommend periodic recalibration in line with your requirements

on accuracy, the type of process in operation and your own experience. The frequency of the need for recalibration depends very much on the specific application, and therefore appropriate intervals cannot be exactly defined here.

5.2.3 Single point calibration

For applications measuring more than 1% Vol. O₂ gas a single point calibration is recommended. In this case the zero point is set manually to 0 nA. For measuring below 1% Vol. O₂ a zero point calibration with e.g. nitrogen gas is recommended.

By carrying out a single point calibration, the slope or offset of the sensor can be set. The calibration medium can be gas with known water-vapor saturation (e.g. water-vapor saturated air).

After the sensor signal has stabilized, the complete measuring system can be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O₂, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual of the transmitter).

Zero-point calibration is recommended for the InPro 6950 i G only when measuring less than 125 ppm (by vol. O₂) in gases.

Attention! Incorrect zero point calibration is a frequent source for measurement error. For correct calibration, we recommend the use of carbon dioxide gas with a level of purity of at least 99.9995 % for measurements in CO₂ applications. For all other applications we recommend to use nitrogen (N₂) gas.

After the sensor signal has stabilized (after 6 ... 12 hours), the sensor can be calibrated through the transmitter to zero % value of the desired measurable variable, e.g. 0 % air, 0.0 % O₂, or 0.0 ppm (see instruction manual for the transmitter).

5.2.4 Dual point calibration

By carrying out a dual point calibration both slope and zero point of the sensor can be set.

Important! In case of a **dual point calibration**, **always start with the zero point calibration** before calibrating the slope.

6 Maintenance

6.1 Inspection of the sensor

6.1.1 Visual inspection

To check your sensor, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Before calibration always examine the membrane foil for signs of damage. The foil must be intact and clean. Dirty membranes should be wiped clean using a soft, moist tissue.

Note: An undulated membrane has no influence on the sensor performance, assuming the membrane is intact.

- The membrane body must be replaced if the sensor has too long a response time, the reading is unstable or subject to drift, and if the sensor cannot be calibrated or the membrane shows sign of mechanical damage.
- Check the cathode area for discoloration, contamination or cracks in the glass. If necessary rinse with demineralized water and clean with a clean soft brush or soft paper tissue.

Attention! Do not use any cleaning agents containing alcohol. This could damage the sensor or lead to fault current.


Attention! The glass body is fragile and sensitive to vibration.

6.1.2 Using the METTLER TOLEDO O₂ sensor master

We recommend to use the METTLER TOLEDO O₂ sensor master to check the status of your sensor as follows:

- Connect the sensor to the O₂ sensor master.

The figure shows the Sensor Master (Order No. 52 200 892) for the analog sensor InPro 6800 G. For the digital ISM sensors use the Sensor Master digital ISM with order no. 52 206 329. For the digital Sensor Master the procedure is similar.

As soon as the sensor is connected to the O₂ sensor master, the polarization function is automatically activated. Please note: if the sensor was disconnected from the transmitter for longer than 5 minutes, the sensor must be polarized first (polarizing time see Section "5.1") to get representative test results.

- **Battery Check:**
Push the switch to the left. If the battery is ok and the O₂ sensor master is operational the green LED lights up. Otherwise, please consult the instruction manual of the O₂ sensor master.
- **Sensor Check:**
For this test the O₂ sensor must be **fully polarized** and the **membrane of the sensor tip must be dry and clean outside**.

Expose the sensor (connected to O₂ sensor master) to air. By pushing the switch to the right to position "Sensor check", the O₂ sensor master checks whether the electrode current for air measurement is within the specified range, i.e. 2500 to 6000 nA for InPro 6950i G.

	Air current	Zero current in % of the air current
6800 G	50 – 110 nA	< 0.1
6850 i G	50 – 110 nA	< 0.1
6900 i G	250 – 500 nA	< 0.03
6950 i G	2500 – 6000 nA	< 0.025

If the green LED lights up the current for air measurement is within the specified range.

If the LED does not light up, you should check the battery of the O₂ sensor master (see instruction manual "Accessories"). If the battery is working, than there is probably a problem with your sensor. You should change the electrolyte and/or the membrane body of your sensor. If after a membrane change the LED still does not light up, this means that there is something wrong with the interior body of the sensor. You should then change it (see Section 6.2).

Important! The Sensor Check function only verifies the correctness of the electrode current for air measurement. In order to be absolutely sure of the functionality of the sensor, the residual signal in an oxygen free medium should also be controlled (see "Section 6.1.3").

6.1.3 Testing the sensor via a transmitter

A periodic zero current measurement (**no zero point calibration!**) is recommended for verification of proper sensor function.

Note: At the time you carry out the zero current measurement, the sensor must be polarized.

Zero current measurement can be done by using zeroing gel (order no. 30 300 435) or nitrogen (N₂) or carbon dioxide (CO₂) calibration gases with a purity of at least 99.995 %, alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 10 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If the measured values are too high, this suggests a depleted electrolyte or a defective membrane. In the first instance replace the electrolyte, and in the second case exchange both the membrane body and the electrolyte accordingly.

If after such procedures the above mentioned values are still not reached, replace the interior body. If this does not solve the problem send the sensor to your local METTLER TOLEDO representative for inspection.

6.1.4 ISM design

The integrated ISM function allows an advanced monitoring of the sensor. The following parameters are stored in the sensor:

- serial no.
- type of sensor
- order no.
- calibration data
- enhanced diagnostic data
- slope
- zero point

After powering up of the sensor the following automatic test procedures are running:

- digital communication
- plausibility of the stored calibration data

6.2 Changing the electrolyte, the membrane body or the interior body

Note: The InPro 6900i G and 6950i G use a special electrolyte. This electrolyte ensures fast response time and together with the guard ring increases the signal stability of the sensor. The electrolyte must be exchanged on a regular basis or when the sensor has been exposed to air for more than 24 hours without a watering cap filled with conditioning solution.

If the membrane and/or the interior body exhibits signs of failure (long response time, increased current in an oxygen-free medium, mechanical damage, etc.) the membrane body and/or the interior body has to be replaced.

Warning! The O₂ electrolyte has a slight alkaline pH value. Contact of electrolyte with mucous membrane or eyes is to be avoided. **Therefore protective gloves and safety glasses have to be worn for the following dismantling works.**

If such contact occurs, the affected area should be well rinsed with water. In the case of accident, or should ever any adverse signs appear, get immediate medical attention.

When changing the electrolyte, the membrane body or the interior body, please observe the following instructions (see also the following illustration):

Attention! Make sure that this maintenance step is carried out in **clean place**.

1. Unscrew the cap sleeve from the sensor shaft and carefully pull it off the sensor.
2. Pull off the membrane body from the interior body. If it is tight-fitted, eject by pushing it with the flat finger tip. Before electrolyte is refilled, the membrane body must be removed from the cap sleeve!
3. Rinse the interior body with demineralized water and carefully dab it dry with a paper tissue.

Note: steps 4 to 7 may only be carried out when changing the interior body.

4. Unscrew the retainer nut of the interior body with an adjustable wrench or with a 9 mm wrench.
5. Remove the interior body by pulling it out of the sensor shaft. If necessary use a plier.


Warning! Do not twist the interior body. Otherwise the connection pins can be damaged.

6. Insert the new interior body in the sensor shaft. Turn the interior body in the shaft until the slit of the interior body is aligned with the pin placed in the shaft.

7. Press the body in the shaft and screw the new retainer nut in place. Only hand tighten the retainer nut using a small wrench. Applying a too high torque (> 2 Nm) may damage the sensor.
8. Examine the O-rings visually for mechanical defects, and replace if necessary.
9. Half-fill the new membrane body with O₂ electrolyte.

Note: The electrolyte bottle is equipped with a special pouring system. To ensure proper functioning, hold the bottle vertically, upside-down.

Note: make sure that all air bubbles are removed from the membrane body. Air bubbles can be removed by carefully tapping on the membrane body.


10. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and have to be removed with a paper tissue.

Important! No electrolyte, sample media or contamination may be present between the membrane body and the cap sleeve. Please check carefully!

11. Carefully slip the cap sleeve over the fitted membrane body, holding the sensor in a vertical position and screw it tight. The cap sleeve must be clean and dry.
12. After each exchange of electrolyte or membrane body, the sensor has to be repolarized and recalibrated before it is ready for use.

Replacement of the interior body

7 Storage

When the sensor is stored without polarization for more than one week, the electrolyte has to be removed.

Always use clean water to clean the sensor membrane from outside.

Important: Never use alcohol or isopropanol or mixtures with these organic solvents.

Prepare the cleaning and conditioning solution as follows: Dissolve one tablet in 40 ml of deionized water and wait 5 minutes for the tablet to be completely dissolved. Fill the protection cap with this solution and place it over the tip of the sensor. This solution has some cleaning properties which will keep the membrane free of microorganisms. In case you do not have any cleaning and conditioning set, you may also use checking gel or deaerated water in the protection cap. Before mounting the sensor in the process, always remove the protection cap and rinse the tip of the sensor with water.

Attention! If the storage period of the sensor without current supply (transmitter, sensor master) exceeds 1 week, the sensor should be stored dry, i.e. without any electrolyte in the membrane body. A sensor being stored dry (without electrolyte in the membrane body) may on no account be connected to the O₂ sensor master or any other polarization module.

8 Product specification

8.1 Certificates

Each sensor is delivered with a set of **3.1 certificates** (complying with EN 10204).

All wetted metal parts (sensor shaft, cap sleeve and membrane body) are identified with a engraved symbol corresponding to the heat number on the paper certificate delivered with the sensor.

Each wetted metal part (sensor shaft, cap sleeve and membrane body) is polished in order to get a surface roughness lower than 0.4 µm (16 µin). This represents a roughness grade number of N5 (according to ISO 1320:1992).

8.2 Specifications

InPro 6800 G / 6850i G / 6900i G / 6950i G	
Measurement principle	Polarographic Clark electrode
Working conditions	
Pressure resistance (measurement)	6800 G: 0.2 ... 9 bar 6850i G: 0.2 ... 9 bar 6900i G: 0.2 ... 9 bar 6950i G: 0.2 ... 9 bar
Mechanical pressure resistance	Max. 12 bar
Temperature limits (mechanical)	- 5 ... 121 °C [23 ... 250 °F]
Operating temperature limits (ambient)	6800 G: 0 ... 70 °C [32 ... 158 °F] 6850i G: 0 ... 70 °C [32 ... 158 °F] 6900i G: 0 ... 70 °C [32 ... 158 °F] 6950i G: 0 ... 40 °C [32 ... 104 °F] (sterilizable)
Specification	
Temperature compensation	Automatic
Cable connection	AK9
O-ring material	Silicone, FDA and USP Class VI approved Kalrez
Membrane material	PTFE / Silicone / PTFE (reinforced with steel mesh) Kalrez
Material sensor body (wetted parts)	316L stainless steel
Surface roughness of wetted metal parts	N5 (R _a = 0.4 µm [16 µin])
Quick disconnect interior body	Standard
Cathode	Pt
Anode	6800 G: Ag 6850i G: Pt 6900i G: Ag 6950i G: Pt
Guard ring	6800 G: - 6850i G: - 6900i G: Pt 6950i G: Pt
Reference	Ag
Dimensions	
Sensor diameter Ø	6800 G: 12 mm 6850i G: 12 mm 6900i G: 12 mm 6950i G: 12 mm
Immersion length (a) for 12 mm sensor	120, 220 mm [4.7, 8.66"]
Performance	
Measurement range in O ₂ Gas	6800 G: 300 ppm to 100% (1 bar) 6850i G: 300 ppm to 100% (1 bar) 6900i G: 50 ppm to 100% (1 bar) 6950i G: 5 ppm to 20% (1 bar)
Accuracy	6800 G: ≤ ± [1 % + 300 ppm] 6850i G: ≤ ± [1 % + 300 ppm] 6900i G: ≤ ± [1 % + 50 ppm] 6950i G: ≤ ± [1 % + 5 ppm]
Response time at 25 °C / 77 °F	90 % of final value in < 30 s
Sensor signal for ambient air at 25 °C / 77 °F	6800 G: 50 ... 110 nA 6850i G: 50 ... 110 nA 6900i G: 250 ... 500 nA 6950i G: 2500 ... 6000 nA
Residual signal in oxygen-free medium	6800 G: < 0.1 % of the signal 6850i G: < 0.1 % of the signal 6900i G: < 0.03 % of the signal 6950i G: < 0.025% of the signal

Certificates

EHEDG, 3A	Yes
3.1 (EN 10204)	Yes
ATEX Certificate	Yes
IECEX Certificate	Yes
UKEX Certificate	Yes
FM Approval	Yes
FDA / USP VI	Yes
Quality Control	Yes

Compatibility

with METTLER TOLEDO transmitters	see "Section 9.2"
with METTLER TOLEDO housings	see "Section 9.6"

9 Ordering information

For more detailed information refer to the technical data sheet. Ask your local distributor.

9.1 Sensors

Part	Order No.
InPro 6800 G / 12 / 120	52 206 425
InPro 6800 G / 12 / 220	52 206 426
InPro 6800 G / 12 / 120 / Ka	52 206 427
InPro 6800 G / 12 / 220 / Ka	52 206 428
InPro 6800 G / 12 / 120 / C22	52 206 429
InPro 6800 G / 12 / 220 / C22	52 206 430
InPro 6850i G / 12 / 120	52 206 431
InPro 6850i G / 12 / 220	52 206 432
InPro 6850i G / 12 / 120 / Ka	52 206 433
InPro 6850i G / 12 / 220 / Ka	52 206 434
InPro 6850i G / 12 / 120 / C22	52 206 435
InPro 6850i G / 12 / 220 / C22	52 206 436
InPro 6900i G / 12 / 120	52 206 437
InPro 6900i G / 12 / 220	52 206 438
InPro 6900i G / 12 / 120 / Ka	52 206 439
InPro 6900i G / 12 / 220 / Ka	52 206 440
InPro 6950i G / 12 / 120	52 206 443
InPro 6950i G / 12 / 220	52 206 444

9.2 Accessoires

Part	Order No.
Sensor Master digital ISM	52 206 329
Sensor Master analog	52 200 892
InPro 6900 electrolyte pack (3 × 5 ml) for InPro 6900 i G	30 298 425
InPro 6950 electrolyte pack (3 × 5 ml) for InPro 6950 i G	30 298 426

9.3 Spare parts

Part	Order No.
Interior body InPro 6800 G	52 206 449
Membrane body T-6800 Gas	52 201 151
Membrane kit T-6800 Gas	52 201 149
Membrane body T-6800 Gas Ka	52 201 158
Membrane kit T-6800 Gas Ka	52 201 159
Membrane body T-6800 Gas C22 Ka	52 201 163
Membrane kit T-6800 Gas C22 Ka	52 201 164
Interior body InPro 6850i G	52 206 450
Membrane body T-6850i G	52 206 453
Membrane kit T-6850i G	52 206 454
Membrane body T-6850i G Ka	52 206 455
Membrane kit T-6850i G Ka	52 206 456
Membrane body T-6850 i G C22 Ka	52 206 457
Membrane kit T-6850i G C22 Ka	52 206 458
Interior body InPro 6900i G	52 206 451
Membrane body T-6900i Gas	52 206 459
Membrane kit T-6900i Gas	52 206 460
Membrane body T-6900i Gas Ka	52 206 461
Membrane kit T-6900i Gas Ka	52 206 462
Interior body InPro 6950i G	52 206 452
Membrane body T-6950i Gas	52 206 465
Membrane kit T-6950i Gas	52 206 466

9.4 Cables

Part	Order No.
AK9 coax cable with K8S connector, 1 m	59 902 167
AK9 coax cable with K8S connector, 3 m	59 902 193
AK9 coax cable with K8S connector, 5 m	59 902 213
AK9 coax cable with K8S connector, 10 m	59 902 230
AK9 coax cable with K8S connector, 20 m	52 300 204
VP cable VP6-ST / 1 m	52 300 111
VP cable VP6-ST / 3 m	52 300 112
VP cable VP6-ST / 5 m	52 300 113
VP cable VP6-ST / 10 m	52 300 114

9.5 Recommended housings

Housing (12 mm Ø)	Order No.
Static housing	Please ask your local METTLER TOLEDO representative.
InFit 761 e CIP	
Retractable housings	
InTrac 777 e	
InTrac 797 e	
Immersion housing	
InDip™ 550	

Note: The housings are available in different versions. Please contact your local METTLER TOLEDO sales organization to get the right ordering information.

10 Theory of the polarographic sensor

10.1 Introduction

Two types of electrodes are employed in analytical work: **potentiometric** and **amperometric** electrodes.

- Potentiometric electrodes develop a voltage generated by the activity of a particular ion. Examples of such electrodes are glass electrodes (like pH electrodes) and most ion-selective electrodes. Their individual potentials cannot be determined. The measurable quantity is the difference of potential between the measuring electrode and an inert reference electrode. The potential of the reference electrode must be constant.


All potentiometric electrodes are subject to Nernst's law and for this reason electrodes and measuring instruments are in most cases interchangeable.

An important requirement of potentiometric measurements is the virtually currentless determination of the electrode voltage. During measurement no chemical reaction occurs and the solution remains in equilibrium.

- In the case of **amperometric electrodes**, such as the **oxygen electrode**, activity measurement is based on a current measurement.

The conventional oxygen electrode consists of a cathode and an anode conductively connected by an electrolyte. A suitable polarization voltage between the anode and the cathode selectively reduces the oxygen at the cathode.

Cathode reaction:

These chemical reactions result in an electric current which is proportional to the oxygen partial pressure (p_{O_2}). The oxygen electrode reduces constantly oxygen. Thereby, the concentration of dissolved oxygen is reduced. By diffusion this oxygen is replaced. The viscosity and flow rate of the solution are therefore important parameters.

The electrode current of an oxygen electrode is determined not only by the oxygen partial pressure but by many other electrode parameters. The electrode currents of different electrode types may differ by several powers of ten. For this reason oxygen electrodes and amplifiers cannot be freely interchanged.

10.2 Principle of the design of an oxygen electrode

There are two main types of oxygen electrodes:

- Electrodes **without** membrane
- Electrodes **with** gas-permeable membrane (Clark Principle)

The membrane electrode according to Clark is today most widely used. As compared to the electrode without membrane it possesses the following advantages:

- Oxygen measurement in gases and solutions
- No mutual contamination of electrode and solution
- No or little dependency on flow



In the case of the Clark electrode, geometrical configuration is very important. In particular, the thickness of the electrolyte film between the cathode and the membrane must be within narrow tolerances so as to ensure good linearity and a low zero current (current in a nitrogen atmosphere).

METTLER TOLEDO oxygen sensors are available in different designs:

Type A, 2-electrodes system, InPro 6800 G

InPro 6800 G for medium and high oxygen concentrations. Cathode and anode/reference. Anode and reference are united in a silver / silver chloride electrode.


The following equilibration reaction takes place at the anode resp. reference:

Type B, 3-electrodes, InPro 6850i G

InPro 6850i G for medium and high oxygen concentrations. The reference corresponds the conventional silver /silver chloride anode. The anode is a platine electrode and separated from the reference.


The following reaction takes place at the anode:

The reference is a silver /silver chloride electrode. After polarization there is a stable equilibrium between electrode surface and electrolyte. No net-reaction takes place.


Reference reaction in equilibrium:

Reaction: No current



Type C, 3-electrodes, InPro 6900i G

For oxygen measurements in the lower ppb range. Here, the anode and reference are united in a silver / silver chloride electrode (as with type A). The sensors are equipped with an additional guard ring around the cathode. Like the cathode with the anode, this forms a closed electric circuit, which prevents that oxygen diffuses from the side to the cathode and falsifies the measurement result.


In all sensors the function of the cathode is identical.

Between cathode and reference a constant voltage of -500 resp. -675 mV (polarization voltage) is applied. Thereby a reduction of O₂ into OH⁻ at the cathode takes place.




The total equation thus is:


Type D, 4-electrodes, InPro 6950i G

The InPro 6950i G is intended for measuring permanent low oxygen level down to 0.1 ppb. In terms of measurement principle it is a combination of the InPro 6850i G and InPro 6900i G. It consists of 4 electrodes. The anode (platinum) and reference (silver / silver chloride) are separated in two electrodes. A guard ring is placed around the cathode. The cathode has the highest diameter from all amperometric oxygen sensors from METTLER TOLEDO.

The electrodes show the following reactions:

Due to the larger surface of the cathode a higher current appears during measurement. This allows for measurement of oxygen levels below 1 ppb.

As a consequence of the reaction at the cathode and the anode, a current flows between cathode and anode (4 electrons per oxygen molecule), which is proportional to the oxygen partial pressure at the cathode.

The amperage of the current flow depends on the surface of the cathode. Typical values are:

	Air current	Zero current in % of the air current
6800 G	50 – 110 nA	< 0.1
6850 i G	50 – 110 nA	< 0.1
6900 i G	250 – 500 nA	< 0.03
6950 i G	2500 – 6000 nA	< 0.025

This current is measured and converted and indicated in the transmitter into an oxygen value. In ISM sensors this computation takes place in the sensor and the transmitters displays this value.

10.3 Parameters determining current

The quantity of oxygen diffused in and the magnitude of the electrode current are influenced by the following parameters:

- Oxygen partial pressure of the solution
- Membrane material and thickness
- Size of cathode
- Polarization voltage
- Temperature
- Flow conditions in the solution


Fick's law gives the mathematical relationship between these parameters:

$$I = k \times D \times a \times A \times \frac{pO_2}{X}$$

- I = Electrode current
- k = Constant
- D = Coefficient of diffusion of O₂ in the membrane
- a = Solubility of O₂ in the membrane material
- A = Cathode surface
- pO₂ = Oxygen partial pressure in the solution
- X = Thickness of gas-permeable membrane

10.4 Polarization voltage

The voltage between the anode and the cathode is so selected that the oxygen is fully (> A, see polarogram) reduced while other gases are unaffected (< D). The ideal voltage for the Pt / Ag / AgCl or Pt / Pt / Ag / AgCl system is between –500 and –750 mV.

The polarization voltage should remain as constant as possible. Besides a constant source of voltage the following prerequisites must be met: The electrical resistance of electrolyte film must not exceed a particular value so as to avoid a voltage drop.

The anode must have a large surface so as to prevent polarization of the anode by the electrode current.

10.5 Temperature

The temperature dependence of the current passing through an oxygen electrode when referred to a constant O₂ partial pressure is determined mainly by the properties of the gas-permeable membrane.

10.6 Dependence on flow

With most oxygen electrodes the electrode current is smaller in stagnant solutions than in agitated ones. In consequence of the oxygen consumption of the electrode, a reduction of oxygen takes place outside the membrane in close proximity to the cathode. The missing oxygen is replaced by diffusion. If the electrode current is strong, the solution cannot fully restore the reduced oxygen by diffusion. This results in an electrode current weaker than would correspond to conditions in the solution. In agitated solutions the oxygen is conveyed to the surface of the membrane not only by diffusion but additionally by the flow (convection). In that case no oxygen impoverishment occurs at the membrane surface.

A high degree of flow dependence occurs mainly with large cathodes, thin and highly permeable membranes, i.e. where electrode currents are large.

The problem of flow dependence is often solved by prescribing a minimum flow rate.

In METTLER TOLEDO InPro 6950 i G electrodes, the thin PTFE membrane determining the electrode current (i.e. the actual measuring signal) is separated from the sample solution by a relatively thick silicone membrane. This latter is highly permeable to oxygen molecules and thus acts as an oxygen reservoir. The diffusion of oxygen out of the sample solution into the silicone membrane is spread over a wide area. Since this results in less oxygen being extracted from the sample solution per unit area, the double PTFE/silicone membrane forms an effective buffer against disturbances due to hydrodynamic flow.

This membrane coupled with the guard ring and the special electrolyte ensures excellent signal stability even when the hydrodynamic flow stops (as on a beer filling line for example).

10.7 Oxygen partial pressure – oxygen concentration

The electrode current is dependent on the partial pressure of oxygen andw the oxygen permeability of the membrane. Conversion of partial pressure into concentration of oxygen depends on the measurement medium (measurement in liquids or gases).

Measurement in liquids

When measuring in liquids, the concentration of oxygen also depends on the solubility of the oxygen in the measurement medium. However, since this is not measured by the sensor current, the concentration of oxygen must be calculated in the transmitter. To do this, Henry's law is applied which states that the concentration of oxygen is proportional to the partial pressure of oxygen (pO_2).

$$CI = pO_2 \times a$$

a = Solubility factor

If "a" is constant, the oxygen concentration can be determined by means of the electrode. This applies at constant temperature and with dilute aqueous solutions such as drinking-water.

The solubility factor is strongly influenced not only by the temperature but also by the composition of the solution:

Medium, saturated with air	Solubility at 20°C (68°F) and 760 mm Hg
Water	9.2 mg O ₂ /l
4 mol/l KCl	2 mg O ₂ /l
50 % Methanol-water	21.9 mg O ₂ /l

Although the solubilities vary widely, the oxygen electrode gives the same reading in all three solutions.

Thus, determination of the oxygen concentration is only possible with constant and known solubility factors "a".

Solubility may be determined by a Winkler titration or the method developed by Käppeli and Fiechter.

Measurement in gases

The concentration of oxygen when measuring in gases is always given as a proportion by volume of the gas composition. Common units are % (by vol.) ppm (by vol.).

They can simply be converted from one unit of measurement to the other.

Example:

The percent by volume of the composition of air is generally known. For example, air contains 20.95 % oxygen. This corresponds to 209,500 ppm (by vol.). (Conversion: ppm value = 10,000 × value in %)

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US-Normen

Titel	Nummer	Ausgabe-Datum
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety "i"	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	ANSI/UL 61010-1:2018	2018

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Titel	Nummer	Ausgabe-Datum
Explosive Atmospheres – Part 0: Equipment – General Requirements.	CAN/CSA-C22.2 No. 60079-0	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety "i".	CAN/CSA-C22.2 No. 60079-11	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	CAN/CSA-C22.2 No. 61010-1	R2017

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Titre	Numéro	Date de publication
Approval Standard for Electrical Equipment for Use in Hazardous (Classified) Locations – General Requirements	FM Class 3600	2022
Approval Standard for Intrinsically Safe Apparatus and Associated Apparatus for Use in Class I, II & III, Division 1, Hazardous (Classified) Locations	FM Class 3610	2021
Approval Standard for Electrical Equipment for Measurement, Control and Laboratory Use	FM Class 3810	2021
Explosive Atmospheres – Part 0: Equipment – General Requirements	ANSI/UL 60079-0:2020	2020
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”	ANSI/UL 60079-11:2018	2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements	ANSI/UL 61010-1:2018	2018

Normes canadiennes

Titre	Numéro	Date de publication
Explosive Atmospheres – Part 0: Equipment – General Requirements.	CAN/CSA-C22.2 No. 60079-0	2019
Explosive Atmospheres – Part 11: Equipment Protection by Intrinsic Safety “i”.	CAN/CSA-C22.2 No. 60079-11	R2018
Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use – Part 1: General Requirements.	CAN/CSA-C22.2 No. 61010-1	R2017

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