General instructions for pH combination electrode Type Xerolyt Electrodes

cs	Návod k obsluze	2
da	Brugsanvisning	4
de	Betriebsanleitung	6
en	Instruction manual	8
es	Manual de instrucciones	10
fr	Instructions d'utilisation	12
fi	Käyttöohje	14
hu	Használati utasítás	16
it	Istruzioni per l'uso	18
ja	取扱説明書	20
ko	사용 매뉴얼	22
nl	Gebruikershandleiding	24
pl	Instrukcja obsługi	26
pt	Manual de instruções	28
ru	Инструкция по эксплуатации	30
sv	Bruksanvisning	32
th	คู่มือคำแนะนำ	34
zh	说明书	36

Type Xerolyt Electrodes
59 905 939 B

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General instructions for Xerolyt-Electrodes

1 Introduction

These general instructions must be readily available to operators of the electrode at all times. Before using the electrode, these instructions must be read and understood by all operating personnel.

Owing to their special characteristics, pH/redox electrodes with a Xerolyt™ reference system are suitable for the following applications:

- Measurements in highly contaminating solutions
- Measurements in emulsions and suspensions
- Measurements in effluents/wastewater
- Measurements in solutions containing proteins
- Measurements in solutions with a high sulfide content
- Measurements in hot alkaline solutions (with Type HA pH glass)
- Installation in measuring and control loops with wide pressure fluctuations.

Xerolyt electrodes are **not** suitable for the following applications:

- Measurements where the electrode is steam-sterilized
- Measurements with large, long-term, temperature variations
- Measurements in solutions with a pH value of $\text{pH} < 2$
- Measurements in solutions containing chlorine

2 Safety remarks

The Xerolyt electrodes are intended solely for the purpose described in §1 above. The electrodes consist partly of glass. If they are to be cleaned or calibrated using acid or alkaline solutions, protective clothing (goggles and gloves) must be worn. In addition, local and in-house-safety regulations must be observed. We recommend that the electrodes only be put into operation in conjunction with original products/accessories from METTLER TOLEDO. The electrodes may be operated and serviced only by personnel familiar with the equipment in question and who have read and understood these instructions.

3 Product description

The following data apply only to the Xerolyt reference system. (With single-rod pH or redox electrodes assemblies, the specification of the indicator electrode should also be taken into consideration.)

The markings on each electrode provide the following information:

XXX-DXK-S8/yyy designation/length mm

xxx = 405 combination electrode

xxx = 363 reference electrode

pH 2 ... 14, 0° ... + 110 °C = pH measuring range/temperature range

Order number: xx xxx xx xx

4 Installation and preparation for use

- a Carefully inspect the electrode on delivery for mechanical damage to the membrane or the shaft. Electrodes showing any defects should be sent back in the original packaging, accompanied by the relative invoice number.
- b Remove the watering cap and rinse electrode with distilled water. After rinsing, a pH electrode unit should be dried by careful dabbing with a paper handkerchief or Kleenex tissue. Do not rub the glass membrane, since this can lead to very sluggish response.
- c Remove any air bubbles present in the membrane chamber of the glass electrode by lightly swinging round in a vertical plane.
- d pH electrode assemblies with dry glass membranes must be conditioned by immersion for several hours in 3 mol/l KCl in the supplied watering cap (use approx. 1 ml).

5 Calibration and measurement

Xerolyt electrodes require 2-point calibration. To calibrate, first remove the watering cap. The electrode is then dipped in succession into two different buffer solutions with given pH reference values, and the transmitter calibrated to these buffer values (for example buffer pH 4.01 and pH 7.00). Please refer to the instruction manual of the transmitter for details.

6 Maintenance

The electrode with polymer electrolyte Xerolyt cannot be re-filled. The electrode must not be stored dry. In the event that the electrode dries out, it is necessary to recondition it before use by immersion in Friscolyt™ (Order no. 51 340 053) for at least 24 hours. The electrode should be stored with fitted watering cap containing 3 mol/l KCl or METTLER TOLEDO Friscolyt. After each working cycle, the electrode tip and the apertures (diaphragms) must be thoroughly cleaned using de-ionized water. Regeneration of the glass membrane of a pH electrode with HF (hydrofluoric acid) is not recommended, since HF will attack the Xerolyt reference system.

7 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

8 Warranty

In the event of manufacturing faults, 12 months warranty from date of production is granted.

Xerolyt and Friscolyt are trademarks of the METTLER TOLEDO Group.

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All rights reserved