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

显示屏	操作	备注
 <p>The screen displays 'CAL' in large letters, with 'MEDIUM AIR' below it. There are several icons in the top right corner, including a triangle with an exclamation mark and a triangle with a circle. A 'CAL' label is at the bottom left.</p>	<p>选择校准。传感器至于空气中，按 enter 键启动设备进入 Hold 状态</p>	<p>“介质 水”和“介质 空气”在校准时已设置。</p>
 <p>The screen displays '050.0' in large digits, with 'REL HUMIDITY' below it. There are icons in the top right corner. A 'CAL' label is at the bottom left.</p>	<p>按箭头键输入相对湿度按 enter 键继续</p>	<p>空气相对湿度深挖默认值：相对湿度 = 50%</p>
 <p>The screen displays '1000' in large digits, with 'BAR' to the right and 'PRESSURE' below it. There are icons in the top right corner. A 'CAL' label is at the bottom left.</p>	<p>按箭头键输入校准压力按 enter 键继续</p>	<p>默认值： 1.000 bar 单位 bar/kpa/PSI</p>
 <p>The screen displays '613' in large digits, with 'nA' to the right. Below it, '120 5' and '273' are shown. There are icons in the top right corner. A 'CAL' label is at the bottom left.</p>	<p>漂移检测：显示：传感器电流 (nA)、设置时间 (s)、温度 (° C/° F) 按 enter 键继续</p>	<p>漂移检测可能要持续几分钟。</p>
 <p>The screen displays '59.3' in large digits, with 'nA' to the right. Below it, 'ZERO' and '00 13 nA' are shown. There are icons in the top right corner. A 'CAL' label is at the bottom left.</p>	<p>显示校准数据（斜率和零点）。按 enter 键继续</p>	
 <p>The screen displays '20.93' in large digits, with 'MEAS' and 'REPE' below it. There are icons in the top right corner. A 'CAL' label is at the bottom left.</p>	<p>调整测量值中的测量显示（此处为：Vo1%）。设备仍处于 HOLD 状态：安装传感器并检查测量是否正确。MEAS 结束校准。REPEAT 允许重复。</p>	<p>校准结束后，输出端仍短时间保持在 HOLD 状态</p>

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