

Reference Manual for

METTLER TOLEDO


Standard Interface Command Set (MT-SICS)

MT-SICS 0 version 2.30

MT-SICS 1 version 2.20

MT-SICS 2 for Moisture Analyzer MJ33 version 2.30

MT-SICS 3 for Moisture Analyzer MJ33 version 1.30



METTLER TOLEDO

Contents

	Overview of all commands	4
1	Introduction	6
2	Basic information on data interchange	9
2.1	Command formats	9
2.2	Response formats	10
2.2.1	Format of the response with weight value	10
2.2.2	Format of the response without weight value	11
2.2.3	Error messages	12
2.2.4	Tips for the programmer	12
3	Commands and responses	13
3.1	Commands and responses MT-SICS level 0	13
3.2	Commands and responses MT-SICS level 1 for Moisture Analyzer MJ33	21
3.3	Commands and responses MT-SICS level 2 for Moisture Analyzer MJ33	23
3.4	Commands and responses MT-SICS level 3 for Moisture Analyzer MJ33	25
4	System Configuration (MJ33 – Computer)	40
5	What if...?	41

Overview of all commands

Commands and responses MT-SICS level 0

		Page
I0	Inquiry of all implemented MT-SICS commands	14
I1	Inquiry of MT-SICS level and MT-SICS versions	15
I2	Inquiry of instrument data	16
I3	Inquiry of SW version and type definition number	16
I4	Inquiry of serial number	16
I5	SW-Identification number	17
S	Send stable weight value	17
SI	Send weight value immediately	18
SIR	Send weight value immediately and repeat	18
Z	Zero	19
ZI	Zero immediately	20
@	Reset	20

Commands and responses MT-SICS level 1 for Moisture Analyzer MJ33

D	Display	21
DW	Weight display (Display show Weight)	22

Commands and responses MT-SICS level 2 for Moisture Analyzer MJ33

DAT	Date	23
PWR	Power On/Off	24
TIM	Time	24

Commands and responses MT-SICS level 3 for Moisture Analyzer MJ33

Control commands

HA01	Reset application / escape	25
HA02	Set factory settings	25
HA03	Switch keypad on/off	25
HA05	Start / end drying	26
HA06	Trigger audio signal	26
HA07	Report instrument status change	27
HA08	Request printer records	28

Status inquiries

HA20	Inquiry of instrument status	28
HA21	Inquiry of heating module position	29
HA22	Inquiry of last balance adjustment	29
HA23	Inquiry of last heating module adjustment	30
HA24	Inquiry of temperature	30
HA25	Inquiry of drying weights	31
HA26	Inquiry of drying data	32
HA27	Inquiry of drying result	33

Instrument settings

HA40	Inquiry / setting of language	34
HA40X	Inquiry / setting of menu parameters	35
HA61	Inquiry / setting of method parameters (part 1)	36
HA62	Inquiry / setting of method parameters (part 2)	38

1 Introduction

To enable you to integrate your Moisture Analyzer in your electronic data system in a simple manner and utilize their capabilities to the full, most Moisture Analyzer functions are available as appropriate commands via the data interface.

Standardization of the commands

All new METTLER TOLEDO balances and Moisture Analyzers support the standardized command set "METTLER TOLEDO Standard Interface Command Set" (MT-SICS), which is divided into 4 levels, depending on the functionality of the equipment:

- MT-SICS level 0 Command set for the simplest balance, e.g. weighing cell
- MT-SICS level 1 Extension of the command set for standard balances, i.e. balances without integrated applications
- MT-SICS level 2 Extension of the command set family
- MT-SICS level 3 Application-specific commands as independent command set, e.g. MT-SICS level 3 for Moisture Analyzers MJ33

A particular distinguishing feature of this concept is that the commands combined in MT-SICS level 0 and 1 are identical for all balances and instruments. Both the simplest weighing balance and a Moisture Analyzer recognize the commands of MT-SICS level 0.

What do the commands of MT-SICS level 0 and 1 offer?

You can use the commands of MT-SICS level 0 and 1 to perform the following operations of the Moisture Analyzer via the interface:

- request weighing results (raw data),
- zero the balance,
- identify MT-SICS implementation (version number),
- identify the Moisture Analyzer (serial number),
- reset the Moisture Analyzer,
- control the display.

The commands of MT-SICS level 2 and 3 for the Moisture Analyzer

All additional higher level functions for the Moisture Analyzer are collected in the commands of MT-SICS level 2 and 3.

When creating your software application, please note that the commands of MT-SICS level 3 apply to your application and can not be supported by every METTLER TOLEDO balance or instrument.

Most commands in this manual are identical for the MJ33 and the other Moisture Analyzers of METTLER TOLEDO.

Additional documentation on data interface

Settings of the interface such as baud rate, number of data bits, parity, handshake protocols and connector pin assignment are described in the operating instructions of the peripheral instrument or cable in question.

How the Moisture Analyzer operates

Your Moisture Analyzer is used to determine the moisture content of virtually any substance. The instrument operates according to the thermogravimetric principle: At the start of the measurement the Moisture Analyzer determines the weight of the sample, the sample is then rapidly heated with the built-in dryer unit and the moisture evaporates. During the drying, the instrument continuously records the weight of the sample and shows the decrease in the moisture. On completion of the drying, the final result of the moisture or dry substance content of your sample is displayed.

Version number of the MT-SICS

Each level of the MT-SICS has its own version number which can be requested with the command I1 from level 0.

This manual describes

MT-SICS level 0, version 2.30

MT-SICS level 1, version 2.20


MT-SICS level 2 for Moisture Analyzer version 2.30

MT-SICS level 3 for Moisture Analyzer version 1.30

You can use the command I1 via the interface to request the MT-SICS level and MT-SICS versions implemented on your moisture analyzer.

Please make sure that the versions implemented on your moisture analyzer agree with those listed above.

During drying the Moisture Analyzer passes through the following instrument statuses:



The following instrument statuses also exist:

- Entry status
- Startup
- Taring
- Weight adjustment
- Temperature adjustment
- Error status

Some functions or commands can be executed only in particular instrument statuses. In the instrument status "Drying", for example, the value in the display can not be overwritten. If a command can not be executed for this reason, the Moisture Analyzer sends an appropriate message.

You will find detailed information on the functions of the Moisture Analyzer in the operating instructions.

2 Basic information on data interchange

Each command received by the Moisture Analyzer via the data interface is acknowledged by a response of the Moisture Analyzer to the transmitter (e.g. computer). Commands and responses are data strings with a fixed format, and will be described in detail in chapter 3.

2.1 Command formats

Commands sent to the Moisture Analyzer comprise one or more characters of the ASCII character set. Here, the following must be noted:

- Commands must be entered in upper case letters (case sensitive!).
- The possible parameters of the command must be separated from one another and from the command name by a space (ASCII 32 dec., in this description represented as \square).
- Each command must be closed by CRLF (ASCII 13 dec., 10 dec.).

The characters $C_{R\text{L}F}$, which can be inputted using the Enter or Return key of most entry keypads, are not listed in this description, but it is essential they be included for communication with the instrument.

Example

Command to Moisture Analyzer which writes Hallo into the display:

$D\square\text{"Hallo"}$

The command terminator $C_{R\text{L}F}$ is not shown

Comment

The quotation marks "" must be inserted in the entry.

2.2 Response formats


All responses sent by the Moisture Analyzer to the transmitter (e.g. computer) to acknowledge the received command have one of the following formats:

- Response with weight value
- Response without weight value
- Error message

2.2.1 Format of the response with weight value

A general description of the response with weight value is the following.

ID	□	Status	□	WeightValue	□	Unit	C _R	L _F
----	---	--------	---	-------------	---	------	----------------	----------------



ID	Response identification
□	Space (ASCII 32 dec.)
Status	Interface status of the Moisture Analyzer, see description of the commands and responses
WeightValue	Weighing result; shown as number with 10 digits, incl. decimal point and sign – directly in front of the first digit if value negative. The weight value appears right-aligned. Preceding zeros are not shown with the exception of the zero to the left of the decimal point.
Unit	Weight unit
C_R	Carriage Return (ASCII 13 dec.)
L_F	Line Feed (ASCII 10 dec.)

Comment

C_RL_F will not be shown in this description.

Example

Response with stable weight value of 0.256 g:

S□S□□□□□□□□0.256□g

2.2.2 Format of the response without weight value

A general description of the response without weight value is the following

ID	␣	Status	␣	Parameters	C _R	L _F
----	---	--------	---	------------	----------------	----------------

1 – x characters
1 character

ID	Response identification
␣	Space (ASCII, 32 dec.)
Status	Interface status of the Moisture Analyzer, see description of the commands and responses
Parameters	Command-dependent response code
C _R	Carriage Return (ASCII 13 dec.)
L _F	Line Feed (ASCII 10 dec.)

Comment

C_RL_F will not be shown in this description.

Example

Response to D␣"HALLO" when HALLO appears unabridged in the display: D␣A.

2.2.3 Error messages

There are three different error messages. The identification always comprises two characters.

ID	C _R	L _F
----	----------------	----------------

ID	Error identification Possible error messages are	
	ES	Syntax error The Moisture Analyzer has not recognized the received command.
	EL	Logical error The Moisture Analyzer can not execute the received command.
C_R		Carriage Return (ASCII 13 dec.)
L_F		Line Feed (ASCII 10 dec.)

Comment

C_RL_F will not be shown in this description.

2.2.4 Tips for the programmer

Command and response

You can improve the dependability of your application software by having your program evaluate the response of the Moisture Analyzer to a command. The response is the acknowledgement that the Moisture Analyzer has received the command.

Reset

To be able to start from a definite condition when establishing the communication between Moisture Analyzer and system, you should send a reset command to the Moisture Analyzer.

When the Moisture Analyzer is switched off, meaningless characters may be received or sent.

Quotation marks ""

Quotation marks included in the command must always be entered.

3 Commands and responses

The Moisture Analyzer receives commands from the system (e.g. computer) and acknowledges the command with an appropriate response.

The following sections contain a detailed description of all commands of the command set in alphabetical order with the associated responses. Commands and responses are always closed with CRLF. These termination characters are not shown in the following description, but they must always be entered with commands or sent with responses.

3.1 Commands and responses MT-SICS level 0

The commands of MT-SICS level 0 are available with even the simplest balances which support the METTLER TOLEDO Standard Interface Command Set.

IO Inquiry of all implemented MT-SICS commands

Command	IO	Send list of all implemented MT-SICS commands
Response	IO␣B␣x1␣"1.Command"	x1 = ident of the MT-SICS level where the 1. Command belongs to.
	IO␣B␣x1␣"2.Command"	
	:	2nd (next) command implemented
	:	
	IO␣A␣x1␣"last Command"	Last command implemented. End of the list
	IO␣I	The list cannot be sent at present as another operation is taking place

Example

Command	IO	Send list of commands
Response	IO␣B␣0␣"IO"	Level 0 command "IO" implemented
	IO␣B␣0␣"I1"	Level 0 command "I1" implemented
	:	:
	IO␣B␣0␣"S"	Level 0 command "S" implemented
	:	:
	:	:
	IO␣A␣3␣"HA403"	Level 3 command "HA403" implemented

Comments

- The **IO** command lists all commands implemented in the present software.
- All level 0 commands are listed in alphabetical order before all commands of level 1 etc. This order corresponds to the order how the commands are described in this manual.

I1 Inquiry of MT-SICS level and MT-SICS versions

Command	I1	Inquiry of MT-SICS level and MT-SICS versions
Response	I1└┬┴┬"x1"└┬"x2"└┬"x3"└┬"x4"└┬"x5"	<p>x1 = 0 Balance with MT-SICS level 0 (simplest balance)</p> <p>x1 = 01 Balance with MT-SICS level 0 and 1 (standard balance)</p> <p>x1 = 012 Balance with MT-SICS level 0, 1 and 2 (standard balance with extensions)</p> <p>x1 = 03 Balance with MT-SICS level 0 and 3 (simplest balance with a special application)</p> <p>x1 = 013 Balance with MT-SICS level 0, 1 and 3 (standard balance with a special application)</p> <p>x1 = 0123 Balance with MT-SICS level 0, 1, 2 and 3 (standard balance with extensions and a special application)</p> <p>x1 = 3 Application device with MT-SICS level 3 (not necessarily a balance)</p> <p>x2 Version of the implemented MT-SICSO commands</p> <p>x3 Version of the implemented MT-SICS1 commands</p> <p>x4 Version of the implemented MT-SICS2 commands</p> <p>x5 Version of the implemented MT-SICS3 commands</p>

Example

Command	I1	Inquiry of MT-SICS level and versions
Response	I1└┬┴┬"3"└┬"2.30"└┬"2.20"└┬"2.30"└┬"1.30"	<p>3 Application device with MT-SICS level 3</p> <p>2.30 Level 0, version V2.30</p> <p>2.20 Level 1, version V2.20</p> <p>2.30 Level 2, version V2.30</p> <p>1.30 Level 3, version V1.30</p>

I2 Inquiry of instrument data

Command I2 Inquiry of instrument data

Response I2\A\ "text" Instrument data as "text"

Example

Command I2 Inquiry of instrument type

Responses I2\A\ "MJ33\Moisture-Analyzer\35.010\g"

I3 Inquiry of SW version and type definition number

Command I3 Inquiry of Moisture Analyzer SW version and type definition number

Response I3\A\ "TEXT" Moisture Analyzer SW version and type definition number as TEXT

Example

Command I3 Inquiry of SW version number(s) and type definition number

Response I3\A\ "1.00\4.10.5.93.43"
1.00 Software version
4.10.5.93.43 Type definition number

I4 Inquiry of serial number

Command I4 Inquiry of serial number

Response I4\A\ "text" Serial number as "text"

Example

Command I4 Inquiry of serial number

Response I4\A\ "0123456789"

I5 SW-Identification number

Command	I5	Inquiry of SW-Identification number.
Responses	I5␣A␣"x"	SW-Identification number as Text. x: SW-Identification number.
	I5␣I	Command understood, not executable at present.

Example

Command	I5	Inquiry of SW-Identification number.
Response	I5␣A␣"12345678A"	SW-Identification number with index.

Comment

- The SW-Identification number is unique for every Software.

S Send stable weight value

Command	s	Send the current stable weight value
Response	S␣S␣WeightValue␣Unit	Current stable weight value
	S␣I	Command not executable (Moisture Analyzer is currently executing another command)
	S␣+	Balance in overload range
	S␣-	Balance in underload range

Example

Command	s	Send a stable weight value
Response	S␣S␣␣␣␣␣␣␣␣␣1.000␣g	The current, stable weight value is 1.000 g

Comment

- Timeout approx. 30 s.

SI	Send weight value immediately
-----------	--------------------------------------

Command	SI	Send the current weight value, irrespective of balance stability
Response	S S S WeightValue Unit	Stable weight value
	S D D WeightValue Unit	Nonstable (dynamic) weight value
	S I	Command not executable (Moisture Analyzer is currently executing another command)
	S L +	Balance in overload range
	S L -	Balance in underload range

Example

Command	SI	Send current weight value
Response	S D D D D D D D D D 2.907 g	The current weight value is unstable (dynamic) and is 2.907 g

SIR	Send weight value immediately and repeat
------------	---

Command	SIR	Send the weight values repeatedly, irrespective of balance stability
Response	S S S WeightValue Unit	Stable weight value
	S D D WeightValue Unit	Nonstable (dynamic) weight value
	S I	Command not executable (Moisture Analyzer is currently executing another command)
	S L +	Balance in overload range
	S L -	Balance in underload range

ZI Zero immediately

Command	ZI	Zero immediately, i.e. stores immediately the current weight value, which can be stable or non stable (dynamic), as zero value.
Response	ZI┐S	Zero setting performed, stable weight value
	ZI┐D	Zero setting performed, non-stable (dynamic) weight value
	ZI┐I	Zero setting not performed (Moisture Analyzer is currently executing another command)
	ZI┐L	Command understood but not executable (e.g. certified version of balance)
	ZI┐+	Upper limit of zero setting range exceeded
	ZI┐-	Lower limit of zero setting range exceeded

Example 1

Command	ZI	Zero immediately
Response	ZI┐S	Zero setting performed, weight value was stable

Example 2

Command	ZI	Zero immediately
Response	ZI┐D	Zero setting performed, weight value was dynamic (non-stable)

@ Reset

Command	@	Resets the interface to the condition found after switching on, but without a zero setting being performed.
Response	I4┐AL"text"	Serial number of the Moisture Analyzer, the Moisture Analyzer is ready for operation.

Example

Command	@	
Response	I4┐AL"1114350697"	Moisture Analyzer is reset, its serial number is 1114350697.

Comment

All current commands are terminated on @ command

3.2 Commands and responses MT-SICS level 1 for Moisture Analyzer MJ33

The commands of MT-SICS level 1 are available with all standard balances which support the METTLER TOLEDO Standard Interface Command Set. With the MJ33 Moisture Analyzer, only the commands D and DW are supported.

D Display

Write into display

Command	<code>D_L"TEXT"</code>	Write TEXT into Moisture Analyzer display
Response	<code>D_LA</code>	TEXT appears unabridged left-aligned in the Moisture Analyzer display.
	<code>D_LR</code>	The end of the text appears in the Moisture Analyzer display, the start is cut off
	<code>D_LI</code>	Command not executable
	<code>D_LL</code>	Command understood, parameter wrong

Example

Command	<code>D_L"HALLO"</code>	Write HALLO into the Moisture Analyzer display
Response	<code>D_LA</code>	The full text HALLO appears in the Moisture Analyzer display

Clear display

Command	<code>D_L" "</code>	Clear Moisture Analyzer display
Response	<code>D_LA</code>	Display cleared

Comments

- A display command can be cleared with the DW or Reset command.
- This command can be executed only in the instrument statuses "basic mode", "ready for taring", "weighing-in" and "ready for start".

DW	Weight display (Display show Weight)
-----------	---

Command	DW	Switch display to weight mode
Response	DW┐A	Display shows the current weight value
	DW┐I	Command not executable

Comment

This command can be executed only in the instrument statuses "basic mode", "ready for taring", "weighing-in" and "ready for start".

3.3 Commands and responses MT-SICS level 2 for Moisture Analyzer MJ33

DAT	Date
------------	-------------

Inquiry of date

Command **DAT** Inquiry of current date of the Moisture Analyzer

Response **DAT␣A␣dd␣mm␣yyyy**
"dd␣mm␣yyyy" represents the date in the format
day␣month␣year

Set date

Command **DAT␣dd␣mm␣yyyy**
Set date in the format "dd␣mm␣yyyy"

Response **DAT␣A** Date has been set
DAT␣L Command not executed as the date format was not
correct
Inquiry of date of the Moisture Analyzer

Example

Command **DAT** Current date of the Moisture Analyzer is
2 April 2000

Response **DAT␣A␣02␣04␣2000**

Comments

- The set date is retained after the reset command "@".
- Admissible years: 1.1.1901 – 31.12.2099.

PWR Power On/Off

Command	PWR \downarrow x	Switch Moisture Analyzer On or Off x = 0 Set Moisture Analyzer to standby mode x = 1 Switch Moisture Analyzer on
Response	PWR \downarrow A	Moisture Analyzer has been switched off successfully
	PWR \downarrow A I4 \downarrow A \downarrow "text"	Moisture Analyzer with the serial number according to text has been switched on successfully (see also I4 command)
	PWR \downarrow L	Command understood, parameter wrong

Comments

- In the standby mode, the interface remains active; but all commands except **PWR**, **HA07**, **HA20** and **@** are answered with EL.
- On switching on, the Moisture Analyzer also sends the serial number (see also **I4** command).
- On switching off, all current commands are terminated.

TIM Time

Inquiry of time

Command	TIM	Send current time of the Moisture Analyzer
Response	TIM \downarrow A \downarrow hh \downarrow mm \downarrow ss	"hh \downarrow mm \downarrow ss" represents the time in the 24-hour format (hours/minutes/seconds)

Set time

Command	TIM \downarrow hh \downarrow mm \downarrow ss	Set time in 24-hour format (hours \downarrow minutes \downarrow seconds)
Response	TIM \downarrow A	Time has been set, clock running
	TIM \downarrow L	Command not executed as the time format is not correct (e.g. 22 \downarrow 67 \downarrow 25)

Example

Command	TIM	Inquiry of time
Response	TIM \downarrow A \downarrow 22 \downarrow 56 \downarrow 11	The current time of the Moisture Analyzer is 22 hours, 56 minutes and 11 seconds

3.4 Commands and responses MT-SICS level 3 for Moisture Analyzer MJ33

All Moisture Analyzer specified commands are combined in MT-SICS level 3 for Moisture Analyzers HR73, HG53, HB43-S and MJ33.

HA01 Reset application / escape

Command	HA01	Reset application / escape
Response	HA01 ␣A	Application reset

Comment

This command has the same effect as the Reset key, see operating instructions of the Moisture Analyzer. It terminates all current commands and activities.

HA02 Set factory settings

Command	HA02	Set factory setting of the menu and method parameters
Response	HA02 ␣A	Menu and method parameters set to factory setting

Comment

All menu parameters are reset to factory settings except RS interface settings and language. This command terminates a drying.

HA03 Switch keypad on/off

Command	HA03 ␣x	x = 0 Keypad of Moisture Analyzer switched off x = 1 Keypad of Moisture Analyzer switched on (factory setting)
Response	HA03 ␣A	Command executed
	HA03 ␣L	Command understood, parameter wrong

HA05 Start / end drying

Command	HA05 □ x	x = 0	End drying, possible only in instrument status "drying"
		x = 1	Start drying, possible only in instrument status "ready for start"
Response	HA05 □ A	Command executed	
	HA05 □ I	Command not executable as the Moisture Analyzer is not in the relevant instrument status	
	HA05 □ L	Command understood, parameter wrong	

Comment

To abort current drying use **HA01** (Reset application).

HA06 Trigger audio signal

Command	HA06	Trigger audio signal, e.g. at end of drying
Response	HA06 □ A	Command executed

HA07 Report instrument status change

Command	HA07 □ x1	Report each internal status change x1 = 0 Switch off x1 = 1 Switch on
Response	HA07 □ A	Command executed
	HA07 □ A □ x1	Status change (see HA20) x1 = 0 "Standby" x1 = 1 "Basic mode" x1 = 2 "Load pan and tare" x1 = 3 "Weighing-in" x1 = 4 "Ready for start" x1 = 5 "Drying" x1 = 6 "End of drying" x1 = 7 "Entry" x1 = 10 "Startup" x1 = 11 "Taring" x1 = 12 "Weight adjustment" x1 = 13 "Temperature adjustment" x1 = 101 "Error 1" x1 = 102 "Error 2" x1 = 10n "Error n"
	HA07 □ L	Parameter wrong (number, value range, ...)
	HA07 □ I	Response always available, hence not possible

Comments

- aborted with the **HA01** command
- see also **HA20** command
- also active in standby

HA08 Request printer records

Command	HA08 □ x1	Request printer records: x1 = 0 Do not send printer records x2 = 1 Send printer records
Response	HA08 □ A	Command executed
	HA08 □ L	Parameter wrong (number, value range, ...)
	HA08 □ I	Response always available, hence not possible

Comments

- The printer records use the 8-bit ASCII IBM table 4.
- Regardless of menu setting (see **HA403**).
- This setting is not stored.
- To reactivate the menu settings, use the "@" command.

HA20 Inquiry of instrument status

Command	HA20	Inquiry of instrument status		
Response	HA20 □ A □ x	x = 0	Status:	"Standby"
		x = 1	Status:	"Basic mode"
		x = 2	Status:	"Ready for taring"
		x = 3	Status:	"Weighing in"
		x = 4	Status:	"Ready for start"
		x = 5	Status:	"Drying"
		x = 6	Status:	"End of drying"
		x = 7	Status:	"Entry"
		x = 10	Status:	"Startup"
		x = 11	Status:	"Taring"
		x = 12	Status:	"Weight adjustment"
		x = 13	Status:	"Temperature adjustment"
		x = 101	Status:	"Error 1"
			
		x = 10n	Status	"Error n", see operating instructions of the Moisture Analyzer

Comments

- With the message **HA20**□**A**□**6** instrument status "End of drying", it is not apparent whether drying was ended correctly or terminated. This is possible only via the command **HA25** – Inquiry of drying weight.
- Also active in standby.

HA21 Inquiry of heating module position

Command	HA21	Inquiry of heating module position	
Response	HA21 ␣ A ␣ x	x = 0	Heating module closed
		x = 1	Heating module open

HA22 Inquiry of last balance adjustment

Command	HA22	Inquiry of last successful balance adjustment	
Response	HA22 ␣ A ␣ x1 ␣ x2 ␣ x3 ␣ x4 ␣ x5 ␣ x6		
		x1	Number of the successful adjustments
		x2	Day of the last successful adjustment
		x3	Month of the last successful adjustment
		x4	Year of the last successful adjustment
		x5	Hour of the last successful adjustment
		x6	Minute of the last successful adjustment

Example

Command	HA22	
Response	HA22 ␣ A ␣ 15 ␣ 02 ␣ 04 ␣ 2000 ␣ 09 ␣ 34	A total of 15 successful balance adjustments have been performed. The last took place on April 02, 2000 at 9.34.

Comments

- The time of the last successful balance adjustment is specified in the 24-hour format.
- Possible years are 1999 ... 2099.
- The counter for the balance adjustments runs to 65535.

HA23 Inquiry of last heating module adjustment

Command	HA23	Inquiry of the last successful heating module adjustments
Response	HA23└A└x1└x2└x3└x4└x5└x6	
	x1	Number of successful adjustments
	x2	Day of the last successful adjustment
	x3	Month of the last successful adjustment
	x4	Year of the last successful adjustment
	x5	Hour of the last successful adjustment
	x6	Minute of the last successful adjustment

Example

Command	HA23	
Response	HA23└A└15└02└04└2000└09└34	A total of 15 successful heating module adjustments have been performed. The last took place on April 02, 2000 at 9.34.

Comments

- The time of the last successful heating module adjustment is specified in the 24-hour format.
- Possible years are 1999...2099.
- The counter for the heating module adjustments runs to 65535.

HA24 Inquiry of temperature

Command	HA24	Inquiry of current temperature
Response	HA24└A└x	Current temperature in °C

Example

Command	HA24	Inquiry of current temperature
Response	HA24└A└105	The temperature is 105°C.

HA25 Inquiry of drying weights

Command **HA25** Inquiry of drying weight of the last or current drying

Response **HA25**␣A␣x1␣x2␣x3␣x4

x1 Drying status

x1 = 0 No drying exists

x1 = 1 Drying running

x1 = 2 Drying ended

x1 = 3 Drying terminated

x2 Wet weight in grams

x3 Current weight or dry weight in grams

x4 Drying time (seconds)

Example 1

Command **HA25** Inquiry of drying weights

Response **HA25**␣A␣2␣12.345␣7.890␣180

Drying has been ended regularly, wet weight
12.345 g, dry weight 7.890 g,
drying time 180 seconds

Example 2

Command **HA25** Inquiry of drying weights

Response **HA25**␣A␣0␣0.000␣0.000␣0

No drying exists, e.g. as the battery was discharged

Comment

Together with the command **HA07** – Report instrument status change – dryings can be shown in parallel on the host.

HA26 Inquiry of drying data

Command **HA26** \square **x1** Inquiry of drying data in configurable display mode
x1 = 0 currently set display mode
x1 = 1 Grams
x1 = 2 DC (dry content)
x1 = 3 MC (moisture content), (factory setting)
x1 = 4 AM (ATRO moisture content)

Response **HA26** \square **A** \square **x1** \square **x2** \square **x3** \square **x4** \square **x5** \square **x6**
x1 Drying status
x1 = 0 No drying exists
x1 = 1 Drying running
x1 = 2 Drying ended
x1 = 3 Drying terminated
x2 Display mode
x2 = 1 Grams
x2 = 2 DC (dry content)
x2 = 3 MC (moisture content), (factory setting)
x2 = 4 AM (ATRO moisture content)
x3 Wet weight in grams
x4 Current weight or dry weight in grams
x5 Actual result in requested display mode
x6 Drying time (seconds)
HA26 \square **L** Command understood, parameter wrong

Example 1

Command **HA26** \square **3** Inquiry of drying data

Response **HA26** \square **A** \square **2** \square **3** \square **4.762** \square **3.066** \square **35.61** \square **497**
Drying has been ended regularly, result requested in
% moisture content, wet weight 4.762 g, dry weight
3.066 g, 35.61% moisture content, drying ended
at 497 seconds

Example 2

Command **HA26**␣**2** Inquiry of drying data

Response **HA26**␣**A**␣**1**␣**2**␣**2.672**␣**2.467**␣**92.33**␣**143**
Drying is running, result requested in % dry content,
wet weight 2.672 g, dry weight 2.467 g, 92.33%
dry content, drying for 143 seconds in progress

Comment

- If a drying is inexistent (e.g. after a RAM LOST), the parameters x3...x6 are set to 0.
- If the measuring results exceed the tolerances for ATRO result display (L-999.99% AM or >999.99% AD) the selected results in x2 = 4 AM or x2 = 5 AD will automatically be transferred in x2 = 3 MC or x2 = 2 DC respectively.

HA27 Inquiry of drying result

Command **HA27**␣**x1** Inquiry of drying data in configurable display mode

x1 = 0 currently set display mode
x1 = 1 Grams
x1 = 2 DC (dry content)
x1 = 3 MC (moisture content), (factory setting)
x1 = 4 AM (ATRO moisture content)

Response **HA27**␣**A**␣**x1**␣**x2**
x1 **Drying status** (always 7 digit number)
x2 **Display mode** (g, %DC, %MC, %AM, %AD)

HA27␣**I** Response not available (drying in progress)

Example

Command **HA27**␣**3** Inquiry of drying result

Response **HA27**␣**A**␣**3**␣**-73.25**␣**MC**
Drying result -73.25% MC

Comment

If the measuring results exceed the tolerances for ATRO result display (L-999.99% AM or > 999.99% AD) the selected results in x1 = 4 AM or x1 = 5 AD will automatically be transferred in x1 = 3 MC or x1 = 2 DC respectively.

HA40 Inquiry / setting of language

Inquiry of language

Command **HA40** Inquiry of language currently set
Response **HA40□A□x** x Language (see below)

Setting language

Command **HA40□x** Set language
x = 1 English
x = 2 German
x = 3 French
x = 4 Italian
x = 5 Spanish
x = 6 Russian
x = 7 Japanese (Nihongo)
x = 8 Portuguese

Response **HA40□A** Language set
HA40□L Command understood, parameter wrong

HA40X Inquiry / setting of menu parameters

All commands on the inquiry / setting of menu parameters work similarly. As an example the inquiry for the startmode setting and the actual setting of the startmode is shown.

Inquiry of menu parameters

Command	HA401	Inquiry of startmode
Response	HA401┐┐A┐x1	
		x1 = 0 Startmode automatic (factory setting)
		x1 = 1 Startmode manual

Example

Command	HA401	Inquiry of startmode
Response	HA401┐┐A┐1	Startmode manual

Setting menu parameters

Command	HA401┐┐x1	Setting of startmode
		x1 = 0 Startmode automatic (factory setting)
		x1 = 1 Startmode manual
Response	HA401┐┐A	Startmode set
	HA401┐┐L	Command understood, parameter wrong

Example

Command	HA401┐┐0	Setting startmode to automatic
	HA401┐┐A	Startmode set

List of menu parameter inquiries / settings

HA401	Startmode	
		x1 = 0 Startmode automatic (factory setting)
		x1 = 1 Startmode manual
HA402	Protection against change in the settings (menu protection)	
		x1 = 0 Menu protection off, changes possible (factory setting)
		x1 = 1 Menu protection on, changes are not possible
HA403	Printer	
		x1 = 0 Printout off
		x1 = 1 Printout on (factory setting)

Comment

After the usage of HA08 the command HA403 has no effect on the printout but only on the setting of printout in the menu

HA61 Inquiry / setting of method parameters (part 1)

Inquiry of display mode, switch-off criteria and temperature profile

Command **HA61□x1** Inquiry of method parameters regarding display mode, switch-off criteria and temperature profile
x1 = 1 Inquiry of parameters of the current method

Response **HA61□A□x1□x2□...□x11**
Current setting of the method parameters
x1...x11 Represent the individual parameters (see below)

Example

Command **HA61□1** Inquiry of current setting of the method parameters

Response **HA61□A□1□3□1□0□1□105□0□0□0□0**
Get the method parameters: Display mode MC, switch-off Auto, temperature 105°C; all other parameters are set to the default values.

Setting display mode, switch-off criteria and temperature profile

Command **HA61□x1□x2□...□x11**
Set method parameters regarding display mode, switch-off criteria and temperature profile
x1...x11 represent the individual parameters (see below)

Response **HA61□A** Method parameter set
HA61□L Command understood, parameter wrong

Parameters	x1	Number of the method Always 1
	x2	Display mode
	x2 = 1	Grams
	x2 = 2	DC (dry content)
	x2 = 3	MC (moisture content), (factory setting)
	x2 = 4	AM (ATRO moisture content)
	x3	Switch-off criterion
	x3 = 1	Switch off Auto
	x3 = 2	Switch off via timer
	x4	Setting the timer in seconds Possible settings 60 - 5940 in steps of 60 s (only if x3 = 2, otherwise always 0)
	x5	Drying program Always 1
	x6	Set temperature in °C Possible settings 50...160 in steps of 1°C
	x7	Ramp time in seconds Always 0
	x8	Temperature of level 1 of step drying, in °C Always 0
	x9	Time of level 1 of step drying, in seconds Always 0
	x10	Temperature of level 2 of step drying, in °C Always 0
	x11	Time of level 2 of step drying, in seconds Always 0

Example

Command `HA61┘1┘3┘1┘0┘1┘105┘0┘0┘0┘0┘0`

Set method parameters: Display mode MC, switch-off Auto, set temperature 105°C; all other parameters are set to the default values.

Comments

- The parameters x7...x11 must always be 0.
- Setting the method parameters terminates a drying.

HA62 Inquiry / setting of method parameters (part 2)

Inquiry of target weight, print interval, method name and code

Command **HA62** \square **x1** Inquiry of method parameters regarding print interval
x1 = 1 Inquiry of parameters of the current method

Response **HA62** \square **A** \square **x1** \square **x2** \square ... \square **x5**
Current setting of the method parameters
x1 = 1 Set parameters of the current method
x1...x5 Represent the individual parameters (see below)

Example

Command **HA62** \square **1** Inquiry of current setting of the method parameters

Response **HA62** \square **A** \square **1** \square **0** \square **4** \square " " " "
Print interval 30 seconds.

Setting print interval

Command **HA62** \square **x1** \square **x2** \square ... \square **x5**
Set method parameters regarding print interval
x1...x5 represent the individual parameters (see below)

Response **HA62** \square **A** Method parameter set
HA62 \square **L** Command understood, parameter wrong

Parameters	x1	Number of the method x1 = 1 current method
	x2	Target weight in grams Possible settings: 0
	x3	Print interval x3 = 1 No print interval set, manual initiation of printout x3 = 4 Printout every 30 seconds x3 = 5 Printout every 60 seconds x3 = 9 Printout every 300 seconds
	x4	Method name " "
	x5	Code " "

Example

Command `HA62105" " " "`
Print interval set to 60 seconds


Response `HA62A` Method parameter set

Comment

Setting the method parameters terminates a drying.

4 System Configuration (MJ33 – Computer)

The MJ33 Moisture Analyzer is equipped with a 9 pin female RS232C Interface connector. It can be connected to a computer using a cable with order number 11101051 (9 pin) or order number 11101052 (25 pin).



The standard configuration of MJ33 is: 2400 baud, 7 bits, even parity, no handshake. These settings may be adjusted in the menu of MJ33. Additionally the printout setting in the menu should be "off". This avoids that printout strings are sent to the computer. Please refer to the Operating instruction.

5 What if...?

Tips from actual practice when the communication between the system (e.g. computer) and the Moisture Analyzer does not function.

Establishing the communication

Test whether the unidirectional operation is working:

Switch the Moisture Analyzer off with the "Off" key and then on again with the "On" key.

The Moisture Analyzer must now send the identification string **I4**, e.g. **I4**
A"0123456789". If "METTLER TOLEDO" is printed the communication functions properly. Change the printout setting in the menu to the "off" position.

If no identification string is received, check the following points.

Connection

For bidirectional communication, at least three connecting lines are needed:

- Data line from the Moisture Analyzer (TxD signal with RS232 interface).
- Data line to the Moisture Analyzer (RxD signal with RS232 interface).
- Signal ground line (SG with RS232 interface).

Make sure that all these connections are in order. Check the connector pin assignment of the connection cables.

Interface parameters

For the transmission to function properly, the settings of the following parameters must match at both the computer and the Moisture Analyzer:

- Baud rate (send/receive rate)
- Number of data bits
- Parity bit

Check the settings at both devices.

Handshake

For control of the transmission, in part separate connection lines are used (CTS/DTR). If these lines are missing or wrongly connected, the computer or Moisture Analyzer can not send or receive data.

Check whether the Moisture Analyzer is prevented from transmitting by handshake lines (CTS or DTR).

Set the parameter "handshake" for the Moisture Analyzer and the peripheral device to "No Handshake" or "none". The handshake lines now have no influence on the communication.

**To protect your METTLER TOLEDO product's future:
METTLER TOLEDO service assures the quality, measuring
accuracy and preservation of value of all METTLER TOLEDO
products for years to come.**

**Please send for full details about our attractive terms
of service.**

Thank you.



Subject to technical changes.