
Operating instructions
Bedienungsanleitung
Mode d'emploi
Instrucciones de manejo

METTLER TOLEDO

Microbalance filter set for filters up to \varnothing 70 mm
Mikrowaagen-Filterset für Filter bis \varnothing 70 mm
Kit pour le pesage de filtres jusqu'à \varnothing 70 mm sur microbalances
Conjunto para filtros de hasta 70 mm \varnothing para microbalanzas

Page 3

English

Seite 10

Deutsch


Page 17

Français



Página 24

Español

1 Overview and standard equipment


Part no.	Description
1	Coupling
2	Positioning nut
3	Sliding base (with brass pivot)
4	Sliding disk (3 pieces, 2 as spares)
5	Sliding bushing
6	Sealing cover
7	Draft shield (with serial number)
8	Filter pan
9	Knurled nut
10	Weighing pan fastener
11	Filter weighing pan ø 70mm or ø 47 mm
12	Glass cover
13	Tweezers for filter weighing


2 Assembly

Lift off the draft shield ring (A) and the weighing pan (B).

Lift off the intermediate shelf (C) and the sleeve (D).


Unscrew the mounting (E), remove the drip pan (F) and lift off the glass draft shield (G).

Remove the sealing cover (H).

Insert the coupling (1) and screw on the positioning nut (2).

Attach the sliding base (3) and sliding bushing (5). Attach the sliding disk (4) to the brass pivot of the sliding base.

Attach and secure the sealing cover (6) to the draft shield (7). Ensure that the small wheels (I) sit in the slots of the draft shield.

Insert the draft shield (7) and the sealing cover (6) and press them down. Turn the draft shield gently until it clicks into position (the sliding bushing (5) must sit in the slot on the underside of the draft shield).

Attach the filter pan (8) and screw on the knurled nut (9). Test the function of the draft shield using the door operating buttons on the balance (open to the left and to the right).

Insert the weighing pan fastener (10). **Important:** The shorter, slanted part (K) at the top of the weighing pan fastener must face the front of the balance!


Attach the large or small filter weighing pan (11).

Put on the glass cover (12).

3 Instructions for use


The filter set is supplied with specially designed **tweezers**. We recommend that only these tweezers be used for filter weighing. The upper tongue is shortened to ensure that the tweezers do not come into contact with the material on the filter. The lower, longer tongue ensures that the filter can be held steadily in the tweezers.

Insert the filter **sideways** through the opening in the draft shield.


We recommend that you wear **gloves**. By doing so you will avoid transferring body heat into the measuring compartment and thereby achieve an accurate result more quickly.

For high-precision measurements, we recommend that you use the **optional draft shield** and the **ionizer**. Both items are available as accessories (chapter 6.2).

Place the filter **directly in the center**, so that it touches both vertical pins on the back cantilever arms of the filter weighing pan.

After placing the filter, **close the draft shield** (operation of the draft shield door is described in the instruction manual for the balance).

You can use the integrated weight or an external weight for testing or **adjusting the balance**.

External cylindrical adjusting weights are available as accessories (chapter 6.2).

Place the external adjusting weight in the appropriate **grooves in the front cantilever arm of the filter weighing pan** (see figure).

You can find instructions on how to make adjustments in the instruction manual for the balance.

Important information about the draft shield:


Due to its rotary movement, the **sliding disk** will become rough over time. This can cause the draft shield to get jammed or stop turning completely. Should you experience this, replace the sliding disk (2 spare disks are included in the delivery).

4 Measuring methods for filter weighing

There are two main methods used for filter weighing, as described below.

Important: Before weighing, all filters must be sufficiently acclimatized in the measuring room!

4.1 Simple method


With the simple method, only one **test filter** (B) is used.

First, the new and uncontaminated filter is weighed (step 1) and then exposed to the source of contamination (steps 2 – 4).

Then, the contaminated filter is reweighed (step 5).

By comparing the two weighing results, the contamination can be calculated.

4.2 Method with reference filter

This method also allows for changes in ambient conditions and the potential contamination of the measuring room. Two filters are used for this, a **reference filter** (A) and a **test filter** (B).

The new reference filter and the test filter are weighed and the "**D1**" difference is calculated. The reference filter (A) is then stored in a filter case in the measuring room (steps 1 – 5).

The **test filter** (B) is exposed to the contamination (steps 6 – 7).

At the end of the contamination period, the reference filter (steps 8 – 10) and the test filter (steps 11 – 12) are reweighed and the "**D2**" difference is calculated.

The **effective contamination "D"** can be calculated by subtraction:

$$D = D2 - D1.$$

Note: The difference values "**D1**" and "**D2**" are easier to calculate if the reference filter is already loaded and the balance already tared. Then the reference filter is removed and the test filter is loaded. The balance now displays the difference in weight between the two filters.

5 Technical data


5.1 Specifications for the microbalances with the filter set

		MX5	UMX2	UMX5 Comparator
Readability		1 μ g	0.1 μ g	0.1 μ g
Maximum load	with standard weighing pan	510.0mg	210.0mg	510.0mg
	with filter set	510.0mg	210.0mg	510.0mg
Taring range	with standard weighing pan	0 ... 5100mg	0 ... 2100mg	0 ... 5100mg
	with filter set	0 ... 5100mg	0 ... 2100mg	0 ... 5100mg
Repeatability (s)	with standard weighing pan (guaranteed)	0.8 μ g	0.25 μ g	0.25 μ g
	with filter set (typical)	1 – 2 μ g	0.5 – 1.5 μ g	0.5 – 1.5 μ g
Stabilization time	with standard weighing pan	8s	10s	10s
	with filter set	20 – 40s	30 – 50s	30 – 50s
Standard weighing pan		\varnothing 27mm	\varnothing 16mm	\varnothing 16mm
Filter weighing pan		\varnothing 47/70mm	\varnothing 47/70mm	\varnothing 47/70mm

All further specifications remain unchanged and can be found in the operating instructions for the balance.

5.2 Dimensions

All dimensions in millimeters (mm).

6 Spare parts and accessories


6.1 Spare parts

Part no.	Description	Order number
11	Small filter weighing pan \varnothing 47 mm	11122137
11	Large filter weighing pan \varnothing 70 mm	11122138
	Pivot for the sliding base (brass) and knurled nut (part no. 9)	11122139
4	Sliding disk	11122132
6	Sealing cover	11122144
12	Glass cover	11122133
5	Sliding bushing	11122107
13	Tweezers	11122131

6.2 Accessories

Description		Order number
XP W5 complete draft shield		11116043
U ionizer small, electrode (Power supply must be ordered separately, see below)		11140161
Power supply for U ionizer electrode		11107766
5 mg E2 cylindrical adjustment weight with certificate		11122181
50mg cylindrical adjustment weight with certificate		11122183
100mg cylindrical adjustment weight with certificate		11122185
200mg cylindrical adjustment weight with certificate		11122187
500mg cylindrical adjustment weight with certificate		11122163
1.0g E2 cylindrical adjustment weight with certificate		11122165
1.5g cylindrical adjustment weight with certificate		11122189
2.0g E2 cylindrical adjustment weight with certificate		11122167
Tweezers for adjustment weights		

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**To protect your METTLER TOLEDO product's future:
METTLER TOLEDO Service assures the quality, measuring accuracy
and preservation of value of all METTLER TOLEDO products for years
to come.**

**Please send for full details about our attractive terms of service.
Thank you.**

P11781160

Subject to technical changes and to the availability of the
accessories supplied with the instruments.

© Mettler-Toledo AG 2007 11781160 Printed in Switzerland 0711/6.33

Mettler-Toledo AG, Laboratory & Weighing Technologies, CH-8606 Greifensee, Switzerland
Phone +41-44-944 22 11, Fax +41-44-944 30 60, Internet: <http://www.mt.com>