English

Operating Instructions

Microbalance Filter Kit for Filters up to \varnothing 70 mm

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


METTLER TOLEDO



1 Overview and Standard Equipment

Part no.	Description
1	Coupling
2	Positioning nut
3	Sliding base (with brass pivot and eccentric bolt)
4	Sliding bushing
5	Draft shield (with serial number)
6	Sealing cover
7	Filter pan
8	Knurled nut
9	Weighing pan holder
10	Adjustment pan
11	Filter weighing pan \varnothing 70mm or \varnothing 47 mm
12	Glass cover
13	Tweezers for filter weighing
14	Screwdriver "Torx" T6 (M2) and T10 (M3)


English


2 Assembly

Lift off the draft shield ring (A) and the weighing pan (B).


Lift off the intermediate shelf (C) and the sleeve (D).

Unscrew the mounting (E), remove the drip pan (F) and lift off the glass draft shield (G).


Remove the sealing cover (H).

Insert the coupling (1) and screw on the positioning nut (2).

Attach the sliding base (3) and sliding bushing (4).

Remove the plastic cap (I) from the eccentric bolt (K). Loosen the eccentric bolt and turn it clockwise until the sliding base (3) is fixed. Fasten the eccentric bolt and attach the plastic cap (I).


Attach and secure the sealing cover (6) to the draft shield (5). Ensure that the small wheels (L) sit in the slots of the draft shield.


Insert the draft shield (5) and the sealing cover (6) and press them down. Turn the draft shield gently until it clicks into position (the sliding bushing (4) must sit in the slot on the underside of the draft shield).

Attach the filter pan (7) and screw on the knurled nut (8).
Test the function of the draft shield using the door operating buttons on the balance (open to the right).


Insert the weighing pan holder (9).
Important: The shorter, slanted part (M) at the top of the weighing pan fastener must face the front of the balance!

Attach the large or small filter weighing pan (11). Put on the glass cover (12).

3 Changing the Side of the Door Opening


The draft shield can be adapted from right side to left side operation.

1 Remove the glass cover (12).
Remove the large or small filter weighing pan (11).

2 Unscrew the knurled nut (8). Remove the filter pan (7).

3 Unscrew the limit stop (A) in the right position. Change it to the left position.

4 Screw on the limit stop (A) in the left position.


5 Attach the filter pan (7) and screw on the knurled nut (8).
Test the function of the draft shield using the door operating buttons on the balance (open to the left).

6 Attach the large or small filter weighing pan (11).
Put on the glass cover (12).

4 Instructions for Use


The filter set is supplied with specially designed **tweezers**. We recommend that only these tweezers be used for filter weighing. The upper tongue is shortened to ensure that the tweezers do not come into contact with the material on the filter. The lower, longer tongue ensures that the filter can be held steadily in the tweezers.

Insert the filter **sideways** through the opening in the draft shield.

We recommend that you wear **gloves**. By doing so you will avoid transferring body heat into the measuring compartment and thereby achieve an accurate result more quickly.


For high-precision measurements, we recommend that you use the **optional draft shield** and the **ionizer**. Both items are available as accessories (chapter 6.2).

Place the filter **directly in the center**, so that it touches both vertical pins on the back cantilever arms of the filter weighing pan.

After placing the filter, **close the draft shield** (operation of the draft shield door is described in the instruction manual for the balance).

You can use the integrated weight or an external weight for testing or **adjusting the balance**. **External adjusting weights** are available as accessories (chapter 7.2).

Place the **cylindrical weight** in the appropriate **grooves in the front cantilever arm of the filter weighing pan** (11) (see figure).

Place the **standard weight** on the adjustment pan (10) (see figure).


You can find instructions on how to make adjustments in the instruction manual for the balance.

5 Measuring Methods for Filter Weighing

There are two main methods used for filter weighing, as described below.

Important: Before weighing, all filters must be sufficiently acclimatized in the measuring room!

5.1 Simple method


With the simple method, only one **test filter (B)** is used.

First, the new and uncontaminated filter is weighed (step 1) and then exposed to the source of contamination (steps 2 – 4).

Then, the contaminated filter is reweighed (step 5).

By comparing the two weighing results, the contamination can be calculated.

5.2 Method with reference filter

This method also allows for changes in ambient conditions and the potential contamination of the measuring room. Two filters are used for this, a **reference filter (A)** and a **test filter (B)**.

The new reference filter and the test filter are weighed and the “**D1**” difference is calculated. The reference filter (A) is then stored in a filter case in the measuring room (steps 1 – 5).

The **test filter (B)** is exposed to the contamination (steps 6 – 7).

At the end of the contamination period, the reference filter (steps 8 – 10) and the test filter (steps 11 – 12) are reweighed and the “**D2**” difference is calculated.

The **effective contamination “D”** can be calculated by subtraction:

$$D = D2 - D1.$$

Note: The difference values “D1” and “D2” are easier to calculate if the reference filter is already loaded and the balance already tared. Then the reference filter is removed and the test filter is loaded. The balance now displays the difference in weight between the two filters.

6 Technical Data

6.1 Specifications for the Microbalances with the Filter Set


	MX5	UMX2	UMX5 Comparator
Readability	1 μ g	0.1 μ g	0.1 μ g
Maximum load with standard weighing pan with filter set	5100 mg 5100 mg	2100 mg 2100 mg	5100 mg 5100 mg
Taring range with standard weighing pan with filter set	0 ... 5100 mg 0 ... 5100 mg	0 ... 2100 mg 0 ... 2100 mg	0 ... 5100 mg 0 ... 5100 mg
Repeatability (s) with standard weighing pan guaranteed (test load) with filter set (typical)	0.8 μ g (5 g) 1 – 2 μ g	0.25 μ g (2 g) 0.5 – 1.5 μ g	0.4 μ g (5 g) 0.5 – 1.5 μ g
Settling time with standard weighing pan with filter set	8 s 20 – 40 s	10 s 30 – 50 s	10 s 30 – 50 s
Standard weighing pan	\varnothing 27 mm	\varnothing 16 mm	\varnothing 16 mm
Filter weighing pan	\varnothing 47/ 70 mm	\varnothing 47/ 70 mm	\varnothing 47/ 70 mm

	XP6	XP2U	XP6 U
Readability	1 μ g	0.1 μ g	0.1 μ g
Maximum load with standard weighing pan with filter set	6100 mg 6100 mg	2100 mg 2100 mg	6100 mg 6100 mg
Taring range with standard weighing pan with filter set	0 ... 6100 mg 0 ... 6100 mg	0 ... 2100 mg 0 ... 2100 mg	0 ... 6100 mg 0 ... 6100 mg
Repeatability (s) with standard weighing pan guaranteed (test load) with filter set (typical)	0.8 μ g (6 g) 1 – 2 μ g	0.25 μ g (2 g) 0.5 – 1.5 μ g	0.4 μ g (6 g) 0.5 – 1.5 μ g
Settling time with standard weighing pan with filter set	8 s 20 – 40 s	10 s 20 – 40 s	10 s 20 – 40 s
Standard weighing pan	\varnothing 27 mm	\varnothing 16 mm	\varnothing 16 mm
Filter weighing pan	\varnothing 47/ 70 mm	\varnothing 47/ 70 mm	\varnothing 47/ 70 mm

All further specifications remain unchanged and can be found in the operating instructions for the balance.

6.2 Dimensions

All dimensions in millimeters (mm).



7 Spare Parts and Accessories


7.1 Spare Parts

Part no.	Description	Order number
11	Small filter weighing pan \varnothing 47 mm	30006009
11	Large filter weighing pan \varnothing 70mm	30006020
10	Adjustment pan	30006024
	Pivot for the sliding base (brass) and knurled nut (part no. 9)	30006022
4	Sliding bushing	11122107
6	Sealing cover	30006023
12	Glass cover	11122133
13	Tweezers	11122131

7.2 Accessories

	Description	Order number
	XP W5 complete draft shield	11116043
	U ionizer small, electrode (Power supply must be ordered separately, see below)	11140161
	Power supply for U ionizer electrode	11107766
	5 mg E2 cylindrical adjustment weight with certificate	11122181
	50 mg E2 cylindrical adjustment weight with certificate	11122183
	100 mg E2 cylindrical adjustment weight with certificate	11122185
	200 mg E2 cylindrical adjustment weight with certificate	11122187
	500 mg E2 cylindrical adjustment weight with certificate	11122163
	1.0 g E2 cylindrical adjustment weight with certificate	11122165
	2.0 g E2 cylindrical adjustment weight with certificate	11122167
	10 mg E2 standard adjustment weight with certificate (wire weight)	00158337
	20 mg E2 standard adjustment weight with certificate (wire weight)	00158347
	50 mg E2 standard adjustment weight with certificate (wire weight)	00158357
	100 mg E2 standard adjustment weight with certificate (wire weight)	00158367
	200 mg E2 standard adjustment weight with certificate (wire weight)	00158377
	1 g E2 standard adjustment weight with certificate (cylindrical weight with knob)	00158397
	2 g E2 standard adjustment weight with certificate (cylindrical weight with knob)	00158407
	5 g E2 standard adjustment weight with certificate (cylindrical weight with knob)	00158417
	Tweezers for adjustment weights	00015900

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© Mettler-Toledo AG 10/2010
Mettler-Toledo AG Laboratory & Weighing Technologies