Operating instructions

LP16-M Thermometer kit 13701

METTLER TOLEDO

Preparing the LP16-M Infrared Dryer

Check if the protective screen of your LP16-M dryer has a keyhole-like perforation according to **Fig. 1**.

If not, an appropriate modification must be made by METTLER TOLEDO-Service.

Important: Note that temperature measurements cannot be performed during a drying cycle

Temperature measurement is intended for the two following applications:

- Measuring the effective temperature of the sample. (The temperature of the sample is a function of color and composition of the substance. The instrument has been factory set, so that a sample of common salt will reach a temperature of 100°C in 20 minutes.)
- Temperature measurement for calibrating your instrument with a sample of common salt (sodium chloride).


See overleaf 

P703131

Subject to technical changes and to the availability of the accessoires supplied with the instruments.


Preparation for temperature measurement

Fill in approx. 6 mm of sample substance in sample pan (the substance must cover the entire bottom surface). Following, slide the plastic clamp over the thermometer, according to **Fig. 2**.

Caution: Do not laterally compress clamp, instead press against base surface and slide into position (slit facing front, so that the scale remains readable).

Now, place pan containing the sample into LP16-M dryer and close the hood. Hold the thermometer by the clamp and carefully insert it into the hole in the screen until contacting the sample. Note the position on the thermometer graduation matching the

screen surface. Retrieve the thermometer and move the clamp until its lower edge reaches the same position. Then, open the hood and smooth over the sample surface. Close the hood and reinsert the thermometer. Now, an impression of a certain depth is formed within the sample (the thermometer should penetrate approx. half way into the sample layer - note also **Fig. 3**). For fine adjustment the clamp should be moved up or down by a few graduations, followed by smoothing the sample surface and rechecking the thermometer penetration. This must be repeated until the desired value is reached.

Caution: Before opening the hood always remove the thermometer (breaking the thermometer causes the release of mercury).

Temperature measurement procedures

The entire thermometer tip must be in close contact with the sample. The balance must indicate a value in excess of 100 digits. (In case of overload or if weighing pan is empty the measuring cycle cannot be started).

Select the desired temperature, set the timer to 0 minutes (drying process continues until -Stop- is pressed) and press the start key. When the temperature ceases to rise at the thermometer, it can be read off.

For instrument calibration a sample of common salt (sodium chloride) is needed. The potentiometer at the right side of the LP16-M dryer (small cover) permits temperature adjustment. Clockwise rotation causes temperature increase; adjustable within $\pm 10^{\circ}\text{C}$.

Summary of important points


- Measurements only possible in a sample (not in plain air)
- The sample must cover the entire bottom surface of the pan.
- Temperature measurement cannot be performed during a drying cycle.
- The thermometer must penetrate half way into the sample layer.
- The temperature is a function of substance.
- The entire thermometer tip must be in close contact with the sample.
- Before opening the hood of the LP16-M dryer always remove the thermometer.
- An LP16-M applications brochure, containing instructions for moisture determination, is available on request.

Thermometer kit includes:

- Thermometer
- Plastic clamp
- Operating instructions

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