METTLER TOLEDO

InPro 7200 Series

**Inductive Conductivity Sensors
Induktive Leitfähigkeitsmesszellen
Capteurs de conductivité inductifs**

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InPro 7200 Series

Instruction manual

Subject to technical changes without prior notice.

© It is strictly forbidden to reprint this instruction manual or any parts thereof. Without the written permission of Mettler-Toledo GmbH, Process Analytics, Industrie Nord, 8902 Urdorf, Switzerland, no section or excerpt whatsoever may be reproduced or with the assistance of electronic systems be edited, duplicated or distributed, in particular in the form of photocopies, photographs, magnetic media or other recording methods. All rights reserved, especially the right of duplication and translation as well as in regard to patent and registration rights.

Mettler-Toledo GmbH, 8606 Greifensee, Switzerland

Contents

1	Introduction.....	5
1.1	General.....	5
1.2	Safety Instructions.....	5
2	Sensor	6
2.1	Sensor Certification Specifications	6
2.2	Sensor Identification	6
2.3	Standard Specifications	7
2.4	Sensor Specifications	7
3	Sensor Applications	11
4	Sensor Installation.....	12
4.1	Installation Notes	12
4.2	Flange Installation	13
4.3	Bushing Installation	14
4.4	Sanitary Installation.....	14
4.5	Immersion Installation	15
4.6	In-Line Installation	15
5	Junction Box Installation.....	16
5.1	Junction Box Wiring.....	16
6	Degaussing and Sensor Resistance Check.....	17
6.1	Degaussing.....	17
6.2	Sensor Resistance Check	17
7	Disposal.....	18
8	Dimension Drawings.....	18
9	Spare Parts	19

1 Introduction

1.1 General

The inductive Sensors InPro 7200 Series are intended to be used with Transmitter Cond Ind 7100 e, Cond \pm 7100/2H, Cond \pm 7100/2(X)H, or M 700 (module Cond Ind 7700) to measure conductivity and concentration of liquids. Sensors are available in a number of materials suitable for a variety of conditions.

The sensors measure an induced current in a loop of solution. Two toroidally wound coils are encapsulated in close proximity within the sensor which is immersed in the solution. An ac signal, applied to one toroidal coil, induces a current in the second coil which is directly proportional to the conductance of the solution.

1.2 Safety Instructions

Dangers

An electrical shock hazard may be present under sensor fault conditions when sensor is located in process fluids at electrical potentials above 30 V rms or 60 V dc. In order **to prevent electrical shock hazard, remove sensor from process or de-energize power to process prior to performing maintenance on the Transmitter Cond Ind 7100 e, Cond \pm 7100/2H, Cond \pm 7100/2(X)H, or M 700.**

Certain units described in this instruction manual contain electrical shock hazard potential. **Maintenance should be performed by qualified personnel only. During maintenance, remove power to avoid severe injury or death.**

When performing maintenance, wear appropriate, **protective clothing including safety goggles**. Escaping chemicals under pressure can cause severe injury including blindness.

Warnings

Use care when connecting and disconnecting high pressure service connections. Use **proper gloves and follow the recommended procedures to avoid severe injury to personnel or damage to equipment.**

When processing hazardous liquids, follow the recommended procedures. Failure to do so could result in injury to personnel and damage to equipment.

Use only METTLER TOLEDO recommended replacement parts. Substitution parts could result in damage to equipment, damage to the process, and/or injury to personnel.

Be careful to avoid touching exposed circuits and components. Potential shock hazards are present.

2 Sensor

2.1 Sensor Certification Specifications

Testing Laboratory, Types of Protection and Area Classification	Conditions of Certification	Suitable Sensor Types
FM Certification I.S. Class 1 Div 1 Group A – G T5 at maximum ambient temperature of 85 °C	Transmitter Cond \bar{I} 7100/2(X)H	InPro 7200 InPro 7201 InPro 7202
ATEX intrinsically safe for Ⓔ II 1 GD EEx ia IIC T3...T6 According to EC-type test certificate KEMA 03 ATEX 1462 X	Transmitter Cond \bar{I} 7100/2(X)H and M 700X	InPro 7200 InPro 7201 InPro 7202

Table 1: Sensor Certification Specifications

Note: These sensors have been designed to meet the electrical safety descriptions noted in the table above. For detailed information contact your METTLER TOLEDO representative.

2.2 Sensor Identification

A sensor can be identified by the model code located on both ends of the cable. Refer to «Sensor Applications» section for sensor application information.


Figure 1: Sensor Identification

Date code

The date code identifies the year and week of manufacture (see Figure 1). In the example 2B0312 XXX serial 03 identifies the year of manufacture as 2003, and 12, the week of manufacture in that year. The additional numbers are the serial number.

Entity parameters as follows:

$$\begin{aligned}
 U_i &= 7.5 \text{ V} & I_i &= 70 \text{ mA} \\
 P_i &= 0.4 \text{ W} & L_i &= 7.8 \text{ mH}
 \end{aligned}$$

The effective internal capacitance C_i is negligibly small.

2.3 Standard Specifications

Wetted Parts

See Table 2 "Sensor Specifications".

Cable

Integral, 6 m (20 ft) multi-screened (multi-shielded) cable with the following jacket material:

InPro 7200; irradiated polyolefin jacket

InPro 7201; PTFE jacket

InPro 7202; irradiated polyolefin jacket

Mounting

All sensor extensions contain an O-ring and truncated 3/4"-14 NPSM pipe thread.

Insertion Mounting: Sensors are used with threaded bushings or flanges that form a process seal against sensor O-ring.

Note: When mounted in-line, the sensors must be installed (centered) in a minimum pipe size to avoid pipe wall effects (DN 80 or 3").

2.4 Sensor Specifications

	InPro 7200	InPro 7201	InPro 7202
Measurement	0.01...2000mS/cm	0.01...2000mS/cm	0.01...2000mS/cm
Cell factor nominal ¹⁾	2.15	2.15	2.15
Transfer ratio ²⁾	48.36	48.36	48.36
Material Sensor O-rings*	Glass-filled PEEK EPR	Glass-filled PEEK EPR	Polypropylene EPR
Process temperature	-5...+120 °C	-5...+200 °C	-5...+80 °C
Max. process press. ³⁾	0...17.5 bar	0...17.5 bar	0...17.5 bar
Temperature probe	NTC 100 kΩ	Pt 100	Pt 100
Explosion protection	FM I.S. Class 1, Div 1 Group A – G T5 at max. ambient temperature of 85 °C ATEX intrinsically safe for II 1 GD EEx ia IIC T3...T6 According to EC-type test certificate KEMA 03 ATEX 1462 X	FM I.S. Class 1, Div 1 Group A – G T5 at max. ambient temperature of 85 °C ATEX intrinsically safe for II 1 GD EEx ia IIC T3...T6 According to EC-type test certificate KEMA 03 ATEX 1462 X	FM I.S. Class 1, Div 1 Group A – G T5 at max. ambient temperature of 85 °C ATEX intrinsically safe for II 1 GD EEx ia IIC T3...T6 According to EC-type test certificate KEMA 03 ATEX 1462 X

Table 2: Sensor Specifications

* other O-ring materials on request (see chapter 9)

- 1) This value may differ for each individual sensor and depends on the installation conditions. The exact value must be determined with a calibration.
- 2) This typical value may differ for each individual sensor. The exact value must be determined with a calibration.
- 3) The combination of high pressure, high temperature and/or aggressive process medium reduces the sensor life span.

Special conditions

The maximum permissible process temperature should be obtained from the following table according to the temperature class:

Temperature class	Maximum process temperature
T6	85 °C
T5	100 °C
T4	135 °C
T3	200 °C

METTLER TOLEDO offers a variety of fittings and process adapters for a simple installation of the inductive sensors into the process.

The measurement range is defined for sensors used with Transmitter Cond Ind 7100 e, CondI 7100/2H, CondI 7100/2(X)H or M 700.

ATEX certificate

(1) **EC-TYPE EXAMINATION CERTIFICATE**

(2) Equipment or protective system intended for use in potentially explosive atmospheres - Directive 94/9/EC

(3) EC-Type Examination Certificate Number: **KEMA 03ATEX1462 X**

(4) Equipment or protective system: **InPro 7200 sensor series**

(5) Manufacturer: **Mettler-Toledo GmbH**

(6) Address: **Im Hackacker 15, 8902 Urdorf, Switzerland**

(7) This equipment or protective system and any acceptable variation thereto is specified in the schedule to this certificate and the documents therein referred to.

(8) KEMA Quality B.V., notified body number 0344 in accordance with Article 9 of the Council Directive 94/9/EC of 23 March 1994, certifies that this equipment or protective system has been found to comply with the Essential Health and Safety Requirements relating to the design and construction of equipment and protective systems intended for use in potentially explosive atmospheres given in Annex II to the directive.

The examination and test results are recorded in confidential report no. 2037022.

(9) Compliance with the Essential Health and Safety Requirements has been assured by compliance with:

EN 50014 : 1997 + A1, A2	EN 50020 : 2002
EN 50281-1-1:1995 + A1	EN 50284 : 1999

(10) If the sign "X" is placed after the certificate number, it indicates that the equipment or protective system is subject to special conditions for safe use specified in the schedule to this certificate.

(11) This EC-Type Examination Certificate relates only to the design, examination and tests of the specified equipment or protective system according to the Directive 94/9/EC. Further requirements of the directive apply to the manufacturing process and supply of this equipment or protective system. These are not covered by this certificate.

(12) The marking of the equipment or protective system shall include the following:

II 1 GD EEx ia IIC T3 ... T6 T 110 °C ... T 260 °C

Amhem, 15 July 2004
KEMA Quality B.V.

C.G. van Es
Certification Manager

* This Certificate may only be reproduced in its entirety and without any change

KEMA Quality B.V.
Ulrichsweg 310, 6812 AR Amhem, The Netherlands
P.O. Box 5185, 6802 ED Amhem, The Netherlands
Telephone +31 26 3 56 20 68, Telex +31 26 3 52 68 00

ACCREDITED BY THE
DUTCH COUNCIL FOR
ACCREDITATION

Page 1/3

(13) **SCHEDULE**

(14) **to EC-Type Examination Certificate KEMA 03ATEX1462 X**

(15) **Description**

The InPro 7200 sensor series are for connection to electrochemical transmitters.

Ambient temperature range -20 °C ... +60 °C.

The relation between ambient temperature, process temperature and temperature class is as shown in the table below:

Maximum ambient temperature (°C)	Maximum process temperature (°C)	Temperature class
60	85	T6
60	100	T5
60	135	T4
60	200	T3

The relation between ambient temperature, process temperature and surface temperature "T" is as shown in the table below:

Maximum ambient temperature (°C)	Maximum process temperature (°C)	Surface temperature "T" (°C)
60	50	110
60	60	120
60	80	140
60	85	145
60	125	185
60	135	195
60	175	235
60	200	260

Electrical data

Sensor circuit in type of explosion protection intrinsic safety EEx ia IIC, only for connection to a certified intrinsically safe circuit, with the following maximum values:

U _i	=	7,5	V
I _i	=	70	mA
P _i	=	0,4	W
C _i	=	0	nF
L _i	=	7,8	mH

Page 2/3

(13) **SCHEDULE**

(14) **to EC-Type Examination Certificate KEMA 03ATEX1462 X**

(16) **Report**
KEMA No. 2037022.

(17) **Special conditions for safe use**

1. For applications in explosive atmospheres caused by air/dust mixtures and where category 1D apparatus is required, the surface temperature "T" is determined for a dust layer with a thickness of 5 mm maximum.
2. For applications in explosive atmospheres caused by gases, vapours or hazes or air/dust mixtures and where category 1 apparatus is required, precautions shall be taken such that electrostatic charges on the non metallic parts of the sensors are avoided.
3. For ambient temperature range and electrical data see (15).

(18) **Essential Health and Safety Requirements**
Covered by the standards listed at (9).

Compliance with EN 50281-1-1 : 1998 has been assured, taking into account draft IEC 61241-0 : 2002 and draft IEC 61241-11 : 2002 as a guide.

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Mettler-Toledo GmbH, Process Analytics

 Waldemar Rauch
 General Manager PO Urdorf

 Christian Zwicky
 Head of Marketing

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3 Sensor Applications

Notes:

- PEEK is a thermoplastic material with excellent strength and chemical resistance properties over a wide range of process temperatures and pressures. METTLER TOLEDO recommends using PEEK preferentially for all applications that qualify.
- PEEK material displays excellent chemical resistance to most aqueous solutions of acids, bases, and salts. It is also excellent for organic solvents such as toluene, ethyl acetate, acetone, gasoline, and carbon tetrachloride. It is not recommended for sulfuric or nitric acid solutions above 70 %, nor is it recommended for Oleum applications.

Sensor Designation	Application
InPro 7200	This sensor is suitable for the majority of all electrodeless conductivity applications. Its compact size enables it to be mounted in a multitude of methods, including insertion (flange, bushing), and immersion. Typical applications include salinity and brine measurements, steel pickling, scrubbing towers, ion exchange regeneration, plating baths, rinse water, caustic metal cleaning, and textile measurements in scouring, mercerizing, and carbonizing baths.
InPro 7201	This sensor is identical in physical size and appearance to the InPro 7200 sensor and may be applied to any of the applications identified which have intermittent or continuously high temperatures. Additional high temperature applications for which this sensor is suited include alumina-to-caustic ratio, boiler blowdown, and Cleaning-In-Place (CIP) measurements in food and related industries.
InPro 7202	This general purpose, sensor may be used for most routine applications involving low (<20 %) concentrations of inorganic acids (hydrochloric, nitric, sulfuric, etc.), bases (caustic, calcium hydroxide, etc.), and salts (sodium chloride, calcium chloride, sodium sulfate, etc.). When doubtful about the effect of high levels of chemicals, temperatures, or abrasion, the InPro 7200 PEEK sensor is a preferred alternative.

Table 3: Suitable Applications

4 Sensor Installation

4.1 Installation Notes

Proper installation of the sensor is important for efficient and accurate operation. For all applications and sensor configurations, mounting arrangements must be located so that:

- Sample at the sensing area is representative of the solution.
- Solution circulates actively and continuously past the sensing area.
- Position and orientation of the sensor does not trap air bubbles within sensing area.
- Deposits of sediment or other foreign material do not accumulate within the sensing area.

If cable is installed in a conduit (recommended), either flexible conduit should be used or some other provision made for removal of sensor from the process.

Also lubricate all O-rings with a thin film of suitable grease.

Notes:

- There are two grooves on the thread diameter of the sensor. The groove closest to the cable is **not** an O-ring groove and should remain unoccupied. The next groove, as shown in Figure 2, **is** for the O-ring.
- Ensure that the sensor received has the O-ring in the proper position.
- Ensure that future O-ring replacement is in the proper position.

Warning: Incorrect O-ring position as previously indicated could result in process leakage, causing injury to personnel.

- There are **two notches** at the cable cap that line up with the bore of the sensor doughnut. These notches are an aid in positioning the sensor in a pipe or vessel.


Figure 2: Sensor overview

4.2 Flange Installation

METTLER TOLEDO flanges are used for permanent installations in pipelines and tanks. Flanges are used with electrodeless conductivity sensors in systems using 2-inch and larger process piping. Flanges come complete with locknut and a set of spacers. Refer to Table 4 for flange specifications, and to Figure 3 for a typical flange installation.

Note: One or two spacers are used depending on flange thickness. The O-ring must be positioned correctly for proper spacing.

Flange size	Part No (Set)	Material*
ANSI 2"	52401933	1.4435
ANSI 2 1/2"	52401934	1.4435
ANSI 3"	52401938	1.4435
ANSI 4"	52401939	1.4435
DN 50 / PN 16	52401940	1.4435
DN 100 / PN 16	52401941	1.4435

Table 4: Flange Specifications. *Flanges made of stainless steel are supplied as standard with Inspection Certificate 3.1B.

Typical Flange Installation


Figure 3: Flange installation

4.3 Bushing Installation

METTLER TOLEDO bushings (locknuts supplied) are used for installations in pipelines and tanks. Refer to Table 5 for bushing specifications and to Figure 4 for a typical bushing installation.

Typical Bushing Installation


Figure 4: Bushing installation

Material	Part No	Thread size	Rated pressure		Max. temperature at rated pressure	
			MPa	psi	°C	°F
1.4435 with Inspection Certificate 3.1B	52401947	G 1 1/2" NPT	1.75	250	200	390
	52401943	R 1 1/2"				
	52401949	G 2" NPT				
	52401945	R 2"				
PVDF	52401948	G 1 1/2" NPT	1.0	150	25	75
	52401944	R 1 1/2"	0.4	60	80	180
	52401950	G 2" NPT	0.2	30	120	250
	52401946	R 2"				

Table 5: Hex-Head Bushing Specifications

4.4 Sanitary Installation

Part No	Description	Material	Certificate
52 401 952	Dairy adapter DN50 (DIN 11851)	DIN 1.4435	3.1B Inspection Certificate
52 401 953	Aseptic adapter DN50 (DIN 11864-1)	DIN 1.4435	3.1B Inspection Certificate

METTLER TOLEDO offers an adapter according to DIN 11851 (DN50) as well as an aseptic adapter DIN 11864-1 (DN50) for the inductive sensors of the InPro 7200 series.

4.5 Immersion Installation

For open tanks or inlet channel the sensor can be installed in an immersion housing InDip 550Ind.


Figure 5: Immersion housing "InDip 550Ind"

4.6 In-Line Installation

For in-line applications within a metal or plastic pipe, the pipe must have a minimum diameter of DN80 (3"), and the sensor should be aligned as closely as possible with the centerline of the pipe.

Larger pipe diameter should be used wherever possible.

Note: Whenever the minimum size pipe diameter is used, the system should be calibrated with the sensor in place and with a solution of known conductivity (within the range of the instrument's calibration).

5 Junction Box Installation

5.1 Junction Box Wiring

1. Remove the junction box cover and loosen the cable connectors.
2. Insert the sensor cable through the appropriate connector and connect the numbered terminals of the sensor cable to the corresponding numbered terminals on the terminal strip.
Note: For maximum accuracy, sensor should not be located more than 30 m (100 ft) from Transmitter Cond Ind 7100 e, Cond I 7100/2H, Cond I 7100/2(X)H or M 700.
3. Insert extension cable assembly through appropriate connector and connect the numbered terminals of extension cable assembly opposite the corresponding number.

Junction box


Figure 6: Junction box "ZU 0074X"


Figure 7: Junction box "ZU 0307"

6 Degaussing and Sensor Resistance Check

6.1 Degaussing

To degauss a sensor, use a degaussing tool such as that used for erasing magnetic tape. With degaussing tool power on, bring the tool close to the sensor and move the tool slowly around the sensor in close proximity. Then slowly move the tool away from the sensor (min. 20 cm). Turn tool power off.

6.2 Sensor Resistance Check

Note: Use a high impedance meter when making these checks because a low impedance meter can cause the probe to become magnetized. Refer to Figure 8 for sensor cable resistances.

Sensor cable resistances


Figure 8: Sensor cable resistances

7 Disposal

The sensors described in this manual usually cannot be repaired. Their disposal is to the responsibility of the user. Sensors contain electronic components (coils, cables) that require special disposal without any health hazard for humans, and with no risk of harm for the environment. Please refer to your local regulations.

8 Dimension Drawings

Order No.	Designation	Ø A	Ø B	Ø C	D	Spacer in front	Spacer at the back
52 401 940	Flange DN50/PN16 3.1B	125	4x18	165	18	1x4.75	1x9.5
52 401 941	Flange DN100/PN16 3.1B	180	8x18	220	20	1x4.75	1x9.5
52 401 933	Flange ANSI 2" 3.1B	120.7	4x19	152.4	12.7	1x9.5	1x9.5
52 401 934	Flange ANSI 2 1/2" 3.1B	139.7	4x19	177.8	14.2	1x9.5	1x9.5
52 401 938	Flange ANSI 3" 3.1B	152.4	4x19	190.5	15.7	1x9.5	1x9.5
52 401 939	Flange ANSI 4" 3.1B	190.5	8x19	228.6	17.5	1x4.75	1x9.5


Figure 9: Flange


Figure 10: Threaded bushing

Order-No.	Designation	A	B	C
52 401 943	Threaded bushing R 1 1/2" 3.1B	R 1 1/2"	SW60	22
52 401 944	Threaded bushing R 1 1/2" PVDF	R 1 1/2"	SW60	22
52 401 945	Threaded bushing R 2" 3.1B	R 2"	SW75	24
52 401 946	Threaded bushing R 2" PVDF	R 2"	SW75	24
52 401 947	Threaded bushing G 1 1/2" NPT 3.1B	G 1 1/2" NPT	SW60	22
52 401 948	Threaded bushing G 1 1/2" NPT PVDF	G 2" NPT	SW60	22
52 401 949	Threaded bushing G 2" 3.1B	G 2" NPT	SW75	24
52 401 950	Threaded bushing G 2" PVDF	G 2" NPT	SW75	24


Figure 11: Sanitary adapter

Order No.	Designation	DIN	Ø A	Ø B	Ø C	D
52 401 952	Dairy adapter	11851	58	92	87	140
52 401 953	Aseptic adapter	11864-1	61	92	87	140


Figure 12: Immersion housing

9 Spare Parts

Order No.	Description	Material	Used for
52 401 942	Set locknut with spacers	1.4435	Flange installation
52 401 980	Set Locknut	1.4435	Bushing installation Sanitary/Aseptic installation Immersion installation
20 303 1112	O-ring	EPDM (FDA)	InPro 7200 Series
20 304 1100	O-ring	Kalrez (FDA)	InPro 7200 Series

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Page 2/3

	KEMA 
(13)	SCHEDULE
(14)	to EC-Type Examination Certificate KEMA 03ATEX1462 X
(16)	Report KEMA No. 2037022.
(17)	Special conditions for safe use 1. For applications in explosive atmospheres caused by air/dust mixtures and where category 1D apparatus is required, the surface temperature "T" is determined for a dust layer with a thickness of 5 mm maximum. 2. For applications in explosive atmospheres caused by gases, vapours or hazes or air/dust mixtures and where category 1 apparatus is required, precautions shall be taken such that electrostatic charges on the non metallic parts of the sensors are avoided. 3. For ambient temperature range and electrical data see (15).
(18)	Essential Health and Safety Requirements Covered by the standards listed at (9). Compliance with EN 50281-1-1 : 1998 has been assured, taking into account draft IEC 61241-0 : 2002 and draft IEC 61241-11 : 2002 as a guide.
(19)	Test documentation 1. EC-Type Examination Certificate KEMA 00ATEX1160 X dated 2. Drawing No. BS811ME, rev. A 06.05.2004 4S200YY, rev. A (5 sheets) 08.08.2003

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- CH** **Mettler-Toledo (Schweiz) AG**, Im Langacher, CH – 8606 Greifensee, Switzerland
- D** **Mettler-Toledo GmbH**, Prozeßanalytik, Ockerweg 3, D – 35396 Gießen, Germany
- USA** **Mettler-Toledo Ingold, Inc.**, 36 Middlesex Turnpike, Bedford, MA 01730 USA