English
Deutsch
Français
Español
Italiano

Operating Instructions **Compact Antistatic Kit Excellence**
Bedienungsanleitung **Kompakt Antistatik-Kit Excellence**
Mode d'emploi **Kit anti-statique compact Excellence**
Manual de instrucciones **Kit antiestático compacto Excellence**
Istruzioni per l'uso **Kit antistatico compatto Excellence**

METTLER TOLEDO

Operating Instructions
Compact Antistatic Kit

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Declaration of Conformity

EC Declaration of Conformity

This declaration of conformity is issued under the sole responsibility of the manufacturer:

Fraser Anti-Static Techniques Ltd
Scotts Business Park, Bampton GB-EX16 9DN

who declares that the following equipment:

Mettler Toledo Compact Ionizer Part 30090337

is in conformity with the relevant European Community harmonisation legislation:

Low Voltage Directive 2006/95/EC

Standard IEC 601010

EMC Directive 2004/108/EC

Standard IEC 61326-1

RoHS2 Directive 2011/65/EU

Signed for and on behalf of the manufacturer

A handwritten signature in black ink, appearing to read 'Robert Fraser'.

Robert Fraser, Director

5th September 2013

CE marking first affixed in September 2013.

1 Introduction

Your Antistatic Kit is an accessory for METTLER TOLEDO analytical balances and is used to discharge electrostatically charged objects to be weighed.

The guidelines stated in the Operating Instructions of the METTLER TOLEDO balance connected to your ionizer also fully apply to your Antistatic Kit.

1.1 Conventions and symbols used in these Operating Instructions

Key and button designations are indicated by a picture or text in square brackets (e.g. [] or [Define]).

- These symbols indicate an instruction:
- ▶ prerequisites
 - 1 steps
 - 2 ...
 - ⇒ results

2 Safety Information

2.1 Definition of signal warnings and symbols

Safety notes are indicated by signal words and warning symbols and contain warnings and information about safety issues. Ignoring safety notes can lead to personal injury, damage to the instrument, malfunctions and erroneous results.

Signal words

- | | |
|----------------|---|
| WARNING | for a hazardous situation with medium risk, possibly resulting in severe injuries or death if not avoided. |
| CAUTION | for a hazardous situation with low risk, resulting in damage to the device or the property or in loss of data or minor or medium injuries if not avoided. |
| Note | (no symbol)
for useful information about the product. |

Warning symbols

General hazard

Electrical shock

Prohibition symbol

Pacemaker

2.2 Product specific safety notes

Your Antistatic Kit is an accessory for METTLER TOLEDO analytical balances and is used to discharge electrostatically charged objects. Use the Antistatic Kit exclusively for this purpose. Any other type of use and operation beyond the limits of technical specifications without written consent from Mettler-Toledo AG, is considered as not intended.

WARNING

Persons with cardiac pacemakers

Persons with cardiac pacemakers must always maintain a safety distance of more than 50 cm from the Ionizer.


CAUTION

Damage of device


For use only in dry interior rooms.

3 Overview

Antistatic Kit

1	Ionizer	2	Indicator light
3	Side door marked with "Conducting Glass")	4	AC adapter of balance
5	Power supply cable of balance	6	Operating instructions
7	Power connection	8	Control cable

Overview Balance back

9	Aux 1 (connection for ionizer control)	10	Aux 2 (connection for ionizer control)
11	Socket for AC adapter		

4 Installation and Putting into Operation

The compact ionizer produces a corona of ionized air. Ionized air provides ions of the opposite polarity to neutralize static charge in the object being weighed.

4.1 Scope of delivery

- Compact ionizer
- Side door marked with "Conducting Glass"
- Operating Instructions

4.2 Assembling the Antistatic Kit

WARNING

Risk of electric shock

- Use only an approved METTLER TOLEDO AC adapter which has an earthed supply. The use of other SELV adapters without an earthed supply could result in dangerous operation.
- The emitter pin is resistively coupled to the high voltage electronics. The maximum current available from it is less than 40 μ A.
- Use only the 3-pin power cable with equipment grounding conductor which was supplied with your balance.
- Only a 3-pin grounded outlet for connecting your balance must be used.
- Only extension cables which meet this relevant standards and also have an equipment grounding conductor may be used.
- Intentional disconnection of the equipment grounding conductor is prohibited.


The ionizer should arrive in good condition, protected by its packaging.

- 1 Inspect the ionizer and all cables for damage caused in transit.
- 2 Report any damage immediately.
- 3 Do not use if there is any damage to the cables, connectors or the ionizer body.
- 4 Guide the cables so that they cannot become damaged or interfere with the measuring process.
- 5 Remove protective cap without squeezing.

4.2.1 Installing Antistatic Kit

- ▶ Balance is switched off.
 - ▶ Balance is disconnected from the power supply.
- 1 Remove the relevant side door (1).
 - 2 Remove the any of the 3 clips, for example the middle one (2).


- Push the compact ionizer in from back to front until it clicks into place. If necessary, gently squeeze the strips at the back.

- 1 Insert the side door (3) marked with "Conducting Glass" (4) at an angle of approx. 30° into the two slots.
- 2 Lift the side door against the balance until it clicks into place.
- 3 Push the side door completely to the front.
- 4 Check that the side door runs easily, otherwise it is not correctly inserted.
- 5 For XPE balances, return the handle to the upright position.
For XSE balances, clip the handle together.


- 1 Plug the power connection (5) into the socket for AC adapter.
- 2 Plug the control cable (6) into Aux 1 or Aux 2.
- 3 Plug the power supply cable (7) into the ionizer.
- 4 Connect the AC adapter (8) to the power supply.

4.2.2 Installing second Antistatic Kit

Note : For the installation of a second Antistatic Kit you will also need a second AC adapter, **see** Accessories (page 11).


- 1 Remove the relevant side door (1).
- 2 Remove the any of the 3 clips, for example the middle one (2).


- Push the compact ionizer in from back to front until it clicks into place. If necessary, gently squeeze the strips at the back.

- 1 Insert the side door (3) marked with "Conducting Glass" (4) at an angle of approx. 30° into the two slots.
- 2 Lift the side door against the balance until it clicks into place.
- 3 Push the side door completely to the front.
- 4 Check that the side door runs easily, otherwise it is not correctly inserted.
- 5 For XPE balances, return the handle to the upright position.
For XSE balances, clip the handle together.

- 1 Plug the control cable (5) into Aux 1 or Aux 2.
- 2 Plug the power supply cable (6) into the ionizer.
- 3 Connect the second AC adapter (7) to the power supply.

4.2.3 XPE balance settings for Antistatic Kit

The arrow buttons can be used to page forward or back to a menu page.

- 1 Switch on the balance [⏻].
- 2 Press [⏮].
- 3 Choose **Smart & ErgoSens** and tap the button [Define].
⇒ A selection window appears.
- 4 Choose **ErgoSens 1 (Aux1)** or **ErgoSens 2 (Aux2)** and tap on the associated button.
⇒ A selection window appears.
- 5 Scroll and select [AntiStatic Kit].
- 6 To save the settings, tap [OK].
- 7 Choose **Ionizer Setup** and tap the button [Define].
- 8 Choose settings appropriate to your application. Recommended settings are active and 0 seconds **After door closed**.
- 9 To save the settings, tap [OK].
⇒ Indicator light turns on and Status-Icon [🔌] appears in the display when the ionizer is working.

4.2.4 XSE balance settings for Antistatic Kit

The arrow buttons can be used to page forward or back to a menu page.

- 1 Switch on the balance [⏻].
- 2 Press [↩].
- 3 Choose **Smart & ErgoSens** and tap the button [Define].
⇒ A selection window appears.
- 4 Choose **ErgoSens 1 (Aux1)** or **ErgoSens 2 (Aux2)** and tap on the associated button.
⇒ A selection window appears.
- 5 Tap in the scroll bar to select [AntiStatic Kit].
- 6 To save the settings, tap [OK].
- 7 Choose **Function Keys** and tap the button [Define].
⇒ A selection window appears.
- 8 Tap [Ionizer].
- 9 To save the settings, tap [OK].
⇒ The function key [Ionizer] is displayed in the application at the bottom of the display.
⇒ Indicator light turns on and Status-Icon [⚡] appears in the display when the Ionizer is working.

5 Operation

CAUTION

Risk of ignition or sparking

The ionizer must NOT be used inside an inert gas atmosphere or with an inert gas purging system. There is a serious risk of sparks from the ionizer if such an environment is used.

The ionizer is controlled by the balance. It should be activated during loading the balance and deactivated during weighing.

Activating with XPE

With XPE balances the ionizer is controlled by the automatic door.

Activating with XSE

- To activate the ionizer, tap function key [Ionizer].
⇒ The indicator light on the ionizer turns on.
- To deactivate the ionizer, tap again function key [Ionizer].
⇒ The indicator light turns off.

If the indicator light doesn't turn on:

- 1 Check cable connections and settings.
- 2 If the indicator light remains off, contact your local METTLER TOLEDO service agency.

Further information, **see** balance Operating Instructions.

6 Maintenance

6.1 Cleaning

WARNING

Risk of electric shock

- The ionizer must be disconnected from the power supply, before cleaning or other maintenance work.
- Use only the power cable from METTLER TOLEDO, if it needs replacing.
- Ensure that no liquid comes into contact with the balance, the terminal, the ionizer or the AC adapter.
- Never open the balance, terminal, ionizer or AC adapter – they contain no components which can be cleaned, repaired or replaced by the user.
- Do not turn on power until the cleaned parts are dry.

Clean the ionizer with a soft brush or with ethyl alcohol and a lint-free cloth.

Quick cleaning

Clean the ionizer once a month.

- Clean the insulation, white PTFE, around the emitter pin (1) and the emitter pin itself.
Note : Avoid touching this area with fingers to prevent contamination.
⇒ PTFE is white again.
- Check that no ethyl alcohol is left on the ionizer.

Cleaning

Clean the ionizer every three months.

- Deinstall the ionizer.
- Clean the insulation, white PTFE, around the emitter pin (1) and the emitter pin itself.
Note : Avoid touching this area with fingers to prevent contamination.
⇒ PTFE is white again.
- Clean the insulation tube (3) between the steel emitter barrel (2) and the body of the ionizer.
- Check that no ethyl alcohol is left on the ionizer.

See Installing Antistatic Kit (page 5).

Note

Please contact your METTLER TOLEDO representative for details of the available service options.

6.2 Disposal

In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment. If you have any questions,

please contact the responsible authority or the distributor from which you purchased this device. Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.

Thank you for your contribution to environmental protection.

7 Technical Data

CAUTION

Use only with a tested AC Adapter with SELV output current.
Ensure correct polarity 

Ionizer

Input voltage:	9 V to 15 V DC
Input current (operating):	560 mA
Weight:	150 g (complete with packaging)

Protection and standards

Overvoltage category:	II
Degree of pollution:	2
Range of application:	For use only in closed interior rooms

Environmental conditions


Height above mean sea level:	up to 4000 m
Ambient temperature:	5-40 °C
Relative air humidity:	Max. 80% up to 31 °C, linearly decreasing to 50% at 40 °C, noncondensing
Atmosphere:	Do not use under inert gas atmosphere

Ozone

All HV ionizers produce a trace of ozone. The METTLER TOLEDO ionizer ozone level is considerably below the international safety limit of 0.1 ppm.


7.1 Dimensions

Dimensions in mm.

8 Accessories and Spare Parts

8.1 Accessories

	Description	Part No.
	AC/DC adapter (without power cable) 100–240 V AC, 0.8 A, 50/60 Hz, 12 V DC 2.5 A	11107909
	Country-specific 3-Pin power cable with grounding conductor.	00088751 30015268 00087920 30047293 00087452 00087925 00089405 00225297 11600569 00087457 11107881 11107880 00088668 00089728

8.2 Spare Parts

	Pos	Description	Part No.
		Side door "Conducting Glass"	11106263

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GWP®

Good Weighing Practice™

GWP® is the global weighing standard, ensuring consistent accuracy of weighing processes, applicable to all equipment from any manufacturer. It helps to:

- Choose the appropriate balance or scale
- Calibrate and operate your weighing equipment with security
- Comply with quality and compliance standards in laboratory and manufacturing

► www.mt.com/GWP

www.mt.com/lab-accessories

For more information

Mettler-Toledo AG, Laboratory Weighing

CH-8606 Greifensee, Switzerland

Tel. +41 (0)44 944 22 11

Fax +41 (0)44 944 30 60

www.mt.com

Subject to technical changes.

© Mettler-Toledo AG 10/2014

30090073B en, de, fr, es, it