**English**

Installation Instructions **QLL Kits for Bottles** for Liquid Dosing Applications

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
**METTLER TOLEDO**

---

Installation Instructions **QLL Kits for Bottles**

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---

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## Table of Contents

<b>1</b>	<b>Compatibility</b>	<b>3</b>
<b>2</b>	<b>Safety Information</b>	<b>3</b>
2.1	Further applicable documents .....	3
2.2	Definitions of signal words and warning symbols .....	3
2.3	Product-specific safety notes .....	3
<b>3</b>	<b>Installation and Putting into Operation</b>	<b>4</b>
3.1	Scope of delivery .....	4
3.2	Attaching the dosing head to the bottle cap .....	5
3.3	Connecting the tubes .....	5
<b>4</b>	<b>Technical Data</b>	<b>7</b>
4.1	Model-specific data .....	7
4.2	Material specifications .....	8
4.2.1	Glossary of material acronyms .....	8
4.2.2	Liquid dosing head QL001 .....	8
4.2.3	Liquid dosing head QL003 .....	9
4.2.4	QLL standard/advanced kit for bottles .....	9
<b>5</b>	<b>Spare parts</b>	<b>9</b>

## 1 Compatibility

	QL001	QL003
Q3 dosing module	✓	✓
QLX3 liquid module	✓	✓
Q2 powder module	✓	–
QLX45 liquid module	✓	–

## 2 Safety Information

### 2.1 Further applicable documents

Search for documents

► [www.mt.com/library](http://www.mt.com/library)

Refer to the Reference Manual (RM) of your balance or dosing module for more information. Manuals are available online or through your METTLER TOLEDO service representative.

► [www.mt.com/XPR-analytical-RM](http://www.mt.com/XPR-analytical-RM)

► [www.mt.com/XPR-automatic](http://www.mt.com/XPR-automatic)

### 2.2 Definitions of signal words and warning symbols

Safety notes contain important information on safety issues. Ignoring the safety notes may lead to personal injury, damage to the instrument, malfunctions and false results. Safety notes are marked with the following signal words and warning symbols:

#### Signal words

<b>DANGER</b>	A hazardous situation with high risk, resulting in death or severe injury if not avoided.
<b>WARNING</b>	A hazardous situation with medium risk, possibly resulting in death or severe injury if not avoided.
<b>CAUTION</b>	A hazardous situation with low risk, resulting in minor or moderate injury if not avoided.
<b>NOTICE</b>	A hazardous situation with low risk, resulting in damage to the instrument, other material damage, malfunctions and erroneous results, or loss of data.

#### Warning symbols

General hazard

Notice

### 2.3 Product-specific safety notes

#### Intended use

This kit for bottles is designed to be used in analytical laboratories by trained staff. The kit for bottles is intended for dosing liquid samples.

Any other type of use and operation beyond the limits of use stated by Mettler-Toledo GmbH without consent from Mettler-Toledo GmbH is considered as not intended.

#### Responsibilities of the instrument owner

The instrument owner is the person holding the legal title to the instrument and who uses the instrument or authorizes any person to use it, or the person who is deemed by law to be the operator of the instrument. The instrument owner is responsible for the safety of all users of the instrument and third parties.

Mettler-Toledo GmbH assumes that the instrument owner trains users to safely use the instrument in their workplace and deal with potential hazards. Mettler-Toledo GmbH assumes that the instrument owner provides the necessary protective gear.

## Safety notes

### CAUTION

#### Injury due to leaking liquids

Wrongly cut tubing can result in leaking connections.

- Cut the tubes with a tube cutter or a sharp knife.

### NOTICE

#### Damage to the instrument or malfunction due to the use of unsuitable parts

- Only use parts from METTLER TOLEDO that are intended to be used with your instrument.

## 3 Installation and Putting into Operation

### 3.1 Scope of delivery

Your package includes one of the following products:

#### QLL standard kit for bottles (GL45)

- QL001 liquid dosing head
- Bottle, 250 ml, pressure resistant
- Bottle cap (GL45) with tube connector and dosing head support
- Liquid tube 3.2/1.6 mm (2.5 m)
- Air tube 4.0/2.4 mm (2 m)
- Tube adapter (KK2P-04H)
- Micro dosing valve tool
- Suction filter
- Various small parts (rings, nuts, etc.)
- Installation Instructions

#### QLL standard kit for small bottles (GL25)

- QL001 liquid dosing head
- Bottle, 25 ml, pressure resistant
- Stand for bottle and dosing head
- Bottle cap (GL25) with tube connector
- Liquid tube 1.6/0.8 mm (2.5 m)
- Air tube 4.0/2.4 mm (2 m)
- Tube adapter (KK2P-04H)
- Micro dosing valve tool
- Various small parts (rings, nuts, etc.)
- Installation Instructions

#### QLL advanced kit for bottles (GL45)


- QL003 liquid dosing head
- Bottle, 250 ml, pressure resistant
- Bottle cap (GL45) with tube connector and dosing head support
- Liquid tube 3.2/1.6 mm (2.5 m)
- Air tube 4.0/2.4 mm (2 m)
- Tube adapter (KK2P-04H)
- Suction filter
- Various small parts (rings, nuts, etc.)
- Installation Instructions

The QLL kits for bottles are compatible with bottles with a GL45 thread. The QLL kits for small bottles are compatible with bottles with a GL25 thread.


### 3.2 Attaching the dosing head to the bottle cap

- 1 Insert the liquid dosing head in the liquid dosing head support.

- 2 To remove the liquid dosing head from the liquid dosing head support, pull the lever (1) towards the dosing head and remove the liquid dosing head (2).


### 3.3 Connecting the tubes

#### Tubes definition

The liquid tube is the thinner tube (1) used for transporting liquid from the bottle to the liquid dosing head. The air tube is the slightly bigger tube (2) used for pumping air into the bottle. When the pressure reaches the target pressure, between 0.3 and 0.5 bar (4.4 to 7.2 psi), the micro dispensing valve in the dosing head opens and liquid can ascend the liquid tube.

- 1 Liquid tube
- 2 Air tube

#### Preparing the liquid tube

#### ⚠ CAUTION

##### Injury due to leaking liquids

- Wrongly cut tubing can result in leaking connections.
- Cut the tubes with a tube cutter or a sharp knife.


- The liquid dosing head is inserted in the liquid dosing head support of the bottle.
- Using a tube cutter or a sharp knife, cut a sufficient amount of tubing. The appropriate length mainly depends on the distance between the balance and the bottle while dosing.  
Recommended length: about 0.9 m

- NOTICE: Leaking liquid due to wrong assembly. Pay attention to the orientation of the sealing ring when threading it on the tube.**

Place the sealing ring (1) on a flat stable surface, e.g., a table or a workbench, with the wider end downwards.

- Take the end of the liquid tube and press it into the sealing ring.
  - ➔ This is the dosing-head end of the tube. The opposite end is the bottle end.


- Thread the two fastening nuts (2), paying attention to the orientation.
- NOTICE: Leaking liquid due to wrong assembly. Pay attention to the orientation of the sealing ring when threading it on the tube.**  
Thread the sealing ring (3) from the bottle end of the tube.
- Slide the sealing ring until the distance (D) is enough for the tube to reach the bottom of the bottle. Recommended distances for typical bottle volumes are listed next.

#### Typical distances (D) between the sealing ring and the bottle-end of the tube

Thread	Bottle volume	Distance (D)
GL45	1 000 ml	220 mm
	500 ml	170 mm
	250 ml	135 mm
GL25	25 ml	90 mm

#### Connecting the liquid tube

- Insert the dosing-head end of the tube (1) in the dosing head (2).
- Tightly fasten the fastening nut to the dosing head.
- Insert the bottle end of the tube through the corresponding hole in the bottle cap (3). The tube should reach the bottom of the bottle.
- If needed, attach the suction filter (4) to the bottle end of the tube.
- Tightly fasten the fastening nut to the bottle cap.
- Screw the cap to the bottle (5).

The suction filter is used to ensure that no particles or impurities are carried through the liquid dosing head. Using the suction filter will extend the lifetime of the dosing head. However, when dosing solutions, molecules of one of the substances might be absorbed by the suction filter, altering the concentration of the solution. The suction filter should only be used when dosing pure solvents.

#### Connecting the air tube

For information on connecting the air tube to the pump and handling the bottle, consult the manual of the pump or dosing system.

## 4 Technical Data

### 4.1 Model-specific data

#### Liquid dosing heads

For further information about the performance of the liquid dosing heads (QL001, QL003) on your **Q3 dosing module** or **QLX3 liquid module**, consult the corresponding Reference Manual (RM).

► [www.mt.com/Q3-RM](http://www.mt.com/Q3-RM)

► [www.mt.com/QLX3-RM](http://www.mt.com/QLX3-RM)

For further information about the performance of the liquid dosing heads (QL001) on your **Q2 powder module** or **QLX45 liquid module**, consult the corresponding Reference Manual (RM).

► [www.mt.com/Powder-Module-RM](http://www.mt.com/Powder-Module-RM)

► [www.mt.com/Liquid-Module-RM](http://www.mt.com/Liquid-Module-RM)

#### Tubing

	Outer diameter	Inner diameter	Recommended length
Liquid tube, for GL45 bottles (dosing head to bottle)	3.2 mm	1.6 mm	0.9 m
Liquid tube, for GL25 bottles (dosing head to bottle)	1.6 mm	0.8 mm	0.9 m
Air tube (pump to bottle)	4.0 mm	2.4 mm	0.7 m
Tube for exhaust air	6 mm	–	–
Tube for external gas	6 mm	–	–

#### Suction filter

Tubing outer diameter	3.2 mm
Filter pore size	10 µm

#### Bottle

Pressure resistance, min. <sup>1)</sup>	1.5 bar
Burst pressure, min.	3 bar
Volume, max.	2 l

<sup>1)</sup> According to DIN EN 1595: Pressure Equipment made from Borosilicate Glass 3.3 – General Rules for Design, Manufacture and Testing

## 4.2 Material specifications

Materials in contact with the substance to be dosed. Please note that traces of all contact material could be transferred to the sample.

### 4.2.1 Glossary of material acronyms

Acronyms defined in ISO 1043: Plastics – Symbols and abbreviated terms.


<b>ETFE</b>	=	Ethylene tetrafluoroethylene
<b>FEP</b>	=	Fluorinated ethylene propylene
<b>FFKM</b>	=	Perfluoroelastomer
<b>PE</b>	=	Polyethylene
<b>PEEK</b>	=	Polyetheretherketone
<b>PMP</b>	=	Polymethylpentene
<b>POM</b>	=	Polyoxymethylene
<b>PP</b>	=	Polypropylene
<b>PP (ESD)</b>	=	Polypropylene (electrostatic dissipative)
<b>PTFE</b>	=	Polytetrafluoroethylene

### 4.2.2 Liquid dosing head QL001

	Designation	Material	
1	Fastening nut	PEEK	
2	Sealing ring	ETFE	
3	Tube connector	Stainless steel 1.4404	
4	O-Ring	FFKM	
5	Solenoid valve	PEEK Sapphire-Ruby Stainless steel 1.4105IL Stainless steel 1.14301 / 1.4306 Stainless steel 1.4305	

### 4.2.3 Liquid dosing head QL003

	Designation	Material
1	Fastening nut	PEEK
2	Sealing ring	ETFE
3	Nozzle	Ceramic Al <sub>2</sub> O <sub>3</sub> > 99.7%
4	Valve block	PTFE
5	Gasket	EPDM
6	Valve	PEEK


### 4.2.4 QLL standard/advanced kit for bottles

#### Note

The QLL kits for bottle include a liquid dosing head, see [Liquid dosing head QL001 ▶ Page 8], [Liquid dosing head QL003 ▶ Page 9].

	Designation	Material
1	Cap insert	PE
2	Bottle	Borosilicate glass
3	Tubing	FEP
4	Suction filter	PP

## 5 Spare parts

Spare parts are available for this product. Please contact your METTLER TOLEDO sales representative or look online for more details.

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© Mettler-Toledo GmbH 11/2021