METTLER TOLEDO

InPro 7100-VP Series

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4-electrode conductivity sensors

Instruction manual

Subject to technical changes without prior notice.

© It is strictly forbidden to reprint this instruction manual or any parts thereof without the written permission of Mettler-Toledo GmbH, Process Analytics, Im Hackacker 15, 8902 Urdorf, Switzerland. No section or excerpt whatsoever may be reproduced or with the assistance of electronic systems be edited, duplicated or distributed, in particular in the form of photocopies, photographs, magnetic media or other recording methods. All rights reserved, especially the right of duplication and translation as well as in regard to patent and registration rights.

InPro, InTrac, InFit, and InDip are trademarks of the METTLER TOLEDO Group. Viton and Kalrez are trademarks of DuPont Performance Elastomers LLC.

Contents

Contents

1	Product description	5
1.1	Introduction	5
1.2	Equipment and scope of delivery	5
1.3	Technical data	7
2	Safety	8
2.1	Application compatibility	8
2.2	Proper utilization	8
2.3	Safety measures	9
2.4	Use in Ex-zones	10
3	Initial start-up	12
3.1	Installation	12
3.2	Electrical connections	13
4	Maintenance and troubleshooting	14
4.1	Conditions of warranty	14
4.2	Maintenance	14
4.3	Troubleshooting	15
5	Removal from operation, storage, disposal	16
5.1	Removal from operation	16
5.2	Storage	16
5.3	Environmental protection	16
6	Dimensional drawings	17
7	Accessories	18

1 Product description

1.1 Introduction

Thank you for buying this InPro™ 7100-VP Series 4-electrode conductivity sensor from METTLER TOLEDO. InPro 7100-VP Series 4-electrode conductivity sensors are intended to be used to measure conductivity in medium to high conductivity water. The sensors are available in a wide selection of process connections to meet every application need.

Markings

The printed markings on each InPro 7100-VP Series 4-electrode conductivity sensor contain the following information:

METTLER TOLEDO	Sensor manufacturer
InPro 710X-VP	Type of sensor
Cell M: XXXXXX	Cell constant
Temp M: XXXXXX	Temperature constant
Serial No: XXXXXXXX	Serial number
Ex-classification	ⓧ II 1/2 G Ex ia IIB T3...T6 Ga/Gb (plastic sensors) ⓧ II 1/2 G Ex ia IIC T3...T6 Ga/Gb (metallic sensors/metallic sensors with plastic surface)
Examination certificate number:	SEV 14 ATEX 0129X CML 22 UKEX 2211X
Notified body:	Eurofins Electric & Electronic Product Testing AG, (No.: 1258)
Approved Body:	Eurofins E&E CML Limited, (AB: 2503)

1.2 Equipment and scope of delivery

Scope of delivery

The InPro 7100-VP Series sensors are delivered ready for use. Each sensor is accompanied by this instruction manual, an individual «Certificate of Quality», a Material Certificate following EN 10204 (3.1B) and a CE Declaration of Conformity Certificate. Other Certificates may be included as specified for the individual product description. Please check that the details given in the Certificate of Quality match the sensor label. For each sensor, the

relevant cell constant has been determined individually in solutions of known conductivity at 18–25 °C (64–77 °F) during the manufacturing process, and the value documented in the accompanying «Certificate of Quality».

4-electrode sensors are so-called diffuse-field sensors. The individually determined cell constant was established under laboratory conditions without external influences.

Depending on the installation configuration, external influences (e.g. from the pipe) can result in a shift in cell constant. If the application does not guarantee a minimum clearance of 30 mm (1.18") around the sensor, it is recommended to carry out recalibration of the sensor in installed position.

Packaging

The packaging consists of cardboard and plastic material.

Keep the packaging for later use during storage or for transportation of the sensor. Should you wish to dispose of the packaging material, please observe your local regulations as well as those data and instructions given in Chapter 5.3 of this manual.

Unpacking and inspection

Please check the sensors immediately during unpacking in order to determine possible damage or missing items. Any irregularities should immediately be reported to your carrier and to your supplier.

Each InPro 7100-VP Series sensor is equipped with an integral temperature probe Pt 1000, IEC 751 Class A.

The nominal cell constant is 0.25 cm⁻¹. The exact value is printed on the sensor label. For each sensor, the cell constant and temperature constant have been measured/established individually and the values have been documented in the accompanying Quality Certificate. All calibrations can be traced back to NIST and/or ASTM Standards.

Both the measuring range and the system measurement accuracy strongly depend upon which type of transmitter is employed.

1.3 Technical data

	InPro 7108-VP/ CPVC	InPro 7108-VP/ PEEK	InPro 7108/25/*-VP	InPro 7108-TC-VP InPro 7108-VAR-VP
Measuring principle	4-electrode sensor	4-electrode sensor	4-electrode sensor	4-electrode sensor
Electrode material	SS 316L	SS 316L or HA-C22	SS 316L or HA-C22	SS 316L
Body material	CPVC	PEEK	PEEK	PEEK
RTD	Built-in Pt1000	Built-in Pt1000	Built-in Pt1000	Built-in Pt1000
Sensor diameter	See drawing	See drawing	See drawing	See drawing
Insertion length	28 mm (1.10")	28 mm (1.10")	40/65 mm (1.57/2.56")	25 mm (0.98")
Max. sensor length	151 mm (5.96")	126.7 mm (4.99")	123/148 mm (4.86/5.84")	105 mm (4.14")
Process connection	- 1" NPT - 1" NPT Conduit	- 1" NPT	- Cap nut G 1 1/4"	- Tri-Clamp 1.5" - Tri-Clamp 2" - Tuchenhagen- Varivent DN 40-DN 125
Cell constant nominal	0.25 cm ⁻¹	0.25 cm ⁻¹	0.25 cm ⁻¹	0.25 cm ⁻¹
Cell constant accuracy	± 1.0%	± 1.0%	± 1.0%	± 1.0%
System accuracy ^a	± 5.0% or better	± 5.0% or better	± 5.0% or better	± 5.0% or better
Working conditions				
Max. pressure at 25 °C (77 °F)	7 bar (100 psi)	17 bar (250 psi)	17 bar (250 psi)	17 bar (250 psi)
Max. pressure at 95 °C (203 °F) ^b	-	7 bar (100 psi)	7 bar (100 psi)	7 bar (100 psi)
Measuring temperature range	-10...80 °C (14...176 °F)	-10...140 °C ^c (14...284 °F)	-10...140 °C ^c (14...284 °F)	-10...140 °C ^c (14...284 °F)
Temperature range (sterilization)	N/A	Sterilizable -10...140 °C ^c (14...284 °F)	Sterilizable -10...140 °C ^c (14...284 °F)	Sterilizable -10...140 °C ^c (14...284 °F)
Temperature accuracy at 25 °C (77 °F)	± 0.25 °C (± 0.5 °F)	± 0.25 °C (± 0.5 °F)	± 0.25 °C (± 0.5 °F)	± 0.25 °C (± 0.5 °F)
Design				
Temperature compensation	Pt1000 IEC class A	Pt1000 IEC class A	Pt1000 IEC class A	Pt1000 IEC class A
Cable connection	VarioPin (IP68)	VarioPin (IP68)	VarioPin (IP68)	VarioPin (IP68)
Wetted parts:				
- Metals	SS 316L	SS 316L or HA-C22	SS 316L or HA-C22	SS 316L
- Plastics	CPVC	PEEK (FDA)	PEEK (FDA)	PEEK (FDA)
- O-rings	N/A	N/A	EPDM (FDA)	N/A
Certificates				
Cell constant	•	•	•	•
Material certificate				
EN 10204 3.1	•	•	•	•
Material specification 2.1	•	•	•	•
ATEX	•	•	•	•

^a The system accuracy depends on the measuring range and the particular transmitter.

^b According to PED (Pressure Equipment Directive, 97/23/EC, Art. 3, Sec. 3).

^c Short term 150 °C (302 °F).

Practical measuring range

Transmitter	InPro 7108	InPro 7100	InPro 7100 i
M100	–	–	0.02–500
M300	0.02–650	0.02–400	–
M300 ISM	–	–	0.02–500
M400 4-w/2-w	0.02–650	0.02–400	0.02–500
M700	0.02–500	0.02–400	–
M800	0.02–650*	0.02–400	0.02–500
System accuracy	±5 %	±5 %	±5 %

All values in mS/cm. * M800 1-channel only.

Subject to technical changes

2 Safety

2.1 Application compatibility

The wetted material parts of the sensor (several different materials come into contact with the sample medium) can under some circumstance be incompatible with the particular composition of the process medium and/or of the operating conditions. Responsibility to verify application compatibility lies wholly with the user.

The compatibility of different types of material are outlined on

<http://www.coleparmer.com/techinfo/chemcomp.asp>.

Mettler-Toledo GmbH, Process Analytics accepts no responsibility whatsoever for the correctness or accuracy of such details.

2.2 Proper utilization

METTLER TOLEDO InPro 7100-VP Series sensors are intended solely for the precise measurement of conductivity in aqueous solutions in industrial applications.

Any other use, or any operation over and above that intended by the manufacturer, are not recommended and can lead to harm or injury to

material/equipment and persons. This is also relevant for applications which do not comply with the technical data labeled on the sensor. For any damage arising from such misuse, the user assumes full and sole responsibility.

2.3 Safety measures

The InPro 7100-VP Series sensors have been manufactured in line with state-of-the-art technology and in accordance with accepted technical safety regulations. Nevertheless, the sensors can still represent a source of risk and danger:

- if the sensors are operated by insufficiently trained personnel,
- if the sensors are employed incorrectly or not as intended by the manufacturer,
- if the sensors are not regularly maintained or serviced.

Local legislation and regulations must be observed at all times. Such stipulations do not form an integral part of this instruction manual.

It is necessary to use protective gear, including gloves, for persons coming in contact with the fluid. It is recommended to use gloves when handling the electrodes of the sensor to limit contamination of the sensor.

The user is responsible for the instruction and training of personnel. Additional copies of the instruction manual can be ordered from your supplier. This instruction manual is an essential element of the sensor equipment and must at all times be readily available to operators directly at the measurement site.

Before the sensor is removed from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature lowered to a safe range. Any escape of hot process fluid under pressure can cause damage to material/equipment or injury to persons.

No modification whatsoever may be carried out on the sensors. Any unauthorized modification or manipulation of the sensors results in immediate expiration of the full scope of warranty granted by the manufacturer.

2.4 Use in Ex-zones

Note!

For intended installation in an Ex-classified area, please observe the following guidelines:

Ex- classification:

II 1/2 G Ex ia IIB T3...T6 Ga/Gb (plastic sensors)

II 1/2 G Ex ia IIC T3...T6 Ga/Gb (metallic sensors/metallic sensors with plastic surface)

Marking and number of the test certificate:

SEV 14 ATEX 0129X

CML 22 UKEX 2211X

Introduction

1. The conductivity sensors made of plastic are according to Directive 2014/34/EU¹⁾ (ATEX 114) Annex I equipment group II, category 1/2G after RL 99/92/EC²⁾ (ATEX 137) in the zones 0/1 or 0/2 or 1/2 and the gas groups IIA and IIB, explosion hazards due to flammable substances in the field of temperature classes T1 to T6, may be used.
2. The conductivity sensors made of metal or metal coated with plastic to Directive 2014/34/EU¹⁾ (ATEX 114) Annex I to equipment group II, category 1/2G after RL 99/92/EC²⁾ (ATEX 137) in zones 0/1 or 0/2 or 1/2 and gas groups IIA, IIB and IIC, explosion hazards due to flammable substances in the field of temperature classes T1 to T6, may be used.
3. For use/installation, the requirements according to EN 60079-14 must be observed.

¹⁾ For UK statutory Requirements SI 2016 No. 1107

²⁾ For UK regulations 7 and 11 of the Dangerous Substances and Explosive Atmospheres Regulations 2002 (DSEAR)

Electrical ratings

Conductivity measuring circuit and temperature sensing circuit in type of protection intrinsic safety Ex ia IIC only for connection to a certified intrinsically safe circuit.

Maximum values:

$$U_i \leq 16 \text{ V}$$

$$I_i \leq 190 \text{ mA}$$

$$P_i \leq 200 \text{ mW}$$

The effective internal inductance and capacitance are negligible.

The values above apply, each as the sum of all the individual circuits of the associated intrinsically safe supply and evaluation unit.

Special conditions for safe use

1. The relationship between the maximum permissible ambient or media temperature and temperature class is shown in the following table:

For $U_i \leq 16 \text{ V}$, $I_i \leq 190 \text{ mA}$, $P_i \leq 200 \text{ mW}$.

Conductivity measuring circuit and temperature measuring circuits:

Temperature class	Maximum ambient or media temperature
T6	51 °C
T5	63 °C
T4	91 °C
T3	143 °C


2. For installation, the system drawings and wiring diagrams of the measurement systems should be observed.
3. The capacitance and inductance of the connecting cable has to be considered.
4. The metal body of the conductivity sensors or the safety weld-in socket or the fitting InFit™ 76*_*_* or InTrac™ 7**_*_* or other appropriate fitting must be electrically connected to the equipotential bonding system.

- The metal body of the conductivity sensors or the safety weld-in socket or the fitting InFit 76*_*** or InTrac 7**_*** or other appropriate fitting is optionally included in the routine pressure test of the system.

3 Initial start-up

3.1 Installation

The sensors are to be mounted in such a way that the sample medium flows directly into the sensor through the orifice at the tip and re-emerges through the vent holes. Any other installation position of the sensors can create the risk of formation of airlocks, or of contamination through deposits of solid matter.

Flow of medium must be across the face of the sensor. Air bubbles and solid matter deposits are to be avoided. A minimum clearance of 30 mm (1.18") between sensor and pipe wall must be maintained (see drawing above).

Vertical mounting (top entry) in a pipe is recommended only if the pipe is full and no air bubbles are able to develop in the flow.

In the event of side mounting (side entry) of the sensor, vertical upward flow of sample process medium must be ensured.

3.2 Electrical connections

All of the InPro 7100-VP Series sensors can be connected to the associated transmitter using the appropriate multiwire cable from METTLER TOLEDO.

Colors and function of the cable strands (supplied separately)

4-electrode sensor	Color
Cnd outer 1	white/blue
Cnd inner 1	white
Cnd inner 2	blue
Cnd outer 2	black
RTD	green
RTD sense	red
RTD ref/GND	Bare shield

Connection to the relative transmitter must be followed according to the wiring diagram (see respective transmitter instruction manual for detail).

METTLER TOLEDO supplies the cables. Please refer to chapter 7 "Accessories" for ordering information

4 Maintenance and troubleshooting

4.1 Conditions of warranty

METTLER TOLEDO guarantees the quality of materials and workmanship within a narrow range of manufacturing tolerances, so that the product purchased is free from any substantial deviations from material and manufacturing quality standards. The warranty is valid for a period of one year from date of delivery (ex-works). If within this warranty period, any repair or replacement should become necessary, and such cause is not due to misuse or incorrect application, please return the sensor, freight pre-paid, to your appropriate METTLER TOLEDO agency. Repair work will be carried out free of charge. Final decision on whether the defect is due to a manufacturing error or to incorrect operation of the sensor by the customer is made at the discretion of the Customer Service department of METTLER TOLEDO. After expiration of the warranty period, faulty sensors will be repaired or replaced on an exchange basis against payment of the costs involved.

4.2 Maintenance

Dirty or contaminated sensors can deliver incorrect measurement values. If fouling is presumed, the sensor is to be removed from operation and the electrodes as well as the insulation between the electrodes cleaned with a soft cloth. Do not clean electrodes with an abrasive paper or pad. Suitable cleaning solutions are mild detergents or strongly diluted acids (< 0.5 % by wt.) such as nitric acid. Never use cleaning agents that are not compatible with the material to be cleaned.

When handling acids, precautionary measures are to be taken.


The sensor must be thoroughly flushed with distilled or distilled water prior to re-installation.

Following cleaning and re-installation into the process, it can take from several minutes up to several hours until the sensor delivers the originally measured value.

4.3 Troubleshooting

Error	Possible cause	Corrective action
No signal on display	Electrical connections either missing or incorrectly in place	Check all connections and associated cabling
	Sensor is not in contact with the sample medium	Check installation for air pockets
No temperature signal	Temperature probe is not connected	Check all connections and associated cabling
	Transmitter does not support Pt1000 RTD	Replacement of transmitter or sensor necessary
Incorrect or unstable measurement reading	Sensor is contaminated/fouled	Clean the sensor (electrodes)
	Sensor is installed too close to pipe wall	Ensure that a minimum distance of 30 mm (1.18") between sensor tip and pipe wall is maintained (see drawing)

Ensure that a minimum distance of 30 mm (1.18") between sensor tip and pipe wall is maintained.

5 Removal from operation, storage, disposal

5.1 Removal from operation

The sensor is only conditionally subject to aging. When used as intended and appropriately maintained and serviced, the lifetime of the sensor can extend to several years.

Before removing the sensor from the process/process adapter, it must be ensured that the process pressure has been reduced to a safe level and the process temperature has lowered to a safe range. Any escape of hot process medium under pressure can cause damage to material/equipment or injury to persons.

After removal from the process the sensor should first be flushed with distilled water.

If the sensor incurs a defect, it can not be repaired and must be disposed under observance of prevailing local regulations.

5.2 Storage

If the sensor is not use, it can be stored dry. However, it has to be reconditioned in the process fluid accordingly prior to renewed operation. This procedure may take several hours to complete.


5.3 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.


6 Dimensional drawings


InPro 7108-VP/CPVC


InPro 7108-VP/PEEK


InPro 7108-25/40-VP


InPro 7108-25/65-VP

InPro 7108-TC-VP

InPro 7108-VAR-VP

	Tri-Clamp
52 002 008	1.5"
52 002 009	2"

	Tuchenhagen-Varivent®
52 002 790	DN40 - DN125

7 Accessories

O-rings

Sensor	Standard O-ring	Complementary O-rings	Order No.
InPro 7108-25/40 InPro 7108-25/65	EPDM FDA (20.29 × 2.62)		59 909 522
		Silicone, Si USP VI (20.29 × 2.62) Viton™ FDA (20.24 × 2.62) Kalrez™, KA USP VI (20.30 × 2.62)	59 905 677 59 908 834 59 909 490

Cable

The maximum cable length for 4-electrode sensors is 15 m (50ft) with transmitter Cond 7050 e, for all other transmitters 60 m (200ft).

Cable length	Order No.
1.5 m (5 ft)	58 080 201
3.0 m (10 ft)	58 080 202
4.6 m (15 ft)	58 080 203
7.6 m (25 ft)	58 080 204
15.2 m (50 ft)	58 080 205
22.9 m (75 ft)	58 080 206
30.5 m (100 ft)	58 080 207

Adaption DN 25 (for InPro 7108-25/40)

Designation	Material	Order No.
Safety weld-in socket G1 ¼" DIN, straight L = 47 mm	DIN 1.4435	52 400 518
Weld-in socket G1 ¼" DIN, straight L = 40 mm	DIN 1.4435*	59 901 127
Weld-in socket G1 ¼" DIN, straight L = 40 mm	PVC	59 901 205
Weld-in socket G1 ¼" DIN, straight L = 40 mm	PVDF	59 901 206
Safety weld-in socket G1 ¼" DIN, inclined L = 40 mm	DIN 1.4435	52 400 462
Weld-in socket G1 ¼" DIN, inclined L = 40 mm	DIN 1.4435*	59 901 124
Weld-in socket G1 ¼" DIN, inclined L = 40 mm	PVC	59 901 546
Weld-in socket G1 ¼" DIN, inclined L = 40 mm	PVDF	59 901 208
Blind plug, straight L = 40 mm, silicone seal	DIN 1.4435	59 901 287

* other materials available on request

Please contact your METTLER TOLEDO representative for more options.

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