
InFlow 751

Immersion Housing

METTLER TOLEDO

English Page 3

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© It is forbidden to reprint this instruction manual in whole or part.

No part of this manual may be reproduced in any form, or modified, copied or distributed using electronic systems, in particular in the form of photocopies, photographs, magnetic or other recordings, without written consent of Mettler-Toledo GmbH, Process Analytics, CH-8902 Urdorf, Switzerland.

InFlow 751 Immersion Housing

Instruction Manual

Contents

	Page
1 Introduction	4
1.1 Conventions	4
2 Important Notes	4
2.1 General	4
2.2 Safety Precautions	4
3 Description of Product and Materials	5
4 Installation and Start up Procedures	5
4.1 Fitting the Sensor Into the Housing	5
4.2 Installation of the Housing	5
4.3 Checking for Correct Installation	6
5 Maintenance	7
6 Product Specification	7
6.1 Supply	7
6.2 Technical Specification	7
7 Spare Parts	7
8 Terms of Warranty	12
9 Housing Dimensions	12
10 Disposal	12

InFlow is a registered trademark of the METTLER TOLEDO Group.
Viton® is a registered trademark of DuPont Performance Elastomers LLC.
All other trademarks are the property of their respective holders.

1 Introduction

These operating instructions describe how to use the housings InFlow® 751.

These housings are state of the art, high-grade engineered products, tested by METTLER TOLEDO.

Nevertheless, improper handling could be dangerous.

1.1 Conventions

This pictogram represents safety and hazard warnings which, if ignored could result in injuries to personnel and/or in material damage.

2 Important Notes

2.1 General

Immediately on receipt, check that the housing is complete and in good condition. Notify your supplier of any damage or deficiency. Please also refer to your supplier for further information on the ordering of spare parts and accessories.

2.2 Safety Precautions

Safety and hazard warnings could, if ignored, result in injuries to personnel and/or material damage.

Important! Work on the housing should be assigned only to trained personnel.

- Observe local regulations concerning the safety of people and property.
- Be sure to read and follow the instructions in this operating leaflet carefully.
- The InFlow 751 housings are designed to contain only METTLER TOLEDO sensors. Any other kind of use could be dangerous and is not permitted.
- The materials used for the housing are described in Section 3 "Description of Product and Materials". Make sure the materials are suitable for the required application.
- To ensure that the housing is correctly installed and maintained, follow the instructions given in this leaflet.
- Incorrect handling of the housing can result in a broken sensor and leakage from the vessel.
- Before doing anything to the installed housing, ensure that the process facility is in a safe condition (release pressure, empty, rinse, vent, purge, etc.).
- Use only clean sensors, housings and sockets. Replace any damaged seals and housing components immediately.
- Before starting up, always check the measuring system. Inspect the housing/sensor assembly and check for leaks from housing and apparatus.
- If ever in doubt, consult your supplier.

3 Description of Product and Materials

The insertion housings InFlow 751 serve as enclosures for METTLER TOLEDO sensors used for pH, redox and dissolved oxygen, conductivity and turbidity measurements, particularly in industrial waste water.

All parts of the housing in contact with the process (wetted parts) are made of PVC, or PVDF according to order specification. O-rings in contact with the process medium are made of Viton®.

4 Installation and Start up Procedures

The housing is delivered completely assembled. To fit the sensor, and to install the housing proceed as follows.

4.1 Fitting the Sensor Into the Housing

When installing 12 mm sensors with PG13.5 ensure that the white (PTFE-) washer is positioned directly beneath the 12 mm sensor head, with the tapered (chamfered) face pointing downwards. Then follows the O-ring.

Screw the sensor with Pg13.5 or NPT1" thread hand-tight into the housing. Put the housing into the T-piece and tighten the cap nut.

Screw the sensor with NPT ¾" thread laterally directly into the T-piece.

Connect the cable socket to the sensor.

For the removal of the sensor, carry out the above procedure in the reverse order.


4.2 Installation of the Housing

The housing can be directly fixed to the pipe by means of a bonding agent (PVC) or by welding (PVDF)


For installation of the InFlow 751 in a pipe, please follow the instructions of the pipe supplier.

Figure 1


Sensor/electrode with Pg13.5 thread mounted in an InFlow 751.

Sensor with a NPT 3/4" thread mounted in an InFlow 751.

Electrode with a NPT 1" thread mounted in an InFlow 751.

4.3 Checking for Correct Installation

Important! Each time before starting up, check the measuring system, inspect the sensor assembly and examine for leaks from housing and apparatus.

Do not begin operation until the measuring system has been checked and any necessary corrective action taken.

5 Maintenance

The sensor, housing and socket must be kept clean. Replace any damaged seals or components immediately.

Details on maintenance of the sensor are contained in the sensor operating instructions.

Attention: Before mounting the housing, it is imperative to clean and grease the O-rings with anti-seize.

Important: It is advisable that all O-rings are replaced periodically.

6 Product Specification

6.1 Supply

A standard delivery comprises the following items:

- Housing Type InFlow 751;
- Operating instructions.

6.2 Technical Specification

Temperature range:	PVC	0...60 °C;
	PVDF	0...100 °C
Pressure range:	PVC	4 bar/45 °C (58 psig/113 °F);
	PVC	1 bar/60 °C (14.5 psig/140 °F);
	PVDF	4 bar/75 °C (58 psig/167 °F);
	PVDF	1 bar/100 °C (14.5 psig/212 °F)
Material:	All wetted parts made of PVC or PVDF according to order specification. Seals made of Viton®.	
Dimensions:	See Chapter 9.	

For further details on sensors, see sensor data sheets or ask the supplier of your housing.

7 Spare Parts

Articles	Order-No.
Electrode holder PVC d=32	52 400 653
Electrode holder PVC d=50	52 400 654
Electrode holder PVDF d=63	52 400 658

8 Terms of Warranty


The housings are of high technical quality and undergo a policy of continuous design review to incorporate the latest advances. Their reliability is ensured by a thorough final inspection prior to leaving our factory.

The warranty is valid for one year from the date of delivery and covers any defects due to faulty materials or manufacture.

Not covered by the warranty are normal wear and tear and any damage caused by improper use (e.g. chemical incompatibility of the materials, etc.).

The warranty extends only to replacement or repair of deficient products, at our discretion. The warranty is void if the customer or others modify in any way the products supplied by us. Defects must be reported to the supplier immediately upon discovery, and, under all circumstances, within the warranty period.

9 Housing Dimensions

10 Disposal

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

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All rights reserved. 52 400 740 E