
InFlow 76X

Flow-Through Housing

METTLER TOLEDO

English: page 3

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Français: page 56

InFlow 76X

Instruction manual

How to use this instruction manual

This instruction manual is an integral part of the METTLER TOLEDO flow-through housing InFlow® 76X and contains notes and instructions that are important for safety and operation.

All persons working on or with the flow-through housings InFlow 76X must have first read and understood the chapters appropriate to the work in hand.

Please read this instruction manual carefully before using the flow-through housing InFlow 76X. Keep this document close to the unit, so that operating personnel may easily be able to refer to it at any time.

Please first read Chapter 1 "Introduction" and Chapter 2 "Safety instructions"!

Proprietary designations

InFlow, InTrac and InFit are registered trademarks of the METTLER TOLEDO Group. All other trademarks are the property of their respective holders.

Term explanation

The generic term InFlow 76X Series used in this instruction manual refers to flow-through housings of the following types:

- **InFlow 761:** 1.4404/316L version
- **InFlow 762:** PVDF version

© It is forbidden to reprint this instruction manual in whole or part.

No part of this manual may be reproduced in any form, or modified, copied or distributed using electronic systems, in particular in the form of photocopies, photographs, magnetic or other readings, without written consent of Mettler-Toledo GmbH, Process Analytics, CH -8902 Urdorf, Switzerland.

All rights reserved, in particular reproduction, translation and patenting/registration.

Contents

1	Introduction	6
2	Safety	7
2.1	Introduction	7
2.2	Declaration of conformity/type examinations	7
2.3	Flow-through housing designations	8
2.4	Intended use	8
2.5	Inappropriate use	8
2.6	Basic principles	9
2.7	Warning notices and symbols	9
2.8	Responsibilities, organizational measures	10
2.8.1	Responsibilities of the operator	10
2.8.2	Responsibilities of the personnel	10
2.8.3	Selection and qualification of personnel – basic duties	11
2.9	Product-specific hazards	11
2.9.1	Removal of electrode/sensor	11
2.9.2	Manipulation and maintenance work on the housings	12
2.9.3	Plastic flow-through housings	12
2.9.4	Installation in pressurized systems	12
2.9.5	Installation in potentially explosive areas (hazardous areas)	13
2.10	Residual hazards	13
2.10.1	Leaky connections	13
2.10.2	Medium residues	14
2.10.3	Heat protection	14
2.10.4	External impacts	14
2.11	Emergency measures	14
2.12	Safety measures	14
2.13	Modifications	15
3	Product and material description	16
3.1	Scope of delivery	16
3.2	Packing	16
3.3	Checking the shipment	16
3.4	Product overview	17
3.5	Recommended combination InFlow 76X and InTrac/InFit	17
4	Installation and startup	18
4.1	Installation of the InFlow 76X	18
4.2	Adaption of a InTrac or InFit housing to the InFlow 76X	19
4.3	Verification of the correct installation	21
5	Maintenance	22
6	Product specifications	23
6.1	Technical specifications	23
7	Accessories/spare parts	24
8	Terms of warranty	24
9	Storage and disposal	25
9.1	Storage	25
9.2	Disposal	25
10	Dimensional drawings InFlow 76X	26
10.1	Flow-through housing with flange connection	26
10.2	Flow-through housing with welded connection	26
10.3	Flow-through housing in combination with InFit and InTrac housing	27
10.4	Product key InFlow 76X	28

1 Introduction

- The flow-through housings InFlow 76X are safe to operate and has been tested by METTLER TOLEDO and dispatched ready for installation.
- Before starting to use the flow-through housing, carefully read this instruction manual: the safety precautions and warnings contained in it must be observed.

In addition to this instruction manual please also note the following:

- all local safety regulations
- all instructions and warning remarks in the publications of the products that are used in conjunction with the flow-through housing (electrodes, sensors, etc.)
- all safety precautions for the plant into which the flow-through housing InFlow 76X will be installed
- all instructions and warnings labeled on the flow-through housing InFlow 76X.

2 Safety

2.1 Introduction

The instruction manual contains the most important information for using the flow-through housings InFlow 76X efficiently and in accordance with regulations. A basic condition for safe handling and operation without malfunctions is the knowledge of these safety instructions and the observance of the further warnings in the instruction manual.

This instruction manual, and in particular the safety regulations, are intended for personnel entrusted with the operation and maintenance of the flow-through housings. It is assumed that these persons are familiar with the equipment in which the flow-through housing is installed. Therefore, before any work is started with the flow-through housing, this instruction manual must be read and understood by those persons involved.

The instruction manual must be stored where it is constantly accessible and available to any person working with the flow-through housing InFlow 76X.

On receipt of the shipment, check immediately:

- the flow-through housing and accessories for any sign of transport damage. Report any damage immediately to the carrier and to your supplier
- the type designation on the flow-through housing body
- for completeness of the supply.

Please notify your supplier immediately if the shipment is incomplete or in any way incorrect (see Section 3.1 "Scope of delivery").

2.2 Declaration of conformity/type examinations

NOTICE

The "Declarations of Conformity and Type Examinations" are dependent on the design and the individual type of the flow-through housing, and have no general validity for the complete InFlow product range.

2.3 Flow-through housing designations


Housing designation as well as part and serial numbers can be noted from the type plate and used for clear identification when communicating with the manufacturer.

NOTICE

The specifications shown on the type plate are dependent on the design and the individual type of flow-through housing and have no general validity for the complete InFlow product range.

Example InFlow 76X label:

2.4 Intended use

The flow-through housings InFlow 76X are designed and intended solely for safe and reliable mounting of the InTrac® and InFit® series sensor housings directly into the process or in a bypass (pipe). An electrode/sensor from METTLER TOLEDO, namely pH and Redox (ORP) combination electrodes as well as oxygen, CO₂, conductivity or turbidity sensor can be installed for measurement. Use the flow through housings only for this purpose.

The following are also part of the stipulations for the correct and appropriate use:

- compliance with the instructions, regulations and information contained in this instruction manual
- adherence to the prescribed inspection and maintenance/servicing intervals
- correct maintenance of the flow-through housings
- operation in compliance with prevailing regulations concerning the environmental and operating conditions as well as with the admissible mounting positions
- observance of local legislation.

2.5 Inappropriate use

Any utilization other than the above mentioned, as well as any utilization that is not consistent with the technical data is taken as being not in conformance with regulations. The operator bears the sole risk for any damage caused by such utilization.

2.6 Basic principles

The flow-through housing InFlow 76X is built in accordance with state-of-the-art technology and recognized technical safety regulations. However, the flow-through housing can be a source of risk and danger:

- if the flow-through housing is operated by insufficiently trained persons
- if the flow-through housing is not used in compliance with regulations and/or stipulations for appropriate use.

The flow-through housings InFlow 76X must be used only in perfect technical condition and for the purpose intended by the manufacturer. It is assumed that the user is fully aware of safety/danger issues and respects the instruction manual and all local safety regulations. Malfunction and damage that can affect safety and function of the flow-through housing must immediately be remedied by the operator or an expert, and notified to the manufacturer in writing.

2.7 Warning notices and symbols

This instruction manual identifies safety instructions by means of the following symbols:

**DANGER**

Warns of a hazard that could lead to extensive material damage, death or grievous bodily injuries.

**CAUTION**

Warns of a potentially hazardous situation that could lead to minor bodily injury and/or material damage.

**NOTICE**

Information relating to technical requirements. Non-observance could lead to defects, inefficient operation and possible loss of production.

2.8 Responsibilities, organizational measures

2.8.1 Responsibilities of the operator

- The operator must restrict permission to work with InFlow 76X flow-through housings to persons who are familiar with the basic requirements of work safety and accident prevention, and who have been instructed in the handling of the flow-through housing. This instruction manual serves as the basic document.
- In addition to the instruction manual there are also generally applicable legal and other binding regulations for work safety and accident prevention as well as for environmental protection, and these must be provided by the operator and instructed to personnel using the flow-through housings.
- The operator/user must periodically check that the personnel is fully aware of regulations on safety and risk prevention.
- Measures must be taken to ensure that the flow-through housings are only operated in a safe and fully functional condition.
- If the flow-through housings are employed in hazardous areas, compliance with prevailing regulations must be ensured.

CAUTION

Before the flow-through housing is put into operation, the operator has to make sure that use of the flow-through housing in conjunction with the other associated equipment and resources is fully authorized.

2.8.2 Responsibilities of the personnel

- All persons operating the flow-through housings must have read Chapters 1 "Introduction" and 2 "Safety" as well as the warning notices in this instruction manual.
- In addition to the instruction manual, generally applicable legal and other binding regulations for work safety and accident prevention must be adhered to.
- Avoid any kind of working that is doubtful from a safety perspective or which exceeds the admissible scope of use.
- Do not use high-pressure cleaning equipment for polymer /plastic components (of the flow-through housing).

NOTICE

Before every startup, the flow-through housing must be checked for:

- damage to the connections, fastenings, etc.
- leakage
- proper functioning
- authorization for use in conjunction with other plant equipment and resources.

DANGER

Defective flow-through housings must neither be installed nor put into operation. Leakage and inappropriate installation may lead to the escape of medium or to pressure surges (explosion), potentially harmful both to persons and to the environment.

2.8.3 Selection and qualification of personnel – basic duties

- Work on or with the flow-through housings may only be carried out by authorized and appropriately trained or instructed personnel. The personnel must have read this instruction manual in advance.
- Clear responsibilities must be established for the personnel entrusted with operation, service, repair, etc. of the flow-through housings.
- It must be ensured that only specifically assigned personnel may operate the flow-through housings.

**DANGER**

Incorrect manipulation or operation of the flow-through housing or non-observance of safety regulations can lead to malfunction of the flow-through housing and to the escape of process medium, thus presenting a potential hazard to the environment, personnel and material.

2.9 Product-specific hazards

2.9.1 Removal of electrode/sensor

**DANGER**

The electrode or sensor must be removed only after the piping/vessel has been depressurized and emptied. Otherwise, the removed electrode/sensor may cause the process medium to escape thus presenting a severe danger to persons, material and the environment. Any toxic or aggressive medium may cause severe poisoning or causticization.

**NOTICE**

Any broken sensor or defective/cut O-ring must be replaced without delay.

**DANGER**

Broken sensors prevent accurate measurements and thus adversely affect process safety.

2.9.2 Manipulation and maintenance work on the housings

NOTICE

Before dismantling a flow-through housing or commencing any maintenance work on it, ensure that the equipment in which the flow-through housing is installed is in a safe condition (depressurized, explosion-proof, emptied, rinsed, vented, etc.).

Flow-through housings may only be stripped down after having been completely dismantled. Manipulation of the sensor and the flow-through housing may only take place after it has been ensured that no process medium can escape through the flow-through housing in the event of incorrect manipulation. For this reason, the complete system must be emptied and vented in advance (safe condition).

It is mandatory to wear personal protective outfit such as protective goggles and clothing. Only maintenance and repair work specified in this instruction manual may be performed on the flow-through housing.

Exclusively use genuine spare parts from METTLER TOLEDO when replacing defective components (see Chapter 7 "Accessories/spare parts").

DANGER

Non-compliance with the prescribed maintenance instructions may endanger personnel and the environment.

2.9.3 Plastic flow-through housings

NOTICE

- Flow-through housings made of plastic do not have the same mechanical resistance as steel flow-through housings.
- Flow-through housings made of plastic require more frequent servicing than steel flow-through housings.
- Defective components may be replaced by an authorized service center only.

DANGER

Plastic flow-through housings may not be used in applications subject to high mechanical stress as this could lead to leakage or breakage of the flow-through housing.

2.9.4 Installation in pressurized systems

NOTICE

The specified maximum temperature and pressure must not be exceeded. The respective specifications depend on the design and type of the flow-through housing and are given on the individual type plates.

**DANGER**

If temperature and pressure limits are exceeded, there is a risk to the integrity of the system, thus presenting a potential threat to human life and to the environment.

**NOTICE**

Ample specifications of maximum permissible temperature and pressure are given in Chapter 6 "Product specifications".

2.9.5 Installation in potentially explosive areas (hazardous areas)

**NOTICE**

- The flow-through housings InFlow 76X have no own ignition source.
- It should be considered to include the METTLER TOLEDO flow-through housing InFlow 76X and the process connections in your recurring pressure test program for the complete plant as a whole.
- The operator must ensure that the housing can be used safely in conjunction with other associated plant resources.
- The flow-through housing and the process connections must be connected to the main potential equalization system of the plant.

**DANGER**

Non-observance of legal regulations concerning use in hazardous areas can endanger human life and the environment.

2.10 Residual hazards

**NOTICE**

Despite all precautionary measures taken, residual hazards still remain.

2.10.1 Leaky connections

- Connections can become loose through the effects of vibration.

**NOTICE**

The connections between the flow-through housing and the process adaptor must be checked regularly by the customer/operator, and kept in full working condition.

**DANGER**

Leaky connections can cause the process medium to escape to the environment, presenting a hazard for persons and the environment.

2.10.2 Medium residues

**DANGER**

When retracting/replacing an electrode/sensor, small quantities of process medium will remain at the electrode/sensor. If the medium is a toxic or environmentally harmful substance, or contains pathogenic germs, then such contamination must be removed and disposed of in accordance with the applicable regulations!

2.10.3 Heat protection

**DANGER**

The flow-through is not equipped with heat protection. During steam-sterilization procedure, the surface of the flow-through housing can reach high temperatures and cause burns.

2.10.4 External impacts

**NOTICE**

Objects falling on the flow-through housing can damage or destroy the unit, or cause leaks etc.

2.11 Emergency measures

**NOTICE**

Always observe and comply with local regulations!

2.12 Safety measures

**NOTICE**

Always observe and comply with local laws and regulations! These are not an integral part of this instruction manual.

**DANGER**

It is mandatory to wear protective equipment such as protective goggles and protective clothing. Aggressive process medium escaping from the system may be hazardous to persons or the environment.

**NOTICE**

The operator is responsible for the instruction of personnel. Additional copies of this instruction manual can be ordered from the equipment supplier. As an integral part of the flow-through housing, this instruction manual must at all times be readily accessible to users at the point of operation of the flow-through housing. The operator must inform the supplier/manufacturer of the insertion housing immediately about any safety-relevant incidents, or observations made, during use of the housing.

**DANGER**

Incorrect manipulation and/or instruction errors can lead to potential hazards for persons and for the environment.

**NOTICE**

Before every startup, the insertion housing must be checked for:

- damage to the connections, fastenings, etc.
- leakage
- defective cables and lines etc.
- authorization for use of the housing in conjunction with the associated plant resources.

**DANGER**

Defective flow-through housings must neither be installed nor put into operation. Leakage and inappropriate installation may lead to the escape of medium and cause a potential threat to life (including the risk of explosion).

2.13 Modifications

**NOTICE**

No attachments or modifications to the flow-through housings are allowed.

**DANGER**

The manufacturer/supplier accepts no responsibility for any damage caused by unauthorized attachments and alterations or for the incorporation of spare parts which are not of METTLER TOLEDO provenance. The risk is borne entirely by the operator.

3 Product and material description

3.1 Scope of delivery

Standard supply of the flow-through housing is made up of the following:

- flow-through housing InFlow 76X
- an instruction manual and a quick set-up guide

Additional with flange version:

- fixing screws
- washers
- nuts.

3.2 Packing

The packing consists of cardboard with protective paddings.

Keep the packing for possible later use, such as for storage or transport of the housing. However, if you wish to dispose of the packing, please observe your local regulations on waste disposal.

3.3 Checking the shipment

When unpacking the shipment, carefully examine for signs of damage. Any damage must be reported to the carrier and your supplier without delay. Check that the shipment meets the delivery papers and your order.

**NOTICE**

Damaged flow-through housings must not be installed or put into operation (see Chapter 2 "Safety").

3.4 Product overview

The InFlow 76X Series serves as a flow-through housing for METTLER-TOLEDO InTrac and InFit housings for the measurement of pH, ORP, dissolved oxygen, conductivity and turbidity, particularly in the field of industrial processes.

The flow-through housings type InFlow 76X are designed for the optimal fitting of the InTrac and InFit housings in a process or bypass pipe.

The attuned design allows for an optimal inflow of the sensor at maximum safety. All parts of the flow-through housing being in contact with the process medium are made of stainless steel 1.4404 (316L) or PVDF. The O-rings which come into contact with the process medium are made of Viton® (only InFlow 762 with Ingold DN25).

The flow-through housings are available in different versions, based on the product key (see Section 10.4 "Product key InFlow 76X").

3.5 Recommended combination InFlow 76X and InTrac/InFit

InFlow	InTrac	InFit
InFlow 761, 1.4404/316L		
InFlow 761 INGOLD DN25	InTrac series housing with recommended insertion length 70 –100 mm; stainless steel, and DN25 Ingold housing adaptation.	InFit 76X/xS/70/4435 / D00/xxxx InFit 76X/xS/70/4435 / D11/xxxx InFit 76X/xS/70/4435 / D10/xxxx
InFlow 761 flange DN50	InTrac series housing with recommended insertion length 70 –100mm; stainless steel, and DN50 flange housing adaptation.	on request
InFlow 762, PVDF		
InFlow 762 INGOLD DN25	InTrac series housing with recommended insertion length 70 –100mm; PVDF, and DN25 Ingold housing adaptation.	InFit 76X/xS/70/PVDF / D00/xxxx
InFlow 762 flange DN50	InTrac series housing with recommended insertion length 70 –100mm; PVDF, and DN50 flange housing adaptation.	InFit 76X/xS/70/PVDF / D04/xxxx

Please contact your local METTLER TOLEDO representative for more options.

4 Installation and startup

The flow-through housing InFlow 76X is delivered as a complete assembled unit and serves for fitting of an InTrac or InFit housing. In order to create a flow-through housing, proceed as follows:

4.1 Installation of the InFlow 76X


General

The flow-through housing can be directly integrated into a pipe. Special attention has to be paid to proper parallel alignment of the pipes (see drawing 1). Excessive pressure, tractive and torsional forces should be absolutely avoided. Pay attention to the appropriate pressure/temperature diagrams.

NOTICE

Attention has to be paid to the permitted installation position of housing and sensor. Appropriate information is available in the corresponding operating instructions.

Drawing 1

Welded connection

Direct welding of the InFlow 76X in the pipe has to be performed according to the instructions of the pipe supplier (pipe diameters see chapter 10 „Dimensional drawings InFlow 76X“).


Flange connection

A medium proof flat gasket has to be inserted between the flanges. The fixing screws have to be tightened uniformly and crosswise. The maximum torque is 10 Nm.

4.2 Adaption of a InTrac or InFit housing to the InFlow 76X


Flange version

Place the InTrac- or InFit housing onto the InFlow 76X. A medium proof flat gasket has to be inserted between housing and mounting flange (see drawing 2). The fixing screws have to be tightened uniformly and crosswise. The maximum torque is 60 Nm.



INGOLD DN25 socket version

Place the housing in the INGOLD socket (see drawing 3) and tighten the union nut by hand.
The torque of the union nut has no influence to the imperviousness of the connection!

Examples for installation:

Examples of possible assembly of an InFlow flow-through housing on a pipe, either as a 180° or 90° version. There is the choice of combination with either an InTrac or an InFit housing.

**NOTICE**

Avoid any installation positions that can create the risk of formation of airlocks, or contamination through deposits of solid matter to a sensor.

4.3 Verification of the correct installation

**NOTICE**

The measuring system has to be verified each time before startup. This includes inspection of the electrode/sensor assembly as well as examination of housing and system for leaks.


Do not commence operation until the measuring system has been checked and any necessary corrective action has been taken.

5 Maintenance

Always keep electrode/sensor, housing and socket clean. Replace defective seals and other components without delay. Information about maintenance of electrodes and sensors is available in the corresponding operating instructions.

**NOTICE**

It is recommended to regularly replace the O-rings (see drawing 4) of the InFlow 762 when using the INGOLD DN25 version (all other versions do not have any O-ring seals).


6 Product specifications

6.1 Technical specifications

Ambient conditions									
Temperature	–40 ... +100 °C (–40 ... +212 °F)								
Max. admissible pressure [PS]/[TS] (linear relation)	<table border="0"> <tr> <td>PVDF:</td> <td>6 bar/80 °C (87 psig/176 °F)</td> </tr> <tr> <td></td> <td>1 bar/140 °C (14.5 psig/284 °F)</td> </tr> <tr> <td colspan="2"><hr/></td> </tr> <tr> <td>1.4404/316L:</td> <td>16 bar/140 °C (232 psig/284 °F)</td> </tr> </table> <p>Please refer to the pressure/temperature diagram on the following page!</p>	PVDF:	6 bar/80 °C (87 psig/176 °F)		1 bar/140 °C (14.5 psig/284 °F)	<hr/>		1.4404/316L:	16 bar/140 °C (232 psig/284 °F)
PVDF:	6 bar/80 °C (87 psig/176 °F)								
	1 bar/140 °C (14.5 psig/284 °F)								
<hr/>									
1.4404/316L:	16 bar/140 °C (232 psig/284 °F)								
Medium-wetted materials	<table border="0"> <tr> <td>InFlow 761:</td> <td>1.4404/316L</td> </tr> <tr> <td>InFlow 762:</td> <td>PVDF (version INGOLD DN25 with Viton™ O-ring)</td> </tr> </table>	InFlow 761:	1.4404/316L	InFlow 762:	PVDF (version INGOLD DN25 with Viton™ O-ring)				
InFlow 761:	1.4404/316L								
InFlow 762:	PVDF (version INGOLD DN25 with Viton™ O-ring)								
Weight	<table border="0"> <tr> <td>InFlow 761:</td> <td>approx. 10 kg (22 lbs)</td> </tr> <tr> <td>InFlow 762:</td> <td>approx. 6 kg (13 lbs)</td> </tr> </table>	InFlow 761:	approx. 10 kg (22 lbs)	InFlow 762:	approx. 6 kg (13 lbs)				
InFlow 761:	approx. 10 kg (22 lbs)								
InFlow 762:	approx. 6 kg (13 lbs)								
Pressure domain	According to PED 97/23/EC, Article 1, Section 2.2: "Pressure" means pressure relative to atmospheric pressure, i.e. gauge pressure. As a consequence, vacuum is designated by a negative value.								
Explosion protection	The body of a flow-through housing type InFlow 76X must be electroconductively connected to the plant potential equalization system (functional ground).								
Certificates/approvals	<p>Depending on housing model, following certificates/approvals are included:</p> <ul style="list-style-type: none"> • CE Declaration of Conformity • Pressure Equipment directive, • Certificate of Conformity according to EN10204 2.1, material certificate according to 3.1B 								

For more detailed information regarding electrodes and sensors please refer to the corresponding data sheets or contact your housing supplier.

Pressure/temperature diagram

7 Accessories/spare parts

Designation	Order No.
O-ring set for InFlow 762 (PVDF) with INGOLD DN25, EPDM	52 403 716
O-ring set for InFlow 762 (PVDF) with INGOLD DN25, Kalrez®	52 403 717
O-ring set for InFlow 762 (PVDF) with INGOLD DN25, Viton®	52 403 715

8 Terms of warranty

The flow-through housings employ state-of-the-art technology and are very reliable. They are permanently revised to meet the latest technical evolutions and leave the factory only after having passed thorough quality testing.

The warranty period of 12 months starts from the date of delivery. The warranty covers all defects due to faulty material or manufacturing.

Normal wear or damage caused by inadequate use (such as chemical incompatibility, etc.) are not covered by the warranty.

In any case, the warranty is limited to the replacement or repair of faulty products, depending on manufacturer's choice. The warranty is void if the customer or any third party apply modifications to products supplied by METTLER TOLEDO. Any defect that may be subject to warranty must be reported to the supplier without delay or, at the latest, before the warranty period ends.

9 Storage and disposal

NOTICE

The safety notes in Chapter 2 "Safety" must be observed. Decommissioning may only be carried out by persons with appropriate training or by skilled technicians.

9.1 Storage

Store the InFlow 76X in a dry place.

9.2 Disposal

It is mandatory that the operator disposes of the device in accordance with local regulations. The operator must deliver the device either to a licensed private or public disposal company, or dispose of it himself in accordance with prevailing regulations. Waste is to be recycled or disposed of without causing any risk to human health, and without using procedures or methods that might harm the environment.

Sorting

Sorting into waste groups takes place when dismantling the device. The groups are listed in the current European Waste Catalogue. This catalog is valid for all kinds of waste, whether intended for disposal or for recycling.

The packaging is made up of the following materials:


- cardboard
- foam plastic.

The flow-through housing is made up of the following materials:


- steel
- PVDF (polyvinylidene fluoride)
- medium-wetted materials (plastics) according to specifications.

10 Dimensional drawings InFlow 76X

10.1 Flow-through housing with flange connection

180° version

90° version


a**InFlow 761, 1.4404/316L**

125 mm/4.92"


InFlow 762, PVDF

150 mm/5.91"

10.2 Flow-through housing with welded connection

180° version

90° version

a**b****c****InFlow 761, 1.4404/316L**

DN25 125 mm/4.92"

33.7 mm/1.33"

2.0 mm/0.08"

DN50 125 mm/4.92"

60.3 mm/2.37"

2.0 mm/0.08"

InFlow 762, PVDF

DN25 150 mm/5.91"

32 mm/1.26"


2.4 mm/0.09"

DN50 150 mm/5.91"

63 mm/2.48"

3.0 mm/0.12"

10.3 Flow-through housing in combination with InFit and InTrac housing

InFlow 76X with InTrac 777 e

InFlow 76X with InFit 761 e

	L	L1
InTrac 776 e	650 mm/25.59"	715 mm/28.15"
InTrac 777 e	470 mm/18.5"	535 mm/21.06"
InFit 764 e	510 mm/20.08"	640 mm/25.2"
InFit 761 e	363 mm/14.29"	493 mm/19.41"

10.4 Product key InFlow 76X

Material *									
1	/	4	4	0	4	1.4404/316L			
2	/	P	V	D	F	PVDF (polyvinylidene fluoride)			
Flow-through direction									
1	8	0	180°						
-	9	0	90°						
Process connection *									
D	2	5	DIN flange DN25 PN16						
D	5	0	DIN flange DN50 PN16						
A	0	1	ANSI flange A150-1"						
A	0	2	ANSI flange A150-2"						
W	2	5	Welding connection DN25 (1")						
W	5	0	Welding connection DN50 (2")						
Housing connection *									
D	0	0	INGOLD DN25						
D	0	4	DIN flange DN50						
Special									
-	Standard								
S	Special								

InFlow 76 / / / / / /

9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28

* others on request

Seals:

The PVDF version with INGOLD DN25 socket is fitted with a medium-wetted O-ring made of Viton®. O-ring sets made of EPDM and Kalrez® are available as accessories.

- Not all configurations are possible.
- Please contact your local METTLER TOLEDO representative for more options.

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All rights reserved. 52 403 711 E