# InSUS Head for Optical Single-Use O<sub>2</sub> Sensors

**Instruction manual  
Bedienungsanleitung  
Instructions d'utilisation**

InSUS Head  
30 283 375 F

**English** **Page** **3**

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# **InSUS Head for Optical Single-Use O<sub>2</sub> Sensors**

## **Instruction manual**

InSUS is a trademark of the METTLER TOLEDO Group.

## Contents

|          |  |           |
|----------|--|-----------|
| <b>1</b> | <b>Introduction</b> .....  | <b>5</b>  |
| <b>2</b> | <b>Important notes</b> .....   | <b>6</b>  |
| 2.1.     | Notes on operating instructions.....                                       | 6         |
| 2.2.     | Intended use.....  | 6         |
| 2.3.     | Safety instructions .....  | 7         |
| 2.4.     | Environmental protection.....  | 7         |
| 2.5.     | Examples of some typical applications .....                                | 8         |
| <b>3</b> | <b>Product description</b> .....   | <b>9</b>  |
| 3.1.     | General information.....   | 9         |
| 3.2.     | Principle .....  | 9         |
| 3.3.     | Scope of delivery.....   | 10        |
| <b>4</b> | <b>Installation</b> .....  | <b>11</b> |
| 4.1.     | Connection.....  | 11        |
| 4.1.1.   | Connection of the single-use sensor to the InSUS head .....                | 11        |
| 4.1.2.   | Digital connection of the InSUS head to a transmitter .....                | 11        |
| 4.1.3.   | Analog connection of the InSUS head to a transmitter or<br>controller..... | 12        |
| 4.1.4.   | Power connection of sensor head .....                                      | 12        |
| <b>5</b> | <b>Operation</b> .....   | <b>13</b> |
| 5.1.     | Configuration.....   | 13        |
| 5.1.1.   | Sensor head detection .....  | 13        |
| 5.1.2.   | Sampling rate.....   | 13        |
| 5.1.3.   | LED mode .....   | 13        |
| 5.2.     | Calibration .....  | 14        |
| 5.2.1.   | Purpose of calibration.....  | 14        |
| 5.2.2.   | Factory calibration.....   | 15        |
| 5.2.3.   | Single point calibration (Slope or Process Calibration) .....              | 15        |
| 5.2.4.   | Process calibration .....  | 15        |
| 5.2.5.   | Dual point calibration with M400 / M800 transmitter.....                   | 16        |
| 5.2.6.   | Calibration when connected with analog signal (via M50)..                  | 17        |
| <b>6</b> | <b>Maintenance</b> .....   | <b>18</b> |
| 6.1.     | Sensor head inspection.....  | 18        |
| 6.1.1.   | Visual inspection.....   | 18        |
| 6.1.2.   | Testing the sensor with the transmitter .....                              | 18        |
| <b>7</b> | <b>Storage</b> .....   | <b>19</b> |
| <b>8</b> | <b>Product specification</b> .....   | <b>20</b> |
| 8.1.     | Specifications .....   | 20        |
| <b>9</b> | <b>Ordering information</b> .....  | <b>21</b> |
| 9.1.     | Sensors .....  | 21        |
| 9.2.     | Accessories .....  | 21        |
| 9.3.     | Recommended transmitters .....   | 21        |

## 1 Introduction

Thank you for buying the **optical oxygen sensor head from METTLER TOLEDO**.

The construction of INGOLD's optical oxygen sensors employs leading edge technology and complies with safety regulations currently in force. Notwithstanding this, improper use could lead to hazards for the user or a third-party, and / or adverse effects on the plant or other equipment.

**Therefore, the operating instructions must be read and understood by the persons involved before work is started with the sensor.**

The instruction manual must always be stored close at hand, in a place accessible to all people working with the sensor.

If you have questions, which are not or insufficiently answered in this instruction manual, please contact your METTLER TOLEDO supplier. They will be glad to assist you.

## 2 Important notes

### 2.1. Notes on operating instructions

These operating instructions contain all the information needed for safe and proper use of the optical sensor head.

The operating instructions are intended for personnel entrusted with the operation and maintenance of the sensors. It is assumed that these persons are familiar with the equipment in which the sensor is installed.

#### Warning notices and symbols

This instruction manual identifies safety instructions and additional information by means of the following symbols:

This symbol draws attention to **safety instructions and warnings of potential danger** which, if neglected, could result in injury to persons and / or damage to property.

This symbol identifies **additional information and instructions** which, if neglected, could lead to defects, inefficient operation and possible loss of production.

### 2.2. Intended use

**METTLER TOLEDO optical O<sub>2</sub> sensors are intended solely for inline measurement of the oxygen partial pressure, as described in this instruction manual.**

Any use of these sensors which differs from or exceeds the scope of use described in this instruction manual will be regarded as inappropriate and incompatible with the intended purpose. The manufacturer / supplier accepts no responsibility whatsoever for any damage resulting from such improper use. The risk is borne entirely by the user / operator.

Other prerequisites for appropriate use include:

- compliance with the instructions, notes and requirements set out in this instruction manual.
- acceptance of responsibility for regular inspection, maintenance and functional testing of all associated components, also including compliance with local operational and plant safety regulations.
- compliance with all information and warnings given in the documentation relating to the products used in conjunction with the sensor (single-use bags, transmitters, etc.).
- observance of all safety regulations governing the equipment in which the sensor is installed.
- correct equipment operation in conformance with the prescribed environmental and operational conditions, and admissible installation positions.
- consultation with Mettler-Toledo Process Analytics in the event of any uncertainties.

## 2.3. Safety instructions

- The plant operator must be fully aware of the potential risks and hazards attached to operation of the particular process or plant. The operator is responsible for correct training of the workforce, for signs and markings indicating sources of possible danger, and for the selection of appropriate, state-of-the-art instrumentation.
- It is essential that personnel involved in the commissioning, operation or maintenance of these sensors or of any of the associated equipment (e.g. single-use bags, transmitters, etc.) be properly trained in the process itself, as well as in the use and handling of the associated equipment. This includes having read and understood this instruction manual.
- The safety of personnel as well as of the plant itself is ultimately the responsibility of the plant operator. This applies in particular in the case of plants operating in hazardous zones.
- The oxygen sensors and associated components have no effect on the process itself and cannot influence it in the sense of any form of control system.
- Maintenance and service intervals and schedules depend on the application conditions, composition of the sample media, plant equipment and significance of the safety control features of the measuring system. Processes vary considerably, so that schedules, where such are specified, can only be regarded as tentative and must in any case be individually established and verified by the plant operator.
- Where specific safeguards such as locks, labels, or redundant measuring systems are necessary, these must be provided by the plant operator.
- A defective sensor head must neither be installed nor put into service.
- Only maintenance work described in this operating instruction may be performed on the sensors.
- When changing faulty components, use only original spare parts obtainable from your METTLER TOLEDO supplier (see spare parts list, "Section 9").
- No modifications to the sensors and the accessories are allowed. The manufacturer accepts no responsibility for damages caused by unauthorised modifications. The risk is borne entirely by the user.

## 2.4. Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

## 2.5. Examples of some typical applications

Below is a list of examples of typical fields of application for the oxygen sensors. This list is not exhaustive.

### **Measurement in liquids:**

- Fermentation
- Bio-Tech

## 3 Product description

### 3.1. General information

The optical single-use oxygen sensor, the InSUS fiber optical cable (with integrated temperature probe), and the InSUS head are used for measurement of oxygen in single-use devices.

The single-use sensors are gamma irradiation sterilizable and are typically integrated in single-use devices such as single-use bags and other similar applications. Please also refer to the documentation provided by your single-use device supplier.

### 3.2. Principle


The optical oxygen sensors are based on an optical detection method, the so called fluorescence quenching. Here is a short summary of the principle. In contrast to the polarographic Clark-electrode, which detects a redox reaction of oxygen at the electrode, the optical method is based on an energy transfer between a chromophore and oxygen.

- A chromophore, embedded in the sensor is illuminated with blue light. This chromophore absorbs the energy and if no oxygen is present emits red fluorescence light with a specific lifetime. This emitted light is being detected by a detector in the sensor head.
- In the presence of oxygen, the chromophore transfers the energy to the oxygen molecule. Oxygen is then able to transfer this energy as heat to the surrounding area and no fluorescence is emitted.
- The total intensity of the fluorescence and the lifetime of the fluorescence is related to the Oxygen partial pressure in the medium.
- To analyze the lifetime of the fluorescence, the excitation light is pulsed with a constant frequency, the emitted light shows the same course but with a time delay to the excitation. This delay is called Phase shift or Phase angle ( $\Phi$ ). The phase shift is dependent on the oxygen level and follows the Stern-Vollmer correlation.
- The sensor head detects this phase shift and calculates the oxygen concentration.
- The oxygen value is digitally transferred from the InSUS Head to the transmitter.

### 3.3. Scope of delivery

Each sensor is supplied fully assembled and factory tested and calibrated for correct function together with:

- a certificate of calibration
- this instruction manual

## 4 Installation

### 4.1. Connection

#### 4.1.1. Connection of the single-use sensor to the InSUS head

Remove the protection cap from the single-use sensor. Connect the single-use sensor and the InSUS head with the corresponding InSUS fiber optical cable. Observe the information sheet that accompanies the InSUS fiber optic cable, then tightly screw both plugs on the sensor and the sensor head.

#### 4.1.2. Digital connection of the InSUS head to a transmitter

Transmitter M400 or M800

The sensor head is connected to the transmitter via a 5 pin data cable. The cables are available from METTLER TOLEDO in different lengths. The data cable ensures a secure connection between the transmitter and the sensor head under harsh industrial conditions. To connect the data cable to the sensor head align the slit of the connector with the pin in the plug, then tightly screw the plug to fasten the two parts.

Transmitter M400

**Note:** For connecting the cable to the terminals of the transmitter, please refer also to the instructions given in the METTLER TOLEDO transmitter manual.

Connect the data cable to the transmitter as described in the tables below.

#### RS485 cable for InSUS Head

|        |          | M400 | M800<br>1-channel | M 8 0 0<br>2/4-channel |
|--------|----------|------|-------------------|------------------------|
| Color  | Function | TB3  | TB3               | TB2                    |
| Brown  | 24 V     | 18   | 7                 | 9                      |
| Black  | GND24V   | 17   | 8                 | 10                     |
| Gray   | Shield   | 13   | 2                 | 12                     |
| Yellow | Shield   | 13   | 5                 | 15                     |
| Blue   | RS485B   | 14   | 3                 | 13                     |
| White  | RS485A   | 15   | 4                 | 14                     |

#### **4.1.3. Analog connection of the InSUS head to a transmitter or controller**

A METTLER TOLEDO M50 signal converter RS485 / nano-Ampere mimics a polarographic Clark sensor. The M50 signal converter is connected to an analog signal input of a transmitter, e.g. M300, or bio-controller. Please refer to the instructions given in the METTLER TOLEDO M50 converter manual.

#### **4.1.4. Power connection of sensor head**

If the sensor is used with a 4-wire transmitter (M400 / M800/M50), the Sensor Head is powered over the transmitter.

Power Supply Specification: 24 VDC; 500 mA

## 5 Operation

### 5.1. Configuration

#### 5.1.1. Sensor head detection

Before installing an optical sensor, please refer to the manual for the transmitter and configure the transmitter for automatic sensor detection. Date and time must be set correctly in the transmitter.

In the case wrong date and time are set, the calibration and setup might be corrupted

#### 5.1.2. Sampling rate

Optical oxygen sensors do not measure permanently. Each measurement cycle has a duration of approx. 1 second. The measurement interval can be set to any value between 1 and 60 seconds. Please choose the appropriate setting. Default setting is 3 seconds which is sufficient for most applications.

#### 5.1.3. LED mode

The measurement can be switched off if the system is not needed.

When the sensor head is not measuring, the sensor LED is off. In this state the sensor head sends a constant measurement value of  $-1\%$  air to the transmitter and the transmitter is set to the "Hold mode". To configure the "Hold mode" please refer to the transmitter manual.

#### **Automatic switch off at high temperature**

If the LED mode is set to "Auto" (default setting) the sensor LED will be switched off as soon as a specific process temperature is reached.

#### **Temperatures**

|            | Maximum operating temperature | Default switch off temperature |
|------------|-------------------------------|--------------------------------|
| InSUS Head | 60 °C / 140 °F                | 60 °C / 140 °F                 |

This limit can be set to an individual value by the user. Using the transmitter (M400 or M800) or with iSense. These settings are also active if the sensor is operated with an M50 signal converter. The switch off temperature should be set at least  $5^{\circ}$  higher than the highest process temperature. For example, if the process temperature is  $37^{\circ}\text{C} / 99^{\circ}\text{F}$ ,  $42^{\circ}\text{C} / 104^{\circ}\text{F}$  should be the minimum set-point. In this situation, as soon as the temperature exceeds  $42^{\circ}\text{C} / 104^{\circ}\text{F}$  the sensor will stop measuring and the LED will be switched off. For the switch on, a hysteresis of  $3^{\circ}$  is implemented, meaning that the sensor (and LED) will be switched on as soon as the temperature drops below  $39^{\circ}\text{C} / 101^{\circ}\text{F}$ .

#### **Manual switch off of the sensor head (M400 / M800 transmitter)**

The sensor can be switched off manually via the transmitter menu (see the transmitter manual) by setting

the LED mode to "off". To restart the measurement, the LED mode needs to be set manually to "on" via the transmitter menu, or via a remote signal (digital input).

### **Remote switch off of the sensor head (M400 / M800 transmitter)**

The M400 transmitter can be set to "Hold" by applying an external digital signal (see the transmitter manual). In this situation the sensor and the sensor LED are switched off. As soon as the "Hold Mode" is off, the optical sensor will continue to measure using the previous settings.

## **5.2. Calibration**

### **5.2.1. Purpose of calibration**

Information about the calibration, you find also in the manual of the transmitter.

Calibration must be performed every time a new single-use sensor is connected to the InSUS head via the InSUS optical fiber cable.

Since the correlation between the measured phase and the oxygen value is not linear, a calibration of an optical sensor must be performed very accurately. Wrong calibrations may significantly reduce the measurement accuracy.

Each oxygen sensor has its own individual phase angle at zero oxygen ( $\phi 0$ ) and hundred percent air saturation ( $\phi 100$ ).

Several methods for calibration are available for the optical oxygen sensors. The highest measurement accuracy is achieved by performing a 2-point calibration with air and a zero gas e.g.  $N_2$  or  $CO_2$  with a purity of at least 99.99 %.

**Note: Please take into account the correct air pressure and humidity. Small deviations in air ( $\pm 3\%$ ) are due to differences in humidity and process pressure settings. The sensor calculates for 100 % humidity if it is set to dissolved oxygen measurement.**

#### **General remarks:**

- **Calibration must be performed after gamma sterilization** of the single-use device.
- **For calibration in gas** the single-use device in which the single-use sensor is integrated **must be filled with air ( $\phi 100$ ) or zero point gas ( $\phi 0$ )**.
- Make sure that the settings for **oxygen saturation** of the calibration is **correct and remains constant** during calibration.
- In the event of calibration in water or sample medium, the **calibration medium must be in equilibrium with the air**. Oxygen exchange between water and air is very slow. Therefore it takes quite long time until water is saturated with atmospheric oxygen.
- **Make sure that all other parameters, such as temperature and pressure, are constant.**

- Calibration always needs accurate pressure and temperature measurement. Only process scaling is independent of those parameters (see chapter 5.2.4).
- Make sure that the correct calibration pressure, humidity and salinity values are set in the transmitter before the calibration is started.
- Please refer also to the transmitter or biocontroller manual for detailed informations

### 5.2.2. Factory calibration

The sensor head is delivered pre calibrated and ready for use.

The factory calibration data are stored in the sensor and can not be changed by the user.

### 5.2.3. Single point calibration (Slope or Process Calibration)

For most applications, a single point calibration should be sufficient.

By carrying out a single point calibration, the factual phase at the desired oxygen value e.g. at hundred percent oxygen (phi 100) of the sensor can be established. The corresponding calibration curve is calculated.

The calibration medium can be either air or a calibration gas with known O<sub>2</sub> concentration or water with a known oxygen concentration.

Before starting the calibration in gas, the correct pressure and the correct humidity have to be set in the transmitter.

**Note: Wrong pressure values are the most common reasons for bad measurement accuracy.**

E.g. 50 mbar difference between the ambient pressure and the value set in the transmitter result in 5% measurement error at air.

For calibration in gas it is important that the temperature reading of the sensor is stable and represents the real gas temperature.

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O<sub>2</sub>, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

### 5.2.4. Process calibration

For detailed information please refer also to the transmitter manual.

Two different routines for process calibration are possible:

- **Process calibration**
- **Process scaling (M400 / M800 transmitter)**

**Process calibration** is performed when a reliable control value is available and process pressure is known. Process pressure is only needed if the system is measuring in saturation (% air or % O<sub>2</sub>) or gas (ppm gas) units. During this calibration the phase values of the calibration curve are adjusted.

**Process scaling** is performed mainly in biopharma applications when the user desires to set the system to an initial value. During this calibration the phase values of the sensor are not adjusted, only the displayed values and the 4 – 20 mA output are rescaled to the desired value.

**Note:** For process calibration the operator can use either the process pressure or the calibration pressure, depending on how the reference value is taken.

After the sensor signal has stabilized, the complete measurement system can be calibrated to the desired variable, e.g. % air, % O<sub>2</sub>, ppm or ppb (see instruction manual for the transmitter).

**Note:** For this type of calibration an accurate reference value and correct pressure settings are essential.

### 5.2.5. Dual point calibration with M400 / M800 transmitter

To receive a maximum accuracy of the measured values over the full measuring range, a dual point calibration is required.

By carrying out a dual point calibration both phase angles at zero oxygen (phi 0) and at hundred percent oxygen (phi 100) of the sensor can be established.

Point 1: Slope correction (with air or other calibration media with known O<sub>2</sub> value)

After the sensor signal has stabilized, the complete measuring system can then be calibrated to the 100 % value of the desired measurable variable, e.g. 100 % air, 20.95 % O<sub>2</sub>, or 8.26 ppm at 25 °C (77 °F) and normal pressure (see instruction manual for the transmitter).

Point 2: Zero point

After the sensor signal has stabilized, the sensor can be calibrated to the 0% value of the desired measurable variable, e.g. 0% air, 0.0 % O<sub>2</sub>, or 0 ppm at 25 °C / 77 °F (see instruction manual for the transmitter).

**Note:** Incorrect zero point calibration is a frequent source of measurement error. For correct calibration, we recommend the use of nitrogen gas or other oxygen-free medium with a level of purity of at least 99.99 %.

### **5.2.6. Calibration when connected with analog signal (via M50)**

As soon as a sensor is connected with an analog communication, only rescaling of the M50 output is available.

The procedures can be performed the same way as with amperometric sensors:

1. Slope correction
2. Offset correction

## 6 Maintenance

 **Note: All maintenance work can be done without any tools.**

### 6.1. Sensor head inspection

#### 6.1.1. Visual inspection

To check your sensor head, we recommend the following procedure:

- The contacts of the connector must be dry. Moisture, corrosion and dirt in the connector can lead to false readings.
- Check the cable for buckling, brittle areas or ruptures.
- Check the contacts of the InSUS fiber optical cable.

 **Attention! Do not use any cleaning agents containing alcohol or any solvents. This could damage the sensor head.**

#### 6.1.2. Testing the sensor with the transmitter

If the measured values differ from the expected value, a air calibration should be performed.

**Appropriate phase values after a correct calibration:**  
Depending on the age of the single-use sensor the phase values typically decrease over time compared with a new single-use sensor (see table).

| New sensor                 |                            | Limit for old sensor |                |
|----------------------------|----------------------------|----------------------|----------------|
| Phi0                       | Phi100                     | Phi0                 | Phi100         |
| $62^{\circ} \pm 3^{\circ}$ | $35^{\circ} \pm 3^{\circ}$ | $< 58^{\circ}$       | $> 38^{\circ}$ |

The phase values of the sensor are stored in the calibration history. The actual phase value can be checked in the "Calibration – Verify" menu.

If after such procedures the above mentioned values are still not reached, the sensor should not be used.

Zero oxygen measurement can be done by using CO<sub>2</sub> or nitrogen (N<sub>2</sub>), alternatively in a sample medium saturated with one of these gases.

After 2 minutes in an oxygen-free sample medium, the reading on the transmitter should drop to below 5 % of the reading in ambient air, and within 10 minutes the value should have dropped to below 1 %.

If after such procedures the above mentioned values are still not reached, the sensor should not be used.

## **7 Storage**

For storage, the sensor head should be clean and dry. The protection caps have to be placed on the sensor and the cable connectors.

## 8 Product specification

### 8.1. Specifications

|  | <b>InSUS Head</b>  |
|--|--|
| <b>Measurement principle</b>   | optical  |
| <b>Power requirements</b>  | $U_{\min} = 19.5 \text{ VDC}$<br>$U_{\max} = \leq 25 \text{ VDC}$<br>$P_{\max} = 0.75 \text{ W} \rightarrow 33 \dots 40 \text{ mA}$                    |
| <b>Design features</b>   |  |
| Temperature compensation   | automatic with built-in RTD in InSUS fiber optical cable   |
| Cable connection (digital)   | data cable 5 pin   |
| <b>Single-use O<sub>2</sub> sensor integrated in single-use device</b> |  |
| Sensing element  | Optical spot   |
| Measuring range  | 0 ... 250 % air  |
| Accuracy under defined laboratory conditions                           | $< 2.5 \%$ for the range 50 – 100 % air after 1-point calibration in 100 % air<br><br>$< 1 \%$ after two point calibration in 100 % air and 0 % oxygen |
| Response time  | $T_{98} < 30 \text{ s}$  |
| Temperature probe  | Pt 1000 class B according to DIN EN 60751, integrated in InSUS optical fiber cable   |
| Recommended wetting time (preconditioning time)                        | 10 minutes   |
| Operating temperature  | 5 °C to 60 °C (41 °F to 140 °F)  |
| Sensor cable connection  | Glass fiber extension, M 16 × 1 thread   |

## 9 Ordering information

For more detailed information ask your local distributor.

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All rights reserved