4000TOCe & M300 Transmitter Quick Start Guide

English	3
Deutsch	17
Français	31
Italiano	45
Español	59
Česky	73
Dansk	87
Nederlands	101
Suomi	115
Magyar	129
Polski	143
Português	157
Svenska	171
Русский	185
中文	201
日本語	215
한국어	229
ไทย	243

4000TOCe & M300 Transmitter
30 415 918 Rev A

METTLER TOLEDO

A graphic element consisting of a series of parallel, slightly curved lines that create a sense of depth and movement, resembling a stylized 'M' or a series of steps, positioned behind the company name.

Content

1	Safety Instructions	4
2	Unpacking and Inspection	5
3	Dimensional Drawings	6
4	4000TOCe Sensor Installation	6
5	Electrical Connection	7
6	Terminal Block (TB) definitions	8
7	4000TOCe Start-up	8
8	4000TOCe Sensor	9
9	Menu Structure	9
10	Operation	10
11	General Setup	10
12	Disposal	12
13	Specifications	12
14	Service and Maintenance	13

Statement of intended use – The 4000TOCe is designed to measure Total Organic Carbon (TOC) concentration in high purity water applications. If this equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment against hazards may be impaired.

Warning!

Please refer to the 4000TOCe Operation Manual for safety instructions regarding installation and start-up. Follow all warnings, cautions, and instructions indicated on and supplied with this product.

1 Safety Instructions

This manual includes safety information with the following designations and formats.

Definition of equipment and documentation symbols and designations

WARNING: POTENTIAL RISK OF PERSONAL INJURY.

CAUTION: possible instrument damage or malfunction.

NOTE: Important operating information.

On the transmitter or in this manual text indicates: Caution and/or other possible hazard including risk of electric shock (refer to accompanying documents).

The following is a list of general safety instructions and warnings. Failure to adhere to these instructions can result in damage to the equipment and/or personal injury to the operator.

- The 4000TOCe and M300TOC transmitter should be installed and operated only by personnel familiar with the transmitter and who are qualified for such work.
- The 4000TOCe and M300TOC transmitter must only be operated under the specified operating conditions.
- Repair or service of the 4000TOCe and M300TOC transmitter must be performed by authorized, trained personnel only.
- With the exception of routine maintenance, cleaning procedures or fuse replacement, as described in this manual, the 4000TOCe and M300TOC transmitter must not be tampered with or altered in any manner.
- METTLER TOLEDO accepts no responsibility for damage caused by unauthorized modifications to the transmitter.
- Follow all warnings, cautions, and instructions indicated on and supplied with this product.
- Install equipment as specified in this instruction manual. Follow appropriate local and national codes.
- Protective covers must be in place unless qualified personnel are performing maintenance.
- If this equipment is used in a manner not specified by the manufacturer, the protection provided by it against hazards may be impaired.

WARNINGS:

- Installation of cable connections and servicing of this product require access to shock hazard voltage levels. Switch or circuit breaker shall be in close proximity to the equipment and within easy reach of the OPERATOR; it shall be marked as the disconnecting device for the equipment.
- Main power and relay contacts wired to separate power source must be disconnected before servicing.
- Main power must employ a switch or circuit breaker as the disconnecting device for the equipment.
- Electrical installation must be in accordance with the National Electrical Code and/or any other applicable national or local codes.

1 Safety Instructions continued

– Safety and performance require that this 4000TOCe instrument be connected and properly grounded through a three-wire power source.

NOTE! RELAY CONTROL ACTION: the M300TOC transmitter relays will always de-energize on loss of power, equivalent to normal state, regardless of relay state setting for powered operation. Configure any control system using these relays with fail-safe logic accordingly.

NOTE! PROCESS UPSETS: Because process and safety conditions may depend on consistent operation of this sensor, provide appropriate means to maintain operation during sensor cleaning, replacement or sensor or instrument calibration.

NOTE: This is a 4-wire-product with an active 4–20 mA analog output. Please do not supply power to Pin1–Pin3 of TB2.

2 Unpacking and Inspection


4000TOCe: Inspect the shipping container. If it is damaged, contact the shipper immediately for instructions. Do not discard the container. If there is no apparent damage, carefully unpack the 4000TOCe Sensor. The box should contain the following items:

- 4000TOCe Sensor
- Installation Kit (see operation manual for list of contents)
- Operation Manual CD-ROM
- Quick Start Guide
- Product Insert Card
- Certificate of Calibration
- Declaration of Conformity

M300TOC: Inspect the shipping container. If it is damaged, contact the shipper immediately for instructions. Do not discard the container. If there is no apparent damage, unpack the container. Be sure all items shown on the packing list are present.

If items are missing or shipping container is damaged, notify METTLER TOLEDO immediately.

3 Dimensional Drawings


4 4000TOCe Sensor Installation

Installation

Mount the 4000TOCe in an appropriate location. A minimum clearance dimension of 18" (45.7 cm) is required on the left side of the sensor for UV lamp change-out. Two plumbing connections are provided on the right side of the enclosure. The lower connection is labeled 'Sample Inlet' and the upper connection 'Sample Outlet'. A third connection, a safety drain, is located on the bottom of the enclosure. Figure 1 shows a typical installed configuration. The installation kit provided includes installation tube assemblies, a sample conditioning coil and a high capacity inlet filter. Attach the high capacity inlet filter to the sample inlet connection as shown in Figure 2.

Figure 1


4 4000TOCe Sensor Installation continued

Connect the 0.125" (3 mm) PTFE tubing to a sample point isolation valve (not supplied). Warning: For inlet pressure above 85 psi (5.9 bar), a pressure regulator (Thornton p/n 58 091 552) is required. Flush sample tubing before installing at the sensor. Close the isolation valve. Attach the other end of the PTFE tubing to the sample conditioning coil then attach the sample conditioning coil to the high capacity inlet filter.

Connect the stainless steel drain tube (Thornton p/n 58 091 553) as shown in Figure 2 and direct outlet to an atmospheric drain.

Figure 2


5 Electrical Connection

WARNING: MAIN POWER MUST EMPLOY A SWITCH OR CIRCUIT BREAKER AS THE DISCONNECTING DEVICE FOR THE EQUIPMENT. THE BREAKER SHOULD BE PROPERLY LABELED AND EASILY ACCESSIBLE FOR SYSTEM ISOLATION BY SERVICE PERSONNEL. THE CIRCUIT BREAKER SHALL HAVE A MAXIMUM OF 20 AMP RATING.

Two fittings are provided on the left side of the 4000TOCe enclosure: the AC Power bulkhead fitting for power cable, and the M300TOC patch cable connection. Feed the power cable through the bulkhead fitting on the side-wall. Line and Neutral AC power connections are positioned on the right side of the printed circuit board. The AC Ground connection is provided on a grounding plate located at the inside-top of the sensor door. See Figure 3 for terminal wiring. Ensure line fuse is properly installed. See product label for specified fuse size.

Figure 3

6 Terminal Block (TB) definitions


M300TOC Transmitter wiring connections

For complete instructions see M300TOC Transmitter operation manual.

Power connections are labeled **-N** for Neutral and **+L** for Line, for 100 to 240 VAC or 20–30 VDC.

¼
DIN

TB2 for ¼ DIN	
1	AO1+
2	AO1-/AO2-
3	AO2+
4	-
5	-
6	-
7	-
8	-
9	-


TB1A for ¼ DIN	
1	NO2
2	COM2
3	NC2
4	-
5	-
4	NO4
5	COM4

TB1B for ¼ DIN	
1	NO1
2	COM1
3	NC1
4	-
5	-
4	NO3
5	COM3

½
DIN

TB2 for ½ DIN	
1	AO1+
2	AO1-/AO2-
3	AO2+
4	-
5	-
6	-
7	-
8	-
9	-

TB1 for ½ DIN			
1	NO1	8	-
2	COM1	9	-
3	NC1	10	-
4	NO2	11	NO3
5	COM2	12	COM3
6	NC2	13	NO4
7	-	14	COM4

TB 3 provides access to channel A signal inputs

4000TOCe Sensor Cable Connection

TOC sensor uses 58 080 27X series cables.

Pin no.	Sensor wire Color	Function
1.	-	not used
2.	-	not used
3.	-	not used
4.	-	not used
5.	-	not used
6.	White	GND
7.	Black	RS485-B
8.	Red	RS485-A
9.	Blue	+5V

7 4000TOCe Start-up


Introduce Sample Flow

Apply power, slowly open the sample point isolation valve (not supplied) to start flow to the sensor. Allow approximately 3-5 minutes for the sensor to fill, and observe flow at the drain tube. Once there is flow, check that there are no leaks inside the sensor enclosure. Allow the sensor to initially flush with sample water (4-24 hours is recommended). Close front cover and connect the M300TOC patch cable to the connector on the bottom left side of enclosure.

8 4000TOCe Sensor

The 4000TOCe Sensor housing provides four LED status lights and a single push button for local UV lamp control. See Figure 4.

Figure 4

Sensor status LED: Illuminates whenever AC power is on and the M300TOC and 4000TOCe are connected.


UV Lamp ON LED: When illuminated continuously, UV Lamp is on. When flashing, indicates sensor is in rinse and auto-balance mode.

UV Lamp ON key: Allows the user to turn UV lamp on or off at the sensor. Sensor Status LED: Illuminates continuously when sensor is powered and the patch cable is plugged in. Not lit if patch cable or power is disconnected.

Error LED: Illuminates continuously when an Error condition exists.

Fault LED: Flashing when a Fault condition exists.

9 Menu Structure

10 Operation

Navigation of data entry fields

Use the ► key to navigate forward or the ◀ key to navigate backwards within the changeable data entry fields of the display.

Entry of data values, selection of data entry options

Use the ▲ key to increase or the ▼ key to decrease a digit. Use the same keys to navigate within a selection of values or options of a data entry field.

Note: Some screens require configuring multiple values via the same data field (ex: configuring multiple setpoints). Be sure to use the ► or ◀ key to return to the primary field and the ▲ or ▼ key to toggle between all configuration options before entering to the next display screen.

Navigation with ↑ in Display

If a ↑ is displayed on the bottom right hand corner of the display, you can use the ► or ◀ key to navigate to it. If you click [Enter] you will navigate backwards through the menu (go back one screen). This can be a very useful option to move back up the menu tree without having to exit into the measuring mode and re-enter the menu.

Exit menu

Note: Exit the menu at any time by pressing the ◀ and ► key simultaneously (escape). The transmitter returns to the Measurement mode.


“Save changes” dialogue

Three options are possible for the “Save changes” dialogue:

- “Yes & Exit”: Save changes and exit to measuring mode
- “Yes & ↑”: Save changes and go back one screen
- “No & Exit”: Don’t save changes and exit to measuring mode

The “Yes & ↑” option is very useful if you want to continue configuring without having to re-enter the menu.

11 General Setup (PATH: Menu/Quick Setup)

The 4000TOCe Sensor is configured from the M300TOC transmitter.

While in Measurement mode press the [MENU] key on the M300TOC to bring up the Menu selection. Select “Quick Setup” and press the [ENTER] key.

Display Convention:

1st line on display → a 3rd line on display → c
2nd line on display → b 4th line on display → d

Select the units of measurement for a and b. Only lines a and b can be configured in Quick setup. Go to the Configuration Menu to configure lines c and d.

11 General Setup (PATH: Menu/Quick Setup) continued

5.0 ppbC
0.16 µScm
Aout1 min= 0.000 ppb
Aout1 max= 100.0 ppb

Analog outputs

By selecting Yes the linear 4–20 mA analog output Aout1 will be set up for measurement a when [ENTER] is pressed. Selecting No means that no analog output is set up.

Aout1 min, Aout1 max are the minimum and maximum measurement values for the 4 and 20 mA values respectively.

5.0 ppbC
0.16 µScm
a Set Point Yes
SPl Type= High

Set Points

After configuring the Analog Output a setpoint can be configured for that output. If No is selected and [ENTER] is pressed then the quick setup is done and the menus are exited without setting up any setpoint.

Selecting Yes and pressing [ENTER] means a setpoint can be configured for measurement a.

Following setpoint Types can be selected:

- High (High value has to be set)
- Low (Low value has to be set)
- Between (High and Low value has to be set)
- Outside (High and Low value has to be set)

- USP (% safety margin below U.S. Pharmacopoeia limits)
- EP PW (% safety margin below European Pharmacopoeia limits for Purified Water)
- EPWFI (%safety margin below European Pharmacopoeia limits for Water for Injection)
- JP Cond (% safety margin below Japanese Pharmacopoeia limits)

5.0 ppbC
0.16 µScm
SPl High = 0.000 ppb

After setting the set point value(s) a Relay (none, 1, 2, 3, 4) can be configured for that setpoint. The Relay delay is set to 10 seconds and the Hysteresis is set to 5%.

UV Lamp

Once flow has been established, the 4000TOCe is ready to begin taking TOC measurements. If sensor activation is desired at this time, set the UV Lamp to On. Allow approximately one minute for readings to begin. Additional time (approximately 4-24 hours) may be necessary for the readings to stabilize as the sample line rinses clean, and the TOC sensor reaches thermal stability.

Auto Start

The 4000TOCe sensor can be configured such that the UV lamp will automatically turn on after recovering from a power loss or fault condition. If automatic recovery is desired, set the Auto Start option to 'Yes'.

It is recommended that Auto Start be set to 'Yes'. The default setting is 'No'.

12 Disposal

When the sensor is finally removed from service, observe all local environmental regulations for proper disposal.

13 Specifications

TOC Sensor

Measurement Range	0.05 - 1000 ppbC ($\mu\text{gC/L}$)
Accuracy	± 0.1 ppbC for TOC < 2.0 ppb (for water quality > 15 M Ω -cm [0.067 $\mu\text{S/cm}$]) ± 0.2 ppbC for TOC > 2.0 ppb and < 10.0 ppb (for water quality > 15 M Ω -cm [0.067 $\mu\text{S/cm}$]) $\pm 5\%$ of measurement for TOC > 10.0 ppb (for water quality 0.5 to 18.2 M Ω -cm [2.0 to 0.055 $\mu\text{S/cm}$])
Repeatability	± 0.05 ppbC < 5 ppb, $\pm 1.0\%$ > 5 ppb
Resolution	0.001 ppbC ($\mu\text{gC/L}$)
Analysis Time	Continuous
Initial response time	< 60 seconds
Limit of Detection	0.025 ppbC

Conductivity Sensor

Conductivity Accuracy	$\pm 2\%$, 0.02-20 $\mu\text{S/cm}$; Constant Sensor*
Cell Constant Accuracy	$\pm 2\%$
Temperature Sensor	Pt1000 RTD, Class A
Temperature Accuracy	$\pm 0.25^\circ\text{C}$

Sample Water Requirements

Temperature	0 to 100 $^\circ\text{C}$ **
Particle Size	<100 micron
Minimum Water Quality	≥ 0.5 M Ω -cm (≤ 2 $\mu\text{S/cm}$), pH < 7.5 ***
Flow Rate	20 mL/min
Pressure	4 to 200 psig (0.3 bar(g)) to 13.6 bar(g)) at sample inlet connection ****

General Specifications

Case Dimensions	11" [280mm] W x 7.4" [188mm] H x 5.25" [133mm] D
Weight	5.0 lb. (2.3 kg)
Enclosure material	Polycarbonate plastic, flame retardant, UV and chemical resistant UL # E75645, Vol. 1, Set 2, CSA #LR 49336
Enclosure rating	NEMA 4X, IP65 Industrial environment
Ambient Temperature/ Humidity rating	5 to 50 $^\circ\text{C}$ / 5 to 80% Humidity, non-condensing
Power requirements	100 - 130VAC or 200 - 240VAC, 50/60 Hz, 25W Maximum
Local Indicators	Four LED lights for Fault, Error, Sensor Status and UV Lamp ON
Ratings/Approvals	CE Compliant, UL and cUL (CSA Standards) listed, Conductivity and temperature sensors traceable to NIST, ASTM D1125 and D5391. Meets ASTM D5173 Standard Test Method for On-Line Monitoring of Carbon Compounds in Water by UV Light Oxidation

13 Specifications continued

Sample Connections

Inlet connection	10-32 female threaded port (6' [2m] FDA compliant PTFE tubing supplied)
Outlet connection	10-32 female threaded port (Fixed right angle 316SS drain tube provided)
Inlet Filter	316SS, inline 60 micron
Wetted parts	316SS/Quartz/PEEK/Titanium/PTFE/EPDM/FFKM
Wall Mount	Standard, mounting tabs provided
Pipe Mount	Optional, with pipe-mount bracket accessory for nominal pipe sizes 1" [2.5 cm]
Maximum Sensor Distance	300ft [91 m]

* Readout in equivalent S/m ranges selectable at M300TOC

** Temperature above 70°C requires Sample Conditioning Coil (included)

*** For power plant cycle chemistry samples, pH may be adjusted by measurement after cation exchange.

**** Process pressure above 85 psig (5.9 bar) requires optional High Pressure Regulator p/n 58 091 552.

Specifications subject to change without notice.

14 Service and Maintenance

The 4000TOCe Sensor is designed to minimize service and maintenance. Therefore, normal operating wear and tear is negligible. This reduces the amount of consumable components as well as the amount of time needed to maintain the sensor. Listed below are instructions on how to perform simple periodic maintenance, which includes UV Lamp change (every 4500 hours of operation), filter replacement (typically every 6 months), and general cleaning.

UV Lamp Replacement


WARNING: UV RADIATION HAZARD

Apply power to UV lamp only when installed in housing in accordance with operation manual. **DO NOT** remove UV lamp from housing unless power is off.

Always protect eyes and skin from direct exposure to UV light.

METTLER TOLEDO Thornton recommends replacement of the UV lamp inside the 4000TOCe Sensor after 4500 hours of operation, or 6 months of continuous use, not to exceed one year. This is a simple procedure that requires only a few minutes to complete. The following steps explain the proper procedure for the UV lamp replacement. Refer to the illustration below.

14 Service and Maintenance continued


CAUTION: Use of a UV lamp other than those provided by METTLER TOLEDO Thornton specifically for use with the 4000TOCe Sensor will affect performance and void the warranty of this product.

1. At the sensor, turn off the UV lamp by pressing the lamp control button (UV lamp ON LED will turn off). If the LED does not turn off, check that the Sensor Key Lock is in the off position in the M300TOC. Refer to Section 10.3.1 Sensor Key Lock function in the operation manual.
2. Once power to the UV lamp is off, open the front cover of the sensor enclosure with the front cover tool.
3. Remove the side cover labeled 'UV LAMP REPLACEMENT COVER' on the left side of the sensor enclosure. Use a large-bladed flat-head screwdriver and turn the cover counter clockwise to loosen and unscrew the cover.
4. Disconnect the power cable to the UV lamp. This connector is located on the backside of the front cover, above the circuit board.
5. Loosen the UV lamp holding screw located on the left side of the oxidation chamber.
6. Slide the cable of the UV lamp through the side opening of the enclosure and gently slide the UV lamp out of the oxidation chamber assembly (stainless steel cylinder). Be careful not to let the UV lamp hit the quartz glass tube inside the chamber.
7. Use the gloves supplied with each replacement bulb. Hold the new lamp from the ends of the lamp. Do not touch the bulb. Slide the new UV lamp into the side opening of the enclosure and into the oxidation chamber opening until it stops. Do not use excessive force to insert the UV lamp as this may cause damage to the lamp or the internal components of the oxidation chamber.
8. Tighten the UV lamp holding screw until snug. Do not over-tighten.
9. Feed the power cable through the side opening of the enclosure. Re-connect it to the power connector on the front door.
10. Close front cover of the sensor and secure fasteners with the front cover tool.
11. Install the UV Lamp replacement cover on the opening on the side of the enclosure.
12. At the M300TOC, press the 'Menu' key, and select Path: Menu/Configure/TOC Setup/Lamp Control. Press the 'Enter' key until 'Lamp Time Rese' is displayed. Select 'Yes', then press enter twice to confirm your selection and save the changes.
13. After the lamp replacement has been completed, a TOC calibration must be performed. In addition, a System Suitability Test is recommended in applications where compliance with Pharmacopeia regulations is necessary.

14 Service and Maintenance continued

High Capacity Inlet Filter Replacement

The 4000TOCe Sensor includes a high capacity filter (shown below) containing a filter element that should be replaced (p/n 58 091 551, package of 2) every 6 months or sooner, depending on water quality conditions. Detailed instructions to replace this filter are included in the replacement package.

REPLACING THE FILTER ELEMENT FOR THE HIGH CAPACITY INLET FILTER

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© Mettler-Toledo Thornton 08/2017. Printed in United States.
Switzerland **Mettler-Toledo (Schweiz) GmbH**
METTLER TOLEDO, GmbH
METTLER TOLEDO Thornton, Inc.