InPro 3300 Non-Glass pH Electrode

Instruction Manual

Contents

cs	Návod k obsluze	2
da	Brugsanvisning	4
de	Betriebsanleitung	6
en	Instruction Manual	8
es	Manual de instrucciones	10
fi	Käyttöohje	12
fr	Instructions d'utilisation	14
hu	Használati utasítás	16
it	Istruzioni per l'uso	18
ja	取扱説明書	20
ko	사용 매뉴얼	22
nl	Gebruikershandleiding	24
pl	Instrukcja obsługi	26
pt	Manual de instruções	28
ru	Инструкция по эксплуатации	30
sv	Bruksanvisning	32
th	คู่มือคำแนะนำ	34
zh	说明书	36

InPro 3300

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
InPro 3300

Non-Glass pH Combination Electrode

Instruction Manual

1. Introduction

METTLER TOLEDO pH sensor type InPro® 3300 is a low-maintenance, autoclavable and pressure-resistant electrode with a glass-free sensing head and thickened reference electrolyte. Please read these instructions carefully before putting the electrode into operation, in order to ensure perfect operation. All wetted parts are made of materials in compliance with the FDA directives (Food and Drug Administration). We recommend that the sensor may only be in operation with the original accessories from METTLER TOLEDO. The sensor may be operated and serviced only by personnel familiar with the unit, and who have read and understood this instruction manual.

2. Safety remarks

If an electrode is to be cleaned or calibrated using acid or alkaline solutions, safety goggles and gloves should be worn.

3. Product designation

The markings on the polymer shaft include the following information:

METTLER TOLEDO	Electrode/sensor manufacturer
InPro 3300/120/Pt xxx	Designation/shaft length in mm/ temperature sensor
Combination pH	Type of electrode
pH 0...14	pH measurement range
0...80 °C (130 °C)	Temperature range for operation/ (sterilization)
Order No. 52 xxx xxx	Article number for ordering

Additionally, the electrode is supplied with a serial number marked on the connector cap to enable individual identification of each electrode produced.

4. Installation and preparation for use

- Do not touch, hold or make any contact with the plug and/or pins of the electrode. Electrostatic discharge can damage the chip.

2. On unpacking, check the electrode for mechanical damage. Report any signs of damage immediately to your METTLER TOLEDO supplier.
3. Remove the watering cap and briefly rinse the electrode with de-ionized water. An already used electrode should first be cleaned by lightly scrubbing the sensing element with a soft-bristled toothbrush and a mild soap solution, followed by rinsing with de-ionized water.
4. Connect the electrode to the pH transmitter, using cable VP6-ST/ISFET. Please refer to the wiring and terminal diagram accompanying the cable.
5. Prior to calibration, soak the electrode in buffer pH 7.00 or pH 4.01 for at least 10 minutes, with the electrode connected to the transmitter.
6. Before putting the electrode into operation, carry out a zero point calibration. Exact instructions for zero point calibration are given in the instruction manual of the associated pH transmitter.

5. Operation

5.1 Calibration of the electrode and the pH transmitter

A 2-point pH calibration is recommended e.g. using buffers pH 7.00 and pH 4.01. Please refer to the instruction manual of the pH transmitter for further details.

5.2 CIP-cleanability (Cleaning-In-Place)

A standard procedure for maintaining hygienic conditions in food processing is Cleaning-In-Place (CIP). Periodic flushing with hot caustic soda (NaOH) will however shorten the lifetime of the electrode. We therefore recommend to retract the electrode from the process and to clean it with hot steam (sterilization procedure).

6. Maintenance

1. Carefully rinse the sensor tip and diaphragm with de-ionized water after each production cycle. Under no circumstances must measuring solution be allowed to dry on these parts!
2. When the sensor is not in operation, store it with sensor tip and diaphragm well submerged into 9823/3M KCl or 9816/Viscolyt™ electrolyte.
3. If the sensor is stored mounted in its housing, the same rule as in 2) above applies, but the storage electrolyte has to be slightly modified, with 2 parts of buffer solution pH 9.2 to 10 parts of electrolyte, this in order to avoid corrosion of housing parts.
4. If a sensor is stored dry for a few days by mistake, let it soak in the normal storage electrolyte for several hours before use.
5. Occasionally check the connector for possible traces of moisture. If necessary, clean well with de-ionized water or alcohol, and afterwards dry carefully.
6. The silicone O-ring (marked No. 4 in the drawing) is to be replaced at least every six months or, if necessary, at a shorter interval.

7. Trouble-shooting

For cleaning, lightly scrub the sensing element of the electrode with a soft toothbrush and water to remove any possible residues. Use mild detergent if necessary. Possible protein precipitation and blockage of the diaphragm can be treated with cleaner 51 340 070.

 **Attention:** HF-containing solutions will destroy the pH sensitive element of the electrode.

8. Disposal

If no official regulations are in force, used or defective electrodes and the packaging can be disposed as ordinary household waste.

9. Warranty

In the event of manufacturing faults, a warranty period of 12 months from date of production is granted.

InPro and Viscolyt are trademarks of the Mettler-Toledo Group.

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02/2016 © Mettler-Toledo GmbH
Analyse Industrielle S.A.S.
- Germany** **Mettler-Toledo GmbH**
- Switzerland** **Mettler-Toledo (Schweiz) GmbH**
Mettler-Toledo GmbH