METTLER TOLEDO

InPro 485X i pH Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction manual	14
es	Manual de instrucciones	18
fr	Instructions d'utilisation	22
fi	Käyttöohje	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

InPro 485Xi pH elektroda

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Obsah

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Instruction Manual

Contents

1	Introduction	15
2	Safety instructions	15
3	Product description	15
4	Installation and commissioning	16
5	Operation: Calibration of the electrode and pH transmitter	16
6	Maintenance.....	17
7	Removal of dirt and residue.....	17
8	Disposal	17
9	Guarantee	17
10	Storage conditions	17

InPro and ISM are registered trademarks of
METTLER TOLEDO Group in Switzerland, the USA,
the European Union and a further five countries.

1 Introduction

METTLER TOLEDO's combination pH electrodes with pNa reference system and temperature sensor offer reliable pH measurement in solutions with a near stable concentration of sodium ions. This electrode offers top performance in harsh brine processes with oxidizing compounds e.g. chlorine.

Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

2 Safety instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis, both protective spectacles and protective gloves should be worn. The sensor can contain up to 0.3 % boric acid by weigh.

Do not touch the electrode at the plug-in contacts, as electrostatic discharges can damage the electronics.

3 Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
InPro™ 485Xi/SG/yyy	Type designation x = membrane glass, i = ISM Digital/SG = ancillary electrode/yyy = rod length in mm
combination pH	Type of electrode
pH 0...14 (InPro 4850i)	pH measurement range
- 10...120 °C (14...248 °F)	Temperature range for operation
Order No. 52 00X XXX	Order number
 SEV 14 ATEX 0168 X	Identification and number of certificate
 IECEx SEV 14.0025X	Identification and number of certificate
 CML 22 UKEX 2108X	Identification and number of certificate
IS/I, II, III/1/ABCDEFG/T6	FM tested
See instruction manual	Instruction (Observe the operating instructions)

In addition, each electrode is allocated a serial number on the red part of the plug-in head to enable identification.

Important! Please take note of the following technical specifications of this pH electrode:

Reference system	Sodium sensitive glass membrane
Min. Na ⁺ concentration	10 mg/L when pH > 7; 100 mg/L when 7 > pH > 2; 1g/L when pH < 2

4 Installation and commissioning

1. While unpacking, check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the protective cap and rinse the electrode for a short time with Buffer pH = 4.01 / Na⁺ 3.9M. After rinsing, the electrode should only be lightly dabbed. Any rubbing of the pH or the sodium-sensitive glass can increase the response time as a result of electrostatic charge.
3. Check the space behind the pH and the sodium sensitive glass for air bubbles and remove any air bubbles by shaking the electrode slightly in a vertical direction.
4. Place the electrode into the housing as described in the instructions for the housing.
5. Connect the electrode and the pH transmitter with the corresponding cable (AK9). Observe the connections scheme that accompanies the cable or the transmitter.
6. ISM™ sensors such as InPro™ 485Xi enable Plug and Measure and expanded diagnostics. For installation, commissioning and operation of the ISM system please refer to instructions for the transmitter, module and cable.

5 Operation: Calibration of the electrode and pH transmitter

Before a calibration, immerse the electrode for 10 minutes in a Buffer pH = 7.00 / Na⁺ 3.9M or Buffer pH = 4.01 / Na⁺ 3.9M buffer solution while still connected to the transmitter. A 2-point calibration is recommended, e.g. Buffer pH = 7.00 / Na⁺ 3.9M and Buffer pH = 4.01 / Na⁺ 3.9M buffer solutions.

For further details please consult the operating instructions for the pH transmitter.

6 Maintenance

1. After each working cycle, the sensitive part can be cleaned with a soft toothbrush and a diluted HCl-solution (0.1 and 1 mol/L), and then rinsed with 20% NaCl solution.
2. If the electrode is not in use it must be stored with the electrode tip immersed in Buffer pH = 4.01 / Na⁺ 3.9M.
3. If an electrode is inadvertently stored in a dry condition for several days, the membrane must be immersed in Buffer pH = 4.01 / Na⁺ 3.9M for at least 6 hours.
4. The plug-in head should occasionally be checked for possible traces of moisture. If necessary, clean it thoroughly with deionized water or alcohol and then dry it off carefully.

7 Removal of dirt and residue

With a soft toothbrush and water any dirt/residue can easily be removed from the electrode. Mild washing agents can also be used.

Warning! Solutions containing hydrofluorides damage the pH sensitive part of the electrode.

8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9 Guarantee

On manufacturing defects, 12 months after delivery.

10 Storage conditions

We recommend that the electrodes are stored at room temperature. Keep the electrode tip in the supplied cap with Buffer pH = 4.01 / Na⁺ 3.9M.

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本表依据SJ/T 11364的规定编制。本产品符合以下标志规范。
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○: 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○: Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x: 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x: Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 08/2022 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 52 004 112 B

Management System
certified according to
ISO 9001 / ISO 14001

www.mt.com/pro

* 5 2 0 0 4 1 1 2 B *