METTLER TOLEDO

InPro 42XX (i) Series pH Electrode

Instruction Manual

cs	Návod k obsluze	2
da	Brugsanvisning	6
de	Betriebsanleitung	10
en	Instruction Manual	14
es	Manual de instrucciones	18
fr	Instructions d'utilisation	22
fi	Käyttöohje	26
hu	Használati utasítás	30
it	Istruzioni per l'uso	34
ja	取扱説明書	38
ko	사용 매뉴얼	42
nl	Gebruikershandleiding	46
pl	Instrukcja obsługi	50
pt	Manual de instruções	54
ru	Инструкция по эксплуатации	58
sv	Bruksanvisning	62
th	คู่มือคำแนะนำ	66
zh	说明书	70

Řada InPro 42XX (i) pH elektroda

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Obsah

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Instruction Manual

Contents

1	Introduction	15
2	Safety instructions	15
3	Product description.....	15
4	Installation and commissioning	16
5	Operation: Calibration of the electrode and pH transmitter.....	16
6	Maintenance	16
7	Removal of sources of malfunction	17
8	Disposal.....	17
9	Guarantee.....	17
10	Storage conditions.....	17

InPro, ISM and Xerolyt are registered trademarks of METTLER TOLEDO Group in Switzerland, the USA, the European Union and a further five countries.

1 Introduction

The METTLER TOLEDO InPro™ 42XX(i) pH-single rod measuring cell is a low maintenance, high-pressure resistant pH electrode with a solid polymer electrolyte (Xerolyt® Extra) together with a silver ions barrier. The diaphragm is of an open aperture type, so that direct contact exists between the medium being measured and the electrolyte.

Please read through these operating instructions carefully before commissioning, in order to ensure trouble-free use.

2 Safety instructions

When cleaning or calibrating the electrode using solutions containing acids or alkalis both protective spectacles and protective gloves should be worn.

Do not touch the electrode at the plug-in contacts, since electrostatic discharges can damage the electronics.

We recommend that you only operate the electrode in combination with original parts from METTLER TOLEDO. Operation and maintenance should be carried out only by trained personnel and staff who have read and understood the operating instructions.

3 Product description

The inscription on each electrode contains the following information:

METTLER TOLEDO	Manufacturer of the electrode
InPro™ 42XXi/SG/yyyy/zzz	Type designation; XX=shaft material, membrane glass; i=ISM Digital/ SG=ancillary electrode/ yyyy=rod length in mm/zzz=type of temperature sensor
combination pH	Type of electrode
pH 0 ... 14	pH measurement range
pH 1 ... 14 (InPro 4281i)	
0 ... 100 °C (130 °C)	Temperature range for operation (sterilization)
Order No. 52 00X XXX	Order number
 SEV14 ATEX 0168 X	Identification and number of certificate
 IECEx SEV 14.0025X	Identification and number of certificate
 CML 22 UKEX 2108X	Identification and number of certificate
IS/I, II, III/1 /ABCDEFG/T6	FM tested
See instruction manual!	Instruction (Observe the operating instructions)

In addition each electrode is allocated a serial number on the red part of the plug-in head to enable identification.

4 Installation and commissioning

1. While unpacking check the electrode for mechanical damage. Please report any damage immediately to your METTLER TOLEDO supplier.
2. Remove the protective cap and rinse the electrode for a short time with deionized water. If the electrode has already been used the sensitive part can be cleaned with a soft toothbrush and a mild soap solution, and then rinsed for a short time with deionized water. After rinsing the electrode should only be dabbed. Any rubbing of the pH sensitive glass can increase the response time as a result of electrostatic charge accumulation.
3. Check the space behind the pH sensitive glass for air bubbles and remove any air bubbles that may be present by shaking the electrode slightly in a vertical direction.
4. Place the electrode into the housing as described in the instructions for the housing.
5. Connect the electrode and the pH transmitter with the corresponding cable (VP6, VP8, AK9). Observe the connections scheme that accompanies the cable or the transmitter.
6. Intelligent ISM™ sensors such as InPro 426Xi enable "Plug and Measure" and expanded diagnostics. For installation, commissioning and operation of the ISM system please refer to instructions for the transmitter, module and cable.

5 Operation: Calibration of the electrode and pH transmitter

Before a calibration immerse the electrode for 10 minutes in a pH 7.00 or pH 4.01 buffer solution whilst connected it to the transmitter. A 2-point calibration is recommended, e.g. pH 7.00 and pH 4.01 buffer solutions. For further details please consult the operating instructions for the pH transmitter.

6 Maintenance

1. After each working cycle the electrode tip and diaphragm have to be carefully rinsed with deionized water. Any remaining of the solution being measured on these parts must be avoided at all costs!
2. If the electrode is not in use it has to be stored with the electrode tip and the diaphragm immersed in a tank of electrolyte (9823/3M KCl or FriscoIyft™ 9848).

3. If the electrode is stored in its housing the advices described under Item 2 apply, although the electrolyte has to be slightly modified for storage purposes (2 parts of buffer solution pH 9.2 to 10 parts of electrolyte) in order to avoid any corrosion of the housing parts.
4. If an electrode is inadvertently stored in a dry condition for several days it must be rinsed for several hours in a standard storage electrolyte.
5. The plug-in head should occasionally be checked for possible traces of moisture. If necessary clean it thoroughly with deionized water or alcohol and then dry it off carefully.

7 Removal of sources of malfunction

With a soft toothbrush and water any residues can easily be removed from the electrode. Mild washing agents can also be used.

8 Environmental protection

Waste electrical products should not be disposed of with household waste. Please recycle where facilities exist. Check with your Local Authority or retailer for recycling advice.

9 Guarantee

On manufacturing defects 12 months after delivery.

10 Storage conditions

We recommend that the electrodes are stored at room temperature.

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 本表依据SJ/T 11364的规定编制。本产品符合以下标志规范：
 Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○ : 表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
 ○ : Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x : 表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
 x : Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

产品中有毒物质的名称及含量 Toxic and hazardous substance name and containment in product						
部件名称 Part Name	有毒有害物质或元素 Toxic and hazardous substances					
	铅 Lead (Pb)	汞 Mercury (Hg)	镉 Cadmium (Cd)	六价铬 Hexavalent chromium (Cr6+)	多溴联苯 Polybrominated biphenyls (PBB)	多溴二苯醚 Polybrominated diphenyl ethers (PBDE)
插头 Elect. Connector	x	○	○	○	○	○
电路板 PCBA	x	○	○	○	○	○

本表依据SJ/T 11364的规定编制。本产品符合以下标志规范。
Table composed in accordance with SJ/T 11364 (CN). This product is bearing the following symbol:

○：表示该有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下
○：Indicates that the content of the hazardous substance in all homogeneous materials of the part is below the limit specified in GB / T 26572

x：表示该有害物质至少在该部件的某一均质材料中的含量超出GB/T 26572规定的限量要求。
x：Indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limits specified in GB/T 26572.

For addresses of METTLER TOLEDO
Market Organizations please go to:
www.mt.com/contacts

METTLER TOLEDO Group
Process Analytics
Local contact: www.mt.com/contacts

Subject to technical changes
© 08/2022 METTLER TOLEDO
All rights reserved
Printed in Switzerland. 52 005 365 B

Management System
certified according to
ISO 9001 / ISO 14001

**UK
CA**₂₅₀₃

CE₁₂₅₈

EAC

www.mt.com/pro

* 5 2 0 0 5 3 6 5 B *