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
Ø 12 mm / Ø 25 mm

- S-96 Type Kit with 4membranes, Order No. 52 200 025
- T-96 Type Kit with 4membranes, Order No. 52 200 024
- T-6800 Gas Type Kit with 4 membr., Order No. 52 201 149
- S-96 Type Kit with 20 membranes, Order No. 52 200 790
- T-96 Type Kit with 20 membranes, Order No. 52 200 791

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
Kathode

Membrankörper

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Instructions for electrolyte and membrane exchange of InPro® 6000 series O₂ sensors

ENGLISH

Depending on your order the membrane kit contains of either T-96 or S-96 type membrane body, which can be identified by the marks T or S on the body. Please make sure you received the right membrane body.

Important notice:

The O₂ electrolyte has an alkaline pH of 13. Contact of electrolyte with the skin, especially mucous membrane or eyes, should be avoided. If such contact occurs, the affected area should be well rinsed with water. Get medical attention if adverse signs appear.


Great care should be exercised when handling the glass inner bodies, since any hairline cracks resulting from knocks, adversely affect sensor performance.

Because contact with the electrolyte is very likely during exchange of electrolyte or membrane body, the use of protective gloves is recommended.

When the membrane begins to exhibit signs of failure (long response time, mechanical damage, increased residual current in oxygen free medium, etc.), it has to be replaced.

When replacing the membrane body and electrolyte, strictly observe the following instructions:

1. Unscrew the cap sleeve from the shaft and carefully pull it from the sensor.
2. Usually the membrane body remains in the cap sleeve and should be ejected by pushing it with the flat finger tip. Otherwise pull the membrane body from the interior body.
Before electrolyte is refilled, the membrane body **must** be removed from the cap sleeve.
3. Clean the interior body with a paper handkerchief or paper tissue.
4. Check all O-rings and the trapezoid silicone sealing for mechanical defects and replace if necessary.
5. Half-fill the membrane body with O₂ electrolyte and make sure that all bubbles are removed.
6. Slip the membrane body over the interior body while holding the sensor in a vertical position. The excess electrolyte will be displaced and must be removed with a paper cloth.
7. Carefully slip the cap sleeve over the fitted membrane body and screw it down.
8. Remove expelled electrolyte with a paper tissue.
9. After each exchange of electrolyte or membrane body the sensor has to be repolarized and recalibrated (see instruction manual InPro® 6000 series O₂ sensors, chapter 2 and 3).

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09/18 © Mettler-Toledo GmbH, Switzerland