

Short Instructions

Security Codes:

Set-up (operation parameters)

S.CD 22

Calibration

C.CD 11

How to enter in set-up menu

Press  , display shows000
S.CD.Press  to scroll to s.CD. 22.

TC

Press  , display showsPress  to select desired function.

How to return to measurement

Press   together, escape function.

How to enter in calibration menu

Press , display shows

000
C.CD.

Press  to scroll to C.CD. 11.

Press , display shows

SETUP HOLD
CAL
Con

Calibration

SETUP HOLD
CAL
Con
6

Enter calibration menu and press  to start calibration.

SETUP HOLD
0.100
CEL
6

Adjust cell constant with   keys and press .

SETUP HOLD
14.3
25.0
μS
6

Put sensor in calibration solution and use   to match measured value ($\mu\text{S}/\text{mS}$) to standard solution, press .

SETUP HOLD
1.01
FCt
6

The display shows the revised factor. Press  to return to measurement.

Set ATC or manual temperature compensation

SETUP HOLD
SEt
°C °F

SETUP HOLD
On
ATC

SETUP HOLD
OFF
ATC

Enter Setup menu with S.CD 22 and press   until display shows "SEt °C °F", press .

Select unit °C or °F with   keys and press .

Select ATC ON or ATC OFF with   keys and press .

If ATC OFF is selected, set manually calibration and process temperature.

Press  to confirm.

Entering set points

SETUP HOLD
SP1

SETUP HOLD
200 μs
SP1

SETUP HOLD
LO
SP1

Enter in Setup menu and scroll with   to show SP1
Press .

Scroll with   to desired value ($\mu\text{s}/\text{ms}$)
Press .

Set relay in LO or HI mode.
Press . Proceed to adjust hysteresis and delay timer.
Continue same with SP2.